**SECRETARY OF STATE**  
STATE OF INDIANA

Connie Lawson  
Secretary of State

**APOSTILLE**

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1. Country: United States Of America
2. This public document had been signed by *Johnathan Liu*
3. Acting in the capacity of notary public in & for *Monroe County*
4. And bears the seal/stamp of notary public in & for the State of Indiana

**CERTIFIED**

5. at Indianapolis, Indiana
6. This *Fourth* day of *October*, 2018
7. by the Secretary of State of Indiana
8. No. A2018-1004102934
9. Seal/Stamp:
10. Signature:

*Connie Lawson*

**Connie Lawson**  
Indiana Secretary of State

*Effective May 1, 2003 all apostilles from the Indiana Secretary of State will have an electronically printed seal.  
This document was revised March 17, 2012.*

COOK INCORPORATED  
750 DANIELS WAY  
BLOOMINGTON, IN 47404 USA  
PHONE: 812.339.2235 TOLL FREE: 800.457.4500  
WWW.COOKMEDICAL.COM

### CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for "Bakri Postpartum Balloon" Submission.

Kris Weathers  
Regulatory Affairs Manager  
Cook Incorporated


State of Indiana )

ss

County of Monroe )

Subscribed and sworn to before me the 3 day of October 2018.

My Commission expires November 30, 2022

  
Johnathan Liu

## INSTRUCTIONS FOR USE

### Bakri Postpartum Balloon

#### INTENDED USE

Used for temporary control or reduction of postpartum bleeding when conservative management is warranted (see Instructions for Use).

**CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or a properly licensed practitioner).**

#### DEVICE DESCRIPTION

##### Table of conversions

| Parameter | Conversion |
|-----------|------------|
| 1/3 Fr    | 1 mm       |
| 1 Inch    | 25,4 mm    |

#### PERFORMANCE CHARACTERISTICS

##### Strength characteristics of the catheter:

Bond strength of components - not less than 15 N;  
Tensile strength - not less than 8.1 MPa;  
Relative elongation at break point - not less than 150%;  
The maximum burst pressure - 81 atm (1200psi)


##### Strength characteristics of the connecting tube:

Bond strength of components - not less than 15 N;

##### Strength characteristics of the one-way stopcock:

Force at break - not less than 50N;


**Bakri Postpartum Balloon (J-SOS-100500, J SOSR-100500)**

J-SOS-100500

**The Bakri Postpartum Balloon catalogue number J-SOS-100500 comes with a stopcock, and Syringe 60 m.**

The balloon catheter Bakri Postpartum Balloon is comprised of a silicone tube with a diameter of 24.0 Fr / 8 mm and a length of 54.0 cm and a silicone balloon with a volume of 500 ml. Near the distal tip of the balloon catheter there are two sideports (for inflation and deflation of the balloon) located opposite each other. The distal tip is rounded to a smooth finish. A one-way stopcock is connected to an inflation adapter, which is located at the proximal end of the catheter. At the proximal end of the catheter there is also sidearm for deflating the balloon.


J-SOSR-100500

The balloon catheter Bakri Postpartum Balloon, catalogue number J-SOSR-100500 comes with a one-way stopcock, a 60 ml syringe, a dual check valve, a connecting tube with components for rapid installation, namely with an adapter and a cap. The connecting tube has a diameter of 14 Fr / 4.66 mm and a length of 180 cm.

The balloon catheter Bakri Postpartum Balloon is comprised of a silicone tube with a diameter of 24.0 Fr / 8 mm and a length of 54.0 cm and a silicone balloon with a volume of 500 ml. Near the distal tip of the balloon catheter there are two sideports (for inflation and deflation of the balloon) located opposite each other. The distal tip is rounded to a smooth finish. A one-way stopcock is connected to an inflation adapter, which is located at the proximal end of the catheter. At the proximal end of the catheter there is also a sidearm for deflating the balloon.

| Product Catalogue Number | Catheter Outer Diameter Fr/mm | Catheter Inner Diameter mm | Catheter Length (cm) | Inflation Volume (ml) | Remarks                            |
|--------------------------|-------------------------------|----------------------------|----------------------|-----------------------|------------------------------------|
| Tolerance                | ±10%                          | (±0,51) mm                 | ±10%                 | ±10%                  |                                    |
| J-SOS-100500             | 24/8,0                        | 6,2                        | 54                   | 500                   | 100% silicone material             |
| J-SOSR-100500            | 24/8,0                        | 6,2                        | 54                   | 500                   | With Rapid Instillation components |

Scheme 1. General View of the Catheter Balloon Bakri Postpartum Balloon.


Parameters of the Catheter Balloon Bakri Postpartum Balloon.

Length: A – 54 cm (± 2,0 cm); B – 7,9 cm (± 2,0 mm); C – 3,2 cm ( ± 2,0 mm); D – 2,2 mm (± 2,0 mm).

Diameter: E – 8,0 mm (± 0,51mm).

Scheme 2. The Proximal Side of the Catheter Balloon Bakri Postpartum Balloon.

Parameters of Proximal Side of the Catheter Balloon Bakri Postpartum Balloon:

Length: A – 5,21cm ( $\pm 0,127$  mm); B – 3,81 cm ( $\pm 0,127$  mm); C – 2,33 cm ( $\pm 0,127$  mm);


Diameter: D – 4,75 mm ( $\pm 0,51$ mm); E – 8,0 mm ( $\pm 0,51$ mm).

Thickness of the catheter shaft: F – 0,0762 mm ( $\pm 0,051$  mm).

Parameters Adapter for Inflation at the Proximal End of the Catheter:


Diameter: D – 5,334mm ( $\pm 0,051$ mm); E – 8,738mm ( $\pm 0,051$ mm);

Length: F – 4,369 mm ( $\pm 0,051$ mm); G – 7,137mm ( $\pm 0,051$ mm); H – 7,925 mm ( $\pm 0,051$ mm).

| "A"                    | "B"                     | "C"                    |
|------------------------|-------------------------|------------------------|
| 1,778mm $\pm 0,051$ mm | 19,863mm $\pm 0,051$ mm | 1,575mm $\pm 0,051$ mm |

**Scheme 3. Distal Side of the Catheter Balloon Bakri Postpartum Balloon.**

Parameters of Balloon:

Length: A – 10,3 cm ( $\pm 0.127$  mm); B – 8,0 cm ( $\pm 0.127$  mm); C – 7,9 cm ( $\pm 2,0$  mm);  
D – 3,2 cm ( $\pm 2,0$  mm).  
E – 1,4 cm ( $\pm 2,0$  mm).  
Diameter: F – 3,0 mm ( $\pm 1,0$  mm); G – 4,0 mm ( $\pm 2,0$  mm); H – 1,1 cm ( $\pm 4,0$  mm);  
I – 2,0 mm ( $\pm 1,0$  mm).

Scheme 4. Dual Check Valve (only for J-SOSR-100500).


Parameters of Dual Check Valve:

Length: A – 41mm; B – 22mm;

Parameters for Rigid Snap-on Marker:

Diameter: A- 5.537- 6.35 mm ( $\pm 0.127$  mm);  
Length: B- 6.35 mm ( $\pm 0.127$  mm).


Scheme 5. Connecting Tube (only for J-SOSR-100500).

Parameters for Connecting Tube:

Length: A - 180 cm ( $\pm 10$  mm).

**Scheme 6. The Tip of Connecting Tube (only for J-SOSR-100500).**

Parameters for The Tip of Connecting Tube:

Length: A - 33,6 mm ( $\pm 0.127$  mm); B - 19,7 mm ( $\pm 0.127$  mm); C - 9,6 mm ( $\pm 0.127$  mm);

D - 7,8 mm ( $\pm 0.127$  mm); E - 1,2 mm ( $\pm 0.127$  mm); F - 6,5 mm ( $\pm 0.127$  mm).

Diameter: G - 17,8 mm ( $\pm 0.127$  mm); H - 7,3 mm ( $\pm 0.127$  mm); I - 5,0 mm ( $\pm 0.127$  mm).


Bevel Angle - 30°

Bevel Angle - 60°

Diameter - Ø3,0 mm

Bevel Angle - 30°

**Scheme 7. Adapter of Connecting Tube (only for J-SOSR-100500).**


Parameters for Adapter of Connecting Tube:

Diameter: A - 1.829 mm ( $\pm 0.127$  mm);

Length: B - 17.145 mm ( $\pm 0,127$ mm).

#### Scheme 8. One-Way Stopcock


#### Parameters for One-Way Stopcock:

Length: A - 40.49 mm ( $\pm 0.127$  mm);  
 Inner Diameter: B - 2.49 mm ( $\pm 0.127$  mm);  
 Inner Diameter: C - 4.3 mm ( $\pm 0.127$  mm);  
 Height of Handle: D - 18.001mm ( $\pm 0.127$  mm).


#### Parameters of Luer-Lock «Female»:

Diameter: A - 8.509mm ( $\pm 0.127$ mm); B - 3.277mm ( $\pm 0.127$ mm); C - 4.749mm ( $\pm 0.127$ mm).  
 Length: D - 7.62mm ( $\pm 0.127$ mm); E - 7.137mm ( $\pm 0.127$ mm); F - 1.397mm ( $\pm 0.127$ mm);  
 G - 7.366 mm ( $\pm 0.127$ mm); H - 18.415 mm ( $\pm 0.127$ mm); I - 4.699mm ( $\pm 0.127$ mm).

#### Parameters of Luer-Lock «Male»:

Diameters: B - 3.937mm ( $\pm 0.127$ mm); C - 4.775mm ( $\pm 0.127$ mm); D - 3.048mm ( $\pm 0.127$ mm).  
 E - 2.515 mm ( $\pm 0.127$ mm).  
 Length: A - 19.761 mm ( $\pm 0.051$ mm).

#### Syringe Parameters:

1. Total scale divisions: major 10 mL, minor 1 mL
2. Scale length: There is a mark every 5 mL from 0 to 60 mL.
3. The total length of the scale to the nominal capacity mark is not less than 75 mm
4. Capacity scale tolerance: ( $\pm 5\%$ )
5. Syringe plunger activation force: 300 KPa ( $\pm 5\%$ )
6. Maximum length that the plunger extends from syringe: 115 mm (this is not a controlled measurement)
7. Tip length: 8.75 mm
8. Diameter of the central hole of the tip- not less than 1.2 mm.
9. Maximum "dead space": 0.20 mL
10. Volume 60 mL ( $\pm 5\%$ )
11. The smallest outer diameter of the syringe tapered tip (cone):  
3,925 - 3,990 mm
12. The largest outer diameter of the syringe tapered tip (cone):  
4,270 - 4,315 mm
13. The length of the outer tapered tip (cone) - not less than 7,500 mm
14. The length of the inner tapered tip (cone) is not less than 7,500 mm
15. External diameter of the lock - 8,0 ( $\pm 0,1$ ) mm
16. The inside diameter of the lock is -7.0 ( $\pm 0.1$ ) mm
17. The maximum internal diameter of the external thread -  
6.73 mm
18. External diameter of the external thread - 7.83 (+ 0 / -0.10) mm
19. Taper of 6%

#### INTENDED USE/INDICATION FOR USE:

Used for temporary control or reduction of postpartum bleeding when conservative management is warranted

#### CONTRAINDICATIONS:

- Arterial bleeding requiring surgical exploration or angiographic embolization
- Cases indicating hysterectomy
- Pregnancy
- Cervical cancer
- Purulent infections in the vagina, cervix, or uterus

- Untreated uterine anomaly
- Disseminated intravascular coagulation
- A surgical site that would prohibit the device from effectively controlling bleeding

#### **POTENTIAL ADVERSE EVENTS:**

- Fitting separation
- Balloon rupture/separation
- Improper drainage

#### **COMPLICATIONS:**

- Tissue damage/trauma from balloon placement
- Uterine rupture from balloon over-inflation
- Concealed bleeding from obstructed drainage
- Bleeding continuing from inadequate tamponade
- Loss of tamponade from balloon rupture

#### **INSTRUCTIONS FOR USE**

**IMPORTANT:** Prior to transvaginal or transabdominal placement of the catheter balloon Bakri Postpartum Balloon, the uterus should be free of all placental fragments, and patient should be evaluated to ensure that there are no lacerations or trauma to the genital tract and that the source of the bleeding is not arterial.

##### **Transvaginal Placement**

1. Determine uterine volume by direct examination or ultrasound examination.
2. Insert the balloon portion of the catheter into the uterus, making certain that the entire balloon is inserted past the cervical canal and internal ostium.
3. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.

##### **Transabdominal Placement, Post-Cesarean Section**

1. Determine uterine volume by direct examination.
2. From above, via access of the cesarean incision, pass the tamponade balloon, inflation adapter first, through the uterus and cervix.  
**NOTE:** Remove the one-way stopcock to aid in placement and reattach prior to filling balloon.
3. Have an assistant pull the shaft of the balloon through the vaginal canal until the dilated balloon base comes into contact with the internal cervical ostium.
4. Close the incision per normal procedure, taking care to avoid puncturing the balloon while suturing.  
**NOTE:** Ensure that all product components are intact and the hysterotomy is securely sutured prior to inflating the balloon. If clinically relevant, the abdomen may

remain open upon inflation of the balloon to closely monitor uterine distention and confirm the hysterotomy closure.

**NOTE:** If clinically relevant, a B-Lynch compression suture may be used in conjunction with the Bakri Postpartum Balloon.

#### **Balloon Inflation.**

##### **With Syringe**

**WARNING:** Always inflate the balloon with liquid. Never inflate with air, carbon dioxide or any other gas.

**WARNING:** The maximum inflation is 500 ml. Do not overinflate the balloon. Over inflation of the balloon may result in the balloon being displaced into the vagina.

**NOTE:** To ensure that the balloon is filled to the desired volume, it is recommended that the predetermined volume of fluid be placed in a separate container, rather than relying on a syringe count to verify the amount of fluid that has been instilled into the balloon.

1. Place an indwelling urinary bladder Foley catheter at this time, if not already in place, to collect and monitor urine output.
2. Using the enclosed syringe, begin filling the balloon to the predetermined volume through the one-way stopcock.
3. Once the balloon has been inflated to the predetermined volume, confirm placement via ultrasound.

**NOTE: See Fig. 1 for proper placement.**

4. If desired, traction can be applied to the balloon shaft. In order to maintain tension, secure the balloon shaft to the patient's leg or attach to a weight, not to exceed 500 grams.

**NOTE:** To prevent displacement of the balloon into the vagina, counterpressure can be applied by packing the vaginal canal with iodine- or antibiotic-soaked vaginal gauze.

5. Connect the drainage port to a fluid collection bag monitor hemostasis.

**NOTE:** To adequately monitor hemostasis, the balloon drainage port and tubing may be flushed clear of clots with sterile isotonic saline.

6. Monitor the patient continuously for signs of increased bleeding and uterine cramping.

#### **Balloon Inflation**

##### **With Rapid Instillation Components**

See Figs. 2-8,

**NOTE:** Ultrasound should be used to confirm proper placement of the balloon once the balloon is inflated to the predetermined volume.


#### **Balloon Removal**

**NOTE:** The timing of balloon removal should be determined by the attending clinician upon evaluation of the patient once bleeding has been controlled and the patient has been stabilized. The balloon may be removed sooner upon the clinician's determination of hemostasis. The maximum indwell time is 24 hours.

1. Remove tension from the balloon shaft.
2. Remove any vaginal packing.


3. Using an appropriate syringe, aspirate the contents of the balloon until fully deflated. The fluid may be removed incrementally to allow periodic observation of the patient.

**NOTE:** In an emergent situation, the catheter shaft may be cut to facilitate more rapid deflation.

**Improper Placement**

patient for signs of bleeding.


**Proper Placement**

ently  
retract  
the  
balloo  
n from  
the  
uterus  
and  
vagina  
l canal  
and  
discar  
d.

4. G

5. M

onitor

## MATERIALS

J-SOS-100500, J-SOSR-100500

| Device Component                | Name of Material                                                                                                                                        | Patient Contacting |
|---------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| <b>Bakri Postpartum Balloon</b> | <b>Shaft:</b><br>Silicone: Dow Corning Q7-4780                                                                                                          | Direct Contact     |
|                                 | <b>Adapter:</b><br>Silicone: Dow Corning Q7-4850                                                                                                        |                    |
|                                 | <b>Adhesive Adapter:</b><br>Nusil Med-2000                                                                                                              |                    |
|                                 | <b>Adapter for Inflation:</b><br>Polyamide 6 № 40340                                                                                                    |                    |
|                                 | <b>Balloon:</b><br>Silicone: LR 3003/20                                                                                                                 |                    |
|                                 | <b>Adhesive of Balloon:</b><br>Nusil Med2-4213                                                                                                          |                    |
|                                 | <b>One-Way Stopcock</b>                                                                                                                                 |                    |
|                                 | <b>Handle:</b><br>High Density Polyethylene:<br>Purell GA7760<br>Blue Colorant: PR-149                                                                  |                    |
|                                 | <b>Body:</b><br>Polycarbonate № 80020                                                                                                                   |                    |
| <b>Syringe, volume 60 ml</b>    | <b>Barrel:</b><br>Polycarbonate № 80020<br><b>Plunger:</b><br>Polystyrene HYPAK PR3ML PSTYR<br><b>Stopper:</b><br>Chlorobutyl Natural Rubber BD STOPPER | Indirect Contact   |

#### J-SOSR-100500

| Device Component       | Name of Material                                                                              | Patient Contacting |
|------------------------|-----------------------------------------------------------------------------------------------|--------------------|
| <b>Connecting Tube</b> | <b>Tubing:</b><br>Vinyl Non-Radiopaque: №11160                                                | Indirect Contact   |
|                        | <b>Non-Vented Spike:</b><br>Acrylonitrile Butadiene Styrene:<br>Estadiene 1640 AS Bianco 1306 |                    |
|                        | <b>Adhesive:</b><br>Loctite 4011                                                              |                    |
|                        | <b>Vented Cap for Spike:</b><br>Low Density Polyethylene : 9002-88-4 Riblene                  |                    |
|                        | <b>Adapter of the Connecting Tube:</b><br>Polyamide 66<br>Red Clariant 00058497               |                    |
|                        | <b>Dual Check Valve:</b><br>Polycarbonate: 103598-77-2                                        |                    |
|                        | <b>Swivel Nut:</b><br>Eastar Copolyester: DN003 Natural                                       |                    |
|                        | <b>Rigid Snap-On Marker:</b>                                                                  |                    |

Polyvinylchloride PVC 7020  
Red Colorant: № RS-7

PACKING INFORMATION.

| Component | Name of Material                                                                                                                                                   | Patient Contacting |
|-----------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| Packaging | Paper packaging part:<br>medical paper TYVEK: 1059B<br><br>Tray: Polietilentereftalat glikol W/GLIDEX<br><br>Film part of the package:<br><br>Polyethylene: Mangar | Indirect Contact   |

PACKING INFORMATION.

Surface density from 50 to 80 ( $\pm 2$ ) g/m<sup>2</sup>

Tensile strength - Tensile strength not less than 0.4 kN/m

Strength of the welded joint - Strength of the weld at least 0.12 MPa

Width of welded joint - Width of adhesive layer not less than 8 mm

**PERFORMANCE CHARACTERISTICS**

Products mentioned in this Technical data summary are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.15 mm.

Products mentioned in this Technical data summary, in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s<sup>2</sup> (10g), and duration of shock acceleration of 16 ms.

Products mentioned in this Technical data summary are stable at temperatures in the range of +32 to + 42 ° C.

Products mentioned in this Technical data summary, in packaging are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C.

Products mentioned in this Technical data summary are stable under sudden changes in temperature.

Products mentioned in this Technical data summary are supplied in sterile packages, do not attempt to clean, disinfect or re-sterilize.

**COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS**

Cook Incorporated manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC.

Cook Incorporated manufactures all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, *Quality System Regulations*.

All of Cook facilities are also certified to ISO 13485 (2003/2012) *Quality Management System*.

### **WARRANTY**

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SHAG".

Russia, 119002, Moscow, Karmanitskiy per., 9.office 501a.

phone: (495) 956-13-09;

Fax: (495) 956-1310,

e-mail: info@schag.ru

Manufacturer's warranty binding on quality products.

The manufacturer guarantees the operability and preservation of the sterility of the devices during the entire shelf life (3 years from the date of manufacture indicated on the package), under the storage conditions.

These devices are intended for single use and are sterile if the package is not opened or damaged. Do not use the device if there is any doubt that it is sterile.

### **DELIVERY INFORMATION**

The products are delivered in separate individual packages with instructions for use. They are placed in corrugated cardboard boxes, on which on the sides or along the perimeter on the agricultural control tape.

### **PACKAGING PARAMETERES**

| <b>RPN</b>    | <b>Height, cm</b> | <b>Width, cm</b> | <b>Weight, g</b> |
|---------------|-------------------|------------------|------------------|
| Tolerance     | ±10%              | ±10%             | ±10%             |
| J-SOS-100500  | 50.8              | 20.3             | 183              |
| J-SOSR-100500 | 50.8              | 20.3             | 216              |

### **SHELF LIFE**

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Bakri Postpartum Balloon have a three (3) year shelf life.

## **STERILIZATION DETAILS**

Cook Incorporated sterilizes devices using a validated ethylene oxide (EtO) method to a SAL (sterility assurance level) of 10<sup>-6</sup>. EtO was the sterilization method selected for these devices as it is effective and does not alter the physical properties of the device. EtO product sterilization is addressed in Cook Incorporated's Quality System Procedure, QSP18 (Sterilization) and Cook Incorporated is certified as compliant with ISO 13485:2003 and EN ISO 11135-1:2007.

## **STORAGE RULES**

The catheter balloon Bakri Postpartum Balloon should be stored in a dark, dry, cool place; avoid extended exposure to light.

## **TRANSPORTATION RULES**

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

## **USAGE CONDITIONS**

The catheter balloon Bakri Postpartum Balloon are used in the patient care medical institutions.

## **POTENTIAL CONSUMERS**

The device is intended for use only by qualified specialists who have been trained in procedures for stopping postpartum uterine bleeding.

## **PACKAGING**

### **J-SOS-100500: is supplied in individual packaging:**

The Catheter Balloon Bakri Postpartum Balloon, Diameter 24 Fr, Length of 54 cm -1p.  
One-Way Stopcock -1p.  
Syringe Volume 60 ml-1p.  
Instruction for Use-1p.

### **J-SOSR-100500: is supplied in individual packaging:**

The Catheter Balloon Bakri Postpartum Balloon, Diameter 24 Fr, Length of 54 cm-1p.  
One-Way Stopcock -1p.  
Syringe Volume 60 ml. -1p.  
Dual Check Valve -1p.  
Connecting Tube, Diameter – 14.0 Fr, length 180 cm-1p.  
Instruction for Use -1p.

## **TRANSPORTATION LABELLING**

The catheter balloon Bakri Postpartum Balloon is supplied in the individual packages which are placed in corrugated boxes with sealed with control tape. Transport boxes show the following information: dimensions of the box, unique bar-code, Cook Medical logo and the following symbols:

- symbol for "Fragile"

- symbol for "temperature limitation"

- symbol for "Keep dry"

### **LABELLING**

The catheter balloon Bakri Postpartum Balloon is supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

- company-manufacturer name and manufacturing address;
- address of the EU representative;
- marking "For specialists only";
- medical device name;
- part number;
- serial number (LOT);
- parameters and general device drawing;
- symbol about CE certification;
- symbol on sterilization;
- symbol on reusability;
- manufacturing date;
- symbol for keep away from sunlight;
- symbol for keep dry;
- bar-code;
- shelf life;
- symbol for "Consult Instruction for Use"
- symbol for caution "Attention, see Instructions for Use"

| Product | Label |
|---------|-------|
|         |       |

**Catheter  
Balloon Bakri  
Postpartum  
Balloon**  
(J-SOS-100500)

**Bakri Postpartum Balloon**  
100% Silicone Contains No Latex

(REF) J-SOS-100500

(REF) G30673

STERILE EO 25 10 2005 25 10 2007 Rx only

|             |              |
|-------------|--------------|
| COOK        | J-SOS-100500 |
| REF: G30673 | LOT: H322245 |
| COOK        | J-SOS-100500 |
| REF: G30673 | LOT: H322245 |
| COOK        | J-SOS-100500 |
| REF: G30673 | LOT: H322245 |

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#### REPAIR/MAINTENANCE


The catheter balloon Bakri Postpartum Balloon are single use only, **do not attempt to repair, carry out maintenance or reuse.**

### DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with Council Directive 93/42/EEC

For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms N2.1.7.2790-10 of December 12, 2010 "Sanitary-epidemiological requirements for the treatment of medical waste". Risk class of medical waste: class B.

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