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Abbott Ireland
Diagnostics Division,
Finisklin Business Park,
Sligo,
Ireland.

Tel. +353 43 3331000
Fax +353 43 3331001

I, Mary Heffernan, Supervisor Project Management, Abbott Ireland Diagnostics Division, Finisklin Business Park, Sligo, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Insert for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: Mary Heffernan
Mary Heffernan,
Supervisor Project Management,
Abbott Ireland Diagnostics Division,

Date: 25 Apr 19

SLIGO CHAMBER OF COMMERCE (INC)
Finisklin,
CERTIFIED CORRECT IN ALL GOOD FAITH
Sligo,
Ireland

Abbott Ireland Diagnostics Division
Finisklin Business Park
Sligo
County Sligo
Ireland

Cortisol

REF 08P3301

G71179R02

S8P330

Alinity i

Cortisol Calibrators

Read Highlighted Changes: Revised March 2018.

NAME

Alinity i Cortisol Calibrators (also referred to as Cortisol Cals)

INTENDED USE

The Alinity i Cortisol Calibrators are for the calibration of the Alinity i analyzer when used for the quantitative determination of cortisol in human serum, plasma or urine.

For additional information, refer to the Alinity i Cortisol reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

The **CAL A** - **CAL F** contain human serum. The **CAL B** - **CAL F** contain purified cortisol. Preservatives: sodium azide and ProClin 950.

The calibrators are at the following concentrations:

Calibrator	Quantity	Cortisol CONC	
		(µg/dL)	(nmol/L)
CAL A	1 x 3.0 mL	0.0	0.0
CAL B	1 x 3.0 mL	3.0	82.8
CAL C	1 x 3.0 mL	5.4	149.0
CAL D	1 x 3.0 mL	10.7	295.2
CAL E	1 x 3.0 mL	25.2	695.3
CAL F	1 x 3.0 mL	59.8	1649.9

MATERIALS REQUIRED BUT NOT PROVIDED

- 04R1001 Alinity ci-series Calibrator/Control Replacement Caps

STANDARDIZATION

The Alinity i Cortisol Calibrators are matched to internal reference standards, which are manufactured gravimetrically using cortisol (USP Reference Standard). The concentration values for these internal reference standards are assigned using Liquid Chromatography Mass Spectrometry/Mass Spectrometry (LC-MS/MS). The LC-MS/MS calibration is verified by BCR 192 and BCR 193 certified reference materials.

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use
- **Rx ONLY**

Safety Precautions

-  **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that this product and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The calibrators contain human-sourced materials that is non-reactive for HBSAg, HIV-1Ag or HIV-1 RNA, anti-HIV-1/HIV-2 and anti-HCV.

The following warnings and precautions apply to: **CAL A** - **CAL F**

WARNING	Contains methylisothiazolones and sodium azide.
H317	May cause an allergic skin reaction.
EUH032	Contact with acids liberates very toxic gas.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- Thaw completely at room temperature (15-30°C) for 45-60 minutes.
- Prior to each use, mix thoroughly by gentle inversion (5-10 times).

STORAGE

- This product is shipped on dry ice.
- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	-10°C or colder	Until expiration date	
Opened	2 to 8°C	90 days after thawing	After each use, immediately return the thawed calibrators to refrigerated storage (2 to 8°C). Store tightly capped with new replacement cap.

The analyzer will track In-use Stability, which is the time the calibrator is outside of refrigerated storage while on the analyzer. The analyzer will not allow the use of the calibrator if the In-use Stability has been exceeded. Maximum In-use Stability can be found in the

Abbott

Assay Parameter Report. For additional information on calibrator In-use Stability, refer to the Alinity ci-series Operations Manual, Section 5.

For additional information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

INSTRUMENT PROCEDURE

- Calibrator lots may be configured using the bar code label on the calibrator carton.
- For information on configuring calibrator data, refer to the Alinity ci-series Operations Manual, Section 2.
- For instructions on ordering and loading calibrators on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

QUALITY CONTROL PROCEDURES

A single sample of each control level must be tested to evaluate the assay calibration. Ensure that assay control values are within the ranges specified in the respective control package insert.

For information on ordering controls, refer to the Alinity ci-series Operations Manual, Section 5.

Once a calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Daily quality control results are outside of statistically-based quality control limits used to monitor and control system performance, as described in the Quality Control Procedures section of the associated reagent package insert.
- If statistically-based quality control limits are not available, then the calibration should not exceed a 30-day limit for recalibration frequency.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

For additional information, refer to the assay reagent package insert and the Alinity ci-series Operations Manual.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, if calibration does not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria, or if controls do not meet the appropriate criteria.

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

ISO 15223 Symbols	
	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Upper limit of temperature
	Use by/Expiration date
	In Vitro Diagnostic Medical Device
	Lot Number
	List Number
	Serial number
Other Symbols	
	After thaw store at
	Calibrator (A,B,C,D,E or F)
	Control Number
	Concentration
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Distributed in the USA by
	Information needed for United States of America only
	Produced for Abbott by
	Product of USA
	For use by or on the order of a physician only (applicable to USA classification only).
	Until first use store at

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

 Abbott Ireland
Diagnostics Division
Finisklin Business Park
Sligo
Ireland
+353-71-9171712

 PRODUCED FOR ABBOTT BY
Fisher Diagnostics,
A Div. of Fisher Scientific Company, LLC
A Part of Thermo Fisher Scientific, Inc.
8365 Valley Pike
Middletown VA 22645 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised March 2018.
©2016, 2018 Abbott Laboratories

 DISTRIBUTED IN THE USA BY
Abbott Laboratories
Abbott Park, IL 60064 USA

ADDENDUM

to the Instruction for Use "Alinity i Cortisol Calibrators" (for Russia)

SHIPPING AND STORAGE CONDITIONS

Recommended shipping conditions for device is - 10°C or colder in single-use thermo containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in autorefrigerators using thermal indicators. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier. Storage of Medical device at a temperature - 10 °C should be carried out in cold rooms or in refrigerators providing a regulated temperature regime with a daily temperature record.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 18 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required

REPAIRS

Not required

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste". Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

ISO 9001, ISO 13485, ISO 13485 (including CMDCAS)

POTENTIAL USER

The Product "Alinity i Cortisol Calibrators" is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals, diagnostic centers, reference laboratories and blood transfusion facilities around the world, whose personnel have been trained in the use of the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are absent

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 8th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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PRODUCED FOR ABBOTT BY

Fisher Diagnostics,
A Div. of Fisher Scientific Company, LLC
A Part of Thermo Fisher Scientific, Inc.
8365 Valley Pike
Middletown VA 22645 USA

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SLIGO CHAMBER OF COMMERCE (INC)

CERTIFIED CORRECT IN ALL GOOD FAITH

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