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Abbott

Abbott Ireland Diagnostics
Division,
Lisnamuck,
Longford
Co. Longford
Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

I, Nora Fitzgerald, Regulatory Affairs Specialist, Abbott Ireland Diagnostics Division, Longford, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: Nora Fitzgerald

Date: 03-DEC-2018

Nora Fitzgerald
Regulatory Affairs Specialist
Abbott Ireland Diagnostics Division
Lisnamuck, Longford
Co Longford
Ireland

SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH
Gerelene Courtney 04/12/18 Authorised Official

Alinity i

Total T₄ Controls

en

Total T₄

REF 07P9510

G71330R03

C7P950

Read Highlighted Changes: Revised February 2018.

NAME

Alinity i Total T₄ Controls (also referred to as Total T₄ Ctrl's)

INTENDED USE

The Alinity i Total T₄ Controls are for the verification of the accuracy and precision of the Alinity i analyzer when used for the quantitative determination of thyroxine (Total T₄) in human serum and plasma.

For additional information, refer to the Alinity i Total T₄ reagent package insert and the Alinity ci-series Operations Manual.

CONTENTS

The **CONTROL L**, **CONTROL M**, and **CONTROL H** contain T₄ prepared in human serum. Preservative: sodium azide.

The following concentration ranges may be used for individual replicate control specifications on the Alinity i analyzer:

Control	Quantity	Total T ₄ CONC		RANGE	
		(µg/dL)	(nmol/L)	(µg/dL)	(nmol/L)
CONTROL L	1 x 8.0 mL	4.2	54.1	2.73 - 5.43	35.14 - 69.88
CONTROL M	1 x 8.0 mL	7.4	95.2	5.44 - 9.92	70.01 - 127.67
CONTROL H	1 x 8.0 mL	14.6	187.9	9.93 - 19.71	127.80 - 253.67

NOTE: The insert ranges for the controls are not lot specific and represent the total range of values which may be generated throughout the life of the product. It is recommended that each laboratory establish its own means and acceptable ranges which should fall within the package insert ranges. Sources of variation that can be expected include:

- Calibration
- Control lot
- Reagent lot
- Calibrator lot
- Instrument

PRECAUTIONS

- **IVD**
- For In Vitro Diagnostic Use
- **Rx ONLY**

Safety Precautions

- **CAUTION:** This product contains human-sourced and/or potentially infectious components. Refer to the CONTENTS section of this package insert. No known test method can offer complete assurance that products derived from human sources or inactivated microorganisms will not transmit infection. Therefore, all human-sourced materials should be considered potentially infectious. It is recommended that this product and human specimens be handled in accordance with the OSHA Standard on Bloodborne Pathogens, Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁻⁴
- The human-sourced material used in the controls is nonreactive for HBsAg, HIV-1 RNA or HIV-1 Ag, anti-HCV, and anti-HIV-1/HIV-2.

The following warnings and precautions apply to: CONTROL L , CONTROL M , and CONTROL H	
Contains sodium azide.	
EUH032	Contact with acids liberates very toxic gas.
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

PREPARATION FOR USE

- This product is liquid ready-to-use.
- This product may be used immediately after removal from 2 to 8°C storage.
- Prior to each use, mix by gentle inversion.

STORAGE

- Do not use past expiration date.

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	
Opened	2 to 8°C	Until expiration date	Store tightly capped. Return to refrigerated storage after use.

INSTRUMENT PROCEDURE

- To obtain the recommended volume requirements for the controls, hold the bottle vertically, and dispense 4 drops of the low control, 4 drops of the medium control, and 4 drops of the high control into each sample cup in the assigned position.
- For instructions on ordering and loading controls on the instrument, refer to the Alinity ci-series Operations Manual, Section 5.

INDICATIONS OF INSTABILITY OR DETERIORATION

Instability or deterioration should be suspected if there are precipitates, visible signs of leakage, turbidity, or if controls do not meet the appropriate package insert and/or Alinity ci-series Operations Manual criteria.

Abbott

BIBLIOGRAPHY

1. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
2. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
3. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
4. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

Key to Symbols

ISO 15223 Symbols	
	Caution
	Consult instructions for use
	Manufacturer
	Temperature limitation
	Use by/Expiration date
	In Vitro Diagnostic Medical Device
	Lot Number
	List Number
Other Symbols	
	Control Number
	Concentration
	Contains Sodium Azide. Contact with acids liberates very toxic gas.
	Control Low, Medium, High (L,M,H)
	Distributed in the USA by
	Information needed for United States of America only
	Product of Ireland
	Range
	For use by or on the order of a physician only (applicable to USA classification only).

Alinity is a trademark of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000

DISTRIBUTED IN THE USA BY

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative or find country-specific contact information on www.abbottdiagnostics.com

Revised February 2018.

©2016, 2018 Abbott Laboratories

ADDENDUM

to the Instruction for Use Alinity i Total T4 Controls (for Russia)

SHIPPING

Recommended shipping conditions for device is ambient. Device is stable until the expiration date indicated on the package. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 22 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements (GOST R 51088-2013).

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2790-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Waste".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, 1272/2008/EC (CLP) as amended, 1907/2006/EC (REACH) as amended, 80/181/EEC, 2009/3/EC, Regulation (EC) 1069/2009, Regulation (EC) 142/2011, CLASS 2013 (as amended), GHS (as amended), WHSR Part 7.1 (as amended), EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641, EN 13612, EN ISO 15223-1, EN 980, EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23640, EN ISO 20776-1, EN 13975, EN ISO 17511, ISO 80000-1

POTENTIAL USER

Professional user

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative

STORAGE

Unopened and opened cartridges / bottles are recommended to be stored at a temperature of 2°C to 8°C in refrigerated chambers or in refrigerators providing a regulated temperature regime.

PRECISION

Not applicable.

TRACEABILITY / STANDARDIZATION

The calibrators are matched to an Abbott internal reference standard. This internal reference standard is manufactured by gravimetric methods using USP reference L-Thyroxine.

MATERIALS REQUIRED BUT NOT PROVIDED

- Alinity i Total T4 assay file
- 07P9520 / 07P9530 Alinity i Total T4 Reagent Kit
- 07P9501 Alinity i Total T4 Calibrators
- Alinity Trigger Solution
- Alinity Pre-Trigger Solution
- Alinity i-series Concentrated Wash Buffer

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- US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
- US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
- World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
- Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.

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