
STATE OF MINNESOTA)
)ss.
COUNTY OF HENNEPIN)

I, Deborah L Shaffer, Notary for the State of Minnesota, hereby declare this 11th
day of October, 2016, that the attached Enteer Re-Entry Catheter IFU MS-
501160-001, is a true and complete copy of the original document.

Signature of Notary Public

Enteer[™]

Re-entry Catheter

CE
0120

English

DEVICE DESCRIPTION

The Enteer[™] Re-entry Catheter facilitates the placement, support, steering and controlling of a guidewire into discrete regions of the peripheral vasculature through the central guidewire lumen or through one of two side-ports (identified by radiopaque markers). The side-ports communicate with the central guidewire lumen and facilitate guidewire steering (at an angle to the central lumen) by allowing the guidewire to exit the Enteer Catheter. The Enteer Catheter contains a small balloon and two radiopaque markerbands on the distal end of a flexible shaft. The radiopaque markers are used for fluoroscopic orientation and to identify guidewire exit ports in the balloon. The distal portion of the catheter has a hydrophilic coating. The Enteer Catheter is compatible with $\geq 0.066"$ (1.68 mm) ID support catheters and may be used with guidewires $\leq 0.018"$ (0.46 mm). All components of the Enteer Catheter are provided sterile and non-pyrogenic. For dimensional information refer to label.

INDICATIONS

The Enteer Catheter is indicated for directing, steering, controlling and supporting a guidewire in order to access discrete regions of the peripheral vasculature.

CONTRAINDICATIONS

The Enteer Catheter is not indicated for use in the cerebral vasculature.

WARNINGS

- Only physicians thoroughly trained in interventional procedures should use the Enteer Catheter.
- The Enteer Catheter is supplied **STERILE** for single use only. Do not reprocess or re-sterilize. Reprocessing or re-sterilizing may increase the risk of patient infection and risk of compromised device performance.
- Do not use the Enteer Catheter if package is open and sterile barrier is broken.
- To reduce the potential for vessel damage, the inflated balloon should only be used in vessels that have a diameter (as measured proximal and distal to the inflation site) equal to or greater than the maximum balloon width.
- Balloon pressure should not exceed the rated burst pressure. Use of a pressure monitoring device is recommended to prevent over pressurization.
- Use only the recommended balloon inflation medium.
- Do not advance or retract the Enteer Catheter unless the balloon is fully deflated under vacuum.
- Pay special attention when directing a guidewire out a side port. Never push, withdraw or torque components that meet resistance. Device kinking/breaking or vessel damage may occur. Fluoroscopy should always be used to aid device manipulation.

PRECAUTIONS

- The Enteer Catheter should be stored in a cool and dry place.
- Use prior to the "Use By" date.
- Before insertion of the catheter, administer appropriate anticoagulant and vasodilator therapy.
- The Enteer Catheter should be handled with care. Prior to use and during the procedure, inspect the packaging and catheter for bends, kinks, or other damage. Discontinue use if the catheter becomes damaged.
- Do not expose to organic solvents.
- Always use wire support for advancement, manipulations, and withdrawal of the Enteer Catheter.
- The Enteer Catheter should only be manipulated under fluoroscopic observation.

POTENTIAL ADVERSE EFFECTS

Potential adverse effects include, but are not limited to, the following:

- Hemorrhage or hematoma
- Vessel trauma that may require further intervention or surgical repair
- Artery spasm
- Embolism
- Stroke
- Neurological deficit
- Drug reactions, allergic reaction to contrast media
- Infection

SUPPLIES

The following supplies are not provided and need to be available and prepped prior to use of the Enteer Catheter:

- Introducer sheath
- Support catheter (guide catheter or guiding sheath)
- Angiographic imaging supplies (i.e. radiopaque contrast, manifold, tubing, etc.)
- Sterile heparinized saline
- Guidewire
- Other supplies as needed to complete the procedure

PRODUCT PREPARATION

1. Inspect the device package prior to opening. Do not use if package is opened or damaged.
2. Remove the Enteer Catheter from the sterile packaging.
3. Remove the balloon protector and stylet.
4. Flush the guidewire lumen with saline.
5. Ensure all device surfaces are exposed to sterile heparinized saline.
6. Prepare a 1:1 mixture of contrast medium and sterile saline.
7. Aspirate 3 to 4 cc (ml) of contrast mixture into a 20 cc (ml) or larger syringe and purge air from the syringe barrel.
8. Attach the syringe to a female luer lock on a 3-way stopcock and fill the stopcock lumens with contrast medium mixture.
9. Attach the 3-way stopcock and syringe to the Enteer Catheter inflation port.
10. Hold the syringe with the nozzle pointing downward and aspirate for 5 seconds.
11. Close the Enteer Catheter inflation port with the 3-way stopcock, remove the syringe, and evacuate air from the barrel.
12. Reconnect the syringe to the 3-way stopcock and purge the 3-way stopcock.
13. Open the stopcock, aspirate until bubbles no longer appear, release the plunger, and disconnect the syringe from the 3-way stopcock.
14. Prepare a balloon inflation device according to the manufacturer's instructions.

DIRECTIONS FOR USE

- Access the vascular system and place the appropriate size introducer sheath per manufacturer's instructions.
- Engage the vessel to be treated with an appropriate guiding catheter per manufacturer's instructions. Connect a hemostatic Y-adapter to the support catheter hub.
- Assure the Enteer Catheter has been prepped per instructions.
- Remove an appropriate guidewire from its package and inspect for damage.
- Advance the guidewire to the target location using fluoroscopic guidance.
- Back-load the Enteer Catheter onto the proximal portion of the guidewire.
- Loosen the knob on the hemostatic Y-adapter.
- Track the Enteer Catheter over the wire and into the proximal support catheter lumen.
- Tighten the hemostatic Y-adapter to create a seal around the Enteer Catheter. Do not over-tighten.

ev3 Inc.,
4600 Nathan Lane North
Plymouth, MN 55442, USA
Tel: +1 763 398 7000
www.ev3.net

Bobby Emann
Associate RA Specialist

10. Advance the Enteer Catheter over the guidewire using fluoroscopic guidance. Position the Enteer Catheter ports (using the radiopaque markers) at the target location.
11. Exchange the guidewire for a shaped guidewire (~45° bend at the distal 1.5 mm segment).
12. Connect the balloon inflation device to the 3-way stopcock and purge the stopcock with contrast medium mixture.
13. Inflate the balloon to 3-4 atm (304-405 kPa).
14. Adjust the fluoroscopic view until the balloon is shown at its minimum width.
15. Steer the shaped guidewire through the desired port by slowly rotating the shaped guidewire immediately proximal of the marker band until the shaped guidewire begins to exit the port.
16. Assess proper guidewire position and direction.
17. If the guidewire position or direction is incorrect, target the other side port by steering the shaped guidewire through the desired port by slowly rotating the shaped guidewire immediately proximal of the marker band until the shaped guidewire begins to exit the port.
18. Advance the shaped guidewire through the port as needed.
19. After completing the procedure, deflate the balloon, withdraw the Enteer Catheter and guidewire. Discard according to hospital standards and national legal requirements.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, ev3 Inc. has no control over the conditions under which this product is used. ev3 Inc. therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ev3 Inc. shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind ev3 Inc. to any representation or warranty with respect to the product. The exclusions and limitations set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

DEFINITION OF SYMBOLS

	Manufacturer		Sterilized using ethylene oxide
	Catalogue Number		Batch Code
	Authorized representative in the European community		For prescription use only
	Keep Dry		Keep away from sunlight
	Use By		Do not reuse
	Do not use if package is damaged		Consult instructions for use
	Rated Burst Pressure		Nominal Pressure
	Telephone		Facsimile
	Package Quantity		

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
STATE OF MINNESOTA)
)ss.
COUNTY OF HENNEPIN)

I, Deborah L Shaffer, Notary for the State of Minnesota, hereby declare this 11th
day of October, 2016, that the attached IFU Viance MS-501158-
001d 19FEB14, is a true and complete copy of the original document.

Signature of Notary Public

Viance™
Crossing Catheter

English

DEVICE DESCRIPTION

The Viance™ Crossing Catheter facilitates the placement and support of a guidewire in discrete regions of the peripheral vasculature through its central guidewire lumen. The Viance Catheter contains a rounded distal tip (1 mm diameter) mounted to a flexible and irremovable proximal shaft. The Viance Catheter is packaged with a dedicated torque device positioned at its proximal end. The distal end of the catheter is hydrophilic coated. The Viance Catheter is compatible with $\geq 0.066"$ (1.68 mm) ID support catheters and may be used with guidewires $\leq 0.014"$ (0.36 mm). All components of the Viance Catheter are provided sterile and non-pyrogenic. For dimensional information refer to label. The Viance Catheter comes in two shaft configurations: flexible and standard stiffness.

INDICATIONS

The Viance Catheter is intended for use with a guidewire to access discrete regions of the peripheral vasculature.

CONTRAINDICATIONS

Do not use with guidewire extension systems with a coupling profile larger than 0.014" (0.36 mm) diameter (i.e. wave pattern coupling mechanism).

WARNINGS

- Only physicians thoroughly trained in interventional procedures should use the Viance Catheter.
- The Viance Catheter is supplied **STERILE** for single use only. Do not reprocess or re-sterilize. Reprocessing or re-sterilizing may increase the risk of patient infection and risk of compromised device performance.
- Do not use the Viance Catheter if package is open and sterile barrier is broken.
- To reduce the potential for vessel damage, the Viance Catheter should only be used in vessels that are ≥ 1 mm in diameter.
- Always use the included torque device during catheter advancement and manipulation especially during device rotation. Failure to use the included torque device may result in catheter failure and may result in patient injury.

PRECAUTIONS

- The Viance Catheter should be stored in a cool and dry place.
- Use prior to the "Use By" date.
- Before insertion, if the Viance Catheter, administer appropriate anticoagulant and vasodilator therapy.
- The Viance Catheter should be handled with care. Prior to use and during the procedure, inspect the packaging and catheter for bends, kinks, or other damage. Discontinue use if the catheter becomes damaged.
- Do not expose to organic solvents.
- The Viance Catheter should only be manipulated under fluoroscopic observation.
- Do not use a syringe smaller than 5 cc when flushing the catheter.
- Do not use an inflation or power assist device when flushing the catheter.

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- Vessel trauma requiring surgical repair or intervention
- Hemorrhage or hematoma
- Artery spasm
- Embolism
- Stroke
- Neurological deficit
- Drug reactions, allergic reaction to contrast media
- Infection

SUPPLIES

The following supplies are not provided with the device and need to be available and prepped prior to use of the Viance Catheter:

- Support catheter (guide catheter or guiding sheath)
- Angiographic imaging supplies (i.e. radiopaque contrast, manifold tubing, etc.)
- Sterile heparinized saline
- Guidewire
- Other supplies as needed to complete the procedure

PRODUCT PREPARATION

1. Inspect the device package prior to opening. Do not use if package is opened or damaged.
2. Remove the Viance Catheter from the sterile package.
3. Flush the guidewire lumen with sterile heparinized saline.
4. Ensure all device surfaces are exposed to sterile heparinized saline.
5. Gently tighten the black knob of the Viance Catheter torque device onto the proximal shaft at a position near the proximal hub.

DIRECTIONS FOR USE

1. Access the vascular system and place the appropriate size introducer sheath per manufacturer's instructions.
2. Engage the vessel to be treated with an appropriate guiding catheter per manufacturer's instructions. Connect a hemostatic Y-adapter to the guiding catheter hub.
3. Assure the Viance Catheter has been prepped per instructions.
4. Retrieve an appropriate guidewire from its package and inspect for damage.
5. Advance the guidewire to the target location using fluoroscopic guidance.
6. Back-load the Viance Catheter onto the proximal portion of the guidewire.

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
STATE OF MINNESOTA)
)ss.
COUNTY OF HENNEPIN)

I, Deborah L Shaffer, Notary for the State of Minnesota,
hereby declare this 11th day of October, 2016, that the
attached IFU MS-501062-001f TrailBlazer is an Original
Document.

Signature of Notary Public

CE
0120

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If resistance is encountered at any time during the procedure, do not force passage. Resistance may be due to device or lumen. Carefully withdraw the support catheter. This catheter should be used by physicians with adequate training in the use of the device. Catheter manipulation should only occur under fluoroscopy.

PRECAUTIONS

Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician. Carefully inspect the package and catheter prior to use to verify no damage occurred during shipment. Do not use catheter if it is damaged or compromised. Do not wipe down the catheter surface with dry gauze. Do not attempt to pass the ev3 device through a smaller sized sheath introducer than indicated on the label. Refer to product labeling. Only use guidewire of the recommended diameter and length. To avoid kinking, advance the catheter slowly, in small increments until the proximal end of the guidewire emerges from the catheter.

The catheter should not be advanced into a vessel having a diameter smaller than the catheter outer diameter. After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and applicable local, state, and federal laws and regulations.

POTENTIAL COMPLICATIONS /

ADVERSE EFFECTS

The following complications may result, but may not be limited to:

- Access site complications
- Arterial dissection
- Arterial spasm
- Arterial thrombosis
- Catheter fracture with tip separation and distal embolization
- Death
- Distal embolization (air, blood clots or plaque)
- Intimal disruption
- Local or systemic infection including sepsis
- Perforation and vessel rupture
- Surgical intervention

HOW SUPPLIED

The ev3 device is supplied STERILE for single use only.

STORAGE

The ev3 device should be stored in cool, dry, and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation or ultraviolet light. Rotate inventory so that device is used prior to the Use By date on the package label.

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- 10 cc syringe filled with sterile heparinized saline
- Three way stop cock
- Appropriately sized exchange guidewire (refer to product labeling)
- Appropriately sized hemostatic introducer sheath (refer to product labeling)

Model	Length	Outer Diameter	Inner Diameter	Weight	Material
ev3-100	100 cm	0.08"	0.04"	1.1 g	Polyurethane
ev3-150	150 cm	0.08"	0.04"	1.6 g	Polyurethane
ev3-200	200 cm	0.08"	0.04"	2.1 g	Polyurethane
ev3-250	250 cm	0.08"	0.04"	2.6 g	Polyurethane
ev3-300	300 cm	0.08"	0.04"	3.1 g	Polyurethane
ev3-350	350 cm	0.08"	0.04"	3.6 g	Polyurethane
ev3-400	400 cm	0.08"	0.04"	4.1 g	Polyurethane
ev3-450	450 cm	0.08"	0.04"	4.6 g	Polyurethane
ev3-500	500 cm	0.08"	0.04"	5.1 g	Polyurethane
ev3-550	550 cm	0.08"	0.04"	5.6 g	Polyurethane
ev3-600	600 cm	0.08"	0.04"	6.1 g	Polyurethane
ev3-650	650 cm	0.08"	0.04"	6.6 g	Polyurethane
ev3-700	700 cm	0.08"	0.04"	7.1 g	Polyurethane
ev3-750	750 cm	0.08"	0.04"	7.6 g	Polyurethane
ev3-800	800 cm	0.08"	0.04"	8.1 g	Polyurethane
ev3-850	850 cm	0.08"	0.04"	8.6 g	Polyurethane
ev3-900	900 cm	0.08"	0.04"	9.1 g	Polyurethane
ev3-950	950 cm	0.08"	0.04"	9.6 g	Polyurethane
ev3-1000	1000 cm	0.08"	0.04"	10.1 g	Polyurethane

ev3-1000 is supplied with contrast/sterile saline mix

DISPOSAL

CAUTION: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, ev3, Inc. has no control over the conditions under which this product is used. ev3, Inc. therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. ev3, Inc. shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind ev3, Inc. to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, waive any mandatory provisions of applicable law. If any part of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law, the entire Disclaimer of Warranty shall not be invalidated. In any competent jurisdiction, the validity of the exclusions and limitations of this Disclaimer of Warranty shall not be affected. No oral agreements or understandings shall be construed and no person has any authority to bind ev3, Inc. to any representation or warranty with respect to the product.

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INDIVIDUAL ACKNOWLEDGMENT

State/Commonwealth of Minnesota }
 County of Hennepin } ss.

On this the 19th day of February, 2018, before me,
Day Month Year

Deborah L. Shaffer, the undersigned Notary Public,

Name of Notary Public

personally appeared

Bobby Omann

Name(s) of Signer(s)

☒ personally known to me – OR –

☐ proved to me on the basis of satisfactory evidence

to be the person(s) whose name(s) is/are subscribed to the within instrument, and acknowledged to me that he/she/they executed the same for the purposes therein stated.

WITNESS my hand and official seal.

Place Notary Seal and/or Stamp Above

[Signature]
Signature of Notary Public

Plymouth, MN

Other Required Information (Printed Name of Notary, Residence, etc.)

OPTIONAL

Although the information in this section is not required by law, it may prove valuable to persons relying on the document and could prevent fraudulent removal and reattachment of this form to another document.

Description of Any Attached Document

Title or Type of Document: IFU, Trailblazer Angled

Document Date: _____ Number of Pages: _____

Signer(s) Other Than Named Above _____

Right Thumbprint
of Signer

Top of thumb here

ev3 Inc.,
4600 Nathan Lane North
Plymouth, MN 55442, USA
Tel: +1 763 398 7000
www.ev3.net

Feb 19, 2018

Belen

Bobby Omann

Medtronic

TrailBlazer™ Angled

Support Catheter

Angled Support Catheter with Hydrophilic Coating

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
Figure 1. TrailBlazer™ angled support catheter

Instructions For Use

Device Description

The TrailBlazer™ angled support catheter is an over-the-wire (OTW), single lumen catheter with a tapered, angled tip. The catheter system is compatible with 0.36 mm (0.014 in), 0.46 mm (0.018 in), and 0.89 mm (0.035 in) guidewires. Depending on the model, the support catheter has a working length of 65 cm, 90 cm, 135 cm, or 150 cm. All models are compatible with a 4 Fr introducer sheath. For specific sizing, see Table 1 for model specifications.

As shown in Figure 1, the lumen begins (2) at the proximal end of the manifold (1) and continues through the catheter (3) to the angled distal tip (4). The lumen is used to pass the catheter over a guidewire. The guidewire-compatible diameter and working length are marked on the strain relief (5). The distal catheter shaft has three radiopaque markers (6) that aid in positioning the catheter. The first 40 cm at the distal end of the catheter (7) are coated with a hydrophilic coating.

Indications For Use

The TrailBlazer™ angled support catheters are percutaneous, single lumen catheters designed for use in the peripheral vascular system. The TrailBlazer™ angled support catheters are intended to guide and support a guidewire during access of the vasculature, to allow for wire exchanges and provide a conduit for the delivery of saline solutions or diagnostic contrast agents.

Warnings

- The TrailBlazer™ catheter should be used by physicians experienced with interventional techniques in the vascular system.
- This device is supplied STERILE for single use only. Do not reprocess or resterilize. Reprocessing and resterilizing could increase the risk of patient infection and risk of compromised device performance.
- If resistance is encountered at any time during the insertion procedure, do not force passage. Forcing passage during the insertion procedure can damage the device or the lumen. If this happens, carefully withdraw the support catheter.
- Do not torque the catheter if it is or seems stuck. Torquing with excessive force can cause damage to the device or the vessel.
- Always use fluoroscopy or other visualization when manipulating the catheter.

Precautions

- Federal (U.S.A.) Law restricts this device to sale by or on the order of a physician.
- Carefully inspect the package and catheter before use to verify that no damage occurred during shipment. Do not use the catheter if it is damaged or compromised.
- Do not wipe down the catheter surface with dry gauze.
- Do not attempt to pass the TrailBlazer™ support catheter through a smaller-sized introducer sheath than indicated on the label. See Table 1 for introducer sheath size specifications.
- Ensure that the guidewire diameter and length are compatible with the TrailBlazer support catheter. See Table 1 for guidewire size specifications.
- To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
- Do not advance the catheter into a vessel with a diameter smaller than the catheter's outer diameter.
- Do not exceed the maximum recommended infusion pressure of 300 psi.
- After use, this product is a potential biohazard. Handle and dispose of the device in accordance with accepted medical practice and applicable local, state, and national laws and regulations.

Potential Complications / Adverse Events

Potential complications or adverse events include, but are not limited to:

- Access site complications
- Arterial dissection
- Arterial spasm
- Arterial thrombosis
- Catheter fracture with tip separation and distal embolization
- Death
- Distal embolization (air, blood clots, or plaque)
- Intimal disruption
- Local or systemic infection including sepsis
- Perforation and vessel rupture
- Stroke/TIA
- Surgical intervention

How Supplied

Warning: This device is supplied STERILE for single use only. Do not reprocess or resterilize.

Reprocessing and resterilizing could increase the risk of patient infection and risk of compromised device performance.

Table 1. Sizing guide

MODEL	ASC-014-090	ASC-014-135	ASC-014-150	ASC-018-090	ASC-018-135	ASC-018-150	ASC-035-065	ASC-035-090	ASC-035-135	ASC-035-150
Maximum guidewire diameter mm / in	0.36 mm / 0.014 in	0.36 mm / 0.014 in	0.36 mm / 0.014 in	0.46 mm / 0.018 in	0.46 mm / 0.018 in	0.46 mm / 0.018 in	0.89 mm / 0.035 in	0.89 mm / 0.035 in	0.89 mm / 0.035 in	0.89 mm / 0.035 in
Working length cm	90	135	150	90	135	150	65	90	135	150
Minimum guidewire length cm	150	180	180	150	180	180	135	160	200	225
Marker band spacing cm	15	15	15	15	15	15	15	15	15	15
Shaft outside diameter mm / in	0.76 mm / 0.030 in	0.76 mm / 0.030 in	0.76 mm / 0.030 in	0.86 mm / 0.034 in	0.86 mm / 0.034 in	0.86 mm / 0.034 in	1.40 mm / 0.055 in	1.40 mm / 0.055 in	1.40 mm / 0.055 in	1.40 mm / 0.055 in
Tip outside diameter mm / in	0.51 mm / 0.020 in	0.51 mm / 0.020 in	0.51 mm / 0.020 in	0.58 mm / 0.023 in	0.58 mm / 0.023 in	0.58 mm / 0.023 in	1.04 mm / 0.041 in	1.04 mm / 0.041 in	1.04 mm / 0.041 in	1.04 mm / 0.041 in
Minimum introducer sheath diameter Fr / mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm	4 Fr / 1.40 mm

Instructions for Use

English

Storage

Store the TrailBlazer™ support catheter in a cool, dry, and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the TrailBlazer support catheter is used before the Use by date on the package label.

Recommended Items

Prepare the following items using sterile technique:

- 10 cc syringe filled with sterile heparinized saline
- Three-way stopcock
- Appropriately sized guidewire (refer to Table 1 for guidewire size specifications)
- Appropriately sized introducer sheath (refer to Table 1 for introducer sheath size specifications)

Preparation

Caution: Carefully inspect the package and catheter before use to verify that no damage occurred during shipment. Do not use the catheter if it is damaged or compromised.

1. Carefully remove the catheter from the inner package.
2. The TrailBlazer support catheter is coated with a hydrophilic coating. Before inserting the catheter, activate the coating by immersing the catheter in normal saline, for approximately 30 to 60 seconds, or wiping down the catheter shaft with a saturated gauze sponge.
Caution: Do not wipe down the catheter surface with dry gauze.
3. To displace air, connect a 10 cc syringe filled with heparinized saline to the luer lock on the manifold and flush out the lumen.

Directions For Use

1. Introduce the TrailBlazer support catheter percutaneously using an appropriately sized introducer sheath.

Caution: Do not attempt to pass the TrailBlazer support catheter through a smaller-sized introducer sheath than indicated on the label. Refer to Table 1 for introducer sheath sizing specifications.

2. Place the prepared catheter over a previously inserted guidewire.
Caution: Ensure that the guidewire diameter and length are compatible with the TrailBlazer support catheter. See Table 1 for guidewire size specifications.
Caution: To avoid kinking, advance the catheter slowly, in small increments, until the proximal end of the guidewire emerges from the catheter.
3. Carefully advance the TrailBlazer support catheter to the desired vasculature, using fluoroscopic imaging to visualize the marker bands.

Warning: If resistance is encountered at any time during the insertion procedure, do not force passage. Forcing passage during the insertion procedure can damage the device or the lumen. If this happens, carefully withdraw the support catheter.

Warning: Do not torque the catheter if it is  seems stuck. Torquing with excessive force can cause damage to the device or the vessel.

Warning: Always use fluoroscopy or other visualization when manipulating the catheter.

Caution: Do not advance the catheter into a vessel with a diameter smaller than the catheter's outer diameter.

4. Carefully remove the TrailBlazer support catheter, while maintaining guidewire position.

Note: To use the catheter for infusion, withdraw the guidewire and attach an infusion device (syringe, power injector, and so forth). Do not exceed 300 psi inlet infusion pressure.

Caution: Do not exceed the maximum recommended infusion pressure of 300 psi.

Table 2. TrailBlazer™ angled support catheter flow rates (ml/second)

Model	Guidewire diameter	Length	Saline		Contrast*	
			150 psi	300 psi	150 psi	300 psi
ASC-014-090	0.36 mm 0.014 in	90 cm	0.9	1.3	0.2	0.5
ASC-014-135	0.36 mm 0.014 in	135 cm	0.8	1.1	0.2	0.3
ASC-014-150	0.36 mm 0.014 in	150 cm	0.8	1.0	0.1	0.3
ASC-018-090	0.46 mm 0.018 in	90 cm	1.6	2.1	0.5	1.1
ASC-018-135	0.46 mm 0.018 in	135 cm	1.3	1.8	0.3	0.8
ASC-018-150	0.46 mm 0.018 in	150 cm	1.2	1.7	0.3	0.7
ASC-035-065	0.89 mm 0.035 in	65 cm	7.8	11.2	6.1	8.5
ASC-035-090	0.89 mm 0.035 in	90 cm	7.1	10.2	5.4	8.1
ASC-035-135	0.89 mm 0.035 in	135 cm	5.9	8.3	4.2	7.3
ASC-035-150	0.89 mm 0.035 in	150 cm	5.5	7.7	3.8	7.0

*A mix of 75% Isovue 370 contrast media and 25% saline equivalent solution


Disposal

Caution: After use, this product is a potential biohazard. Handle and dispose of the device in accordance with accepted medical practice and applicable local, state, and national laws and regulations.

Warranty Disclaimer

Although this product has been manufactured under carefully controlled conditions, Medtronic has no control over the conditions under which this product is used. Medtronic therefore disclaims all warranties, both express and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. Medtronic shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind Medtronic to any representation or warranty with respect to the product.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held to be illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

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©2013 Covidien.
Enteer is a trademark of a Covidien company.
Viance is a trademark of a Covidien company.