Registered No. 2018 - 4251

NOTARIAL CERTIFICATE

MOON JUNG HWAN NOTARY OFFICE

20, Hwangsaeul-ro 312beon-gil, Bundang-gu,
Seongnam-si, Gyeonggi-do Korea

210mm X 297mm

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
101501

User manual

Medical LCD monitor Model name: IPH-MED 2703, including component

«APPROVE»

COJE Co. Ltd., Rep. of Korea

(Position)

President

(Name)

Young Se Oh Young Oh

(Signature)

Young Se Oh

Stamp

Medical LCD monitor, model name: IPH-MED 2703

450

Name of medical device and component to be delivered

Medical LCD monitor model: IPH-MED 2703, component:

- 1) Medical LCD monitor model : IPH-MED 2703 : 1
- 2) Power cable – 1
- 3) DVI cable – 1
- 4) User manual – 1
- 5) Monitor base – 1 (if necessary)

Medical device manufacturer

COJE Co., Ltd (COJE Co., Ltd);
47, Gyeongdae-ro 17-gil, Buk-gu, Daegu, Republic of Korea
Tel : +82-70-4659-2139; +82-53-255-7029;
Email: nahee.ryu@COJEdisplays.com

Purpose of medical device

Medical LCD monitor model, IPH-MED 2703 (hereinafter called 'Monitor') equipped with components is used with an endoscope, light source, camera control device, and ultrasonic center for video diagnosis

Terms of use for medical device

A monitor is used only by medical institutions for the purpose of treatment and prevention.
A monitor is designed to be used by a doctor or under his supervision.

Medical LCD monitor, model name: IPH-MED 2703

2

Ys

Application scope

A monitor is used within a medical imaging system along with an endoscope, light source, video information center, camera control device and ultrasonic center.

Aftereffect

The aftereffect directly related to this product has not been known up to now.

Side effect

There is no side effect occurred when storing, transporting and using the product according to safety regulations specified in the user manual.

Risk in use of medical device

A risk management evaluation was carried out and documentation performed in a risk management report and category of this plan including the followings.

1. Risk analysis for each product
2. Characteristics affecting purpose and safety of product
3. Potential hazard
4. Potential clinical hazards and serious aftereffects
5. Analysis for the proportion of clinical risk and usability

The type and aftereffect analysis required to evaluate the risk of a medical device has been developed in order to define usability for the purpose.

Hazard related to use, structure and procedure was defined within the scope of procedure. Hazard of risk level not allowed was reduced to the allowable level through change of composition, test performance and mutually accepted actions.

Now, the manufacturer uses seriousness evaluation procedure for clinical hazard which is classified from minor risk to major risk by stage. When the risk level not allowed is found, actions for reducing risk are required.

All hazards and risks found based on analysis were reduced to allowable levels.


Requirements for risk management shall comply with major requirement of 93/42/EEC directive.

450

Safety regulations

The relevant symbol is used to notify the possibility of life loss or death and property loss to customers

Use accompanied power cable to connect to a outlet complying with your nation's standard.

There shall be no scratch over the power cable. Damaged cable can cause a fire or an electric shock. When pulling out a cord from an outlet, use a plug instead of a wire.

Do not place a monitor or other heavy things above the power cable. A fire or a electric shock may occur due to cable damage.

Do not place other thing in the main body. The thing in the main body can cause damage in the product or an accident.

Y50

Do not place a monitor under rain. Keep away from humidity or dust in order to prevent a fire or an electric shock.

Do not place a monitor in a very bright place or a place to which direct sunlight or other light source is lying on. The light source can be reflected to the monitor screen. Place a monitor below the eye level.

Do not place a thing over the monitor. Penetration fluid can cause a fire or an electric shock.


Since the screen is easy to break, do not give impact onto it or make a scratch.

Do not touch an electric plug with wet hands.

Description of medical device

1. Packaging :

Monitor

User manual



Power cable

DVI cable

A monitor base is delivered as main component if necessary.

2. Cable connection

1) Power and input/output

DVI (input), VGA(input), SDI (input), S-VIDEO (input), C-VIDEO (input), R/Pb (input), G/Y (input), B/Pr (input), CSYNC/H (input), VSYNC (input)

DVI/VGA (output), RS-232CAC, SDI (output), S-Video (output), C-Video (output), R/Pb (output), G/Y (output), B/Pr (output), HSYNC/CSYNC (output), VSYNC (output)

Power input part, power button

2) Cable connection


1) Connect a power cable to AC power terminal.

450

- 2) Use a suitable cable to connect the monitor to the equipment.
- 3) First, turn on the monitor and connected equipment in order.

3. OSD menu (screen menu) and introduction

Shortcut key

- (SOURCE) Change the source.
- (INFO) Show the source, resolution, frequency/size and gamma value information.
- (OVERSCAN) Activate/deactivate the overscan function.
- (BACKLIGHT) Activate/deactivate the backlight function.
- PIP (Picture In Picture) Turn on and off the function.
- PBP (Picture By Picture) Turn on and off the function.
- (LED/indicator) Turn on/off control button indicator.
- (POWER) Turn on/off the power.

- Scroll the roller for "Introduction" and value "setting".
- Press this for location selection.

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Technical description of medical device

1. Physical specification

Dimension (height x width x length, mm)	653.35 x 402.6 x 86.34 mm (± 2 mm)
Distance between fixed holes	VESA 75 for fixing 4 holes 75 mm (± 1 mm) VESA 100 for fixing 4 holes - 100 mm (± 1 mm)
Weight (net/gross)	8.3 kg / 10.8 kg ($\pm 10\%$)

2. Technical specification

Panel

Size (diagonal)	685.8 mm (27") ± 1 mm
Mattress type	IPS
Backlight type	LED
Resolution	1920 x 1080
Number of pixel	2,073,600
Pixel pitch	0.311 x 0.311 mm
Response time	25(± 5) ms (black > white > black)
Color rendering	1670M color
Brightness	Normal : 650 cd/m ² (+5%), maximum : 700 cd/m ² ($\pm 5\%$)
Contrast ratio	1 000:1 ($\pm 10\%$)
Field of view	178° (horizontal)/178° (vertical) ($\pm 1^\circ$)
Adjustment range :	0 ~ 100
Brightness	0 ~ 100
Contrast	0 ~ 100
Horizontal location level	0 ~ 100
Vertical location level	
Horizontal synchronization	128 kHz
Vertical synchronization	60 Hz

Y2

Input/output

Signal type supported	DVI, D-SUB, composite video signal, S-Video, HD-SDI (compatible with SMPTE259M-C292M), RGB, RGB composite signal, analog type RGB (BNCx5)
Input terminal	DVI, DVI-I, D-SUB, composite video signal, S-Video, HD-SDI (compatible with SMPTE259M-C292M), RGB, RGB composite signal, analog type RGB (BNCx5)
Output terminal	DVI, RGB/Component (BNC x 5), DVI, HD-SDI, S-Video, composite video signal, serial port : RS-232C
Signal voltage : DVI, DVI-I signal D-SUB signal Composite video signal, S-Video signal HD-SDI RGB signal, RGB composite signal	0.5 ~ 1.5 Vp.p. 0.7 ~ 1.5 Vp.p. 0.5 ~ 1.5 Vp.p. 0.5 ~ 1.5 Vp.p. 0.7 ~ 1.5 Vp.p.

Power supply

Power	Alternating current : 100V ~ 240V, Alternating current voltage : $\pm 10\%$, 50/60Hz
Power consumption	Maximum 70W

Structure

User control	OSD control using control button and scroll wheel
--------------	---

Power cable

Length	1500 mm ($\pm 2\%$)
Weight	140 g ($\pm 10\%$)

DVI cable

Length	1500 mm ($\pm 2\%$)
Weight	188 g ($\pm 10\%$)

Monitor base

Height adjustment	0-110 mm (± 5 mm)
Rotation angle	90° (horizontal direction or vertical direction)
Angle of inclination	-3°(down)/15°(up) ($\pm 1^\circ$)
Rotation angle	-20°(left)/20°(right) ($\pm 1^\circ$)

450

Dimension	208 x 208 x 100 mm ($\pm 5\%$)
Weight	360 g ($\pm 10\%$)

**Information on whether medicine, animal material
and/or material coming from human body in a
medical device exist or not**

There is no medicine, animal material and/or material coming from human body in the product.

1010.211

450

Use of medical device

1. Before monitor operating

Inadequacy in use of medical device :

Check the monitor operating status before every use. In any doubt for operating status of monitor products, the corresponding product is considered as inadequate to use and contact the manufacturer or representative in full authority to check the use for the purpose is possible. If a monitor is damaged, do not use the monitor and send it to the repair center if necessary.

If there is a mixture of various medical devices, it may not be technically safe. Therefore, check the compatibility with the manufacturer when connecting it to an external device.

2. Monitor setting

Adjust all monitor settings with the main menu.

2.1. Image

Set image quality including brightness, contrast and definition.

450

2.1.1. Brightness

- Press "Left" button (◀) to reduce brightness.
- Press "Right" button (▶) to increase brightness.

2.1.2. Contrast

- Press "Left" button (◀) to reduce contrast.
- Press "Right" button (▶) to increase contrast.

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213 . Definition

Increase contour definition to set image definition.

- Press "Left" button (◀) to reduce definition.
- Press "Right" button (▶) to increase definition.

214 . Color

Set the color characteristics of a monitor.

- Press "Left" button (◀) to reduce color saturation.
- Press "Right" button (▶) to increase color saturation.

450

215 . Tone

Set a monitor tone.


- Press "Left" button (◀) or "Right" button (▶) to adjust the tone.

2.1.1. Backlight (dimming/brightness)

Set the brightness of backlight in the back liquid crystal.


- Press "Left" button (◀) to reduce backlight brightness.
- Press "Right" button (▶) to increase backlight brightness.

2.2. Color and gamma setting

Manually set these palette parameters as "Green, Red and Blue" other than brightness, contrast and color temperature to save them in the customer profile.


- Press "Left" button (◀) to reduce the value.
- Press "Right" button (▶) to increase the value.

2.3. Image adjustment

Select image mode as entire screen mode or the original resolution.

- Entire screen mode : the image shows up in the entire screen irrespective of resolution.
- Common mode : the image shows up as the original resolution.


If input signal comes in from DVI, D-SUB 15, RGB or SOG, the following menu items can be used additionally :

Y50

- Open/close the overscan.
- Horizontal location : adjust the horizontal location of the image.
- Vertical location : adjust the vertical location of the image.
- Auto - adjustment : adjust the image location automatically.

2.4. Setting transparency, latency time and screen menu location

Adjust transparency, OFF time interval and screen menu location.

If input signal comes in from DVI, D-SUB 15, RGB ~~15/15~~, the following menu items can be used additionally :

- Phase : Delete visual defect in horizontal direction in order to improve image quality.
- Frequency : Delete visual defect in vertical direction in order to improve image quality.

241 . Transparency

450

Set the transparency of screen menu window.

242 . Screen menu latency time

Set the interval time (10/20/30 seconds) at which the screen menu is turned off when there is no keyboard input.

243 . Screen menu location

Set the screen menu location within the display.

480

2.5. PIP (Picture In Picture)/POP (Picture By Picture) function

The monitor additionally shows a video signal coming from additional source. If you use this menu, you can change PIP and POP setting.

251 . Window configuration

Since display area is divided into lots of sections, it is easy to manage various input video signals.

Y50

252 . Video signal source of main window

Set the first video signal source.

253 . Video signal source of the second window

Set the second video signal source.

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254 . Video signal source of main window

Set the first video signal source.

255 . Video signal source of the second window


Set the second video signal source.

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256 . Window switching

Switch between the first video signal source and the second video signal source.

257 . PIP location

Set the location of PIP window.

2.6. Service setting

Change the various additional settings.

2.6.1. Resetting of factory setting (initialization)

Initialize all parameters of the monitor before basic setting.


2.6.2. Working time

Display total working time of backlight.


2.6.3. Working time

Display total working time of backlight.

2.6.4. Display power control


If you do not use a monitor for a certain period of time, it enters sleep mode for saving power.

730

2.6.5. Screen menu blocking

Block the screen menu in order to prevent setting change, which occurs accidentally.

Life of medical device

Life is 5 years.

It is recommended that the device be disposed of in accordance with the manufacturer's requirements after the end of its life.

Environment information

1. Disposal information

It is required that natural resources are mined and used in order to produce products. A product may include materials harmful to health and environment. To prevent alien substance from penetrating into the surrounding environment and reserve natural resources, it is recommended that monitors be sent to electrical and electronic waste collection centers when the monitors have reached end of life. In this case, most of the monitor materials whose service life passed are recycled or processed.

The symbols marked with X in the trash mean that the product's disposal through the electrical and electronic waste collection centers is recommended after use. If you want additional information about the collection, processing and disposal center, contact your relevant city government's waste disposal service. Or if you want additional information about the environment and disposal of our product, please contact your local representative.

Monitors belong to mechanically hazardous wastes and must be disposed of in accordance with SanPiN 2.1.7.2790-10 "Sanitary requirements for medical waste disposal", A grade.

Y 10

2. Requirement for environment protection

When using, transporting and storing monitors, it does not adversely affect the human body or the environment.

3. Cleaning

The screen is made of thin glass with a surface treated with plastic and when it drops and given a shock or under physical influence, it can be damaged. Do not clean the front panel of a monitor using the materials containing ketone (such as acetone), ethyl alcohol, toluene, ethyl acid, methyl or chloride since it brings about the damage of the front panel.

Panel : Clean the panel using a soft wool or non-woven cloth. Use a diluted soft glass cleaning solution.

Main body : Immerse a soft cloth in a soft solution made by mixing detergent and water and use it for cleaning. And then, use a cloth wet with only water to clean it, and wipe it dry with a dry towel.

4. Double packaging

Do not discard boxes and packing materials. It is recommended that the box be stored well for exchange or repair for future use when returning the products. The box is the most suitable for monitor transportation and exchange and it is not easy to find the same packaging.

Safety information

General recommendations

Read the safety regulations and the recommendations for use before you use the monitor.

Keep the documents containing safety regulations and precautions for use for future use. Observe all precautions in the user manual.

Observe all precautions for use.

It is possible to earth surely when connected to only ~~the equipment~~ marked with «Hospital Only» or «Hospital Grade».

	Warning Do not open as it can cause electric shock.	
Warning: Do not open the lid (or hood section) in order to reduce the hazard of electric shock. There is no part which users can repair inside the monitor.		

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Transportation, storage and terms of use

Monitors can be transported in accordance with transportation regulations by freight type through the use of closed transportation means of all types. When transporting the freights, the protection against humidity and unexpected movement is required. When transporting and storing the monitors, it is recommend that the monitors be stored in the original manufacturer's package.

The following product environment conditions are recommended :

Monitor use condition is temperature of 0°C - $+50^{\circ}\text{C}$, relative humidity of air of 10% - 80%, and atmospheric pressure of 700 gPa - 1060 gPa. The monitor stands physical action under the vibration load having frequency of 10 - 55 Hz and movement amplitude of 0.15 mm.

Store the monitor under the influence of climatic factor such as temperature of -20°C - $+60^{\circ}\text{C}$, relative humidity of air of 0% - 90%, and atmospheric pressure of 500 gPa - 1060 gPa.

When transporting monitors, make sure they are under the influence of climatic factor such as temperature of -20°C - $+60^{\circ}\text{C}$, relative humidity of air of 0% - 90%, and atmospheric pressure of 500 gPa - 1060 gPa and maintain frequency of 10 - 55 Hz, vibration load of 0.35 mm for movement amplitude, and impact load of 10g for impact acceleration of 16 ms.

Packaging medical device

Parts and components of a monitor are packaged with polyethylene film and a foam insert is contained in corrugated box. The box cover (flap) made of corrugated board is glued with a polyethylene tape having an adhesive layer.

Y50

Warranty obligations

COJE Co. Ltd. warrants the absence of defects in a product including manufacture materials and quality during 12 months from the date at which a product is delivered to a customer.

Regarding warranty COJE Co. Ltd. hold the responsibility only to defects occurred under the condition that the product is used properly.

COJE Co. Ltd. does not hold any responsibility for defects occurred due to improper storage, transportation, fake product, structure change, use not suitable for the purpose or over-use, use of product or part of other manufacturer.

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Medical device labeling

1. Safety notice

This manual specifies the important information for safety using the following graphic symbols. Read the symbols carefully.

	Power button. Press the button when turning on or off a monitor (standby mode).
	Protective earth (earth).
	Fragile Handle it carefully.
	Alternating current
	Series number
	Indicates the manufacturer
	Indicates a production date
	Follow the user manual
	Dispose monitors according to European directive (WEEE) regarding electric and electronic wastes.
	This product meets the requirement regarding European directive 2002/95/EC «Guideline for use limits on electric and electronic equipment of some dangerous materials».
	This product meets the requirements of Federal Communications Commission (FCC).
	Underwriters Laboratories certified the safety of the corresponding product.
	This product meets RoHS guideline.
	This product holds a certificate of China Compulsory Certification (CCC).
	This product meets the technical specification of Industry Canada (IC) standard.
	This product operates at a height of 2000m above sea level.
	This product complies with European directive.
	Temperature limit
	Humidity limit

Y52

Avoid contact with moisture.

Keep away from heating or radiation.

Avoid direct sunlight.

Shipping packing symbol

Packaging materials are made of processed materials.

Above the product

When stacking the products, please do not stack them in more than five stories.

An electronic device can affect or limit monitor operating, or emit electromagnetic waves which interfere with the normal operation. Install the device at a place where this kind of effect is not applied.

Product characteristics may vary depending on a region. Refer to your local user manual for technical characteristics of the equipment.

150

22

EMC, Manufacturer's notification for radio frequency

Manufacturer's recommendations and notifications – electromagnetic Emission

Emission check	Compliance for standard	Electromagnetic Emissions – recommendations
Radio frequency emissions	Group I	Radio frequency energy is used only for indoor equipment. Therefore, radio frequency emission of this system is so small and does not interfere with operating of nearby electronic equipment.
Radio frequency emissions	Grade B	The equipment can be used in all spaces connected to public low voltage electrical network providing the power to residential buildings.
Harmonic emissions	Grade D	
Voltage fluctuation/1 time emission	Suitable	

Manufacturer's recommendations and notifications – protection against electromagnetic interference

Examination of protection level	Examination level according to EN/IEC 60601 standard	Compliance level	Electromagnetic emissions – recommendations
Electrostatic discharge	± 6 kW (upon contact) ± 8 kW (via air)	± 6 kW (upon contact) ± 8 kW (via air)	The floor shall be made of timber, concrete or ceramic tile. If the floor is made of synthetic materials, the relative humidity shall be more than 30 %.
Drop and surge : - 4	± 2 kW (for power circuit) ± 1 kW (for input and output circuit)	± 2 kW (for power circuit) Absence of input and output circuit	Electrical network shall meet typical requirement of a medical building or clinic.
Rapid voltage rise	± 1 kW (differential connection) ± 2 kW (in-phase connection)	± 1 kW (differential connection) ± 2 kW (in-phase connection)	Electrical network shall meet typical requirement of a building or clinic.

Y50

Examination of protection level	Examination level according to EN/IEC 60601 standard	Compliance level	Electromagnetic Emissions – recommendations
Voltage drop in input circuit, instant outage and voltage fluctuation	$< 5\% U_T$ $(> 95\% \text{ drop } U_T) \text{ during } 0.5 \text{ cycle}$ $40\% U_T$ $(60\% \text{ drop } U_T) \text{ during } 5 \text{ cycles}$ $70\% U_T$ $(30\% \text{ drop } U_T) \text{ during } 25 \text{ cycles}$ $< 5\% U_T$ $(> 95\% \text{ drop } U_T) \text{ for } 5 \text{ seconds}$	$< 5\% U_T$ $(> 95\% \text{ drop } U_T) \text{ during } 0.5 \text{ cycles}$ $40\% U_T$ $(60\% \text{ drop } U_T) \text{ during } 5 \text{ cycles}$ $70\% U_T$ $(30\% \text{ drop } U_T) \text{ during } 25 \text{ cycles}$ $< 5\% U_T$ $(> 95\% \text{ drop } U_T) \text{ for } 5 \text{ seconds}$	Electrical network shall meet typical requirement of a building or clinic. If you have to disconnect the power and carry out the equipment work, it is recommended that you connect uninterruptible power supply or a battery.
Magnetic field of electrical network frequency (50/60 Hz)	3 A/m	3 A/m	Magnetic field level of electrical network frequency shall meet general characteristics of typical building or clinic.
Radio frequency penetration energy	3 V (effective value) 150 kHz - 80 MHz	3 V	The gap between portable wireless communication means and component (including cable) shall be more than geographical gap calculated according to formula of transmitter frequency. Recommended geographical gap : $d = 1,2\sqrt{P}$ (less than 80 MHz)

Y 20

Examination of protection level	Examination level according to EN/IEC 60601 standard	Compliance level	Electromagnetic emissions – recommendations
Radio frequency emission energy	3 V/m 80 MHz – 2.5 GHz	3 V/m	$d = 1,2\sqrt{P}$ (80 MHz - 800 MHz) $d = 2,3\sqrt{P}$ (800 MHz - 2,5 GHz) Where, P – maximum output of transmitter notified by transmitter manufacturer, (W), d – recommended geographical gap (m) The electric field strength (a) of fixed radio frequency emission source is measured when surveying the space, and it shall be less than allowable numerical value (b) of each frequency range. Electronic interference may occur next to the equipment with the following symbols: 

Note. U_1 – voltage of power network before applying a test level

Note 1. In case of 80 and 800 MHz frequency, the requirements for the higher frequency applies.

Note 2. These recommendations are not suitable for all cases. The diffusion of electromagnetic waves are influenced by building elements, objects and people.

a) For example, it is impossible to theoretically predict electric field strength of fixed radio frequency emission sources with enough accuracy for wireless phone (wireless and cell phone) base station, wireless transmitter base station, AM and FM radio broadcasting station and TV broadcasting station. Space survey is required in order to determine electromagnetic emissions level. If the electric field strength measured in a place where the equipment is used exceeds the allowable level recorded as above, observe that an error occurs in the equipment work. For example, if an error occurs such as change in equipment location or movement, an additional action is required.

b) The electric field strength shall be less than 3 V/m in the frequency range of 150 kHz - 80 MHz.

Recommended space gap value between portable (movable) radio frequency communication means and the monitor

Recommended space gap value between portable (movable) radio frequency communication means and the monitor			
If the interference level of a monitor can be managed, the monitor is designed to be used in the electromagnetic environment. If the minimum space gap between portable (movable) radio frequency communication means and the monitor is warranted considering the maximum output of the communication means as recommended as below, a purchaser or a user can be out of the influence of electromagnetic interference for a monitor.			
Rated maximum output of transmitter, P , W	Space gap in accordance with transmitter frequency, d , (m)		
	$d = \left[\frac{35}{V_1} \right] \sqrt{P}$ 150 kHz - 80 MHz per one line	$d = \left[\frac{15}{E_1} \right] \sqrt{P}$ 80 - 800 MHz as per one line	$d = \left[\frac{25}{E_1} \right] \sqrt{P}$ 800 MHz - 2.5 GHz per one line
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
In case of a transmitter of rated maximum output not specified in the table, in order to specify recommended space gap, d , suggest the rated maximum output, P (W) specified in the transmitter manufacturer's document.			
Note 1 The larger electric field strength value is applied to 80 and 800 MHz frequency. 2 The above statement does not apply to all cases. Electromagnetic wave diffusion is influenced by absorption or reflection caused by structure, a thing and people. 3 In case of a transmitter of rated maximum output not specified in the table, in order to specify recommended space gap, d suggest the rated maximum output, P (W) specified in the transmitter manufacturer's document.			

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450

List of international standard met by medical device

EN 1041:2008 «Information supplied by the manufacturer of medical devices»

ISO 13485:2003/Cor.1:2009 «Medical devices. Quality management system. Requirements for the purpose of adjustment»


ISO 14971:2007 «Medical devices. Application of risk management to medical devices»

ISO 15223-1:2012 «Medical devices. Symbols used in supplied information when labeling medical devices. Part 1. Overall requirements/Medical devices. Symbols to be used with medical device labels, labeling and information to be supplied Part 1: General requirements»

EN 60601-1-2:2014 «Electric medical device. Part 1-2. Overall requirement for major safety and main functional characteristics. Parallel standard. Electromagnetic compatibility. Requirements and test»

EN 60601-1:2006 «Electric equipment medical part 1 : overall requirement for major safety and significant productivity»

Directive 93/42/EEC «Medical device, unit and equipment»

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Technical support

Manufacturer

Contact technical support department as below to obtain the answer for technical questions.

+82-70-4659-2139;

+82-53-255-7029

If additional help is requested, contact the manufacturer's local representative.

Manufacturer's representative in full authority

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COJE Co. Ltd., Rep. of Korea-----
CEO/President , Young Su Oh-----
personally appeared before me and
admitted his(her) subscription to the
attached -----
User Manual (Medical LCD monitor)---

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This is hereby attested on this
31st day of Oct. 2018 at this office.

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Belong to Suwon District Prosecutors`
Office

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20, Hwangsaoul-ro 312beon-gil, Bundang-gu,

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Seongnam-si, Gyeonggi-do Korea

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Signature of the Notary Public

MOON JUNG HWAN

This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
November 09, 2013 under Law No.157

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