K0423 / 2019

Direction for use

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Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals

legal manufacturer

«VDW GmbH», Bayerwaldstraße 15, 81737 München, Germany,

Tel: +49 89 62734-0

Regulatory Affairs Manager

18.Feb. 2019

«VDW GmbH»

Kamilla Kwiatkowski Stamp VDW

VDW GmbH
81737 MÜNCHEN, BAYERWALDSTR. 15
TEL. 089 / 6 27 34-0. FAX 089 / 6 27 34 - 100

Direction for use

1.Product subtitle:

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals

2.Product release:

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals of different sizes:

1. Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals in blisters of 6 instruments, size 12,5, length 21mm
2. Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals in blisters of 6 instruments, size 12,5, length 25mm
3. Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals in blisters of 6 instruments, size 12,5, length 31mm
4. Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals in blisters of 4 instruments, size 12,5, length 21mm
5. Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals in blisters of 4 instruments, size 12,5, length 25mm
6. Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals in blisters of 4 instruments, size 12,5, length 31mm

3.Indications for use

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals is a glide path instrument used in reciprocating motion to prepare the root canal system before the shaping with a rotary or a reciprocating instrument.

4.Description

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals will be used mechanically in reciprocating motion with endodontic motors with the RECIPROC settings. The glide path will be obtained with

one single mechanized file in each cases an initial manual scouting should be used.

Following the intended use of this device, the main interaction of endodontic R – PILOT Instruments of nitinol (NiTi) for preparation of root canals with the human body will be a limited contact with the internal part of the root canals that have to be treated: dentin and residual pulp. Residual pulp is defined as the remaining pulp after manual scouting using hand files. It can be assumed that contact with skin or mucosa may happen unintentionally during treatment


Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals are intended as single use (i.e. single patient) devices. In order to prevent sterilization, an ABS (acrylonitrile butadiene styrene) ring, located on the handle, will expand in case of sterilization. The dimension of this ring after sterilization makes the connection of the file with the contra-angle impossible. Therefore, the manufacturer can insure that the files will not be reprocessed

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals used in reciprocating motion to prepare the root canal system before the shaping with a rotary or a reciprocating instrument.



These instruments are to be used only in a clinical or hospital environment by qualified users (general practitioners (dentist) as well as Endo Specialists), following Good Dental Practices.

All patients can be treated with the product.

Photo. Example of the product

Carton box with approximate position of Russian sticker (Front view)

	Carton box (back view)
	Blister with approximate positioning of Russian sticker

5.Application: for dental use only

6.Composition of the working length:

Working length of endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals is made of nitinol (NiTi) type (M – Wire)

Shaft of made of zinc and copper (basic components)

Stopper is made of silicone

Non – autoclavable plastic ring (ABS) is made of acrylonitrile butadiene styrene

Stopper and non – autoclavable plastic ring (ABS) have brown color in order to distinguish endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals from other instruments due to its interim size (12,5)

7. Technical characteristics

Size: 12,5

Working length: 21, 25, 31 mm

Taper: 4%

Shaft diameter: 2,35 mm

Depth marks, tapering and cylindrical part

	Working length 21 mm	Working length 25 mm	Working length 31 mm
Number of depth marks	3	4	5
Depth marks placement			
Depth mark 1	18 +/-0,10 mm measuring from cutting tip (D1 side) to the middle of the depth mark 1	18 +/-0,10 mm measuring from cutting tip (D1 side) to the middle of the depth mark 1	18 +/-0,10 mm measuring from cutting tip (D1 side) to the middle of the depth mark 1
Depth mark 2	1 +/-0,05 mm from the middle of the depth mark 1 till the middle of the depth mark 2	1 +/-0,05 mm from the middle of the depth mark 1 till the middle of the depth mark 2	1 +/-0,05 mm from the middle of the depth mark 1 till the middle of the depth mark 2
Depth mark 3	1 +/-0,05 mm from the middle of the depth mark 2 till the middle of the depth mark 3	1 +/-0,05 mm from the middle of the depth mark 2 till the middle of the depth mark 3	1 +/-0,05 mm from the middle of the depth mark 2 till the middle of the depth mark 3
Depth mark 4	Not applicable	2 +/-0,05 mm from the middle of the depth mark 3 till the middle of the depth mark 4	2 +/-0,05 mm from the middle width of the depth mark 3 till the middle of the depth mark 4
Depth mark 5	Not applicable	Not applicable	2 +/-0,05 mm from the middle of the depth mark 4 till the middle of the depth mark 5
Tapering part	17 +/-0,5 mm	17 +/-0,5 mm	17 +/-0,5 mm
Cylindrical part	4 +/-0,5 mm	8 +/-0,5 mm	14 +/-0,5 mm

8. Indications for Use

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals instrument is a Glide Path instrument used in reciprocating motion to prepare the root canal system before the shaping with a rotary or a reciprocating instrument.

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals instrument may only be used in reciprocating motion with a designated drive system

with original VDW RECIPROCR settings, e.g. VDW.CONNECT Drive or VDW.GOLDR RECIPROC. The use of reciprocating drive systems without original VDW RECIPROCR settings can lead to instrument fracture and misuse.

9. Contraindications

Similar to all machine driven root canal instruments, endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals should not be used in a root canal with an abrupt apical curvature in the apical region which is managed with pre-curved hand files.

10. ⚠ Warnings

None known

11. ▽ Precautions

The instruments may only be used by qualified persons in dentists' offices or clinics.

Use a rubber dam.

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals instruments are sterile, single use instruments, designed to be used in no more than one molar. The instrument cannot be re-sterilized.

Examine the instrument during treatment for signs of wear such as untwisting, bending or damaged cutting edges. In these cases, instruments should be discarded.

Clean the flutes from debris in an interim stand every three pecks and irrigate the canal according to the appropriate irrigation protocol. (INTERIM-STAND with REF No. 495)

Only light pressure should be applied on endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals instruments. If resistance is encountered or if the instrument does not advance, clean the instrument, irrigate the canal according to the appropriate irrigation protocol and check the canal for possible blockage with a C-PILOT file ISO size 08 at minimum or similar hand file. Then use the endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals instrument briefly with a lateral brushing motion.

12. Adverse Events

During the application of equivalent devices, it was found out that there is a possibility of instrument breakage. All risks associated with this case were considered and analyzed in the Risk Management Report. In conclusion, due to low risk rate reported for equivalent devices, Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals can be considered to be safe when used for dental root canal preparation, following Precautions stated in 11. Precautions of this Direction for use.

13. Step – by – step Instructions

The application of endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals is recommended for the following cases:

- Initially to render the existing natural canal path smooth prior to the use of shaping instruments
- During root canal preparation with RECIPROC whenever advancement with RECIPROC becomes difficult (e.g. RECIPROC R25 or RECIPROC blue R25)

Prior to the use of shaping instruments

1. Establish a straight-line coronal and radicular access to the root canal entrance.
2. Use a C-PILOT file ISO size 08 at minimum or similar hand file to ensure patency and determine working length with the help of a preoperative radiograph [mandatory] and/or an apex locator.
3. Afterwards use endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals as described in sections a-f until working length has been reached
 - a) Place irrigant in the access cavity of the root canal.
 - b) Introduce endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals free without binding into the canal. Press the ON/OFF button or foot pedal of the motor when the instrument is at the root canal orifice.
 - c) Move the instrument in a slow in-and-out pecking motion. The amplitude of the in-and-out movements should not exceed 3 mm. Only very light pressure should be applied. It is also possible to brush laterally. The instrument will advance easily in the canal. One in-and-out movement = one peck.
 - d) Clean the instrument in the interim stand after 3 pecks.

e) Irrigate the canal.

f) Repeat steps c-e until full working length has been reached

4. Continue your shaping with the VDW file system of your choice.

14. During root canal preparation with RECIPROC

In some canals the RECIPROC instrument may stop from advancing in the canal or its advancement may become difficult.

1. Use a C-PILOT file ISO size 08 at minimum or similar hand file to determine working length with the help of a preoperative radiograph [mandatory] and/or an apex locator.
2. Afterwards use Endodontic R – PILOT instrument of nitinol (NiTi)) For Preparation of root canals as described in section 8. 3a-f until working length has been reached.

15. Disinfection, cleaning and sterilization and disposal

Endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals is intended for single-use only. Disinfection, cleaning and sterilization procedures are not recommended. Re-use of endodontic R – PILOT instrument of nitinol (NiTi) for preparation of root canals can increase the risk of cross contamination and breakage.

16. Conditions of storage and shelf life

Store at temperature 0°C - +30°C

Humidity 50 – 60%

Away from direct sunlight

Shelf life is 5 years

17. Transportation conditions

Product should be transported in compliance with regulations valid for each type of transport at temperature - 20°C - +60°C, humidity 50 – 60%, away from direct sunlight

18. Environmental protection requirements

The product does not have negative impact for the environment during its life cycle.

19. Application conditions

Instrument only be used by qualified persons in dentists' offices or clinics at a temperature of 0 ° C to + 30 ° C, at a relative humidity of 50 to 60%, away from sunlight. The products are properly used inside the oral cavity, with the presence of biological fluids at the nominal temperatures as from 32 ° to 42 ° C

20. Maintenance and repair

Product is intended for single use only. Product is not a subject for reuse and therefore also not for repair

21. Symbols used

	Batch code
	Use-by date
	Catalog number
	4 instruments in each blister
	6 instruments in each blister
	Do not use if packaging is damaged
 RECIPROC	RECIPROC trade mark
	Do not re-use
	M – Wire nickel – titanium alloy sign
	CE marking
	Sterilized using irradiation
	Manufacturer
Rx only For Dental Use only	For dental use only

23. Waste

Disposal in accordance with sanitary rules and regulations of the Russian Federation SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for handling medical waste". Risk class of medical waste is B.

24. Warranty

VDW products are guaranteed to be free of defects in materials and craftsmanship to the original purchaser.

VDW does not take any responsibility and does not pay compensation for possible damage or accidents caused by:

- Use of product in contravention of its Directions for Use;
- Use of product in contravention of its Indications for Use;
- Use of product in conjunction with a reciprocating drive system (i.e., endodontic motor) without the manufacturer's product-specific settings;
- Accident, abuse, neglect, or normal wear and tear.

In no event shall VDW be liable for any special, direct, indirect, exemplary, punitive or consequential damages or costs. For product quality issues, please contact VDW's Authorized Representative:

In Russian Federation:

Dentsply Sirona LLC

115432 Moscow, prospect Andropova, dom 18, korpus 6, offices 1 39 40

+7 (495) 725 1087

cis@dentsply.ru

25. List of international applied standards:

1	EN 1:2001/AC:2006 556-	Sterilization of medical devices - Requirements for medical devices to be designated "STERILE" - Part 1: Requirements for terminally sterilized medical devices
2	EN 980:2008	Symbols for use in the labelling of medical devices
3	EN 1041:2008	Information supplied by the manufacturer of medical devices
4	EN ISO 15223-1: 2012	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
5	EN ISO 21531:2009	Dentistry - Graphical symbols for dental instruments
6	EN ISO 10993-1:2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
7	EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
8	EN ISO 11137-3: 2006	Sterilization of health care products - Radiation - Part 3: Guidance on dosimetric aspects of development, validation and routine control
9	EN ISO 11137-2:2015	Sterilization of health care products - Radiation - Part 2: Establishing the sterilization dose
10	EN ISO 11137-1: 2015	Sterilization of health care products - Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
11	EN ISO 11607-1:2009 + A1:2014	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
12	EN ISO 11607-2:2006 + A1:2014	Packaging for terminally sterilized medical devices - Part 2: Validation requirements for forming, sealing and assembly processes
13	EN ISO 13485: :2012/AC:2012	Medical devices - Quality management systems - Requirements for regulatory purposes
14	EN ISO 14971:2012	Medical devices - Application of risk management to medical devices
15	IEC 62366:2007	Medical devices - Application of usability engineering to medical devices
16	ISO 1797-2:1992	Dental rotary instruments - Shanks - Part 2: Shanks made of plastics
17	ISO 1797-1:2011	Dentistry - shanks for rotary instruments - Part 1: Shanks made of metals
18	EN ISO 3630-1:2008	Dentistry - Root-canal instruments - Part 1: General requirements and test methods

26. Manufacturer

«VDW GmbH», Bayerwaldstraße 15, 81737 München, Germany,

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info@vdw-dental.com • www.vdw-dental.com

27. Authorized Representative in the territory of Russian Federation:

Dentsply Sirona LLC

115432 Moscow, prospect Andropova, dom 18, korpus 6, offices 1 39 40

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