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APOSTILLE

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1. Country : United States of America
2. This public document has been signed by **CYNTHIA KELSEY**
3. acting in the capacity of **NOTARY PUBLIC**
4. bears the seal/stamp **CYNTHIA KELSEY , NOTARY PUBLIC, BERKS COUNTY, COMMONWEALTH OF PENNSYLVANIA**

Certified

5. at Harrisburg, Pennsylvania
6. The 1st day of June, 2016
7. by Pedro A. Cortés, Secretary of the Commonwealth of Pennsylvania
8. No: 201614705
9. Seal/Stamp
10. Signature

Pedro A. Cortés
Pedro A. Cortés

This Apostille only certifies the authenticity of the signature and the capacity of the person who has signed the public document, and, where appropriate, the identity of the seal or stamp which the public document bears.


This Apostille does not certify the content of the document for which it was issued.

This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

TRUTHFUL AND ACCURATE STATEMENT


To Whom It May Concern:

As a responsible official of Angiotech Puerto Rico, Inc, located at Road 459 Km 0.6, Montana Industrial Park, Aguadilla, Puerto Rico 00603, U.S.A., I hereby certify the information in the attached document to be, to the best of my knowledge, truthful and accurate.


Susan Hauck, Documentation Services
Manager
Surgical Specialties Corporation

Sworn to and subscribed before me this 12th day of May 2016

Cynthia Kelsey
Notary Public

Instruction for Use
Stratafix Spiral PDO (Polydioxanone) Knotless Tissue Closure Device
(Bidirectional)

Stratafix Spiral PDO (Polydioxanone) Knotless Tissue Control Device (hereinafter – the Suture) consists of suture material with barbs and a needle at each end. Barbs allow for tissue approximation without the need to tie surgical knots.

The Stratafix Spiral PDO sutures are made from dyed (violet) polyester, poly(p-dioxanone). The pigment for the dye is D&C Violet No. 2. It is stated that polydioxanone has been found to be nonantigenic, nonpyrogenic and elicits only slight tissue reaction during absorption.

While the formation of barbs in the Stratafix Spiral PDO reduces the tensile strength relative to non – barbed suture material of the same size, tying of knots in non – barbed suture materials also reduce their effective strength. For this reason, the strength of the Stratafix Spiral PDO can be compared to USP knot strength of non – barbed sutures. Additionally, USP designations for diameter are used to describe the Stratafix Spiral PDO prior to barbing, except for minor variation in suture diameter with a maximum overage of 0.1 mm.

ACTIONS

The effectiveness of absorbable suture *in vivo* is characterized by two essential parameters: an ability to retain tensile strength and absorption rate (mass loss rate). The Stratafix Spiral PDO formula allows for minimization of variability of these parameters and wound support during the critical wound healing period and further prolonged healing period.

Findings of studies related to implantation in animals with use of the Stratafix Spiral PDO sutures demonstrate that sutures of a size 3-0 or greater retain over 80% of their original strength after two weeks and four weeks post implantation.

After six weeks following the implantation, sutures retain nearly 40-70% of their tensile strength. The sutures of 4-0 size retain nearly 67% of their tensile strength after two weeks and nearly 50% after four weeks, and nearly 37% after a six week period.

Findings of tests with implantation show that the Stratafix Spiral PDO suture polymer subsequent resorption is minimal for up to approximately 120 days, and is essentially complete within 180 days (six months).

APPLICATIONS.

This medical product is applied only for therapeutic purposes.

PRESCRIPTION

Assignment for compound (cross-linking) of soft tissues.

INDICATIONS

Stratafix Spiral PDO is indicated for use in soft tissue approximation where use of absorbable sutures is appropriate.

CONTRAINDICATIONS

Stratafix Spiral PDO is not to be used where prolonged (beyond six weeks) approximation of tissues under stress is required and is not to be used in conjunction with or for fixation of prosthetic devices (e.g. heart valves or synthetic grafts), that are non-absorbable in nature.

WARNINGS

Do not resterilize. Discard opened, unused Stratafix Spiral PDO and associated surgical needles.

Doctor should be familiar with surgical procedure and techniques involving absorbable suture before employing Stratafix Spiral PDO for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material

used. Physicians should consider the *in vivo* performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

The safety and effectiveness of Stratafix Spiral PDO has not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissues, osseous tissues, tendinous tissue, ophthalmic surgery or for use in microsurgery, therefore this product should not be used for these purposes.

As with any foreign body prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable suture, Stratafix Spiral PDO may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable suture material, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or destination or which may require additional support.

PRECAUTIONS

Stratafix Spiral PDO device contains bidirectionally oriented barbs to anchor tissues and does not require knots to approximate opposing edges of a wound. Tying knots on the barbed section of the material will damage the barbs and potentially reduce the suture tensile strength and barb effectiveness. For the bidirectional forces to be created and for the device to function properly, both sides of the Stratafix Spiral PDO device must be engaged in the tissue. Additionally, when

completing placement, an additional backstitch or bite of tissue lateral to the end of the incision is required to lock the device in place.

Avoid contacting the Stratafix Spiral PDO and associated needles with other materials (e.g. surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction on the needle to disengage it from the barbs.

Care should be taken to avoid damage when handling. Avoid crushing or crimping the suture material with surgical instruments, such as needle holders and forceps. Do not attempt to remove memory in the polymer by running fingers down the suture material as this can damage the barbs.

Infections, erythema, foreign body reactions, transient inflammatory, reactions and in rare instances wound dehiscence are typical or foreseeable risks associated with any suture and hence are also potential complications associated with Stratafix Spiral PDO.

When using Stratafix Spiral PDO device comprised of dyed PDO subcutaneously, the device should be placed as deeply as possible in order to minimize erythema and induration normally associated with absorption. Acceptable surgical practice should be followed with respect to drainage and Closure of infected wounds.

To avoid damaging needle points and swage areas, grasp the needle in an area one-third ($1/3$) to one-half ($1/2$) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "sharps" containers.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds.

ADVERSE REACTIONS

Adverse effects associated with the use of this device may include wound dehiscence, failure to provide adequate wound support in closer of the sites where expansion, stretching or distention occur, failure to provide adequate wound support in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, suture extrusion and delayed absorption in tissue with poor blood supply, calculus formation in urinary and biliary tracts when prolonged contact with salt solutions occurs, and transitory local infection at the wound site. Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood borne pathogens.

Due to prolonged suture absorption, some irritation and bleeding may occur.

APPLICATION

Use as required per surgical procedure. Deeper layers of the wound are to be closed in such a fashion as to take tension off of the dermal edges.

To deploy Stratafix Spiral PDO device, the first end should be pulled through the tissue until the center transition zone has reached the tissue. Estimation of the center can be aided by taking a single bite of tissue, then aligning the two needles together. Taking one arm of the Stratafix Spiral PDO device, at least two loose bites of tissue should be taken. This is then repeated with the other arm of the Stratafix Spiral PDO device. Once at least four bites, or passes, through the tissue have been completed, each of the strands can then be grasped and the tissue can be

approximated to the desired tension. Bites, or passes through the tissue, can be taken to approximate the wound. Care should be taken to utilize Stratafix Spiral PDO device on the barbed segments only. Do not attempt to approximate wounds using the non-barbed segments, as the barbs are required for successful wound approximation with Stratafix Spiral PDO device.

For Stratafix Spiral PDO device with straight needles, Stratafix Spiral PDO device can be deployed at right angles to the wound. Starting at one edge of the wound, the Keith needle is deployed at right angles to the wound along the desired tissue plane and the needle is advanced through the tissue until the desired amount of tissue has been engaged, then the needle can exit at the desired location remote from the wound edge. The same procedure is repeated with the second arm of the Stratafix Spiral PDO device. After successful deployment of both arms of the Stratafix Spiral PDO device and needle exits, the center transition zone of the Stratafix Spiral PDO device should be roughly in the center of the open wound. The ends of the Stratafix Spiral PDO device where they exit remote to the wound can now be grasped and the tissue can be pushed together to facilitate the soft tissue approximation to the desired tension.

To remove the entire Stratafix Spiral PDO device comprised of dyed PDO: Cut the Stratafix Spiral PDO device at the midline, between the opposing barbed segments, and then pull the lateral ends to remove the device.

STERILITY

The Stratafix Spiral PDO device is sterilized by ethylene oxide gas. Do not restirilize. Do not use if package is opened or damage. Discard opened unused sutures. Do not use after expiration date.

HOW SUPPLIED

The Stratafix Spiral PDO is available sterile, dyed violet in various lengths and with a USP diameter size 2 through 4-0 (metric size 5 through 1,5). Each Stratafix Spiral PDO is double armed with various needles sizes.

MANUFACTURER'S WARRANTY

Angiotech Puerto Rico, Inc. undertakes a warranty obligation to end users in relation to the medical device, and ensures that it has no defects and complies with the stated manufacturer's characteristics within the entire shelf life of the device, provided that transportation, storage and operation conditions are observed.

DISPOSAL

Dispose medical institution in accordance with the requirements 2.1.7.2790-10 SanPiN "Sanitary-of the requirements for the management of medical waste." Products after the intended use are potentially infectious wastes and uchi lysed as waste Class B products with expired shelf life, and packed in damage to refer to medical waste class A. Shelf life - 5 years.

MANUFACTURER

Angiotech Puerto Rico, Inc.

Road 459 Km 0.6 Montana Industrial Park Aguadilla, 00603 Puerto Rico, USA

REPRESENTATIVE IN RUSSIAN FEDERATION

"Johnson & Johnson" LLC

121614, Krylatskaya str., 17 bld. 2, Moscow, Russia

Phone: +7 (495) 588 77 77, Fax: +7 (495) 588 78 78.

SYMBOL USED FOR LABELING

	Device Item Number
	Lot Number
	Manufacturer
	Best before
	Sterilized by ethylene oxide
	The CE Label including the SGS (authorized body) identification number, the Product complies with the essential requirements of the European Directive 93/42/EEC on Medical Devices
	Qualified persons only
	Do not reuse product
	Do not resterilize product
	Temperature Range for Storage
	Temperature Range for Transportation

	Attention, address shipping documents
	Do not use if package is damaged

Supplied:

Stratafix Spiral PDO (Polydioxanone) Knotless Tissue Closure Device, in configurations – 12 шт. (each thread in an individual sterile packaging);

Unstruction for using – 1 pcs;

Secondary package – 1 pcs

REQUIREMENTS FOR MEDICAL DEVICE

Наименование	Значение
Suture Material Structure	Monofilament
Suture Material Color	Violet
Coating	Uncoated
Suture Material Size Metric Size (USP) Diameter	M1(5/0) – 0,100–0,149 мм M1.5(4/0) – 0,150–0,199 мм M2(3/0) – 0,200–0,249 мм M3(2/0) – 0,300–0,349 мм M3.5(0) – 0,350–0,399 мм M4(1) – 0,400–0,449 мм

Наименование	Значение
Simple Knot Breaking Strength, at least (thread for the specified value)	M1.5(4/0) – 6,8 H M2(3/0) – 9,5 H M3(2/0) – 21,0 H M3.5(0) – 26,8 H M4(1) – 39,0 H
Barbed Segment Length	7 x 7 cm; 10 x 10 cm; 14 x 14 cm; 24 x 24 cm; 30 x 30 cm; 36 x 36 cm permissible variation $\pm 2\%$
Nominal Length of Suture Material	7x7 cm – nom. length 21.2 cm; 10 x 10 cm – nom. length 22.8 cm; 14x14 cm – nom. length 35.2 cm; 24x24 cm – nom. length 55.2 cm; 30x30 cm – nom. length 67.2 cm; 36 x 36 cm – nom. length 79.2 cm; permissible variation $\pm 2\%$
Atraumatic Needle Suture Holding Strength, at least (thread for the specified value)	M1(5/0) – 1,7 H M1.5(4/0) – 2,3 H M2(3/0) – 4,5 H M3(2/0) – 9,0 H M3.5(0) – 11,0 H M4(1) – 15,0 H
Distance between Barbs	0.61 or 0.95 mm
Barb Length	0.34 or 0.58 mm
Barb Height	0.082(± 0.031) mm

Наименование	Значение
Barb Angle	15°(±5°)
Needle Curvature	Needle Type: CTX, CT-3, SH, MH, CT-1, MO-4 – 1/2 curved needle Needle Type: OS-8, FS, CP-2, OS-6, FS-2, TAPERCUT V-4, TAPERCUT V-26, PS-2 prime – 3/8 curved needle
Needle Point Sectional Shape	Needle Type: CTX, CT-3, SH, MH, CT-1, MO-4 – round cross-section Needle Type: OS-8, FS, CP-2, OS-6, FS-2, TAPERCUT V-4, TAPERCUT V-26, PS-2 prime – triangular cross-section
Needle Hardness	HV 4900-6475 N/mm ²
Needle Corrosion Resistance	corrosion resistant
Needle Attachment Segment Diameter	M1(5/0) – 0.35 mm M1.5(4/0) – 0.40 mm M2(3/0) – 0.50 mm M3(2/0) – 0.50 mm M3.5(0) – 0.60 mm M4(1) – 0.70 mm
Needle Resilience Under test force of 5N, residual deformation, not exceeding	M1(5/0) – 0.60 mm M1.5(4/0) – 0.60 mm M2(3/0) – 0.40 mm M3(2/0) – 0.40 mm M3.5(0) – 0.30 mm

Наименование	Значение
	M4(1) – 0.30 mm
Needle Surface	shiny, without cracks, shells, scratches and burrs
Needle Surface Roughness, not exceeding	at the bar and point – Ra 0.32 μ m; at suture material attachment segment Rz 20 μ m
Applied Part Bluntness Radius, not exceeding	0.03 mm
Sharpness of Needle Lancinating Part, not exceeding	0.025 mm
Triangular Cutting Point Needle Width, not exceeding	along the length from 2 to 5 diameters of the shank – 0.025 mm
Curved Needle Curvature Radius	CTX – 48 mm CT-3 – 22 mm SH – 26 mm MH – 36 mm CT-1 – 36 mm MO-4 – 36 mm OS-8 – 40 mm FS – 26 mm CP-2 – 26 mm OS-6 – 36 mm FS-2 – 19 mm TAPERCUT V-4 – 17 mm

Наименование	Значение
	TAPERCUT V-26 – 26 mm PS-2 prime– 19 mm permissible variation ± 1 mm
Complete rated life should be at least	25 passes
Complete average life should be at least	50 passes
Sterilization Method	sterilized by ethylene oxide

The materials of which the product is made:

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APOSTILLE

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1. Country : United States of America
2. This public document has been signed by **CYNTHIA KELSEY**
3. acting in the capacity of NOTARY PUBLIC

bears the seal/stamp CYNTHIA KELSEY , NOTARY PUBLIC, BERKS COUNTY, COMMONWEALTH OF PENNSYLVANIA

Certified

at Harrisburg, Pennsylvania

6. The 1st day of June, 2016

by Pedro A. Cortés, Secretary of the Commonwealth of Pennsylvania

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


Cynthia Kelsey
Notary Public

Instruction for Use
Stratafix Spiral PDO (Polydioxanone) Knotless Tissue Closure Device
(Unidirectional)

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
The Stratafix Spiral PDO sutures are made from dyed (violet) polyester, poly(p-dioxanone). The pigment for the dye is D&C Violet No. 2. It is stated that polydioxanone has been found to be nonantigenic, nonprogenic and elicits only slight tissue reaction during absorption.

While the formation of bards in the Stratafix Spiral PDO reduces the tensile strength relative to non – barded suture material of the same size, tying of knots in non – barded suture materials also reduce their effective strength. For this reason, the strength of the Stratafix Spiral PDO can be compared to USP knot strength of non – barded sutures. Additionally, USP designations for diameter are used to describe the Stratafix Spiral PDO prior to barding, except for minor variation in suture diameter with a maximum overage of 0.1 mm.

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APPLICATIONS.

This medical product is applied only for therapeutic and diagnostic purposes.

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The loop of the STRATAFIX™ Spiral PDO Device for wound closure should be deployed by users familiar with surgical procedure, techniques and tissues. Physicians should consider the quantity and quality of tissue in which the loop will be anchored. The use of this suture may be inappropriate in highly vascularized tissue

or fragile tissue that cannot withstand potential further tightening/constriction within the loop.

Doctor should be familiar with surgical procedure and techniques involving absorbable suture before employing Stratafix Spiral PDO for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. Physicians should consider the *in vivo* performance (under ACTIONS section) when selecting a suture for use in patients. The use of this suture may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

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PRECAUTIONS

Stratafix Spiral PDO contains a loop end and barbs to anchor tissues and does not require knots to approximate opposing edges of a wound. Tying knots on the

barbed section of the material will damage the barbs and potentially reduce the suture tensile strength and barb effectiveness. For the device to function properly, both sides of the Stratafix Spiral PDO must first be anchored in the tissue by the deployment of the around robust tissue, and then with the other segment subsequently engaged in the tissue by the barbs. Additionally, when completing placement of the barbed segment an additional backstitch or bite of tissue lateral to the end to the incision is required to lock Stratafix Spiral PDO device in place.

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Reshaping needles may cause them to lose strength and be less resistant to banding and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in designated biohazard "sharps" containers.

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APPLICATION

Use the required per surgical procedure. Secure the non – barbed loop portion to a robust tissue. Start by taking a bite of the designated anchoring tissue and then pass the needle through the loop and pull in order to tighten loop around entrapped tissue. Once the Stratafix Spiral PDO is anchored into the tissue, take at least two loose bites of this tissue to be approximated. Once at least two bites, or passes through the tissue have been completed, the strand can then be grasped and the tissue can be approximated to the desired tension. Bites, or passes through the tissue, can be taken in a continuous suturing technique to approximate the wound. Care should be taken to utilize the Stratafix Spiral PDO on the barbed segments only. Do not attempt to approximate wound using the non – barbed segments, as the barbs are required for successful wound approximation with the Stratafix Spiral PDO.

To remove the entire Stratafix Spiral PDO: Cut the Stratafix Spiral PDO at the end anchored with the loop and pull the distal and to remove the device.

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HOW SUPPLIED

The Stratafix Spiral PDO is available sterile, dyed violet in various lengths and with a USP diameter size 0 through 3-0 (metric size 3.5 through 2). Each Stratafix Spiral PDO is single armed with various needles sizes.

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Phone: +7 (495) 588 77 77, Fax: +7 (495) 588 78 78.

SYMBOL USED FOR LABELING

	Device Item Number
	Lot Number
	Manufacturer
	Best before
	Sterilized by ethylene oxide
	The CE Label including the SGS (authorized body) identification number, the Product complies with the essential requirements of the European Directive 93/42/EEC on Medical Devices
	Qualified persons only
	Do not reuse product
	Do not resterilize product
	Temperature Range for Storage
	Temperature Range for Transportation
	Attention, address shipping documents

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Наименование	Значение
Nominal Length of Suture Material	20 cm – nom. length 44.7 cm; 30 cm – nom. length 54.7 cm; 45 cm – nom. length 69.7 cm; permissible variation $\pm 2\%$
Atraumatic Needle Suture Holding Strength, at least (thread for the specified value)	M1(5/0) – 1,7 H M1.5(4/0) – 2,3 H M2(3/0) – 4,5 H M3(2/0) – 9,0 H M3.5(0) – 11,0 H M4(1) – 15,0 H
Distance between Barbs	0.61 or 0.95 mm
Barb Length	0.34 or 0.58 mm
Barb Height	0.082(± 0.031) mm
Barb Angle	15°($\pm 5^\circ$)
Needle Curvature	Needle Type: CTX, CT-3, SH, MH, CT-1, MO-4 – 1/2 curved needle Needle Type: OS-8, FS, CP-2, OS-6, FS-2, TAPERCUT V-4, TAPERCUT V-26, PS-2 prime – 3/8 curved needle
Needle Point Sectional Shape	Needle Type: CTX, CT-3, SH, MH, CT-1, MO-4 – round cross-section Needle Type: OS-8, FS, CP-2, OS-6, FS-2, TAPERCUT V-4, TAPERCUT V-26, PS-2 prime – triangular cross-section
Needle Hardness	HV 4900-6475 N/mm ²
Needle Corrosion Resistance	corrosion resistant
Needle Attachment Segment Diameter	M1(5/0) – 0.35 mm M1.5(4/0) – 0.40 mm M2(3/0) – 0.50 mm

Наименование	Значение
	M3(2/0) – 0.50 mm M3.5(0) – 0.60 mm M4(1) – 0.70 mm
Needle Resilience Under test force of 5N, residual deformation, not exceeding	M1(5/0) – 0.60 mm M1.5(4/0) – 0.60 mm M2(3/0) – 0.40 mm M3(2/0) – 0.40 mm M3.5(0) – 0.30 mm M4(1) – 0.30 mm
Needle Surface	shiny, without cracks, shells, scratches and burrs
Needle Surface Roughness, not exceeding	at the bar and point – Ra 0.32 μm; at suture material attachment segment Rz 20 μm
Applied Part Bluntness Radius, not exceeding	0.03 mm
Sharpness of Needle Lancinating Part, not exceeding	0.025 mm
Triangular Cutting Point Needle Width, not exceeding	along the length from 2 to 5 diameters of the shank – 0.025 mm
Curved Needle Curvature Radius	CTX – 48 mm CT-3 – 22 mm SH – 26 mm MH – 36 mm CT-1 – 36 mm MO-4 – 36 mm OS-8 – 40 mm FS – 26 mm CP-2 – 26 mm OS-6 – 36 mm

Наименование	Значение
	FS-2 – 19 mm TAPERCUT V-4 – 17 mm TAPERCUT V-26 – 26 mm PS-2 prime– 19 mm permissible variation $\pm 1\%$
Complete rated life should be at least	25 passes
Complete average life should be at least	50 passes
Sterilization Method	sterilized by ethylene oxide
Curved Needle Curvature Radius	CTX – 48 mm CT-3 – 22 mm SH – 26 mm MH – 36 mm CT-1 – 36 mm MO-4 – 36 mm OS-8 – 40 mm FS – 26 mm CP-2 – 26 mm OS-6 – 36 mm FS-2 – 19 mm TAPERCUT V-4 – 17 mm TAPERCUT V-26 – 26 mm PS-2 prime– 19 mm permissible variation $\pm 1\text{mm}$

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