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Registered No. 20 *18* — 2235

NOTARIAL CERTIFICATE

Busan International Law & Notary Offices.

(Geoje-dong, Busan Judicial Offices Town Bldg)

#205, 28, Beobwon-ro, Yeonje-gu,

Busan, Korea

Instruction For Use

Attn: Federal Service of Health Care and Social Development Control of Russia

We, **OSSTEM IMPLANT Co., Ltd.**, hereby declare that the **instruction for use** provided to Federal Service of Health Care and Social Development Control of Russia for Imported Medical Product Recordation is complete, consistent and authentic.

Signature: _____

Name: Eom Tae Kwan

Title: President

OSSTEM IMPLANT Co., Ltd.

66-16, Bansong-ro 513beon-gil, Haeundae-gu, Busan, 48002, Korea

1. Designation of the medical device

1.1. Medical Device Name

Tools for surgical installation of dental implants (hereinafter - instruments)

1.2. Intended Use

Instruments are designed for surgical placement of Osstem implant attachments

1.3. Indications for Use

Use in accordance with general implant surgical guides. For more details please surgical manual, prosthetic manual and catalog.

1.4. Contraindications

Osstem implants and also tools for their installation should not be used in cases where the remaining areolar bone is diminished as such to provide adequate width or height to surround the implant. Lack of osseointegration or subsequent implant failure may occur in cases where there is insufficient bone available, and poor bone quality, also poor oral hygiene, heavy smoking, tobacco abuse, or medical conditions such as blood disorders or uncontrolled diabetes may also cause failure of the implant.

Contraindications include the following, but are not limited to:

- Patients with hemophilia or difficulties related to bone or wound treatment
- Patients with uncontrollable diabetes, heavy smoker or alcoholic
- Patients whose immunity system is inactive due to chemical therapy or radiation therapy
- Patients with oral infection or inflammation (improper oral hygiene, bruxism)
- Patients with untreatable occlusion/joint disorder, insufficient dental arch space
- Any patient who is not suitable for an operation

1.5. Possible Adverse Effects of Medical Device Use

A few problems may occur after the operation (loss of implant stability, damage of prosthesis, etc.). Deficient quality and quantity of the remaining bone, infection, allergic reaction, inferior oral hygiene or uncooperativeness of patient, implant mobility, partial deterioration of tissue, and improper position or arrangement of implants may cause the above mentioned problems.

1.6. Possibility of using Medical Device for different category of patients

The use of this product is possible according to the prescription of the doctor and is not possible if there are contra-indications. The product can be used for all ages, gender and ethnic groups.

1.7. Possible Effects of Medical Device on Ability to Drive and Use Machines

The possibility of the impact of using the product on the ability to drive vehicles and mechanisms is not revealed.

1.8. Operating Environment

Environmental requirements during operation

Temperature	From + 10 °C to + 40 °C
Humidity	From 0% to 100%
Pressure	From 700 hPa to 1060 hPa

1.9. Conditions of Use

The customer who uses these things is dentist who had been educated officially.

2. Data for Medical Device Development and Manufacturing

2.1. Name

Osstem Implant Co., Ltd.

2.2. Address of the developer and manufacturer

66-16, Bansong-ro 513beon-gil, Haeundae-gu, Busan, Korea 48002

2.3. Address of the manufacture sites

Osstem Implant Co., Ltd., Haeundae plant, 66-16, Bansong-ro 513beon-gil, Haeundae-gu, Busan, Korea 48002

2.4. The authorized representative of the manufacturer

«OSSTEM» LLC

115432, Moscow, Prospekt Andropova, house 18, building 7, floor 8, office 1

3. Medical Device Classification

Class depending on the potential risk of a medical device, in accordance with the nomenclature classification	I
---	---

Type of contact with the body	Surgically invasive transient use devices (< 60 minutes of patient contact).
Sterilization	Non-sterile

4. Principles of Operation of the Medical Device

Tools for surgical installation of dental implants are for Surgical placement of Osstem implant fixtures. Multi Abutment Outer Driver is to fix and separate when using Multi Abutment. Locator Torque Driver is to use for the tightening torque to Locator Abutment. Angled Torque Driver is to use for the tightening torque to Angled Abutment. Torque Wrench is used for fixture grafting or screw tightening. Its neck is bended automatically when apply too much force beyond allowable. Multi Tool is used for combination or disconnecting Male Cap and Male. Tools will be used as aids in the Implant procedure.

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5.2. Possible Combinations with Other Equipment

Medical devices are used with each abutment:

- Multi Abutment Outer Driver, Multi Abutment Machine Driver _ Multi Abutment
- Locator Torque Driver _ Locator Abutment, Stud Abutment
- Angled Torque Driver _ Angled Driver
- The other products are used with screws and abutment of OSSTEM.

5.3. Materials of the medical device

	Code	Type and grade of steel of Steel
Locator Torque Driver	TWLDSK	Stainless Steel X15TN (C - 0.40%, Mn - 2.00%, Mo – 2,0%, Cr - 15.50%, N - 0,2%, V – 0.30%, Fe – 79,6%)
	TWLDLK	
Angled Torque Driver	ATD12S	
	ATD12L	
Multi Abutment Outer Driver	MAOD	Stainless Steel CarTech® TrimRite® (C 0.15-0.30%, Mn 1.0%, P 0.04%, S 0.03%, Si 1.0%, Cr 13.50-15.00%, Ni 0.25-1.00, Mo 0.40-1.00%, Fe 80.63-83.63%)
Torque Wrench	TW30	
Torque Wrench Bar Type	TW30B	
Multi Tool	PTCMT	Stainless Steel SUS316L (Fe 61.4-72%, C<0.03%, Cr 16-18.5%, Ni 10-14%, Mo 2-3%, Mn <2%, Si<1%, P <0.045%, S <0.03%) TiN coating


6. Basic Parameters and Performance Data of Medical Device

Contact your authorized manufacturer's representative or website www.osstem.ru for more information*

*- see Annex 1 for information

7. Labeling

7.1. Draft marking in Russian

7.2. Explanation of symbols

Symbol	Definition
	Caution. Indicates there is a possible danger to the user. In the manual, the symbol indicates information the user must read.
	Reorder number
	Lot number
	Manufacturer

Symbol	Definition
	Date of Manufacture
	Expire date
	Mark of conformity of voluntary certification

8. Packaging

Devices have two options for delivery : Carton Box and Pack Case

One or three unit of the Angled Torque Driver (ATD12S, ATD12L, ATD12S3S, ATD12L3S), One or unit of the Locator Torque Driver (TWLDSK, TWLDLK) or Multi Abutment Outer Driver (MAOD) packaged in Inner blister. Inner blister is placed in the Carton Box with the instruction for use.

One unit of the Torque Wrench(TW30), Torque Wrench Bar Type(TW30B) and Multi Tool(PTCMT) packaged in Pack Case with the instruction for use.

9. Transportation and Storage Conditions

	Environmental conditions for storage	Environmental conditions for shipping
Temperature	From + 1°C to + 30 °C	From -18 °C to + 60 °C
Humidity	From 20% to 70%	From 10% to 80%
Pressure	From 700 hPa to 1060 hPa	From 500 hPa to 1060 hPa

10. Service Life

Drivers are 50 times of use, and the others tools don't have expire date.

11. Safety Requirements, Precautions and Warnings

Warnings *(The warnings are directly linked to the implants, which are placed with tools)*

For the safe and effective use of Osstem implants, it is strongly suggested that specialized training be undertaken since the surgical techniques required to place dental implants are highly specialized and complex procedures. Improper patient selection and technique can contribute to implant failure and/or loss of supporting bone. Osstem implants are intended for use only in the indicated applications. Dental implants must not be altered in any way. The use of electro-surgical instruments or lasers around metallic implants and their abutments is not recommended due to the risk of electric shock and/or burns. Implant mobility, bone loss, or chronic infection may indicate implant failure. If the implant becomes contaminated by the patient's bodily fluids in any way, the implant cannot be used in any other patient.

Cautions

- The driver can also be used up to 50 times. The recommended number of uses for the Angled Torque Driver is 10, and the recommended torque is 30Ncm or below.
- Separate all used tools immediately after surgery and after cleaning and drying, then store at room temperature.
- Do not leave the instruments in a place vulnerable to contamination.
- Be sure to sterilize surgical tools in the autoclave (132°C for 15 minutes) before and after each new surgical procedure.
- Do not use hydrogen peroxide as a disinfectant or cleaner (it may cause damage or discoloration of the TiN coating, laser marking and color coding).
- This product comes with a 1 year warranty from the date of purchase.

12. Environment Protection

12.1. Environment protection requirements

Medical device during use, transportation and storage has no negative effects on humans and the environment.

12.2. Medical Device Recycling or Disposal

After use the device may be potential biohazard. During disposing of medical device, it is required to provide the mandatory implementation of all the necessary procedures in according to local laws and hospital practices.

13. Methods and Detergents for Disinfection and Pre-cleaning

These products don't be guided certain detergents and method.

Before using tools are sterilized in autoclave of 132°C for 15 minutes and after using, clean products with water.

14. Medical and Animal Substance

Medical Device does not contain any substance which may be considered a medicinal or animal origin substance.

15. Sterility

Tools for surgical installation of dental implants are supplied non-sterile.

Autoclave sterilization process at a temperature of 132°C for 15 minutes before and after each new surgical procedure.

16. The list of documents/standards applicable to the medical device

Standard Number	Standard Name
EN 1041:2008	Information supplied by the manufacturer of medical devices
EN 1641:2009	Dentistry - Medical devices for dentistry - Materials
EN 1642:2011	Dentistry - Medical devices for dentistry - Dental implants
EN 980:2008	Symbols for use in the labelling of medical devices
EN ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity
EN ISO 13485:2012	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 14971:2012	Medical devices – Application of risk management to medical devices
ASTM F899-12b	Standard Specification for Wrought Stainless Steels for Surgical Instruments

17. Installation Requirements

Not applicable.

18. Maintenance Requirements

Not applicable.

19. Warranty

This product comes with a 1 year warranty from the date of purchase.

20. Complaints

«OSSTEM» LLC

115432, Moscow, Prospekt Andropova, house 18, building 7, floor 8, office 1


www.osstem.ru

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E-mail: info@osstem.ru

Annex 1. Basic Parameters and Performance Data of Medical Device


(The client receives this information on request when contacting an authorized representative of the manufacturer or on the website www.osstem.ru)

Multi Abutment Outer Driver is to fix and separate when using Multi Abutment.

The tolerance is $\pm 5\%$, unless otherwise noted

Name of characteristic	Value
A(± 0.1), mm	18.5
B(± 0.1), mm	8.0
Weight, g	2.57
Outer diameter, mm	4.5(+0.05,-0)
Inner diameter, mm	3.31(+0.01,-0)


Locator Torque Driver is to use for the tightening torque to Locator Abutment.

The tolerance is $\pm 5\%$, unless otherwise noted

TWLDSK	
Name of characteristic	Value
A(± 0.1), mm	16.5
B(± 0.05), mm	3.175
C(± 0.05), mm	1.78
Rotation torque, N*m	30
Weight, g	1.26

TWLDLK	
Name of characteristic	Value
A(± 0.1), mm	21.5
B(± 0.05), mm	3.175
C(± 0.05), mm	1.78
Rotation torque, N*m	30
Weight, g	1.56

Angled Torque Driver is to use for the tightening torque to Angled Abutment.

The tolerance is $\pm 5\%$, unless otherwise noted

ATD12S

Name of characteristic	Value
A(± 0.1), mm	13
B(+0, -0.02), mm	1.23
Length of working part, mm	1.64(± 0.05)
Diameter C, mm	5(+0, -0.1)
Rotation torque, N*m	30
Weight, g	0.89

ATD12L

Name of characteristic	Value
A(± 0.1), mm	20
B(+0, -0.02), mm	1.23
Length of working part, mm	1.64(± 0.05)
Diameter C, mm	5(+0, -0.1)
Rotation torque, N*m	30
Weight, g	0.96

Torque Wrench (TW30)


Torque Wrench is used for fixture grafting or screw tightening. Its neck is bended automatically when apply too much force beyond allowable.

The tolerance is $\pm 5\%$, unless otherwise noted

Name of characteristic	Value
D1(± 0.1), mm	11.2
D2(± 0.1), mm	9.0
Inner diameter of the head, mm	6.52
Rotation torque, N*m	10/20/30
Weight, g	26.50

Torque Wrench Bar Type (TW30B)


Torque Wrench is used for fixture grafting or screw tightening. Its neck is bended automatically when apply too much force beyond allowable.

The tolerance is $\pm 5\%$, unless otherwise noted

Name of characteristic	Value
A(± 0.2), mm	90.4
B(± 0.1), mm	11.2
C (± 0.2), mm	68.2
Inner diameter of the head, mm	3.26
Rotation torque, N*m	10/20/30
Weight, g	28.05

Multi Tool (PTCMT)


Multi Tool is used for combination or disconnecting Male Cap and Male.

The tolerance is $\pm 5\%$, unless otherwise noted

Name of characteristic	Value
A(± 0.1), mm	7.9
B(± 0.2), mm	82.08
Length of working part, mm	0.7
Diameter of working part, mm	2.15(+0,-0.05)
Weight, g	15.69

	Code	Hardness (Minimum)	Tensile Strength (Minimum)	Surface roughness
Multi Abutment Outer Driver	MAOD	255HB	621MPa	No more than 0.99 μm
Locator Torque Driver	TWLDSK	373HB	1344MPa	No more than 0.99 μm
	TWLDLK			No more than 0.99 μm
Angled Torque Driver	ATD12S	373HB	1344MPa	No more than 0.99 μm
	ATD12L			
Torque Wrench	TW30	255HB	621MPa	No more than 0.99 μm
Torque Wrench Bar Type	TW30B	255HB	621MPa	No more than 0.99 μm
Multi Tool	PTCMT	100HB	737.04MPa	No more than 0.99 μm

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Lee, Kyung Min attorney-in-fact of
OSSTEM IMPLANT Co., Ltd.
President Eom, Tae Kwan
appeared before me and admitted said
principal's subscription to the attached
Instruction For Use.

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This is hereby attested on this 29th
day of October, 2018 at this office.

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Busan International Law & Notary Offices
Belong to Busan District Public Prosecutors' Office
(Geoje-dong, Busan Judicial Offices Town Bldg)
#205, 28, Beobwon-ro, Yeonje-gu, Busan, Korea

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
Park Yongsoo

(Signature of the Notary Public)

Park Yong Soo

This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public since
Dec. 31, 1982 under Law No. 3594

Numbered, sewed and sealed 11 sheets in all

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