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Warning!

For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Caution: Federal law restricts this device to sale by or on the order of a physician.

1. Indications for Use

The Hemostasis Clip System is intended for the following use:

- endoscopic marking.
- hemostasis for
 - mucosal / sub-mucosal defects < 3 cm.
 - bleeding ulcers.
 - polyps < 1.5 cm diameter.
 - diverticula in the colon.
- closure of GI tract luminal perforations < 20 mm that can be treated conservatively.

2. Contraindications

The following contraindications can be found in the medical literature with regards to the Hemostasis Clip System:

- Mucosal/submucosal defects greater than 3 cm.
- Polyps greater than 1.5 cm diameter.
- The patient with poor general condition cannot tolerate endoscopy.
- The patient has narrow upper digestive tract where endoscope cannot pass through.
- The patient has serious coagulation disorders and hemorrhagic diseases.
- The patient is allergic to the device and the drugs used in the operation.
- The patient who is not suitable to use the product according to the diagnosis.
- The patient or the families are uncooperative.

With regard to individual endoscopic indications for instrument application, the qualified user is ultimately responsible to assess contraindications and immanent risks of instrument application vs. its benefits for each specific medical situation.

3. Potential Complications

- Inflammation of tissue, perforation, bleeding or mucosal damage for the patient;
- Infection, septicemia, etc;
- Complications which are not currently known or observed may be present.

4. Warning!

Rx Only

In accordance with the laws and regulations, the product is intended for use by trained professional doctor only.

The user of this equipment should be thoroughly trained in the applicable procedure. Furthermore, failure to read and thoroughly understand the contents of this instruction manual may result in serious injury to the patient and/or user. It is essential to follow the instructions enclosed, as well as instructions found in other instruction manuals, infection control requirements and current hospital protocols for usage, cleaning and sterilization. Possible injuries related to endoscopic procedures may include perforation, electrical burns and shock, hemorrhage, infection, explosion, etc.

! Failure to follow these instructions may result in malfunction and/or damage of the instrument or endoscope!

! Do not use this device for any purpose other than the stated indications.

- ! Operation of this instrument is based on the assumption that open surgery is possible as an emergency measure if the clip cannot be detached from the instrument or if any other unexpected circumstance takes place.
- ! It might be impossible to stop bleeding depending on the hemorrhage situation because the clip performance for hemostasis is limited. Prepare more than one hemostasis device. Select the most appropriately hemostatic device to different hemorrhage situations appropriately. Some devices may be used together for best result.
- ! Re-bleeding may occur on the clipping site, depending on the local condition. Check the patient for any re-bleeding after the operation as appropriate.
- ! As the clips are made of stainless-steel, do not use them on a patient who is severely allergic to metals. Otherwise, allergic symptoms may occur. This device does not contain latex.

5. Caution!

- Patient should be informed and expressed his/her acceptance for the details of the operation and all the potential risks and complications, which may lead to injury, illness or death of the patients.
- Do not use this Instrument when hemostasis cannot be verified visually within the endoscopic field of view after application.

6. Storage

Store the device in a cool, dark and dry place, and do not expose the package to organic solvent, ionizing radiation or ultraviolet radiation. The sterilization is valid for 3 years.

7. General Safety Instruction

- Wearing protective garments (i.e., gloves, protective eyewear, gowns, etc.) will assist in the protection against risk of cross-contamination and injury.
- Do not use an instrument if it does not function properly. If the device has any physical defect(s), replace the instrument with a new and/or functioning instrument.
- Check the endoscope's working channel diameter prior to introducing the instrument into the endoscope. To ensure the appropriate fit of the instrument the diameter of the endoscope's working channel needs to be at least "0.2 mm" larger than the instrument's outer diameter.
- ! Never use products which do not comply with the technical specifications!
- Do not force an instrument down the working channel of the endoscope. This could result in endoscope and/or instrument damage.
- Do not under any circumstances make any modifications to the instrument or components described in this manual.

8. Inspection and Testing

- Note any damage/open seal of the packaging.
- ! If damage or an open seal is detected, sterilization may be compromised and the instrument should be discarded and replaced.
- Check date of sterilization.
- ! Never use products with non-valid sterilization date.
- Check the instrument prior to each use.
- ! Should any irregularities occur, dispose the instrument and replace it by a new one.
- ! Using a damaged instrument with physical or functional abnormalities may result in endoscope damage, patient and/or user injury.
- ! The product is intended for single use only and one patient only. Do not reuse the device, nor share it by two or more patients.

9. Notes

- Always observe the endoscopic image during operation. Make sure that the device is present on the image and the operation is normal. If the clip falls off from the endoscope before release, take it out with foreign body forceps.
- Please make sure that the visual field of endoscope is clear before use.
- After operation, follow the doctor's advice not to eat food too early. Improper diet, forcible

defecation, cough and other activities increasing abdominal pressure will result in forcible falling off of the clip and bleeding afterwards.

- Limited studies indicate that lesions located in the esophagus and the lesser curvature of the stomach may be difficult to treat with linear-viewing endoscope.
- Limited studies indicate that the treatment of esophageal varices may require clipping in combination with a sclerosing agent.
- Limited studies indicate that clipping hard or severely fibrotic lesions to achieve hemostasis may be more difficult.
- Limited studies have shown that the number of clips required for hemostasis may vary depending upon the anatomical site, histology, lesion type and patient condition and history. A sufficient quantity of clips should be prepared in consideration of all of these factors prior to the procedure.
- Limited studies indicate that the hemostasis clips remain in place for an average of 9.4 days; re-bleeding may occur if the clips detach within 24 hours.
- Limited studies indicate that the use of clips in the presence of bacterial contamination may potentiate or prolong infection.

10. Application

- Remove yellow protective sleeve from the distal end.
- Carefully insert device into working channel ensuring that the clip (1) is closed.
! The endoscope must remain as straight as possible when inserting the instrument.
- Pass the clip (1) through the endoscope channel towards targeted site.
- Open the clip (1) by pushing the fingerslide (2) forward.
- The clip (1) can be rotated clockwise or anti-clockwise by turning the handle (3) until correct position is reached.
- Advance instrument until contact is made with the targeted site.
- When satisfied with the clip (1) positioning, close the clip (1) against the tissue by pulling the fingerslide (2) back until tactile resistance is felt.
- The clip (1) position may now be assessed prior to deployment.
- If the clip (1) is not in desired position, the clip (1) may be re-opened and repositioned.
! Do not continue pulling back the fingerslide (2) further beyond the tactile resistance until you are ready to deploy the clip (1), otherwise you may not be able to re-open the clip (1). If you hear or feel a click, the clip (1) cannot be reopened.
- To deploy the clip (1), continue pulling back the fingerslide (2) beyond the tactile resistance point. You will hear or feel a click as the clip (1) detaches.
- Do not again push forward the fingerslide (2)! The clip (1) is already released.
- Carefully retract the instrument from the endoscope.
- Endoscope must remain as straight as possible when withdrawing the device.
- After use, dispose the entire instrument according to legal requirements and your facility's infection control protocols.

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The Spirit of Care.

Medi-Globe GmbH, Medi-Globe-Str. 1-5, 83101 Achenmuhle - Germany

Medi-Globe Offices and Sales Representations

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URNr.
UT

№ 0152 /2019

I hereby certify that the signature on this document was acknowledged to be his signature by

Mr. Ralph Rothe,
born on 24th of June 1967,
of Rupertistraße 54,
83278 Traunstein, Germany
personally known to me

I am not capable of the Russian language and didn't examine that part of the document for that reason.

Prien a. Chiemsee, this 29th of January 2019

v. Daumiller, notary a.D.,
as officially announced representative
of Bertrand Koller, notary public in
Prien a. Chiemsee, Germany

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