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
**Instructions for use
of medical device**

**Gel for Vaginal Microflora Normalization
Multi-Gyn® FloraFem, 5 ml, 5 pieces per pack**

**Manufactured by BioClin B.V., Delftechpark 55, 2628 XJ, Delft,
Netherlands)**

I confirm the accuracy, correctness and authenticity of this text document	
Surname, name and position of the person authorized to sign such documents <u>Mr. M.J. BAST</u> <u>Operational Director</u>	Signature  Date <u>21-2-2019</u>

BIOCLIN
BIO-ACTIVE SELF CARE

INSTRUCTION FOR USE of a medical device

Name

Multi-Gyn® FloraFem, gel for normalization of vaginal flora, 5 ml, 5 pieces per pack (hereinafter referred to as "Multi-Gyn® FloraFem")

Purpose

Multi-Gyn® FloraFem is designed for prevention of vaginal microflora and pH level disorders, with the aim of its normalization and natural maintenance due to moisture.

Description

Non-sterile, semitranslucent, colorless with yellowish tone, low viscosity gel with a slight smell of aloe.

Composition:

Active component: Galactoarabinan Polyglucuronic Acid Crosspolymer*

*2QR – complex, patented bio-active polysaccharides

Excipient components: lactic acid, xanthan gum, fructooligosaccharides, lactitol (lactitol monohydrate), hydroxypropyl methylcellulose, sodium hydroxide, glycerin, levulinic acid, p-anisic acid, caprylyl glycol, water.

Characteristics and properties

Multi-Gyn® FloraFem is a gel for intravaginal application, the effectiveness of which is due to the combination of following properties: barrier (anti-adhesive) activity of 2QR, blocking the mechanism of adhesion of harmful microorganisms to vaginal epithelium by correcting or maintaining natural acidic environment in vagina. These properties provide an effective mechanism aimed at eliminating symptoms: such as itching, irritation, unpleasant odor and excessive discharge. The mechanism by which Multi-Gyn® FloraFem eliminates these vaginal problems, is based on neutralization of harmful microorganisms.

Fructooligosaccharides and lactitol (disaccharide analogue of lactulose) are used in the product as nutrients that stimulate lactobacilli growth. Xanthan gum and hydroxyl propyl methyl cellulose are natural bioadhesive thickeners that contribute to film formation. Combination of fragrances of natural origin (water, glycerin, levulinic acid, p-anisic acid) with a mild inherent odor and biostabilizing properties at the same time. Caprylyl glycol is used as a humectant agent and has antimicrobial properties. Sodium hydroxide and lactic acid are used as a pH buffer system that corrects and maintains stable physiological acid pH value (3.7-4.1) after gel application.

Multi-Gyn® FloraFem is a topical gel and is not absorbed into body. All components of the product are included in the register of cosmetic components of the European Commission. As ingredients of the medical device, these substances are used to ensure its use for moisture and do not have any pharmacological effect.

Indications for use

The gel is used by women of all ages for prevention of vaginal microflora and pH disorders:

- Itching, irritation, unpleasant odor and discharge (including caused by vaginal candidiasis),
- disruption/imbalance of the vaginal microflora disorders associated with menstruation, sexual activity, swimming, taking antibiotics, health problems;
- (recurrent) distribution of vaginal yeast microorganisms (of the genus Candida).

Contraindications

Hypersensitivity to the components of the device. Before the first use, apply a small amount of the gel on the inside of the forearm. It is not advisable to use Multi-Gyn® FloraFem in case of itching or reddening.

Dosage and administration

Intravaginal. The gel is applied directly to vaginal mucous membrane using a special form applicator, which is a continuation of the tube

Recommended doses:

Frequency and duration of application is set individually depending on experience or according to recommendations of a gynecologist.


To prevent vaginal diseases (e.g. recurrent candidiasis), maintain vaginal lactobacillus flora and protect against vaginal diseases use 1 tube of Multi-Gyn® FloraFem with 1 dose of gel every 3 days.

To relieve itching, burning and irritation in vagina, eliminate unpleasant odor and discharge use 1 tube of Multi-Gyn® FloraFem with 1 dose of gel whenever you need to eliminate these symptoms.

To eliminate symptoms of vaginal candidiasis use 1 tube of Multi-Gyn® FloraFem with 1 dose of gel before bedtime, for 5 days.

Multi-Gyn® FloraFem should not be used for more than 30 days, it allows vaginal flora to self-regulate without additional support. In case of persistent or recurrent symptoms consult your physician.

Method of administration into the vagina: See the drawings below.

1. Shake the tube contents against the tip of the tube

2. Break off the tip of the tube by twisting.

3. Carefully insert the thin end of the applicator (of the tube) into the vagina

4. Squeeze slightly the tube to deposit the gel into the vagina.

5. Remove the tube, still squeezing

Note: The tube is designed to dispense the recommended amount of gel.

Preliminary disinfection or sterilization of the tube before gel application is not required.

Application of the product does not require special skills and training of consumers and/or medical workers.

Side effect

No side effects known.

Warnings and precautions

For intravaginal application only!

Please be aware! Immediately after gel application you may feel warmth and easy tingling. It is normal and completely disappears in a few minutes.

Since Multi-Gyn® FloraFem has optimal acidity (low pH) for vagina, and sperm has alkaline reaction (high pH), its use can reduce sperm activity. Therefore, it is not recommended to use Multi-Gyn® FloraFem during ovulation if you are trying to get pregnant.

Multi-Gyn® FloraFem is not a contraceptive, does not treat or prevent sexually transmitted diseases.

Interaction with other drugs

Multi-Gyn® FloraFem is not a drug and has no known interactions with drugs.

Pregnancy and lactation

No contraindications for use of Multi-Gyn® FloraFem during pregnancy and lactation are established.

Presence of a medicinal product for medical use or materials of animal and (or) human origin in the medical device

Multi-Gin® FloraFem does not contain medicines for medical use, materials of animal and (or) human origin.

All components of the product are included in the registry of cosmetic components of the European Commission. In the composition of the medical device, these components do not possess a pharmacological action, but perform a mechanical function (barrier, anti-adhesive) and moisture.

Multi-Gyn® FloraFem contains a Lactitol monohydrate, which according to data of the State Register of Drugs of the Ministry of Health of the Russian Federation is registered as a pharmaceutical substance (28.09.2011 ΦС-000158 Manufacturer by Danisco USA Inc.). In this composition, this component has nutritional functions for the growth of lactobacilli (lactitol is a disaccharide analogue of lactulose).

Technical specifications

A tube with an applicator for intravaginal application, which is a continuation of the tube, is made of white plastic (white polyethylene resin Dowlex 2552E, certified for food and medical devices production) and sealed.

Weight of 1 tube is 7.9 - 8.3 g, total length of the tube with an applicator is 111.0 ± 0.7 mm. Weight of the product (5 single-dose tubes with instructions for use in a cardboard pack) is 52.3 - 52.7 g.

Multi-Gyn® FloraFem does not comprise parts (units), accessories that are not medical, but intended for use in combination with it.

Multi-Gyn® FloraFem is produced in accordance with the requirements of national and international standards: Council Directive 93/42/EEC, EN ISO 13485, and EN ISO 14971.

In terms of biological effects on the patient, the materials used meet the requirements of the EN ISO 10993 series of standards. The gel is biologically safe and resistant to biological fluids and body tissue secretions with which it comes into contact during use.

Labeling and packaging

Labeling of the device contains the following information (symbols on the primary and secondary packaging):

Symbol	Explanation (symbol title)	Symbol	Explanation (symbol title)
	Manufacturer		Read the instruction carefully before use
	Temperature limits		Do not reuse
 0482	Conformity to European standards and identification number of the competent certification authority.		

Additional advertising information is allowed on the product packaging.

Presentation

Vaginal gel, 5 ml. 5 single-dose, plastic tubes and the Instruction for Use in a cardboard pack.

Shelf life

3 years. Do not use after the expiration date.

Medical device maintenance and repairs

The device is for single use only. The device is not subject to maintenance and repairs.

Manufacturer's warranty

The manufacturer guarantees the preservation of the product's properties provided that device transportation parameters and storage and operation conditions are maintained over the shelf life of the product.

Storage conditions

Keep out of reach of children. Store in a place protected from light at a temperature of +10 to +25°C and relative humidity of 60%-80%. Do not freeze. The device is for single use only. Do not store an opened device.

Operation conditions

The device is intended for individual use at home and/or in healthcare institutions. To prevent contamination, the tip of a tube must not come into contact with any other surface. If the primary package of the device is visibly damaged, it not advisable to use the product. The recommended operation conditions of the device are the nominal temperature of +32 to +42°C and relative humidity of 40%-80%.

Transportation conditions

All types of enclosed transport in accordance with the rules for carriage of goods. The medical device can be transported at a temperature of +10 to +25°C and relative humidity of 60%-80%. Do not freeze.

Disposal and destruction

In the home setting, any unused (expired) device or used device should be discarded with household waste.

In medical and/or pharmaceutical institutions, any unused (expired) device or used device is discarded in accordance with the requirements/rules of the purchasing country. At the territory of the Russian Federation, the medical device belongs to the hazard class A (epidemiologically safe waste) in accordance with SanPiN 2.1.7.2790-10 and is subject to disposal with household waste.

Pharmacy purchasing terms: Over the counter.


Manufacturer:

BioClin B.V., Delftechpark 55, 2628 XJ, Delft, Netherlands,
Tel. +31 (15) 2561999, e-mail: info@bio-active.com

Authorized representative of the manufacturer in the Russian Federation:

AO "Pharm-Sintez", 2nd Kabelnaya Street, 2, bld. 46, Moscow, 111024.
Tel. +7(495)796-94-33, e-mail: info@pharm-sintez.ru

Date of revision of the Instruction: November 2018.

I, mr Klaas Johannes van den Dool, civil-law notary practising in Delft, have authenticated the above signature of **M.J. Bast**, whose personal details are stated on the identity document of which a copy has been attached hereto. The person who acts in this matter has not been informed by the undersigned about the content and the consequences.

The statement of the undersigned only serves to establish the identity of the person who acts in this matter.

The authentication of the signature of the person who acts in this matter does not imply any assessment of his authority and/or of the capacity in which he acts.

Any liability that does not pertain to the identity of the person appearing is explicitly excluded.

Signed in Delft on January 8th, 2019.

APOSTILLE

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1. Country: THE NETHERLANDS
This public document
2. has been signed by mr. K.J. van den Dool
3. acting in the capacity of notary at Delft
4. bears the seal/stamp of aforesaid notary

Certified

5. in Den Haag
6. on 08-01-2019
7. by the registrar of the district court of Den Haag
8. no. 19 207

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Page reserved for the authorities responsible for issuing the passport.

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Apeldoorn

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M/M

3 delum van elgite / dele of issen / dele de deliverance
DECEMBER 2016

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