APPROVED BY

Silvana BROWN
Sr Director, International Regulatory Affairs
BALT EXTRUSION SAS

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M. BROWN

(non exclusively to certify the above signature)

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INSTRUCTIONS FOR USE

CRISTAL

Sterile Disposable Balloon Catheters for Peripheral Angioplasty and Valvuloplasty

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These Instructions for Use have been prepared to provide information for registering a medical device in the Russian Federation. The information is provided as required in accordance with the Rules for State Registration of Medical Devices (Decree of the Government of the Russian Federation No. 1416 dated December 27, 2012). The contents of the Instructions for Use comply with the requirements of Order of the Ministry of Healthcare of the Russian Federation No. 11н On Approval of Requirements to the Contents of Technical and Operational Documentation of the Manufacturer (Producer) of the Medical Device dated January 19, 2017.

The device is developed, manufactured and sold in accordance with the requirements of EN ISO 13485:2012 + AC : 2012 (Certificate No 513975 MP2012, date of issue February 6, 2018, valid until March 31, 2019, notified body for medical devices: DQS Medizinprodukte GmbH) and is labelled with a CE mark (EC certificate number: 513975 MR2 annex II excluding section 4, date of issue May 12, 2018, date of expiry July 31, 2022, EC certificate number: 537614 MRA annex II.4, date of issue February 1, 2018, date of expiry January 31, 2023 notified body: DQS Medizinprodukte GmbH).

Device name:

CRISTAL Sterile Disposable Balloon Catheters for Peripheral Angioplasty and Valvuloplasty, design options:

- CBVP8x20 / 70;
 - CBVP10x20 / 70;
 - CBVP12x20 / 70;
 - CBVP15x25 / 70;
 - CBP8x3035 / 110;
 - CBP10x30 / 110;
 - CBP12x40 / 110;
 - CBV12x30 / 110;
 - CBV15x30 / 110;
 - CBV15x40 / 110;
 - CBV16x40 / 110;
 - CBV18x30 / 110;
 - CBV18x40 / 110;
 - CBV18x45 / 110;
 - CBV20x30 / 110;
 - CBV20x40 / 110;
 - CBV20x45 / 110;
 - CBV23x45 / 110;
 - CBV25x40 / 110;
 - CBV25x50 / 110;
 - CBV28x50 / 110;
 - CBV30x60 / 110;
 - CBV35x60 / 110;
 - CBV40x40 / 110.
- Instructions for use (supplied with each device) – 1 pc.

Device classification in accordance with the Nomenclature Classifier of Medical Devices approved by Order of the Ministry of Healthcare of the Russian Federation No. 4н dated June 6, 2012 – 136170;

OKPD2 – 32.50.13.110 – Syringes, needles, catheters, cannulas and similar instruments;

Class depending on the potential risk degree of medical application in accordance with the Nomenclature Classifier of Medical Devices approved by Order of the Ministry of Healthcare of the Russian Federation No. 4н dated June 6, 2012 – 3;

Type of contact with the human body – Short-term contact (less than 24 hours) with the internal environment and body tissues, blood.

Conditions of use: Medical treatment facilities and healthcare centers, only when prescribed by a medical specialist.

2. Full name of the medical device manufacturer

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3. Intended use of the device and principles of operation:

3.1. Intended use of the device: The device is intended for peripheral percutaneous transluminal angioplasty (PTA) and valvuloplasty.

4. Device functional specifications:

CRISTAL sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, are coaxial double-lumen catheters with a balloon at the distal end. The outer lumen is intended for inflating and deflating the dilatation balloon, and the inner lumen – for the guidewire allowing guiding the catheter and facilitating its passing through stenosis to be dilated. The balloon is designed so that its diameter and length correspond to the set pressure. One or two radiopaque markers allow for the right position of the catheter's balloon under angiographic control. Standard Luer lock hubs ensure connection with the material.

The device is available for delivery in several modifications. They vary only in balloon diameters, balloon lengths and total catheter lengths. Specific recommendations are intended for angioplasty (Peripheral CRISTAL BALLOON CBP) and valvuloplasty (valvuloplasty CRISTAL BALLOON CBV and petit valvuloplasty CRISTAL BALLOON CBVP).

Each catheter is placed on a stand. The balloon is protected by a special sheath. Then the device is packaged in a removable pouch and a carton containing 2 catheters and the instructions for use.

The CRISTAL BALLOON catheters are disposable devices and are sterilized with ethylene oxide. Pursuant to Rule 7 of the European Directive 93/42/EEC these are Class III devices. In fact, these products are surgically invasive devices intended for short-term use (less than 24 hours) in contact with the central circulatory system (circulating blood).

These devices shall meet the requirements with regard to biocompatibility of materials, compatibility with the introducer sheaths being used, navigation ability and visibility when performing interventional radiology procedures, mechanical stability and balloon inflation. Preclinical evaluation results are given in Chapter 17 of the Technical File. A gap analysis was performed based on the biocompatibility data and showed that additional biocompatibility testing was required to ensure compliance with the essential requirements of Directive 93/42/EEC Rule 7.

5. Risks of using the medical device, contraindications, expected and predictable side effects:

5.1 Indications:

CRISTAL sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty are intended for peripheral percutaneous transluminal angioplasty (PTA) and valvuloplasty.

5.1 Contraindications:

Application of Cristal Balloon catheters is contraindicated in the following cases:

- if a guidewire cannot pass the lesion;
- coronary and intracerebral arteries.

5.3. Potential complications and side effects

Complications which can result from balloon-assisted dilatation include, but are not limited to:

- allergic reaction (to the device, contrast medium, or medicinal products);
- hypotension;
- arrhythmia;
- arteriovenous fistula;
- embolism (air, the device, an atherosclerotic plaque, etc.);
- haematoma;
- bleeding, including at the puncture site;
- sepsis/infection;
- pseudoaneurysm;
- infection;
- thromboembolic episodes;
- vascular wall damage (dissection, perforation, rupture);
- vessel occlusion;
- vasospasm;
- recurrent expanded vessel stenosis;
- death.

5.4 Precautions

- Do not use if the package is opened or damaged. The devices are sterile if the package is undamaged.
- Cristal Balloon catheters are intended for single use. Do not reuse. Reuse of the device may result in a high risk of contamination to the patient and performance degradation. Do not resterilize.
- Store in a dry and dark place at room temperature. Do not use past the expiration date.
- These devices should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology and/or cardiology.
- Never move an intravascular device if resistance is encountered until the causes of resistance have been determined by means of angiographic control. Movement against resistance may result in damage to the catheter or injury to the patient.
- Use a suitable pressure gage to monitor pressure during the procedure.
- Observe the recommended inflation pressure specified on the device label.
- Never use air or gas to inflate the balloon.
- Follow the instructions for use provided with the agents used in the catheter.

5.5. Information on potential consumers: The device should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology. Conditions of use: Medical treatment facilities and healthcare centers, only when prescribed by a medical specialist.

6. Technical specifications of the medical device:


Figure 1: Appearance of CRISTAL BALLOON Sterile Disposable Balloon Catheters for Peripheral Angioplasty and Valvuloplasty, design options,

Figure 2: Appearance of CRISTAL BALLOON Sterile Disposable Balloon Catheters for Peripheral Angioplasty and Valvuloplasty, design options, in the terminal sterilization packaging

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Filled balloon diameter
Balloon length
Filling/discharge hole
Body outer diameter
Gold connection ring
Total length
Angiographic hole

There are several device models. They vary in balloon diameters, balloon lengths and total catheter lengths. Some models are designed for angioplasty (CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, CBP series) and valvuloplasty (CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, CBV and CBVP series).

All models are provided in the table below:

Table 1: Technical specifications of the device

Catheters for valvuloplasty – Petit (CBVP):

References	Balloon ($\varnothing \times L$ mm)	Length of the catheter (cm)	Guide wire
CBVP8x20/70	8,0 x 20	70	.035''
CBVP10x20/70	10,0 x 20	70	.035''
CBVP12x20/70	12,0 x 20	70	.035''
CBVP15x25/70	15,0 x 25	70	.035''

Catheters for peripheral angioplasty (CBP):

References	Balloon ($\varnothing \times l$ mm)	Length of the catheter (cm)	Guide wire
CBP8x3035/110	8,0 x 30	110	.035''
CBP10x30/110	10,0 x 30	110	.035''
CBP12x40/110	12,0 x 40	110	.035''

Catheters for valvuloplasty (CBV)

References	Balloon ($\varnothing \times l$ mm)	Length of the catheter (cm)	Guide wire
CBV12x30/110	12,0 x 30	110	.038''
CBV15x30/110	15,0 x 30	110	.038''
CBV15x40/110	15,0 x 40	110	.038''
CBV16x40/110	16,0 x 40	110	.038''
CBV18x30/110	18,0 x 30	110	.038''
CBV18x40/110	18,0 x 40	110	.038''

References	Balloon ($\varnothing \times l$ mm)	Length of the catheter (cm)	Guide wire
CBV18x45/110	18,0 x 45	110	.038''
CBV20x30/110	20,0 x 30	110	.038''
CBV20x40/110	20,0 x 40	110	.038''
CBV20x45/110	20,0 x 45	110	.038''
CBV23x45/110	23,0 x 45	110	.038''
CBV25x40/110	25,0 x 40	110	.038''
CBV25x50/110	25,0 x 50	110	.038''
CBV28x50/110	28,0 x 50	110	.038''
CBV30x60/110	30,0 x 60	110	.038''
CBV35x60/110	35,0 x 60	110	.038''
CBV40x40/110	40,0 x 40	110	.038''

Table 2.1: Technical specifications of the device

	CBV1 8x20/ 70	CBVP 10x20 /70	CBVP 12x20 /70	CBVP 15x25 /70	CBP 8x30x35/ 110	CBP 10x30/ 110	CBP 12x40/ 110	CBV 12x30/ 110	CBV 15x30/ 110	CBV 15x40/ 110	CBV 16x40/ 110	CBV 18x30/ 110	CBV 18x40/ 110	CBV 18x45/ 110
Total length, mm	700±10	700±10	700±10	700±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10
Effective length, mm	635±10	635±10	635±10	635±10	1020±10	1020±10	1020±10	1020±10	1020±10	1020±10	1020±10	1020±10	1020±10	1020±10
Distal part diameter and length, mm	0.90x1.25	0.90x1.25	1.00x1.40	1.00x1.40	0.90x1.35	0.90x1.35	1.00x1.40	1.00x1.40	1.00x1.40	1.00x1.40	1.00x1.40	1.00x1.40	1.00x1.40	1.00x1.40
Balloon overall dimensions (max. diameter x length), mm	8x20	10x20	12x20	15x25	8x30	10x30	12x40	12x30	15x30	15x40	16x40	18x30	18x40	18x45
Balloon maximum volume (hyper-compatible balloon, without tolerances), mL	2	3	4	8	2,5	3,5	5	5	9	10	12	12	14	15
Inner diameter (main lumen), mm	0.90	0.90	1.00	1.00	0.90	0.90	1.00	1.00	1.00	1.00	1.00	1.00	1.00	1.00
Dual lumen thickness/width, mm	0.35	0.35	0.40	0.60	0.35	0.35	0.40	0.40	0.60	0.60	0.60	0.60	0.60	0.60
Weight, g	From CBVP8X20/70 to CBV 40X40/110 minimum weight 4.5 g, maximum weight 15 g													

Table 2.2: Technical specifications of the device

	CBV 20x30/110	CBV 20x40/110	CBV 20x45/110	CBV 23x45/110	CBV 25x40/110	CBV 25x50/110	CBV 28x50/110	CBV 30x60/110	CBV 35x60/110	CBV 40x40/110
Total length, mm	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10	1100±10
Effective length, mm	1020±10	1020±10	1020±10	1030±10	1030±10	1030±10	1030±10	1030±10	1030±10	1030±10
Distal part diameter and length, mm	1.00x1.40	1.00x1.40	1.00x1.40	1.10x1.70	1.10x1.70	1.10x1.70	1.10x1.70	1.10x1.70	1.10x1.70	1.10x1.70
Balloon overall dimensions (max. diameter x length), mm	20x30	20x40	20x45	23x45	25x40	25x50	28x50	30x60	35x60	40x40

Balloon maximum volume (hyper-compatible balloon, without tolerances), ml.	14	17	19	25	29	34	43	56	73	68
Inner diameter (main lumen), mm	1.00	1.00	1.00	1.10	1.10	1.10	1.10	1.10	1.10	1.10
Dual lumen thickness/width, mm	0.60	0.60	0.60	0.80	0.80	0.80	0.80	0.80	0.80	0.80
Weight, g	From CBVF8X20/70 to CBV 40X40/110 minimum weight 4.5 g, maximum weight 15 g									

Table 2.3 Technical specifications of the device

Surface condition (clean, smooth, free from scratches, etc.)	Yes
Min. tear force, N	3
Maximum tear strength of the catheter, N	≥ 25
Maximum tear strength of the hub, N	> 35

Unspecified acceptable values: ±1%

Device appearance: The device surface has no foreign particles, bubbles or stains. The surface is smooth and free from technological and surface defects, thus minimizing vessel trauma during use. The catheter is biologically safe and compatible with human body fluids. The device is not implantable and is removed from the human body after use.

7. Description of accessories and medical devices used with the reported device:

CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, can be used with Indeflator Syringes applied for balloon inflation/deflation and compatible guidewires for device passing through vessels.

The diameter of the largest guidewire used with the CRISTAL BALLOON Sterile Disposable Balloon Catheter for Peripheral Angioplasty and Valvuloplasty, design options, (if any), is provided in the table below.

Table 3. Size options of guidewires used with CRISTAL BALLOON Sterile Disposable Balloon Catheters for Peripheral Angioplasty and Valvuloplasty, design options.

Part No.	To be used with a guidewire
CBVP8x20/70	0.89 mm
CBVP10x20/70	0.89 mm
CBVP12x20/70	0.89 mm
CBVP15x25/70	0.89 mm
CBP8x3035/110	0.89 mm
CBP10x30/110	0.89 mm
CBP12x40/110	0.89 mm
CBV12x30/110	0.97 mm
CBV15x30/110	0.97 mm
CBV15x40/110	0.97 mm
CBV16x40/110	0.97 mm
CBV18x30/110	0.97 mm
CBV18x40/110	0.97 mm
CBV18x45/110	0.97 mm
CBV20x30/110	0.97 mm
CBV20x40/110	0.97 mm
CBV20x45/110	0.97 mm
CBV23x45/110	0.97 mm
CBV25x40/110	0.97 mm
CBV25x50/110	0.97 mm
CBV28x50/110	0.97 mm
CBV30x60/110	0.97 mm
CBV35x60/110	0.97 mm
CBV40x40/110	0.97 mm

Note: Unspecified tolerances are ±1%.

8. Information on the use of a medicinal product, materials of animal or human origin in medical device:

Not applicable. CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, do not contain any medicinal products or pharmaceutical substances or materials of animal or human origin.

9. Information on the procedure for installation, assembly, setting, and calibration required to commission the medical device:

METHOD OF APPLICATION

9.1 Catheter preparation

- Inspect the package thoroughly: it should not be opened or damaged. CRISTAL balloon catheter for peripheral angioplasty and valvuloplasty should be placed on a special insert inside its package.
- Open the package and remove the catheter using an aseptic technique. Inspect the catheter visually upon removal to ensure absence of any damages.
- Remove a sterile protective cap carefully from the distal part of the catheter on which the balloon is located. Remove the mandrel from the catheter.
- Avoid bending and kinking of the catheter while performing the above actions.

9.2 Catheter washing

- Fill a syringe having sufficient capacity with a sterile radiopaque agent (diluted with sterile saline solution, if necessary).
- Never use air or any gas to fill the balloon during washing and particularly during the procedure.
- Fill the balloon with radiopaque agent after attaching the syringe to the side port of the balloon. Inflate and deflate the balloon by means of a one-way stopcock. While holding the balloon in an inclined position, tap at it gently to remove air bubbles. With the syringe above the balloon the end of which is directed downward, evacuate the radiopaque agent and air. Close the stopcock between the syringe and Luer lock connection, then remove the syringe and evacuate air from it.
- Repeat the procedure until no air bubble remains in the balloon during inflation and in the syringe during deflation. The balloon of the catheter should be completely filled with radiopaque agent.

9.3 Balloon testing

- Evacuate the air completely from the dilatation balloon by rinsing it with radiopaque agent following the manufacturer's instructions contained in the catheter's IFU enclosed. Attach the valve Y-shaped connector, put the introducer selected to perform the procedure into the central lumen of the Cristal Balloon catheter and rinse this lumen with the sterile saline solution.

- To ensure proper functioning of the dilatation balloon, inflate it to the recommended pressure using the appropriate pressure gauge, then check the balloon. Empty the balloon completely and wind it around the tube during deflation to minimize its profile, if necessary.

9.4 Recommendations for choosing a guidewire and a balloon for angiography

A guidewire used should allow for rinsing the inner catheter lumen during the procedure. According to the general PTA recommendations, the guidewire must remain inside the catheter during the procedure not being removed from it. In case of significant catheter kinking during its advancement in some vessels (e.g. at the level of aortic bifurcation in contralateral femoral approach) the remaining lumen at the kink level may be insufficient to perform angiographic control. In such cases, that can be envisaged in advance, it is recommended to use a guidewire slightly smaller in diameter than the recommended maximum diameter, or a moveable-core guidewire that allows infusion through its own lumen. Moreover, the removal of the guidewire poses a risk during the procedure, since it must remain in place to provide access to the site of angioplasty.

The size and strength of the balloon should be selected in accordance with the nature of the lesion to be treated. It is recommended to use the balloons, the strength of which is designed for higher pressure, if stenosis is in the calcified zone. It is recommended to control the pressure being used. This control aims at preventing excessive inflation of the balloon, in which it may lose compliance with the selected nominal size, and then burst or disconnect from the catheter. The consequences of such cases can be dangerous to the patient's health.

Never use air or any gas to fill the balloon during washing and particularly during the procedure. The use of gas may result in the following:

- gas compressibility creates technical difficulties at inflation and deflation due to compression and expansion of gases;
- there is a danger of gas embolism in case of gas leakage from the balloon; in such case the volume of gas is much larger than the volume of the balloon filled with a fluid;
- the use of a fluid containing radiopaque agent allows for the angiographic control.

For the above reasons, it is necessary to completely remove residual gas from the balloon system.

9.5 Use

- Set the introducer and guiding catheter (in case the latter is used) according to the standard method.
- Rinse the inner catheter lumen with the sterile saline solution.
- Following complete deflation of the balloon, gently introduce the CRISTAL balloon catheter equipped with a guidewire for peripheral angioplasty and valvuloplasty through a haemostatic valve. If resistance is met, suspend the insertion of the catheter and check it for the entire length. Clamp the valve around the CRISTAL Balloon Catheter for peripheral angioplasty and valvuloplasty to prevent reverse flow through the guiding catheter, but not too tight in order not to impede the catheter advancement and cause compression that may

interfere with inflation and deflation of the balloon. Purge air from the valve. Rinse guiding catheter with the sterile saline solution.

- Advance the CRISTAL Balloon Catheter for peripheral angioplasty and valvuloplasty the end of the guiding catheter avoiding the catheter spinning around the guidewire.
- Advance the guidewire under angiographic control by directing it at the arterial bifurcation sites, until it crosses the stenotic area to be dilated. Advance the CRISTAL Balloon Catheter for peripheral angioplasty and valvuloplasty over the wire until crossing stenotic site. Set the balloon in this site using the radiopaque marker (-s). Advancement of the balloon through the stenotic site and its setting shall be carried out gradually, the actions can take a significant amount of time.
- Fill the balloon with radiopaque agent applying the appropriate pressure. Never use air or gas to inflate the balloon. The pressure shall in no case reach the level at which the balloon can burst, it should be below or equal to the maximum recommended pressure. The pressure should be monitored by an appropriate pressure gauge.
- The patient's condition should be monitored by appropriate means, depending on the exposed area. Under angiographic control allowing to evaluate the dilatation process, deflate the balloon completely when dilatation is achieved, or if the patient's condition requires it. Keep the guidewire in place.
- Following dilatation, check the results and make sure there are no complications introducing radiopaque agent into the guiding catheter.
- After angioplasty, remove the deflated balloon from the target area. Keep the guidewire in place passing through the stenosis within 10 minutes after angioplasty. After making sure that there are no procedure-related lesions in the dilated area, remove the guidewire from the CRISTAL Balloon Catheter for peripheral angioplasty and valvuloplasty from the guiding catheter. Then remove the guiding catheter. The valve introducer is usually left in place until the patient's general and hemodynamic status stabilizes. Provide hemostasis at the puncture site.
- In case of any difficulties when removing the balloon (balloons with larger valvuloplasty catheters), make sure that the balloon is empty; inflate the balloon if necessary and deflate again for a better shaping.

PRECAUTIONS FOR USE

- Do not use if the package is opened or damaged. The devices are sterile if the package is undamaged.
- CRISTAL BALLOON Sterile Disposable Balloon Catheters for Peripheral Angioplasty and Valvuloplasty are intended for single use. Do not reuse. Reuse of the device may result in a high risk of contamination to the patient and performance degradation. Do not resterilize.
- Store in a dry and dark place at room temperature. Do not use past the expiration date.
- These devices should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology and/or cardiology.

- Never move an intravascular device if resistance is met until the causes of resistance have been determined by means of angiographic control. Movement against resistance may result in damage to the catheter or injury to the patient.
- Use a suitable pressure gage to monitor pressure during the procedure.
- Observe the recommended inflation pressure specified on the device label.
- Never use air or gas to inflate the balloon.
- Follow the instructions for use provided with the agents used in the catheter.

9.1 Principle of operation:

CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, are coaxial double-lumen catheters with a balloon at the distal end. The outer lumen is intended for inflating and deflating the dilatation balloon, and the inner lumen – for the guidewire allowing to guide the catheter and facilitating its passing through stenosis to be dilated. The balloon is designed so that its diameter and length correspond to the set pressure. One or two radiopaque markers allow for the right position of the catheter's balloon under angiographic control. Standard Luer lock hubs ensure connection with the material.

The device is available for delivery in several modifications. They vary only in balloon diameters, balloon lengths and total catheter lengths. Specific recommendations are intended for angioplasty (Peripheral CRISTAL BALLOON CBP) and valvuloplasty (valvuloplasty CRISTAL BALLOON CBV and petit valvuloplasty CRISTAL BALLOON CBVP).

Each catheter is placed on a stand. The balloon is protected by a special sheath. Then the device is packaged in a removable pouch and a carton containing 2 catheters and the instructions for use.

The CRISTAL BALLOON catheters are disposable devices and are sterilized with ethylene oxide. Pursuant to Rule 7 of the European Directive 93/42/EEC these are Class III devices. In fact, these products are surgically invasive devices intended for short-term use (less than 24 hours) in contact with the central circulatory system (circulating blood).

These devices shall meet the requirements with regard to biocompatibility of materials, compatibility with the introducer sheaths being used, navigation ability and visibility when performing interventional radiology procedures, mechanical stability and balloon inflation. Preclinical evaluation results are given in Chapter 17 of the Technical File. A gap analysis was performed based on the biocompatibility data, and it showed that additional biocompatibility testing was required to ensure compliance with the essential requirements of Directive 93/42/EEC Rule 7.

10. Requirements to the premisses intended for installation of the medical device:

Ambient temperature from +20 °C to +25 °C, relative humidity from 35 % to 65 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa, away from direct sunlight. The device is temporarily operated inside the human body at a temperature ranging from +36°C to +42°C. The device is operated inside the human body for less than 24 hours.

11. Information for verifying the correctness of installation of the medical device and possibility of safe operation/use, including:

a) **Scope and frequency of maintenance, including cleaning and disinfection of the medical device:** CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, are not subject to maintenance, cleaning or disinfection.

b) **The list of information, keys and access passwords required for assembly and maintenance of the medical device provided by the manufacturer:** CRISTAL BALLOON sterile disposable balloon catheter for peripheral angioplasty and valvuloplasty, design options, do not contain software, information, keys or access passwords required for assembly and are not subject to maintenance.

c) **List of consumables, as well as procedure for use and replacement thereof:** CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, do not contain any consumables.

d) **Necessity of the medical device calibration to ensure appropriate and safe operation throughout the service life:** Not applicable. CRISTAL BALLOON sterile disposable balloon catheter for peripheral angioplasty and valvuloplasty, design options, do not require calibration to ensure an appropriate and safe operation throughout the service life.

e) **Methods of mitigating risks related to the installation, calibration or maintenance of the medical device:**

To reduce the risks related to the placement of CRISTAL BALLOON sterile disposable balloon catheter for peripheral angioplasty and valvuloplasty, design options, read the Instructions for Use for this device carefully as well as the Instructions for Use for a device used with the reported device. For additional information please refer to Paragraph 9 Information on the procedure for installation, assembly, setting, and calibration required to commission the medical device.

The device is intended for single use and does not require maintenance. Calibration is also not required; all the necessary settings are preset by the manufacturer.

f) **Information on assembly, adjustment, setting, calibration and other actions required to commission the medical device and for correct use (operation) thereof:** Not applicable. CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, do not require assembly, adjustment, setting, calibration and other actions necessary for commissioning the medical device and its correct use (operation).

g) **Information on the list of the main specifications for the use of the medical device during transportation and storage conditions:**

Transportation conditions: Ambient temperature from -10 °C to +35 °C, relative humidity from 35 % to 80 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa, away from direct sunlight.

Storage conditions: Store the device in its original package in a dry and clean place, away from heat sources and direct sunlight, at controlled temperature. The shelf life is indicated on each device label. Ambient temperature from +15 °C to +27 °C, relative humidity from 35 % to 80 % (non-condensing), atmospheric pressure from 500 hPa to 1,060 hPa.

12. **Information on the sterile condition of the medical device, sterilization method and procedure in case of breach of sterile packaging:** Each CRISTAL BALLOON sterile disposable balloon catheter for peripheral angioplasty and valvuloplasty, design options, is supplied as

sterile disposable device. CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, are packed in a terminally sterilized package and sterilized with ethylene oxide. The device should not be used and should be disposed of immediately if the package sterility is compromised. Reuse of the device may result in a high risk of contamination to the patient and performance degradation.

13. Information on processing of the medical device for reuse: The device is supplied sterile and shall not be subject to processing, disinfection or re-sterilization. After use, the device shall not be reused and shall be disposed of immediately. Reuse of the device may result in a high risk of contamination to the patient and performance degradation.

14. Information required identifying medical devices to obtain a safe combination and information on the known restrictions of concurrent use of medical devices:

Paragraph 7 Description of accessories and medical devices used with the reported device contains information on the devices which may be used with CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options. Please refer to the instructions for use for each BALT medical device for all the necessary information on the use and known limitations.

15. Information on the nature, type and intensity and distribution of radiation during the use of the medical device:

Not applicable. CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, do not require radio-frequency energy and therefore do not create electromagnetic field and interference induced thereby.

16. Information on the precautions taken if:

A) The medical device is defective, malfunctions or features deviations of operation when these issues may affect safety of the medical device, including visually detectable defects:

Please refer to Paragraph 9 for the information on identifying potential failures, malfunctions or deviations in operation that may affect safety of the medical device, including those identified by means of visual inspection for the possibility of their elimination.

B) The medical device is exposed to external factors related to the use of the medical device with other medical devices:

Refer to the instructions for use of the devices being used with CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, (Paragraph 7 of the Instructions for Use);

b) A risk of electromagnetic interference created by the medical device affecting other medical equipment: Not applicable. CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, do not require radio-frequency energy and therefore do not create electromagnetic field and interference induced thereby.

17. Warnings and precautions taken by the consumers when using the medical device containing a human medicinal product, carcinogenic, mutagenic or toxic materials, possible extraction or washout whereof results in sensitization, allergic reaction or negatively affects fertility: Not applicable. CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design options, do not contain any medicinal products or pharmaceutical substances or materials of animal or human origin, carcinogenic, mutagenic or

toxic materials, possible extraction or washout whereof results in sensitization, allergic reaction negatively affects fertility.

18. Warnings and precautions taken by the consumer when disposing of the medical device accessories and consumables used therewith:

Dispose of the device as Class B waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10. After physical decontamination and waste appearance change excluding the possibility of reuse, Class B waste can be accumulated, temporarily stored, transported, destroyed and buried together with Class A waste (epidemiologically safe waste with the composition similar to that of municipal solid waste). The package of decontaminated Class B medical waste should have a label indicating that the waste has been decontaminated.

19. Information on the circumstances requiring consultation of the consumer with a health professional:

The device should only be used by medical specialists trained in the field of interventional neuroradiology and/or interventional radiology. Conditions of use: Medical treatment facilities, healthcare centers, only when prescribed by a medical specialist. In view of the foregoing, it follows that in any case, the consumer should consult with a health care professional to determine the necessity to perform the procedure, and as a consequence, the use of CRISTAL BALLOON sterile disposable balloon catheters for peripheral angioplasty and valvuloplasty, design optimized prior to undergoing the cerebral embolization.

20. Information on the first issue or the latest revision of the operating instructions:

Not applicable. The information on the first issue and the latest revision of the operating instructions is missing.

21. Procedure and conditions for the disposal or destruction of the medical device: Dispose of the device as Class B waste (epidemiologically hazardous waste) in accordance with SanPiN 2.1.7.2790-10. After physical decontamination and waste appearance change excluding the possibility of reuse, Class B waste can be accumulated, temporarily stored, transported, destroyed and buried together with Class A waste (epidemiologically safe waste with the composition similar to that of municipal solid waste). The package of decontaminated Class B medical waste should have a label indicating that the waste has been decontaminated.

22. Warranty

Shelf life – 5 years. Sterilization validity – 5 years.

Address for complaints:

Center for Perinatal Medicine, Limited Liability Company

Legal/actual address: 20/4 Tovarishchesky per., 109004, Moscow

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