
Medtronic Inc. 110 Medtronic Parkway Minneapolis, MN 55432 USA Tel: 1-800-633-8766	 CINDI A. BIGELBACH Notary Public State of Minnesota My Commission Expires January 31, 2022 <i>Cindi A. Bigelbach</i> 3/5/18	 Medtronic  Calero Gloria Regulatory Affairs Specialist Medtronic, Inc. Date <i>March 5, 2018</i>
--	---	--

Addition to the technical manual

1. Name of medical device	The defibrillation electrode of the Sprint Quattro Secure 6947M of various versions
2.1. Prescription of medical device	Electrode Sprint Quattro Secure 6947M is designed for stimulation, sensing, cardioversion and defibrillation
2.2. Operational principles of medical device	Device Description Sprint Quattro Secure® Model 6947M  Principle of Operation The leads feature an extendable and retractable helix electrode, silicone insulation, and parallel conductors. The four electrodes of the

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	lead are the helix, ring, RV coil, and SVC coil described as follows: 1. The helix electrode is electrically connected to the connector pin of the DF4 pin connector. 2. The ring electrode is electrically connected to the most proximal connector ring of the DF4 connector. 3. The RV coil electrode is electrically connected to the center ring of the DF4 connector. 4. The SVC coil electrode is electrically connected to the most distal ring of the DF4 connector. The DF4-LLHH (two low power and two high power contacts) four-pole high voltage (HV) inline connector on the lead facilitates device connection during implant. The DF4 connector pin has a color band indicator that may be used to visually confirm proper connection to the device. The RV and SVC coils deliver cardioversion and defibrillation therapies. An AccuRead analyzer cable interface tool (AccuRead tool) is in place on the lead as shipped to facilitate accurate electrical measurements during implant. The helix electrode can be actively fixed into the endocardium. The helix electrode can be extended or retracted by rotating the DF4 connector pin with one of the purple fixation tools included in the package. The purple tools are specifically sized for the DF4 pin. The basic functionality of the Sprint Quattro Secure Family of leads included in this evaluation is to sense, pace, and to either defibrillate or cardiovert arrhythmias in the ventricle. Consequently, these leads include both high power and low power electrodes and conductors. The functions of the low power electrodes and conductors are to receive ventricular electrical activity, conduct it to the ICD ventricular sense amplifier and to deliver pacing pulses from the ICD to the ventricle. The tip electrode (cathode) is designed to efficiently deliver pacing pulses for right ventricular pacing. Typically, pacing and sensing occur between the helix and ring electrodes. However, devices can be programmed to sense between the helix and RV coil electrode in some devices. The pace-sense ring has a higher surface area than the tip to improve pacing efficiency. The high power electrodes and conductors are designed to efficiently and reliably deliver high energy cardioversion/defibrillation pulses to the heart.
3. Indications and Contraindications for medical device	MRI conditions for use A complete SureScan defibrillation system is required for use in the MR environment. A complete SureScan defibrillation system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan defibrillation system components displays the MR Conditional symbol. Warning: Do not scan a patient without first programming MRI SureScan to On. Scanning the patient without programming MRI SureScan to On may result in patient harm or damage to the SureScan defibrillation system. Note: MRI SureScan cannot be programmed to On if the device is recommended for replacement. Patients and their implanted systems must be screened to meet the following requirements: Cardiology requirements: • no implanted lead extenders, lead adaptors, or abandoned leads • for patients who have multiple MR Conditional devices present, the MR labeling conditions for all implants are satisfied • no broken leads or leads with intermittent electrical contact as confirmed by lead impedance history • a SureScan defibrillation system implanted in the left or right pectoral region • the SureScan device is operating within the projected service life • pacing capture thresholds of ≤ 2.0 volts (V) at a pulse width of 0.4 milliseconds (ms) • no diaphragmatic stimulation at a pacing output of 5.0 V and at a pulse width of 1.0 ms in patients whose device will be programmed

to an asynchronous pacing mode when MRI SureScan is programmed to On.

Note: Patients experiencing atrial fibrillation may be scanned if all other pre-MRI scan requirements are satisfied.

Radiology requirements

Scanning patients who have multiple MR Conditional devices present is acceptable when the MR labeling conditions for all implants are satisfied.

The safety and reliability of the SureScan defibrillation system has been evaluated for scanning patients using MRI equipment that has the following operating characteristics:

Scanner type	Horizontal field, cylindrical bore, clinical system
Scanner characteristics	• Static magnetic field of 1.5 T operating at a frequency of 64 MHz, maximum spatial gradient of ≤ 20 T/m (2000 gauss / cm) • Gradient systems with maximum gradient slew rate performance per axis of ≤ 200 Tesla per meter per second (T/m/s)
Scanner operation	MRI radio frequency (RF) power – Normal Operating Mode. The whole body averaged specific absorption rate (SAR) must be ≤ 2.0 watts per kilogram (W/kg). The head SAR must be ≤ 3.2 W/kg.

Patient monitoring and rescue requirements


Continuous patient monitoring is required while MRI SureScan is programmed to On.

An external defibrillator must be immediately available while MRI SureScan is programmed to On.


Potential adverse events that may occur in the MRI environment

The SureScan defibrillation system has been designed to minimize the potential adverse events that may cause patient harm. The following potential adverse events may occur in the MR environment:


- lead electrode heating and tissue damage resulting in loss of sensing or capture or both
- spontaneous tachyarrhythmia occurring during the scan that is not detected and treated because tachyarrhythmia detection is suspended while MRI SureScan is programmed to On
- potential for VT/VF induction when the patient is programmed to an asynchronous pacing mode during MRI SureScan
- device heating resulting in tissue damage in the implant pocket or patient discomfort or both
- MR-induced stimulation on leads resulting in continuous capture, VT/VF, hemodynamic collapse, or all three
- damage to the device or leads causing the system to fail to detect or treat irregular heartbeats, or causing the system to treat the patient's condition incorrectly
- damage to the functionality or mechanical integrity of the device resulting in the inability of the device to communicate with the programmer

	• movement or vibration of the device or leads resulting in dislodgment
4. Information on potential consumers of medical device	This medical device is intended for use in clinical practice by qualified physicians and medical specialists.
5. Description of the main functional parts of medical device	Photographic image of the Sprint Quattro Secure MRI SureScan 6947M.  The image shows the packaging and components of the Medtronic Sprint Quattro Secure MRI SureScan 6947M. The main box is white with black and red text, featuring the Medtronic logo and product name. It also includes a diagram of the device and various safety warnings. In front of the box, there is a clear plastic bag containing the device, which is a long, thin, flexible lead. Several other items are visible on the surface, including a small white box, a syringe, and some papers. The background is a plain, light-colored surface.

Photographic image lead with a radiopaque anchoring sleeve, stylet, and AccuRead tool


The scheme lead with a radiopaque anchoring sleeve, stylet, and AccuRead tool:

Feature	6947M MRI SURESCAN	
	55 cm	62 cm
Lead length		
Length, mm	551,2	619,8
Diameter distal anchoring sleeve, mm		5,5
Diameter lead, mm		2,5
Mass, g	8	9

The scheme radiopaque anchoring sleeve:

Feature	Value
Length, mm	20,3
Outer diameter, mm	5
Inside diameter, mm	3
Mass, g	0,1

The Sprint Quattro Secure MRI SURESCAN lead lengths are listed in table:

Device	Model number(s)	Lead length
Sprint Quattro Secure MRI SURESCAN	6947M	55 cm
Sprint Quattro Secure MRI SURESCAN	6947M	62 cm

Specification

Feature	6947M MRI SURESCAN
General	
Polarity	Quadripolar
Sensing Polarity	True Bipolar
Fixation	Active Helix
Connector	DF4
MR Conditional	Yes
Lead Diameters	
Body	2,8 mm (8.6 Fr)
Coils	2,8 mm (8.6 Fr)
Tip	2,8 mm (8.6 Fr)
Active Fixation	Helix
Fixation Depth	1.8 mm
Helix Diameter	1,45 mm (0.057 in)
Tip Electrode Wire Diameter	0,25 mm (0.010 in)
Helix	1.4 mm
Helix length (extended)	1.8 mm
Introducer Size	without guide wire: 3 mm (9,0 French) with guide wire: 3.7 mm (11.0 French)
Pacing Details	
Tip to Ring Spacing	8 mm
Tip/Cathode Surface Area	5.7 mm ²
Ring/Anode Surface Area	25.2 mm ²
High Voltage Details	
Electrode Configuration	Dual Coil

Defibrillation Coil	RV and SVC
Coil Treatment	Silicone Backfill
Inter-electrode Spacing	
Tip to Ring	8 mm
Tip to RV Coil	12 mm
Tip to SVC	180 mm
Electrode Surface Area	
Tip	5.7 mm ²
Ring	25.2 mm ²
RV Coil	614 mm ²
SVC Coil	860 mm ²
Electrode Shadow Area	
RV Coil	506 mm ²
SVC Coil	709 mm ²
Electrode Lengths	
RV Coil	5.7 cm
SVC Coil	8 cm
Connector	
High Voltage	DF4
Pace / Sence	DF4

Sprint Quattro Secure MRI SureScan 6947M

D.C. Resistance, (ohms)						
Lead Length	tip to m-4 pin a-j (ohms)	ring to connector ring b-i ohms max	distal RV to connector c- h ohms max	proximal RV to connector d- h ohms max	distal SVC to connector e-g ohms max	proximal SVC to connector f-g ohms max
55 cm	24,7±4,9	1,2	1,38	1,27	1.00	0.83
62 cm	27,9±5,6	1,4	1,53	1,42	1.15	0.98

Conductor torsional strength: a minimum torque of 1.977 $\mu\text{N}\cdot\text{m}$ (.28 inch-ounce).


Helix assembly torsional strength: a minimum torque of 1.271 $\mu\text{N}\cdot\text{m}$ (.18 inch-ounce).

Torque requirement: 2.118 $\mu\text{N}\cdot\text{m}$ (.3 inch-ounce of torque)


Connector torsional strength: a minimum torque of 1.977 $\mu\text{N}\cdot\text{m}$ (0.28 inch-ounce).

Maximum number of rotations to extend or retract the helix electrode:

Lead length	Number of rotations
	6947M
55 cm	18


6. Description of all accessories, medical devices or non-medical devices, intended for use with the medical device submitted for registration (if applicable)	62 cm		20	
	All unspecified dimensional tolerances are assumed to be ± 5%			
	Accessories:			
				
	Stylet guide	Fixation tool	Vein lifter	
				
	Slit anchoring sleeve		AccuRead analyzer cable interface tool (ACI)	
			Stylet	
	<u>Schemes of accessories:</u>			
	The scheme stylet:			
				
Feature		6947M MRI SURESCAN		
Lead length		55 cm	62 cm	
Working part length, mm		572	641,1	
Working part diameter, mm		0,4		
Diameter of the stylet handle, mm		7		
Mass, g		0,6	0,7	

The scheme purple stylet guide:


Characteristics	Value
Length, mm	11,4
Diameter of the wide part, mm	7,8
The diameter of the narrow part, mm	3,5
Mass, g	0,2

The scheme slit anchoring sleeve:


Characteristics	Value
Length, mm	21,8
Outer diameter, mm	$4,6 \pm 0,08$
Inside diameter, mm	$2,9 \pm 0,08$
Width, mm	8,7
Mass, g	0,2

The scheme AccuRead tool:


Characteristics	Value
Length, mm	25
Width, mm	12,1
Height, mm	8,9
Working part diameter, mm	1,7
Mass, g	1,1

The scheme vein lifter:

Characteristics	Value
Length, mm	63
Width, mm	9,5
Height, mm	10
Mass, g	1,2

The scheme purple fixation tools:

Characteristics	Value
Length, mm	49
Width, mm	37
Height, mm	3,9
Mass, g	3,1

All unspecified dimensional tolerances are assumed to be $\pm 5\%$.

7. List and description of medical device materials with direct or indirect contact with the patient's body (human body)

The SPRINT QUATTRO SECURE MRI SURESCAN 6947M lead is composed of the following materials and components that have the potential for direct and/or indirect patient body tissue/fluid contact.

Table of Tissue Contacting Materials:

Name of the product or part and accessories	Materials used in the manufacture	Type of contact with the body
Electrode implantable quadropole ventricular for defibrillation Sprint Quattro Secure MRI SureScan 6947M, dimensions: 55 cm, 62 cm in composition	Top coating of the electrode: Polyurethane Seal zone: Polyetheretherketone Distal sleeve, Proximal sleeve, Four-pole lead connector, endcap, MCRD cap: Polyurethane Helix: Platinum / Iridium Indicator ring, Drive shaft, Electrode ring, Helix electrode: Platinum / Iridium Axle connection, Proximal connector sleeve, Distal connector sleeve: Silastic silicone rubber Slit anchoring sleeve: Titanium of the 1st class Binding agent for steroid: Silicone rubber Steroid: Dexamethasone acetate Sodium phosphate dexamethasone Total amount not more than 1.0 mg Anchoring sleeve: Silicone rubber White dye: titanium dioxide Adhesive: - Silicone Adhesive (polydimethylsiloxane) - Polyurethane adhesive - Polyurethane adhesive	Permanent contact (more than 30 days) with the internal environment and body tissues

	Vein lifter	Acetal copolymer yellow	Short-term (less than 24 hours) contact with tissues and blood
	Stylets	Stilet material - stainless steel Stiletto Stick Material -Polycarbonate Stylette stain dye - purple and gray dyes White ink	No contact with the patient
	AccuRead tool	AccuRead cable interface analyzer tool (Analyzer Cable Interface) - acetal	No contact with the patient
	Purple fixation tools	Acetal Purple dye	No contact with the patient
	Purple stylet guide	Acetal Purple dye	No contact with the patient
8. The possibility and ways of integrating with other medical devices (device compatibility) (if applicable)	Consult the instructions for use for guidelines in this section to select implantable device that are compatible with the lead. The appropriate techniques for implanting the device may vary according to physician preference and the patient's anatomy or physical condition. Consult the product literature packaged with the implantable device for instructions on proper lead connections verification of sensing, pacing, cardioversion, and defibrillation efficacy. A complete SureScan defibrillation system is required for use in the MR environment. A complete SureScan defibrillation system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan defibrillation system components displays the MR Conditional symbol. Before performing an MRI scan, refer to the MRI Technical Manual for MRI-specific warnings and precautions. Do not use any implantable device with this lead without first verifying lead and connector compatibility. Device compatibility Although the lead conforms to the International Connector Standard for DF4, do not attempt to use the lead with any device other than a commercially available implantable defibrillator system with which it has been tested and demonstrated to be safe and effective. The potential adverse consequences of using such a combination may include, but are not limited to, undersensing cardiac activity and failure to deliver necessary therapy. Do not use another manufacturer's implantable device without demonstrated compatibility with Medtronic leads. If an implantable device is not compatible with a Medtronic lead, the result may be undersensing of cardiac activity, failure to deliver necessary therapy, or a leaking or intermittent electrical connection.		

9. Labelling and Packaging of medical device

Example Shelf Box Label

Medtronic

Sprint Quattro Secure MRI™ SureScan™ 6947M - 62 cm

en Lead / fr Sonde / de Elektrode / es Cable / nl Geleidsdraad / it Electrocatetere / no Ledning / sv Elektrod / el Ηλεκτρόδιο / da Elektrode / pt Electrodo / hu Vezeték / pl Elektroda / cs Elektroda / ru Электрод / sk Elektroda / tr Lead / ro Derivade / bg Електрод / sr Elektroda / it Elektroda / hr Odvod / et Elektrodijune / lt Elektrodas / lv Elektrods

Package contents: 1x Lead + Accessories + Product documentation

62 cm

DF4-LLHH

Series → Pace

Defibrillation

Steroid eluting

Retractable screen

MR Conditional

Transcatheter lead with two defibrillation electrodes

Lead introducer 3.0 mm (9.0 French)

Lead introducer with guide wire 3.7 mm (11.0 French)

REF 6947M62

SN TDK234203V

2019-01-27 YYMMDD

2017-02 YYMMYY

PIN 6947M00071

Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432 USA
Manufactured in
Vladivostok, Russia
www.medtronic.com

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

The following provides a diagram of the Standard Lead Package (Standard Lead Package) packaging configuration.

Primary Packaging:

Feature	6947M MRI SURESCAN	
	55 cm	62 cm
Lead length		
Primary packaging mass, g	94	102
Packaging dimensions (Length x Width x	205,6x180x24	

Height), mm

Secondary Packaging:

Feature	6947M MRI SURESCAN	
Lead length	55 cm	62 cm
Secondary packaging mass, g	394	402
Packaging dimensions (Length x Width x Height), mm	225x190x35	

The tolerance is $\pm 5\%$, unless otherwise specified

Sterile and Shelf Box Packaging Materials

The following tables provide the materials used for both the primary and secondary packaging:

Table 4: Sterile Package - Primary Packaging

	Package Part	Material
Sterile packaging	Tray	Polyethyleneterephthalateglycol (PETG)
	Tray coverage	High-density polyethylene Tyvek
	Adhesive layer	Adhesive

Table 5: Shelf Box - Secondary Packaging

Package Part	Material
Shelf Box	Paperboard
Shrink Wrap/ Tear Tape	Polyethylene / PVC (Polyvinyl chloride)

Information on the Russian-language sticker

For sale on the territory of the Russian Federation, the final packaging is accompanied by an additional sticker.

Products come with a sticker that can contain the following information:

- The name of the medical device in Russian;
- Name of the manufacturer's factory;
- Country of origin;
- UPN code (unique product code);
- CFN code (model number);
- Date and number of the Registration Certificate;
- Mark of conformity in Russia;
- Warning and informative signs and inscriptions (if applicable).

Sample sticker:

	Name of medical device
---	------------------------

	<table><tr><td></td><td>Manufacturer:</td><td>Code UPN</td></tr><tr><td></td><td>Country of origin:</td><td>Code CFN</td></tr><tr><td></td><td colspan="2">Manufacturer representative in Russia:</td></tr><tr><td></td><td colspan="2">Registration Certificate № ____ of ____</td></tr><tr><td></td><td>Date of production / Series / Expiration date: see packaging</td><td>See Instructions See symbols on the package</td></tr><tr><td></td><td colspan="2">Not toxic Pyrogen-free</td></tr></table>		Manufacturer:	Code UPN		Country of origin:	Code CFN		Manufacturer representative in Russia:			Registration Certificate № ____ of ____			Date of production / Series / Expiration date: see packaging	See Instructions See symbols on the package		Not toxic Pyrogen-free	
	Manufacturer:	Code UPN																	
	Country of origin:	Code CFN																	
	Manufacturer representative in Russia:																		
	Registration Certificate № ____ of ____																		
	Date of production / Series / Expiration date: see packaging	See Instructions See symbols on the package																	
	Not toxic Pyrogen-free																		
10. List of materials of animal and (or) human origin, indicating information on their biocompatibility and safety, on selection of sources (donors), sampling, processing, storage and handling of these materials (if applicable)	Not applicable.																		
11. References to previous modifications of medical device or similar modifications of medical devices	Not applicable																		
12. Information in accordance with the data of the state register of medicinal products for medical use about those contained in the medical device	For the steroid and the coating material of the helix: Registered in Russia: LSR-002453/08 from 03.04.2008 dexamethasone acetate dexamethasone sodium phosphate FS-000453 12.12.2012 Sanofi Chimie, Le Bourg, 63480 Vertolaye, France Manufacturing process The drug containing component of steroid-eluting devices (the MCRD) is manufactured by combining drug DXP with medical grade silicone. Silicone, mixed with steroid, serves as a binder that holds but does not chemically interact with the steroid allowing steroid to																		

	be released upon exposure to fluid. MCRDs containing silicone are manufactured using one of three types of medical grade silicone. Different silicones are used to provide manufacturing properties appropriate for each MCRD design. The general steps in forming a silicone MCRD are as follows: • Silicone adhesive or elastomer is blended with steroid. • Silicone-steroid mixture is transferred to mold cavities for curing. • After curing and cooling, molded pieces are removed from the mold and trimmed (if required) to final dimensions for each type of MCRD. 12:13																																																
13.Sterilization method, Sterilization Validation, Packaging Validation for Sterile devices	Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This device is for single use only and is not intended to be resterilized. Model 6947M leads are provided sterile and are sterilized using the Ethylene Oxide (EO) method. The 30-minute 100% EO sterilization process is an “all-in-one” chamber process that is conducted internally by Medtronic at a variety of manufacturing facility locations, including MVC (Villalba) and MSO (Singapore) facility locations that model 6947M leads is manufactured and sterilized at. Sterilant residual comply with the requirements specified in ISO 10993-7. Routine sterilization parameters are listed in Table 6 below. Table 6: Routine Sterilization Parameter Specifications <table><tr><th>Cycle Phase</th><th>Sterilization Parameter</th><th>Acceptance Criteria</th></tr><tr><td rowspan="9">Conditioning</td><td>Vacuum Rate</td><td>7-13 kPa/minute</td></tr><tr><td>Pressure at the beginning of leak check</td><td>14.0-20.0 kPa</td></tr><tr><td>Leak Check Duration</td><td>15 ±2 minutes</td></tr><tr><td>Leak Check Vacuum Loss</td><td><4 kPa</td></tr><tr><td>Time</td><td>30 - 35 minutes</td></tr><tr><td>Number of Humidity Pulses</td><td>4</td></tr><tr><td>Chamber Temperature№</td><td>45 - 55°C</td></tr><tr><td>Chamber Relative Humidity№</td><td>55 - 85 %RH</td></tr><tr><td>Pressure, just prior to gas injection</td><td>8-14 kPa</td></tr><tr><td rowspan="6">Exposure</td><td>Gas Inlet Temperature</td><td>15 - 75°C</td></tr><tr><td>Exposure Time</td><td>27 - 33 minutes</td></tr><tr><td>Chamber Temperature & Load Probel</td><td>45 - 55°C</td></tr><tr><td>Pressure</td><td>55 - 75 kPa</td></tr><tr><td>Gas ConcentrationI</td><td>> 720 mg/L</td></tr><tr><td>EO Gas Cartridge Weight</td><td>127 ± 5 grams</td></tr><tr><td rowspan="2">Post Vacuum</td><td>Vacuum Rate</td><td>1.0 - 2.0 kPa/minute</td></tr><tr><td>Pressure</td><td>20.0-30.0 kPa</td></tr><tr><td>Air Flush</td><td>Time</td><td>> 25 minutes</td></tr><tr><td rowspan="2">Aeration</td><td>Temperature</td><td>44 - 56°C</td></tr><tr><td>Time</td><td>8-16 hours (1X, 2X, 3X, or 4X)</td></tr></table>	Cycle Phase	Sterilization Parameter	Acceptance Criteria	Conditioning	Vacuum Rate	7-13 kPa/minute	Pressure at the beginning of leak check	14.0-20.0 kPa	Leak Check Duration	15 ±2 minutes	Leak Check Vacuum Loss	<4 kPa	Time	30 - 35 minutes	Number of Humidity Pulses	4	Chamber Temperature№	45 - 55°C	Chamber Relative Humidity№	55 - 85 %RH	Pressure, just prior to gas injection	8-14 kPa	Exposure	Gas Inlet Temperature	15 - 75°C	Exposure Time	27 - 33 minutes	Chamber Temperature & Load Probel	45 - 55°C	Pressure	55 - 75 kPa	Gas ConcentrationI	> 720 mg/L	EO Gas Cartridge Weight	127 ± 5 grams	Post Vacuum	Vacuum Rate	1.0 - 2.0 kPa/minute	Pressure	20.0-30.0 kPa	Air Flush	Time	> 25 minutes	Aeration	Temperature	44 - 56°C	Time	8-16 hours (1X, 2X, 3X, or 4X)
Cycle Phase	Sterilization Parameter	Acceptance Criteria																																															
Conditioning	Vacuum Rate	7-13 kPa/minute																																															
	Pressure at the beginning of leak check	14.0-20.0 kPa																																															
	Leak Check Duration	15 ±2 minutes																																															
	Leak Check Vacuum Loss	<4 kPa																																															
	Time	30 - 35 minutes																																															
	Number of Humidity Pulses	4																																															
	Chamber Temperature№	45 - 55°C																																															
	Chamber Relative Humidity№	55 - 85 %RH																																															
	Pressure, just prior to gas injection	8-14 kPa																																															
Exposure	Gas Inlet Temperature	15 - 75°C																																															
	Exposure Time	27 - 33 minutes																																															
	Chamber Temperature & Load Probel	45 - 55°C																																															
	Pressure	55 - 75 kPa																																															
	Gas ConcentrationI	> 720 mg/L																																															
	EO Gas Cartridge Weight	127 ± 5 grams																																															
Post Vacuum	Vacuum Rate	1.0 - 2.0 kPa/minute																																															
	Pressure	20.0-30.0 kPa																																															
Air Flush	Time	> 25 minutes																																															
Aeration	Temperature	44 - 56°C																																															
	Time	8-16 hours (1X, 2X, 3X, or 4X)																																															

			sterilization)
	Fans	Current (conditioning, exposure, aeration)	0.030 - 0.300 amps
	¹ Relative humidity & temperature requirements must be met at the end of conditioning - just prior to gas inject. ² During the first five minutes of exposure (directly after the completion of gas injection), temperature and gas concentration spikes or drops may exceed the process parameters while the chamber is equilibrating. Has been successfully qualified into the MVC and MSO 30-minute 100% EO sterilization process per manufacturing document POD003074 or any of Medtronic's equivalent 30-minute 100% EO sterilization processes. A Sterility Assurance Level (SAL) in excess of 10^{-6} is achieved when the model 6947M Lead product family is sterilized at MVC in load configuration "D" (or equivalent load configuration) per manufacturing document POD003074 (or equivalent process). Model 6947M Lead product family requires a minimum of 8 hours of forced heated aeration following 1X, 2X, 3X, and 4X sterilization exposures. The maximum aeration allowed is 16 hours of forced heated aeration following each sterilization exposure. All internal and external requirements have been met. Model 6947M Lead product family is qualified for up to 4X sterilization in the 30-minute 100% EO sterilization process.		
14. Requirements for maintenance and repair of a medical device	Not applicable. There are no maintenance requirements.		
15. Procedure and conditions for the disposal or destruction of a medical device	Consider the following information related to device disposal: - Abandoned leads should be capped to avoid transmitting electrical signals. - Follow local hospital and clinic packaging and shipping instructions of biohazardous materials. - Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. - Please return all removed or unused leads, or lead sections, to Medtronic for analysis and disposal. Note: Disposal of explanted devices or leads is subject to local, state, and federal regulations. After use, this product may be a potential biohazard. Unused, expired devices are classified as medical waste class «A» Catheters extracted from the patient's body are medical waste Class B.		
16. Information about manufacturer of medical device (name, location address, mailing address, phone number, address of place of production)	Legal Manufacturer: Medtronic, Inc. 710 Medtronic Parkway Minneapolis, MN 55432 USA Design Site: Medtronic, Inc. 8200 Coral Sea Street NE Mounds View, MN 55112 USA Manufacturing Facilities:		

	Medtronic Puerto Rico Operations Co., Villalba Road 149, Km. 56.3, Call Box 6001, Villalba, PR 00766, USA Medtronic Singapore Operations Pte. Ltd. 49 Changi South Avenue 2, Nasaco Tech Centre, Singapore 486056, Singapore Russian Authorized Representative: «Medtronic» LLC 123317, Moscow, Presnenskaya Naberezhnaya., 10 Tel.: (495) 580-73-77. Fax: (495) 580-73-78. www.medtronic.ru Russian Legal Name for Certificate: «Medtronic» LLC 123317, Moscow, Presnenskaya Naberezhnaya., 10 Tel.: (495) 580-73-77, Fax: (495) 580-73-78. www.medtronic.ru
17. Class risk depending on the potential risk of a medical product, in accordance with the nomenclature	Classification of medical device: class risk 3 B - products whose failure reduces the effectiveness or delays the treatment and diagnostic process in non-critical situations, or increases the burden on medical or maintenance personnel; Device of permanent contact - single use, the contact of which exceeds 30 days, long contact with blood and tissue, body fluids. After use, this product may be a potential biohazard; Unused, expired devices are classified as medical waste class «A»; Catheters extracted from the patient's body are medical waste Class B; Class of electrical safety, laser safety, and other types of hazards: not applicable; The device is sterile, not re-sterilized.
18.Type of contact with the body	Device of permanent contact – single use, the contact of which exceeds 30 days. Long term contact with blood and tissue, body fluids.
19.Methods and apparatus for disinfection and sterilization cleaning	Not Applicable – single use only
20.Shelf life	Shelf life 24 months (2 years) from the date of sterile packaging
21.Durability in the human body (if applicable)	More than 30 days (the doctor makes the decision independently)
22.Operating conditions	In a medical and preventative treatment facility. The device should be used only by or under the supervision of a trained physician.

		Service conditions	
	Operating temperature	from + 32°C (0/ +2°C) to + 42°C (-2/0°C)	
	Humidity rate	100% (-5/0%)	
	Atmospheric pressure	from 700 (0/+3,5) hPa to 1060 (-3,5/0) hPa	
23.Conditions of transportation and storage	Handling and storage/transportation instructions		
	Carefully observe these guidelines when handling or storing the device.		
	• Device handling		
	Before opening the sterile package tray, visually check for any signs of damage that might invalidate the sterility of the package contents.		
	• If the package is damaged		
	The device packaging consists of an outer tray and inner tray. Do not use the device or accessories if the outer packaging tray is wet, punctured, opened, or damaged. Return the device to Medtronic because the integrity of the sterile packaging or the device functionality may be compromised. This device is not intended to be resterilized.		
	• Device temperature		
	Do not store device above 40°C ± 0,5°C		
	• “Use by” date		
	Use device by the "Use by" date printed on the label. Resterilizing a device does not extend the “Use by” date.		
There are no other special conditions related to the transport and storage of this product.			
		Conditions of transportation	Conditions of storage
	Temperature	From – 18 °C to + 60 °C	no more 40 ± 0,5°C
	Atmospheric pressure	From 700 hPa (0/+3,5) to 1060 hPa (-3,5/0)	From 700 hPa (0/+3,5) to 1060 hPa (-3,5/0)
	Humidity rate	not more than 60%	not more than 60%
24.Installation requirements	Not applicable. There are no installation requirements.		
25.List of information provided by the manufacturer (manufacturer) of the medical device, keys, access passwords, programs required for installation, commissioning, operation and maintenance of the	Not applicable.		

medical device																			
26. Warranty	Disclaimer of warranty and general warning Although the Medtronic lead, hereafter referred to as “Product”, has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. Warnings contained in the Product labeling provide more detailed information and are considered an integral part of this DISCLAIMER OF WARRANTY. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect, or failure of the Product, whether the claim is based on warranty, contract, tort, or otherwise. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in any conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid. General warning Medtronic implantable leads are implanted in the extremely hostile environment of the human body. Leads are necessarily very small in diameter and must still be very flexible, which unavoidably reduces their potential performance or longevity. Leads may fail to function for a variety of causes, including, but not limited to: Medical complications, body rejection phenomena, allergic reaction, fibrotic tissue, or failure of leads by breakage or by breach of their insulation covering. In addition, despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, leads may be easily damaged before, during, or after insertion by improper handling or other intervening acts. Consequently, no representation or warranty is made that failure or cessation of function of leads will not occur; or that the body will not react adversely to the implantation of leads; or that medical complications (including perforation of the heart) will not follow implantation of leads; or that the lead will, in all cases, restore adequate cardiac function. Warranty period of storage: 2 years.																		
27. The list of national and international regulatory documents/standards, which matches medical device	<table> <tr> <th data-bbox="537 965 828 997">Standard</th><th data-bbox="828 965 2184 997">Name of the standard</th></tr> <tr> <td data-bbox="537 997 828 1061">EN 556-1</td><td data-bbox="828 997 2184 1061">Sterilization of medical devices – Requirements for medical devices to be labeled «Sterile» - Part 1: Requirements for terminally sterilized medical devices</td></tr> <tr> <td data-bbox="537 1061 828 1093">EN 980</td><td data-bbox="828 1061 2184 1093">Graphical symbols for use in the labelling of medical devices</td></tr> <tr> <td data-bbox="537 1093 828 1125">EN 1041</td><td data-bbox="828 1093 2184 1125">Information supplied by the manufacturer with medical devices</td></tr> <tr> <td data-bbox="537 1125 828 1204">EN 45502-1</td><td data-bbox="828 1125 2184 1204">Active implantable medical devices – Part 1: General requirements for safety, marking and for information to be provided by the manufacturer</td></tr> <tr> <td data-bbox="537 1204 828 1268">EN 45502-2-2</td><td data-bbox="828 1204 2184 1268">Active implantable medical devices- Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)</td></tr> <tr> <td data-bbox="537 1268 828 1300">EN 62366</td><td data-bbox="828 1268 2184 1300">Medical devices – Application of usability engineering to medical devices</td></tr> <tr> <td data-bbox="537 1300 828 1332">EN ISO 14971</td><td data-bbox="828 1300 2184 1332">Medical devices - Application of risk management to medical devices</td></tr> <tr> <td data-bbox="537 1332 828 1412">EN ISO 11607-1</td><td data-bbox="828 1332 2184 1412">Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems</td></tr> </table>	Standard	Name of the standard	EN 556-1	Sterilization of medical devices – Requirements for medical devices to be labeled «Sterile» - Part 1: Requirements for terminally sterilized medical devices	EN 980	Graphical symbols for use in the labelling of medical devices	EN 1041	Information supplied by the manufacturer with medical devices	EN 45502-1	Active implantable medical devices – Part 1: General requirements for safety, marking and for information to be provided by the manufacturer	EN 45502-2-2	Active implantable medical devices- Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	EN 62366	Medical devices – Application of usability engineering to medical devices	EN ISO 14971	Medical devices - Application of risk management to medical devices	EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
Standard	Name of the standard																		
EN 556-1	Sterilization of medical devices – Requirements for medical devices to be labeled «Sterile» - Part 1: Requirements for terminally sterilized medical devices																		
EN 980	Graphical symbols for use in the labelling of medical devices																		
EN 1041	Information supplied by the manufacturer with medical devices																		
EN 45502-1	Active implantable medical devices – Part 1: General requirements for safety, marking and for information to be provided by the manufacturer																		
EN 45502-2-2	Active implantable medical devices- Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)																		
EN 62366	Medical devices – Application of usability engineering to medical devices																		
EN ISO 14971	Medical devices - Application of risk management to medical devices																		
EN ISO 11607-1	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems																		


	EN ISO 11135-1	Sterilization of health care products - Ethylene oxide – Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices
	EN ISO 10993-1	Biological evaluation of medical devices. Part 1. Evaluation and testing
	ISO 10993-3	Biological evaluation of medical devices. Part 3. Tests for genotoxicity, carcinogenicity and reproductive toxicity
	ISO 10993-4	Biological evaluation of medical devices. Part 4. Selection of tests for interactions with blood
	ISO 10993-5	Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity
	ISO 10993-6	Biological evaluation of medical devices. Part 6. Tests for local effects after implantation
	ISO 10993-7	Biological evaluation of medical devices. Part 7. Ethylene oxide sterilization residuals
	ISO 10993-10	Biological evaluation of medical devices. Part 10. Tests for irritation and delayed-type hypersensitivity.
	ISO 10993-11	Biological evaluation of medical devices. Part 11. Tests for systemic toxicity.
	ISO/TS 10974	Assessment of the safety of magnetic resonance imaging for patients with an active implantable medical device
	ISO 27186	Active implantable medical devices - Four pole connector system for implantable cardiac rhythm management devices - Dimensional and test Requirements
28.Customer Service	For more information, contact your local Medtronic representative, or contact Medtronic Inc. as listed below: Medtronic Inc., USA Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative in Russia: Medtronic, LLC, 123317, Moscow, Presnenskaya Naberezhnaya., 10 Tel.: (495) 580-73-77, Fax: (495) 580-73-78. www.medtronic.ru	

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Sprint Quattro Secure® Model 6947M

Sprint Quattro Secure Model 6947M
SVC Defibrillation Coil

RV Defibrillation Coil

Ring Electrode

Helix Electrode

Anchoring Sleeve

DF4 Connector

SVC Defibrillation Coil Connector

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