B. BRAUN
SHARING EXPERTISE

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Sylvie PLUMET
Catherine BOISMENU
(seen exclusively to certify the above signature)

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B. Braun Medical

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Etablissement de CHASSENEUIL-DU-POITOU

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86360 CHASSENEUIL-DU-POITOU

TO THE
**FEDERAL SERVICE ON
SURVEILLANCE IN HEALTHCARE
(ROSZDRAVNADZOR)**
Department of Registration of Medical Devices

Chasseneuil, le

10 MAI 2019

Covering letter

visé 0286

With this letter we

B. BRAUN MEDICAL, 26 rue Armengaud, 92210 Saint-Cloud, France, as the responsible manufacturer of the medical device

Peripherally inserted central catheters Celsite® PICC-Cel and accessories

Herewith declare that we provided the IFU with Annex 1 with updated information about legal address of manufacturer of the medical device in English language:

1. IFU for *Peripherally inserted central catheters Celsite® PICC-Cel and accessories*

The Annex 1 to the IFU contains next items*:

- 1.1. Name of the medical device;
- 1.2. Name of manufacturer and manufacturing site;
- 1.3. Operating conditions of the catheter Celsite® PICC-Cel in the internal environment of the human body;
- 1.4. Conditions of transportation and storage;
- 1.5. Delivery set and packaging data;
- 1.6. The information about the symbol the conformity Declaration of conformity in the Russian Federation;
- 1.7. The information about the symbol of non-pyrogenic;
- 1.8. Shelf life for the medical device;
- 1.9. Warranty/Complaints;
- 1.10. Procedure for disposal or destruction of a medical device;
- 1.11. Indications for use/ Intended use;
- 1.12. System contents.

*- the description of data, see at Annex 1 to the present IFU.

for the registration purposes of *Peripherally inserted central catheters Celsite® PICC-Cel and accessories* and submission of the listed above document in the Federal Service on Surveillance in Healthcare (Roszdravnadzor) to the Department of Medical Devices Registration.

For and on behalf of

B. BRAUN MEDICAL

Sylvie PLUMET
CoE Chasseneuil, Director

B.Braun Medical

Catherine BOISMENU
CoE Chasseneuil, Quality & Environmental
Management System and Regulatory Affairs, Manager
B.Braun Medical

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
Tél. +33(0)1 41 10 53 00
Fax +33(0)1 41 10 53 99

Instruction for use

Peripherally inserted central catheters Celsite® PICC-Cel and accessories

Manufacturer: B. BRAUN MEDICAL, 26 rue Armengaud, 92210 Saint-Cloud, France.

Manufacturing site: B. BRAUN MEDICAL, 30 avenue des Temps Modernes, 86360 CHASSENEUIL-DU-POITOU, France

2019

System contents

1F = 0.33 mm

1 inch = 25.4 mm

21G/0.8 mm

 PICC single lumen	 PICC double lumen	 ECG adaptor
 Thin wall vein puncture needle 0.8 mm x 0.55 mm x 70 mm	 Straight guidewire 0.018 in x 450 mm	 Straight guidewire 0.018 in x 1300 mm
 Caresite luer access device	 Safety scalpel E11	 Peelable sheath
 Omnifix syringe 10 mL Luer slip	 Measuring tape 1300 mm (+/- 5 mm)	 GRIP-LOK CS securement device

Refer to unit package label for actual content list

I- Description of the medical device

Celsite® PICC-Cel are peripherally inserted central venous catheters (PICC) made of polyurethane (PUR) with a reverse taper design. The catheters are radio-opaque and are marked every cm.

The sideport adaptor is used for priming the catheter with saline prior and during insertion and then removed after the catheter is inserted.

II- Indications

Catheter Celsite® PICC-Cel may be used when a central line is needed for repeated short to long-term intravenous therapy.

Catheter Celsite® PICC-Cel are indicated for intravenous therapies which include drugs such as chemotherapy, antibiotics, antivirals, blood sampling, transfusion, or nutrition.

They may also be used for Contrast Enhanced Computed Tomography (CECT) using high pressure injection (see "XII Recommendations for High Flow rate/High Pressure Injection) and for Central Venous Pressure (CVP) monitoring.

The catheter Celsite® PICC-Cel may be used for up to 3 months. Time period of catheter stay is agreed with treating physician.

The auxiliary products included in the kit are designed to facilitate the installation of the catheter and to ensure the functioning of the catheter during its use.

III - Contraindications

- Known or suspected infection of the puncture site, bacteraemia or septicaemia.
- Known allergy to the materials contained in the kit (main raw materials: Polyurethane, Carbon steel, Stainless steel, Polyvinyl Chloride, Styrene butadiene polymers, Silicone, Polycarbonate, Polypropylene, Isoprene, Polytetrafluoroethylene, Polyether block amide, Polyethylene, Nitinol).
- If the patient's anatomy does not allow the insertion or navigation of the catheter to the chosen site (body size, e.g. obese) or if the patient has had previous radiotherapy to the insertion site.
- Previous venous thrombosis or vascular surgical procedures at the prospective placement site.
- Severe coagulopathies.
- Implantation in agitated patients (e.g. mental disorder).
- Dermatitis, cellulitis or other local tissue factors e.g. over burn sites, which could impact on implantation or maintenance.

IV- Potential complications

Before implantation of the Celsite® PICC-Cel these and other well documented potential complications and their emergency treatment should be considered.

• Accidental removal • Brachial plexus injury • Bleeding • Cardiac arrhythmia • Cardiac tamponade • Catheter rupture, migration or fragmentation • Catheter occlusion • Catheter malposition /retraction • Catheter/air/thrombo-embolism • Catheter erosion through the skin	• Dysrhythmias • Endocarditis • Exit Site infection/necrosis • Fibrin sheath formation • Haematoma • Intolerance/reaction to implant device • Inadvertent arterial puncture • Laceration/Perforation of vessels or viscus • Leakage • Local or general infection/ sepsis/inflammation	• Mediastinal injury • Myocardial erosion • Nerve injury or irritation • Phlebitis • Risks associated with local or general anaesthesia • Risks associated with iodine based contrast media • Thrombophlebitis • Venous thrombosis • Ventricular thrombosis • Vessel erosion
---	--	---

V-MR Information

The catheter Celsite® PICC-Cel is MR safe.

Warning: Verify that the side port and stylet are removed before the MR examination.

VI- Warnings

Please note that the following warnings contain critical information for the patient safety. Failure to follow this Instruction For Use could result in serious injury to the patient:

- The product is sterilised using Ethylene Oxide and is sterile and non-pyrogenic in unopened, undamaged, individual packaging. For single use only.
 - The device and its accessories are not reusable nor designed to be reused. Any re-use would definitely compromise the performance and safety of the device. Components may not be cleaned or re-sterilized.
 - Destroy system components after implantation and explantation.
 - Check that the sterile packaging is undamaged before opening. Confirm the date of sterile expiration on the labeling. Never implant a device which has exceeded its sterile expiration date.
- To maintain sterility, all normal aseptic room sterile procedures should be followed.
- Always verify the integrity of the device before placement.
 - Do not open the tray until ready for use and ensure that all handling, maintenance, use and removal is under strict aseptic conditions.
 - Read instructions for use carefully before using this device. This product should only be implanted, accessed or removed by personnel experienced in the technical and clinical aspects of PICC lines.
 - Special attention should be paid to insertion of any intravenous device in patients who have a pacemaker or other intracardiac device inserted. Insertion of other devices in patients who have a PICC inserted should be carefully considered.
 - Assess the patient for heparin sensitivity. Heparin-induced Thrombocytopenia has been reported after the use of heparin flush solution.
 - Extreme care must be exercised during approaches through tortuous or twisted anatomies which can cause catheter kinking and make catheter insertion difficult or impossible. Always stop pushing on the catheter if significant resistance is encountered.
 - To prevent shearing of the guidewire never withdraw the guidewire through the vein puncture needle. Remove the guidewire and the dilator together. Do not remove the guidewire through the dilator as this may result in unravelling of the guidewire tip.
 - Do not insert or withdraw the guidewire forcibly from any component.
- The guidewire tip could break or unravel.

- When trimming, the catheter should be cut at right angles to reduce the risk of vessel damage.
- Never attempt to cut the stylet.
- The catheter must be secured in place to minimise the risk of catheter breakage and subsequent tip movement and embolisation.
- Never use incompatible drugs together or in sequence: rinse the catheter lumen(s) of all residues with at least 10 mL of sodium chloride (NaCl 0.9%) between each medication and at the end of treatment.
- Avoid repeated and prolonged exposure to alcohol, antiseptics or drugs containing alcohol, acetone and ointments containing polyethylene glycol to avoid any catheter deterioration. Check the ingredients of all preparations used for disinfection or instillation.
- Ensure that all caps and connections prior to, and between manipulations, are correctly secured.
- Clamping the catheter repeatedly in the same location could damage the extension tubing.
- Patient movement can cause catheter tip movement.
- After insertion the patient should be informed about activities or other restrictions which should be observed during the implantation period e.g. avoid excessive arm movements.

VII- Special attention

VII- 1 The Patient

- Check the condition of the skin over the puncture/insertion site.
- The PICC line diameter should be adapted to the diameter of the vein at the insertion site to minimise the risk of thrombosis. It is recommended that the catheter diameter should not exceed 1/3rd of the vessel diameter.

VII- 2 Identification of the Insertion Vein

- Patient diagnosis, age and size of patient, unusual anatomical variables and anticipated catheter dwell time should be taken into account prior to placement.
- Use direct visualisation examinations to assess vein integrity, blood aspirate colour is not always a reliable indicator of venous access. The use of ultra-sound guidance to identify a suitable vein, and during puncture of the vein is strongly recommended.
- PICCs are typically inserted into basilic, cephalic or medial/lateral brachial veins. The basilic vein above the antecubital fossa is the preferred site.

VII- 3 Catheter and Catheter Position

- The catheter tip should always be placed in the Superior Vena Cava at the entrance to the right atrium.
- Ensure that there is no kinking or looping of the catheter along the vein.
- A radiograph should be taken to confirm catheter tip position. Tip position should be monitored routinely according to local protocol.
- Always verify that the catheter lumen(s) is (are) functional by aspirating 2 mL of blood into a syringe and injecting 5 mL of sodium chloride (NaCl 0.9%) before attempting to start an infusion.
- Stop treatment immediately if any pain or swelling is noted or if blood return is absent.
- Do not use sharp instruments near the catheter, or scissors to remove dressing. Properly dispose of sharps in sharps container in accordance with your local protocols.
- Use only the in-line clamps provided.

VII- 4 Occlusion

- If resistance is encountered by aspiration or injection, the lumen of the catheter may be partially or completely occluded. Attempt to inject a few mL of sodium chloride (NaCl 0.9%) into the catheter. If resistance to injection continues or if swelling occurs along the catheter, device malfunction should be suspected.
- In case of obstruction, never try to clear the blockage using a small syringe or fluid under high pressure, this carries the risk of catheter fracture and migration.

VII- 5 Infection

In case of infection, if appropriate, treat with antibiotic drugs according to local protocols. If this fails, or is not appropriate, the catheter should be removed.

VII-6 Safety devices

Celsite® PICC-Cel contain a Safety vein puncture needle and a Safety scalpel designed to reduce the risks of accidental sharps injury. Care must be taken to avoid sharps injuries in accordance with local Infection Control standards and avoid the risk of exposure to contaminated blood.

Handle and dispose of all needles and scalpels according to local protocols.

Use of the Safety vein puncture needle

Access the vein with the needle in the normal manner. After withdrawal of the Safety vein puncture needle activate the safety clip. To do this, grasp the clear clip housing, which houses the metal safety clip, with your free hand.

Advance the clear clip housing forward along the Safety vein puncture needle towards the needle bevel in one continuous motion. Advance the clip housing over the needle bevel activating the safety mechanism. The clip will remain on the needle bevel. The clip housing will separate from the clip and needle. Dispose of the needle and the clear clip housing immediately into sharps container.

Do not attempt to re-insert the needle into the patient or to re-cap the needle with the protective guard.

Use of the Safety Scalpel

To activate the Safety scalpel press down on the green button and slide the blade forwards until the stop point. The blade should now be held in place. To activate the safety mechanism, press down again on the green button, retract the blade, resistance will be felt towards the end of travel. To lock the blade, continue to maintain pressure on the green button and retract the blade against the resistance. The blade will then be locked inside the housing. Safety activation is confirmed by the protrusion of a green point from the housing. The scalpel can then be disposed of in a sharps container.

VIII- Implantation techniques

VIII-1 General considerations

- Ensure that use of the PICC line is indicated taking into consideration the patients needs and proposed use of the PICC line (e.g. patient life-style, duration of use, etc).
- Always ensure that the system is patent before using.

- It is recommended that sodium chloride (NaCl 0.9%), with or without heparin according to local protocols, be used to flush the catheter and stylet before and during implantation
- When trimming the catheter, do not cut the stylet.
- Never use force to remove the stylet. Resistance can damage the catheter. Take care not to damage, tear or perforate the catheter when using the stylet or a scalpel.
- During implantation ensure that the catheter is not damaged by unguarded forceps, suture needles or other sharp instruments.
- Never close the clamp on the catheter stylet.
- It is recommended that the catheter be positioned above the antecubital fossa to reduce the risk of phlebitis and catheter fracture.
- To minimise blood loss and reduce the risk of air embolism during insertion of the guidewire into the needle or of the catheter into the peelable shell, the patient should perform the Valsalva manoeuvre or a syringe or end cap may be used on the hub, or a thumb placed over the hub opening.
- During dual lumen catheter insertion, ensure that the lumen which does not contain the stylet is clamped.
- The length of catheter inserted into the patient should be noted in the patients file.

VIII- 2 References with a short guidewire

- Prepare the placement site according to local aseptic procedures.
- Determine the optimal catheter length using either anatomical landmarks and/or the tape measure included in the kit.
- Unscrew the torque device holding the stylet then pull back the stylet from the catheter, trim the catheter to the desired length, ensuring that the stylet is not cut or damaged. Advance the stylet to the distal end of the catheter and re-secure the torque device.
- Flush the catheter and stylet with either normal sodium chloride (NaCl 0.9%) or heparinised saline.
- Position the arm at 90° to the patient's body; place a tourniquet above the intended insertion site to dilate the veins.
- Select the vein (recommended veins: basilic, cephalic, medial or lateral brachial). The use of ultrasound guidance for vein location is strongly recommended to avoid accidental arterial or nerve puncture or repeated punctures.
- Insert the vein puncture needle, attached to the syringe into the chosen vein; verify the position by ultrasound guidance or by observing blood reflux.
- Remove the syringe.
- Insert the guidewire through the needle 15-20 cm into the vein. Release the tourniquet.
- Fluoroscopic control is recommended when advancing the guidewire past the axilla.
- If necessary, the puncture site can be enlarged using the scalpel.
- Leaving the guidewire in place, withdraw the needle.
- Thread the assembled peelable shell over the guidewire using a twisting motion to pass through the skin planes.
- Remove the guidewire and dilator from the vein.
- Insert the catheter into the peelable shell.
- Slowly advance the catheter to the desired position. Check catheter tip position using fluoroscopy.
- Remove the stylet from the catheter, the catheter can be stabilised by maintaining pressure on the catheter proximal to the peelable shell.
- Remove the peelable shell by pushing down on the 2 handles of the shell, then peel the shell from the catheter ensuring that the catheter is not dislodged from the vessel. The '0' mark on the catheter should be at the insertion site.
- Connect the syringe and aspirate to ensure free blood reflux. A pulsatile flush with sodium chloride (NaCl 0.9%) or heparinised saline should be performed.
- Secure the catheter with the catheter securement device and cover with a transparent dressing. A tape may also be placed around the catheter luer connector for additional security.

- u) Remove the side port and connect the Caresite(s) device to the catheter luer connector.
- v) The final catheter tip position should be checked radiographically and recorded.

VIII- 3 References with a long guidewire

- a) Prepare the placement site according to local aseptic procedures.
- b) Flush the catheter and stylet with either normal sodium chloride (NaCl 0.9%) or heparinised saline.
- c) Position the arm at 90° to the patient's body; place a tourniquet above the intended insertion site to dilate the veins.
- d) Select the vein (recommended veins: basilic, cephalic, medial or lateral brachial). The use of ultrasound guidance for vein location is strongly recommended to avoid accidental arterial or nerve puncture or repeated punctures.
- e) Insert the vein puncture needle, attached to a syringe, into the chosen vein; verify the position by ultrasound guidance or by observing blood reflux.
- f) Remove the syringe.
- g) Insert the guidewire through the needle 15-20 cm into the vein. Release the tourniquet.
- h) Using fluoroscopic control, advance the guidewire to the desired catheter tip location (at the entrance to the right atrium).
- i) Measure the length of guidewire outside the patient (L2) (from the patient skin to the guidewire extremity). Subtract this measurement from the total length of the guidewire (L1) (130 cm) to obtain the correct catheter length (L3).
Example: L1 = 130 cm, L2 = 90 cm, L1 - L2 = L3 = 40 cm, measure 40 cm from the 'O' mark on the catheter to have the correct catheter length.
- j) Unscrew the torque device holding the stylet then pull back the stylet from the catheter, trim the catheter to the desired length, ensuring that the stylet is not cut or damaged. Advance the stylet to the distal end of the catheter and re-secure the torque device.
Leaving the guidewire in place, withdraw the needle.
- k) Thread the assembled peelable shell over the guidewire using a twisting motion to pass through the skin planes.
- l) Remove the guidewire and dilator from the vein.
- m) Insert the catheter into the peelable shell.
- n) Slowly advance the catheter to the desired position. Check catheter tip position using fluoroscopy.
- o) Remove the stylet from the catheter, the catheter can be stabilised by maintaining pressure on the catheter proximal to the peelable shell.
- p) Follow steps r) to v) **References with a short guidewire.**

VIII-4 Over-The-Wire (OTW) Implantation Technique (with long guidewire)

- a) Follow steps a) to m) **References with a short guidewire.**
- b) Then remove the dilator from the shell.
- c) Remove the stylet and side port from the catheter.
- d) Pass the catheter over the guidewire to the desired position at the junction of the Superior Vena Cava and the right atrium. Check catheter tip position using fluoroscopy.
- e) Remove the guidewire from the catheter, the catheter can be stabilised by maintaining pressure on the catheter proximal to the peelable shell.
- f) Follow steps r) to t) **References with a short guidewire.**
- g) Then connect the Caresite (s) device to the catheter luer connector.
- h) The final catheter tip position should be checked radiographically and recorded.

IX- ECG assisted implantation technique

This technique can only be used with the ECG references which have specially adapted accessories.

IX-1 Warnings

- To protect against the risks of current leakage, only use Class 1 ECG monitor CE marked, with floating input (Type CF) and connect additional potential equalization (earth) on the monitor.
- The Celsite® PICC-Cel with ECG technique should not be placed if any of the following is applicable to the patient:
 - The presence of a cardiac pacemaker.
 - Known cardiac abnormalities, or heart rhythm disturbances.
 - Protection against the risks of electrostatic discharge:
 - Use only in rooms protected from electrostatic discharge (conductive flooring).
 - Theatre shoes should be made of a recommended polymer material.
 - Do not allow the guidewire to come into contact with external conductive parts.
 - Before use, check that all connections are secure.
 - Do not autoclave.
 - Do not immerse in liquid.
 - Repairs must be carried out only by the manufacturer.
- If no ECG change is noted during placement, catheter tip position has to be confirmed with fluoroscopy as catheter tip malposition is suspected. A post-operative radiograph should be taken to confirm catheter position.

IX- 2 Implantation Technique with Saline Solutions

- a) Place ECG pads on the patient's chest (see XV-Implantation images Fig.1) and prepare the placement site according to local aseptic procedures.
- b) Determine the optimal catheter length using either anatomical landmarks and/or the tape measure included in the kit.
- c) Remove the stylet and side-port from the catheter and trim the catheter to the desired length.
- d) Follow steps d) to m) **VIII – 2 Implantation Techniques, References with a short guidewire.**
- e) Then attach the luer lock of the ECG adaptor to the catheter hub (see XV-Implantation images Fig.2),
- f) The syringe filled with saline solution (NaCl) 0.9%, is connected to the ECG adapter, ensure there is no air in the system. It is recommended that a minimum of 10 mL of saline is used,
- g) Connect the plug end of the cable to the central socket of the Certodyn® universal adapter (see XV-Implantation images Fig.3) and switch to internal monitoring mode (see XV-Implantation images Fig.4).
- h) Remove the guidewire and dilator from the vein.

Note: If inserting a dual lumen catheter, the extension leg not being used for the saline column should be clamped.

- i) Insert the catheter into the shell.
- j) Advance the catheter towards the heart.
- k) As the catheter approaches the right atrium slowly start to inject the saline solution (NaCl 0.9%) and look for the position where the P wave starts to elevate. Maximal P wave height is reached and maintained when the catheter enters into the right atrium. Note that the amplitude of the P wave will decrease if the catheter is advanced down into the inferior vena cava. After identifying the area where the P wave begins to develop its maximal amplitude (see XV-Implantation images Fig.5) (which corresponds anatomically to the junction between superior vena cava and the right atrium) advance the catheter a further 2 cm. This is the final position of the catheter tip. As the procedure is performed with the patient in the supine position this catheter position allows for the 2-3 cm cranial movement of the catheter which occurs when the patient is upright.

- j) Remove the peelable shell by pushing down on the 2 handles of the shell, then peel the shell from the catheter ensuring that the catheter is not dislodged from the vessel. The '0' mark on the catheter should be at the insertion site.
- m) Aspirate to ensure free blood reflux and flush with sodium chloride (NaCl 0.9%) or heparinised saline.
- n) Secure the catheter with the catheter securement device and cover with a transparent dressing. A tape may also be placed around the catheter luer connector for additional security.
- o) Connect the Caresite(s) device to the catheter luer connector.

Note: A post-operative radiograph should be taken to confirm catheter position if there are any doubts about the correct position of the catheter.

X- Catheter securement device: GRIP-LOK CS

- a) Select the area for the placement of the GRIP-LOK CS when the catheter is already inserted and prepare the skin according to the local hospital protocol for dressing application. Skin must be clean and dry.
- b) Open the top flap and slide the GRIP-LOK CS under the catheter hub centering it in the exposed adhesive area.
- c) Pull the paper backing from one side of the GRIP-LOK CS, then the other, to secure in desired position on the skin.
- d) Remove the interior liner to expose the adhesive area.
- e) Secure top flap over the catheter. Apply firm pressure over the adhesive and around.

Note: use of an alcohol swab will ease removal of securement device from the skin.

Note: Re-use of this device may change its mechanical or biological features and may cause device failure, allergic reactions or bacterial infections.

Note: replace securement device if soiled or saturated in fluid.

Warning: Always use the GRIP-LOK CS system to secure the catheter. Do not suture the catheter directly on the arm to secure the device.

XI- Maintenance of the PICC lines

XI-1 Flushing and patency

- Rinsing of the catheter lumen(s) is essential.
- Rinse the catheter with at least 10 mL of sodium chloride (NaCl 0.9%) before treatment, between administration of different drugs and after each use. If heparinised saline is used, the system should be rinsed first with sodium chloride (NaCl 0.9%) alone, as some drugs may react with heparin and result in blockage of the catheter due to the formation of precipitates.
- Follow local protocols for flushing frequency and heparin concentration, if used.

- Special care should be taken with drugs which have a higher risk of precipitation, with anticoagulant agents, after blood sampling or transfusion to reduce the risk of catheter occlusion.
- Failure to maintain the system may result in occlusion of the catheter.

XI-2 Dressing

- Check under the dressing for accumulation of blood, fluid or moisture in the first 24 hours.
- An occlusive dressing should cover the insertion site at all times.
- Check the dressing regularly for cleanliness, change the dressing every seven days or according to local protocols or if the dressing becomes soiled, wet or non-occlusive.
- Check catheter position during the dressing changes by measuring the external length of the catheter. The external length must remain the same from the insertion until the retrieval.
- Catheter tip position, placement and patency should be checked regularly.
- The Caresite device must be changed at least every 7 days or according to local protocols.

XI-3 Catheter removal

- Remove the dressing.
- Hold the catheter near the insertion site and pull slowly. Do not use excessive force. Do not pull against resistance.
- If resistance is felt, apply a warm compress, wait 20-30 minutes and attempt removal again.
- Verify that the length of catheter removed corresponds to the length inserted (noted in patient's file).
- Place a dressing over the exit site to avoid bleeding after catheter removal.
- Handle and dispose of the removed device in accordance with local Infection Control standards to avoid the risk of exposure to contaminated blood.

XII- Recommendations for high flow rate/high pressure injection

Warnings:

- Always verify that the catheter is functional by aspirating 2 mL of blood into a syringe and injecting a 5 mL of sodium chloride (NaCl 0.9%) into the catheter before attempting to start an injection.
- Do not exceed the recommended pressure and flow rate as system failure may occur.
- Contrast media should be warmed to 37°C (98.6°F) according to drug manufacturers recommendations. Failure to follow this recommendation will result in up to 50% lower flow rates and/or catheter or injection system failure.
- Depending on the technical characteristics of the injection system, the target flow rate might not be attained.
- Flush the catheter with at least 10 mL of sodium chloride (NaCl 0.9%) before and after using the catheter for CECT, followed by usual rinsing procedures.

XIII - Central venous pressure (CVP) monitoring

Central venous pressure (CVP) monitoring should always be used in conjunction with other patient assessment metrics when evaluating cardiac function. Prior to conduct central venous pressure monitoring:

- Ensure proper positioning of the catheter tip
- Flush catheter with sodium chloride (NaCl 0.9%)
- Ensure the pressure transducer is at the level of the right atrium

It is recommended that a continuous infusion of sodium chloride (NaCl 0.9%) (3 mL/hr) is maintained through the catheter while measuring CVP to improve the accuracy of CVP results.

Use your institutions protocols for central venous pressure monitoring procedures.

XIV - TECHNICAL CHARACTERISTICS OF THE Celsite® PICC-Cel

Catheter diameter / N° lumen(s) / Length (mm)	Guidewire length (mm)	Catheter material	Guidewire material	Gravity flow rate NaCL 0.9% (per lumen) (mL/min)	High pressure injection Contrast medium viscosity 11.4 mPa.s (Cp)		Priming volume (ml)
					Max pressure (psi/bar)	Max recommended flow rate	
4F / SL / 510	450	PUR	Stainless steel	≥ 20	300 / 20.7	5	0,6
4F / SL / 510	1300	PUR	Stainless steel	≥ 20	300 / 20.7	5	0,6
4F / SL / 510	450	PUR	Nitinol	≥ 20	300 / 20.7	5	0,6
4F / SL / 510	1300	PUR	Nitinol	≥ 20	300 / 20.7	5	0,6
5F / SL / 610	450	PUR	Stainless steel	≥ 20	300 / 20.7	5	0,7
5F / SL / 610	1300	PUR	Stainless steel	≥ 20	300 / 20.7	5	0,7
5F / SL / 610	450	PUR	Nitinol	≥ 20	300 / 20.7	5	0,7
5F / SL / 610	1300	PUR	Nitinol	≥ 20	300 / 20.7	5	0,7
5F / DL* / 460	450	PUR	Stainless steel	≥ 8	300 / 20.7	5	0,5
5F / DL* / 460	1300	PUR	Stainless steel	≥ 8	300 / 20.7	5	0,5
5F / DL* / 460	450	PUR	Nitinol	≥ 8	300 / 20.7	5	0,5
5F / DL* / 460	1300	PUR	Nitinol	≥ 8	300 / 20.7	5	0,5
5F / DL / 560	450	PUR	Stainless steel	≥ 8	300 / 20.7	4	0,5
5F / DL / 560	1300	PUR	Stainless steel	≥ 8	300 / 20.7	4	0,5
5F / DL / 560	450	PUR	Nitinol	≥ 8	300 / 20.7	4	0,5
5F / DL / 560	1300	PUR	Nitinol	≥ 8	300 / 20.7	4	0,5
6F / DL / 610	450	PUR	Stainless steel	≥ 8	300 / 20.7	5	0,6
6F / DL / 610	1300	PUR	Stainless steel	≥ 8	300 / 20.7	5	0,6
6F / DL / 610	450	PUR	Nitinol	≥ 8	300 / 20.7	5	0,6
6F / DL / 610	1300	PUR	Nitinol	≥ 8	300 / 20.7	5	0,6

SL = single lumen, DL = double lumen

1F = 0.33 mm, 0.018 in = 0.46 mm

Maximum recommended flow rates have been determined according to ISO10555-1:2013, Sterile, single-use intravascular catheters – Part 1: General requirements.

* Length = 460 mm. 5 ml/sec may only be reached if the catheter is trimmed to a length of 460 mm or less. If the catheter length is 560 mm, the maximum flow rate achievable is 4ml/sec.

Catheter diameter / N° lumen(s)	Catheter length (mm)	Peelable sheath size (F)	Peelable sheath ID/length (mm)	Dilator ID/OD/length (mm)
4F/SL	510	4.5	1.5/70	0.4/1.6/86.8
5F/SL	610	5.5	1.8/70	0.4/1.9/86.8
5F/DL	560	5.5	1.8/70	0.4/1.9/86.8
5F/DL	460	5.5	1.8/70	0.4/1.9/86.8
6F/DL	610	6.5	2.1/70	0.4/2.2/86.8

XV - Implantation images

FIG. 1

FIG. 2


FIG. 3


FIG. 4




FIG. 5

XVI – Symbols glossary

 Do not re-use	 Do not resterilize	 Keep dry	
 Power injectable catheter	 Keep away from sunlight	 Latex free	
 Batch code	 Caution, consult accompanying documents	 Do not use if package is damaged	
 Use-by date	 Sterilized using ethylene oxide	 Manufacturer	 Date of manufacture
 Catalogue number	 MR safe	 DEHP-free	

№	Name of additional item	Description
1.1	Name of the medical device	Name of the medical device: "Peripherally inserted central venous catheter Celsite® PICC-Cel in the set"
1.2	Name of manufacturer and manufacturing site	Manufacturer: B. BRAUN MEDICAL, 26 rue Armengaud, 92210 Saint-Cloud, France, Manufacturing site: B. BRAUN MEDICAL, 30 avenue des Temps Modernes, 86360 CHASSENEUIL-DU-POITOU, France.
1.3	Operating conditions of the medical device in the internal environment of the human body	Temperature, (C°): From 32 C° to 42 C°; Humidity, (%): 100 %.
1.4	Conditions of transportation and storage	Finished products are stored and transported at the temperature from -20 ° C to + 50 ° C and relative humidity up to 80%. Individual and group packages of the medical device must be stored in a dry place protected from direct sunlight.
1.5	Delivery set and packaging data	The delivery set for medical device, includes the following: 1. <u>Secondary packaging:</u> Cardboard corrugated box - 1pc. Cardboard corrugated box contains: 2. Instruction for use - 1 pc. 3. Patient Card Celsite® PICC-Cel - 1 pc. 4. Patient Label Celsite® PICC-Cel - 1 × 6 pcs. 5. <u>Sterile Individual packaging:</u> 5.1 Outer packaging: "Breathable" bag - 1 pc. Material of Outer packaging: • Seven-layer nylon transparent film; • High density polyethylene; • Heat-insulating breathable coating: High-density polyethylene (PEHD). 5.2 Inner packaging: Covering Tyvek + Tray - 1 pc. Material of Inner packaging: • Covering Tyvek: High density polyethylene; • Tray: Polyethylene terephthalate. Sterile Inner packaging contains: 6. Medical device: "Peripherally inserted central venous catheter Celsite® PICC-Cel in the set" - 1 pc.
1.6	The information about the symbol the conformity Declaration of conformity in the Russian Federation	The symbol the conformity Declaration of conformity in the RF: 
1.7	The information about the symbol of non-pyrogenic	The symbol of non-pyrogenic : 
1.8	Shelf life for the medical device	Shelf life for the medical device consists: <u>3 years.</u>
1.9	Warranty/Complaints	The company B. BRAUN MEDICAL, France guarantees that the medical device: "Peripherally inserted central venous catheter Celsite® PICC-Cel in the set" meets the requirements of the technical documentation, subject to conditions, transporting and storage.

0287 10 MAI 2019

		Company B. BRAUN MEDICAL, France is not responsible for the use of medical device for other purposes. Claims for civil liability resulting from improper handling or use in combination with other medical devices or accessories will be rejected. On the circulation issues in the territory of the Russian Federation of the medical device: "Peripherally inserted central venous catheter Celsite® PICC-Cel in the set", contact to the authorized organization (importer) in the Russian Federation: Limited Liability Company "B. Braun Medical", Russia, 191040, St. Petersburg, ul. Pushkinskaya, 10 Tel / fax.: 8 (812) 320-40-04, E-mail: office.spb.ru@bbbraun.com.
1.10	Procedure for disposal or destruction of a medical device	Open unused or damaged packaging can not be used, they should be disposed of. Disposal of the medical device after use is carried out in accordance with the requirements of SanPiN 2.1.7.2790-2010: • Unused medical device and unused packaging must be disposed in accordance with the class A. After use, the medical device must be disposed in accordance with the class B. If you have any questions regarding disposal of a medical device, you should contact the authorized organization (importer) in the Russian Federation.
1.11	Indications for use/ Intended use	The Celsite® PICC-Cel catheter is designed to be inserted into the peripheral vein and moved to a large central vein to create short-term or long-term intravascular access. The catheter allows multiple access to the vascular system for up to 30 days or more in order to: • carrying out multiple infusion therapy, including the introduction of such drugs as: antibiotics, chemotherapeutic, antiviral drugs; • taking blood samples; • carrying out transfusion; • introduction of parenteral nutrition; • Conducting Contrast Computed Tomography using high-pressure infusion; • control of central venous pressure (CVP); The auxiliary products included in the kit are designed to facilitate the installation of the catheter and to ensure the functioning of the catheter during its use. <u>Duration of use:</u> The catheter Celsite® PICC-Cel may be used for up to 3 months. Time period of catheter stay is agreed with treating physician. Intended patient population: adults and pediatric (taking into consideration the patients anatomy needs). <u>Potential consumer:</u> This product should only be implanted, accessed or removed by personnel experienced in the technical and clinical aspects of PICC lines.
1.12	System contents	The sets with ECG adaptors are not include in the current registration process (registration dossier) in Russian Federation.

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For and on behalf of
B. BRAUN MEDICAL

Sylvie PLUMET
CoE Chasseneuil, Director
B.Braun Medical

Catherine BOISMENU
CoE Chasseneuil, Quality & Environmental
Management System and Regulatory Affairs, Manager
B.Braun Medical

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