Medtronic

Abre™

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501181-001 Rev. F

Nicole A. Weulander  
Regulatory Affairs Specialist  
Date: 13 AUG 2018

 ООО ЦИПМИ  [INFO@CIRMI.RU](mailto:INFO@CIRMI.RU)  [WWW.GOSZDRAVNADZOR.RU](http://WWW.GOSZDRAVNADZOR.RU)

[illegible]


Figure 1. Abre™ venous self-expanding stent

### Device description

The Abre™ venous self-expanding stent system (Abre™ system) is a venous self-expanding nitinol stent system. The Abre™ stent is intended for permanent implant and comes pre-mounted on a 9 Fr (3.1 mm) delivery system inserted over a 0.89 mm (0.035 in) wire.

- The Abre™ self-expanding stent, as shown in Figure 1, is cut from a nickel titanium alloy (nitinol) tube in an open lattice design (1), with integral nitinol markers at the trailing end (2) and the leading end (3) of the stent. Upon deployment, the Abre stent exerts an outward force to establish patency.
- The Abre™ delivery system, as shown in Figure 1, has a triaxial shaft design, comprised of an inner shaft assembly (4), a retractable sheath (silver, 5), an isolation sheath (blue, 6), and an ergonomic handle (7). The inner shaft assembly terminates with a flexible catheter tip (9) and originates at the luer hub (10). The deployment handle (7) has a removable locking pin (11), thumbwheel (12), and luer hub (10).
- The radiopaque markers, at the trailing end and the leading end of the stent, guide stent positioning at the target lesion before Abre stent deployment. There is a single luer hub on the deployment handle. A thumbwheel on the deployment handle rotates to pull back the retractable sheath (silver). A locking pin prevents the stent from being deployed before use, and must be removed to actuate the thumbwheel. The Abre stent is fully deployed when the radiopaque marker (8) in the retractable sheath (silver) reaches beyond the integral nitinol markers on the trailing end of the stent. A uniquely designed isolation sheath (blue, attached to the deployment handle) improves control and accuracy of stent delivery.

### Indications for use

The Abre™ venous self-expanding stent system is intended for use in the iliofemoral veins for the treatment of symptomatic venous outflow obstruction.

### Contraindications

- Do not use the Abre system with patients in whom anticoagulant or antiplatelet therapy is contraindicated.
- Do not use the Abre system with patients with known hypersensitivity to nickel titanium (nitinol).

### Warnings

- The device is provided sterile, for single use only. Do not reprocess or re-sterilize. Reprocessing or re-sterilizing could increase the risk of patient infection and the risk of compromised device performance.
- If unusually high resistance is encountered when advancing the Abre delivery system over the guidewire, assess the cause of the resistance before proceeding.
- If high resistance is felt when initially rotating the thumbwheel, do not force deployment. Carefully withdraw the system and do not use it.

### Precautions

- Use of the Abre stent in vessel beds other than the iliofemoral vein region has not been studied.
- Abre stent system implant procedures should only be performed by physicians that have experience with interventional techniques in the vascular system.
- Appropriate stent size selection is crucial. Use Table 1 for guidance in selecting the appropriate stent size. Proper size selection reduces stent migration and ensures appropriate stent apposition to the vessel wall.
- Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Do not use

the stent system if it is damaged or compromised.

- Carefully remove the system from the tray and tube without kinking the system.
- Do not use the stent system if the Abre stent is partially deployed before starting the procedure.
- Gaining access too close to the intended treatment site can lead to difficulty with stent deployment.
- Always use an introducer sheath during the implant procedure to protect both the vessel and puncture site.
- Ensure that the locking pin remains in locked position until the target site is reached.
- If the retractable sheath (silver) is exposed beyond the hemostatic valve, friction between the retractable sheath and the hemostatic valve can impair deployment. To avoid friction between the retractable sheath (silver) and the hemostatic valve, pull the introducer sheath back so that the hemostatic valve covers the end of the isolation sheath.
- Do not reposition the Abre stent after establishing apposition against the vessel wall.
- The Abre delivery system is not designed for recapturing the stent.
- The Abre stent is not designed to be lengthened or shortened from its nominal length. Excessive stent lengthening or shortening can increase the risk of stent fracture.
- If resistance is felt when withdrawing the system, do not force withdrawal.
- Use care when crossing a deployed Abre stent with any adjunct device to avoid stent dislodgement or damage to the adjunct device.

### Potential adverse events

- Access failure
- Access site infection
- Allergic reaction to contrast medium or procedure medications
- Allergic reaction to nitinol or other device materials
- Arrhythmia
- AV fistula
- Bleeding
- Bruising
- Death
- Device breakage
- Device maldeployment
- Edema
- Fever
- Hematoma
- Hypotension, nausea, or other vasovagal response
- Infection
- Myocardial infarction
- Pain
- Pseudoaneurysm
- Pulmonary embolism
- Renal insufficiency/renal failure (new or worsening)
- Sepsis
- Stent fracture
- Stent malapposition
- Stent malposition
- Stent migration
- Stroke/paradoxical embolism/transient ischemic attack/intracerebral hemorrhage
- Tissue necrosis
- Transfusion reaction following blood transfusion for treatment

- of major bleeding
- Vessel damage, including perforation or rupture
- Venous occlusion/thrombosis, within or outside of stented segment

#### Stent size selection

Considering the estimated anatomic vessel diameter, use Table 1 to select the Abre™ stent diameter size. Choose a stent length that extends beyond both ends of the target lesion, with at least 1 cm on each side of the lesion to reduce the risk of restenosis.

Table 1. Sizing guide

| Stent diameter (mm) | Estimated anatomic vessel diameter (mm) | Stent length (mm)         |
|---------------------|-----------------------------------------|---------------------------|
| 10                  | 7.5–9.5                                 | 40, 60, 80, 100, 120, 150 |
| 12                  | 9.5–11.5                                | 60, 80, 100, 120, 150     |
| 14                  | 11.5–13.5                               | 60, 80, 100, 120, 150     |
| 16                  | 13.5–15.5                               | 60, 80, 100, 120, 150     |
| 18                  | 15.5–17.5                               | 60, 80, 100, 120, 150     |
| 20                  | 17.5–19.0                               | 60, 80, 100, 120, 150     |

**Caution:** Appropriate stent size selection is crucial. Use Table 1 for guidance in selecting the appropriate stent size. Proper size selection reduces stent migration and ensures appropriate stent apposition to the vessel wall.

#### How supplied

**Warning:** The device is provided sterile, for single use only. Do not reprocess or resterilize. Reprocessing or resterilizing could increase the risk of patient infection and the risk of compromised device performance.

#### Directions for use

##### Required items for implant procedure

- 5–10 cc syringe filled with saline
- 0.89 mm (0.035 in) guidewire
- 9 Fr minimum hemostatic introducer sheath
- Balloon catheter

##### Preparation

1. Open the shelf box and remove the contents. Carefully inspect the pouch, looking for damage to the sterile barrier, then carefully peel open the pouch and extract the device by the handle.  
**Caution:** Carefully inspect the sterile package and device before use to verify that no damage occurred during shipment. Do not use the stent system if it is damaged or compromised.  
**Caution:** Carefully remove the delivery system from the tray and tube without kinking the system.
2. Attach a syringe filled with saline to the luer lock injection hub. Inject the saline through the guidewire lumen until it comes out of the catheter tip.
3. Examine the catheter to ensure that the Abre stent is fully contained within the retractable sheath (silver).  
**Caution:** Do not use the stent system if the Abre stent is partially deployed before starting the procedure.

##### Stent deployment

1. Insert the introducer sheath and the guidewire.
  - a. Gain access at the appropriate site, using a compatible introducer sheath with a hemostatic valve.  
**Caution:** Gaining access too close to the intended treatment site can lead to difficulty with stent deployment.
  - b. Insert a guidewire of appropriate length across the target lesion via the introducer sheath.  
**Caution:** Always use an introducer sheath during the implant procedure to protect both the vessel and puncture site.
2. Predilate the lesion with a balloon catheter using standard techniques. Remove the balloon from the patient while maintaining access with the guidewire.
3. Introduce the Abre™ stent delivery system over the guidewire through the hemostatic valve and introducer sheath.

**Warning:** If unusually high resistance is encountered, when advancing the Abre™ delivery system over the guidewire, assess the cause of the resistance before proceeding.

**Caution:** Ensure that the locking pin remains in locked position until the target site is reached.

4. Advance the Abre delivery system until the leading edge of the stent is beyond the target lesion.
5. Ensure that the isolation sheath (blue) remains inside the introducer sheath valve.  
**Caution:** If the retractable sheath (silver) is exposed beyond the hemostatic valve, friction between the retractable sheath and the hemostatic valve can impair deployment. To avoid friction between the retractable sheath (silver) and the hemostatic valve, pull the introducer sheath back so that the hemostatic valve covers the end of the isolation sheath.
6. Remove the locking pin by holding the device handle with one hand and gently pulling the locking pin with the other hand.
7. While holding the device with one hand, use the other hand to gently support and position the isolation sheath at the hemostatic valve (see Figure 2). Do not constrict the retractable sheath during deployment.

Figure 2. Hand placement for stent deployment

8. Initiate Abre stent deployment by slowly rotating the thumbwheel in the direction indicated by the arrow on the handle.
9. After the leading end of the stent has emerged from the retractable sheath, but before achieving vessel apposition, reposition the Abre stent as needed.

**Warning:** If high resistance is felt when initially rotating the thumbwheel, do not force deployment. Carefully withdraw the system and do not use it.

**Caution:** Do not reposition the Abre stent after establishing apposition against the vessel wall.

**Caution:** The Abre delivery system is not designed for recapturing the stent.

10. Continue to rotate the thumbwheel until the stent is fully deployed.

**Caution:** The Abre stent is not designed to be lengthened or shortened from its nominal length. Excessive stent lengthening or shortening can increase the risk of stent fracture.

**Note:** If placement of sequential stents is necessary, ensure that there is sufficient overlap.

##### Post stent deployment

1. Withdraw the entire delivery system from the patient.  
**Caution:** If resistance is felt when withdrawing the system, do not force withdrawal.
2. Perform post-deployment balloon dilation as needed, using an appropriately sized balloon catheter with conventional dilation techniques.  
**Caution:** Use care when crossing a deployed Abre stent with any adjunct device to avoid dislodging the stent or damaging the adjunct device.
3. Confirm that the Abre stent is fully expanded, then remove the balloon catheter from the patient.  
**Note:** Always visualize the stent by using imaging techniques to verify full vessel wall apposition.
4. Remove the guidewire and introducer sheath from the patient and discard the delivery system, guidewire, and introducer sheath.

##### MRI safety information

###### MR Conditional

Nonclinical testing demonstrated that the Abre stent in single and overlapped conditions is MR Conditional for stents up to 150 mm. A patient with this device can be scanned safely, immediately after stent placement, under the following conditions:

- Static magnetic field of 1.5 Tesla or 3.0 Tesla
- Maximum spatial gradient magnetic field of 4000 Gauss/cm or

less (40 T/m).

- Maximum MR system-reported, whole-body-averaged specific absorption rate (SAR) of 2.0 W/kg (normal operating mode).

MR image quality may be compromised if the area of interest is in the exact location or close to the position of the Abre stent.

#### **MRI-related temperature rise**

Under the scan conditions as defined in the previous section, "MR Conditional," the Abre stent is expected to produce a maximum temperature rise less than or equal to 5.2 °C after 15 minutes of continuous scanning (per pulse sequence). The effect of temperature rise in the MRI environment for stents with fractured struts is not known.


It is recommended that patients register conditions under which the implant can be scanned safely with the MedicAlert Foundation ([www.medicalert.org](http://www.medicalert.org)) or equivalent organization.

#### **Artifact information**

In non-clinical testing, the maximum artifact size as seen on the gradient echo pulse sequence at 3.0 Tesla extends approximately 5 mm, relative to the size and shape of the Abre™ stent. The lumen of the stent can be visualized using the T1-weighted spin echo pulse sequence and the T1-weighted gradient echo pulse sequence at 3.0 Tesla.

#### **Warranty**

This product is subject to standard Medtronic warranty terms. For standard terms and conditions of supply for Australia and New Zealand, go to [www.medtronic.com](http://www.medtronic.com).

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[illegible][illegible]

To avoid this problem, the *dhfr* gene was removed from the vector and replaced with the *dhfr* gene from the *Shigella flexneri* T3000 strain. The resulting vector was designated pT3000-*dhfr* (Fig. 1a). The other three vectors, pT3000-*dhfr* (Fig. 1b), pT3000-*dhfr* (Fig. 1c), and pT3000-*dhfr* (Fig. 1d), were constructed by inserting the *dhfr* gene from the *Shigella flexneri* T3000 strain into the *dhfr* gene from the *Shigella flexneri* T3000 strain. The resulting vectors were designated pT3000-*dhfr* (Fig. 1e), pT3000-*dhfr* (Fig. 1f), and pT3000-*dhfr* (Fig. 1g). The other three vectors, pT3000-*dhfr* (Fig. 1h), pT3000-*dhfr* (Fig. 1i), and pT3000-*dhfr* (Fig. 1j), were constructed by inserting the *dhfr* gene from the *Shigella flexneri* T3000 strain into the *dhfr* gene from the *Shigella flexneri* T3000 strain. The resulting vectors were designated pT3000-*dhfr* (Fig. 1k), pT3000-*dhfr* (Fig. 1l), and pT3000-*dhfr* (Fig. 1m).

do so, start with a new copy of the text and then delete the text that you don't want. In this case, you would delete the text that says "I'm happy." To delete text, click the text that you want to delete, and then press the **Delete** key on the keyboard. To delete a whole paragraph, click the paragraph, and then press the **Delete** key. To delete a whole page, click the page, and then press the **Delete** key. To delete a whole document, click the document, and then press the **Delete** key.

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- To proceed, first let's discuss the question, "What is the *best* method?" Many factors are involved, and you may choose to use a different method depending on the nature of your data. For instance, if you are analyzing time series data, you may want to use a method that accounts for autocorrelation. Or, if you are analyzing data with a large number of variables, you may want to use a method that can handle high-dimensional data.
- For the purposes of this discussion, let's assume that you are analyzing time series data. In this case, there are two main methods: the autocorrelation function (ACF) and the partial autocorrelation function (PACF). The ACF measures the correlation between a time series and its lagged values, while the PACF measures the correlation between a time series and its lagged values, controlling for the values of the time series at all other lags.
- For example, if you are analyzing a time series of monthly sales data, you might use the ACF to determine if there is a seasonal pattern. If the ACF shows a significant peak at a lag of 12 months, this would suggest a seasonal pattern. Alternatively, you might use the PACF to determine if there is a linear trend. If the PACF shows a significant peak at a lag of 1 month, this would suggest a linear trend.
- Of course, there are many other methods available, and the choice of method will depend on the specific characteristics of your data. It is always a good idea to consult with a statistician or data analyst to help you choose the best method for your analysis.

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
I Nicole Weulander, with the title of Regulatory Affairs Specialist, International hereby swear (or affirm) that the attached Appendix to the Instruction for use for Abre is a true and accurate document.

Country: Russia

  
Nicole Weulander

21 FEB 2019


*A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document, to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document*

STATE OF MINNESOTA )  
 )ss.  
COUNTY OF HENNEPIN )

Subscribed and sworn to (or affirmed) before me this 21<sup>st</sup> day of February, 2019 by Nicole Weulander, personally appeared who is personally known to me –or– proved to me through satisfactory evidence to be the person whose name is subscribed to the within instrument, and acknowledged to me that he/she/they executed the same for the purposes therein stated.

WITNESS my hand and official seal

  
Deborah L. Shaffer, Notary Public

Description of Attached Document(s):

Appendix to the Instruction for use for Abre

**Medtronic, Inc.**  
**710 Medtronic Parkway**  
**Minneapolis, MN 55432**  
**www.medtronic.com**

**Appendix to the Instruction for use**

**Venous Self-expanding Stent System Abre, model:**

AB9G10040090, AB9G10060090, AB9G10080090, AB9G10100090,  
AB9G10120090, AB9G10150090, AB9G12060090, AB9G12080090,  
AB9G12100090, AB9G12120090, AB9G12150090, AB9G14060090,  
AB9G14080090, AB9G14100090, AB9G14120090, AB9G14150090,  
AB9G16060090, AB9G16080090, AB9G16100090, AB9G16120090,  
AB9G16150090, AB9G18060090, AB9G18080090, AB9G18100090,  
AB9G18120090, AB9G18150090, AB9G20060090, AB9G20080090,  
AB9G20100090, AB9G20120090, AB9G20150090

- 1.1. With regards to the manufacturer of the medical device - full name and, if available, abbreviated name, including company name, legal form of the legal entity, location address or last name, first name and (if available) patronymic, details of the identity document, address of residence of the individual entrepreneur, as well as phone numbers and (if available) e-mail address of the legal entity or individual entrepreneur;**

**Legal Manufacturer**

Medtronic, Inc.  
710 Medtronic Parkway, Minneapolis, MN 55432,  
USA  
[www.medtronic.com](http://www.medtronic.com),  
Tel: 763-514-4000

**Addresses of manufacturing facility of medical device:**

Medtronic, Inc., 4600 Nathan Ln N, Plymouth, MN  
55442 USA

**Design Facility**


Medtronic, Inc., 2300 Berkshire Lane N Suite 5,  
Plymouth, MN 55441 USA

- 1.2. With respect to authorized representative of the manufacturer of the medical device - full and, if available, abbreviated name, including company name, legal form of the legal entity, address (location), telephone numbers and (if available) e-mail address of the legal entity**

**Authorized Representative of Manufacturer:**

«Medtronic» LLC  
123317, Moscow, Presnenskaya Naberezhnaya., 10  
**Tel.:** (495) 580-73-77, **Fax:** (495) 580-73-78, [www.medtronic.ru](http://www.medtronic.ru)

  
Nicole A. Weulander  
Sr. Regulatory Affairs Specialist  
Date: 21 FEB 2019

Available upon request

«Medtronic» LLC

123317, Moscow, Presnenskaya Naberezhnaya., 10

Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru

**Holder of Registration Certificate in Russia**

“Medtronic” LLC, 123317, Russia, Moscow, Presnenskaya Naberezhnaya, 10

Tel: (495) 580-73-77, Fax: (495) 580-73-78

www.medtronic.ru

**2. Technical characteristics of the medical device**

The Abre system consists of a stent and stent delivery system designed specifically for implantation in the peripheral venous system. The proposed system consists of a flexible self-expanding stent made of a nickel-titanium alloy (nitinol) provided in multiple lengths and diameters. All stent sizes are available on an over-the-wire stent delivery system consisting of a 9F catheter and a working length of 90cm.

Table 1 below provides the stent sizes of the Abre system.

Table 1. Abre Stent Sizes

|                     |    | Stent Length (mm) |    |    |     |     |     | Nitinol Tube Wall Thickness |
|---------------------|----|-------------------|----|----|-----|-----|-----|-----------------------------|
|                     |    | 40                | 60 | 80 | 100 | 120 | 150 |                             |
| Stent Diameter (mm) | 10 | x                 | x  | x  | x   | x   | x   | 0.018"                      |
|                     | 12 |                   | x  | x  | x   | x   | x   |                             |
|                     | 14 |                   | x  | x  | x   | x   | x   |                             |
|                     | 16 |                   | x  | x  | x   | x   | x   | 0.028"                      |
|                     | 18 |                   | x  | x  | x   | x   | x   |                             |
|                     | 20 |                   | x  | x  | x   | x   | x   |                             |

0.018" = 0.4572 mm

0.028" = 0.7112 mm

The Abre venous self-expanding stent system (Abre system) is a vascular self-expanding nitinol stent system. The Abre stent is intended for permanent implant and comes premounted on a 9 Fr (3.1 mm) delivery system inserted over a 0.89 mm (0.035 in) wire.

- The Abre self-expanding stent, as shown in Figure 1, is cut from a nickel titanium alloy (nitinol) tube in an open lattice design (1), with integral nitinol markers at the trailing end (2) and the leading end (3) of the stent. Upon deployment, the Abre stent exerts an outward force to establish patency.

- The Abre delivery system, as shown in Figure 1, has a triaxial shaft design, comprised of an inner shaft assembly (4), a retractable sheath (silver, 5), an isolation sheath (blue, 6), and an ergonomic handle (7). The inner shaft assembly terminates with a flexible catheter tip (9) and originates at the luer hub (10).

The deployment handle (7) has a removable locking pin (11), thumbwheel (12), and luer hub (10).

- The radiopaque markers, at the trailing end and the leading end of the stent, guide stent positioning at the target lesion before Abre stent deployment. There is a single luer hub on the deployment handle. A thumbwheel on the deployment handle rotates to pull back the retractable sheath (silver). A locking pin prevents the stent from being deployed before use, and must be removed to actuate the thumbwheel. The Abre stent is fully deployed when the radiopaque marker (8) in the retractable sheath (silver) reaches beyond the integral nitinol markers on the trailing end of the stent. A uniquely designed isolation sheath (blue, attached to the deployment handle) improves control and accuracy of stent delivery.

Available upon request

Appendix to the Instruction for use

DxTerity Diagnostic Angiographic Catheters



Figure 1a. Abre venous self-expanding stent

Tolerances (if not mentioned on picture)

.x  $\pm 0.2$  mm

.xx  $\pm 0.13$  mm

.xxx  $\pm 0.013$  mm

Figure 1b. retractable sheath location with dimensions

Available upon request

Appendix to the Instruction for use

 DxTerity Diagnostic Angiographic Catheters

| STENT Ø<br>(mm) | STENT<br>LENGTH<br>(mm) | A cm | B cm<br>REF |
|-----------------|-------------------------|------|-------------|
| 10              | 40                      | 12.2 | 27.2        |
|                 | 60                      | 12.2 | 27.2        |
| 12              | 40                      | 12.2 | 27.2        |
|                 | 60                      | 12.2 | 27.2        |
|                 | 100                     | 12.2 | 27.2        |
| 14              | 60                      | 12.2 | 27.2        |
|                 | 100                     | 12.2 | 27.2        |
| 16              | 60                      | 12.2 | 27.2        |
|                 | 100                     | 12.2 | 27.2        |
|                 | 150                     | 17.0 | 37.0        |
| 18              | 60                      | 12.2 | 27.2        |
|                 | 100                     | 12.2 | 27.2        |
|                 | 150                     | 17.0 | 37.0        |
| 20              | 60                      | 12.2 | 27.2        |
|                 | 100                     | 12.2 | 27.2        |
| 10              | 100                     | 12.2 | 27.2        |
|                 | 150                     | 17.0 | 37.0        |
| 12              | 150                     | 17.0 | 37.0        |
| 14              | 40                      | 12.2 | 27.2        |
|                 | 150                     | 17.0 | 37.0        |
| 20              | 150                     | 17.0 | 37.0        |

|    |     |      |      |
|----|-----|------|------|
| 10 | 80  | 12.2 | 27.2 |
|    | 120 | 17.0 | 37.0 |
| 12 | 80  | 12.2 | 27.2 |
|    | 120 | 17.0 | 37.0 |
| 14 | 80  | 12.2 | 27.2 |
|    | 120 | 17.0 | 37.0 |
| 16 | 80  | 12.2 | 27.2 |
|    | 120 | 17.0 | 37.0 |
| 18 | 80  | 12.2 | 27.2 |
|    | 120 | 17.0 | 37.0 |
| 20 | 80  | 12.2 | 27.2 |
|    | 120 | 17.0 | 37.0 |

Figure 1c. Table with parameter ratios regarding Figure 1b.



Figure 1d. Abre Left Handle Assembly, location of logo "Medtronic"

Available upon request


DxTerity Diagnostic Angiographic Catheters

.xxx  $\pm 0.013$  mm


| A cm | ITEM NO. 1 |
|------|------------|
| 56.0 | 313032-001 |
| 66.0 | 313032-002 |

A- overall length of the isolation sheath (the differences between items is that there are two isolation sheaths depending size of the stent)

**Figure 1e.** Isolation sheath drawing with dimensions and tolerances

**Figure 1f.** Abre Stent Radiopaque markers location,

.x  $\pm 0.2$  mm  
 .xx  $\pm 0.13$  mm  
 .xxx  $\pm 0.013$  mm

**Figure 1g.** Drawing of Removable locking pin with main parameters and tolerances.

Stent parameter:

1. Stent inner diameter (ID) after deployment must be:

Available upon request

Appendix to the Instruction for use

 DxTerity Diagnostic Angiographic Catheters

| Ref | Inch            | mm             |
|-----|-----------------|----------------|
| 10  | 0.394<br>±0.016 | 10.00<br>±0.40 |
| 12  | 0.472<br>±0.019 | 12.00<br>±0.48 |
| 14  | 0.551<br>±0.022 | 14.00<br>±0.56 |
| 16  | 0.630<br>±0.025 | 16.00<br>±0.64 |
| 18  | 0.709<br>±0.028 | 18.00<br>±0.72 |
| 20  | 0.787<br>±0.031 | 20.00<br>±0.80 |

2. The length of the stent in catheter must be:

| Ref | Inch            | mm             |
|-----|-----------------|----------------|
| 40  | 1.575<br>±0.157 | 40.0<br>±4.0   |
| 60  | 2.362<br>±0.236 | 60.0<br>±6.0   |
| 80  | 3.150<br>±0.315 | 80.0<br>±8.0   |
| 100 | 3.937<br>±0.394 | 100.0<br>±10.0 |
| 120 | 4.724<br>±0.472 | 120.0<br>±12.0 |
| 150 | 5.906<br>±0.591 | 150.0<br>±15.0 |

**3. Information on the list of the main characteristics for the operation (application) of the medical device, transportation and storage conditions (for example, temperature and humidity, lighting and other characteristics)**

*Operation condition:*

Do not use if package is opened or damaged. Use at controlled room temperature.

The product is intended for use by trained medical personnel in the operating room: at room temperature (18 ° C) not higher than 25 ° C, relative humidity of 75%, a pressure of from 86.0 to 106.0 kPa. During the procedure the product is inside of the human body and is operated at a body temperature from 32 to 42 ° C

*Storage condition:*

Catheter should be stored in cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, or ultraviolet light. Rotate inventory so that the device is used prior to the «Used by» date on the package label

Temperature: from 10 °C till +40 °C

Relative Humidity: from 30% till 75%

Pressure: from 585 till 1060 kPa

*Transportation condition:*

Care should be taken when transporting the devices. Packed by the manufacturer devices are transported by all modes of transport in closed vehicles in accordance with the rules of carriage of goods, operating in each mode of transport.

Temperature: from -18 °C till +55 °C

Relative Humidity: from 30% till 85%

Available upon request

Pressure: from 585 till 1060 kPa

**4. List of national standards applied by the manufacturer**

Venous Self-expanding Stent System Abre is compliant with active national standards mentioned in the Declaration of Conformity.

**5. Procedure and conditions for the disposal or destruction of a medical device.**

After use, this product may be a potential biohazard (hazard class B). Handle and dispose of all such device is accordance with accepted medical practice and applicable local, regional and national laws and regulations.

**6. List and description of medical device materials**

| Device          | Material                              | Component                                                                                      |                     |  |
|-----------------|---------------------------------------|------------------------------------------------------------------------------------------------|---------------------|--|
| Sub-System      |                                       |                                                                                                |                     |  |
| Stent           | Nitinol                               | Stent                                                                                          |                     |  |
|                 | Silicone lubricant                    | Lubricant                                                                                      |                     |  |
| Delivery System | Polytetrafluoro ethylene              | PTFE liner                                                                                     |                     |  |
|                 | Nylon                                 | Nylon tube (proximal)                                                                          |                     |  |
|                 | Nylon                                 | Nylon tube (distal)                                                                            |                     |  |
|                 | HDPE                                  | Isolation sheath                                                                               |                     |  |
|                 | Platinum-irridium<br>Pt 10%<br>Ir 90% | Radiopaque marker                                                                              |                     |  |
|                 | PEBAX                                 | Flexible catheter tip                                                                          |                     |  |
|                 | High Density Polyethylene             | Stop tube                                                                                      | Inner shaft         |  |
|                 |                                       | Tip Tube                                                                                       |                     |  |
|                 | High Density Polyethylene             | Retractable sheath                                                                             |                     |  |
|                 | RTP                                   | Handle Left & Right, Thumbwheel with colorant (colorant mark or colorant manufacture) Overmold |                     |  |
|                 | Ink                                   | Logo “Medtronic”<br>Left Handle Ink                                                            |                     |  |
|                 | Pebax                                 | Tip                                                                                            |                     |  |
|                 | PC/ABS                                | Removing locking tip                                                                           |                     |  |
|                 | Polycarbonate                         | Proximal luer                                                                                  |                     |  |
|                 |                                       | thermoplastic elastomer                                                                        | Strain Relief       |  |
|                 |                                       | Ink                                                                                            | Strain Relief (Ink) |  |

**7. Labelling and Packaging of medical device**

The Abre packaging system consists of a thermoformed tray, protective tube, sterile barrier pouch, and folding carton. The delivery system is protected by the tray system.

The tray system is comprised of an ABS tube, into which all catheter lengths are loaded into, and is then snapped into the tray. The handle is protected by snapping it into a pocket in the tray. A breathable Nylon/Tyvek pouch forms the sterile barrier for the packaging system. An IFU and patient card are included within the carton. The sealed pouch assembly is then loaded into a folded carton for shipment and storage.

|                            |                             |
|----------------------------|-----------------------------|
| Primary packages dimension |                             |
| Abre L x W $\pm 5\%$       | 6.13 cm (W) x 56 cm (L)     |
| Primary packages weight    |                             |
| Abre g $\pm 5\%$           | 164.3 grams                 |
| Second packages dimension  |                             |
| $\pm 5\%$                  | 5.56 (W) x 1.25 (H) x 55(L) |

Available upon request

Appendix to the Instruction for use

 DxTerity Diagnostic Angiographic Catheters

|                                                                                                                                                                      |                                                                                                                                                          |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| Second packages weight                                                                                                                                               |                                                                                                                                                          |
| ± 5%                                                                                                                                                                 | 0.7 kg                                                                                                                                                   |
| Transport packages dimension                                                                                                                                         |                                                                                                                                                          |
| Abre, L x W x T, ± 5%                                                                                                                                                | 22 x 18 x 57 ½                                                                                                                                           |
| Transport packages weight                                                                                                                                            |                                                                                                                                                          |
| Abre, kg± 5%                                                                                                                                                         | 33.6 kg + 3.2 kg (the weight of shipper) = 36.8 kg                                                                                                       |
| <b>8. Information about insertion, relocation, withdraw of Abre Stent system</b>                                                                                     |                                                                                                                                                          |
| Information about insertion, relocation, withdraw of Abre Stent system:                                                                                              |                                                                                                                                                          |
| If high resistance is felt when initially rotating the thumbwheel, do not force deployment. Carefully withdraw the system and do not use it.                         |                                                                                                                                                          |
| Do not reposition the Abre stent after establishing apposition against the vessel wall. Caution: The Abre delivery system is not designed for recapturing the stent. |                                                                                                                                                          |
| Track,<br>Deploy, and<br>Withdrawal                                                                                                                                  | All units must be able to be inserted and tracked over a guidewire in a sheath in a representative path without pre-deployment.                          |
|                                                                                                                                                                      | All units must successfully deploy while placed in a representative anatomical model                                                                     |
|                                                                                                                                                                      | Post deployment, all units must be able to be removed from sheath over a guidewire in a representative path with delivery system remaining fully intact. |
| <b>9. Shelf life</b>                                                                                                                                                 |                                                                                                                                                          |
| 3 years                                                                                                                                                              |                                                                                                                                                          |
| <b>10. Sterilization method</b>                                                                                                                                      |                                                                                                                                                          |
| 1. This device is single-use.                                                                                                                                        |                                                                                                                                                          |
| 2. Device is sterile and non-pyrogenic.                                                                                                                              |                                                                                                                                                          |
| 3. Sterilized by EtO.                                                                                                                                                |                                                                                                                                                          |
| 4. Sterility Assurance Level (SAL)- 10 <sup>-6</sup>                                                                                                                 |                                                                                                                                                          |
| <b>11. Warranty</b>                                                                                                                                                  |                                                                                                                                                          |
| Guaranteed service life: not applicable, as the medical device is disposable.                                                                                        |                                                                                                                                                          |
| Guaranteed storage life: 3 years.                                                                                                                                    |                                                                                                                                                          |
| This product is subject to standard Medtronic warranty terms.                                                                                                        |                                                                                                                                                          |

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