
24

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APPROVED

abia DT disposable tips

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DC.001.01.3-IU

1. Product information

Product name: Disposable non-sterile conductive tips for use on robotic systems.
The product is manufactured in two variants: abia DT, 300 µl and abia DT, 1.100 µl.
This IFU covers the following variants of the IVD product: abia DT, 300 µl and abia DT, 1.100 µl.
Product codes are listed in Table 1.

Product codes

Table 1

Medical device code	Medical device code of manufacturing site	Title	Tips per tray, pcs	Tips per package, pcs
DC.001.01.3	49000-0000	abia DT, 300 µl	96	960
DC.001.11.3	49001-0000	abia DT, 1.100 µl	96	960

2. Intended use

For reagent dosage in compatible IVD robotic systems by the trained health care professionals in clinical laboratory conditions.

3. Indications/Contraindications

For reagent dosage in compatible IVD robotic systems. For in vitro diagnostic use.
For professional use.
For use by the trained health care professionals in clinical laboratory conditions.
There are no specific contraindications for use.

4. Warnings and precautions

Non-sterile IVD product.
Disposable device. Do not re-use.
Do not use the product without labelling or in case labelling/package is damaged.
Do not use after expiration date.
Do not use robotic systems that are not recommended by the manufacturer.
Follow general laboratory precautions, wear protective clothing and protective equipment (gloves).
Low level of personal and environmental risk.

5. Principal technical characteristics

Table 2

Potential risk class: 1	abia DT, 300 µl	abia DT, 1.100 µl
Colour	black	black
Sterility	non-sterile	non-sterile
Conductivity of blank tip	6,20 – 6,30 kΩ	6,20 – 6,30 kΩ
Working volume, µl	300	1.100
Maximum volume, µl	300	1.100
Coaxiality, mm	< 1	< 2
External diameter, mm	7.80	7.80
Internal diameter, mm	0.50	0.80
Lenght, mm	55.0	97.0
Weight, g	0.40	0.88

6. Disposal

Medical device must be disposed as potentially infectious medical waste. Acc to Sanitary Regulations SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for the management of medical waste", this type of waste is classified as Class B (epidemiologically hazardous waste).

The system for collection, temporary storage and transportation of medical waste should be organized in accordance with SanPiN 2.1.7.2790-10 "Sanitary and epidemiological requirements for the management of medical waste".

Collection, temporary storage and disposal of waste should be carried out in accordance with the scheme for handling of medical waste established in the medical and / or pharmaceutical organization.

Warning: Single use medical device! Do not reuse! Dispose off in a container for medical waste! Dispose of used medical device together with household waste is prohibited.

After the expiration date, unused devices can be disposed of as a non-hazardous medical waste of Class A (epidemiologically safe waste approximated to the composition of solid household waste).

Follow SanPiN 2.1.7.2790-10, as well as local / national / federal requirements for medical waste disposal and follow general precautions.

7. The order of application

Open the box with the tips. Put the required quantity of trays with the tips into the specialized compartment of a compatible robotic system.

Further actions are carried out inside used robotic system and described in appropriate Instruction/ Manual of use for a compatible robotic system.

Conditions of Use: at 5-30 °C, humidity of 40-60 %. Atmosphere pressure is not specified.

Compatible robotic IVD systems are listed in Table 3.

4

Compatible robotic systems

Table 3

Robotic system	Producer
GEMINI System	STRATEC Biomedical AG
EVOLIS™ System	Bio-Rad Laboratories, Inc.
BEP 2000 Advance System	Siemens Healthcare GmbH
ETI-Max 3000	DiaSorin S.p.A.
Elisys Quattro	Human GmbH
Immunomat™ TWINSsystem	Institut Virion/Serion GmbH
EUROIMMUN Analyzer I	EUROIMMUN AG

8. Interaction with other medical devices

There are no adverse effects when used in combination with other medical devices, when using as intended in accordance with the instructions for use.

9. Shelf life. Storage and transportation

Shelf life is 5 years. Please check the manufacturing date before use. To be stored at 5-30°C, humidity of 40-60 %. Atmosphere pressure is not specified.

Can be transported by all types of enclosed transport at 5-30°C, humidity of 40-60 %.

Atmosphere pressure is not specified.

10. Requirements for maintenance and repair

There are no parts, which need service, in the product and it is not subject to maintenance and repair.

11. Guarantee

This product is designed, manufactured, tested and packaged in accordance with all applicable requirements.

Manufacturer guarantees conformity of the product to the requirements of regulatory and technical documentation. Quality and safety of the product is guaranteed within established shelf life in accordance with storage/transportation conditions. Do not use the product without labelling or in case labelling/package is damaged.

Complaints should be sent to the authorized representative of the manufacturer in the Russian Federation.

Responsible manufacturer:
 AB Diagnostic Systems GmbH
 Sportfliegerstraße 4
 12487 Berlin, Germany
 Phone: +49 30 208 987 160
 Fax: +49 30 208 987 199
 Email: info@ab-ds.de
 Website: www.ab-ds.de

Authorized representative of the manufacturer in the Russian Federation:
 Research and production company "Diagnostic Systems", Ltd.
 (RPC "Diagnostic Systems", Ltd., Russian Federation)
 47 Komintern Street, 603014, Nizhny Novgorod, Russia
 The address to send complaints:
 22 Yablonevaya Street, 603093, Nizhny Novgorod, Russia
 Phone: +7 (831) 434-86-83, 8 (800) 555-03-00
 Email: info@npods.ru, help-ds@npods.ru
 Website: www.npods.ru

12. Key to symbols used

Manufacturer	
For In vitro diagnostic use	
Catalog number	
Batch code	
Expiry date	
Storage temperature limitation	
Do not use if package is damaged	
Quantity (in a box)	
Non-sterile	
Date of manufacture	
Can not be reused	
Consult Instructions for use	
CE marking	

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