800

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Calibrator	Con. and Pos. of F Activity (Isozyme P) (Isozyme Q)	Con. and Pos. of M isozyme (Isozyme Q)
	(concentration) (position)	(concentration) (position)
Calibrator (1)	Blank concentration [S1] Concentration value of calibrator F [S2]	Blank concentration [S1] 0 0
Calibrator (2)	0 [S3] (Isozyme M Calibrator) 0 [S4] (Isozyme N Calibrator)	0 [S3] (Isozyme M Calibrator) 0 [S4] (Isozyme N Calibrator)
Calibrator (3)		
Calibrator (4)		

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Offline Status	Temperature:	Current User: 001 (Management)	2008-8-8 11:10:04
----------------	--------------	--------------------------------	----------------------

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 Sample Register	 Test Result	 Reagent Info	 Calibration	 QC	 System Setup	 Management	 Maintenance	 Help
--	--	---	--	---	---	---	--	---

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Stand-by...

Temperature: 37.0 C
Current User: 001 (Management)
6/16/2009 08:35 Tuesday

Sample Register
Test Result
Reagent Info.
Calibration
QC
System Setup
Management
Maintenance
Help

Search Result

Sample search
Item search

Test Date From ☒ 5/ 5/2009
Name
Case No.
To ☒ 6/16/2009
Sample No.
Bed No.
Barcode
Dept.
Doctor

Stat	Test Date	Sample Type	Case No.	Patient Name	Sex	Age	Unit	Dept.	Bad N
☐	6/11/2009	Serum							
☒	6/11/2009	Serum							
☐	6/11/2009	Serum							
☐	6/11/2009	Serum							
☐	6/11/2009	Serum							
☐	6/11/2009	Serum							
☐	6/11/2009	Serum							
☒	6/12/2009	Serum					Y		
☒	6/12/2009	Serum		1			Y		
☒	6/12/2009	Serum							
☒	6/12/2009	Serum							
☒	6/12/2009	Serum							
☐	6/12/2009	Serum		1.D8			Y		
☐	6/12/2009	Serum							
☐	6/12/2009	Serum							

Test Item	Test Result	Unit	Reference Range	Hint	Rerun Results	Modify
TP	74.3	g/l				

Record number conform to condition: 215
Print list
Edit Result
Search
Print
Close

Start Analysis
Stop Samples
Stop
Sys Monitor
Alarm Info.
User Logout
Exit
Sleep
Off-line
S 800 Auto Chemistry Analyser

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Save failure, please input the concentration of standard liquid!

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Start Analysis
Stop Samples
Stop
Sys Monitor
Alarm Info.
User Logout
Exit
Sleep
On-line
S 800 Auto Chemistr, Analyser

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Select function key by single mouse click.

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Stand-by... Temperature: 37.0 C Current User: 001 (Management) 6/16/2009 08:36 Tuesday

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System Setup

Chemistry Parameters Profile Item Calculated Item Cross Contamination Report Format ISE Setup Other Setup Manual Item Host Communication

Modify Profile Item

No.: Profile Item Name:

☒ ALT	☐ CRE
☒ AST	☐ Mg
☐ TP	☐ UA
☒ ALB	☐ TC
☒ P	☐ TG
☐ TBIL	☐ HDL-C
☐ DBIL	☐ LDL-C
☐ TBA	☐ GLU
☐ CHE	☐ FMN
☐ ALP	☐ LDH
☐ GGT	☐ HBDH
☐ GLDH	☐ CK
☐ LAP	☐ CK-MB
☐ ICDH	☐ Ca
☐ AMY	☐ Cl
☐ BUN	

No.	Profile Name	Profile Item

Add Delete

Close

Start Analysis Stop Sampling Stop Sys. Monitor Alarm Info User Logout Exit Sleep Off-line S800 Auto-Chemistry Analyser

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Help

QC

QC Register
Individual QC
Monthly QC

Execute	Pos.	QC Name	QC Lot No	Sample Type	QC Item	Target Mean	Target SD
☐	C1	LANDOX-1	456123	Serum	ISE-Na+	123	12
☐	C1	LANDOX-1	456123	Serum	ISE-K+	12	2
☐	C1	LANDOX-1	456123	Serum	ISE-Cl-	10	2
☐	C1	LANDOX-1	456123	Serum	ALT	10	2
☐	C1	LANDOX-1	1	Serum	ALB	15	5
☒	C1	LANDOX-1	1	Serum	TP	15	5

QC Name:
Pos.: C1
QC Lot No:
Sample Type:
QC Test:
Target Mean:
Target SD:

Add
Delete

QC Name Setup
QC Regulation

Close

Start Analysis
Stop Sampling
Stop
Sys Monitor
Alarm Info
User Logout
Exit
Sleep
Off-line
5800 Auto-Chemistry Analyser

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Stand by...

Temperature: 37.0 C
Current User: 001 (Management)
6/16/2009 08:41 Tuesday

Sample Register

Test Result

Reagent Info

Calibration

QC

System Setup

Management

Maintenance

Help

Maintenance item list

Instrument reset
Water tank
Light quantity checkup
Cell blank test
Syringe exhaust
Rinsing detergent pipeline exhaust
Rinsing reaction cuvette
Rinsing ISE
Rinsing reaction cuvette and ISE
Rinsing incubation bath
Sample probe vertical checkup
Sample probe horizontal checkup
Reagent probe vertical checkup
Reagent probe horizontal checkup
Stirring mechanism horizontal checkup
Mechanism operation checkup
Barcode reader checkup
ISE checkup
ISE rinsing reagent pipeline
ISE pipeline exhaust
Auto Effluent pipeline rinsing
Manual Effluent pipeline rinsing

Cell Blank Check

Cup No.	N	340	380	405	450	480	505	546	570	600	660	700	750
1	7	14270	13671	13495	13114	12944	12822	15515	12471	12352	12173	12034	1189
2	7	-91	-83	-73	-55	-46	-40	-34	-32	-28	-24	-26	-24
3	7	30	5	-1	1	-1	-4	-15	-22	-32	-39	-47	-54
4	7	107	70	60	48	40	34	12	2	-9	-23	-37	-55
5	7	41	8	1	3	-1	-7	-28	-36	-50	-65	-76	-86
6	7	55	25	19	12	6	3	-16	-23	-32	-44	-54	-66
7	7	48	19	13	12	10	8	-4	-10	-16	-23	-36	-47
8	7	-26	-20	-14	3	11	14	11	8	1	1	-7	-14
9	7	-1	-2	2	14	19	19	13	9	3	-1	-9	-15
10	7	2	3	11	27	33	35	25	21	15	7	-3	-13
11	7	46	23	19	28	28	25	14	6	-4	-11	-23	-32
12	7	52	30	27	32	35	36	30	27	19	14	5	-5
13	7	-14	-24	-21	-7	-4	-4	-12	-15	-25	-29	-35	-44
14	7	13	3	6	16	17	16	10	6	1	-1	-6	-9
15	7	91	46	39	36	28	22	1	-7	-19	-33	-44	-57
16	7	-4	-7	0	20	23	21	13	7	-2	-9	-13	-14
17	7	54	7	-2	-8	-15	-19	-45	-53	-63	-82	-94	-111

Print
Execute
End Maintenance
Close

Start Analysis
Stop Sampling
Stop
Sys. Monitor
Alarm Info.
User Logout
Exit
Sleep
Off-line

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Reagent residual volume scan begin...

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Reagent residual volume scanning... Temperature: 37.0 C Current User: 001 (Management) 2009.7.13 08:28

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Residual Reagent Vol. (ml)
 Reference liquid Internal standard liquid Diluent liquid
 Residual: 498.94 997.95 498.98
 Save

Calibration Parameters

STD	Conc.	Pos.	Sample Vol.
(1)	120	S4	15
(2)	160	S5	15
(3)	132	S6	15

Reference Range
☐ Specific Value:
 Age - Male - Female
 0 - 0 Years 0 0 0 0
 0 - 0 Years 0 0 0 0
 0 - 100 Years 0 0 0 0
☒ Default value: 20 - 30
 Save
 Close

Start Analysis Stop Sampling Stop Sys. Monitor Alarm Info User Logout Exit Sleep Offline S800 Auto-Chemistry Analyser

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Calibration type	Standard solution volume	Object	Applied assay mode	Application example
Blank calibration	Blank solution	renew reagent blank value	All calibration methods	K coefficient method, when standard solution test is omitted
Span calibration	One point except of blank solution	Renew k value	Two point linear, multi-points linear	Recheck standard solution point 1
2 points calibration	Blank solution and span points	Renew reagent blank value and K value	Two point linear, multi-points linear	Linear 2 points method finish standard curve. Multi-point method, omitted the quantity of calibration solution
Full points calibration	All concentration registered	Renew standard curve by all points	Multi-point linear, isozyme Q, isozyme P, nonlinear working curve calibration	Multipoint standard curve, isozyme method

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Offline

S800 Auto-Chemistry Analyser


Calibration application registry successfully.

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Reaction cuvette rinsing... Temperature: 37.0 C Current User: 001 (Management) 6/16/2009 09:08 Tuesday

Sample Register Test Result Reagent Info Calibration QC System Setup Management Maintenance Help

Sample Analysis

Sample Info:
 Sample No.: 2 Disk/Pos.: 1 2 Bar Code: Repeat: 1
 Sample Type: Serum Sample Vol.: Normal Cup Type: Tube Test Date: 6/16/2009
☐ Sample blank ☐ Stat ☐ Dilute

Item Info:

☐ ALT	☐ BUN	☐ Cl
☐ AST	☐ CRE	
☐ TP	☐ Mg	
☐ ALB	☐ UA	
☐ P	☐ TC	
☐ TBIL	☐ TG	
☐ DBIL	☐ HDL-C	
☐ TBA	☐ LDL-C	
☐ CHE	☐ GLU	
☐ ALP	☐ FMN	
☐ GGT	☐ LDH	
☐ GLDH	☐ HBDH	
☐ LAP	☐ CK	
☐ ICDH	☐ CK-MB	
☐ AMY	☐ Ca	

Profiling

Stat Disk Pos. S.No. Name
 1 1 1 Jack Jones

Previous Next Register Sample Patient Info.
 Emergency Test Delete Samples Batch Register Close

Start Analysis Stop Sampling Stop Sys. Monitor Alarm Info User Logout Exit Sleep Off-line
 F2 F3 F4 F5 F6 F7 F8

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ISE is rinsing... Temperature: 37.0 C Current User: 001 (Management) 2009.7.13 08:39

Sample Register Test Result Reagent Info Calibration QC System Setup Management Maintenance Help

System Setup

Chemistry Parameters Profile Item Calculated Item Cross Contamination Report Format ISE Setup Other Setup Manual Item Host Communication Reagent Setup

Analyze Parameters Calibration Parameters Range Parameters

Test Item: ALT Test Full Name: Decimal Digit: 0 Unit: U/L ☐ Sample Blank

Assay: Rate A Test Time: 10 Point: 21 31 0 0 Control Interval: 0

Main Wave: 340 Second Wave: 405 Instrument Factor(Y = aX + b) a = 1 b = 0 ☒ Always dilution

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Stand by... Temperature: 37.0 C Current User: 001 (Management) 2009-7-13 08:32

Sample Register Test Result Reagent Info Calibration QC System Setup Management Maintenance Help

Sample Analysis

Sample Info: Sample No.: 2 Diet/Reagent: 1 x 1 2 Bar Code: Report: 1

Sample Type: Serum ☐ Sample blank

Item Info:

☒ ALT	☐ BUN	☐
☐ AST	☐ CRE	☐
☐ TP	☐ Mg	☐
☐ ALB	☐ UA	☐
☐ P	☐ TC	☐
☐ TBIL	☐ TG	☐
☐ DBIL	☐ HDL-C	☐
☐ TBA	☐ LDL-C	☐
☐ CHE	☐ GLU	☐
☐ ALP	☐ FMN	☐
☐ GGT	☐ LDH	☐
☐ GLDH	☐ HBDH	☐
☐ LAP	☐ CK	☐
☐ ICDH	☐ CK-MB	☐
☐ AMY	☐ Ca	☐

Current Test Date: 2009-7-13

Sample No. Range: (From) 1 (To) 1

Disk: 1 Pos.: 1

Sample Type: Serum Sample Vol.: Normal Cup Type: Tube

Profiling

Item Info:

☒ ALT	☐ LAP
☐ AST	☐ ICDH
☐ TP	☐ AMY
☐ ALB	☐ BUN
☐ P	☐ CRE
☐ TBIL	☐ Mg
☐ DBIL	☐ UA
☐ TBA	☐ TC
☐ CHE	☐ TG
☐ ALP	☐ HDL-C
☐ GGT	☐ LDL-C
☐ GLDH	☐ GLU

Register Close

Previous Next Register Sample Patient Info. Emergency Test Delete Samples Batch Register Close

Start Analysis Stop Sampling Stop Sys Monitor Alarm Info User Logout Exit Sleep Off-line S800 Auto-Chemistry Analyser

Batch registry failure, existed sample number is contained in selected sample number range.

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Нажмите "Patient Info." в меню "sample register":

Stand-by... Temperature: 37.0 C Current User: 001 (Management) 6/16/2009 09:12 Tuesday

Sample Register Test Result Reagent Info Calibration QC System Setup Management Maintenance Help

Sample Analysis

Patient Info:

Sample No.: 1 Case No.: 12 Patient Type: In hospital

Patient Name: Jack Jones Age: 21 Y Sex: F

Dept.: Doctor: Bed No.: 02

Checker: 001 Auditor: 001 Sending Date: 6/16/2009

Diagnosis: Remark:

Previous Next Register Patient Sample Info

Emergency Test Delete Samples Batch Register Close

Stat Disk Pos. S.No. Name

1	1	1	Jack Jones
---	---	---	------------

Start Analysis Stop Sampling Stop Sys. Monitor Alarm Info. User Logout Exit Sleep Off-line S 800 Auto-Chemistry Analyser

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Help

Rerun item setup

Sample No.	Test Item	Sample Pos.	Rerun reason
1	ISE-K+	11	

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


Stand by...

Temperature: 37.0 C
Current User: 001 (Management)
2009.8.25 09:47
Tuesday

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☒ The results within current day
☐ The results within three days

Send To Host
Delete Result
Superadd Result
Audit
Batch Audit
Reaction Curve
Comm Log
Search Result
Rerun
Preview
Batch Print
Close

Records of the total number of samples: 7

Start Analysis F2	Stop Sampling F3	Stop F4	Sys Monitor F5	Alarm Info. F6	User Logout F7	Exit F8	Sleep	Off-line	S 800 Auto-Chemistry Analyser
----------------------	---------------------	------------	-------------------	-------------------	-------------------	------------	-------	----------	-------------------------------

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Stand by... Temperature: 37.0 C Current User: 001 (Management) 2009-8-25 09:48 Tuesday

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Records of the total number of samples: 7

☒ The results within current day
☐ The results within three days

Send To Host Delete Result Superadd Result Audit Batch Audit Reaction Curve
 Comm Log Search Result Rerun Preview Batch Print Close

Start Analysis F2 Stop Sampling F3 Stop F4 Sys Monitor F5 Alarm Info F6 User Logout F7 Exit F8 Sleep Off-line F9 SB00 Auto-Chemistry Analyser

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Calibration Type	Standard solution volume	Calibration result	Calibration method	Application example
Blank calibration	Blank solution	renew reagent blank value	All points calibration method	K coefficient method, when standard solution test is omitted.
Span calibration	One point except blank solution.	Renew k value	Two points linear and multi-points linear Calibration Logit- Log3P Logit- Log4P	Recheck standard solution point 1
2 points calibration	Blank solution and span points	Renew reagent blank value and K value	Linear two points and multi-points linear calibration Logit- Log3P Logit- Log4P	Linear 2 points method finish standard curve. Multi-point method omitted the quantity of calibration solution.
Full points calibration	All registered calibrator	Renew standard curve by all points	Linear multi-point Calibration isozyme Q isozyme P nonlinear working curve calibration	Multipoint standard curve isozyme method

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Reagent Probe	From Test	From Type	To Test	To Type	Detergent Pos.	Vol.

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Stand by... Temperature: 37.0 C Current User: 001 (Management) 2009.7.13 08:41

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Reagent Probe Reaction cuvette Sample Probe

Test Item	R1 Pos.	R1 Vol.	R2 Pos.	R2 Vol.

Test Item: ALT

R1

R1 Vol.:

R2

R2 Vol.:

Add Delete Close

Start Analysis End Analysis Start Analysis Alarm Info End Analysis Sleep Off line S 800 Auto-Chemistry Analyser

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Нажмите "Sample Probe" для установки "sample cross contamination obviation":

Stand by... Temperature: 37.0 C Current User: 001 (Management) 2009 7 13 08:42

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System Setup

Chemistry Parameters | Profile Item | Calculated Item | Cross Contamination | Report Format | ISE Setup | Other Setup | Manual Item | Host Communication | Reagent Setup

Reagent Probe | Reaction cuvette | Sample Probe

Test Item	Pos.

Test Item:

Pos: W1

Add Delete Close

Start Analysis F2 | Stop Compiling F3 | Stop F4 | Alarm Info F6 | Sleep | Online | S 800 Auto-Chemistry Analyser

Please input correct reagent sampling volume:20-350ul.

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Stand by... Temperature: 37.0 C Current User: 001 (Management) 2009 7 13 08:42

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☒ Print Report Endnote

First line of under title

This report only effective for sample.

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Close

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Offline status		Temperature: 37.0 C		Current User: 001 (Management)		2009 7 13 08:43	
Sample Register	Test Result	Reagent Info	Calibration	QC	System Setup	Management	Maintenance
System Setup							
Chemistry Parameters	Profile Item	Calculated Item	Cross Contamination	Report Format	ISE Setup	Other Setup	Manual Item
Barcode setup ☐ Check Barcode device ☐ Scan barcode before test sample ☐ Scan reagent open area barcode ☐ User reagent barcode a parameter ISE Setup ☒ Execute ISE Time waken setup ☒ Execute time waken Waken date: 2009- 7-13 Waken time: 8:24:59 Reagent Alarm Reagent Residual Vol.(L): 3 ml Reagent Residual Vol.(S): 1 ml ISE Reagent Alarm Reference liquid: 50 ml Internal standard liquid: 50 ml Diluent liquid: 50 ml Save Close							
Display	Print Sample	Stop	System Alarm	Alarm On	User Logout	Exit	Back
F1	F2	F3	F4	F5	F6	F7	F8
							On-line
S 800 Auto-Chemistry Analyser							

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Offline status			Temperature: 37.0 C		Current User: 001 (Management)		2009 7 13 08:43																	
System Setup																								
Chemistry Parameters	Profile Item	Calculated Item	Cross Contamination	Report Format	ISE Setup	Other Setup	Manual Item	Host Communication																
<table style="width: 100%;"> <tr> <td style="width: 20%;"> <table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 10%;">Test Item</th> <th style="width: 90%;">Test Full Name</th> </tr> <tr> <td>▶</td> <td></td> </tr> </table> </td> <td style="width: 20%;"> Test Item: </td> <td style="width: 20%;"> Test Full Name: </td> <td style="width: 20%;"> Unit: ▼ </td> <td style="width: 20%;"> Reference: - </td> <td style="width: 20%;"></td> </tr> <tr> <td colspan="3" style="text-align: center;"> Add </td> <td colspan="3" style="text-align: center;"> Delete </td> </tr> </table>									<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 10%;">Test Item</th> <th style="width: 90%;">Test Full Name</th> </tr> <tr> <td>▶</td> <td></td> </tr> </table>	Test Item	Test Full Name	▶		Test Item:	Test Full Name:	Unit: ▼	Reference: -		Add			Delete		
<table border="1" style="width: 100%; border-collapse: collapse;"> <tr> <th style="width: 10%;">Test Item</th> <th style="width: 90%;">Test Full Name</th> </tr> <tr> <td>▶</td> <td></td> </tr> </table>	Test Item	Test Full Name	▶		Test Item:	Test Full Name:	Unit: ▼	Reference: -																
Test Item	Test Full Name																							
▶																								
Add			Delete																					
								Close																
Stop/Cancel	Stop/Cancel	Stop	Stop/Cancel	Stop/Cancel	User Logout	Exit	Stop	On-line																
F2	F3	F4	F5	F6	F7	F8		S 800 Auto-Chem'istry Analyser																

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Offline status				Temperature: 37.0 C		Current User: 001 (Management)		2009.7.13 08:43	
System Setup									
Chemistry Parameters		Profile Item		Calculated Item		Cross Contamination		Report Format	
ISE Setup		Other Setup		Manual Item		Host Communication		Reagent Setup	

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Offline status				Temperature: 37.0 C		Current User: 001 (Management)		2009 7 13 08:44	
System Setup									
Chemistry Parameters Profile Item Calculated Item Cross Contamination Report Format ISE Setup Other Setup Manual Item Host Communication Reagent Setup									
ISE Setup Reference liquid barcode: Volume: 334.53 ml OK Internal standard liquid barcode: Volume: 613.14 ml OK Second line of under title: Volume: 1035.75 ml OK Alkaline detergent Barcode: Volume: 1000.0 ml OK Anti-bacterial phosphor-free detergent Barcode: Volume: 300.0 ml OK									
									Close
Reagent Info Reagent Barcode Unit Sys. Admin. Admin Info User Logout Exit Print On-line S 800 Auto-Chemistry Analyser									

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Mnemonic:

Add Delete Close

User Logout Exit On-line S 800 Auto-Chemistry Analyser

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Offline Status				Temperature:		Current User: 001 (Management)		2009-6-17 13:33 Wednesday							
Management															
User Info Hospital Info Other Info Workload Statistics Database Maintenance System Log															
Dept Doctor <table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 15%;">Dept. No.</th> <th style="width: 35%;">Dept</th> <th style="width: 50%;">Mnemonic</th> </tr> </thead> <tbody> <tr> <td>12</td> <td>q</td> <td>q</td> </tr> </tbody> </table> Dept. No.: Dept.: Mnemonic: Add Delete Close										Dept. No.	Dept	Mnemonic	12	q	q
Dept. No.	Dept	Mnemonic													
12	q	q													
Stat Analysis F2 Stop Sampling F3 Stop F4 Stop Monitor F5 Alarm Info F6 User Logout F7 Exit F8 Stop On-line S 800 Auto-Chemistry Analyser															

Single click table to delete the record.

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Offline Status
Temperature:
Current User: 001 (Management)
2009 6 17 13:34

Sample Register
 Test Result
 Reagent Info
 Calibration
 QC
 System Setup
 Management
 Maintenance
 Help

Management

User Info.
Hospital Info.
Other Info.
Workload Statistics
Database Maintenance
System Log

Dept. Doctor

Doctor No.	Name	Dept.	Mnemonic
12	zdd	2	we

Doctor No.:

Name:

Dept:

Mnemonic:

Add
Delete
Close

Sample Register
Test Result
Reagent Info
Calibration
QC
System Setup
Management
Maintenance
Help

Start Analysis
Data Management
Print
Data Monitor
Alarm Info
User Logout
Exit

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Offline Status				Temperature:	Current User: 001 (Management)		2009.6.17 13:34 Wednesday	
Management								
User Info.		Hospital Info.		Other Info.		Workload Statistics		Database Maintenance
Sample Type		Patient Type		Clinic Diagnosis		Remark		Unit

No.	Patient Type	Mnemonic
1	In hospital	l
2	outpatient	o
3	medicare	m

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Offline Status
Temperature:
Current User: 001 (Management)
2009.6.17 13:35 Wednesday

Sample Register
 Test Result
 Reagent Info
 Calibration
 QC
 System Setup
 Management
 Maintenance
 Help

Management

User Info
Hospital Info
Other Info
Workload Statistics
Database Maintenance
System Log

Sample Type
Patient Type
Clinic Diagnosis
Remark
Unit

No.	Clinic Diagnosis	Mnemonic
10	input	qw

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
Offline Status				Temperature:		Current User: 001 (Management)		2009.6.17 13:35 Wednesday							
Management															
User Info Hospital Info Other Info Workload Statistics Database Maintenance System Log															
Sample Type Patient Type Clinic Diagnosis Remark Unit															
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>No.</th> <th>Remark</th> <th>Mnemonic</th> </tr> </thead> <tbody> <tr> <td>001</td> <td>exma</td> <td>122</td> </tr> </tbody> </table>					No.	Remark	Mnemonic	001	exma	122	No.: Remark: Mnemonic:				
No.	Remark	Mnemonic													
001	exma	122													
					Add Delete Close										
Start Analysis F2	Start Sampling F3	Print F4	Print Manual F5	Alarm Info F6	User Logout F7	Exit F8	Start	On-line	S 800 Auto Chemistry Analyser						

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Offline Status			Temperature:		Current User: 001 (Management)		2009.6.17 13:35 Wednesday		
Management									
User Info		Hospital Info		Other Info		Workload Statistics		Database Maintenance System Log	
Sample Type		Patient Type		Clinic Diagnosis		Remark		Unit	
Unit g/dl mg/dl ug/dl mg/l umol/l mEq/l U/L g/l mol/l mmol/l		Unit		mg/l					
						Add		Delete	
Close									
User Logout		Exit		On-line		S 800 Auto-Chemistry Analyser			
F7		F8							

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The screenshot displays the 'Database Maintenance' window of the S 800 Auto-Chemistry Analyser. At the top, a status bar shows 'Offline Status', 'Temperature:', 'Current User: 001 (Management)', and the date/time '2009.6.17 13:36 Wednesday'. Below this is a row of icons for 'Sample Register', 'Test Result', 'Reagent Info.', 'Calibration', 'QC', 'System Setup', 'Management', 'Maintenance', and 'Help'. The 'Management' menu is open, showing sub-options: 'User Info.', 'Hospital Info.', 'Other Info.', 'Workload Statistics', 'Database Maintenance' (selected), and 'System Log'. Inside the 'Database Maintenance' window, there are two radio buttons: 'Database Backup' (selected) and 'Database Recovery'. Below them, a text box says 'Start after select the backup storage folder, the default path is the BACKUP folder under the application program.' A text field contains the path 'C:\Program Files\RuSoftware of CS-400 auto-Chemistry Analyser\backup\'. There are 'Select Path' and 'Backup' buttons. A note at the bottom states 'Database backup and recovery should be executed in off-line status.' A 'Close' button is in the bottom right. At the very bottom, a row of function keys (F2-F8) and 'Stop' and 'On-line' buttons are visible, along with the text 'Record number conform to condition: 1'.

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Database Maintenance

System Log

Date	Operator	Event
2009-6-17 13:20:33	001	Successfully login.

Type

☒ Login
 ☐ Operation
 ☐ Maintenance
 ☐ Alarm

Date

From

2009- 6-17

To

2009- 6-17

Search

Delete

Close

Test Results

Sample Register

Stop

Log Monitor

Alarm Info

User Logout

Exit

Stop

On-line

S 800 Auto-Chemistry Analyser

F2

F3

F4

F5

F6

F7

F8

Record number conform to condition: 1

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Code	Level	Name	Describe	Countermeasure
0-1	Warning	Absorbance over	The absorbance to be used for calculation after cell blank correction exceeds 3.3ABS	1.Check if the reagent is collected or positioned in the right way. 2.Check if the incubation bath is contaminated 3.Check if there is the impurity interfuse. 4.Check if the reaction cuvette is blurred. 5.Check if there is an obstacle on the optical path of photometer.
0-2	Warning	Reagent blank horizontal checkup over.	The mean value of two tests with standard solution 1 is larger than the setup value of cell blank horizontal check.	Check if reagent is still valid.
0-3	Warning	Divergence checkup over	The difference between two times test standard solution absorbance is larger than deviation permissible absorbance.	1.Check if the value of deviation permissible absorbance is correct 2.Check if sample syringe is leakage and if the incubation bath is contaminated. 3.Check if there is a cross contamination.
0-4	Warning	K factor check value over	The difference between calculate K factor and previous one is more than plus or minus 20%	1.Check if K factor is proper. 2.Check if the standard solution is deteriorate 3.Check if reagent is deteriorate.
0-5	Warning	Technical limit over	The measured value is out of the technical limit range of test item.	1.Dilution the sample and rerun the test. 2.Check setup value.
0-6	Warning	Linearity abnormal	Calculation in the specified photometric range, and there is a different exceeding the linearity limit value.	1.Check if there is impurity in the sample 2.Check if there is impurity and bubble in the reagent. 3.Dilution the sample and rerun the test. 4.Check Stirring mechanism.
0-7	Warning	Sensitivity abnormal	Sensitivity is checked for linear(2-6 points), non-linear or Isozyme-P calibration. This error is indicated if the difference in mean absorbance between standard solution 1 and standard solution (N) is smaller than the sensitivity limit (input value).	1.Call up main screen, check if sensitivity limit value is properly setup. 2.Check if reagent is deteriorate. 3.Check the input concentration value of standard solution is properly setup. 4.Check if standard solution deteriorate.

Code	Level	Name	Describe	Countermeasure
0-22	Warning	Detergent placed error at 45 position on R1 disk	Detergent should be placed at 45 position on R1 disk or barcode stick error at 45 position on R1 disk.	1. Check the barcode stick information. 2. Scan barcode again after check.
0-23	Warning	R1 disk barcode blank out	1. When reagent residual volume is less than 2 barcode is used after blank out.	1. Check the barcode stick information. 2. Scan barcode again after check.
0-24	Warning	R1 disk barcode invalid	1. Barcode cleanup failure 2. Barcode reagent name do not occur in "Barcode item setup"	1. Check barcode stick information. 2. Scan barcode again after check. 3. Corresponding item is added in barcode item setup.
0-25	Warning	R2 disk barcode repeat	R2 disk barcode repeat	1. Check the barcode stick information. 2. Scan barcode again after check.
0-26	Warning	Reagent R1 and R4 is found on R2 disk	R1 and R4 detergent is found on R1 or R4	1. Check the barcode stick information. 2. Scan barcode again after check.
0-27	Warning	Detergent placed error at 45 position on R2 disk.	Detergent should be placed at 45 position on R2 disk Barcode stick error at 45 position on R2 disk.	1. Check the barcode stick information. 2. Scan barcode again after check.
0-28	Warning	R2 disk barcode blank out	1. When reagent residual volume is less than 2 barcode is used after blank out.	1. Check the barcode stick information. 2. Scan barcode again after check.
0-29	Warning	R2 disk barcode invalid	1. Barcode cleanup failure 2. Barcode reagent name do not occur in "Barcode item setup"	1. Check barcode stick information. 2. Scan barcode again after check. 3. Corresponding item is added in barcode item setup.

Code	Level	Name	Describe	Countermeasure
0-8	Warning	Prozone check value over	In the Antigen addition assay and reaction assay, the prozone check value (PC value) exceeds the specified upper/lower limit.	1. Dilution the sample and return the test. 2. Check if the reagent is collected or positioned in the right way. 3. Check if the specified upper/lower limit is properly set up.
0-9	Warning	Absorbency reaction limit over	In Rate A assay and Rate B assay, The check value of substrate is larger than the setup limit value.	1. Test after dilute the sample. 2. Check if the reagent is collected or positioned in the right way. 3. Check if the check value is properly setup.
0-10	Warning	Calculation disable	1. The denominator becomes zero in calculation. 2. An overflow occurs in logarithmic or exponential calculation. 3. The absorbance of test item is less than the absorbance of standard liquid 1.	1. Check if there is a logical error in calculation formulas. 2. Check the concentration setup of standard liquid when multi-point calibration. 3. Check the absorbance value of photometry point of the test.
0-11	Warning	Calculation disable	In multi-point calibration, the absorbance of standard liquid do not increase when concentration increase.	1. Check if calibration liquid is positioned in the right way. 2. Check if absorbance value of different concentrations increase. 3. Please return the calibration.
0-12	Warning	ISE calculated error	In calculation the denominator became zero	Please check if ISE reagent and sample are enough and in right position
0-13	Warning	ISE request calibration	New ISE calibration is required if overtime calibration occurs in ISE calibration	Execute ISE calibration
0-14	Warning	Concentrated waste liquid pipeline need to Run	Concentrated waste liquid pipeline need to Run	Call up the system maintenance menu, carry out the program "Auto Running Concentrated waste liquid pipeline"
0-20	Warning	Barcode repeat	Barcode repeat on R1 disk	1. Check the barcode stick information. 2. Scan barcode again after check.
0-21	Warning	Reagent R2 and R3 is found on R1 disk.	Reagent R2 and R3 is found on R1 disk	1. Check the barcode stick information. 2. Scan barcode again after check.

Code	Level	Name	Describe	Countermeasure
0-31	Warning	ISE exceed the lower limit of measuring range	ISE QC exceed the lower limit of measuring range	Please check if Controls volume and reagent volume is enough. If the position is right
0-32	Warning	ISE exceed the upper limit of measuring range	ISE QC exceed the upper limit of measuring range	Please check if Controls volume and reagent volume is enough. If the position is right
0-33	Warning	ISE exceed the lower limit of measuring range	ISE QC exceed the lower limit of measuring range	Please check if Controls volume and reagent volume is enough. If the position is right
0-34	Warning	ISE exceed the upper limit of measuring range	ISE QC exceed the upper:lower limit of measuring range	Please check if Controls volume and reagent volume is enough. If the position is right
0-35	Warning	ISE exceed the measuring range	ISE QC exceed the upper limit of measuring range	Please check if sample volume and reagent volume is enough. If the position is right
0-254	Stop	Communication abnormal	Communication when instruction send to host computer. Abnormal occur when communicate with host computer and communication halt	1. Check if serial port line is well connected. 2. Restart the instrument and software and then carry out on-line instruction.
0-255	Stop	Communication abnormal	The R1 stirring mechanism does not reach the upper dead point in ascending motion (at the reaction cuvette side)	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-1	Stop	R1 Stirrer abnormal	The R1 stirring mechanism does not reach the upper dead point in ascending motion	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-2	Stop	R1 Stirrer mechanism abnormal	The R1 stirring mechanism does not reach the upper dead point in ascending motion	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-3	Stop	R1 Stirrer abnormal	The R1 stirring mechanism does not leave the upper dead point in descending motion	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-4	Stop	R1 Stirrer abnormal	The R1 stirring mechanism does not reach the running bath side	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-5	Stop	R1 Stirrer abnormal	The R1 stirring mechanism does not reach the reaction cuvette position in movement to the reaction cuvette side	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-6	Stop	R1 Stirrer abnormal	In resetting, the R1 stirring mechanism does not reach the running bath position	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-7	Stop	R1 Stirrer mechanism abnormal	In resetting, the R1 stirring mechanism does not leave the running bath position	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-8	Stop	R1 Stirrer abnormal	In rotation of R1 stirring mechanism, it is not set at the upper dead point	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-9	Stop	R1 Stirrer abnormal	The R1 stirring mechanism has ascend /descend error	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department
1-10	Stop	R1 Stirrer abnormal	The R1 stirring mechanism has sway error	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department

Code	Level	Name	Describe	Countermeasure
2-1	Stop	R2 Stirrer mechanism abnormal	The R2 stirring mechanism does not reach the upper dead point in ascending motion (at the rinsing bath side)	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-2	Stop	R2 Stirrer mechanism abnormal	The R2 stirring mechanism does not reach the upper dead point in ascending motion (at reaction cuvette side)	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-3	Stop	R2 Stirrer mechanism abnormal	The R2 stirring mechanism does not leave the upper dead point in descending motion.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-4	Stop	R2 Stirrer mechanism abnormal	The R2 stirring mechanism does not reach the rinsing bath position in movement to the rinsing bath side	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-5	Stop	R2 Stirrer mechanism abnormal	The R2 stirring mechanism does not reach the reaction cuvette position in movement to the reaction cuvette side	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-6	Stop	R2 Stirrer mechanism abnormal	In resetting, the R2 stirring mechanism does not reach the rinsing bath position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-7	Stop	R2 Stirrer mechanism abnormal	In resetting, the R2 stirring mechanism does not leave the rinsing bath position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
2-8	Stop	R2 Stirrer mechanism abnormal	In rotation of R2 stirring mechanism, it is not set at the upper dead point.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-9	Stop	R2 Stirrer mechanism abnormal	The R2 stirring mechanism has ascend /descend error	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
2-10	Stop	R2 Stirrer mechanism abnormal	The R2 stirring mechanism has sway error	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
3-1	Stop	Rinsing mechanism abnormal	The rinse mechanism does not reach the upper dead point in ascending motion	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
3-2	Stop	Rinsing mechanism abnormal	The rinse mechanism dose not leave the upper dead point in the descending motion.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
4-1	Stop	Reaction disk abnormal	The reaction disk cannot detect the stop position.	1.This trouble is liable to occur after the reaction disk is washed. In this case, thoroughly wipe water droplets off the bottom of the reaction disk. 2.Check if water droplets adhere to the light coupler and code wheel located below the reaction disk. 3.If trouble can not be restored, contact the service department.
4-2	Stop	Reaction disk abnormal	The reaction disk does not stop at the specified position	1.This trouble is liable to occur after the reaction disk is washed. In this case, thoroughly wipe water droplets off the bottom of the reaction disk. 2.Check if water droplets adhere to the light coupler and code wheel located below the reaction disk. 3.If trouble can not be restored, contact the service department.

Code	Level	Name	Describe	Countermeasure
4-3	Stop	Reaction disk abnormal	In resetting, the reaction disk cannot detect the home position	1. This trouble is liable to occur after the reaction disk is washed. In this case, thoroughly wipe water droplets off the bottom of the reaction disk. 2. Check if water droplets adhere to the light coupler and code wheel located below the reaction disk. 3. If trouble can not be restored, contact the service department.
4-4	Stop	Reaction disk abnormal	In resetting, the first reaction cuvette on the reaction disk does not stop at the specified position	1. This trouble is liable to occur after the reaction disk is washed. In this case, thoroughly wipe water droplets off the bottom of the reaction disk. 2. Check if water droplets adhere to the light coupler and code wheel located below the reaction disk. 3. If trouble can not be restored, contact the service department.
4-5	Stop	Reaction disk abnormal	In rotation of the reaction disk any of the sample probe, reagent probe, starting mechanism and rise mechanism is not at the upper dead point (at the reaction cuvette side) Warning: any other alarm may occur simultaneously	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.
4-6	Warning	CS-alkaline detergent shortage	CS-alkaline detergent place in front of instrument is shortage.	Please add CS-alkaline detergent in detergent bottle.
5-1	S. Stop	Sample probe abnormal	The sample probe does not reach the upper dead point in descending motion (at other than the reaction cuvette side)	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.
5-2	Stop	Sample probe abnormal	The sample probe does not reach the upper dead point in descending motion (at the reaction cuvette side)	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.
5-3	Warning	Sample size is insufficient in sample cup	sample size is insufficient in sample cup	1. Check if there is sample in sample cup 2. Check if sample cup place in the correct position 3. Check if sample cup is bend

Code	Level	Name	Describe	Countermeasure
5-4	Stop	Sample probe abnormal	The sample probe goes down abnormally in descending motion (at the reaction cuvette side)	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department. 3. Check if the sample cup is placed defective 4. Call up the system maintenance menu, carry out "sample probe horizontal checkup" 5. Check if the sample cup is placed defective
5-5	S. Stop	Sample probe abnormal	The sample probe does not leave the upper dead point in descending motion (at the sample cup side)	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.
5-6	Stop	Sample probe abnormal	The sample probe does not leave the upper dead point in descending motion (at the reaction cuvette side)	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.
5-7	S. Stop	Sample probe abnormal	Sample probe stay in liquid level detect effective status.	Please contact service department when this error occur repetitively.
5-8	S. Stop	Sample probe abnormal	The sample probe cannot detect the cell position in rotation to the reaction cuvette side	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.
5-9	S. Stop	Sample probe abnormal	The sample probe does not leave the reaction cuvette position in rotation from position	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.
5-10	S. Stop	Sample probe abnormal	Sample probe descend abnormally when discharge ISE test sample	1. Call up the mechanism check screen and carry out the check program 2. Various malfunctions or trouble can not be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
5-11	Warning	Sample probe abnormal	Abnormal touch occur (at the sample cup side) when sample probe descend.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
5-12	S. Stop	Sample probe abnormal	In rotation of reagent-1 probe, it is not set at the upper dead point.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
5-13	Warning	Sample cup is not placed	Sample probe cannot find sample cup	Check if there is sample cup on sample disk.
5-14	S. Stop	Sample probe abnormal	Sample probe cannot reach sample liquid when descend.	1.Check if there is obstacle on sample disk. 2.Various malfunctions or trouble can not be restored. Contact the service department.
5-15	Stop	Sample probe abnormal	Sample probe mechanism ascend/ descend abnormal	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
5-16	Stop	Sample probe abnormal	Sample probe mechanism sway abnormal	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
5-17	Warning	Sample probe abnormal	Detergent is not detected when sample descending	Check if detergent is placed at W1, W2, W3.
6-1	S. Stop	Sample disk abnormal	The sample disk does not detect the specified position on the outer track	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-2	S. Stop	Sample disk abnormal	The sample disk does not stop in the specified position on the outer track	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
6-3	S. Stop	Sample disk abnormal	The sample disk cannot detect the stop position on the intermediate track.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-4	S. Stop	Sample disk abnormal	The sample disk cannot stop in the specific position on the intermediate track	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-5	S. Stop	Sample disk abnormal	The sample disk cannot detect the stop position on the inner track.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-6	S. Stop	Sample disk abnormal	The sample disk cannot stop in the specific position on the inner track	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-7	Stop	Sample disk abnormal	In resetting, the sample disk cannot detect the home position.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-8	Stop	Sample disk abnormal	In resetting, the sample disk does not leave the home position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-9	Stop	Sample disk abnormal	In resetting the sample disk does not stop at sample 1 position on the outer track.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
6-10	Warning	sample barcode device abnormal	Sample barcode reader cannot be found.	1.Check if barcode device is well connected. 2.Trouble cannot be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
7-1	S. Stop	Sample syringe abnormal	The sample syringe is not at the upper dead point.	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
7-2	S. Stop	Sample syringe abnormal	The sample syringe does not leave the upper dead point.	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
8-1	Stop	R1 probe abnormal	The R1 probe does not reach the upper dead point in ascending motion	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
8-2	Warning	R1 probe abnormal	The R1 probe goes down abnormally in descending motion	1. Check if there is still reagent in the bottle. 2. Check if the reagent bottle is open. 3. Check if the reagent disk cover is put in its proper place. 4. Check if the reagent-1 probe is bent. 5. Check if reagent bottle is lean placed.
8-3	Stop	R1 probe abnormal	The R1 probe does not leave the upper dead point in descending motion	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
8-4	Warning	R1 probe abnormal	R1 probe stay in liquid level detect effective status.	Please contact service department when this error occur repetitively.
8-5	Stop	R1 probe abnormal	The R1 probe cannot detect the cell position in rotation to the reaction cuvette side.	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
8-6	Stop	R1 probe abnormal	The R1 probe cannot leave the cell position in rotation to the reaction cuvette side.	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
8-7	Stop	R1 probe abnormal	In rotation of R1 probe, it is not set at the upper dead point.	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
8-8	Warning	R1 probe abnormal	R1 probe cannot find reagent bottle	Check if there is reagent bottle on R1 reagent disk.
8-9	Stop	R1 probe abnormal	R1 probe cannot reach reagent liquid level when descend	1. Check if there is obstacle on R1 reagent disk. 2. trouble cannot be restored, contact the service department.
8-10	Stop	R1 probe abnormal	R1 probe mechanism ascend/ descend abnormal	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
8-11	Stop	R1 probe abnormal	R1 probe mechanism sway abnormal	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
9-1	Stop	R2 probe abnormal	The R2 probe does not reach the upper dead point in ascending motion	1. Call up the system maintenance menu and carry out the check program "reset". 2. Various malfunctions or trouble can not be restored. Contact the service department.
9-2	Warning	R2 probe abnormal	The R2 probe goes down abnormally in descending motion	1. Check if there is still reagent in the bottle. 2. Check if the reagent bottle is open. 3. Check if the reagent disk cover is put in its proper place. 4. Check if the reagent-1 probe is bent. 5. Check if reagent bottle is lean placed.
9-3	Stop	R2 probe abnormal	The R2 probe does not leave the upper dead point in descending motion	1. Call up the mechanism check screen and carry out the check program "mechanism check". 2. Various malfunctions or trouble can not be restored. Contact the service department.
9-4	Warning	R2 probe abnormal	R2 probe stay in liquid level detect effective status.	Contact the service department when this error occurs repetitively.

Code	Level	Name	Describe	Countermeasure
9-5	Stop	R2 probe abnormal	The R2 probe cannot detect the cell position in rotation to the reaction cuvette side.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
9-6	Stop	R2 probe abnormal	The R2 probe cannot leave the cell position in rotation to the reaction cuvette side.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
9-7	Stop	R2 probe abnormal	In rotation of R2 probe, it is not set at the upper dead point	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
9-8	Warning	R2 probe abnormal	R2 probe cannot find reagent bottle	Check if there is reagent bottle on R2 reagent disk.
9-9	Stop	R2 probe abnormal	R2 probe cannot reach reagent liquid level when descend.	1.Check if there is obstacle on R2 reagent disk. 2.Trouble cannot be restored, contact the service departmen..
9-10	Stop	R2 probe abnormal	R2 probe mechanism ascend/descend abnormal	Various malfunctions or trouble can not be restored. Contact the service department.
9-11	Stop	R2 probe abnormal	R2 probe mechanism ascend/descend abnormal	Various malfunctions or trouble can not be restored. Contact the service department.
10-1	Stop	R1 disk abnormal	The R1 disk cannot detect the stop position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
10-2	Stop	R1 disk abnormal	The R1 disk does not stop at the specified position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
10-3	Stop	R1 disk abnormal	The R1 disk cannot detect the home position.	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
10-4	Warning	Bar code device of R1 reagent disk abnormal	Barcode reader of R1 disk cannot be find	1.Check if barcode device is well connected. 2.Various malfunctions or trouble can not be restored. Contact the service department.
10-5	Stop	R1 reagent number invalidity	R1 reagent number invalidity.	Please contact the service department.
11-1	Stop	R2 disk abnormal	The R2 disk cannot detect the stop position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
11-2	Stop	R2 disk abnormal	The R2 disk does not stop at the specified position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
11-3	Stop	R2 disk abnormal	The R2 disk cannot detect the home position	1.Call up the mechanism check screen and carry out the check program "mechanism check". 2.Various malfunctions or trouble can not be restored. Contact the service department.
11-4	Warning	Bar code device of R2 reagent disk abnormal	Barcode reader of R2 disk cannot be find	1.Check if barcode device is well connected. 2.Various malfunctions or trouble can not be restored. Contact the service department.
11-5	Stop	R2 reagent number invalidity	R2 reagent number invalidity.	Please contact the service department.
14-1	Stop	R1 Syringe abnormal	The R1 syringe is not at the upper dead point	When this error occurs repetitively, please contact service department

Code	Level	Name	Describe	Countermeasure
14-2	Stop	R1 Syringe abnormal	The R1 syringe does not leave from the upper dead point	When this error occurs repetitively, please contact service department
15-1	Stop	R2 Syringe abnormal	The R2 syringe is not at the upper dead point.	When this error occurs repetitively, please contact service department
15-2	Stop	R2 Syringe abnormal	The R2 syringe does not leave from the upper dead point	When this error occurs repetitively, please contact service department
16-1	Warning	Electrolyte SIP Sipper abnormal	SIP sipper does not reach the lower dead point (when reset and operate)	1.Call up the system maintenance menu,first carry out the program "reset",then go on "mechanism check" 2.Various malfunctions or trouble can not be restored. Contact the service department.
16-2	Warning	Electrolyte SIP Sipper abnormal	SIP sipper does not leave the lower dead point	1.Call up the system maintenance menu,first carry out the program "reset",then go on "mechanism check" 2.Various malfunctions or trouble can not be restored. Contact the service department.
17-1	Warning	Electrolyte vacuum Sipper abnormal	The vacuum nozzle does not reach the lower dead point	1.Call up the system maintenance menu,first carry out the program "reset",then go on "mechanism check" 2.Various malfunctions or trouble can not be restored. Contact the service department.
17-2	Warning	Electrolyte vacuum Sipper abnormal	The vacuum nozzle does not leave the lower dead point.	1.Call up the system maintenance menu,first carry out the program "reset",then go on "mechanism check" 2.Various malfunctions or trouble can not be restored. Contact the service department.
18-1	Warning	Electrolyte sample Syringe abnormal	SIP syringe does not at the upper dead point.	Contact the service department.
18-2	Warning	Electrolyte sample Syringe abnormal	SIP syringe does not leave the upper dead point.	Contact the service department.
18-3	Warning	Electrolyte sample Syringe abnormal	The diluent syringe is not at the upper dead point.	Contact the service department.

Code	Level	Name	Describe	Countermeasure
18-4	Warning	Electrolyte sample Syringe abnormal	The diluent syringe is not leave the upper dead point.	Contact the service department.
18-5	Warning	Electrolyte sample Syringe abnormal	The internal standard solution syringe is not at the upper dead point.	Contact the service department.
18-6	Warning	Electrolyte sample Syringe abnormal	The internal standard solution syringe is not leave the upper dead point.	Contact the service department.
19-1	Warning	ISE stop	The ISE function stops due to an alarm.Warning: Displayed if re-start is given in sample stop status.	1.Call up the system maintenance menu,first carry out the program "reset",then go on "mechanism check" 2.Various malfunctions or trouble can not be restored. Contact the service department.
20-1	Warning	Incubation bath water temp abnormal	The water temperature in the incubation bath is over 45C	1.Check if cooling fan is operate normally. 2.Check if cooling fan disproof cover is clogged by dirt 3.Trouble can not be restored, contact the service department.
20-2	Warning	Incubation bath water temp abnormal	The water temperature in the incubation bath is over the range of 37 plus or minus 0.5C (check in operating)	1.Check if the room temperature is within a range of 15-32C. 2.Check if cooling fan is operated normally. 3.Check if cooling fan disproof cover is clogged by dirt. 4.Check if the incubation bath water is circulation. 5.Trouble can not be restored, contact the service department.
21-1	Warning	Incubator water level error	The water level in the incubation bath is too low	1.Clean liquid level sensor probe of incubation bath 2.Check if there are water in purified water machine and pipeline 3. Please restart the instrument. 4.If the trouble can't be restored, contact the service department.
22-1	Warning	Incubation bath water exchange	A period of 24 hours has passed since exchanging water in the incubation bath.	Carry out "Incubator water exchange" or restart the instrument.
23-1	Stop	Water level too low	The water level in the water tank is too low	The distilled water tank is not supplied with water. Check the water pressure, water code, and other water supply systems.

Code	Level	Name	Describe	Countermeasure
28-1	Warning	Waste water in Vacuum Tank	The vacuum tank contains some waste water	Empty the water in vacuum tank according to the (User manual)
31-1	Warning	Reagent shortage(R1)	Reagent residual volume in R1 or R4 reagent bottle on R1 disk is less than the setup value of reagent residual volume in system setup interface	1 Check if reagent residual volume on the reagent information screen is lower than its setup value 2 Check if the setup value of reagent residual volume is reasonable
32-1	Warning	Reagent shortage(R2)	Reagent residual volume in R2 or R3 reagent bottle on R2 disk is less than the setup value of reagent residual volume in system setup interface	1 Check if reagent residual volume on the reagent information screen is lower than its setup value 2 Check if the setup value of reagent residual volume is reasonable
33-1	Warning	Detergent Low(R1)	The Detergent volume on No 4/5 position of R1 disk is less than 3ml	Call for the reagent information screen, check if detergent is enough
34-1	Warning	Detergent Low(R2)	The Detergent volume on No 4/5 position of R2 disk is less than 3ml	Call for the reagent information screen, check if detergent is enough
35-1	Warning	Reagent information setup error(R1)	1 Reagent name setup error, reagent residual test time cannot be checked 2 reagent type setup error, reagent residual test time cannot be checked 3 the item is not set in analyze parameter	1 Check barcode stick information, scan the barcode again after check. 2 register the reagent information manually
36-1	Warning	Reagent information setup error(R2)	1 Reagent name setup error, reagent residual test time cannot be checked 2 reagent type setup error, reagent residual test time cannot be checked 3 the item is not set in analyze parameter	1 Check barcode stick information, scan the barcode again after check. 2 register the reagent information manually
37-1	Warning	ISE Reference Liquid shortage	The remaining volume of Reference Liquid is shortage	1 Add ISE internal standard liquid 2 Call up system maintenance screen, carry out ISE internal standard liquid pipetime washing. 3 Carry out ISE calibration
38-1	Warning	ISE internal standard liquid shortage	Internal standard liquid shortage	
39-1	Warning	ISE Diluent liquid shortage	Diluent liquid shortage.	Please check if ISE reagent volume is enough
41-1	Warning	Sample flushed	Sample registered flushed sampling, stop to adding	
51-1	Warning	Automatic calibration disable	Turning calibration is not carried out when reagent residual volume cannot complete 10 tests	1 Call up the system maintenance screen, first carry out the program "reset", then go on "mechanism check". 2 Various malfunctions or trouble can not be restored. Contact the service department.
60-1	Warning	ISE LEVEL Error	The mean electromotive force of internal standard solution (EAV) is out of the following range at 3 of 5 measuring points(internal standard solution)Na-90.0mV-EAV<-10mV	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
60-2	Warning	ISE LEVEL Error	The mean electromotive force of internal standard solution (EAV) is out of the following range at 3 of 5 measuring points(internal standard solution)K-90.0mV-EAV<-10mV	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
60-3	Warning	ISE LEVEL Error	The mean electromotive force of internal standard solution (EAV) is out of the following range at 3 of 5 measuring points(internal standard solution)CI-80.0mV-EAV<-160mV	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
61-1	Warning	ISE Noise Error	The difference between the maximum and minimum electromotive forces of internal standard solution (EIV) is within the following range at 3 of 5 measuring points (internal standard solution)Na-0.7mV-EIV<-2-FIV(4)	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".

Code	Level	Name	Describe	Countermeasure
61-2	Warning	ISE Noise Error	The difference between the maximum and minimum electromotive forces of internal standard solution (FIV) is within the following range at 3 of 5 measuring points (internal standard solution, sample) K: $1.0\text{mV} \leq \text{FIV}(2) - \text{FIV}(4)$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
61-3	Warning	ISE Noise Error	The difference between the maximum and minimum electromotive forces of internal standard solution (FIV) is within the following range at 3 of 5 measuring points (internal standard solution, sample) Cl: $0.8\text{mV} \leq \text{FIV}(2) - \text{FIV}(4)$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
62-1	Warning	ISE Prepare Error	1 Upon calibration, the slope level is within the following range. 2 The electrode response is degraded (the crossover rate A is as follows) Na: $(1)32.0\text{mV} \leq \text{SLOPE} \leq 37.0\text{mV}$ or $8.1\text{mV} \leq \text{SLOPE}$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
62-2	Warning	ISE Prepare Error	1 Upon calibration, the slope level is within the following range. 2 The electrode response is degraded (the crossover rate A is as follows) K: $(1)32.0\text{mV} \leq \text{SLOPE} \leq 37.0\text{mV}$ or $68.1\text{mV} \leq \text{SLOPE}$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not easy for the user to remove the cause of trouble, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
62-3	Warning	ISE Slope Error	1 Upon calibration, the slope level is within the following range. 2 The electrode response is degraded (the crossover rate A is as follows) Cl: $(1)30.0\text{mV} \leq \text{SLOPE} \leq 25.0\text{mV}$ or $68.1\text{mV} \leq \text{SLOPE}$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not a single malfunction, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".

Code	Level	Name	Describe	Countermeasure
63-1	Warning	ISE Slope Error	1 Upon calibration, the slope level is within the following range. 2 The electrode response is degraded (the crossover rate A is as follows) Na: $(1)\text{SLOPE} < 32.0\text{mV}$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not a single malfunction, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
63-2	Warning	ISE Slope Error	1 Upon calibration, the slope level is within the following range. 2 The electrode response is degraded (the crossover rate A is as follows) K: $(1)\text{SLOPE} < 32.0\text{mV}$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not a single malfunction, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
63-3	Warning	ISE Slope Error	1 Upon calibration, the slope level is within the following range. 2 The electrode response is degraded (the crossover rate A is as follows) Cl: $(1)\text{SLOPE} < 25.0\text{mV}$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not a single malfunction, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
64-1	Warning	ISE concentration error	The concentration of internal standard solution C(1S) is within the following range: Na: $\text{C}(1\text{S}) < 120.0\text{mmol/l}$ or $160.0\text{mmol/l} < \text{C}(1\text{S})$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not a single malfunction, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
64-2	Warning	ISE concentration error	The concentration of internal standard solution C(1S) is within the following range: K: $\text{C}(1\text{S}) < 3.0\text{mmol/l}$ or $7.0\text{mmol/l} < \text{C}(1\text{S})$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not a single malfunction, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
64-3	Warning	ISE concentration error	The concentration of internal standard solution C(1S) is within the following range: Cl: $\text{C}(1\text{S}) < 80.0\text{mmol/l}$ or $120.0\text{mmol/l} < \text{C}(1\text{S})$	1 Carry out "ISE Check" on the Mechanisms Check screen. 2 If it is not a single malfunction, contact the service department. 3 For details refer to "electrolyte device maintenance" in "User manual".
65-1	Warning	ISE require calibration	Carry out ISE calibration after ISE maintenance (ISE running, rinsing all)	Carry out ISE calibration.

Code	Level	Name	Describe	Countermeasure
143-10	Warning	R1 liquid level detect failure	when adding detergent	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-11	Warning	R2 liquid level detect failure	when adding detergent	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-12	Warning	Incubation bath liquid level detect failure	Incubation bath liquid level detect failure	1. Check if liquid level sensor of incubation bath is clean and if there are water in incubation bath 2. Restart the instrument 3. Various malfunctions or trouble can not be restored. Contact the service department.
143-13	Warning	Incubation bath liquid path error	Add water overtime in incubation bath	1. Check if liquid level sensor of incubation bath is clean and if there are water in incubation bath 2. Restart the instrument 3. Various malfunctions or trouble can not be restored. Contact the service department.
143-14	Warning	Barcode scanning overtime	Barcode scanning overtime	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-15	Warning	Degas overtime	Degas overtime	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-16	Stop	Reaction initialize failure	Reaction initialize failure	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-17	Stop	Reaction disk stop failure	Reaction disk stop failure	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-18	Stop	Sample probe block	Sample probe is block	To remove the block please according to 12.4.1 in "user manual"

Code	Level	Name	Describe	Countermeasure
143-1	Stop	Time sync failure	Time sync instruction transmitting failure	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-2	Warning	Add water overtime	Water tank liquid path system error, add water overtime	1. Check the supplying pipe whether an exists in check the height of water tank whether it is more than 1.5 meters (distance from ground) if it will be used 2. Check if water supplier, pipeline and filter is normal 3. Various malfunctions or trouble can not be restored. Contact the service department.
143-3	Warning	AD board reset failure	Error occurs when AD board reset	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-4	Warning	Reaction disk reset failure	Error occurs when reaction disk resetting and reaction reset failure	1. Check if there is waterdrop on the code wheel 2. Call up the system maintenance menu and carry out the check program "reset" 3. Various malfunctions or trouble can not be restored. Contact the service department.
143-5	Warning	Sample disk reset failure	Error occurs when sample disk resetting and sample disk reset failure	1. Check if there is substance in sample container 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-6	Warning	R1 disk reset failure	Error occurs when R1 disk resetting and R1 disk reset failure	1. Check if there is substance in reagent R1 container 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-7	Warning	R2 disk reset failure	Error occurs when R2 disk resetting and R2 disk reset failure	1. Check if there is substance in reagent R2 container 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-8	Warning	Water discharge failure	Incubation bath discharge failure	1. Check if incubation filter and diluent waste pipeline is clog 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-9	Warning	Adding detergent overtime	Detergent is not added completely in specified time by reagent probe 1 and reagent 2	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
143-31	Warning	Vacuum pump error	Vacuum pump negative pressure too low	Please contact the maintenance department.
143-32	Warning	Barcode scanning overture	Barcode scanning overture when testing	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-33	Stop	ISE resetting failure	Error occur in ISE resetting process, and the resetting do not accomplish.	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-34	Warning	ISE check operation overture	ISE check operation overture.	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-35	Warning	ISE pipeline running overture	ISE pipeline running overture.	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-36	Warning	ISE testing error	Error occur when ISE testing, ISE testing stop.	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-37	Stop	Gear wheel pump error	Gear wheel pump pressure too low	Please contact the maintenance department.
143-38	Warning	R1 disk cover open	R1 disk cover is open.	
143-39	Warning	R1 disk cover close	R1 disk cover is close, instrument will automatically carry out reagent horizontal scanning.	
143-40	Warning	R2 disk cover open	R2 disk cover is open.	
143-41	Warning	R2 disk cover close	R2 disk cover is close, instrument will automatically carry out reagent horizontal scanning.	

Code	Level	Name	Describe	Countermeasure
143-19	S. Stop	Add sample failure	Last add sample failure	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-20	S. Stop	Add R1 failure	Last add R1 failure.	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-21	S. Stop	Add R2 failure	Last add R2 failure	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-22	S. Stop	Add R3 failure	Last add R3 failure	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-23	S. Stop	Add R4 failure	Last add R4 failure	
143-24	S. Stop	R1 stirring failure	Last R1 stirrer failure	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-25	S. Stop	R2 stirring failure	Last R2 stirrer failure	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-26	S. Stop	R3 stirring failure	Last R3 stirrer failure	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-27	S. Stop	R4 stirring failure	Last R4 stirrer failure	After testing upon added sample, call up system maintenance menu, carry out "reset"
143-28	Warning	Waste bottle full	The waste bottle is full	Empty the waste bottle.
143-29	Stop	Floot switch error	Error on float switch, high liquid level float switch could detect the signal but the low liquid level float switch cannot detect the signal.	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.
143-30	Warning	Reagent horizontal scanning overture.	Reagent horizontal scanning overture.	1. Call up the system maintenance menu and carry out the check program "reset" 2. Various malfunctions or trouble can not be restored. Contact the service department.

Code	Level	Name	Describe	Countermeasure
143-42	Stop	Mapping information failure	Transmit reagent mapping information failure and add reagent may failure too	1 Call up the system maintenance menu and carry out the check program "reset". 2 Various malfunctions or trouble can not be restored. Contact the service department
143-43	Warning	Cooling overtime	When water tank temperature exceed 36.5 degree, induce cold water into water tank. and water tank temperature cannot reduce to 35.5 degree in 1 minute.	1 Room temperature is too high or water supply equipment error occur. 2 Please contact the maintenance department
143-44	Stop	Instrument module error	Error occurs between instrument modules.	1 Power off and restart the instrument. 2 Various malfunctions or trouble cannot be restored. Contact the service department.
143-45	Stop	Contaminated cup occur	Continually 5 contaminated cup occur	Carry out cell blank test to judge cup states. If cell blank value is abnormal, please change reaction curve/cup. If cell blank value is normal please contact maintenance department.
143-46	Stop	AD data error	AD data error	In resetting, if trouble cannot remove, or other malfunctions exist, please contact maintenance department.
143-47	Warning	Test ISE error	ISE test internal standard liquid failure	In resetting, if trouble cannot remove, or other malfunctions exist, please contact maintenance staff.
143-48	Warning	Read edition number overtime	Read edition number overtime	1 Call up the system maintenance menu and carry out the check program "reset". 2 Various malfunctions or trouble can not be restored. Contact the service department.
143-49	Warning	Illegal test item	Illegal test number and reagent volume	Please contact the service department
144-1	Warning	Cooling system abnormal	Cooling time abnormal	1 Check if reagent disk lid is covered. Check if environmental temperature is conform to instrument requirement. 2 If trouble cannot be removed or other malfunction exist, please contact the maintenance department
144-2	Warning	Cooling system abnormal	Cooling circuit abnormal	Please contact the maintenance department.
144-3	Warning	Cooling system abnormal	Cooling chip abnormal	Please contact the maintenance department.

Code	Level	Name	Describe	Countermeasure
144-4	Warning	Cooling system abnormal	Cooling liquid level abnormal	Please add water in cooling system water take according to 12.4.8 in the "user manual".
144-5	Warning	Cooling system abnormal	Cooling status abnormal	Please contact the maintenance department
144-6	Warning	Ambience temperature over	Ambience temperature is over 15-32 degree	Please contact the service department.
144-7	Warning	Cooling status cannot be continued	(1). Cooling system is not connected. (2). Cooling system malfunction.	Please contact maintenance staff
144-8	Warning	Cooling fan malfunction	Cooling fan running abnormal	Please contact maintenance staff
145-1	Warning	Path 1 cooling chip error	Path 1 cooling circuit is less than 5A	Please contact the maintenance department.
145-2	Warning	Path 2 cooling chip error	Path 2 cooling circuit is less than 5A	Please contact the maintenance department.
145-3	Warning	Path 3 cooling chip error	Path 3 cooling circuit is less than 5A	Please contact the maintenance department.
145-4	Warning	Path 4 cooling chip error	Path 4 cooling circuit is less than 5A	Please contact the maintenance department.
145-5	Warning	Path 5 cooling chip error	Path 5 cooling circuit is less than 5A	Please contact the maintenance department.
145-6	Warning	Path 6 cooling chip error	Path 6 cooling circuit is less than 5A	Please contact the maintenance department.
146-1	Warning	Path 1 AD collector error	Path 1 AD collector value has exceed normal range.	Please contact the maintenance department.
146-2	Warning	Path 2 AD collector error	Path 2 AD collector test value has exceed normal range.	Please contact the maintenance department.

Code	Level	Name	Describe	Countermeasure
146-3	Warning	Path 3 AD collector error	Path 3 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-4	Warning	Path 4 AD collector error	Path 4 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-5	Warning	Path 5 AD collector error	Path 5 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-6	Warning	Path 6 AD collector error	Path 6 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-7	Warning	Path 7 AD collector error	Path 7 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-8	Warning	Path 8 AD collector error	Path 8 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-9	Warning	Path 9 AD collector error	Path 9 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-10	Warning	Path 10 AD collector error	Path 10 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-11	Warning	Path 11 AD collector error	Path 11 AD collector test value has exceed normal range.	Please contact the maintenance department.
146-12	Warning	Path 12 AD collector error	Path 12 AD collector test value has exceed normal range.	Please contact the maintenance department.

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Item name	Assay mode	Time	Measurable wavelength (nm)		Photo metric point	Sample volume	Reagent volume				Calibration method	Calibration point	Span point	Reaction direction	Deviation rate check	Discrete check	Sensitivity check	Blank horizontal check
			Main wavelength	Sub wavelength			R1	R2	R3	R4								
ALT	Rate A	10	340	405	21-31	30	240	0	60	0	2-point linearity	2	2	Negative	3.3	0.05	0	0.8-2.5
AST	Rate A	10	340	405	21-31	30	240	0	60	0	2-point linearity	2	2	Negative	3.3	0.05	0	0.8-2.5
ALP	Rate A	10	405	505	21-26	4	200	0	50	0	2-point linearity	2	2	Positive	3.3	0.05	0	0-1.0
GGT	Rate A	10	405	505	21-31	25	200	0	50	0	2-point linearity	2	2	Positive	3.3	0.05	0	0-1.2
TBA	2point rate	10	405	505	20-26	4	270	0	90	0	2-point linearity	2	2	Positive	3.3	0.05	0	0-1.2
TBIL	1point assay	10	546	660	31	20	0	0	200	0	2-point linearity	2	2	Positive	3.3	0.05	0	0-0.5
DBIL	1 point assay	10	546	660	31	20	0	0	200	0	2-point linearity	2	2	Positive	3.3	0.05	0	-0.1-0.5
TP	1 point assay	10	546	700	31	5	250	0	0	0	2-point linearity	2	2	Positive	3.3	0.05	0	-0.5-0.5
ALB	1 point assay	5	500	700	10	2	300	0	0	0	2-point linearity	2	2	Positive	3.3	0.05	0	0-0.5
LAP	Rate A	10	405	505	21-31	15	240	0	60	0	2-point linearity	2	2	Positive	3.3	0.05	0	0-1.2
SCHE	Rate A	10	405	505	20-26	3	240	0	60	0	2-point linearity	2	2	Positive	3.3	0.05	0	0-1.0
ICDH	Rate A	10	340	405	21-31	15	240	0	60	0	2-point linearity	2	2	Positive	3.3	0.05	0	-0.1-1.0

Item name	Assay mode	Time	Measurable wavelength (nm)		Photo point	Sample volume	Reagent volume				Calibration method	Calibration point	Span point	Reaction direction	Deviation rate check	Discrete check	Sensitivity check	Blank horizontal check
			Main wavelength	Sub wavelength			R1	R2	R3	R4								
CK-B	Rate A	10	340	405	25-31	10	200	0	50	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
CK	Rate A	10	340	405	25-31	5	200	0	50	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
HBDM	Rate A	10	340	405	21-31	6	240	0	60	0	2-point linearity	2	2	2	Negative	3.3	0.05	0
LDH	Rate A	10	340	405	21-31	5	240	0	60	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
UA	2point assay	10	546	700	16-31	4	200	0	50	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
P	1point assay	5	340	405	10	4	200	0	0	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
Ca-CPC	2point assay	10	570	660	15-31	6	150	0	150	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
CI	1point assay	5	505	660	15	3	300	0	0	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
Ca-ARS	1point assay	5	660	750	15	3	300	0	0	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
Mg	1point assay	4	546	750	10	3	300	0	0	0	2-point linearity	2	2	2	Positive	3.3	0.05	0

Item name	Assay mode	Time	Measurable wavelength (nm)		Photo point	Sample volume	Reagent volume				Calibration method	Calibration point	Span point	Reaction direction	Deviation rate check	Discrete check	Sensitivity check	Blank horizontal check
			Main wavelength	Sub wavelength			R1	R2	R3	R4								
GLDH	Rate A	10	340	405	21-31	15	240	0	60	0	2-point linearity	2	2	2	Negative	3.3	0.05	0
AMY	Rate A	10	405	505	21-26	5	160	0	40	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
BUN	2point Rate	10	340	405	20-26	3	240	0	60	0	2-point linearity	2	2	2	Negative	3.3	0.05	0
CRE-E	2point Rate	10	546	700	10-31	8	240	0	60	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
CRE	2point rate	10	505	660	18-23	20	150	0	150	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
GLU	2point assay	15	505	660	16-49	2	240	0	60	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
GLU (HK)	2point assay	10	340	405	15-31	2	240	0	60	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
FMN	2point rate	10	546	700	10-31	15	300	0	0	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
TC	2point assay	10	505	660	10-31	3	240	0	60	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
TG	2point assay	10	546	700	10-31	3	240	0	60	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
HDL-C	2point assay	10	546	660	16-31	3	225	0	75	0	2-point linearity	2	2	2	Positive	3.3	0.05	0
LDL-C	2point assay	10	546	660	10-31	4	300	0	100	0	2-point linearity	2	2	2	Positive	3.3	0.05	0

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