MEDRAD® Intego Administration Sets Наборы для введения MEDRAD® Intego

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©2017 Bayer. All Rights Reserved. Reproduction of this manual is strictly prohibited without express written consent of Bayer Medical Care Inc.

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Bayer, the Bayer Cross, MEDRAD, Intego, and MEDRAD Intego may be registered trademarks of Bayer in the US and other countries.

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UNITED STATES

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

EUROPE

Bayer Medical Care B.V.
Horsterweg 24
6199 AC Maastricht Airport
The Netherlands
Phone: +31 (0) 43-3585601
Fax: +31 (0) 43-3656598

© 2007-2018 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, and MEDRAD Intego are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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document completely before installing or using the MEDRAD Intego PET Infusion System Administration Sets. Understanding this information will assist the user in safely operating the MEDRAD Intego System.

Intended Use: The MEDRAD Intego PET Infusion System is intended to deliver accurate doses of ^{18}F Fluorodeoxyglucose (FDG) and ^{18}F Sodium Fluoride (NaF) radiopharmaceuticals and commonly used flushing solutions to patients during molecular imaging (nuclear medicine) diagnostic procedures. The MEDRAD Intego PET Infusion System is also intended to provide effective radiation shielding to medical personnel from Fluorine-18 (^{18}F) radiation exposure during nuclear medicine diagnostic procedures.

Indications for Use: The MEDRAD Intego PET Infusion System is indicated for the administration of Fluorodeoxyglucose (FDG), ^{18}F Sodium Fluoride (NaF), and commonly used flushing solutions to patients during

the delivery of NaF with software version CR10.1 and later. Please contact the Bayer representative for additional information regarding administration of ^{18}F NaF.

Restricted Sale: The United States Food and Drug Administration (FDA) restricts sale of this device to physicians or those with written authorization from a physician.

Warnings

For devices labeled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

This product contains DEHP, a chemical known to the state of California to cause birth defects or other reproductive harm. Please visit www.medrad.com/DEHP for more information regarding DEHP in our products.

Biological contamination can result from reusing the Patient Administration Set or failure to follow sterile technique. Do NOT reuse or attempt to resterilize. Properly discard disposable items after use. If there is any possibility that contamination may have occurred during set-up or use, disassemble and set-up a new sterile product.

Severe operator injury can occur if care is not taken when opening or closing the lids over the shielded chamber or the Patient Administration Set compartment. Be certain that all fingers are clear before opening or closing the shielded chamber lid or the Patient Administration Set compartment lid.

Only the fluid path of the set is sterile and non pyrogenic. Do not use in a sterile or aseptic area without proper precautions.

Air embolism or improper injection dosage may occur if saline line is damaged. Make certain not to pinch the saline tube within the saline pump as this may damage the administration set.

Biological contamination may occur if aseptic techniques are not followed. Follow facility guidelines to disinfect the vial septum.

Radiation Exposure Hazard: Injector may contain radioactive materials. Do not reach under cart while system contains activity. Do not remove access panels while system contains activity.

Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Puncture hazard: The SAS contains sharps. To prevent injury, take care when handling the SAS to control the sharps. Dispose of the sharps per individual facility guidelines.

DO NOT use the Source Administration Set with more than one multi-dose vial. Using the Source Administration Set with more than one multi-dose vial could result in injury to patient or operator.

Radiation exposure hazard. Using the lid of the Shielded Chamber as a work space may expose the operator to excess radiation. Do not use the lid to the Shielded Chamber as a work space while there are radioactive materials present within the system. Using the lid as a work space or placing items on the lid while it is open may damage the lid, hampering its ability to provide an adequate radiation shield.

Radiation exposure hazard. Radiopharmaceuticals emit gamma radiation. Operator should exercise extreme caution when handling radiopharmaceuticals. Doses may exceed maximum permissible dose if appropriate precautions are not properly exercised. Follow all individual facility guidelines for the handling and disposal of radioactive materials. DO NOT open the Shielded Chamber without first taking radiation precautions as outlined in facility procedures.

Warnings

DO NOT use any radiopharmaceutical products other than ^{18}F Fluorodeoxyglucose (FDG) or ^{18}F Sodium Fluoride (NaF) with the MEDRAD Intego System. Use of radiopharmaceuticals other than ^{18}F FDG or ^{18}F NaF may damage the MEDRAD Intego System and may cause inaccurate dose measurement and/or excessive radioactive exposure.

Improper injection dosage or radiation spillage may result if the Shielded Chamber lid is closed onto the tubing set. Ensure that the tubing set and all administration set components are properly installed before closing the Shielded Chamber lid.

Patient infection may result from the use of non-sterile components. Maintain sterility of all disposable components using aseptic techniques.

INSTALLING THE SOURCE ADMINISTRATION SET (SAS)

Always install a new SAS when installing a new multi-dose vial. Install the multi-dose vial and the SAS in the Shielded Chamber using the following procedure:

1. Release the latch securing the Chamber Lid and slide the lid towards the front of the cart to expose the Shielded Chamber.
2. Open the disposable tray and remove the SAS from the packaging.
3. Place the saline bag onto the hanger on the rear of the cart near the handle and insert spike at the end of the SAS into the saline bag.
4. Insert the Saline Tube into the Saline Pump so that the connector is between the pump and the saline bag, with the connector close to the pump.
5. Insert the tubing into the hold-down clip and route the Saline Tube along the SAS track, starting at the Saline Pump and ending at the Vial Shield Compartment where the Vial Shield will be installed.

NOTE: Make certain that the Saline Tube is centered on the pump rollers. The tubing must not fall below the rollers. When closing the pump door, make certain that the Saline Tube is centered in the pump door so that the tubing is not pinched.

6. Insert the Sharps Cartridge into the clip on the shield lid. The tube exiting the Sharps Cartridge should be towards the RP Pump.

NOTE: Do not remove the protective sheath from the needles of the Sharps Cartridge.

7. Insert Confluence into the SAS Confluence holder.
8. Verify dose calibrator is empty before installing coil.
9. Insert the SAS Coil into the Dose Calibrator so that the loose ends of the tubing are towards the top.

NOTE: Be sure that the SAS Coil is placed at the bottom of the Dose Calibrator.

10. Unroll the Waste Bag and insert it into the shielded Waste Storage.
11. Insert the T-connector into the retainer.

12. Install the section of tubing between the T-connector and the Waste Bag into the waste pinch valve.

NOTE: The pinch valve is opened by pushing on the button on the top of the valve. To install the tubing, slide the tubing into the valve while pushing down on the button.

13. Install the pinch valve cover on the waste pinch valve.
14. Install the section of tubing between the T-connector and the Swabbable Valve into the PAS pinch valve.

NOTE: The pinch valve is opened by pushing on the button on the top of the valve. To install the tubing, slide the tubing into the valve while pushing down on the button.

15. Install the pinch valve cover on the PAS pinch valve.
16. Insert the free end of the SAS (exiting the PAS pinch valve) into the Air Detector.

17. Using the Vial Shield Compartment handle, lift the shield lid straight up and then rotate the lid counter-clockwise.
18. The shield lid is hinged. Lift the lid to expose the Sharps Cartridge Clip.

19. Lower the shield lid over the Sharps Cartridge.
20. Insert the Vial Shield containing a full vial of the radiopharmaceutical indicated for use into the Vial Shield Compartment. DO NOT remove the vial access cap of the Vial Shield. Leave the carrying handle on the Vial Shield.

NOTE: The Vial Shield is keyed so that it can only be properly inserted into the Vial Shield Compartment with the Vial Shield Compartment posts resting in the keyed areas located next to the flat surfaces of the Vial Shield.

21. Remove the cap from the Vial Shield by pushing down on the carrying handle, rotating the vial access cap counter-clockwise, and then lifting the cap straight off.

22. Follow aseptic techniques as outlined in individual facility guidelines to disinfect the vial septum.

23. Insert the sharps into the multi-dose vial by pivoting the shield lid back over the well and pressing the needles straight down into the vial.

24. Insert the tubing from the Sharps Cartridge into the RP pump. Ensure that the tube collar rests immediately outside the RP pump, on the side closest to the Vial Shield Compartment. Close the pump lid.

25. Verify that all of the SAS is lying flat within the channel, that the Vial Shield Compartment handle is down and that the door of the RP pump is closed with the tubing completely installed. Close the chamber lid.

NOTE: Follow facility procedures to ensure that the radiopharmaceutical indicated for use is properly secured within the MEDRAD Intego PET Infusion System.

NOTE: The MEDRAD Intego PET Infusion System is designed to be used with standard locks. Bayer does not supply a lock with the system.

26. Enter multi-dose vial assay information by following the procedure outlined in the Entering Multi-dose Vial Assay Information section of the MEDRAD Intego Operation Manual.

INSTALLING THE PATIENT ADMINISTRATION SET (PAS)

Always install a new PAS for each patient. Install the PAS using the following procedure:

1. Insert the PAS into the PAS area of the Shielded

2. Clean the Swabbable Valve at the distal end of the SAS following aseptic techniques described in individual facility guidelines.
3. Attach the PAS to the Swabbable Valve at the distal end of the SAS.
4. Depending on which side of the patient the MEDRAD Intego PET Infusion System is located, route the PAS out either the left or right side of the shielding.
5. Ensure that the PAS will not be pinched, and close the shield cover.

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