A. GIACOMINI 7th June 2017
Regulatory Affairs, Association Diderot/Giacomini

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




EN

COMPOSITION

SURGIDERM® 24XP	Hyaluronic acid gel: Phosphate buffer pH 7.2 q.s. One syringe contains 0.8 mL.	24 mg 1 mL
SURGIDERM® 30	Hyaluronic acid gel: Phosphate buffer pH 7.2 q.s. One syringe contains 0.8 mL.	24 mg 1 mL
SURGIDERM® 30XP	Hyaluronic acid gel: Phosphate buffer pH 7.2 q.s. One syringe contains 0.8 mL.	24 mg 1 mL

DESCRIPTION

SURGIDERM® 24XP / 30 and 30XP are sterile pyrogen-free physiological solutions of cross-linked hyaluronic acid which is not of animal origin. The gel is presented in a graduated, pre-filled, disposable syringe. Each box contains an instruction leaflet and a set of labels in order to ensure traceability and:

- For **SURGIDERM® 24XP**, two 0.8mL syringes and four single-use 30G1/2" sterile needles;
- For **SURGIDERM® 30**, two 0.8mL syringes and four single-use 27G1/2" sterile needles;
- For **SURGIDERM® 30XP**, two 0.8mL syringes and four single-use 27G1/2" sterile needles.

STERILISATION

The contents of the **SURGIDERM® 24XP / 30 and 30XP** syringes is sterilised by moist heat. The 30G1/2" and 27G1/2" needles are sterilised by radiation.

INDICATIONS

SURGIDERM® 24XP	Injectable implant used for filling any medium depression of the skin via mid-dermis injection, as well as for definition and pouting of the lips.
SURGIDERM® 30	Injectable implant used for filling any deep depression of the skin via deep dermis injection, as well as for lip enhancement and cheekbone augmentation.
SURGIDERM® 30XP	Injectable implant used for filling mid and/or deep depression of the skin via mid and/or deep dermis injection, as well as for volume increase, definition and pouting of the lips.

CONTRA-INDICATIONS

SURGIDERM® 24XP	Do not inject in the eyelids. The application in the bags under the eyes is reserved to specialists specifically trained in this technique and having a sound knowledge of the physiology for this particular area.
SURGIDERM® 30	Do not inject in the periorbital area (eyelids, bags under the eyes, eye wrinkles) and glabellar region.
SURGIDERM® 30XP	Do not inject in the periorbital area (eyelids, eye wrinkles) and glabellar region (forehead). The application in the bags under the eyes is reserved to specialists specifically trained in this technique and having a sound knowledge of the physiology for this particular area.

In addition,

- Do not inject into the blood vessels (intravascular). Intravascular injection may lead to embolization, occlusion of the vessels, ischemia or infarction.
- Do not overcorrect.
- **SURGIDERM® 24XP / 30 and 30XP** must not be used in:
 - Patients who tend to develop hypertrophic scarring;
 - Patients with known hypersensitivity to hyaluronic acid and/or to gram positive bacterial proteins as hyaluronic acid is produced by *Streptococcus* type bacteria;
 - Women who are pregnant or breastfeeding;
 - Children.
- **SURGIDERM® 24XP / 30 and 30XP** must not be used in areas presenting cutaneous inflammatory and/or infectious processes (acne, herpes, etc.).
- **SURGIDERM® 24XP / 30 and 30XP** should not be used simultaneously with laser treatment, deep chemical peels or dermabrasion. For surface peels, it is recommended not to inject the product if the inflammatory reaction generated is significant.

PRECAUTIONS FOR USE

- **SURGIDERM® 24XP / 30 and 30XP** are indicated only for intra-dermal injections and in lips.
- **SURGIDERM® 24XP / 30 and 30XP** are not intended for use in breast augmentation/reconstruction.
- As a matter of general principle, injection of a medical device is associated with a risk of infection. Standard precautions associated with injectable materials shall be followed.
- There is no available clinical data about injection of **SURGIDERM® 24XP / 30 and 30XP** into an area which has already been treated with a non-ALLERGAN dermal filler.
- It is recommended not to inject into a site which has been treated with a permanent implant.
- No clinical data is available regarding the efficiency and tolerance of **SURGIDERM® 24XP / 30 and 30XP** injections in patients having a history of, or currently suffering from, autoimmune disease or autoimmune deficiency or being under immunosuppressive therapy. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the disease and its corresponding treatment, and shall also ensure the specific monitoring of these patients. In particular, it is recommended that these patients undergo a preliminary skin testing for hypersensitivity, and to refrain from injecting the product if the disease is active.
- There is no available clinical data concerning the tolerance of the **SURGIDERM® 24XP / 30 and 30XP** injection in patients presenting a history of severe and/or multiple allergies. The medical practitioner shall therefore decide on the indication on a case-by-case basis, according to the nature of the allergy, and shall also ensure the specific monitoring of these at-risk patients. In particular, the decision may be taken to propose a skin testing for hypersensitivity or suitable preventive treatment prior to any injection. In case of history of anaphylactic shock, it is recommended not to inject the product.
- Patients showing a history of streptococcal disease (recurrent sore throats, acute rheumatic fever) shall be subjected to a skin testing for hypersensitivity before any injection is administered. In the event of acute rheumatic fever with heart complications, it is recommended not to inject the product.
- Patients on anti-coagulation medication or using substances that can prolong bleeding (warfarin, acetylsalicylic acid, nonsteroidal anti-inflammatory drugs, or other substances known to increase coagulation time such as herbal supplements with ginkgo or ginkgo biloba, etc.) must be warned of the potential increased risks of bleeding and haematomas during injection.
- There is no data available regarding the safety of injecting greater amount than 20 mL of ALLERGAN dermal fillers per 60 kg (130 lbs) body mass per year.
- Please recommend that the patient not use any makeup during the 12 hours following the injection treatment and that any extended exposure to the sun, UV rays and temperatures below 0°C be avoided, as well as any saunas or hammam sessions during the two weeks following the injection treatment.
- The composition of this product is compatible with fields used for magnetic resonance imaging.

INCOMPATIBILITIES

Hyaluronic acid is known to be incompatible with quaternary ammonium salts such as benzalkonium chloride. **SURGIDERM® 24XP / 30 and 30XP** should therefore never be placed in contact with these substances or with medical-surgical instrumentation which has been treated with this type of substance.

UNDESIRABLE EFFECTS

The patients must be informed that there are potential side effects associated with implantation of this product, which may occur immediately or may be delayed. These include, but are not limited to:

- Inflammatory reactions (redness, oedema, erythema, etc.) which may be associated with itching and/or pain on pressure and/or paresthesia, occurring after the injection. These reactions may last for a week.
- Haematomas.
- Induration or nodules at the injection site.
- Staining or discoloration of the injection site might be observed, especially when HA dermal filler is injected too superficially and/or in thin skin (Tyndall effect).
- Poor effect or weak filling effect.
- Rare but serious adverse events associated with intravascular injection of dermal fillers in the face and tissue compression have been reported and include temporary or permanent vision impairment, blindness, cerebral ischemia or cerebral hemorrhage, leading to stroke, skin necrosis and damage to underlying structures. Immediately stop the injection if a patient exhibits any of the following symptoms, including changes in the vision, signs of stroke, blanching of the skin or unusual pain during or shortly after the procedure. Patients should receive prompt medical attention and possibly evaluation by an appropriate medical practitioner specialist should an intravascular injection occur. Abscesses, granuloma and immediate or delayed hypersensitivity after hyaluronic acid injections have also been reported. It is therefore advisable to take these potential risks into account.
- Patients must report inflammatory reactions which persist for more than one week, or any other side effect which develops, to their medical practitioner as soon as possible. The medical practitioner should use an appropriate treatment.
- Any other undesirable side effects associated with injection of **SURGIDERM® 24XP / 30** and **30XP** must be reported to the distributor and/or to the manufacturer.

METHOD OF USE - POSOLOGY

- This product is designed to be injected into the dermis by an authorized medical practitioner in accordance with local applicable regulation. In order to minimize the risks of potential complications and as precision is essential to a successful treatment, the product should be only used by medical practitioners who have appropriate training, experience and who are knowledgeable about the anatomy at and around the site of injection. The nappage technique can also be used with **SURGIDERM® 30XP**.
- **SURGIDERM® 24XP / 30** and **30XP** are to be used as supplied. Modification or use of the product outside the Directions for Use may adversely impact the sterility, homogeneity and performance of the product and it can therefore no longer be assured.
- Prior to treatment, medical practitioners shall inform their patients about the product's indications, contra-indications, incompatibilities and potential undesirable effects/risks associated with dermal fillers injection and ensure that patients are aware of signs and symptoms of potential complications.
- The area to be treated should be disinfected thoroughly prior to the injection.
- Remove tip cap by pulling it straight off the syringe as shown in fig. 1. Then firmly push the needle provided in the box (fig. 2) into the syringe, screwing it gently clockwise. Twist once more until it is fully locked and has the needle cap in the position shown in fig. 3. If the needle cap is positioned as shown in fig. 4, it is incorrectly attached.
- Next, remove the protective cap by holding the body of the syringe in one hand, the protective cap in the other, as shown in fig. 5, and pulling the two hands in opposite directions.
- Prior to injecting, depress the plunger rod until the product flows out of the needle. Inject slowly and apply the least amount of pressure necessary.
- If the needle is blocked, do not increase the pressure on the plunger rod. Instead, stop the injection and replace the needle.
- Failure to comply with these precautions could cause a disengagement of the needle and/or product leakage at luer-lock level and/or increase the risk of vascular compromise.
- After needle insertion and before injection, it is recommended to withdraw slightly the plunger to aspirate and verify the needle is not intravascular.
- If immediate blanching occurs at any time during the injection, the injection should be stopped and appropriate action taken such as massaging the area until its return to a normal color.

- The degree and duration of the correction depend on the character of the defect treated, the tissue stress at the implant site, the depth of the implant in the tissue and the injection technique. The amount injected will depend on the areas which are to be corrected based on the experience of the medical practitioner.
- Do not overcorrect as injection of an excessive volume can be at the origin of some side effects such as tissue necrosis and oedema.
- A touch up (for achieving optimal correction) and/or a repeat (for maintaining optimal correction) treatment of **SURGIDERM® 24XP / 30** and **30XP** might be required.
- It is recommended to wait until side effects are resolved (with a minimal interval of 2 weeks) between two injections.
- It is important to massage the area treated after the injection in order to ensure that the substance has been uniformly distributed.

WARNINGS

- Check the expiry date on the product label.
- In the event that the content of a syringe shows signs of separation and/or appears cloudy, do not use the syringe.
- Do not re-use. Sterility of this device cannot be guaranteed if the device is re-used.
- Do not re-sterilise.
- For the needles (CE 0123 TSK Laboratory, Japan):
EC-Representative: Emergo Europe
Molenstraat 15
2513 BH The Hague (NL)
- Used needles must be thrown away in the appropriate containers. Do the same for the syringes. Please consult the current applicable directives to ensure their correct elimination.
- Never try to straighten a bent needle; throw it away and replace it.

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