Instructions for Use

ANKYLOS® Prosthetics devices for Intraoral or Extraoral Welding

Bensheim, September 19th 2019

Dr. Aline Frare Mocruha
Manager Regulatory Affairs

**Dentsply
Sirona**
Implants
Dentsply Implants
Manufacturing GmbH
Rodenbacher Chaussee 4
63457 Hanau, Germany

**Dentsply
Sirona**
Implants

For Instructions for Use
and symbols glossary refer
to
ifu.dentsplysirona.com

For the list of products covered by this IFU please refer to the ANKYLOS® WeldOne Concept brochure, catalogue or contact your local dealer. For representative product design please refer to the brochure ANKYLOS® WeldOne Concept.

Intended use

ANKYLOS® Prosthetics devices for intraoral or extraoral welding are designed to provide the framework for a prosthetic restoration, in order to restore the chewing function of the patient. Prosthetics devices are intended to produce an intraoral or extraoral-welded titanium frame for immediate and delayed loading treatment.

Indications

ANKYLOS® Welding Sleeve for Balance Base Abutment narrow

Fabrication of an intraorally or extraorally welded titanium framework for immediate and late restorations on ANKYLOS® Balance Base Abutments narrow exclusively by use of the WeldOne welding unit and the titanium wires designed for that purpose

ANKYLOS® Temporary Welding Abutment

Fabrication of an intraorally or extraorally welded titanium framework for immediate temporary restorations on ANKYLOS® Implants exclusively by use of the WeldOne welding unit and the titanium wires designed for that purpose

ANKYLOS® Retaining Screw for Temporary Welding Abutment

Fixation of the ANKYLOS® Temporary Welding Abutments on the ANKYLOS® implants

ANKYLOS® Welding Cap for Standard Abutments a or b and ANKYLOS® Welding Cap for Standard Abutments a or b lateral screw

Fabrication of an intraorally or extraorally welded titanium framework for immediate and late restorations on ANKYLOS® Standard C/ Abutments exclusively by use of the WeldOne welding unit and the titanium wires designed for that purpose

ANKYLOS® Parallel Gauge for Standard Abutment a or b

Alignment aid to verify the insertion direction of Standard abutments during abutment placement

ANKYLOS® Welding Cap for SynCone 5°

Fabrication of an intraorally or extraorally welded titanium framework for immediate and late restorations on ANKYLOS® SynCone Abutments 5° exclusively by use of the WeldOne welding unit and the titanium wires designed for that purpose

FRIADENT® Titanium Wires

Fabrication of an intraorally or extraorally welded titanium framework for prosthetic restorations on implants from Dentsply Implants exclusively by use of the WeldOne welding unit and in conjunction with the abutments designed for that purpose

Contraindications

- In the case of patients with pacemakers or other active implantable devices, intraoral welding is not indicated and must not be performed under any circumstances.
- Wearers of pacemakers or other active implantable devices must not perform intraoral and extraoral welding.
- The devices are not suitable for single-tooth restorations.

Warnings

- These instructions for use must be read prior to using ANKYLOS® devices for intraoral or extraoral welding.
- ANKYLOS® devices for intraoral or extraoral welding may only be used for their intended purpose according to the general rules for dental/surgical treatment, occupational safety and accident prevention.
- ANKYLOS® devices for intraoral or extraoral welding must only be used together with the ANKYLOS® implant system and the WeldOne welding unit.
- The welding of these devices intraorally may only be performed by the dentist; dental technicians may also perform welding extraorally.
- The user must be familiar with dental surgery and prosthetics, including diagnostics and preoperative planning and/or laboratory procedures.
- Dentsply Implants Manufacturing GmbH is not liable for damage resulting from use outside the intended use of the product.
- If the indication and use are unclear, treatment should be suspended until these points have been clarified. The following instructions are not adequate to allow inexperienced clinicians to administer professional implant dentistry.

- Intraoral and extraoral welding may only be performed subject to relevant training in order to ensure that the treatment procedure runs smoothly!

CAUTION: ANKYLOS® Balance Base Abutments and ANKYLOS® Standard C/ Abutments are generally not allowed to be ground.

CAUTION: The system has not been evaluated for safety and compatibility in the Magnetic Resonance (MR) environment. The system has not been tested for heating or migration in the MR environment.

When using prosthetic devices, general complications related to prosthetic restoration could arise, such as:

- Devices used in the patient's mouth have been aspirated or swallowed.
- The screw of the abutment has fractured for example due to excessive torque.
- Loosening of the abutments/ superstructure for example as a consequence of inadequate torque being applied on the screw
- Deformation of parts when using non-system devices as well as devices which are not indicated for the intraoral and extraoral welding concept.

Precautions

CAUTION: Wearers of pacemakers or other active implantable devices must not perform intraoral and extraoral welding.

CAUTION: Using flammable anaesthetic gases, nitrous oxide or oxygen, as well as alcohol in the oral cavity of the patient before and during the welding procedure is prohibited.

The following precautions are to be met prior to or during treatment:

- Prior to each procedure it must be ensured that all necessary devices, instruments and materials are complete, available in the required quantities and fully operational.
- All products intended for single use must not be reused. Failure to observe this precaution leads to a risk of infections and misfitting devices. The associated complications can cause fractures through to implant loss!
- Position the patient such that the danger of aspiration of devices is minimised.
- All devices used in the patient's mouth must be secured to prevent aspiration and swallowing.
- Always wear protective clothing for your own safety.
- Make sure that the patient and user are wearing eye protection during the welding procedure.
- Ensure that all the devices used during welding are free of grease and other residue.
- It is imperative to follow the welding parameter settings according to the table in the instructions for use and in the user manual.
- Before the application of resins with suitable tools, like polymerisation sleeves, make sure that no plastic gets between the implant, abutment or welding abutment.
- If, due to unfavourable anatomical conditions, instruments do not fit or cannot be used for other reasons, the course of treatment planned with them must not be continued and alternatives must be sought.
- Ensure that no threads are trapped when placing the SynCone Cap.
- Observe specified torques.
- Ensure adequate extraction when grinding devices.

Adverse Reactions

No adverse reactions are known

Application

NOTE: Detailed instructions on surgery, prosthetics and laboratory methods can be found in the respective ANKYLOS® manuals.

Follow the instructions supplied when using the ANKYLOS® Balance Base Abutment System, the ANKYLOS® Standard Abutment System and the ANKYLOS® SynCone 5° devices.

FRIADENT® Titanium wires

Select the thickness of the FRIADENT® titanium wire according to the application:

- Ø 1.2 mm for additional retention
- Ø 1.5 mm for inter-implant distances of up to 8 mm
- Ø 2.0 mm for inter-implant distances of between 8 mm and 15 mm

 **Manufacturer:**
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.

Not all products are available in all countries.

Status of information:
IFU 2024-en_RU Rev. 003 A 2019-09

**Dentsply
Sirona**
Implants

Instructions for Use

Shorten the titanium wire to the required length. The titanium wires should be bent using standard dental bending instruments so that these lie against the welding sleeves without tension or gaps. Afterwards, the welding procedure is carried out.

Welding procedure

In order to obtain an appropriate welding joint with the welding unit provided by Dentsply Implants, ANKYLOS® devices for intraoral or extraoral welding have to be used.

Observe the following instructions before starting the welding procedure: Shorten the titanium wire for intraoral or extraoral welding accordingly and place on the ANKYLOS® welding devices to be welded without tension or gaps.

Align the titanium wire in a vestibular-oral direction and adjust to the prosthetic conditions.

IMPORTANT: Follow the correct welding sequence!

Start welding the ANKYLOS® welding devices with the titanium wire at the most accessible position. Continue welding on the adjacent abutment.

IMPORTANT: It is imperative to follow the welding parameter settings according to the table in the instructions for use and in the user manual.

The welding parameters of the welding unit are to be set as follows based on the diameter of the titanium wire used:

FRIADENT Titanium wire Diameter	Welding parameter setting on welding unit
1.2 mm	Low = 60 %, single pulse
1.5 mm	Medium = 75 %, double pulse
2.0 mm	High = 100 %, double pulse

Follow these important instructions during the welding procedure:

After welding, hold the welding clamp in position for about 3 seconds or, optionally, spray the welding area with water during welding. Carry out a visual inspection of the weld seam after each welding procedure.

IMPORTANT: Should the titanium wire separate from the adjacent welding abutments during welding and no longer sit without tension or gaps remove the previously welded framework and bend the titanium wire so that it is joined again.

NOTE: Parallel Gauges are not suitable for welding.

Follow these important instructions after completion of the titanium framework:

- Check the titanium framework in the patient's mouth again to make sure that it is free of tension (e.g. Sheffield test) before incorporating it into a plastic prosthesis. If the titanium framework does not sit without tension or gaps, it must be detached and prepared again.
- The manufactured titanium framework can be furnished with other retention elements if necessary. Shorten the titanium wire accordingly and weld extraorally with the titanium framework.
- Further processing of the titanium framework is carried out according to the generally accepted dental regulations.

IMPORTANT: Before the application of resins with suitable tools, like polymerisation sleeves, make sure that no plastic gets between the implant, abutment or welding abutment. This prevents the ingress of resin at the sulcus/gingiva and the soft tissue is protected.

ANKYLOS® Welding Sleeve for Balance Base Abutment narrow

CAUTION: ANKYLOS® Balance Base Abutments and ANKYLOS® Standard C/ Abutments are generally not allowed to be ground.

Insert and screw the ANKYLOS® Balance Base Abutments narrow or the ANKYLOS® Balance Base Abutments angled into the implant according to the instructions. Position the ANKYLOS® Welding Sleeve for Balance Base Abutment narrow and screw it hand-tight with the Fixation Screw, occlusal M1.6 mm Hex, extra long provided. The Insert for Prosthetic Ratchet Screwdriver 1.0 mm Hex 15 Ncm is used for screw retention.

Shorten the titanium wire accordingly and place on the ANKYLOS® Welding Sleeves to be welded so as to be free of tension and gaps. Position the titanium wire and the ANKYLOS® Welding Sleeve centrally between the two electrodes of the welding unit's welding clamp. The clamping force of the welding clamp secures the devices in position. Support the welding clamp with the palm of your hand during welding and do not open it under any circumstances.

After welding all ANKYLOS® Welding Sleeves with the titanium wire, loosen the Fixation Screw, occlusal M1.6 mm Hex, extra long and remove the framework from the ANKYLOS® Balance Base Abutments. The ANKYLOS® Welding Sleeve can be shortened if necessary. If the Fixation Screw, occlusal M1.6 mm Hex, extra long is too long, either the long or short (blue anodised) Fixation Screw, occlusal M1.6 mm Hex can be used.

Use polymerisation sleeves if necessary. Place the titanium framework in the patient's mouth on the ANKYLOS® Balance Base Abutments and screw hand-tight using the ANKYLOS® Fixation Screws, occlusal M1.6 mm Hex provided at the most accessible points. Furnish the finished plastic prosthesis with cold- or light-cured resin and place on the framework. After the plastic has hardened, undo the fixation screws and remove the prosthesis, then prepare according to the generally accepted dental regulations.

Insertion of the prosthetic restoration using the ANKYLOS® Welding Sleeve for Balance Base Abutment narrow

Check the secure fit of the ANKYLOS® Balance Base Abutments narrow or angled. Then place the prosthetic restoration on the ANKYLOS® Balance Base Abutments narrow or angled. Check the passive fit of the prosthetic restoration again (e.g. Sheffield test) and screw onto the Balance Base Abutments using the Insert for Prosthetic Ratchet Screwdriver 1.0 mm Hex 15 Ncm applying a torque of 15 Ncm and check the occlusion.

ANKYLOS® Temporary Welding Abutment and ANKYLOS® Retaining Screw for Temporary Welding Abutment

Insert the ANKYLOS® Implants as parallel as possible and remove the premounted Placement Head from the implant with the Insert for Prosthetic Ratchet Screwdriver 1.0 mm Hex 15 Ncm. Insert the ANKYLOS® Temporary Welding Abutment in the implant and screw hand-tight with the ANKYLOS® Retaining Screw for Temporary Welding Abutments provided. Insert the polymerisation sleeves. Position the titanium wire and the ANKYLOS® Temporary Welding Abutments centrally between the two electrodes of the welding unit's welding clamp. The clamping force of the welding clamp secures the devices in position. Support the welding clamp with the palm of your hand during welding and do not open it under any circumstances.

Weld all ANKYLOS® Temporary Welding Abutments with the titanium wire, loosen the retaining screws and remove the framework from the implants. The ANKYLOS® Temporary Welding Abutment can be shortened occlusally only up to the height of the screw head.

Insert the prepared titanium framework in the patient's mouth in the implants and screw hand-tight with the retaining screw provided at the most accessible location. Furnish the finished plastic prosthesis with cold- or light-cured resin and place on the framework. After the plastic has hardened, undo the retaining screw and remove the prosthesis from the patient's mouth, then prepare according to the generally accepted dental regulations.

Insertion of the prosthetic restoration using the ANKYLOS® Temporary Welding Abutment

Insert the prosthetic restoration in the implants. Check the passive fit of the prosthetic restoration again (e.g. Sheffield test) and screw into the implant using the Insert for Prosthetic Ratchet Screwdriver 1.0 mm Hex 15 Ncm and check the occlusion.

ANKYLOS® Welding Cap for Standard Abutments a or b and ANKYLOS® Welding Cap for Standard Abutments a or b lateral screw

Inserting Standard Abutment a or b, straight

Select the abutment according to the prosthetic situation. Pick up the abutment with the Insert for Prosthetic Ratchet for a- or b-Abutments key and screw it in the implant with a torque of 25 Ncm. Check the joint angle of insertion for the abutments with the ANKYLOS® Parallel Gauge for Standard Abutment a or b. Following insertion, seal the inner thread of the abutment with a Headless Screw.

Inserting Standard Abutment a or b, angled

Select the abutment according to the prosthetic situation. Position the angled abutment in the implant. Adjust the joint angle of insertion for the abutments with the ANKYLOS® Parallel Gauge for Standard Abutment a or b. Tighten the central clamping screw using the Insert for Prosthetic Ratchet Screwdriver 1.0 mm Hex 15 Ncm.

NOTE: Excessive torque may loosen the threaded sleeve.

Place the ANKYLOS® Welding Cap for Standard Abutment a or b on the respective abutments, make sure it is securely fitted (friction fit) and when using the Welding Cap for lateral connection, screw tight using a lateral screw and a torque of 10 Ncm. Position the titanium wire and only the pin of the ANKYLOS® Welding Cap for Standard Abutment a or b centrally between the two electrodes of the welding unit's welding clamp. The clamping force of the welding clamp secures the devices in position. Support the welding clamp with the palm of your hand during welding and do not open it under any circumstances.

Weld all ANKYLOS® Welding Caps for Standard Abutments a or b with the titanium wire and remove the framework from the Standard Abutments. The pin of the ANKYLOS® Welding Cap for Standard Abutment a or b can be shortened.

Place the finished titanium framework in the patient's mouth on the Standard Abutments and check the passive fit (e.g. Sheffield test).

NOTE: When using the ANKYLOS® Welding Cap for Standard Abutment a or b, insert the corresponding ANKYLOS® Retaining Screw laterally and use the Hexagon Screwdriver 1.0 mm to screw it hand-tight to the ANKYLOS® Standard Abutment. Furnish the finished plastic prosthesis with cold- or light-cured resin

Subject to change without prior notice.

Not all products are available in all countries.

Status of information:
IFU 2024-en_RU Rev. 003 A 2019-09

 **Manufacturer:**
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

 **Dentsply
Sirona
Implants**

Instructions for Use

and place on the framework. For lateral screw connections, make absolutely sure that adequate access is ground free in the plastic prosthesis. After the plastic has hardened, loosen the lateral screw and remove the plastic prosthesis along with the polymerised titanium framework. Then prepare according to the generally accepted dental regulations.

Insertion of the prosthetic restoration using the ANKYLOS® Welding Cap for Standard Abutment a or b
Make sure the ANKYLOS® Standard Abutments are securely fitted and place the prosthetic restoration on the ANKYLOS® Standard Abutments. When cementing the prosthetic restoration, seal the screw access above the central clamping screw of the ANKYLOS® Standard Abutments with gutta percha, cement, cotton pellets, etc. Check to ensure the prosthesis is free of tension. Follow the manufacturer's instructions. Remove all cement residue. Retain laterally using the Insert for Prosthetic Ratchet Screwdriver 1.0 mm Hex, short 10 Ncm, applying the recommended torque of 10 Ncm. After final insertion of the prosthetic restoration, adjust the occlusion.

ANKYLOS® Welding Cap for SynCone 5°

IMPORTANT: Ensure that no threads are trapped when placing the SynCone Cap in the case of immediate loading.

Position the ANKYLOS® Welding Cap for SynCone 5° on the corresponding parallelly aligned SynCone 5° Abutments and check that it is securely fitted (friction fit). Make sure that the SynCone Caps are correctly fitted, otherwise the prosthesis can either be too loose or too tight.
Position the titanium wire and only the pin of the ANKYLOS® Welding Cap for SynCone 5° centrally between the two electrodes of the welding unit's welding clamp. The clamping force of the welding clamp secures the devices in position. Support the welding clamp with the palm of your hand during welding and do not open it under any circumstances. After welding, hold the welding clamp in position for about 3 seconds or, optionally, spray the welding area with water during welding.
Weld all the ANKYLOS® Welding Caps for SynCone 5° with the titanium wire and remove the framework from the ANKYLOS® SynCone 5° Abutments. The pin of the ANKYLOS® Welding Cap for SynCone 5° can be shortened.
Place the SynCone Polymerisation Sleeves on the SynCone Abutment up to below the abutment equator. This prevents the ingress of resin at the sulcus/gingiva and the soft tissue is protected.
Place the finished titanium framework in the patient's mouth securely on the ANKYLOS® SynCone 5° Abutments. Furnish the finished plastic prosthesis with cold- or light-cured resin and place on the framework. Check the passive fit (e.g. Sheffield test). After the plastic has hardened, remove the prosthesis from the patient's mouth. Then prepare according to the generally accepted dental regulations.

Insertion of the prosthetic restoration using the ANKYLOS® Welding Caps for SynCone 5°
Check that the ANKYLOS® SynCone 5° Abutments are fitted securely.
Use the Hexagon Screwdriver 1.0 mm applying the recommended torque of 15 Ncm.
Finally, place the prosthetic restoration on the ANKYLOS® SynCone 5° Abutments.
After final insertion of the prosthetic restoration, adjust the occlusion.

Delivery – Storage – Sterilisation

ANKYLOS® devices for intraoral or extraoral welding are supplied **non-sterile** and, with the exception of the Parallel Gauge, are intended for **single use** only. Do not use ANKYLOS® devices for intraoral or extraoral welding after the expiration date or if the blister is damaged. Store the devices at room temperature in the original packaging under conditions normal for dental practices. Do not expose to direct sunlight. For type, size and, if applicable, angulation see product label.

IMPORTANT: Prior to clinical use, all devices and instruments delivered non-sterile must be cleaned, disinfected and sterilized according to a validated method. Instruments that can be disassembled must be disassembled for cleaning and may be assembled again for sterilization only in a dry condition.

Description	Sterilisation/Disinfection Reuse	Material	Torque
ANKYLOS® Welding Sleeve for Balance Base Abutment narrow	Sterilisation ②	Ti Grade 4	n/a
ANKYLOS® Temporary Welding Abutment	Sterilisation ②	Ti Grade 4	n/a
ANKYLOS® Retaining Screw for Temporary Welding Abutment	Sterilisation ②	Ti6Al4V-ELI	15 Ncm
ANKYLOS® Welding Cap for Standard Abutments a or b	Sterilisation ②	Ti Grade 4	n/a
ANKYLOS® Parallel Gauge for Standard Abutments a or b	Sterilisation Yes	Ti6Al4V-ELI	n/a
ANKYLOS® Welding Cap for SynCone 5°	Sterilisation ②	Ti Grade 2	n/a

FRIADENT® Titanium Wires	Sterilisation	②	Ti Grade 2 Ti Grade 4
② = do not reuse			

Manual pre-cleaning
Coarse soiling must be removed immediately after use (soft nylon brush until visibly clean).

Ultrasonic cleaning (optional for very heavy soiling)
The action times and concentrations specified by the manufacturer of the cleaning agent must be observed.

Mechanical processing
Manual pre-cleaning is necessary to ensure efficient mechanical cleaning. Suitable washer disinfectors as well as cleaning and neutralization agents are to be used for mechanical cleaning. The instructions from the manufacturer of the washer disinfectant must be followed. Cleaning and neutralization agents are to be dosed and used in accordance with the manufacturer's instructions. A cleaning program with thermal disinfection (A₀ ≥ 3000) is recommended. Either demineralized water or water which satisfies this level of purity is recommended for cleaning.

- The mechanical processing method was validated with the following materials and devices:
- Washer disinfectant, type G 7836 CD, Miele & Cie. KG, Gütersloh
 - Cleaning program VARIO-TD
Pre-rinsing with cold water 4 min
Cleaning 10 min from 45–55 °C
Neutralization 6 min
Intermediate rinsing with cold water 3 min
Thermal disinfection 5 min at 90 °C (A₀ ≥ 3000)
Drying 30 min at 105 °C
 - Cleaning agent neodisher® MediClean (0.5 %)/neutralization agent neodisher® Z (0.1 %), Dr. Weigert, Hamburg.

For detailed information on manual and mechanical processing, please refer to the Dentsply Implants Instructions for Sterilization – Instrument Care.

- Sterilization**
Sterilization is to be performed corresponding to the following scheme:
- Preparation for sterilization:** Place devices in a sterilization pouch which is legally marketed for use with the recommended sterilization parameters. Every sterilization package must have a sterilization indicator and sterilization date.
 - Sterilization:**

Method	Cycle	Temperature	Exposure Time*	Dry Time
Steam	Dynamic air removal (prevacuum)	132 °C 135 °C	4 Min. 3 Min.	20 Min.
Steam	Gravity displacement	121 °C	30 Min.	20 Min.

* Minimum exposure times, the operating times are longer and may vary depending on the device.

- Storage:** Store the sterilized devices dry and dust-free at room temperature.

Transportation, Storage and Operating Conditions

Products should be transported and carried carefully in a manner that corresponds to the special transport precautions for each medical device's nature.

Store the devices in their original packaging in a dry place at room temperature. Do not expose to direct sunlight.

- The permissible conditions are:
- Transportation: T -35 °C ... +60 °C, Humidity 10–80 %
 - Storage: T +18 °C ... +25 °C, Humidity 30–60 %
 - Operating: T +32 °C ... +42 °C, Humidity 30–80 %

Terms of Use

Dental clinics and laboratories.
May only be used by dentists/ physicians who have had the appropriate education and training.

Influence on the Environment

The product does not affect the environment during operation.

Disposal of Products

For disposal of the packaging and devices, comply with the currently applicable national waste disposal regulations in your country.
Unused products and packaging material are medical waste hazard class A. Used products are medical waste hazard class B.

Shelf Life

Devices shelf life is 10 years in their respective original packaging. Expiration date indicated on the labeling.

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

Subject to change without prior notice.

Not all products are available in all countries.

Status of information:
IFU 2024-en_RU Rev. 003 A 2019-09

Instructions for Use

Implant System Warranty

For full implant systems warranty information, refer to the detailed terms and conditions in the Implant Systems Warranty leaflet.


Authorised Representative in Russian Federation

Dentsply Sirona LLC
115432, Moscow, Andropov av., 18, bld. 6, offices1, 39, 40, Russia
Tel. +7 (495) 725 10 87

Copyright and Trademarks

All rights reserved. No part of these instructions may be reproduced, processed with electronic systems, duplicated or distributed in any other manner without the express written permission of the manufacturer.
ANKYLOS®, FRIADENT® are registered trademarks of Dentsply Implants Manufacturing GmbH.

Explanation of Symbols

Symbol	Title
	Manufacturer
	Date of manufacture
	Use-by date
	Batch Code
	Catalog Number
	Non-Sterile
	Keep away from sunlight
	Keep dry
	Do not reuse
	Consult Instructions for Use
	European conformity (CE) mark
	Prescription only: U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.
	Russian certification symbol according to GOST

Subject to change without prior notice.

Not all products are available in all countries.

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XiVE® Prosthetics devices for Intraoral or Extraoral Welding

For the list of products covered by this IFU please refer to the XiVE® WeldOne Concept brochure, catalogue or contact your local dealer. For representative product photo please refer to the brochure XiVE® WeldOne Concept.

Intended use

XiVE® Prosthetics devices for intraoral or extraoral welding are designed to provide prosthetic devices, in order to restore the chewing function of the patient. Prosthetics devices are intended to produce an intraoral or extraoral-welded titanium frame for immediate and delayed loading treatment.

Indications

XiVE® MP Welding Sleeve

Fabrication of an intraorally or extraorally welded titanium framework for immediate and late restorations on MP Abutments exclusively by use of the WeldOne welding unit and the titanium wires designed for that purpose

XiVE® Retaining Screw for Welding Sleeve

Fixation of the XiVE® MP Welding Sleeve on MP Abutments

XiVE® Temporary Welding Abutment

Fabrication of an intraorally or extraorally welded titanium framework for immediate temporary restorations on XiVE® Implants exclusively by use of the WeldOne welding unit and the titanium wires designed for that purpose

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Contraindications

- In the case of patients with pacemakers or other active implantable devices, intraoral welding is not indicated and must not be performed under any circumstances.
- Wearers of pacemakers or other active implantable devices must not perform intraoral and extraoral welding.
- The devices are **not** suitable for single-tooth restorations.

Warnings

- These instructions for use must be read prior to using XiVE® devices for intraoral or extraoral welding.
- XiVE® devices for intraoral or extraoral welding may only be used for their intended purpose according to the general rules for dental/surgical treatment, occupational safety and accident prevention.
- XiVE® devices for intraoral or extraoral welding must only be used together with the XiVE® implant system and the WeldOne welding unit.
- The welding of these devices intraorally may only be performed by the dentist; dental technicians may also perform welding extraorally.
- The user must be familiar with dental surgery and prosthetics, including diagnostics and preoperative planning and/or laboratory procedures.
- Dentsply Implants Manufacturing GmbH is not liable for damage resulting from use outside the intended use of the product.
- If the indication and use are unclear, treatment should be suspended until these points have been clarified.
- The following instructions are not adequate to allow inexperienced clinicians to administer professional implant dentistry.
- Intraoral and extraoral welding may only be performed subject to relevant training in order to ensure that the treatment procedure runs smoothly!**

CAUTION: Compatibility of the XiVE® MP Abutments angled

The abutment D3.8-D5.5, anodised in yellow, is compatible with XiVE® S Implants with a diameter of D3.8, D4.5 and D5.5. The use of incompatible devices may cause fitting problems, reduction of fatigue strength or damage to the implant.

CAUTION: MP Abutments are generally not allowed to be ground.

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- Components used in the patient's mouth have been aspirated or swallowed.
- The screw of the abutment has fractured, for example, due to excessive torque.

- Loosening of the abutments/ superstructure for example as a consequence of inadequate torque being applied on the screw.
- Deformation of parts when using non-system devices as well as devices which are not indicated for the intraoral and extraoral welding concept

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- Prior to each procedure it must be ensured that all necessary devices, instruments and materials are complete, available in the required quantities and fully operational.
- All products intended for single-use must not be reused. Failure to observe this precaution leads to a risk of infections and misfitting devices. The associated complications can cause fractures through to implant loss!
- Position the patient such that the danger of aspiration of devices is minimised.
- All devices used in the patient's mouth must be secured to prevent aspiration and swallowing.
- Always wear protective clothing for your own safety.
- Make sure that the patient and user are wearing eye protection during the welding procedure.
- Ensure that all the devices used during welding are free of grease and other residue.
- It is **imperative** to follow the welding parameter settings according to the table in the instructions for use and in the user manual.
- Before the application of resins with suitable tools, like polymerisation sleeves, make sure that no plastic gets between the implant, abutment or welding abutment.
- If, due to unfavourable anatomical conditions, instruments do not fit or cannot be used for other reasons, the course of treatment planned with them must not be continued and alternatives must be sought.
- Observe specified torques.
- Ensure adequate extraction when grinding devices.

Adverse Reactions

No adverse reactions are known

Application

NOTE: Detailed instructions on surgery, prosthetics and laboratory methods using can be found in the respective XiVE® manuals.

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Ø 1.2 mm for additional retention

Ø 1.5 mm for inter-implant distances of up to 8 mm

Ø 2.0 mm for inter-implant distances of between 8 mm and 15 mm

Shorten the titanium wire to the required length with cutting pliers. The titanium wires should be bent using standard dental bending instruments so that these lie against the welding sleeves without tension or gaps. Afterwards, the welding procedure is carried out.

Welding procedure

In order to obtain an appropriate welding joint with the WeldOne welding unit, the XiVE® devices for intraoral or extraoral welding have to be used.

Observe the following instructions before starting the welding procedure: Shorten the titanium wire for intraoral or extraoral welding accordingly and place on the XiVE® welding devices to be welded without tension or gaps.

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Start welding the XiVE® Welding Sleeve devices with the titanium wire at the most accessible position. Continue welding on the adjacent abutment.

Manufacturer:
Dentsply Implants Manufacturing GmbH
Rodenbacher Chaussee 4, 63457 Hanau, Germany
www.dentsplyimplants.com
www.dentsplysirona.com

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Check the titanium framework in the patient's mouth again to make sure that it is free of tension (e.g. Sheffield test) before incorporating it into a plastic prosthesis. If the titanium framework does not sit without tension or gaps, it must be detached and prepared again.

The manufactured titanium framework can be furnished with other retention elements. Shorten the titanium wire and weld extraorally with the titanium framework. Further processing of the titanium framework is carried out according to the generally accepted dental regulations.

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XiVE® MP Welding Sleeve and XiVE® Retaining Screw for Welding Sleeve

NOTE: The XiVE® MP Welding Sleeve and the respective XiVE® MP Retaining Screw for Welding Sleeve can only be used on MP Abutments of the right diameter.

Following the instructions for MP Abutments straight and angled, insert the XiVE® and FRIADENT® MP Abutments in the implant and tighten to a torque of 24 Ncm. Place the XiVE® MP Welding Sleeve onto the MP Abutment and screw hand-tight with the XiVE® MP Retaining Screw for Welding Sleeve. Use the Ratchet Hex Driver 1.22 mm to secure the screw. Shorten the titanium wire accordingly and place on the XiVE® MP Welding Sleeves to be welded so as to be free of tension and gaps. Position the titanium wire and the XiVE® MP Welding Sleeve centrally between the two electrodes of the welding unit's welding clamp. The clamping force of the welding clamp secures the devices in position. Support the welding clamp with the palm of your hand during welding and do not open it under any circumstances. After welding, hold the welding clamp in position for about 3 seconds or, optionally, spray the welding area with water during welding. Carry out a visual inspection of the weld seam after each welding procedure.

Weld all XiVE® MP Welding Sleeves with the titanium wire, release the XiVE® MP Retaining Screws for Welding Sleeves and remove the framework from the MP Abutments. The XiVE® MP Welding Sleeve can be shortened. If, in such cases, the XiVE® MP Retaining Screw for Welding Sleeve is too long, the FRIADENT® MP Screw for castable sleeve can be used.

Before incorporating the titanium framework into a plastic prosthesis, place it in the patient's mouth again to ensure that it is free of tension (e.g. Sheffield test). If the titanium framework does not sit without tension or gaps, it must be detached and prepared again.

Insert the polymerisation sleeves. Place the finished titanium framework in the patient's mouth on the MP Abutments and screw hand-tight at the most easily accessible point using a FRIADENT® MP Screw for castable sleeve. Furnish the finished plastic prosthesis with cold- or light-cured resin and place on the titanium framework. After the plastic has hardened, undo the retaining screw and remove the prosthesis, then prepare according to the generally accepted dental regulations.

Insertion of the prosthetic restoration using the XiVE® MP Welding Sleeve
Recheck the secure seating of the MP Insert Screw or, with angled MP Abutments, the insert head with a torque of 24 Ncm, then place the prosthetic restoration on the MP Abutments. Check the passive seating of the prosthetic restoration (e.g. Sheffield Test), tighten the prosthetic screws with the Ratchet Hex Driver 1.22 mm and a torque of 24 Ncm and check the occlusion.

XiVE® Temporary Welding Abutment

After inserting the XiVE® S Implants as parallel as possible, remove the premounted TempBase abutment from the implant with a Hex Driver 0.9 mm. Place the XiVE® Temporary Welding Abutment in the implant and screw hand-tight using the retaining screw supplied and the Hex Driver 1.22 mm.

Position the titanium wire and the XiVE® Temporary Welding Abutments centrally between the two electrodes of the welding unit's welding clamp. The spring of the welding clamp secures the devices in position. Support the welding clamp with the palm of your hand during welding and do not open it under any circumstances. After welding, hold the welding clamp in position for about 3 seconds or, optionally, spray the welding area with water during welding. Carry out a visual inspection of the weld seam after each welding procedure. The XiVE® Temporary Welding Abutment can be shortened occlusally only up to the height of the screw head.

Before incorporating the titanium framework into a plastic prosthesis, place it in the patient's mouth again to ensure that it is free of tension (e.g. Sheffield test). If the titanium framework does not sit without tension or gaps, it must be detached and prepared again.

Insert the polymerisation sleeves. Insert the finished titanium framework in the patient's mouth in the implants and screw hand-tight with the retaining screw provided at the most accessible location. Furnish the finished plastic prosthesis with cold- or light-cured resin and place on the framework. After the plastic has hardened, undo the retaining screw and prepare the plastic prosthesis according to the generally accepted dental regulations.

Insertion of the prosthetic restoration using the XiVE® Temporary Welding Abutment

Insert the prosthetic restoration in the implants. Make sure that the abutment is correctly inserted in the implant and that no soft tissue is wedged between the abutment and implant. Check the passive seating of the prosthetic restoration, tighten the prosthetic screws with the Ratchet Hex Driver 1.22 mm and torque of 24 Ncm and check the occlusion.

Delivery – Storage – Sterilisation

XiVE® devices for intraoral or extraoral welding are supplied **non-sterile** and are intended for **single use** only. Do not use XiVE® devices for intraoral or extraoral welding after the expiration date or if the blister is damaged. Store the devices at room temperature in the original packaging under conditions normal for dental practices. Do not expose to direct sunlight. For type, size and, if applicable, angulation see product label.

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(2) = do not reuse

Manual pre-cleaning

Coarse soiling must be removed immediately after use (soft nylon brush until visibly clean).

Ultrasonic cleaning (optional for very heavy soiling)

The action times and concentrations specified by the manufacturer of the cleaning agent must be observed.

Mechanical processing

Manual pre-cleaning is necessary to ensure efficient mechanical cleaning.

Suitable washer disinfectors as well as cleaning and neutralization agents are to be used for mechanical cleaning. The instructions from the manufacturer of the washer disinfectant must be followed. Cleaning and neutralization agents are to be dosed and used in accordance with the manufacturer's instructions. A cleaning program with thermal disinfection ($A_0 \geq 3000$) is recommended. Either demineralized water or water which satisfies this level of purity is recommended for cleaning.

The mechanical processing method was validated with the following materials and devices:

- Washer disinfectant, type G 7836 CD, Miele & Cie. KG, Gütersloh
- Cleaning program VARIO-TD
Pre-rinsing with cold water 4 min
Cleaning 10 min from 45–55 °C
Neutralization 6 min
Intermediate rinsing with cold water 3 min

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Sterilization

Sterilization is to be performed corresponding to the following scheme:

1. **Preparation for sterilization:** Place devices in a sterilization pouch which is legally marketed for use with the recommended sterilization parameters. Every sterilization package must have a sterilization indicator and sterilization date.

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3. **Storage:** Store the sterilized devices dry and dust-free at room temperature.

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Terms of Use

Dental clinics and laboratories

May only be used by dentists/ physicians who have had the appropriate education and training.

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Disposal of Products

For disposal of the packaging and devices, comply with the currently applicable national waste disposal regulations in your country.

Unused products and packaging material are medical waste hazard class A.

Used products are medical waste hazard class B.

Shelf Life

Devices shelf life is 10 years in their respective original packaging.

Expiration date indicated on the labeling.

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Dentsply Sirona LLC

115432, Moscow, Andropov av., 18, bld. 6, offices1, 39, 40, Russia

Tel. +7 (495) 725 10 87

Copyright and Trademarks

All rights reserved. No part of these instructions may be reproduced, processed with electronic systems, duplicated or distributed in any other manner without the express written permission of the manufacturer.

XiVE® and FRIADENT® are registered trademarks of Dentsply Implants Manufacturing GmbH.

Explanation of Symbols

Symbol

Title

Manufacturer

Date of manufacture

Use-by date

Batch Code

Catalog Number

Non-Sterile

Keep away from sunlight


Keep dry

Do not reuse

Consult Instructions for Use


European conformity (CE) mark

Prescription only:

U.S. Federal law restricts this device to sale by or on the order of a licensed dentist or physician.

Russian certification symbol according to GOST

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

Facilitating
stable and
durable
chairside
solutions

Ankylos®


The WeldOne® concept

 **Dentsply
Sirona**
Implants


Courtesy Dr. Marco Degidi, Bologna, Italy

With the help of WeldOne Welding Unit and special welding components, a titanium wire is permanently connected onto the abutment.

Immediate reinforcement of the temporary by a titanium framework.

Specially designed Ankylos welding components.

What makes
WeldOne®
for Ankylos®
so special:

WeldOne® – worldwide unique

The first complete evidence based welding concept

The WeldOne concept is worldwide unique as Dentsply Sirona Implants offers the first complete, proven and sophisticated concept for intraoral welding based on resistance spot welding principle.

With the help of the WeldOne Welding Unit and special welding components, a titanium wire is permanently connected to the abutment with precisely the correct force to achieve an optimum weld¹. An electric current leads to fusion of the titanium material in the contact points between the titanium wire and the abutment in milliseconds.

The permanent joint between the titanium wire and the abutment enables an immediate and stable titanium framework to be made which will provide the optimum support for temporary or durable restorations.

This concept is exclusively designed for dental implants made by Dentsply Sirona Implants.

Different procedures of fabrication can be selected for this unique concept of prosthetic restoration:

- Temporary restorations on implant level
- Temporary or durable restorations on abutment level

The WeldOne concept provides a stable, passively fitting framework for temporary or durable prostheses for immediate restorations on the same day of surgery.

It is a fast and economical solution to deliver partial- and full-arch restorations – suitable for immediate or late loading.

- **Meets the patients' demand for immediate restorations^{1,2}**

Surgical procedure and prosthetic delivery on the same day of surgery.

- **Significant time and cost reduction¹**

Standardized protocol with specially designed components for a reduction of appointments to a minimum and the opportunity to reduce total treatment costs.

- **Less risks of fractures on temporaries^{4,5}**

Immediate reinforcement of the temporary by a titanium framework.


- **Safely and predictably for immediate loading treatment^{1,3,5}**

As a consequence of reduced mobility of implants due to immediate rigid splinting a risk of implant failures during healing phases can be minimized.

- **No additional impression-taking of implants necessary^{1,4}**

Simultaneous delivery of a titanium framework on implant abutments at one time. Advantage of the one abutment at one time philosophy in case of abutment level restorations.

1 | Degidi, M., D. Nardi, and A. Piattelli, Immediate rehabilitation of the edentulous mandible with a definitive prosthesis supported by an intraorally welded titanium bar, in *Int J Oral Maxillofac Implants*, Yes, Editor. 2009, 01 until 03 years. p. 342-7. 2 | Avvanzo, P., et al., Use of intraoral welding to stabilize dental implants in augmented sites for immediate provisionalization: a case report. *J Oral Implantol*, 2012. 38(1): p. 33-41. 3 | Degidi, M., G. Dapriile, and A. Piattelli, Implants inserted with low insertion torque values for intraoral welded full-arch prosthesis: 1-year follow-up, in *Clin Implant Dent Relat Res*, Yes, Editor. 2012, 00 until 01 year. p. e39-40. 4 | Degidi, M., et al., Syncrystallization: a technique for temporization of immediately loaded implants with metal-reinforced acrylic resin restorations, in *Clin Implant Dent Relat Res*, Yes, Editor. 2006, 00 until 01 year. p. 123-34. 5 | DENTSPLY Implants internal fatigue, strain and welding tests on file

A close-up portrait of an elderly man with white hair, smiling gently. He is wearing a dark blue suit jacket over a light blue button-down shirt. The background is a soft-focus outdoor scene with green foliage and sunlight filtering through the leaves.

“Vision
becomes
reality”

Dr. Marco Degidi
Bologna, Italy

Placement of the Ankylos temporary welding abutments

Intraorally welded titanium framework

Completed restoration

Courtesy Dr. Marco Degidi,
Bologna, Italy

Intra-oral welding

Known for decades – now brought to perfection

The inventors of this concept are Dr. Marco Degidi and his dental technician Gianluca Sighinolfi both from Bologna, Italy. Their excellent knowledge, craftsmanship and teamwork are the basis for the successful creation of this unique stable and durable chairside concept. The teamwork of the surgeon and the dental technician is key for success.

Although the WeldOne concept is cutting-new, an impressive amount of publications in peer reviewed journals is already available (see below).

The WeldOne concept is an advanced protocol addressed to surgeons and prosthodontists that already have experience with immediate loading treatment.

Special training courses tailored to meet the indication requirements and the experience level are highly recommended. Please see last page of this brochure for course offers or contact your local representative.

Dr. Marco Degidi Bologna, Italy

Active member of the Academy of Osseointegration, Diplomate of International College of Oral Implantology, Private Practice in Bologna, restricted to implant surgery & prosthetics, Adjunct Professor, Postgraduate Program in Clinic Implantology and Biomaterials, University of Chieti, Adjunct Professor, Postgraduate Master in Clinic Implantology, University of Bologna, Speaker in many national and international seminars, congresses and courses; many national and international publications.




Gianluca Sighinolfi, Bologna, Italy

In 1989 he was awarded his diploma from the Villaggio del Fanciullo in Bologna, Italy. He immediately began work in an important dental lab in Bologna, later to become a partner. He has had extensive specialist training and conducted courses both in Italy and internationally. Since 1996, Gianluca has collaborated with Dr Marco Degidi and Prof Piattelli on immediate loading research, contributing to the development of new surgery protocols and numerous articles on this subject that have been published in leading international reviews and journals.

Step-by-step

The WeldOne® concept for Ankylos®


Specially designed Ankylos
Temporary Welding
Abutments

Ankylos
Welding Cap
for Standard
Abutments

Ankylos
Welding
Sleeve for
Balance Base
Abutment
narrow

Ankylos
Welding
Cap for
SynCone 5°


Bending, alignment
and intraoral welding of
titanium wire

Extraorally welded
retentions are added
to the framework

Finalized and
opaque coated titanium
framework

Completed
restoration

The WeldOne® concept

Specially created WeldOne
Welding Unit with touch screen
and easy menu guidance


It is recommended to use the following custom-built tools manufactured by USTOMED. This handful of simple tools allows you to create a titanium framework with a few simple steps. To get detailed information, please contact USTOMED directly: www.ustomed.de


The WeldOne®
concept for
Ankylos® -
designed function


Sharp Curver
bends the bar into
more acute angles
for single abutments.


Soft Curver
bends the titanium
bar along the
curve of the arch.



Holder
grasps and holds
the titanium
bar firmly in
place during the
welding process.

**Cutter with
tungsten carbide
inserts**
easily and
efficiently cuts
off excess
material.

Specially designed Ankylos
welding components

Three sizes of titanium wires to en-
sure the stability of the framework
according to the right application
and respective indication.

Products

WeldOne Welding Unit

Reduced
scale

Titanium Wire
D 1.2 / L 100
Grade 2
(10 pcs)

Titanium Wire
D 1.5 / L 150
Grade 2
(10 pcs)

Titanium Wire
D 2.0 / L 150
Grade 2
(10 pcs)

Order no.	47 - 4710	47 - 3012	47 - 3013	47 - 3014
-----------	-----------	-----------	-----------	-----------

Ankylos Welding Cap for Standard Abutment

Scale 1.2:1

a / 4.0

a / 4.0 lateral screw

a / 6.0

a / 6.0 lateral screw

Order no.	3105 6091	3105 6092	3105 6093	3105 6094
-----------	-----------	-----------	-----------	-----------

Scale 1.2:1

b / 4.0

b / 4.0 lateral screw

b / 6.0


b / 6.0 lateral screw

Order no.	3105 6095	3105 6096	3105 6097	3105 6098
-----------	-----------	-----------	-----------	-----------

Ankylos Parallel Gauge for Standard Abutment

Scale 1.2:1

a / 4.0

a / 6.0

b / 4.0

b / 6.0

Order no.	3103 3605	3103 3606	3103 3607	3103 3608
-----------	-----------	-----------	-----------	-----------

Ankylos Welding Sleeve for BalanceBase Abutment narrow

Scale 1.2:1

Ankylos Temporary Welding Abutment

GH 1.5

GH 3.0

GH 4.5

Ankylos Welding Cap for SynCone 5°

Order no.	3105 6219*	3102 2661*	3102 2662*	3102 2663*	3105 6085
-----------	------------	------------	------------	------------	-----------

* incl. retaining screw

Expert training

The WeldOne® concept for Ankylos®


Please contact your local sales representative to learn more about WeldOne product availability and specifications, trainings and congress presences.

- Practical and theoretical approach
- Hands-on supported

To gain the maximum benefit from the course, it is recommended that your lab technician attends the course as well.

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DENTSPLY Implants
Manufacturing GmbH
Rodenbacher Chaussee 4
63457 Hanau
Germany

About Dentsply Sirona Implants

Dentsply Sirona Implants offers comprehensive solutions for all phases of implant therapy, including Ankylos[®], Astra Tech Implant System[®] and Xive[®] implant lines, digital technologies, such as Atlantis[®] patient-specific solutions and Simplant[®] guided surgery, Symbios[®] regenerative solutions, and professional and business development programs, such as STEPPS[™]. Dentsply Sirona Implants creates value for dental professionals and allows for predictable and lasting implant treatment outcomes, resulting in enhanced quality of life for patients.

About Dentsply Sirona

Dentsply Sirona is the world's largest manufacturer of professional dental products and technologies, with a 130-year history of innovation and service to the dental industry and patients worldwide. Dentsply Sirona develops, manufactures, and markets a comprehensive solutions offering including dental and oral health products as well as other consumable medical devices under a strong portfolio of world class brands. As The Dental Solutions Company[™], Dentsply Sirona's products provide innovative, high-quality and effective solutions to advance patient care and deliver better, safer and faster dentistry. Dentsply Sirona's global headquarters is located in York, Pennsylvania, and the international headquarters is based in Salzburg, Austria. The company's shares are listed in the United States on NASDAQ under the symbol XRAY.

Visit www.dentsplysirona.com for more information about Dentsply Sirona and its products.

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
Specially designed Xive welding components.

What makes
WeldOne®
for Xive® so
special:

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With the help of the WeldOne Welding Unit and special welding components, a titanium wire is permanently connected to the abutment with precisely the correct force to achieve an optimum weld⁵. An electric current leads to fusion of the titanium material in the contact points between the titanium wire and the abutment in milliseconds.

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

“Vision
becomes
reality”

Dr. Marco Degidi
Bologna, Italy

Placement of the Xive temporary welding abutments

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Additional product information:

Dr. Marco Degidi Bologna, Italy

Active member of the Academy of Osseointegration, Diplomate of International College of Oral Implantology. Private Practice in Bologna, restricted to implant surgery & prosthetics. Adjunct Professor, Postgraduate Program in Clinic Implantology and Biomaterials, University of Chieti. Adjunct Professor, Postgraduate Master in Clinic Implantology, University of Bologna. Speaker in many national and international seminars, congresses and courses; many national and international publications.

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

Step-by-step

The WeldOne® concept for Xive®


Specially designed Xive
Temporary Welding
Abutments

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


The WeldOne® concept for Xive® – designed function


Sharp Curver
bends the bar into more acute angles for single abutments.


Soft Curver
bends the titanium bar along the curve of the arch.


Holder
grasps and holds the titanium bar firmly in place during the welding process.

Cutter with tungsten carbide inserts
easily and efficiently cuts off excess material.

Specially designed Xive welding components

Three sizes of titanium wires to ensure the stability of the framework according to the right application and respective indication.

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

DENTSPLY Implants
Manufacturing GmbH
Rodenbacher Chaussee 4
63457 Hanau
Germany

About Dentsply Sirona Implants

Dentsply Sirona Implants offers comprehensive solutions for all phases of implant therapy, including Ankylos[®], Astra Tech Implant System[®] and Xive[®] implant lines, digital technologies, such as Atlantis[®] patient-specific solutions and Siplant[®] guided surgery, Symbios[®] regenerative solutions, and professional and business development programs, such as STEPPS[™]. Dentsply Sirona Implants creates value for dental professionals and allows for predictable and lasting implant treatment outcomes, resulting in enhanced quality of life for patients.

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