
78815

WeldOne™ Welding Unit

Follow instructions for use!

 Manufacturer:

CE
1936

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Distributor:

**Dentsply Implants
Manufacturing GmbH**

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Phone: +49 (0) 621
4302-010 Fax: +49 (0) 621
4302-011

www.dentsplyimplants.com

Authorized Representative in the Russian
Federation:

Dentsplat Sirona LLC
115432, Russia, Moscow
Andropov Avenue, 18, building. 6,
offices 1 39 40

Tel.: +7 (495) 725 10 87
E-Mail: cis@dentsplysirona.com
www.dentsplysirona.com
www.dentsplycis.com

Subject to change without notice.
Not all products are available in all countries.
For further information please contact your DENTSPLY
Implants representative.

Status of information: 2018-06 Rev. 005

Table of Contents

Intended use	4
Contraindications	4
Warnings	4
Safety precautions	5
Adverse reactions	6
Conventions used in the manual	7
Explanation of symbols	7
Product Labeling	7
Type Plate	7
Scope of delivery	8
Inspection on delivery	8
Technical data	8
Installation	9
Installation site requirements	9
EMC requirements	9
Electrical installation requirements	9
Place of installation	9
Parts identification	10
Connections on the welding unit	10
Initial Setup	11
Setup	11
Power on	11
User interface	12
Operation	13
Power on	13
Start page	13
Main Menu	13
Selecting the display language	13
Standard Welding Program "SMART"	14
Starting a welding cycle in the "SMART" mode	14
Welding programs	14
Creating a welding program	15
Welding parameter settings	15
Troubleshooting	16
Error messages	16
Recoverable errors	16

Non-recoverable errors	16
Correct use of the welding clamp	17
Positioning of the welding clamp	17
Cleaning, disinfection and sterilization	18
Cleaning and disinfection of the welding unit	18
WeldOne™ housing	18
Clamp connecting cables	18
Reprocessing of the welding clamp	18
Manual pre-cleaning	18
Ultrasonic pre-cleaning	18
Automatic cleaning	18
Preparation for sterilization	19
Sterilization	19
Storage	19
Maintenance	19
Routine maintenance	19
Periodic maintenance	20
Replacing the copper electrodes	20
Replacing the fuses	21
Decommissioning	21
Change of location	21
Disposal at the end of the unit's lifetime	22
Technical Support	22
Warranty and repair	23
Warranty	23
Repair	23
Appendix	24
Electromagnetic Compatibility (EMC)	24
Tables 1, 2, 4, 6 provided for by standard 60601-1/2 (2010)	25

User's Manual

Intended use

WeldOne™ is an electric resistance-pressure spot welding unit designed for the intraoral and extraoral welding technique in dental implantology. The welding unit is intended for welding titanium wires to abutments/caps on dental implants made by DENTSPLY Implants directly inside the oral cavity of the patient.

Contraindications

- The welding unit must not be used neither on patients nor by operators who are carriers of pacemakers, implantable defibrillators or other electronic life support devices.
- The welding clamp parts that come into direct contact with the patient are made from brass, copper and trivalent chromium. Even if these materials are usually well tolerated, it cannot be ruled out that some people may be allergic to them.

Warnings

Read this User's Manual carefully prior to using the welding unit for the first time to fully understand the operation and the function of the welding unit in order to prevent insufficient information from being the cause of an inadequate and unsafe use of the welding unit.

Always keep this User's Manual near the welding unit so that it is readily accessible when needed.

The welding unit may be used only in conjunction with the welding components by DENTSPLY Implants and on implants made by DENTSPLY Implants. The use of components from other suppliers may affect the safety and the result of the welding.

The welding unit may be used only by dentists or maxillofacial surgeons knowledgeable of the dental technology in a dentist's surgery under conditions normal for dental practices. Intraoral welding may be performed by the dentist only; extraoral welding may also be performed by dental technicians.

The following instructions are not adequate to allow inexperienced clinicians to perform the intraoral welding technique. Intraoral or extraoral welding with the welding unit may only be performed after successful completion of adequate training. This is especially necessary to ensure that the treatment goes smoothly.

Importers, distributors, dealers and agents are required to train their interlocutors (buyers and non-buyers) on correct

use of the welding unit. NEVER allow the welding unit to be used by unauthorized personnel.

The welding unit must be installed and operated strictly adhering to the instructions provided in this manual paying special attention to the information related to safety.

The welding unit may be used only for its intended purpose and according to the general rules for dental/surgical treatment, occupational safety and accident prevention.

Use the equipment only for the planned operations. Any use outside the intended use is inappropriate and, as a consequence, dangerous and therefore is prohibited. The manufacturer of the welding unit is not responsible for the safety, reliability and performance effects caused by inappropriate use of the welding unit.

Installation

- The welding unit is not designed for outdoor use.
- The welding unit must not be placed near sources of heat such as radiators or heaters.
- The welding unit is compliant with the EMC requirements only when used with the original accessories and spare parts. The use of other accessories/spare parts can lead to an increased emission of electromagnetic interference or to reduced resistance against electromagnetic interference.
- The welding unit requires particular precautions regarding EMI and must be installed and put into service in compliance with the information on EMI provided in the appendix to this manual.
- Portable and mobile radio communication devices might affect operation of the welding unit.
- Follow the instructions given in this user's manual to correctly assemble and install the welding unit and its accessories.
- ALWAYS connect the welding unit to an electrical system that is suitable for rooms used for medical purposes and that is compliant with IEC 60364 standard (Electrical Installation for buildings).
- ALWAYS check that the residual-current device of your facility is operational.
- Always check that the mains cable is in good condition.
- Check that the electrical mains power supply meets the requirements listed in the technical data of the welding

unit, particularly regarding power supply voltage, deliverable power and mains frequency.

- ALWAYS cover the mains cable with appropriate devices so that it is secured to the floor, and ensure that it is positioned in a way that it will not become a tripping hazard.
- Do not put the welding unit in a position that makes it difficult to operate the disconnection device (power supply cable plug). The power supply socket must be nearby and easy to access.

Operation

- **WARNING:** DO NOT use the welding unit in the presence of an anesthetic mixture flammable with air, oxygen or nitrous oxide.
- **WARNING:** Make sure that all flammable solutions normally used for disinfecting have evaporated completely before turning the welding unit on since any sparks generated by the discharge during welding may cause fire and explosion.
- **WARNING:** Potentially flammable anesthetics or gases such as oxygen and nitrous oxide that can cause fire on other materials, such as cotton if exposed to sparks, for example, must NOT be used.
- Do not use portable and mobile HF communication equipment (e.g. mobile phones) during operation. These may affect medical electrical equipment.
- Stop using the welding unit immediately in case of malfunctioning or abnormal noises.
- NEVER get distracted when the welding unit is in operation.
- NEVER disconnect the clamp connecting cable from the welding unit while the welding clamp has contact with the patient.

Welding parameter settings

- When using an excess power factor for a very thin wire, there is a risk of completely melting the metal and generating spits of molten metal in the patient's mouth.

Maintenance

- Routinely check that the accessories - especially the electric cables, the welding clamp and the copper electrodes for the clamp - are in good condition.

- Follow the maintenance instructions of the welding unit provided in this user's manual to prevent damage over time due to wear or negligence.
- **WARNING:** NEVER remove, damage or modify the protection and safety devices; if the mains fuses are replaced, use only those shown on the rating plate fastened near the power socket.

Cleaning, sterilization

- Do not use any abrasive substances or solvents for cleaning the housing.
- Do not use any liquid or spray directly on the panel or near the electrical sockets.
- The welding unit is not sterilizable.
- The welding clamp supplied with the welding unit is non-sterile upon delivery and therefore must be sterilized and processed before use.

Safety precautions

The following precautions are to be met prior to the treatment:

- **WARNING:** Users with pacemakers or other active implantable devices must not perform intraoral welding.
- **WARNING:** Using flammable anesthetic gases, nitrous oxide or oxygen, as well as alcohol in the oral cavity of the patient before or during the welding procedure is prohibited.
- Always wear protective clothing and protective goggles for your own safety.
- Make sure that the patient wears eye protection during the welding procedure. Position the patient such that the danger of aspiration of components is minimized. All components used in the patient's mouth must be secured to prevent aspiration and swallowing.
- Position the patient in such a way so that when in contact with the welding clamp, he/she does not come into contact with accessible grounded metallic parts (of the welding unit or other nearby devices).
- For working under sterile conditions it is recommended to cover the clamp connecting cables with a disposable sterile sleeve.
- ONLY use the welding unit if it is fully operational.
- Ensure that all necessary components and materials are available in the required quantities, fully operational,

User's Manual

properly connected and fastened before beginning a new weld.

- It is strongly recommended to test the chosen parameter settings extraorally in advance to verify whether they give good welding results, and to validate them.
- Before starting the welding phase the patient must be informed about the steps of the treatment to prevent any imprevident reaction of the patient.
- Ensure that third parties not involved in the treatment are not near the welding unit when in use.

The following precautions are to be met during the treatment:

- **CAUTION:** Make sure that the mucous membranes, soft tissues, suture thread or other materials are not trapped between the electrodes.
- **IMPORTANT:** It is **imperative** to observe the welding parameter settings provided in this User's Manual.
- The welding unit does not require connection with other devices, so prevent its metallic parts (including the applied parts) from coming into contact with other devices.
- **CAUTION:** When the user is in contact with the patient, he/she must avoid touching the accessible grounded metallic parts of the welding unit or of any other nearby devices.
- **CAUTION:** Extreme heat is generated during welding. Therefore avoid contact between the welded metal joints and the patient's tongue or the soft tissues in general. It is imperative to hold the welding clamp in position for a few seconds after the welding in order to facilitate thermal dissipation of heat through the tips.
- **CAUTION:** Correct use of the welding clamp is fundamental to ensure the safety level necessary for the patient and user.
- **CAUTION:** Avoid slipping of the welding clamp during the welding process, and **NEVER** apply manual pressure on the two arms of the clamp during the welding operation. Sparks can be caused by a failed contact or incorrect contact of the clamp ends or by clamp openings during the welding cycle.

Adverse reactions

- The welding unit generates very high pulses of current during normal operation which leads to a considerable emission of pulsed electromagnetic fields. Therefore, the welding unit must not be used in the presence of patients and operators who are carriers of pacemakers, implantable defibrillators or other electronic life support devices. When reliable data is missing **DO NOT** use the welding unit.
- The welding clamp is made from brass, copper and trivalent chromium. The copper electrodes are made from copper. Even if these materials are usually well tolerated, it cannot be ruled out that some people may be allergic to them.

The parts of the welding unit that come into direct contact with the patient are limited and the application time is very short. Whether or not it is advisable to use this instrument, after having made a careful evaluation, is left to the dentist's discretion for the patient to be treated.

Conventions used in the manual

The warning triangle symbol is used to warn against operating procedures that could result in an accident and/or personal injury and to highlight cautions against handling procedures that can cause damage to the welding unit and its accessories.

Pay utmost attention to the text in this box as it contains important information regarding patient and operator safety and the correct use of the welding unit.

Additional information and useful explanations are marked with an information symbol.

Includes additional information and further explanations for the user.

Explanation of symbols

The following symbols and abbreviations are used on the welding unit, product labels, and in the User's Manual:

Reference number

Serial number

Production lot

Year of manufacture

Manufacturer

1936

CE mark of conformity with number of Notified Body

Type BF Applied Part

Follow instructions for use

Consult instructions for use

WEEE (Waste of Electrical and Electronic Equipment) symbol

IPX0

Degree of protection against the penetration of solid bodies, powders and liquids inside the housing

Electric fuse

Product Labeling

Rating Plate

The following labels and symbols are fixed to the unit in a permanent and legible manner.

The rating plate is attached to the rear of the housing.

- Name of equipment
- Manufacturer's name and address
- Country of origin
- Reference number
- Year of manufacture
- Serial number
- Mains voltage
- Frequency
- Protection grade against the penetration of solid bodies, powders and liquids inside the housing
- Maximum power consumption
- Fuses
- Degree of protection against electrical shock

User's Manual

Scope of delivery

The standard delivery includes the following:

Product designation	DENTSPLY Implants Ref.
WeldOne™ welding unit	47-4710

The following accessories are supplied with the welding unit:

Type	DENTSPLY Implants Ref.
Welding clamp	47-4711
Set of replacement electrodes	47-4712
Clamp connecting cable	47-4713
Pedal switch	47-4714
Mains cable	47-4715
Electrode replacement tool	47-4716
Clamp retainer	47-4717
Fuses	47-4718

Spare parts can be reordered from DENTSPLY Implants as well.

The shipment includes in addition:

- User's Manual in five languages
- Unpacking instructions including list of components
- A safety sign to be placed near the instrument
- Repair Request form of the manufacturer

The name and address of the manufacturer is indicated on the shipping label on the packaging.

Inspection on delivery

Upon delivery of the shipment it is important that the following checks be made:

- Check the transport box for external damage.
- Check the contents of the package against the enclosed list of components and correspondence with the relevant codes.
- Check the welding unit and the accessory components for visible damage.

For all the problems related to transport and delivery, the laws currently in force apply.

We strongly recommend to keep all the packaging material. This ensures a safe transport for the welding unit and the accessories.
Use only the original packaging for return shipments.

Technical data

Power supply voltage:	220-240 V AC
Mains frequency:	50/60 Hz
Power consumption:	150 W, 8 W in standby
Fuses:	2x T3 15A
Medical device class in accordance with Directive 93/42/EEC:	IIb
Protection class:	I
Insulation class:	Type BF Applied Part 
Operating conditions:	
Ambient temperature:	+10 °C – +30 °C
Relative air humidity:	0% – 90%
Atmospheric pressure:	70 – 110 kPa
Storage and transport conditions:	
Temperature range:	-30 °C – +80 °C
Relative air humidity:	0% – 90%
Dimensions (w x h x d):	380 x 310 x 120 mm
Weight:	8.4 kg ±2%

The welding unit is CE marked in accordance with the applicable European Directives 93/42/EC and Amending Directive 2007/47/EC

Installation

Observe the instructions given below in order to prevent malfunctions caused by external influence.

Installation site requirements

For the installation requirements regarding temperatures, air humidity and atmospheric pressure see **Technical data**.

- The welding unit is not designed for outdoor use.
- The welding unit is not suitable for use in the presence of anesthetic mixtures flammable with air, oxygen or nitrous oxide, which constitute an explosion hazard!
- The welding unit must not be placed near sources of heat such as radiators or heaters.
- Do not put the welding unit in a position that makes it difficult to operate the disconnection device (mains cable plug). The mains socket must be nearby and easy to access.

EMC requirements

- The welding unit is compliant with the EMC requirements only when used with the original accessories and spare parts. The use of other accessories/spare parts can lead to an increased emission of electromagnetic interference or to reduced resistance against electromagnetic interference.
- The welding unit requires particular precautions regarding EMI (electromagnetic interference) and must be installed and put into service in compliance with the information on EMI provided in the appendix to this manual.
- Portable and mobile radio communication devices might affect operation of the welding unit.
- The welding unit must not be in contact with other electronic devices and must be positioned at a safe distance from the devices with which it could have electromagnetic interferences (see the tables with the recommended distances).

Electrical installation requirements

- Connect the welding unit to an electrical system that is suitable for rooms used for medical purposes and that is compliant with IEC 60364 standard (Electrical Installation for buildings).
- Ensure that the electrical system to which the welding unit is connected is compliant with the current standards and is adequately sized to the characteristics of the welding unit.
- Check that the residual-current device of your facility is operational.
- Check that the electrical mains power supply meets the requirements listed in the technical data of the welding unit, particularly regarding power supply voltage, deliverable power and mains frequency.


Place of installation

The welding unit, including its accessories, weighs about 8.6 kg. Take appropriate precautions when handling it.

- Place the welding unit on a table or on an appropriate support of the dentist's unit with the adequate stability and load bearing capacity to ensure the safe operation of the unit.
- The housing of the welding unit has an IPX0 protection rating, therefore avoid positioning it where accidental drops of liquid on it are possible.
- ALWAYS cover the mains cable with appropriate devices so that it is secured to the floor, and ensure that it is positioned in a way that it will not become a tripping hazard.

User's Manual

Parts identification


Product designation

- 1 Welding clamp with copper electrodes
- 2 Electrode replacement tool
- 3 Clamp retainer
- 4 Clamp connecting cable
- 5 On/off switch
- 6 Fuse holder
- 7 Mains socket
- 8 Connecting socket for pedal switch
- 9 Pedal switch
- 10 Holder for welding clamp

Function/Purpose

- For welding titanium wires to prosthetic components on dental implants
- To remove used electrodes from the clamp and replace by new ones
- To hold the clamp open during sterilization
- For connecting the welding clamp to the welding unit
- Switching the welding unit on/off
- For storing the fuses
- For connecting the mains cable
- For connecting the pedal switch to the welding unit
- When enabled initiates the welding cycle by pressing the pedal switch.
- For holding the welding clamp when not in use.

Connections on the welding unit


10



IFU 4700 Rev. 05 – 2018-06

IFU 47

Setup

Insert the cylindrical connectors of the welding clamp connecting cable into the socket on the left side of the welding unit and turn the locking ring clockwise until fully coupled.

Do not press any key on the touch panel until the initial diagnostic procedure is completed and the installed software is displayed. Otherwise the automatic calibration of the keys will be disturbed.


					W	e	i	d	O	n	e				
					S	O	F	T	W	A	R	E			
					V	e	r	X	X						

User's Manual

User interface


The user interface consists of four status LEDs, an LCD, and the keys.

The four LEDs [POWER ON] [READY] [WELD] and [FAULT] indicate the status of the various functions.

LED	Status	Function/Purpose
1 [POWER ON] BLUE	on	Welding unit turned on and operational
2 [READY] GREEN	flashing	Ready for starting a welding cycle by activating the pedal switch
3 [WELD] YELLOW	on	Welding cycle in progress
4 [FAULT] RED	on	Indicates a malfunction

The display is an alphanumeric four-line LCD, each line with 20 single characters. When the welding unit is on, the display is always illuminated.

[▲] [▼]
Press to scroll the menu and to increase or decrease the data values.

[MENU]
Press to return to the previous menu.

[ENTER]
Press to confirm the selected value or the menu line to proceed to the corresponding submenu.

[WELD]
Press to enable welding, which will start when the pedal switch is pressed.

Operation

Power on

Switch on the welding unit by pressing the ON/OFF switch.

Do not press any key on the panel until the initial self-test is completed and the installed software is displayed. Otherwise the automatic calibration of the keys will be disturbed.

Start page

					W	e	l	d	O	n	e				
					S	O	F	T	W	A	R	E			
					V	e	r		X		X				

The installed software version is indicated on the LCD for a few seconds after the initial self-test is completed.

All LEDs [READY] [WELD] [FAULT] [WELD] are off.

The pedal switch and all keys are still disabled

After a few seconds the following three warnings are displayed:

				O	X	Y	G	E	N		O	F	F	!	
			N	O		P	A	C	E	M	A	K	E	R	!
		P	R	O	T	E	C	T	I	V	E		G	O	G
				P	R	E	S	S		M	E	N	U		

OXYGEN OFF

Do not use the welding unit in the presence of oxygen or any potentially explosive substance.

NO PACEMAKER

Make sure the patient and people near the welding unit/clamp connecting cable are not carriers of implanted electronic devices, that might be sensitive to the magnetic impulses generated during the welding phase.

PROTECTIVE GOGGLES

Ensure that both the patient and the user are wearing eye protection during the welding procedure.

Main Menu

To enter the main menu, press the [MENU] key.

S	M	A	R	T		P	R	O	G	R	A	M	S		
P	R	O	G	R	A	M	S								
S	E	R	V	I	C	E									
L	A	N	G	U	A	G	E								

In the menu you can select energy level setting in the SMART PROGRAM mode, change and/or select a welding program, verify/change operating parameters or language selection by pressing the [▲] and [▼] keys.

Select from the following settings:

LANGUAGE:

Select the display language.

SMART PROGRAMS:

Select the energy level for welding from three modes:

HIGH = 100% of energy with two current pulses
MEDIUM = 75% of energy with two current pulses
LOW = 60% of energy with a single current pulse

PROGRAMS:

Modify the energy level between 30% and 100% in steps of 5%.

Select one or two current pulses

SERVICE:

NOTE: Can only be accessed by authorized service and maintenance personnel. This menu is accessible only to service personnel authorised by the manufacturer.

Selecting the display language

- Select the LANGUAGE from the main menu by using the [▲] and [▼] keys.
- Confirm your selection by pressing the [ENTER] key.

The available languages are displayed.

D	E	U	T	S	C	H									
E	N	G	L	I	S	H									
I	T	A	L	I	A	N	O								
F	R	A	N	C	A	I	S								

NOTE: Scroll to see all the languages. Available languages: English, German, French, Italian, Spanish.

- Return to the previous menu with the [MENU] key.

User's Manual

Standard Welding Program "SMART"

- Select SMART PROGRAMS from the main menu by using the [▲] and [▼] keys.
- Confirm by pressing the [ENTER] key.
- Select the welding energy level from the three modes displayed.

S	M	A	R	T	S	E	T	T	I	N	G	H	I	G	H		
S	M	A	R	T	S	E	T	T	I	N	G	M	E	D	I	U	M
S	M	A	R	T	S	E	T	T	I	N	G	L	O	W			

HIGH = 100% of energy with two current pulses

MEDIUM = 75% of energy with two current pulses

LOW = 60% of energy with a single current pulse

- Use the [▲] and [▼] keys to select the energy level.
- Confirm your selection by pressing the [ENTER] key.

The selected energy level is displayed:

S	M	A	R	T	S	E	T	T	I	N	G			H	I	G	H
P	R	E	S	S	W	E	L	D									

- Return to the previous menu with the [MENU] key, or
- enable a welding cycle by pressing the [WELD] key.

Starting a welding cycle in the "SMART" mode

S	M	A	R	T	S	E	T	T	I	N	G			H	I	G	H
P	R	E	S	S	W	E	L	D									

- Enable the welding cycle by pressing the [WELD] key.

LED status:

READY	flashing (for eight seconds)
WELD	off
FAULT	off
WELD KEY	off

- To start welding, press the pedal switch within eight seconds.

NOTE: Three seconds after welding the next welding cycle will be enabled (once the READY LED is flashing).

- Press the pedal switch within eight seconds without pressing the [WELD] key again.

LED status:

READY	off
WELD	flashing
FAULT	off
WELD KEY	off

NOTE: After eight seconds, the READY LED switches off and the welding cycle is disabled.

- For the next welding press the [WELD] key again, and press the pedal switch within eight seconds.

NOTE: The display automatically returns to the previous menu when the welding phase is completed.

NOTE: If a malfunction occurs, the welding unit goes in stand-by-mode, showing the malfunction type on LCD (refer to the chapter ERROR MESSAGES).

Welding programs

- Select PROGRAMS from the main menu by pressing the [▲] and [▼] keys.
- Confirm by pressing the [ENTER] key.

S	E	L	E	C	T	P	R	O	G	R	A	M					
C	R	E	A	T	E	P	R	O	G	R	A	M					

- Choose SELECT PROGRAM by using the [▲] and [▼] keys.

NOTE: You can only choose a program you have already created before.

- Confirm your choice by pressing the [ENTER] key.

		S	E	L	E	C	T	P	R	O	G	R	A	M			
P	R	O	G	R	A	M									0	0	
N	O			O	F		C	Y	C	L	E					0	2
E	N	E	R	G	Y		L	E	V	E	L				0	3	0 %

- Select the number of the program (between 00 and 99) by using the [▲] and [▼] keys.
- Confirm the selected program number by pressing the [ENTER] key.

Welds can be made by selecting a one of the default (SMART) programs or by setting the values as required and by saving them as a new program.

- 

In both cases it is strongly recommended to test the chosen parameter settings extraorally in advance to verify whether they give good welding results, and to validate them.

- 

When using an excess power factor for a very thin wire, there is a risk of completely melting the metal and generating spits of molten metal in the patient's mouth.

- IMPORTANT:** The welding parameters of the welding unit must be set as follows based on the diameter of the titanium wire used.

Titanium wire diameter	Welding parameter setting on the welding unit		
1.2 mm	LOW	60%	single pulse
1.5 mm	MEDIUM	75%	double pulse
2.0 mm	HIGH	100%	double pulse

15

Troubleshooting

Error messages

There are two malfunction subsets: recoverable and non-recoverable errors.

A numerical code (#1-#8) is associated with every error, which is indicated in the display. At the same time, the [FAULT] lamp flashes.

A [FAULT] lamp that stays constantly on, while the welding unit is blocked, signals that independent hardware has put the welding unit in the safety mode because of a software failure, thus preventing any welding cycle from starting.

In this case, switch the welding unit off and on again.

Recoverable errors

Error Code	Possible root cause	Remedy
Clamp open #1	User opens the clamp during a welding process.	Press the [MENU] key as indicated on the display and start a new welding cycle.
Clamp open #2	User opens the clamp before starting a welding cycle and keeps the clamp open during welding.	Press the [MENU] key as indicated on the display and start a new welding cycle.
Pedal disconnected #5	Pedal switch defective or disconnected	Check if the plug of the pedal switch is properly connected to the socket. Correct if necessary and press the [MENU] key as indicated on the display.
Over temperature #4	Overheating of electronic switch due to abnormal, very fast welding sequences or a capacitor fault	Wait for electronic switch to cool, then press the [MENU] key and start a new welding cycle.

Non-recoverable errors

Error Code	Possible root cause	Remedy
Incorrect charge #3	Mosfet short circuit Hardware failure which can only occur during the first part of the welding phase. Electronic switch for capacitor charging may be defective.	Call the Technical Service.
Pedal defective #6	Pedal switch defective pedal or short circuit input	Call the Technical Service.
Incorrect charge #7	Time-out in charging phase	Call the Technical Service.

Correct use of the welding clamp


The welding clamp, supplied with the welding unit, is non-sterile upon delivery. It must be sterilized and processed before use.

IMPORTANT: Correct handling of the welding clamp is imperative to ensure the safety level necessary for the patient and operator.

- Make sure that the patient is positioned in such a way so that when in contact with the applied part (welding clamp), he/she does not come into contact with accessible grounded metallic parts of the unit or other nearby devices.
- NEVER disconnect the cable of the clamp from the welding unit while the clamp has contact with the patient.
- When the operator is in contact with the patient, he must avoid touching accessible grounded metallic parts of the welding unit or of any other nearby devices.
- The welding unit does not have functions requiring connection with other devices, so prevent its metallic parts including the applied parts from coming into contact with other devices.

A handle can be pulled out on the front of the welding unit, where the operator can hang up the welding clamp during treatment. The handle can be pushed back when not needed.

Positioning of the welding clamp

Make sure that the mucous membranes, soft tissues, suture thread or other materials are not trapped between the electrodes.

- Position the welding components centrally between the two electrodes of the welding clamp.

The clamping force of the welding clamp secures the components in position. The pressure exerted by the clamp electrodes on the welding point must be constant.

The pressure exerted by a spring that holds the clamp closed and prevents the formation of voids between the electrodes and the metal.

CAUTION: NEVER apply manual pressure on the two arms of the clamp during the welding operation. Imperfect contact on the joints to be welded can generate the formation of sparks during welding.

- Verify the correct positioning and stability before starting a welding cycle.
- Support the welding clamp with the palm of your hand during welding and DO NOT apply manual pressure on the two arms of the clamp during the welding operation under any circumstances.
- Do not move nor remove the clamp from its position until the welding is completed.

A complete welding cycle takes between 0.5 and 5 seconds and is indicated by the flashing **WELD** lamp.

Cleaning, disinfection and sterilization

All components are supplied non-sterile!
The welding unit cannot be sterilized. Only the welding clamp can be sterilized and has to be processed prior to first use.

IMPORTANT: All components that might come into contact with the patient during the treatment, in particular the clamp and the electrodes must be cleaned, disinfected and sterilized before use on the patient using a validated procedure.

Cleaning and disinfection of the welding unit

WeldOne™ housing

- No regular cleaning intervals envisaged.
- Use a dry, clean cloth to clean the welding unit for routine cleaning. Use a damp cloth without detergents for more thorough cleaning.
- Do not use any abrasive substances or solvents for cleaning the housing.
- Do not use any liquid or spray directly on the panel or near the electrical sockets.

Clamp connecting cables

- Clean and disinfect by wiping with disinfectant.
- For working under sterile conditions it is recommended to cover the clamp connecting cables with a disposable sterile sleeve.

Reprocessing of the welding clamp

The reprocessing of the welding clamp was validated with the following equipment and materials:


- Washer-disinfector G 7836 CD, Miele & Cie. KG, Gutersloh, Germany
- Cleaning program VARIO-TD (with thermal disinfection)
- Cleaning agent Neodisher MediClean Forte, neutralization agent Neodisher Z, Dr. Weigert, Hamburg, Germany

- **Ultrasonic bath:** Make sure that the components are fully immersed so that all inner and outer surfaces are wetted.
- For automated cleaning use appropriate washer- disinfectors and cleaning agents recommended and prescribed as by the manufacturer of the washer- disinfecter.
- Observe the instructions of the manufacturer of the cleaning and neutralization agents regarding concentration, exposure time and, if applicable, temperature.

For cleaning and sterilization use the clamp retainer to hold the arms of the welding clamp open so that the contact surfaces of the welding electrodes are exposed to the cleaning agent and thermal disinfection.

Manual pre-cleaning

Remove severe contamination (e.g. blood, saliva) from the welding clamp immediately after use.

- Immerse the welding clamp with the cleaning agent and remove residues thoroughly with a soft nylon brush until visually clean.

Ultrasonic pre-cleaning

- Sonicate the welding clamp in an ultrasonic bath filled with cleaning agent for at least 5 minutes.

Automatic cleaning

- Place the welding clamp in the cleaning basket of the washer-disinfector.
- Use the VARIO-TD cleaning program that is capable of removing blood at a temperature of 45 - 55 °C.
- Final rinsing should take place with deionized water, or water that corresponds to that degree of purity.

Preparation for sterilization

- Wrap the welding clamp in a sterilization pouch (single or double) appropriate for steam sterilization sealed with strips, and/or a sterilization container.
- Every sterilization pouch must have a sterilization indicator and the sterilization date.

Sterilization

- Sterilize in the autoclave with fractionated pre-vacuum at 132 °C for 3 minutes (minimum holding times; the operating times are longer and may vary depending on the device).

Storage

Store the sterilized components dry and dust-free at room temperature.

Reprocessing and sterilization of the welding clamp was validated according to the following standards:

- DIN EN ISO 17664 Sterilization of medical devices - Information to be provided by the manufacturer for the processing of resterilizable medical devices
- Guideline of the German Association for Hospital Hygiene (DGKH), German Society for Sterile Supply (DGSV) and the Workgroup for Reprocessing of Instrumentation (AKI) on the validation and routine control of machine cleaning and thermal disinfection processes for medical devices and on basic principles of equipment selection
- DIN EN ISO 15883-1 Washer-disinfectors – Part 1: General requirements, terms and definitions and tests

Maintenance**General maintenance instructions**

The welding unit does not contain any user-serviceable parts. Therefore there is no need for the operator to get access to the components inside the housing.

Only service personnel authorized by the manufacturer may open the housing to get access to the internal components of the welding unit for maintenance and repair work

If an unauthorized person repairs or modifies the welding unit without authorization, he assumes total responsibility for any damage to the unit or injury to the patient, the operator, or third parties.

Routine maintenance

In order to maintain the safety level for the patient and operator, check the state and wear of the components at regular intervals.

Perform routine maintenance at least every 12 months.

1. Disconnect the power supply cable and check the state to be sure there are no signs of wear or breakage
2. Disconnect the mains cable of the pedal switch and check that there are no signs of wear or breakage.
3. Check the conditions of the clamp connection cables with both the unit and the clamp. The connection must be stable and free of slack.
4. Check the conditions of the cables to make sure there are no signs of wear or breakage and that the protective part made of insulating material shows no breakage
5. Check the condition of the copper electrodes, which must be replaced after 300 welding applications or after one year of use.
6. The recommended replacement of the discharge capacitor is after 2-3 years.

Replace the entire accessory when a part is damaged.

Contact the service center for any aspect that should raise doubts concerning safety and reliability of the welding unit.

Maintenance by the user is restricted to the following:

- Replacement of the copper electrodes
- Replacement of fuses

User's Manual

Periodic maintenance

CAUTION: The efficiency of the copper electrodes is essential for safety and must be kept under control.

The clamp is shipped with the electrodes already inserted and in full contact.

When the contact flatness deteriorates, it is necessary to use abrasive paper interposed between them, and if this proves to be insufficient for restoring adequate contact, it is necessary to replace the electrodes using the replacement tool supplied with the welding unit.

- The copper electrodes must be replaced after 300 welding applications or after one year of use.
- The recommended replacement of the discharge capacitor is after 2-3 years.

Replacing the copper electrodes


The welding unit is delivered with a set of replacement electrodes.

WARNING: Switch off the welding unit and unplug mains cable before replacing the copper electrodes.

To replace the electrodes proceed as follows:

1. Switch off the welding unit and unplug the mains cable.
2. Turn the locking rings of the welding clamp connecting cable counterclockwise until the clamp is disconnected from the clamp connecting cable.
3. Keep the clamp open by pressing the clamp arms together.

4. Unscrew the electrodes by hand or use the electrode replacement tool one by one and remove them
5. Insert the new electrodes by hand and tighten them firmly with the electrode replacement tool.

CAUTION: After replacing the copper electrodes make sure that they are securely fixed to prevent swallowing or aspiration of the copper electrodes by the patient.


Replacing the fuses


- **WARNING:** Before exchanging the fuses, switch off the welding unit and unplug the mains cable.
- When exchanging the fuses, only use replacement fuses with the same specification. For correct technical specifications, please refer to the chapter *Technical data*.

The two mains fuses are located inside the fuse drawer. To replace the fuses proceed as follows:

- Switch off the welding unit.

Unplug the mains cable from the socket and from the welding unit plug

Place the tips of the two screwdrivers into the cut-outs on both sides of the fuse drawer.

Push very carefully, using the screwdrivers as levers, to remove the fuse drawer from the housing.

Remove the old fuses and insert the new fuses of the same technical specification into the fuse holders.

Reinsert the fuse drawer into the housing and push slightly to lock it in place.

- Reinsert the main cable into the welding unit plug and into the socket. Switch on the welding unit.

Decommissioning

Change of location

The handling and the movement of the device must be performed taking appropriate precautionary measures.

The manufacturer assumes no liability in case of non-fulfillment or non-observance of the instructions given below.

Take the following precautions whenever it is necessary to move the welding unit:

- Make sure that the welding unit is switched OFF by checking the ON/OFF switch on the right side of the welding unit.
- Disconnect all cables connected to the welding unit.
- The welding unit may be transported in the original packaging from the manufacturer only.
- When not used, store the welding unit in a dry place that fulfills the storage conditions specified at *Technical data*.
- Adhere to the instructions at *Initial setup*.

Disposal at the end of the unit's lifetime

The unit contains electric and electronic components that may prove harmful for the environment and health unless properly disposed of when the unit is dismantled at the end of its lifetime.

The following symbol is shown on the type plate on the unit:

Implementation of Directives 2002/96/EC and 2003/108/EC regarding reducing the use of hazardous substances in electric and electronic equipment, and waste disposal:

- The symbol of the crossed-out dumpster placed on the unit means that the product must be collected separately from the other waste at the end of its lifetime. Therefore, this unit cannot be disposed of like normal urban waste.
- The user wishing to dispose of this used or dismantled unit can contact one of the WEEE (Waste Electrical and Electronic Equipment) collection centers in his area of residence that is authorized and set up by the public administration, or deliver the unit to his dealer at the time he takes delivery of a new one.
- Appropriate separate collection for the subsequent sending of the dismantled unit to recycling, processing and eco-friendly disposals contributes to preventing possible negative impacts on the environment and health, and helps materials forming the unit to be reused and recycled.
- Abusive disposal of the product by the holder involves application of the administrative sanctions provided for by current regulations.
- Please refer to the above-mentioned regulations for more detailed information.

The shelf life and life cycle of the Weld One welding machine and its accessories is 10 years, with the exception of the copper electrodes whose life cycle is 300 welding cycles or one year of operation.

The recommended replacement of the discharge capacitor is after 2-3 years.

Technical Support

Studioemme S.a.s is determined to provide to their customers the best after-sales support possible. The after-sales service aims at solving problems in the shortest time possible and at avoiding unnecessary returns for repair.

The **Service Hotline** is available at:

Email: aftersales@studioemme.it

Phone: +39 (0) 51 6415239

Address: STUDIOEMME S.a.s.

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40132 Bologna, Italy

If the problem cannot be solved, instructions will be given for return shipment and repair.

Warranty and repair

Warranty

Studioemme S.a.s. warrants that the welding unit and its accessories are free of defects at the time the customer purchases the product, and therefore any malfunctioning that should arise during the 24 months following the date of purchase will be covered by warranty.

The warranty covers:

1. The costs of the components proved to be defective
2. The labor necessary for their replacement and testing of the functions
3. Shipping cost for the unit from the customer to Studioemme and back again.
4. Free-of charge loaner unit for the time of the repair.

The warranty **does not** cover:

1. Damage caused by improper use
2. Incorrect maintenance or repair attempts carried out by personnel not authorized by Studioemme S.a.s
3. Non-observance of the instructions in the User's Manual
4. Replacement of components subject to wear listed in the User's Manual
5. Damage due to transport; the customer must always check that the goods delivered are free of damage.
6. Any compensation for direct or indirect personal injuries or material damage of any kind.

Studioemme S.a.s. does not consider itself liable, nor does it provide for any type of indemnity to the customer for everything that results from the non-use of the welding unit for the time elapsing between the moment the defect arises and the return of the welding unit after repair.

Repair

Any defect detected by the customer shall be reported to Studioemme S.a.s by using the Repair Request Form included in the shipment.

Studioemme S.a.s undertakes to check the request and to reply on whether the warranty covers the repair and how to proceed with the repair.

Studioemme S.a.s. will send the customer a repair authorization in writing including the shipping address for the defective welding unit.

Appendix

Electromagnetic Compatibility (EMC)


The welding unit complies with current electromagnetic compatibility and electrical safety regulations but, like for all electronic devices, special precautions must be taken with regard to electromagnetic interference (EMI) and this is why it must be installed and put into service in compliance with the EM information of this document.

WARNING: Portable radio communication devices may affect unit operation.

The welding must not be used in contact with or in the proximity of other electronic devices

The cables connected to the unit have the following arrangement, characteristics and lengths:

Do not use cables and accessories other than those specified, and use only those supplied by the manufacturer.

The use of other cables may cause undesirable malfunctioning and behaviors such as increased emission, greater susceptibility to interference and therefore reduced electrical safety.

Use

Tables

The We

The use

Table

Guide
The W
The u
Emis
RF er
CISP
RF er
CISP
Harm
IEC 6
Flicke
IEC 6

Tables 1, 2, 4, 6 provided for by standard 60601-1/2 (2010)

The WeldOne™ welding unit is designed to operate in the electromagnetic environment specified in the tables below. The user must ensure that it is used in conformity.

Table 1

Guide and manufacturer's declaration - electromagnetic emissions		
The WeldOne™ medical unit is designed to operate in the electromagnetic environment specified below. The user of the unit should ensure that it is used in this environment.		
Emission test	Conformity	Electromagnetic environment
RF emissions CISPR 11	Group 1	The welding unit uses RF energy only for its internal operation. As a result, its RF emissions are very low and in all likelihood do not cause any interference in electronic devices nearby.
RF emissions CISPR 11	Class B	The welding unit is suitable for use in all environments, including domestic environments and those directly connected to a low voltage public power source that supplies power to buildings used for domestic purposes.
Harmonic emission IEC 61000-3-2	Class A	
Flicker emissions IEC 61000-3-3	Compliant	The welding unit is suitable for use in all environments, including domestic environments and those directly connected to a low voltage public power source that supplies power to buildings used for domestic purposes

Table 2

Immunity test	IEC 60601 test level	Level of conformity	Electromagnetic environment
Electrostatic discharge (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kv air	± 6 kV contact ± 8 kv air	Floors must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity must be at least 30%
Transients / sequence of electric pulses IEC 61000-4-4	± 2 kV for power supply lines ± 1 kv for input/output lines	± 2kV for power supply lines	The grade of the mains voltage should be that of a typical commercial or hospital environment.
Surges IEC 61000-4-5	± 2 kV common mode ± 1 kV differential mode	± 2 kV common mode ± 1 kV differential mode	The grade of the mains voltage should be that of a typical commercial or hospital environment.
Voltage flicker, brief cut-off and voltage variations on the power supply input lines IEC61000-4-11	< 5% U_t (95% of hole) for 0.5 cycles 40% U_t (60% of hole) for 5 cycles 70% U_t (30% of hole) for 25 cycles < 5% U_t (95% of hole) for 5 seconds	< 5% U_t (95% of hole) for 0.5 cycles 40% U_t (60% of hole) for 5 cycles 70% U_t (30% of hole) for 25 cycles < 5% U_t (95% of hole) for 5 seconds	The grade of the mains voltage should be that of a typical commercial or hospital environment. It is recommended to power the device with a UPS or batteries if the welding unit user requires continuous operation even during a black-out.
Magnetic field at mains frequency (50/60 Hz) IEC6100-4-8	3 A/m	3 A/m	The magnetic fields at mains frequency must have levels characteristic of a typical commercial or hospital environment.

Note: U_t is the mains AC power supply before the test level application.

Table 4

<i>Immunity test</i>	<i>IEC 60601 test level</i>	<i>Level of conformity</i>	<i>Electromagnetic environment</i> <i>A device for mobile RF communication should not be closer to the WeldOne™ welding unit (including cables) than what is specified as calculated distance according to the equations below at the different frequencies</i>
conducted RF IEC 61000-4-6	3 Vrms 150Khz – 80Mhz	3 Vrms	Recommended distance between a transmitter and the device: $d = 1,167 \cdot \sqrt{P}$ $d = 1,167 \cdot \sqrt{P}$ 80Mhz – 800Mhz $d = 2,333 \cdot \sqrt{P}$ 800Mhz – 2.5Ghz where P is the maximum nominal output power of the transmitter in W according to the transmitter manufacturer and d is the separation distance in meters. An electromagnetic interference may be present near a device marked with the symbol: 
irradiated RF IEC 61000-4-3	3V/m 80Mhz at 2.5Ghz	3 V/m	

Note: This guideline is not applicable in all situations. Electromagnetic propagation depends on the absorption and reflection of structures, objects and persons.

Table 6

The WeldOne™ welding unit is designed to operate in an electromagnetic environment where the interferences irradiated at RF are under control. The operators that use the WeldOne™ welding unit can contribute to preventing electromagnetic interferences by ensuring a minimum distance between both mobile and portable RF communication devices (transmitters) and the device as recommended below in connection to the maximum output power of the radio communication devices.

Maximum output power of the Transmitter specified (W)	Separation distance at the frequency of the transmitter (m)		
	From 150 kHz to 80MHz $d=1.17*\sqrt{P}$	From 80MHz to 800MHz $d=1.17*\sqrt{P}$	From 800Mhz to 2.5GHz $d=233*\sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.370	0.370	0.740
1	1.17	1.17	2.33
10	3.70	3.70	7.40
100	11.7	11.7	23.3

Note: For transmitters with maximum power not in the above list, the separation distance d in (m) can be calculated using the equation applicable to the transmitter frequency, where P represents the maximum output power of the transmitter in watts, in agreement with the indications of the transmitter manufacturer.

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