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Johnson & Johnson International, Belgium

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Associate Director, Regulatory Affairs

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
М.П. / Stamp

Instruction for use

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles

2022

Page 1 of 6

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles

DESCRIPTION

Monocryl Plus Antibacterial sterile synthetic absorbable surgical suture with atraumatic needles is a sterile, synthetic, absorbable, monofilament surgical suture manufactured from a copolymer of glycolide and ϵ -caprolactone. The empirical molecular formula of the polymer is $(C_2H_2O_2)_m (C_6H_{10}O_2)_n$. Poliglecaprone 25 copolymer has no antigenic or pyrogenic properties and elicits only a slight tissue reaction during absorption. Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles contains Irgacare[®] MP (triclosan), a broad-spectrum antibacterial agent, at no more than 2360 $\mu\text{g}/\text{m}$.

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles is available dyed (with D&C Violet No. 2, color index 60725) and undyed. Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles may have various gauge sizes and lengths and is supplied with stainless steel needles of various types and sizes. The needles may be attached permanently or as control release (CR) ones, enabling the needles to be pulled off without cutting the suture off. Please see catalogue for detailed description of product range.

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles complies with all the requirements of the standards of European Pharmacopoeia for Sterile Synthetic Absorbable Monofilament Sutures and United States Pharmacopoeia for Absorbable Surgical Sutures (except for a slight oversize in diameter).

CONDITIONS FOR USE

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles is intended for professional use in medical and healthcare institutions.

TECHNICAL CHARACTERISTICS

Needle color- Silver/Black

Part of Device	Material
Monocryl Plus suture	copolymer of glycolide and ϵ -caprolactone (Poliglecaprone 25 copolymer) D&C violet No.2 (Colour Index Number 60725) – only for dyed sutures Irgacare MP (triclosan)
Suture needles	420 Stainless Steel 4310 Stainless Steel ETHALLOY
Suture needle coating	Silicone, CEREBRUS, MULTIPASS

INDICATIONS

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles is indicated for soft tissue approximation and/or ligation, where an absorbable material is indicated.

APPLICATION

Sutures should be selected and implanted considering the patient's condition, the surgeon's experience, the surgical technique and the wound size.

PERFORMANCE AND ACTION

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles elicits a minimal initial inflammatory reaction in tissues with gradual encapsulation of the suture by newly formed fibrous connective tissue. Progressive loss of tensile strength and eventual absorption of Monocryl Plus sterile synthetic absorbable surgical

suture with atraumatic needles occurs by means of hydrolysis, where polymer degrades to adipic acid, which is subsequently absorbed and metabolized in the body. Absorption begins as a loss of tensile strength followed by a loss of mass. Below are results of implantation studies of the suture in rats:

DYED SUTURE		UNDYED SUTURE	
After implantation	Approximate % of original strength remaining	After implantation	Approximate % of original strength remaining
7 days	60%	7 days	50%
14 days	30%	14 days	20%
Initial tensile strength is essentially lost after 28 days of implantation.		Initial tensile strength is essentially lost after 21 days of implantation.	

Absorption is essentially complete between 91–119 days.

Using in vitro studies, Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles inhibits colonization of the suture by *Staphylococcus aureus*, *Staphylococcus epidermidis*, *Escherichia coli*, *Klebsiella pneumoniae* and *Enterobacter cloacae* microorganisms, as well as methicillin-resistant *S. Aureus* and *S. Epidermidis*. Clinical significance of these data has not been studied.

CONTRAINDICATIONS

Since Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles (dyed and undyed) is absorbable, it should not be used where prolonged approximation of tissues under stress is required. In particular, undyed Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles is contraindicated for fascial tissues closure.

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles should not be used in patients with known allergic reactions to Irgacare† MP (triclosan).

WARNINGS/PRECAUTIONS/INTERACTIONS

The safety and effectiveness of Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles have not been established in the following areas: neurosurgery, cardiovascular surgery, microsurgery, and ophthalmic surgery.

Users should be familiar with surgical procedures and techniques involving absorbable sutures before employing Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles for wound closure, as risk of wound dehiscence may vary with the site of application and the suture material used. When selecting a suture, surgeons should consider its interactions with the body (see PERFORMANCE AND ACTION section). This suture is not recommended for use in elderly, malnourished or debilitated patients, or in conditions which may delay wound healing.

As with any foreign body, prolonged contact of any suture with salt solutions, such as those found in the urinary or biliary tracts, may result in calculus formation. Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles suture may act transiently as a foreign body.

In case of wound contamination or infection, accepted surgical practice should be followed.

When working with absorbable sutures, the use of supplemental nonabsorbable sutures should be considered by the surgeon for closure of sites which may undergo expansion, stretching or distension, or that require additional support.

Skin sutures which remain in place longer than 7 days may cause localized irritation and should be cut off and removed. Subcuticular sutures should be placed as deeply as possible to minimize the erythema and induration, which are usually associated with the absorption process.

Under certain conditions, notably orthopedic procedures, immobilization of joints by external support may be employed at the discretion of the surgeon. Care should be taken in the use of absorbable sutures in tissues with poor blood supply as suture rejection and delayed absorption may occur.

When handling this or any other sutures, damage to the suture should be avoided. Avoid deformation of suture structure with surgical instruments, such as forceps or needle holders.

Adequate knot security requires the standard surgical technique of flat and square ties with additional throws if indicated by surgical indications and the experience of the surgeon. Additional throws may be particularly appropriate when using monofilament sutures.

Care should be taken to avoid damage when handling surgical needles. The needle should be grasped in an area one-third (1/3) to one-half (1/2) of the distance from the attachment end to the point.

Grasping close to the point could impair the penetration performance and cause fracture of the needle. Grasping the needle at the attachment end may cause needle deformation or breakage. Reshaping a needle may cause it to lose strength, which may result in the needle's lower resistance to deformation and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in "Sharps" container.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or cross-contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users.

ADVERSE REACTIONS

Adverse reactions associated with the use of synthetic absorbable sutures include wound dehiscence, failure to provide adequate wound tissue support in closure of the sites where expansion, stretching, or distension occur, failure to provide adequate tissue support in elderly, malnourished or debilitated patients, as well as in conditions which may delay wound healing, minimal acute inflammatory tissue reaction, transitory local irritation or transitory inflammatory reaction to a foreign body at the wound site, suture extrusion, delayed absorption in tissues with poor blood supply, calculi formation in urinary and biliary tracts when prolonged contact with salt solutions such as urine and bile occurs, as well as allergic reaction to Irgacare[®] MP (triclosan). Like any other foreign body, Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles may potentiate an existing infection. Broken needles and excessive traction when pulling the needle may result in needles or fragments remaining in the wound, or may require extended or additional surgeries, or lead to residual foreign bodies. Inadvertent needle-sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

STERILITY

Monocryl Plus sterile synthetic absorbable surgical suture with atraumatic needles is sterilized by ethylene oxide. Do not resterilize. Do not use if package is opened or damaged. Discard opened, unused sutures.

TRANSPORTATION AND STORAGE

The product is transported by any mode of transportation in covered vehicles with protection against atmospheric precipitation, in accordance with the regulations in force for the mode of transportation.

While transported, the product should not be subjected to long-term exposure to heat sources, direct sunlight.

Recommended storage conditions: -60°C to +30°C.

Recommended transportation conditions: at temperature from -60°C to + 50 °C for no more than 288 hours.

EXPIRATION DATE

The expiration date of the product is indicated on the product packaging and is 24 months. Do not use after expiration date.

WARRANTY

The manufacturer guarantees to the end user absence of defects and conformity to product specifications until the

expiration date under the conditions of adherence to the storage, transportation and package integrity requirements.

DISPOSAL

After use, sutures are a potential biohazard due to contamination with blood and/or other biological substances. Disposal should be carried out according to established medical regulations, as well as the current legislation of the country of disposal.


In the Russian Federation, disposal regulations for B-class waste are applied.

HOW SUPPLIED

Monocryl Plus suture is made in various sizes. Monocryl Plus sutures are wound into Plastic Trays with paper lid. Every suture is packed individually into a sterile foil package (primary package). Sutures in individual primary packages are then placed into a cardboard box (secondary package), which may contain 12 or 36 items. Boxes are wrapped with film, using tear tape to facilitate opening. Booklets, dispenser boxes and instruction leaflets contain printed information.

SYMBOLS USED ON LABELING

 Catalogue number

 Sterilized using ethylene oxide

 CE-mark and identification number of notified body. The product meets the essential requirements of Medical Device Directive 93/42/EEC.

 Batch code

 Do not reuse

 Do not resterilize

 Use-by date

 Manufacturer

 Caution

 Number of Units

 Non-pyrogenic

 Do not use if package is damaged

 See Instructions for use

Date of manufacture

† Registered Trademark of BASF Group

Johnson & Johnson International, Belgium
Johnson & Johnson International
c/o European Logistics Centre
Leonardo Da Vincilaan, 15
BE 1831 Diegem
Belgium
1-877-ETHICON
+1-513-337-6928

Authorized representative of manufacturer:

Johnson & Johnson LLC,
121614, Moscow, Krylatskaya St.,
17/2
Phone: 8 (495) 580 77 77,
Fax: 8 (495) 580 78 78

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