The State of Maryland

Office of the Secretary of State

This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille does not certify the content of the document for which it is issued.

Apostille

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1. Country: United States of America
This public document
2. has been signed by **Robert P. Duckworth**
3. acting in the capacity of **Clerk of the Circuit Court for Anne Arundel County**
4. bears the seal/stamp of the **Circuit Court for Anne Arundel County**

Certified

5. at Annapolis, Maryland
6. the **30th day of October, 2017**
7. by The Secretary of State of Maryland
8. **No. 439454**
9. Seal

10. Signature

Secretary of State

State of Maryland, Anne Arundel County, Sct.

I, **ROBERT P. DUCKWORTH**, Clerk of the Circuit Court for Anne Arundel County, Maryland, a court of record, do hereby certify that **Denise Mandel** was a commissioned/appointed and qualified Notary Public commencing on the 28th day of February, 2014.

In Testimony Whereof, I have hereunto set my hand and affixed the seal of the court October 30, 2017

Robert P. Duckworth

Robert P. Duckworth, Clerk
Circuit Court for Anne Arundel County, Maryland

Affidavit

State of Maryland)
)SS:
County of Anne Arundel)

Before me, the undersigned notary public, this day, personally,
appeared to me known or proven, who being duly sworn according to law,
deposes the following:

Sundara Chintaluri

I, _____, hereby certify that the attached
document is

- ☒ an original document
☐ a true copy of an original document

to the best of my knowledge and belief.

Signed

Subscribed and sworn to before me this

30 day of October, 2017

Notary Public

DENISE MANDEL
Notary Public-Maryland
Anne Arundel County
My Commission Expires
April 02, 2018

25 October 2017

Bayer Medical Care Inc

1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.

(412) 767-2400

www.radiology.bayer.com

To whom it may concern,

To be submitted to the Competent Health Authorities in Federal Republic of Russia

We, Bayer Medical Care Inc., USA (1 Bayer Drive Indianola Pennsylvania 15051 USA) hereby state that attached document (Instruction for use) (Disposable Hand Controller MEDRAD VFlow: VF HC) is true.

Sincerely,

Brenda Fanala
Principal Regulatory Affairs Specialist

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©2017 Bayer. All Rights Reserved. Reproduction of this manual is strictly prohibited without express written consent of Bayer Medical Care Inc.

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Bayer, the Bayer Cross, MEDRAD, Mark 7 Arterion, MEDRAD Mark 7 Arterion, VFlow, and MEDRAD VFlow may be registered trademarks of Bayer in the US and other countries.

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©2013-2014, 2017 Bayer. All Rights Reserved. Reproduction of this manual is strictly prohibited without express written consent of Bayer Medical Care Inc.

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Bayer, the Bayer Cross, MEDRAD Mark 7 Arterion and MEDRAD VFlow may be registered trademarks of Bayer in the US and other countries.

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UNITED STATES

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
Phone: +1 (412) 767-2400
+1 (800) 633-7231
Fax: +1 (412) 767-4120

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




Introduction: Read the information contained in this brochure. Understanding the information will assist you in using the **MEDRAD VFlow Sterile Hand Controller** properly.

Important Safety Notice: The MEDRAD VFlow Sterile Hand Controller is intended to be used by medical professionals with adequate training and experience in the operation of the MEDRAD® Mark 7 Arterion Injection System and angiographic procedures and techniques.

Intended Use: The MEDRAD VFlow Sterile Hand Controller is specifically intended for use with the MEDRAD Arterion Injection System for angiographic studies.

Contraindications: None known.

Restricted Sale: Rx Only.

Environmental Contamination Hazard - Serious patient or worker injury or death may result.

- Do not use if sterile package is opened or damaged. Patient or operator injury may result if package is opened or damaged, or if damaged components are used. Visually inspect contents and package before each use.

Cross Contamination Hazard - Serious patient and/or worker injury or death may result.

- For devices labelled for single use, please note: This product is intended for single use only. Do not resterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

Environmental Contamination Hazard - Minor or moderate patient and/or worker injury may result.

- Properly discard disposables after use, in accordance with hospital hazard waste disposal procedures. Do not reuse disposables.

NOTICE


Electro-Mechanical Hazard - Equipment Damage may result.

- Do not use tools to over tighten connections or to assist in the removal of disposables.

NOTE: Refer to the MEDRAD Mark 7 Arterion injector operation manual for a full set of instructions, warnings and precautions.

The MEDRAD VFlow Hand Controller is a sterile device intended for single patient use. The Hand Controller works for Fixed Flow Rate and Variable Flow Rate injections. When in the Variable Flow Rate injection mode, the flow rate increases incrementally as the Hand Controller plunger (A) is depressed, and decreases as the Hand Controller plunger is released. In the Fixed Flow Rate injection mode, the Hand Controller acts as a start switch, and release of the device stops the injection.


The Hand Controller button (B) is non-functional and will beep from the Injector head and Display Control Unit when depressed.

Hand Controller

Hand Controller Installation

1. Ensure that the injector is in the disarmed state.
2. Open the hand controller pouch using sterile technique.
3. Using sterile technique, remove the hand controller from the pouch.
4. Plug the hand controller into the underside of the Injector Head (see drawing below). Listen for an audible "click" to confirm proper connection.

NOTE: If the hand controller is damaged or does not work properly, discontinue use and discard.

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