
BAUSCH + LOMB

See better. Live better.

1400 North Goodman Street
Rochester, NY 14609
585.338.6000
www.bausch.com

Certified True Copy Letter

I, Isabelle B. Lefebvre, hereby swear (or affirm) that the attached Instruction for Use of Medical Device -IOL injector INJ100 is a true and complete photocopy of the document.


Isabelle B. Lefebvre, MSc.RA
Vice President, Regulatory Affairs
Bausch & Lomb Incorporated

09 January 2019
Date

State of New Jersey

County of Union

Subscribed and sworn to (or affirmed) before me on this 9th day of January
2019, by Isabelle B. Lefebvre, proved to me on the basis of satisfactory
evidence to be the person(s) who appeared before me.

Signature of Notary Public

INSTRUCTION FOR USE OF MEDICAL DEVICE

IOL Injector, INJ100

Design facility:

Bausch & Lomb Incorporated, USA

1400 North Goodman Street Rochester, NY 14609, USA

Manufacturer:

Bausch & Lomb Incorporated, USA

1400 North Goodman Street Rochester, NY 14609, USA

Authorized Representative of Manufacturer:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

The legal entity addressed to whom the registration certificate can be issued:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

Addresses of manufacturing facility of medical device:

1. 1400 North Goodman Street Rochester, NY 14609, USA
2. Medice AG, Dornierstrasse 11, CH – 9423 Altenrhein, Switzerland

DESCRIPTION AND THE DESIGN OF MEDICAL DEVICE

IOL Injector, INJ100 (*hereinafter, IOL is intraocular lens*) is indicated for folding and injection of Bausch + Lomb IOLs (enVista™, model MX60) or other Bausch + Lomb IOLs approved for use with this injector. The inserter consists of a syringe-shaped body and tip with a plunger and an integrated cartridge. The INJ100 is a sterile, disposable, single-use device.

The construction (design) of the medical device is introduced below:


Fig. 1.

How supplied:

Each **IOL Injector, INJ100** is packaged in individual blister package, then injectors INJ100 in individual packages are packed into the carton secondary (group) package of 10 pieces.

Each carton secondary (group) package with injectors INJ100 also includes:

1. Instruction for use of medical product "IOL Injector, INJ100" – 1 piece;
2. Sticker "IOL Injector, INJ100" into the patient's chart – 10 pieces.

INDICATIONS FOR USE

The Bausch + Lomb IOL Injector is indicated for folding and injection of Bausch + Lomb IOLs (enVista™, model MX60) or other Bausch + Lomb IOLs approved for use with this injector. The inserter consists of a syringe-shaped body and tip with a plunger and cartridge. The INJ100 is a sterile, disposable, single-use device.

The target patient populations include those patients who need an IOL placed during cataract replacement surgery.

The Bausch+Lomb™ IOL Injector lens injection system is a single use, disposable lens delivery system that supports the delivery of the following enVista and Akreos intraocular lens models:

CONTRAINDICATIONS

There are no contraindications for the INJ100.

ADVERSE EVENTS

1. IOL haptic and/or optic damage.
2. Anaphylactic Shock.
3. Microbial Infection.
4. Posterior Capsule Rupture.
5. Corneal Endothelium Damage.
6. Zonules, Torn, Weakened.
7. Elevated Intraocular Pressure.
8. Enlarged Incision Size/ Inadvertent Tissue Damage.
9. Surgical Procedure Impact, Delay (<30 min).

Adverse reaction reporting

Adverse reactions and/or potentially sight-threatening complications that maybe reasonably regarded as Akreos Single Use Insertion Device related and that were not previously expected in nature, severity or degree of incidence should be reported to your customer service affiliate.

MATERIALS USED IN MEDICAL DEVICE CONSTRUCTION AND ITS COMPONENTS

The designation of the medical device or component of it	Material	Patient contact	Sterilization type
IOL Injector, INJ100	Injector body: Denka CL-301 clear ABS; Plunger body: Denka CL-301 blue ABS; Plunger tip: Wacker 8304 blue silicone, 40 durometer; Cartridge loading chamber: Total PPR10222 clear polypropylene with 0.3% GMS;	Short-term direct contact with inner environment of organism (upto 3 – 5 minutes)	Sterile (sterilized using ethylene oxide), SAL 10 ⁻⁶

	Cartridge tip: Arkema PEBAX 7033 clear polyamide; Coating (applied to tip only): Medicoat A: polyurethane-based, water-soluble hydrophilic coating. Покрытие Medicoat A coating does not have antibacterial effect (confirmed by the letter from the manufacturer, the substance of coating is highly confidential information (confirmed by the letter from the manufacturer). The quantity of Medicoat A coating: thickness – 3 mkm, volume – 2.016 mm³ Spring: 302 stainless steel: Fe: base C: 0.15 % Max Cr: 17 – 19 % Ni: 8 – 10 % Si: 1.0 % Max Mo: 1.0 % Max Mn: 2.0 % Max P: 0.045 % Max S: 0.03 % Max)		
Blister package	Tray: Eastar Copolyester 6763; Lid: Tyvek ®; Typing on the lid: ink (mark and substance is highly confidential information (confirmed by the letter from the manufacturer); The connection of the tray and the lid: The lid is heat sealed to the tray	Short-term direct contact with the skin of the medical personnel in the case of working without gloves.	Sterile (sterilized using ethylene oxide), SAL 10⁻⁶

THE MAIN GUARANTEED CHARACTERISTICS OF MEDICAL PRODUCT

All the dimensions are introduced on the original manufacturing drawings:




Fig. 2. Original manufacturing drawing of the cartridge tip of IOL Injector, INJ100 (dimensions are introduced in millimeters).

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Fig. 2. Original manufacturing drawing of the cartridge loading chamber of IOL Injector, INJ100 (dimensions are introduced in millimeters).

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Fig. 3. Original manufacturing drawing of the Spring of IOL Injector, INJ100 (dimensions are introduced in millimeters).


Fig. 4. Original manufacturing drawing of the Plunger tip of IOL Injector, INJ100 (dimensions are introduced in millimeters).


Fig. 5. Original manufacturing drawing of the Plunger of IOL Injector, INJ100 (dimensions are introduced in millimeters).


Fig. 6. Original manufacturing drawing of the Injector body of IOL Injector, INJ100 (dimensions are introduced in millimeters).

Shelf life of IOL Injector, INJ100 is 24 (twenty four) months from the sterilization date.

INSTRUCTIONS FOR USE

1. Open the blister tray and transfer the sterile Bausch + Lomb IOL Injector into the sterile field.
2. Hold the injector in such a manner that the moveable and fixed cartridge wings can be guided with the index finger and thumb. Open the loading chamber (Fig. 8).

Fig. 8.

3. Entering from the loading chamber side, apply Bausch + Lomb viscoelastic directly into the conical tip and then apply two thin lines into the lateral grooves within the loading chamber (Fig. 9).


Fig. 9.

4. Advance the plunger tip to the outer edge of the of the cartridge wings (Fig. 10).


Fig. 10.

5. Position the lens anterior side up (reverse “S” orientation) in the middle of the loading chamber (Fig. 10). With non-serrated forceps, use slight pressure to push the lens and haptics down to ensure they are properly seated under the grooves (Fig. 11).


Fig. 11.

6. Slightly close the cartridge wings to hold the lens in place and then advance the plunger such that the haptics are compressed. The compression is correct when the haptic is pointing toward, but not touching, the optic (Fig. 12).


Fig. 12.

7. Close the cartridge wings together until the “click-lock” mechanism engages.
8. Push the lens into the conical tip by advancing the plunger until the spring just contacts the outer edge of the cartridge wings (Fig. 13). Pull the plunger back slightly to visually confirm that the lens remains in the conical tip and then push forward again. The lens is now loaded and ready for injection


Fig. 13.

9. With the conical tip bevel facing down, inject the lens by applying continuous pressure on the plunger until the lens is fully expressed from the tip. Clockwise rotation of the injector will compensate for any lens rotation.
10. Dispose of the injector after use in accordance with applicable regulations regarding the disposal of biohazardous materials.

PRECAUTIONS AND WARNINGS

1. The Bausch + Lomb IOL Injector is intended for single use only. Reuse and/ or resterilization may compromise device performance, which could cause serious harm to the patient's health and safety.
 2. Store the injector between +5°C and +30°C (+41°F and +86°F).
 3. Folding and compression of the lens should be done just prior to insertion and delivery.
 4. Failure to follow the instructions for use may result in patient injury.
- * Lenses and viscoelastic shown in images have been colored for visualization. Actual products will be clear.

SURGICAL PROCEDURE

The appropriate surgical techniques are the responsibility of the individual surgeon. Surgeons must determine the suitability of any particular procedure based on their medical training and experience.

For any questions regarding the use of the Bausch + Lomb IOL Injector, please contact Bausch + Lomb Customer Service at 1-800-338-2020 or the Authorized Representative of Bausch + Lomb in Russia:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia
 Tel.: +7 495 510 2879, fax: +7 495 510 2879
Office.RU@valeant.com

OPERATING CONDITIONS

The device must only be used at adequate facilities by a healthcare professional after training on this device.

Ambient Operating conditions

Allowable temperature (calculated), °C	Humidity (Relative), %	The velocity of air motion, meters per second
+21 – +24 (+21)	30-60	not more than 0,2

STORAGE AND TRANSPORTATION CONDITIONS

Store and transport sterile, sealed package with the Injector for introduction of intraocular lenses, INJ100 at a temperature +5 – +30 °C (+41 – +86 °F).

Humidity (Relative): 30-60 %

Contents of a package are sterile unless package is opened or damaged. Do not resterilize. Do not reuse. The device is for Single use only.

WARRANTY AND LIMITATIONS OF LIABILITY

Bausch & Lomb Incorporated warrants that the IOL Injector, INJ100, when delivered, will conform to the manufacturer's then current version of the published specifications for such IOL Injector, INJ100 in all material respects and shall be free from defects in material or workmanship for a period equal to the IOL Injector, INJ100 expiration date.

Bausch & Lomb Incorporated warrants that the IOL Injector, INJ100, when delivered, will conform to the manufacturer's then current version of the published specifications for such IOL Injector, INJ100 in all material respects and shall be free from defects in material or workmanship for a period equal to the IOL Injector, INJ100 expiration date.

Bausch & Lomb Incorporated excludes all other warranties, whether expressed, implied or by operation of law or otherwise, including but not limited to any implied warranties of merchantability or fitness for a particular purpose. Bausch & Lomb Incorporated shall not be liable for any incidental, consequential or exemplary loss, damage or expense, directly or indirectly arising from the use of the IOL Injector, INJ100 even if Bausch & Lomb Incorporated has been advised of the possibility of such loss, damage or expense.

Return Policy:

Bausch + Lomb will, within the return period as specified on your invoice, from the date of invoice, accept return of this product for a full refund less any handling and shipping charges incurred by Bausch + Lomb.

Customer must call their local Bausch + Lomb customer service representative to request a Return Good Authorization prior to expiration of the return period. It is the Customer's responsibility to properly pack all items being returned.

Please call +1-800-338-2020 for the authorized product return and detailed information obtaining.

The Authorized Representative of Bausch + Lomb in Russia:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

Tel.: +7 495 510 2879, fax: +7 495 510 2879

Office.RU@valeant.com

ADVERSE REACTION REPORTING

Adverse reactions and/or potentially sight-threatening complications that maybe reasonably regarded as IOL Injector, INJ100 related and that were not previously expected in nature, severity or degree of incidence should be reported to your customer service affiliate.

The Authorized Representative of Bausch + Lomb in Russia:

LLC VALEANT

Shabolovka St. 31, Bld. 5, Moscow, 115162, Russia

Tel.: +7 495 510 2879, fax: +7 495 510 2879

Office.RU@valeant.com

SEVICE INFORMATION

Bausch & Lomb Incorporated warrants that the IOL Injector, INJ100, when delivered, will conform to the manufacturer's current version of the published specifications for such IOL Injector, INJ100 in all material respects and shall be free from defects in material or quality for a period equal to the IOL Injector, INJ100 shelf life.

REQUIREMENTS FOR THE RECYCLING (THE PROTECTION OF THE ENVIRONMENT IN THE APPLICATION OF MEDICAL PRODUCTS)

IOL Injector, INJ100 after its using is classified as a class B of waste (according Russian Legislation) and should be disposed according to SanPin 2.1.2790-10 ("Sanitary rules and norms #2.1.7.2790-10").

PACKAGING AND LABELING INFORMATION

On the package and label of IOL Injector, INJ100 the following symbols are used:


Fig. 14. The symbols used on the package and label.

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Note: Optimum sealing conditions are highly dependent upon the materials being sealed, the equipment and production rates.

Marketinginformation

Marketing information

Produktfamilie:

Product family

- Eastar™ Copolyester 6763 ist die Handelsmarke von Eastman Chemical Company.
- Eastar™ Copolyester 6763 is the registered trademark of Eastman Chemical Company

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The information contained herein is based on best knowledge and experience of MEDIPACK AG. Users should consider these data only as a supplement to other information. A legally obligatory warranty of certain information as well as a suitability for concrete application cannot be derived from this data sheet.

Determination of the specific suitability of this product for individual applications is the sole responsibility of the purchaser.

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Part 1 (ISO 11607-1): Requirements for materials, sterile barrier systems and packaging systems

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1400 North Goodman Street
Rochester, NY 14609
585.338.6000
www.bausch.com

Certified True Copy Letter

I, Sharon A. Tonetta, hereby swear (or affirm) that the attached Material Safety Data Sheets are true and complete photocopy of the document.

Sharon A. Tonetta, PhD

Sharon A. Tonetta, Ph.D.
Vice President, Regulatory Affairs
Bausch & Lomb Incorporated

6 Dec 2017

Date

State of New Jersey

County of Union

Subscribed and sworn to (or affirmed) before me on this 6th day of December 2017, by Sharon A. Tonetta, Ph.D., proved to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Rita E. Moore

Signature of Notary Public

MATERIAL SAFETY DATA SHEET

According to Regulation (EC) No 1907/2006

EDITION: 20/11/2009

1. IDENTIFICATION OF THE SUBSTANCE / PREPARATION AND OF THE COMPANY / UNDERTAKING

Product name: UN 5068 HELLBLAU TR

Use of the preparation: Colouring agent for the plastic industry.

Company: COLOR-SERVICE GmbH & Co.KG
Offenbacher Landstr. 107-109
D-63512 Hainburg

Product safety phone No. : +49-(0)6182-954-241
Telefax: +49-(0)6182-954-111
e-Mail: k.wald@color-service.de

Emergency telephone No. : +49-(0)6182-954-241, only during the office hours.

2. HAZARDS IDENTIFICATION

No special hazards.
According to the Directive 1999/45/EC the preparation has not to be classified as dangerous.

3. COMPOSITION / INFORMATION ON INGREDIENTS

Chemical characterisation: Dispersion of pigments in an universal carrier (resin ester - polymer).
Hazardous ingredients: None which have to be labelled according to the EC-Directive 67/548/EC.

Ingredients with EC-workplace exposure limits:
None

4. FIRST AID MEASURES

Medical attention: Call a physician immediately
Inhalation: Move to fresh air. If you feel unwell, seek medical advice.
Skin contact: Take off all contaminated clothing immediately. Wash off with soap and plenty of water.
Eye contact: Rinse immediately with plenty of water and seek medical advice.
Ingestion: Rinse mouth. Consult a physician.

5. FIRE-FIGHTING-MEASURES

Suitable extinguishing media: Foam, water mist, CO₂, dry powder
Unsuitable extinguishing media:
Special exposure hazards: Do not release chemically contaminated water into drains, soil or surface water. Sufficient measures must be taken to retain the water used for extinguishing. Dispose of contaminated water and soil according to local regulations.
Main combustion products: Oxides of carbon, other toxic gases and/or vapours are also possible.
Special equipment for firefighters: Wear self-contained breathing apparatus.

SAFETY DATA SHEET according to Regulation (EC) No 1907/2006 EDITION: 20/11/2009

6. ACCIDENTAL RELEASE MEASURES

Personal precautions:	Do not ingest or breath in. Avoid contact with skin, eyes and clothing.
Environmental precautions:	Prevent contamination of soil, drains and surface waters. Do not allow release into sewerage system.
Measures in case of spillage:	Take up mechanically and collect in suitable container for disposal.

7. HANDLING AND STORAGE

Special fire and explosion precautions:	Avoid dust formation. Take precautionary measures against static discharges.
Storage conditions:	Keep tightly closed in a dry and cool place.
Technical protective measures:	Ensure good ventilation and local exhaust. Heat only in areas with appropriate exhaust ventilation.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Exposure limit values:	None
Occupational exposure controls:	Design of appropriate work processes and engineering controls, the use of adequate equipment and materials, the application of collective protection measures at source, where exposure cannot be prevented by other means the use of individual protection measures, such as personal protection equipment.

PERSONAL PROTECTIVE MEASURES

Respiratory protection:	Respirator must be worn if exposed to dust.
Hand protection:	Suitable gloves of nitrile rubber.
Eye protection:	Suitable goggles.
Skin and body protection:	Handle in accordance with good industrial hygiene and safety practice.
Workplace hygiene:	Keep away from food and drink. Follow good industrial hygiene praxis. Before breaks and the end of work wash hands and face.

Environmental exposure controls:	See paragraph 6. and 7., no further measures are necessary.
----------------------------------	---

9. PHYSICAL AND CHEMICAL PROPERTIES

Form:	Pellets
Colour:	Blue
Carrier:	Resin ester - polymer
Melting range:	65-85°C
Odour:	Weak
pH-Value:	n.a.
Boiling range:	Not measured.
Flashpoint:	> 200°C
Flammability (solid/gas):	Combustible
Oxidising properties:	No data available.
Vapour pressure:	n.a.
Density:	appr. 1,15 g/cm³ at 20 °C
Solubility in solvents:	Soluble in some organic solvents.
Water solubility:	Insoluble
Partition coefficient log Pow:	No data available.
Viscosity:	n.a.
Evaporation rate:	n.a.
Vapour density:	n.a.

10. STABILITY AND REACTIVITY

Thermal stability:	Decomposition temperature > 280°C
Thermal sensitivity:	No data available.
Shock sensitivity:	No data available.
Materials to avoid:	Oxidising agents.
Hazardous decomposition products:	Thermal decomposition or burning may release oxides of carbon. Other toxic gases and/or vapours are also possible.

SAFETY DATA SHEET

according to Regulation (EC) No 1907/2006 EDITION: 20/11/2009

11. TOXICOLOGICAL INFORMATION

Acute oral toxicity:	LD50: > 2000 mg/kg	Species: Rat
Acute dermal toxicity:	No data available..	
Acute inhalation toxicity:	No data available.	
Skin irritation acute:	No data available.	
Eye irritation acute:	No data available.	

12. ECOLOGICAL INFORMATION

ECOLOGY TOXICITY

Fish toxicity:	No data available.
Aquatic Invertebrates:	No data available.
Water plants toxicology:	No data available.
Bacteria toxicity:	No data available.
Algae toxicity:	No data available.
Mobility:	No data available.

PERSISTENCY and BIODEGRADABILITY

Persistence:	No data available.
Biodegradability in aqueous systems:	Basic pigments and polymer are barely water soluble. The preparation has not been tested.
Bioaccumulative-potential:	No data available.
Results of PTB-assessment:	No data available.
Other adverse effects:	No data available.
Remarks:	The preparation is not hazardous for the environment., release into the environment has to be avoided.

13. DISPOSAL CONSIDERATIONS

Waste from residues/unused products:	Burning or disposal in accordance with local regulations for chemical waste.
Contaminated packaging:	Uncontaminated emptied packagings may be treated like household waste.
Reutilisation:	Where possible recycling is preferred to disposal or incineration.

14. TRANSPORT INFORMATION

Transport railway and road:	RID/ADR-Class: free	
Transport air:	UN-No.: free	
	IATA-Class: free	
	UN-No.: free	
Transport sea:	IMDG-Class: free	Marine Pollutant: No
	UN-Nr.: free	

15. REGULATORY INFORMATION

Declaration according to EC Directive:	Classification not required.
Evaluation of the product safety:	Product safety evaluations of chemicals in this preparation have not been done.

16. OTHER INFORMATION

Use restrictions:	Not recommended for the direct coloration of food. Not recommended for the direct coloration of articles that come into contact with internal organs or body fluids.
Supporting regulations which apply:	EC-Directive 1999/45/EC with the amendment 2006/8/EC, EC-Directive 67/548/EC with the amendment 2006/121/EC, REACH-Regulation (EC) No 1907/2006.
Recommended limitation of use:	Only for industrial applications.
R-phrases referring to paragraph 2. and 3. :	R-phrases are not regulated.

SAFETY DATA SHEET

according to Regulation (EC) No 1907/2006

EDITION: 20/11/2009

Changes to last edition:

Adjustment according to the REACH – REGULATION (EC) No 1907/2006.

Data sheet issued by:

Dr. Klemens Wald, Tel. +49(0)6182-954-241

The information and recommendation in this publication are to the best of our knowledge, information and belief accurate at the date of publication. Nothing herein is to be construed as a warranty, express or otherwise. In all cases, it is the responsibility of the user to determine the applicability of such information or the suitability of any products for their own particular purpose.

MATERIAL SAFETY DATA SHEET

File NO.: 65280-001

DENKA TRANSPARENT ABS Natural Color

MSDS DATE : 2005/8/8

SECTION 1: PRODUCT AND COMPANY IDENTIFICATION

RPRODUCT NAME : DENKA TRANSPARENT ABS (NATURAL COLOR)
GRADE NAME : CL, CL-301, CU-301, CL-301-Q, CL-20, CL-430

MANUFACTURER : DENKI KAGAKU KOGYO KABUSIKI KAISYA
DIVISION : FUNCTIONAL PLASTIC DIVISION
ADDRESS : NIHONBASHI MITSUI TOWER, 1-1 NIHONBASHI-MUROMACHI
2-CHOME, CHUO-KU, TOKYO 103-8338 JAPAN

TELEPHONE : JPN(03)5290-5414
FAX PHONE : JPN(03)5290-5503

EMERGENCY CONTACT: JPN0436-26-3232 Denki Kagaku Kogyo Chiba Plant

PRINTED DATE : 2005/8/8

SECTION 2: COMPOSITION/INFORMATION ON INGREDIENTS

COMPOSITION : Methylmethacrylate-Acrylonitrile-Butadiene-Styrene-Butylacrylate Copolymer

CHEMICAL FORMULA :

Copolymer	Additive
More than 98wt%	Less than 2wt%
$\left[\begin{array}{c} \text{CH}_2 - \text{C}(\text{CH}_3) \\ \\ \text{COOCH}_3 \end{array} \right]_L \left[\begin{array}{c} \text{CH}_2 - \text{CH} \\ \\ \text{CN} \end{array} \right]_M \left[\text{CH}_2 - \text{CH} = \text{CH} - \text{CH}_2 \right]_N \left[\begin{array}{c} \text{CH}_2 - \text{CH} \\ \\ \text{C}_6\text{H}_5 \end{array} \right]_O \left[\begin{array}{c} \text{CH}_2 - \text{CH} \\ \\ \text{COOC}_4\text{H}_9 \end{array} \right]_P$	

CAS.No. : 25101-28-4/107592-06-3

SECTION3: HAZARDS IDENTIFICATION

ROUTES OF ENTRY : Inhalation, ingestion, eye and skin contact
POTENTIAL HEALTH EFFECTS : No information
CARCINOGENIC EFFECTS : No information
MUTAGENIC EFFECTS : No information
TERATOGENIC EFFECTS : No information
EFFECTS ON THE REPRODUCTIVE SYSTEM: No information

SECTION4: FIRST AID MEASURES

EYES : Remove contact lenses stick to eye balls. Immediately flush the affected eyes well with clean water. Do not rub the eyes to prevent corneas injuries. Seek medical attention if irritation persists.

MATERIAL SAFETY DATA SHEET

File NO.: 65280-001

DENKA TRANSPARENT ABS Natural Color

MSDS DATE : 2005/8/8

- SKIN : Wash well with water. If condenses of gas generated from the substance molten at high temperature, wash well with soapy water. In case of molten plastic contact with skin, cool immediately with water and seek medical attention.
- INGESTION : Not likely to ingest. If a large amount is swallowed, seek medical attention.
- INHALATION : Not likely to inhale pellets due to their physical form. If inhale gases and/or fumes from the molten resin at high temperature, leave the place immediately and breathe fresh air. If coughing and/or difficult breath occurs, seek immediate medical attention.

SECTION5: FIRE FIGHTING MEASURES

- FIRE FIGHTING : This product generates intense heat, dense black smoke, carbon dioxide, carbon monoxide, carbon, hydrogen cyanide and so on. Wear fire protective closing and breathing apparatus for fire fighting.
- EXTINGUISHING MEDIA: Water pouring, water spray and regular fire extinguishers can be used.
- AUTOIGNITABILITY : 465°C(ASTM-D 1929-77)
- EXPLOSION LIMIT : Not available
- Flash Point : Not available
- SPECIAL FIRE FIGHTING PROCIDURES: None
- HAZARDOUS DECOMPOSITION PRODUCTS: The resin generates decomposed gases.

SECTION6: ACCIDENTAL RELEASE MEASURES

If spilled on road or floor, sweep and collect materials to prevent slipping and falling. Released materials must be collected for environmental protection.

SECTION7: HANDLING AND STORAGE

- HANDLING : Gases generated during processing the material may cause irritation to the skin and the respiratory organs. If a large amount of the gases is inhaled, it may cause nausea and headache. Refrain from inhaling them. Prevent dust accumulation from cutting and grinding the materials, which may cause a dust explosion by electrostatic charges or electric spark.
- STORAGE : Store under no direct sunlight and in a dry place avoiding excessive heat and sources of ignition. Prevent electrostatic charges which may cause explosion.

SECTION8: EXPOSURE CONTROLES/PERSONAL PROTECTION

- ENGINEERING MEASURES : Local ventilation system is required for better working environment as gases are generated during the high temperature processing of the materials.
- RESPIRATORY PROTECTION : Wear respirator during machinery processing or sanding of resin which produces dust. When working under high dense gases and/or fumes circumstances, wear organic gas respirator.

MATERIAL SAFETY DATA SHEET

File NO.: 65280-001

DENKA TRANSPARENT ABS Natural Color

MSDS DATE : 2005/8/8

EYE PROTECTION : Wear protective glasses during machinery processing or sanding.

SKIN PROTECTION : Wear long sleeve clothing when handling molten resin.

OTHER PROTECTIVE EQUIPMENT: Wear heat-resistant gloves when handling molten resin.

SECTION9: PHYSICAL AND CHEMICAL PROPERTIES

APPEARANCE : Pellets

ODOR : Odorless

PHYSICAL STATE : Solid

pH : Not available

BOILING POINT : Not available

MELTING POINT : No definite temperature. Softening in the broad range of temperature between 130-150°C.

SPECIFIC GRAVITY : 1.09-1.11

SOLUBILITY IN WATER: Water - Insoluble

SECTION10: STABILITY AND REACTIVITY

FLAMMABILITY : No information

SOLUBILITY : Water - Insoluble

OXIDIZABILITY : Not occur under proper storage condition

SELF-REACTIVITY : There is no self-reactivity at room temperature.

EXPRESSIVENESS : At high temperature(250-400°C) decomposed gases are generated from resin decomposition. In this case cool the molten resin with water as soon as possible.

DUST EXPLOSIVES : Explosive

STABILITY AND REACTIVITY : None under ordinary condition

SECTION11: TOXICOLOGICAL INFORMATION

SKIN CORROSIVENESS: None

IRRITABILITY : Physical irritation is observed(SKIN/EYE)

ACUTE TOXICITY : Oral LD₅₀(rat) > 5g/kg(Assumed value)

SECTION12: ECOLOGICAL INFORMATION

No data available. Release may cause environmental hazard.

SECTION13: DISPOSAL CONSIDERATIONS

Comply with national and/or local regulations.

MATERIAL SAFETY DATA SHEET

File NO.: 65280-001

DENKA TRANSPARENT ABS Natural Color

MSDS DATE : 2005/8/8

SECTION14: TRANSPORT INFORMATION

Avoid moisture and rough handling which may damage its package. If bag is broken and the materials are spilled, sweep and collect them immediately to prevent slipping and falling.

SECTION15: REGULATORY INFORMATION

Comply with national and local regulations.

SECTION16: OTHER INFORMATION

This information is based on currently available technical materials and data, and is subject to change due to new experience and knowledge. All guideline is for ordinary use and not for special application or handling, which requires proper safety measures. This publication is made to provide information only and does not guarantee anything.

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

1400 North Goodman Street
Rochester, NY 14609
585.338.6000
www.bausch.com

Certified True Copy Letter

We, Bausch & Lomb incorporated, hereby swear (or affirm) that the attached Supplement No. 1 to the Instruction for Use on the Injector for Intraocular Lens insertion, INJIOO are true and complete photocopies of the documents on file at Bausch & Lomb Incorporated, 1400 North Goodman Street, Rochester, NY 14609, USA.


Isabelle B. Lefebvre, MSc.RA
Vice President, Regulatory Affairs

08 February 2020
Date

State of New Jersey
County of Union

Subscribed and sworn to (or affirmed) before me on this 8th day of February, 2020,
by Isabelle B. Lefebvre proved to me on the basis of satisfactory evidence to be the
person(s) who appeared before me.

Signature of Notary Public

Supplement No. 1 to the Instruction for Use of Medical Device IOL Injector, INJ100.

Change of the authorized representative of the manufacturer.

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