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0064245

SECRETARY OF STATE
STATE OF INDIANA

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Connie Lawson
Secretary of State


APOSTILLE

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1. Country: United States Of America
2. This public document has been signed by *Janice E. Mercado-Aguila*
3. Acting in the capacity of notary public in & for *Monroe County*
4. And bears the seal/stamp of notary public in & for the State of Indiana

CERTIFIED

5. at Indianapolis, Indiana
6. This *Eighth* day of *May*, 2019
7. by the Secretary of State of Indiana
8. No. A2019-0508121437
9. Seal/Stamp:
10. Signature:

Connie Lawson

Connie Lawson
Indiana Secretary of State

*Effective May 1, 2003 all apostilles from the Indiana Secretary of State will have an electronically printed seal.
This document was revised March 17, 2012.*

COOK INCORPORATED
750 DANIELS WAY
BLOOMINGTON, IN 47404 U.S.A.
PHONE 812.339.2235 TOLL FREE 800.457.4500
WWW.COOKMEDICAL.COM


CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instructions for Use for the "Peritoneal Lavage Set" Submission.

Kris Weathers
Manager, Regulatory Affairs
Cook Incorporated

24 May 2019
Date

United States of America)
State of Indiana) ss:
County of Monroe)

Subscribed and sworn to before me this 2 day of May month of 2019 year.

My commission expires on December 21, 2019

Janice E. Mercado-Aguila, Notary Public
Commission #631710

INSTRUCTIONS FOR USE

Peritoneal Lavage Set

INTENDED USE

Peritoneal Lavage Set is used to infuse and withdraw rinsing solutions in the peritoneal cavity. It is performed when intra-abdominal bleeding secondary to trauma is suspected.

DEVICE DESCRIPTION

Table of conversions

Parameter	Conversion
1/3 Fr	1 mm
1 Inch	25.4 mm

PERFORMANCE CHARACTERISTICS

Strength characteristics of the catheter:

Bond strength of components - not less than 15 N;

Tensile strength - not less than 8.1 MPa;

Relative elongation at break point - not less than 150%;


Strength characteristics of the wire guide:

Bond strength of wire guide core with sheath and sheath with safety wire - not less than 10 N;

Strength characteristics of needle:

Bond strength of hub and needle - not less than 69 N;

Peritoneal Lavage Set C-PLS-901

The Peritoneal Lavage Set is comprised of a peritoneal catheter, percutaneous entry access needle and wire guide. The peritoneal lavage catheter is manufactured as a 9.0 French/3.0 mm (± 0.003 inch/ ± 0.076 mm). This catheter is manufactured in a length of 20 centimeters (± 5 mm) with an angled distal curve configuration characterized by 90 sideports as a .035 inch/0.889 mm diameter and 0.121 inch/3.07 mm between each sideport, all located within 8 centimeters of the distal tip. The sideports are placed in four rows along the distal tip, with the first and second rows containing 22 sideports and the third and fourth row containing 23 sideports. Each row of sideports is offset against the previous row. The peritoneal lavage catheter is assembled with a cap and adapter on the proximal end.

A wire guide is a component of the Peritoneal Lavage Set. This wire guide is manufactured as a 0.038 inch/0.97 mm ($+0.0001$ inch/ -0.0003 inch ($+0.003$ mm/ -0.008 mm) diameter, with a length measurement of 40 centimeters ($+ 5$ cm/ $- 3$ cm). This wire guide is a straight wire, without a distal curve radius.

A percutaneous entry access needle is a component of the Peritoneal Lavage Set. This access needle is manufactured as an 18 gage/1.25 mm (± 0.05 mm) needle measuring 7 centimeters (± 3 mm) in length and designed with a lancet bevel, bevel angle is 60 degrees. The needle has a hub on the proximal end and is compatible with a .038 inch/0.97 mm diameter wire guide.


Component	Catheter		Wire Guide		Entry Access Needle	
	Diameter	Length	Diameter	Length	Diameter	Length
C-PLS-901	9.0 French/3.0 mm (± 0.003 inch/ ± 0.076 mm)	20 centimeters (± 5 mm)	0.038 inch/0.97 mm ($+0.0001$ inch/ -0.0003 inch ($+0.003$ mm/ -0.008 mm)	40 centimeters ($+ 5$ cm/ $- 3$ cm)	18 gage /1.25 mm (± 0.002 inch/ 0.05 mm)	7 centimeters (± 3 mm)

Illustration of the Peritoneal Lavage Set

Parameters of Catheter


Lengths: A – 20 cm (± 5 mm); B – 6 mm-8mm.

Parameters of Catheter Adapter:

Lengths: A – .735 inch /18.67 mm ($\pm .010$ inch / .25 mm); B – .281 inch /7.14 mm ($\pm .005$ inch / .13 mm); C – .312 дюйма/7.92 mm ($\pm .005$ inch / .13 mm).


Diameters: D – .344 inch /8.77 mm ($\pm .005$ inch / .13 mm); E – .210 inch /5.33 mm ($\pm .005$ inch / .13 mm); F – .118 inch /3.0 mm ($\pm .002$ inch / .05 mm).

Parameters of Catheter Cap:

Lengths: A – .41 inch /10.41 mm ($\pm .01$ inch / .25 mm); B – .200 inch /5.08 mm ($\pm .005$ inch / .13 mm).


Diameters: C – .350 inch /8.89 mm ($\pm .005$ inch / .13 mm); D – .120 inch /3.05 mm ($\pm .002$ inch / .05 mm); E – .168 inch /4.27 mm ($\pm .005$ inch / .13 mm); F – .330 inch /8.38 mm ($\pm .005$ inch / .13 mm).

Parameters of Needle Hub:

Lengths: A – .945 inch /24.0 mm ($\pm .005$ inch / .13 mm); B – .055 inch /1.40 mm ($\pm .005$ inch / .13 mm); C – .030 inch / .76 mm ($\pm .005$ inch / .13 mm).

Diameters: D – .207 inch /5.26 mm - .235 inch /5.97 mm.

INTENDED USE:

Peritoneal Lavage Set is used to infuse and withdraw rinsing solutions in the peritoneal cavity. It is performed when intra-abdominal bleeding secondary to trauma is suspected.

CONTRAINDICATIONS:

- Morbid obesity
- Previous abdominal surgery

- Severe cirrhosis
- Blood clotting disorders
- Pelvic Fractures
- Pregnancy past first trimester

POTENTIAL ADVERSE EVENTS / COMPLICATIONS:

- Infection
- Injury to area of catheter insertion
- Bladder puncture
- Inadequate return of infused fluid
- Infusion into abdominal wall

INSTRUCTIONS FOR USE

1. Under local anesthesia (sold separately), using the Seldinger technique, insert the access needle into the peritoneal cavity.
2. Pass a wire guide through the needle and direct toward the right or left pelvic gutter.
3. Withdraw the needle over the wire guide, leaving the wire in place.
4. Pass the catheter over the wire into the peritoneal cavity, gently twisting or rotating the catheter if needed.
5. Remove the wire guide.
6. Infuse saline (sold separately) into the peritoneal cavity.
7. Withdraw the saline through the catheter.
8. If blood is observed in the aspirated fluid, there is bleeding and patient will need a surgical procedure.
9. If no blood is observed, the lavage is completed and catheter may be removed.

Note: this is a Short-term medical product (continuous use for no more than 30 days).

MATERIALS

Device	Material	Patient Contacting
Catheter	Tubing : Polytetrafluoroethylene Radiopaque № 10050 Black colorant: Ferrite Black № 68186 Bismuth Trioxide № 1304 Catheter Adapter: Polyamide 6 № 40340 Markings: Ink № 70720- Black Catheter Adapter Cap: Polyamide 6 № 40340	Direct Contact

	Adhesive: Loctite 4011 Ethyl 2-cyanoacrylate 60-100% Thickener 5-10%	
Wire Guide	Wire Guide stem: Stainless Steel AISI 304 C-Max 0,08 Cr-18 - 20 Fe-66,345 - 74 Mn-Max 2 Ni-8 - 10,5 P-Max 0,045 S-Max 0,03 Cu-Max 1 Wire Guide coil: Stainless Steel AISI 304 C-Max 0,08 Cr-18 - 20 Fe-66,345 - 74 Mn-Max 2 Ni-8 - 10,5 P-Max 0,045 S-Max 0,03 Cu-Max 1 Solder: Silver № 1618 Silver 30-60% Copper 10-30% Zinc 10-30% Potassium Bifluoride 10-30% Tin 3-7% Potassium Fluoroborate 3-7% Potassium Pentaborate Tetrahydrate 3-7% Potassium Tetraborate Tetrahydrate 1-5% Potassium Silicofluoride 0.1-1%	Direct Contact
Entry Access Needle	Cannula: Stainless Steel AISI 304 C-Max 0,08 Cr-18 - 20 Fe-66,345 - 74 Mn-Max 2 Ni-8 - 10,5 P-Max 0,045 S-Max 0,03 Cu-Max 1 Hub: Polycarbonate № 80020 Pink Colorant: CI Solvent Red 207 Adhesive: Dymax 1193-M-SV05	Direct Contact

PACKAGING

Component	Material name	Patient contacting
Packaging	Backing –TYVEK: 1073B/109P. Film part– Polyester Oliver-Tolas (Mylar).	Non contact.

PACKAGING INFORMATION

Surface density from 50 to 80 (± 2)g/m²

Tensile strength - Tensile strength not less than 0.4 kN/m

Strength of the welded joint - Strength of the weld at least 0.12 MPa

Width of welded joint - Width of adhesive layer not less than 8 mm

This device contains only primary packaging.

PERFORMANCE CHARACTERISTICS

Products mentioned in this Technical data summary are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.15 mm.

Products mentioned in this Technical data summary, in packaging are stable under vibration loads in the range 10-55 Hz, amplitude displacement of 0.35 mm and shock loads at peak acceleration of 100 m / s² (10g), and duration of shock acceleration of 16 ms.

Products mentioned in this Technical data summary are stable at temperatures in the range of +32 to + 42 ° C.

Products mentioned in this Technical data summary, in packaging are stable at a transportation and storage temperature in the in the range of -50 to + 50 ° C.

Products mentioned in this Technical data summary are stable under sudden changes in temperature.

Products mentioned in this Technical data summary are supplied in sterile packages, do not attempt to clean, disinfect or re-sterilize.

WARRANTY

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SHAG".

Russia, 119002, Moscow, Karmanitskiy per., 9.office 501a.

phone: (495) 956-13-09;

Fax: (495) 956-1310,

e-mail: info@schag.ru

Manufacturer's warranty binding on quality products.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Cook Incorporated manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC.

Cook Incorporated manufactures all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820. *Quality System Regulations*.

All of Cook facilities are also certified to ISO 13485 (2003/2012) *Quality Management System*.

SHELF LIFE

The "use by" date as denoted by symbols on the labeling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

Peritoneal Lavage Sets have a shelf-life of 3 years.

STERILIZATION DETAILS

Peritoneal Lavage Set is supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10^{-6} . These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile.

STORAGE RULES

Peritoneal Lavage Set should be stored in a dark, dry, cool place; avoid extended exposure to light.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation.

USAGE CONDITIONS

Peritoneal Lavage Set is used in the patient care medical institutions.

POTENTIAL CONSUMERS

This product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

REPAIR/MAINTENANCE

Peritoneal Lavage Set is intended for single use only. Do not attempt to repair or carry out maintenance.

SUPPLY DETAILS

Peritoneal Lavage Set is supplied in individual packaging together with instruction for use. Individual packages are put in the corrugated board boxes. Boxes are being taped with the control tape to cover the sides or along the perimeter.

SUPPLY DETAILS

Peritoneal Lavage Set packaged with: Catheter - 1p, Wire Guide – 1p, Entry Access Needle - 1p., instruction for use – 1p.

PACKAGING PARAMETERES

RPN	Heigth, cm	Width, cm	Weight, g
Tolerance	±10%	±10%	±10%
C-PLS-901	40.4	15.2	25.1

TRANSPORT LABELING

Individual packages are packed in corrugated cardboard boxes. Boxes are sealed with a control tape with a call on the side walls or along the perimeter. Transportation marking contains information about the dimensions of the box, a unique barcode, the logo of the company "Cook Medical", as well as the symbols:

- a symbol for "Fragile, handle with care"

- the symbol for "Temperature Limit"



- the symbol for "Keep Dry"

LABELING


Peritoneal Lavage Set is supplied in the individual packaging, which should show no signs of damage or opening. Upon detecting the bends, breaks and other defects the device cannot be used.

The label of the device shows the following device parameters:

Symbol	Parameters
	Manufacturer
	Batch code
	Catalogue number
	Date of manufacture
	Use-by date
	Do not re-use
	Sterilized using ethylene oxide
	Keep dry
	Keep away from sunlight
	Authorized representative in the European Community

Additional symbols on the label:

- marking "For specialists only";
- medical device name;
- parameters and general device drawing;
- symbol about CE certification;
- bar-code;


Product	Label
Peritoneal Lavage Set C-PLS-901	 The label for the Peritoneal Lavage Set C-PLS-901 includes the following information: - Product Name: Peritoneal Lavage Set - Reference: REF C-PLS-901 - Lot Number: REF G05267 - Size: 9.0Fr / 20cm - Diagram of the lavage set - STERILE EO - LOT 8324743 - Expiry Date: 2020-10-25 - Manufacture Date: 2017-10-25 - Icons for disposal and handling - COOK MEDICAL logo - Barcode with REF C-PLS-901 and LOT 8324743 - COOK MEDICAL logo - EC REP logo - CE 0088


DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with Council Directive 93/42/EEC

For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms N2.1.7.2790-10 of December 12, 2010 "Sanitary-epidemiological requirements for the treatment of medical waste". Risk class of medical waste: class B.

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