Slide-On™ Sheath for Sensory Testing for Flexible ENT Scopes & ENT-1000 Flexible ENT Scopes

Gain Slide-On™ pour l'étude des seuils de sensibilité pour Endoscopes ORL flexibles & endoscopes ORL flexibles ENT-1000 / Guaina Slide-On™ per test sensoriali per Endoscopi flessibili ENT ed ENT-1000 per otorinolaringoiatria / Slide-On™ Schaft zum sensorischen Testen für Flexible ENT Endoskope und flexible ENT Endoskope vom Typ ENT-1000 / Vaina Slide-On™ de evaluación sensorial para Endoscopios flexibles ENT y ENT-1000 / Slide-On™ Schede voor zintuiglijk testen voor flexibele KNO-scopen en KNO-1000 flexibele KNO-scopen / Slide-On™-sheath til sensoriske test til Flexible ENT-skoper & ENT-1000 Flexible ENT-skoper / Slide-On™-suojus sensorista testaamista varten käytettäväksi taipuisien ENT- ja ENT-1000-endoskooppien kanssa / Slide-On™-holje för sensorisk testning för flexibla ENT-endoskop & flexibla ENT-1000-endoskop / Bainha para Testes Sensoriais Slide-On™ para Nasofaringoscópios Flexíveis ENT e Nasofaringoscópios Flexíveis ENT-1000 / Οηκάρι Slide-On™ για Αισθητήριο Ελεγχό για Εύκαμπτα Ενδοσκόπια ENT & Ευκαμπτα Ενδοσκόπια ENT-1000 / Wsuwana osłonka Slide-On™ do badania sensorycznego dla elastycznych wzorników ENT i ENT-1000 / Pouzdro Slide-On™ pro senzorične testovani s použitím ohebných endoskopů ENT a ENT-1000 / Slide-On™ hüvely szenzorális teszteshez Flexible ENT Scopes és ENT-1000 Flexible ENT Scopes eszközökkel való használatához / Duyusal Test için Slide-On™ Kılıf Esnek KBB Skopları ve ENT-1000 Esnek KBB Skopları / Slide-On™ Hylse for sensorisk test for Fleksible ENT-omfang og ENT-1000 fleksible ENT-omfang / včexon Slide-On™ для сенсорного исследования для гибких ЛОР-эндоскопов и гибких ЛОР-эндоскопов ENT-1000 / Slide-On™ obloga za senzorno ispitivanje za savitljive ENT endoskope i ENT-1000 savitljive ENT endoskope / Ovojnica Slide-On™ za senzorno testiranje za fleksibilne ORL endoskope i fleksibilne ORL endoskope ENT-1000 Flexible

INSTRUCTIONS FOR USE / MODE D'EMPLOI / ISTRUZIONI PER L'USO / BEDIENUNGSANLEITUNG / INSTRUCCIONES DE USO / GEBRUIKSAANWIJZINGEN / BRUGSANVISNING / KÄYTTÖOHJEET / BRUKSANVISNING / INSTRUÇÕES DE UTILIZAÇÃO / ΟΔΗΓΙΕΣ ΧΡΗΣΗΣ / INSTRUKCJA OBSŁUGI / NAVOD K POUŽITÍ / HASZNÁLATI UTASÍTÁS / KULLANIM TALIMATLARI / BRUKSANVISNING / ИНСТРУКЦИЯ ПО ПРИМЕНЕНИЮ / UPUTSTVO ZA UPOTREBU / UPUTE ZA UPORABU

Rx Only

Figure 1.A / Figure 1.A / Figura 1.A / Abbildung 1.A / Figura 1.A / Afbeelding 1.A / Figur 1.A / Kuva 1.A / Figur 1.A / Figura 1.A / Σχήμα 1.A / Rysunek 1.A / Obrazek 1.A / 1.A .abra / Şekil 1.A / Figur 1.A / рис.1.A / Slika 1.A / Slika 1.A


Figure 1.B / Figure 1.B / Figura 1.B / Abbildung 1.B / Figura 1.B / Afbeelding 1.B / Figur 1.B / Kuva 1.B / Figur 1.B / Figura 1.B / Σχήμα 1.B / Rysunek 1.B / Obrazek 1.B / 1.B .abra / Şekil 1.B / Figur 1.B / рис.1.B / Slika 1.B / Slika 1.B



Figure 2.A / Figure 2.A / Figura 2.A / Abbildung 2.A / Figura 2.A / Afbeelding 2.A / Figur 2.B / Kuva 2.A / Figur 2.A / Figura 2.A / Σχήμα 2.A / Rysunek 2.A / Obrazek 2.A / 2.A .abra / Şekil 2.A / Figur 2.A / рис. 2.A / Slika 2.A / Slika 2.A

Figure 2.B / Figure 2.B / Figura 2.B / Abbildung 2.B / Figura 2.B / Afbeelding 2.B / Figur 2.B / Kuva 2.B / Figur 2.B / Figura 2.B / Σχήμα 2.B / Rysunek 2.B / Obrazek 2.B / 2.B .abra / Şekil 2.B / Figur 2.B / рис. 2.B / Slika 2.B / Slika 2.B

EN**DEVICE DESCRIPTION**

Warning: Please read this entire manual carefully before using the Endosheath™ System.

The Slide-On™ Sheath for Sensory Testing is a latex-free, sterile, disposable protective covering with an air delivery channel for use with flexible Nasopharyngoscopes. The Sheath slides on and off the nasopharyngoscope for easy installation and removal. The optically clear window at the tip of the Sheath is designed to fit snugly against the endoscope's lens maintaining complete optical integrity. An air channel is provided for the delivery of air via the AP-4000 Air Pulse Sensory Stimulator System. The Sensory Sheath, for use with the AP-4000, is intended to enable a standard flexible Nasopharyngoscope to provide Sensory Testing and a fast and effective method of reprocessing your nasopharyngoscope. The sheath also ensures that the Insertion tube is covered with a sterile, disposable sheath for each procedure.

The user of this equipment should be thoroughly trained in the techniques of nasopharyngoscopy and with sensory testing procedures. The equipment is specifically designed for the applicable procedures and should not be used for any purposes other than those for which it has been designed.

The contents of this manual describe the recommended procedures for preparing and using the Medtronic Slide-On™ Sensory Sheath for use with the AP-4000 Sensory Stimulator system. It also contains pertinent information on the proper care and handling of the endoscope and Sheath during operation.

CAUTION: Sensory Sheath type (as described on product label) must match your endoscope make and model for the system to function properly. Damage to your endoscope or the Sheath may occur if an attempt is made to install the incorrect Sheath model.

NASOPHARYNGOSCOPE		SENSORY SHEATH CATAOG NUMBER	SENSORY TESTING AIR SOURCE
Vision Sciences™ Endoscopes	ENT-1000	445101	AP-4000
	ENT-2000/3000	445401	
	E-F100	445301	
	ENT-4000	445601	
	ENT-5000	445601	
Olympus Models	ENF-P4/GP	444401	
	ENF-P3/P2	444301	
Pentax Models	FNL-13S	442401	
	VNL-1330	442401	
	FNL-10S/10BS/10RBS/10P2/ 10RP3	445401	
Machida Models	ENT-3L	445301	
	ENT-4L	442401	
Karl Storz Model	11101 RP	446301	
	11101 RP1	445401	
Welch Allyn Model	RL-150	445401	
Smith & Nephew/ Gyrus Model	Duraview OL-1	445301	
Medtronic	Sharpsite	445401	
	ACMI	445401	
Gyrus	800-600 3.6mm	445401	
Optim Medical	ENTity	444401	
Xion Model	EV-NE	442401	
	EV-NC	444401	
Richard Wolf Model	7224.001/7224.001S02	445401	

INDICATIONS FOR USE

The EndoSheath™ System provides a sterile, disposable protective covering for the scope to be used during flexible endoscopic examination of the upper airway, vocal cords and/or nasal passages. The System may also be used in conjunction with the AP-4000 Air Pulse Sensory Stimulator to elicit Laryngeal Closure Reflex (Swallow) and to measure the sensory discrimination threshold at which the reflex occurs in the area of the Upper Airway innervated by the Superior Laryngeal Nerve

Design Features

The disposable sheath features a specially designed air channel for Sensory Testing.

Equipment

EN

Equipment required to use the Sensory Sheath includes:

- Flexible Nasopharyngo-Laryngoscope
- Slide-On™ Sheath for Flexible Endoscopic Evaluation of Swallowing with Sensory Testing
- AP-4000 Air Pulse Sensory Stimulator System Set (Vision Sciences™ model or Pentax model)

WARNINGS

- After the sheath is installed, do not touch the scope control body or eyepiece with the hand you use to manipulate the insertion tube. If you must rotate the diopter adjustment ring after the sheath is installed, do so using a gauze pad or replace your glove prior to performing the procedure. Failure to follow instructions may cause patient injury.
- Never attempt to advance any accessory instrument or object into the luer lock connector or air channel.
- Check product upon receipt to ensure damage did not occur during shipment.

- Carefully inspect the package and Sheath before use for evidence of damage during transit or storage. If the Sheath package has been damaged or otherwise compromised, please return to Medtronic.
- Do not use the Sheath on an endoscope unless the endoscope model number is designated on the Sheath label. It is extremely important that the proper sheath be used with an endoscope it was designed to fit.
- Do not use damaged for out-of-specification endoscopes with Sheaths. Please refer to the Insertion Tube Template Card (Lit # 890583) to identify damage and specification for your scope.
- The Sheath is for SINGLE USE ONLY and should be properly disposed of following use. Do not attempt to reuse the Sheath.
- The Sheath is provided STERILE and should not be resterilized. Resterilization may damage the integrity of the product.
- DO NOT place liquids inside the Sheath prior to or during use.
- DO NOT expose Sheath to Alcohol, Alcohol-based liquids, or any cleaning agents prior to or during use. The ENT Sheath is NOT compatible with alcohol-based lubricants.
- Use with the AP-4000 Air Pulse Sensory Stimulator air source ONLY.
- Refer to the AP-4000 Air Pulse Sensory Stimulator Owners Manual for proper set-up and operation.
- Appropriate protective gear should be worn during nasopharyngoscopy procedures and when handling contaminated sheaths, endoscopes and accessory equipment.
- Install and remove the Sheath according to the instructions in this manual. Failure to follow instructions may lead to Sheath or Scope damage.

SLIDE-ON™ SENSORY SHEATH DESCRIPTION

1. Proximal Connector: The Proximal Connector interfaces with the endoscope at the controlbody-strain relief junction to secure the sheath to the scope.
2. Luer Lock Connector: The standard luer lock connector provides a means of securing the AP-4000 air supply tube to the sheath for sensory testing.
3. Inner Sheath Material: The inner sheath material is a tough thermoplastic elastomeric covering the insertion tube and provides a sterile, disposable protective barrier for the scope to perform various diagnostic and therapeutic procedures.
4. Outer Sheath Material: The outer sheath material is a tough thermoplastic elastomeric covering with a larger inside diameter than the inner sheath. The diameter difference creates an annular gap between the inner and outer sheath for the passage of air. The annular gap is sealed at the proximal end with open communication to the luer lock connector and is sealed at the distal end, with an open port to allow pulsed or continuous air to exit the sheath.
5. Sheath Window: A patented, optically clear window that covers the distal tip of the endoscope for visualization during sensory testing.
6. Protector Tube: The Protector Tube protects Sheath from external damage during shipping, storage and sheath loading.

SENSORY SHEATH LOADING

warning

- Do not attempt to install a damaged Sheath onto your endoscope.
- It is extremely important to insert the endoscope carefully to avoid excessive bending of the insertion tube

Caution

- Ensure the insertion tube is clean and dry prior to loading the sheath.
 - DO NOT use force when installing or removing the Sheath. The endoscope should move freely inside the Sheath without force during loading and unloading procedures. Should ANY resistance be experienced DO NOT CONTINUE. Stop and refer to the Troubleshooting section of this manual for further instruction. Medtronic, Inc. will not be responsible for damage to endoscopes caused by excessive force.
1. Before inserting the endoscope, ensure that the insertion tube is clean and dry. If the insertion tube is tacky, apply a light coating of cornstarch to the entire length. Ensure the optics are free of debris by wiping the lens with alcohol and a lint-free cloth.

2. Peel back the top portion of the package, and grasp the Sheath by the Proximal Connector.
3. Hold the endoscope Control Body with the articulation lever facing towards you, (the user), and articulate to straighten the bending section to the neutral position.
4. Position the sheath so that the air channel and luer lock connector is directly opposite the articulation lever, (12:00 position as referenced by the user). *Reference figure 1.A*
5. Gently slide the Sheath onto the endoscope without twisting until the Proximal Connector is positioned on the strain relief of the endoscope or interface if using a Vision Sciences ENT-2000 Endoscope. Avoid touching the Sheath window to prevent image distortion. Be sure not to twist the sheath while loading so that the working channel remains as straight as possible. *Reference figure 1.B*
6. After the sheath is installed, do not touch the scope control body or eyepiece with the hand you use to manipulate the insertion tube. If you must rotate the diopter adjustment ring after the sheath is installed, do so using a gauze pad or replace your glove prior to performing the procedure.

SENSORY SHEATH - AP-4000 CONNECTION AND USE

NOTE: Secure attachment of Air Supply Tubing is essential to ensure optimal performance during Sensory Testing. Always verify a secure fit between tubing and connector prior to performing the procedure.

1. Remove the Protector Tube without twisting
2. Articulate the scope back and forth several times to fully seat the window.
3. If the window was not fully seated during installation or additional articulations, carefully grasp the sides of the window and gently push towards the insertion tube, being careful not to touch the window surface.
4. Connect the loaded sheath to the AP-4000 by securing the Air Supply Tubing luer connector to the luer of the Proximal Connector. The tubing is secured by threading on to the luer of the connector. Ensure that the tube fitting is threading properly on to the luer of the connector and verify firm attachment.
5. Connect the loaded sheath to the AP-4000 by securing the Air Supply Tubing luer connector to the luer lock connector of the proximal connector. The tubing is secured by threading on to the luer lock of the connector. Ensure that the tube fitting is threading properly on to the luer lock of the connector and verify firm attachment.
6. Continue to follow the instructions found in the AP-4000 Owner's Manual for System Operation, including, Pre-Use Checklist, Self Test, Calibration, Positioning the Endoscope, and Delivering the Air Pulse. Heed all additional instructions and information found in the Owner's Manual including Cleaning and Disinfection of the AP-4000 system.
7. Articulate the Sheath as required during the procedure

SENSORY SHEATH REMOVAL

CAUTION: Should difficulty be experienced in removing the endoscope from the sheath, do not force the endoscope out of the Sheath. Please stop and call Medtronic Customer Service Department at 800-874-5797 or 904-296-9600 EST 8:00am to 6:00pm Monday through Friday. Medtronic, Inc. is not responsible for damage to equipment caused by excessive force.

1. Once the endoscope procedure is completed and the endoscope is withdrawn from the patient, the air supply tube can be disconnected from the Sensory Sheath and the AP-4000.
2. Remove the Sheath from the endoscope in an appropriate area to prevent the spread of pathogenic material.
3. Activate the Articulation Lever to straighten the bending section to the neutral position
4. Grasp the Sheath Proximal Connector with the other hand and gently rotate and pull to disengage it from the scope strain relief or interface by 1-2 cm.
5. If the Sheath buckles or folds over on itself, slide the Proximal Connector back toward the scope strain relief, or interface, far enough to eliminate wrinkling.

Reference figure 2.A

6. Gently grasp the middle of the Sheath and begin to slide it toward the distal end of the instrument in order to disengage the clear plastic window from the tip of the scope. (If buckling occurs by mistake, do not push any further. Gently return the connector to the interface to eliminate the buckling and proceed as above.) Once the clear plastic window is disengaged from the tip of the scope, grasp the plastic window portion of the Sheath and slide the Sheath off the scope. If the Sheath does not slide freely, a small amount of alcohol may be used to separate the Sheath from the endoscope insertion tube. Hold the endoscope vertically and place the alcohol inside the Proximal Connector using a syringe without needle, eye dropper, pipette, etc.

Reference figure 2.B

7. Place endoscope on a clean surface.
8. Carefully discard Sheath per hospital policy.

EN

RECOMMENDED CLEANING & DISINFECTION PROCEDURE:

NOTE: Because the possibility exists that a Sheath could come into contact with contaminated surfaces, the user should develop and follow a prophylactic routine which includes exercising care when handling the sheathed and unsheathed scope, and a cleaning routine which includes the following steps after each clinical procedure:

After sheath removal, inspect the insertion tube and distal portion of the scope to ensure these areas are dry. If there is no moisture on the insertion tube, use the following cleaning and disinfection procedure.

CLEANING PROCEDURE

- Clean the scope by gently washing all external surfaces with an appropriate instrument detergent.
- Rinse all surfaces thoroughly with clean lukewarm water.
- Dry the scope thoroughly after rinsing.

DISINFECTION PROCEDURE

- Wipe the entire fiberscope with gauze soaked in 70% alcohol
- Ensure that all external surfaces of the scope are thoroughly dried prior to installing another sheath or storing the scope.

After sheath removal from the scope, if there is any moisture present following inspection of the insertion tube this is an indication that there may have been a leak during the endoscopy procedure, providing the endoscope was dry when the sheath was installed. If a leak is suspected, the scope must be disinfected (high-level) or sterilized following the manufacturer's instructions.

NOTE: Leak test the nasopharyngoscope according to manufacturer's instructions, prior to submerging the scope in any fluid.

RECOMMENDED CLEANING & DISINFECTION PROCEDURE:

NOTE: Because the possibility exists that a Sheath could come into contact with contaminated surfaces, the user should develop and follow a prophylactic routine which includes exercising care when handling the sheathed and unsheathed scope, and a cleaning routine which includes the following steps after each clinical procedure:

After sheath removal, inspect the insertion tube and distal portion of the scope to ensure these areas are dry. If there is no moisture on the insertion tube, use the following cleaning and disinfection procedure.

Cleaning Procedure

- Clean the scope by gently washing all external surfaces with an appropriate instrument detergent.
- Rinse all surfaces thoroughly with clean lukewarm water.
- Dry the scope thoroughly after rinsing.

Disinfection Procedure

- Wipe the entire fiberscope with gauze soaked in 70% alcohol
- Ensure that all external surfaces of the scope are thoroughly dried prior to installing another sheath or storing the scope.

After sheath removal from the scope, if there is any moisture present following inspection of the insertion tube this is an indication that there may have been a leak during the endoscopy procedure, providing the endoscope was dry when the sheath was installed. If a leak is suspected, the scope must be disinfected (high-level) or sterilized following the manufacturer's instructions.

NOTE: Leak test the nasopharyngoscope according to manufacturer's instructions, prior to submerging the scope in any fluid.

TROUBLESHOOTING

PROBLEM	PROBABLE CAUSE	ACTION
Difficulty inserting endoscope into Sheath	Incorrect Sheath type. Endoscope damaged resulting in excessive diameter; (e.g. rolled over angle section cover). Axis of endoscope not aligned with Sheath causing friction. Endoscope very sticky or tacky. Bending section of the endoscope not straight. Sheath is torn or punctured.	Ensure Sheath type matches your Endoscope Model Repair endoscope. Examine with Insertion Tube Evaluation Template. Align endoscope. Clean endoscope insertion tube with alcohol wipe and allow to dry. Apply a light coating of cornstarch to the entire length of the insertion tube. Straighten bending section. Remove and replace with new Sheath. Contact Customer Service.
Glare evident when viewing through endoscope	Sheath Window not in contact with tip of endoscope due to: endoscope not fully inserted. OR Damaged (enlarged) endoscope tip.	Ensure endoscope is fully inserted when loading. Repair endoscope.

PROBLEM	PROBABLE CAUSE	ACTION
Difficulty REMOVING endoscope from Sheath	The Sheath has folded over on itself or become caught under the Sheath Proximal Connector. OR The endoscopes angle section cover has become wrinkled or rolled over under the Sheath causing difficulty. OR Sheath loaded for extended period of time causing sticking to endoscope. OR Sheath is torn or punctured. Insertion tube is tacky.	1) Slide the proximal connector back toward the scope strain relief, or interface, far enough to eliminate wrinkling. 2) Gently grasp the insertion tube at the 10 cm mark (white line) and slowly rotate until Sheath separates from the scope. 3) If necessary, when the distal end of the Sheath is off the scope, it may be grasped to assist in removal. 4) If the Sheath does not slide freely, a small amount of alcohol may be used to lubricate the endoscope insertion tube. Hold the endoscope vertically and place the alcohol inside the Proximal Connector using a syringe without needle, eye dropper, pipette, etc. NOTE: Never use alcohol or lubricants during sheath installation. 5) Contact Customer Service. Apply a light coating of cornstarch to the entire length of insertion tube, prior to next usage.
Issues with Air Pulse Delivery	Air Supply Tube is not attached at both ends or is not fully secured Numerous possible issues with the AP-4000 unit or the air supply connections Damaged sheath	Ensure the Air Supply Tube is properly connected to the Sensory Sheath and to the AP-4000 Air Out connection by rotating the luer connector clockwise until fully seated. Refer to the Preparation Section of the AP-4000 Owner's Manual Refer to all appropriate sections of the AP-4000 Owner's Manual Do not use a sheath or Sensory instrumentation that is damaged. Failure to observe this precaution may result in malfunction during the procedure and serious patient injury and/or damage to the System.

RETURNS AND/OR REPAIRS

Contact Medtronic Xomed Customer Service at 800-874-5797 to obtain a Return Goods Authorization number (RGA#) and shipping information. Please have the original invoice number or purchase order number available to assist in verifying warranty information.

The RGA# should be prominently displayed on the box and included on all paperwork enclosed with the return and/or repair. All products returned should be safely packed in protective wrapping. Customer must supply the Purchase Order number; the correct shipping and billing address; and either a completed Repair Order Form or a statement of the problem or reason for return.

CUSTOMER SERVICE INFORMATION

For further information regarding the use of this product or to report any problems, please contact Medtronic Xomed using the appropriate information provided on the blue and white contact information card packaged with each device; or contact your local distributor

U.S. Customer Service Contact Information

800-874-5797 or
904-296-9600
Monday - Friday:
8:00 AM – 6:00 PM E.S.T.
medtronicENT.com

Medtronic Xomed, Inc., Technical Support:

800-872-9877
904-296-6448 (FAX)
Monday - Friday
8:00 AM – 5:00 PM E.S.T.
medtronicENT.com

LIMITED WARRANTY

A. This LIMITED WARRANTY provides assurance for the customer who purchases a Medtronic Xomed Product (hereinafter the "Product") that should the Product fail to function to Medtronic Xomed's published specifications during the term of this LIMITED WARRANTY (one year from the date of shipment for new Product, 90 days from date of shipment for refurbished or used Product), Medtronic Xomed will either replace, repair, or issue a credit (adjusted to reflect the age of the Product) for the Product or any portion thereof. This LIMITED WARRANTY is extended only to the buyer purchasing the Product directly from Medtronic Xomed or from its affiliate or its authorized distributor or representative.

B. To qualify for this LIMITED WARRANTY, the following conditions must be met:

- (1) The Product must be used on or before its "Use By" or "Use Before" date, if applicable.
- (2) The Product must be used in accordance with its labeling and may not be altered or subjected to misuse, abuse, accident or improper handling.
- (3) Medtronic Xomed must be notified in writing within thirty (30) days following discovery of a defect.

FR

(4) The Product must be returned to Medtronic Xomed within thirty (30) days of Medtronic Xomed receiving notice as provided for in (3) above.

(5) Upon examination of the Product by Medtronic Xomed, Medtronic Xomed shall have determined that: (i) the Product was not repaired or altered by anyone other than Medtronic Xomed or its authorized representative, (ii) the Product was not operated under conditions other than normal use, and (iii) the prescribed periodic maintenance and services have been performed on the Product.

C. This LIMITED WARRANTY is limited to its express terms. THIS LIMITED WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED WHETHER STATUTORY OR OTHERWISE, INCLUDING ANY IMPLIED WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. In no event shall Medtronic Xomed be liable for any consequential, incidental, prospective or other similar damage resulting from a defect, failure, or malfunction of the Product, whether a claim for such damage is based upon the warranty, contract, negligence or otherwise.

D. The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. Users may benefit from statutory warranty rights under legislation governing the sale of consumer goods. If any part or term of this LIMITED WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable, or in conflict with applicable law, the validity of the remaining portion of the LIMITED WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this LIMITED WARRANTY did not contain the particular part or term held to be invalid

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Medtronic

Medtronic

Medtronic Xomed, Inc.
6743 Southpoint Drive North
Jacksonville, Florida 32216
USA

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Custodian:

Alexandra Oliver

Alexandra Oliver

Regulatory Affairs Specialist

Medtronic Xomed, Inc.

01/10/2019
Date

Appendix to IFU

Slide-On EndoSheath System for Flexible ENT Scopes

ATTESTATION


STATE OF FLORIDA
DUVAL COUNTY

On this 10th day of JANUARY, 2019, I attest that the preceding or attached documents are true, exact, complete, and unaltered photocopies made by me of the following, Appendix to IFU, issued by Medtronic Xomed, Inc., presented to me by the document's custodian, Alexandra Oliver, and, to the best of my knowledge, that the photocopied document is neither a public record nor a publicly recordable document, certified copies of which are available from an official source other than a notary public.

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

The Channel Sheaths, as shown in the Figure, consists of a translucent, highly flexible polyurethane sleeve with a transparent thermoplastic polymer (PETG) distal window. The window is attached to the sheath using UV curable adhesive. A PVC proximal fitting provides a means of attaching the sheath to the flexible endoscope. The inner diameter of the proximal fitting stretches slightly to provide an interference fit over the appropriate part of the scope. The Channel Sheaths incorporate an internal working channel into the disposable sheath. This channel has been designed to accept tools that are less than 2.1mm in diameter. The working channel to the basic sheath is made from polyethylene material except for a portion near the distal end to allow for flexibility during articulation while maintaining a low-friction surface for the passage of biopsy instruments. This section that allows passage of instruments is made out of a flexible, high density, porous PTFE material. The sheaths also have a PVC Proximal scope connector. UV cure adhesive is used for all sheath bonds.

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1. Proximal Connector: The Proximal Connector interfaces with the endoscope at the controlbody-strain relief junction to secure the sheath to the scope.

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2. Luer Lock Connector: The standard luer lock connector provides a means of securing the AP-4000 air supply tube to the sheath for sensory testing.

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© 2014 Medtronic, Inc.
EndoSheath and Slide-On are registered trademarks of Vision Sciences, Inc.
Pentax is a registered trademark of Asahi Kogaku Kogyo Kabushiki Kaisha.
Olympus is a registered trademark of Olympus America, Inc.