
## Approval

Frank Zanow, CEO


Date: 12.04.2019

visor2 User Guide

# visor2

TMS Neuronavigation System

## User Guide

*This document is official since 2017-05-08 for visor2, version 2.3.3*

*Document revision: 27*

**This document is applicable to users in the Russian Federation**

eemagine Medical Imaging Solutions GmbH  
Gubener Straße 47  
D-10243 Berlin - Germany  
T: +49 3029048404 F: +49 3029048405  
IBAN: DE73 1204 0000 0458 8885 00  
Registration: Amtsgericht Charlottenburg  
Number: HRB 71287 B  
VAT: DE201608185

Approved: DRN-PDO-1250

## Manufacturer

### **eemagine Medical Imaging Solutions GmbH**

Gubener Str. 47

D-10243 Berlin, Germany

E-Mail: [info@eemagine.com](mailto:info@eemagine.com)

## The authorized representative of the manufacturer in Russia

### **Neurosoft Limited Liability Company**

#### **Neurosoft LLC**

5, Voronin str., Ivanovo,

153032, Russia

E-mail: [info@neurosoft.com](mailto:info@neurosoft.com)

## Important notices:

The visor2 software is a CE class IIa device according to the MDD Annex IX.

Copyright © 2011-2017 eemagine Medical Imaging Solutions GmbH. All rights reserved.

### **Warning**

*No part of this document may be copied or reproduced without the express permission of the authors.*

# TABLE OF CONTENTS

|          |                                                               |           |
|----------|---------------------------------------------------------------|-----------|
| <b>1</b> | <b>General System Information</b>                             | <b>6</b>  |
| 1.1      | Intended Use                                                  | 6         |
| 1.2      | Contraindications and potential side effects                  | 6         |
| 1.3      | visor2 System components                                      | 7         |
| 1.3.1    | visor2 software                                               | 7         |
| 1.3.2    | 3D Tracking subsystem                                         | 7         |
| 1.3.3    | amplifier subsystem for electrophysiological signal recording | 7         |
| 1.3.4    | Tracking tools                                                | 7         |
| 1.3.5    | Trigger box                                                   | 8         |
| 1.3.6    | Isolation transformer                                         | 8         |
| 1.3.7    | Visor2 compartment                                            | 8         |
| 1.3.8    | cart for transportation                                       | 8         |
| 1.3.9    | EMG Sensor cables with snap-on contact                        | 8         |
| 1.4      | Compatible Components                                         | 9         |
| 1.4.1    | Compatibility of TMS device                                   | 9         |
| 1.4.2    | Compatible electrodes and Sensor cables                       | 9         |
| 1.4.3    | Compatibility of computer, monitor, keyboard                  | 9         |
| 1.4.4    | Compatibility of other extensions                             | 10        |
| 1.5      | System Overview with all Mechanical Components                | 11        |
| 1.5.1    | The tracking tools                                            | 12        |
| 1.5.2    | The visor2 Cart                                               | 14        |
| 1.5.3    | Isolation transformer                                         | 14        |
| 1.6      | Safety Warnings                                               | 14        |
| 1.7      | Safety Measures                                               | 15        |
| 1.7.1    | General                                                       | 15        |
| 1.7.2    | Amplifier for electrophysiological signal recording           | 15        |
| 1.7.3    | Optical Tracking System                                       | 16        |
| 1.8      | Operating Environment Requirements                            | 17        |
| 1.9      | Disposing of waste equipment                                  | 17        |
| 1.10     | Contact Information                                           | 17        |
| <b>2</b> | <b>General information for visor2 software</b>                | <b>18</b> |
| 2.1      | Intended Use                                                  | 18        |
| 2.2      | Contraindications                                             | 18        |
| 2.3      | Classification                                                | 18        |
| 2.4      | Features                                                      | 18        |
| 2.5      | Limitations                                                   | 19        |
| 2.6      | Safety Warnings                                               | 19        |
| 2.7      | Safety Measures                                               | 20        |
| 2.7.1    | Special Recommendations for recording                         | 20        |
| 2.8      | Patient Data Management                                       | 20        |
| 2.9      | Application Workflows                                         | 21        |
| <b>3</b> | <b>Setup Patient</b>                                          | <b>23</b> |
| 3.1      | Creating a new patient entry                                  | 23        |

|       |                                                           |    |
|-------|-----------------------------------------------------------|----|
| 3.2   | Edit existing patient <b>data</b> .....                   | 24 |
| 3.2.1 | Defining <b>new targets</b> .....                         | 24 |
| 3.2.2 | Provide a <b>new MRI</b> .....                            | 24 |
| 3.3   | <b>Load MRI</b> .....                                     | 24 |
| 3.3.1 | Standard <b>MRI</b> .....                                 | 24 |
| 3.3.2 | Patient-specific <b>MRI</b> .....                         | 25 |
| 3.3.3 | <b>MRI Requirements</b> .....                             | 25 |
| 3.4   | <b>Specify MRI Fiducials</b> .....                        | 26 |
| 3.4.1 | Talairach <b>markers</b> .....                            | 27 |
| 3.5   | <b>Create Model</b> .....                                 | 28 |
| 3.5.1 | Manual Brain Segmentation .....                           | 29 |
| 3.5.2 | Automatic Brain Segmentation .....                        | 32 |
| 3.6   | <b>Define Targets</b> .....                               | 33 |
| 3.6.1 | Single <b>targets</b> .....                               | 33 |
| 3.6.2 | <b>Target grid</b> .....                                  | 34 |
| 3.7   | <b>File structure</b> .....                               | 35 |
| 4     | <b>TMS Settings</b> .....                                 | 36 |
| 4.1   | <b>Coil selection</b> .....                               | 36 |
| 4.2   | Before the first calibration <b>of a new coil</b> .....   | 37 |
| 4.3   | <b>Calibration procedure</b> .....                        | 37 |
| 4.4   | Step 1: Calibrate Tool Faces .....                        | 37 |
| 4.5   | Step 2 (Pointer): Record calibration points .....         | 38 |
| 4.6   | Step 2 (Calibration Board): Calibrate Control Point ..... | 39 |
| 4.7   | Step 3 (Calibration Board): Calibrate Coil Center .....   | 40 |
| 5     | <b>EMG / EEG Recording</b> .....                          | 42 |
| 5.1   | <b>Recording Settings</b> .....                           | 42 |
| 5.2   | <b>Amplifier Calibration</b> .....                        | 43 |
| 6     | <b>Neuronavigation</b> .....                              | 46 |
| 6.1   | Choose coil and <b>stimulator</b> .....                   | 46 |
| 6.2   | <b>Validate coil</b> .....                                | 47 |
| 6.3   | <b>Select Patient</b> .....                               | 48 |
| 6.3.1 | Read stimuli from previous <b>session</b> .....           | 48 |
| 6.3.2 | Read target markers .....                                 | 49 |
| 6.3.3 | Session with motor mapping (EMG) .....                    | 50 |
| 6.4   | <b>Define Coordinate System</b> .....                     | 50 |
| 6.5   | <b>Head Shape Points Generation</b> .....                 | 51 |
| 6.6   | <b>NeuroNavigation</b> .....                              | 53 |
| 6.6.1 | Neuronavigation main <b>workflow</b> .....                | 53 |
| 6.6.2 | Neuronavigation with EMG Recording .....                  | 54 |
| 6.6.3 | The induced electric <b>field</b> .....                   | 56 |
| 6.6.4 | Reproduce stimuli and Stimulate at targets .....          | 57 |
| 6.6.5 | <b>rTMS</b> .....                                         | 61 |
| 6.6.6 | Online EMG <b>Settings</b> .....                          | 62 |
| 6.6.7 | <b>Replace Coil</b> .....                                 | 63 |
| 6.6.8 | <b>Validate</b> .....                                     | 63 |

|                                                                                         |                                                   |     |
|-----------------------------------------------------------------------------------------|---------------------------------------------------|-----|
| 6.6.9                                                                                   | Finding the motor hot spot .....                  | 65  |
| 6.6.10                                                                                  | Finding the hot spot <b>of</b> more muscles ..... | 65  |
| 6.6.11                                                                                  | Stimulating at the found hot spot .....           | 67  |
| 6.7                                                                                     | Terminating the <b>session</b> .....              | 68  |
| 6.7.1                                                                                   | Export to DICOM .....                             | 68  |
| 7                                                                                       | Offline <b>Review</b> .....                       | 70  |
| 7.1                                                                                     | Review <b>Session</b> .....                       | 70  |
| 7.1.1                                                                                   | Load Session Data .....                           | 70  |
| 7.1.2                                                                                   | Review Session .....                              | 71  |
| 7.1.3                                                                                   | Export to DICOM .....                             | 73  |
| 7.2                                                                                     | Review Exported <b>DICOM</b> .....                | 74  |
| Appendices                                                                              | .....                                             | 76  |
| Appendix 1: Montages                                                                    | .....                                             | 76  |
| The Montage Editor                                                                      | .....                                             | 76  |
| Appendix 2: MRI coordinate systems                                                      | .....                                             | 79  |
| Appendix 3: Licensing                                                                   | .....                                             | 86  |
| License types                                                                           | .....                                             | 86  |
| Installation of a single license                                                        | .....                                             | 86  |
| Trial <b>license</b>                                                                    | .....                                             | 86  |
| Licensing support                                                                       | .....                                             | 87  |
| Appendix 4: Supported data formats                                                      | .....                                             | 87  |
| Appendix 5: Warranty, Support and Reclamations                                          | .....                                             | 88  |
| Appendix 6: Shortcuts and mouse wheel functions                                         | .....                                             | 90  |
| Appendix 7: System components and guidance related to EMC                               | .....                                             | 91  |
| Appendix 8: Measured, Displayed and Essential Parameters                                | .....                                             | 96  |
| Appendix 9: Explanation <b>of</b> symbols as used in this manual and/or on system parts | .....                                             | 98  |
| Appendix 10: Maintenance of the system and disinfection                                 | .....                                             | 99  |
| Appendix 11: Packing, transportation and storage                                        | .....                                             | 100 |

# 1 GENERAL SYSTEM INFORMATION

This document serves as the User Guide for the **visor2 system** including the **visor2 software**. It contains important information about the operation of the system, its specification and limitations. Especially, it contains several safety warnings that are **to** be obeyed at any time.

This document is intended to introduce you to the basic use of the system. Example cases will help you to start a TMS neuro navigation session, set markers, targets, create and navigate through the brain. You will also learn how to perform more advanced procedures such as offline motor mapping. About two hours are required to become familiar with the main features of the product, following the steps as described in this document.

**visor2 system** consists of the **visor2 software**, selected components and accessories as required to operate the system.

This chapter will focus on the aspects relevant for the system in below, describing the components, main system features and related safety warnings. The next chapter, and most of the rest of the manual, will focus on the use of the software, where necessary describing the interaction with other parts of the system.

## 1.1 INTENDED USE

The **visor2 system** is used to aid 3D navigation for positioning of Transcranial Magnetic Stimulation (TMS) coils over preselected brain regions in clinical procedures (including mapping and speech mapping) or scientific research. The compatible amplifier is used for measuring electrophysiological signals. The system is used with compatible TMS stimulator.

The systems can be used in different health care facilities including clinics, hospitals, health centers, ambulance centers, specialized medical facilities including prevention centers, medicine centers for emergency, military and medical expertise centers, neurosurgery clinics, and research labs of research facilities.

The **visor2 system** must be operated by qualified medical personnel trained to operate on it and knowledgeable of application in neurophysiology, transcranial magnetic stimulation and electromyography (if motor cortex mapping is used). Before attempting to use this device, the operator must read and understand the contents of the user manual and other accompanying documents. The system must not be used during surgeries.

The **visor2 system** does not have age restrictions when used for its intended purpose, but age restrictions can be prescribed by the use of a magnetic stimulator. Before use, it is necessary to read in details the information given in the technical manual of the magnetic stimulator used together with the **visor2 system**.

## 1.2 CONTRAINDICATIONS AND POTENTIAL SIDE EFFECTS

**visor2** is a system aiming at supporting you in positioning a TMS coil over anatomically or functionally interesting areas. The system optionally records and processes electrophysiological signals to guide the user in the coil placement. **visor2** is a 3D navigation system, and, as such, has no contraindications. The TMS stimulator you are using has contraindications.

**Warning:** Before performing transcranial magnetic stimulation, read carefully the complete list of contraindications and possible side effects in the documentation of your TMS stimulator.

### 1.3 VISOR2 SYSTEM COMPONENTS

The following components are part of the **visor2 system**:

#### 1.3.1 VISOR2 SOFTWARE

The core component **of** visor2 system is the **visor2 software**

|                                 |                                         |
|---------------------------------|-----------------------------------------|
| <b>Product name(s):</b>         | visor2 software (on USB stick)          |
| <b>Product code (eemagine):</b> | 666601.01 (USB stick)                   |
| <b>Product name(s):</b>         | visor2 software license (version 2.3.3) |
| <b>Product code (eemagine):</b> | LN-114, LN-131                          |

#### 1.3.2 3D TRACKING SUBSYSTEM

The following optical tracking systems are fully supported for combination with the **visor2 software** in **visor2 system**

|                                 |                                                                                                                    |
|---------------------------------|--------------------------------------------------------------------------------------------------------------------|
| <b>Product name(s):</b>         | NDI Polarís VICRA Position Sensor, M00066.C14                                                                      |
| <b>Product code (eemagine):</b> | M00150 (camera sensor), M000261 (USB converter), M000262 (power supply) or equivalent (compliant with IEC 60601-1) |

#### 1.3.3 AMPLIFIER SUBSYSTEM FOR ELECTROPHYSIOLOGICAL SIGNAL RECORDING

The following amplifiers are supported for combination with the **visor2 software** in **visor2 system**

|                                 |                                                                                          |
|---------------------------------|------------------------------------------------------------------------------------------|
| <b>Product name(s):</b>         | eego amplifier, EE-410.A                                                                 |
| <b>Product code (eemagine):</b> | EE-410 (eego amplifier), M01160 (USB cable), LE-210.01 (USB stick with software package) |

#### 1.3.4 TRACKING TOOLS

The following tracking tools are fully supported for combination with the **visor2 software** in **visor2 system**

|                                 |              |
|---------------------------------|--------------|
| <b>Product name(s):</b>         | Pointer tool |
| <b>Product code (eemagine):</b> | NT-115       |

|                                 |                                                                                  |
|---------------------------------|----------------------------------------------------------------------------------|
| <b>Product name(s):</b>         | Reference tracker tool, NT-116                                                   |
| <b>Product code (eemagine):</b> | NT-103 (head reference tool),<br>NT-144 (clip, mounting tool for head reference) |

|                                 |                                            |
|---------------------------------|--------------------------------------------|
| <b>Product name(s):</b>         | Coil calibration board standard for visor2 |
| <b>Product code (eemagine):</b> | NT-104                                     |
| <b>Product name(s):</b>         | Coil calibration board for Neurosoft       |
| <b>Product code (eemagine):</b> | NT-142                                     |

|                                 |                                                         |
|---------------------------------|---------------------------------------------------------|
| <b>Product name(s):</b>         | Coil tracker set for Neurosoft, NT-135                  |
| <b>Product code (eemagine):</b> | NT-112 (coil tracker tool), NT-141 (coil tracker mount) |

Especially relevant are the following tracking tools:

Calibration board. The board is used for determining the position and orientation of the coil's tracking tool relative to the coil itself. This needs to be done after fixation of the tracking tool. Note that for some coils, special calibration boards are provided. Note that Mag&More coils can also be calibrated without a calibration board. Those coils have position markers directly on the coil and no calibration board is required.

- Reference tracking tool on patient's head: **NT-103/NT-116**. The reference tool allows tracking head movements and makes sure that measured positions remain relative to the patient's head during sessions of electrode digitization or TMS stimulation.
- Pointer tracking tool **NT-115** to mark nasion-ear-positions and other head shape points.

### 1.3.5 TRIGGER BOX

The following device is used together with the vIsor2 software in vIsor2 system.

|                          |                            |
|--------------------------|----------------------------|
| Product name(s):         | Trigger box set            |
| Product code (eemagine): | NH-134                     |
| Product name(s):         | BNC to MCX Trigger Adapter |
| Product code (eemagine): | XS-294                     |

### 1.3.6 ISOLATION TRANSFORMER

The following device is used in the vIsor2 system for isolation:

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| Product name(s):         | Isolation Transformer                             |
| Product code (eemagine): | M01646 or equivalent (compliant with IEC 60601-1) |

### 1.3.7 VISOR2 COMPARTMENT

The following item is used in the vIsor2 system to assemble smaller system components in a defined way and **for mounting to a cart solution:**

|                          |                              |
|--------------------------|------------------------------|
| Product name(s):         | vIsor2 Isolation compartment |
| Product code (eemagine): | M01592 or equivalent         |

### 1.3.8 CART FOR TRANSPORTATION

The vIsor2 system is based on a mobile cart:

|                          |                                                                       |
|--------------------------|-----------------------------------------------------------------------|
| Product name(s):         | vIsor2 cart set                                                       |
| Product code (eemagine): | M01596 or equivalent (compliant with relevant clauses in IEC 60601-1) |

### 1.3.9 EMG SENSOR CABLES WITH SNAP-ON CONTACT

The vIsor2 system includes patient lead cables connecting to the amplifier:

|                          |                                                                                                  |
|--------------------------|--------------------------------------------------------------------------------------------------|
| Product name(s):         | 2 Channels Bipolar Adapter with Shielded Flat Snap Electrodes for Bipolar EE-410 eego amplifiers |
| Product code (eemagine): | XE-292                                                                                           |

## 1.4 COMPATIBLE COMPONENTS

### 1.4.1 COMPATIBILITY OF TMS DEVICE

The vIsor2 system fully supports TMS stimulator devices in combination with TMS coils. Compliant systems are required to meet IEC 60601-1 и IEC 60601-1-2.

### 1.4.2 COMPATIBLE ELECTRODES AND SENSOR CABLES

The vIsor2 system can be used in combination with disposable adhesive Ag/AgCl electrodes for medical application registered in RA according to the established procedure:

|                                 |                                                                        |
|---------------------------------|------------------------------------------------------------------------|
| <b>Product name(s):</b>         | ECG/EMG electrode                                                      |
| <b>Product code (eemagine):</b> | XE-164 or equivalent (adhesive electrodes - flat-snap binder, Ag/AgCl) |

|                                 |                                                                          |
|---------------------------------|--------------------------------------------------------------------------|
| <b>Product name(s):</b>         | vIsor2 starter kit                                                       |
| <b>Product code (eemagine):</b> | M00064 or equivalent (passive sphere markers + ECG electrodes compliant) |

### 1.4.3 COMPATIBILITY OF COMPUTER, MONITOR, KEYBOARD

The vIsor2 system can be used in combination with the following devices.

Desktop PC-based or notebook PC-based system:

|                                 |                                                                                                                                                                                                                                                                                                        |
|---------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product name(s):</b>         | bluechip BUSINESSline S3300 PC                                                                                                                                                                                                                                                                         |
| <b>Product code (eemagine):</b> | XS-611, or equivalent (compliant with IEC 60950, IEC 62311), with Intel Core Duo 3GHz processor or better, 2 GB RAM, Video card: DirectX 9.0 Pixel/Vertex shader 2.0 or better, VGA interface port, 6x USB 2.0 ports, MS Windows 7, 8.1 or 10 OS in "Pro" 32 or 64 bit mode, dimensions 425x326x95 mm) |

#### Monitor:

|                                 |                         |
|---------------------------------|-------------------------|
| <b>Product name(s):</b>         | Monitor (24") TFT (LED) |
| <b>Product code (eemagine):</b> | M00345                  |

#### Keyboard:

|                                 |                                                                    |
|---------------------------------|--------------------------------------------------------------------|
| <b>Product name(s):</b>         | keyboard with Touchpad, USB                                        |
| <b>Product code (eemagine):</b> | M00634, or equivalent (compliant with IEC 60950), IP22 and higher) |

Alternatively, other computers, displays and keyboards can be used when they meet the requirements defined in *Appendix 7: System components*.

## 1.4.4 COMPATIBILITY OF OTHER EXTENSIONS

**Remote control:**

|                                 |                                                            |
|---------------------------------|------------------------------------------------------------|
| <b>Product name(s):</b>         | Wireless remote control                                    |
| <b>Product code (eemagine):</b> | XH-193 or equivalent (compliant with IEC 60950, IEC 62479) |

|                                 |                                                                                                                                                                                                                                          |
|---------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product name(s):</b>         | Webcam                                                                                                                                                                                                                                   |
| <b>Product code (eemagine):</b> | M00263 (webcam), and USB 2.0 extension 2m (M00084), or equivalent (compliant with IEC 60601-1 and IEC 60601-1-2, compatible with MS Windows 7, 8, or 10, minimum resolution 720p, min. 24 frames/s, autofocus and auto image adjustment) |

|                                 |                                                                                                                                                                                                                           |
|---------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product name(s):</b>         | Microphone                                                                                                                                                                                                                |
| <b>Product code (eemagine):</b> | M00604, and Rode VC1 Stereo MiniInke 3,5mm (M00605), Extension cable 3m (M00610), Stimulation Headphones (XS-110), or equivalent (compliant with IEC 60950, IEC 62311, frequency range 60Hz – 18kHz, dynamic range >80dB) |

## 1.5 SYSTEM OVERVIEW WITH ALL MECHANICAL COMPONENTS

The visor2 system is shown in the image below:



FIGURE 1: IMAGE OF THE VISOR2 SYSTEM COMPONENTS COMBINED WITH A MONITOR AND KEYBOARD.

Figure 2 below shows the example of the electrical and functional block diagram of the visor2 system (system components shown in blue) with compatible parts like a cart, computer and stimulator (shown in grey).

**FIGURE 2: VISOR2 SYSTEM OVERVIEW AND COMPONENTS**

Make sure that your TMS Stimulator is always connected to Trigger IN1 (BNC IN1 on NH-134)!

### 1.5.1 THE TRACKING TOOLS

With the term *tracking tool* we mean any object than can be tracked by the infrared tracking camera. Tracking tools are provided with sphere markers, which reflect Infrared light emitted by IR-LEDs, attached to the position sensor. The geometrical configuration of a tracking tool is described in a “.rom” configuration files. The tracking tools provided with visor2 are illustrated below.


The available trackers are listed in a separate document "List of trackers" which can be sent on request sent to [info@neurosoft.com](mailto:info@neurosoft.com).

The **pointer tool** is used during the patient registration. You will be asked to digitize with the pointer the positions of different anatomical landmarks on the patient's head. The **coil tracker** has to be rigidly attached to the coil, and allows tracking the position and orientation of the coil during a navigated TMS session. The **reference tool** has to be fixed to the patient's head. In order to allow the patient to move during the TMS session, all positions are expressed relative to the position of the reference tool. The **calibration board** is used when calibrating a coil. After having calibrated a coil, the calibration board will not be necessary anymore during the TMS session.

- It is very important **that** the reference tool does not move **with** respect to the patient's head!
- Please make sure that you always use a calibration board compatible with your coil! The program will present you with a list of boards that can be used with your coil. Never use a different one!


Pointer tool NT-115

Coil tracker. The lower part is specific to your coil model. The wing part with the reflective spheres is universal: NT-112 for the first coil NT-143 for a second coil.

Reference tool NT-103. Note that it comes together with a clip for fixating it with a headband or on an EEG cap.

Calibration board NT-104 (standard). For a wide range of coils, we also offer tight fit calibration boards that allow an even easier calibration. Some coils do not require a calibration board but have markers directly on the coil body for registration of the coil tracker. For details ask your local distributor!

FIGURE 3: MAIN VISOR2 TRACKING TOOLS

### 1.5.2 THE VISOR2 CART

The system can optionally be provided **with** a cart, which facilitates the transport. The cart has two possible configurations:

- **Operation mode.** During operation, the wheels of **the** cart should be locked and the position sensor rigidly fixated. The position sensor should be positioned in such a way that the patient's head is always in the field **of** view and that there are no obstacles which might obstruct the line of sight.
- **Transportation mode.** Before moving the cart, please remove the position sensor, fold down the sensor mount arm and unlock the wheels. The cart must be maneuvered by holding the pull handle.

### 1.5.3 ISOLATION TRANSFORMER

Make sure to check the details provided in the user manual of the isolation transformer before putting the system into operation especially after transportation. **The transformer supports both 110V and 230V but it requires manual adjustment going from one to the other (among other measurers fuses need to be adjusted).** Technical trained personal is required for the adjustment. Contact [help@neurosoft.com](mailto:help@neurosoft.com) if in doubt about the correct setting.

## 1.6 SAFETY WARNINGS

- visor2 is to be used only by medical experts, trained on the use of the system, and aware of the visor2 Indications of use and limitations.
- Use only tracking tools supplied with the system. Use with utmost care; prevent accidental falling tools, as this may result in irreversible damage/breaks! Never use damaged tools!
- Use the pointer tool very carefully when digitizing the positions of the anatomical landmarks, especially preventing the pointed tip to cause harm. When placing the pointer with your hand, use the other to steady the patient's head.
- Handle all tracker tools very carefully. Avoid touching the reflecting sphere markers. Check for the integrity of all tools before each session with visor2. When replacing the sphere markers, press them completely to the bottom.
- Make sure you always use the correct calibration board for coil calibration. Please refer to your distributor for a full list of compatible hardware.
- Maneuver the visor2 cart carefully. Your full attention is required when moving the system, which may only be performed in TRANSPORTATION MODE. In order to prevent tip-over accidents, use your two hands, and make use of the handle; never push the cart sideways (SEE ALSO page 14).
- The use of the eego amplifiers EE-xxx as a part of Visor2 system is for surface EMG recordings only. When the visor2 system is not in use, please make sure to switch it off.
- **Make sure the mains voltage setting is correct based on requirements in your territory and adjust the setting for the isolation transformer to 230 V or 110 V.**
- The system needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in Appendix 7: System components ; the use of accessories other than specified may result in increased emissions or decreased immunity and may result in improper operation.
- The visor2 system should not be used with other equipment which is not a part of the system.

- The use **of** portable and mobile RF communications equipment can adversely affect the overall performance; do not use an operating cellular phone within 30 cm (12 inches) of the system parts, cables and electrodes to avoid excessive noise on the signals. See also Appendix 7: System components
- During EMC testing the essential performance was monitored based on an overall accuracy of coil position relative to the (artificial) head with an upper limit of the position mismatch (measured with a special test tool and a dedicated test workflow) of 10mm, and timing jitter maximum 500ms between TMS onset trigger and camera position stream, and maximum 5ms between TMS onset trigger and EMG stream.
- There is no special servicing procedure required to guarantee EMC over the lifetime of the device.

## 1.7 SAFETY MEASURES

### 1.7.1 GENERAL

- Do not use the system in an operating room.
- **vIsor2** must not be used as sole basis for the planning of surgical procedures.
- The **vIsor2 system** must not be used as a diagnostic tool for inspecting 3D medical images.
- The **vIsor2 system** is not certified for the recording of EEG.
- The system is not to be immersed in any liquid. If any liquids or moisture penetrates the system or a part thereof, remove the plug from the wall socket and have the system checked by the manufacturer.
- Do not store any liquids in and/or on the cart. Before starting to use the system, verify that there are no liquids stored in/on the cart.
- Make sure that the wall socket is well earthed.
- Inspect the power cord on a regular basis for damages. Do not operate the device with a damaged power cord.
- Do not plug in equipment or power cord extensions to the medical isolation transformer when provided as part of the **vIsor2 system**. Do not change the power connections. Do not block access to its main power switch. Also, see Appendix 7: System components for further regulatory information.
- Setup and switch on the equipment before the patient comes in; do not touch any electrical contact/connector simultaneously with the patient.
- Make sure the isolation transformer is correctly set to 110 V or 230 V depending on the mains voltage in your territory.
- Take care when arranging cables to avoid risks of entanglement or strangulation.
- **WARNING:** do not modify this equipment without authorization of the manufacturer.
- **ELECTRICAL SHOCK HAZARD:** Do not connect electrode inputs to earth ground. Connecting an earth ground might result in electrocution.
- Conductive parts of **ELECTRODES** and associated connectors for **APPLIED PARTS**, including the **NEUTRAL/GROUND/REFERENCE ELECTRODE**, should not contact other conductive parts including earth.
- **ELECTRICAL SHOCK HAZARD:** Servicing of the equipment is to be done by approved technicians only.

### 1.7.2 AMPLIFIER FOR ELECTROPHYSIOLOGICAL SIGNAL RECORDING

The amplifier comes with a dedicated user manual, in which you will find relevant safety warnings and precautionary measures.

#### 1.7.2.1 RECOMMENDATIONS FOR RECORDING

- Before the first use for EMG recording, each amplifier must be calibrated to the **vIsor2 system**.
- Check the signal quality of the EMG before the first stimulation in an EMG-TMS-session

- It is recommended to prevent overload on the CPU resources during EMG acquisition, as may be caused by network activity (e.g., internet browser, or e-mail program), or virus scanners.
- If possible, disconnect the PC from a network during EMG acquisition.
- If you must use a virus scanner, please check if you can disable it during recording, and verify its settings for any scheduled periodic activity.
- After installation and verification of the proper recording procedures, it is best to use the same USB port always for the connection of the amplifier/USB interface.
- Disable the Windows selective suspend feature for energy saving for every USB port and all visor2 devices connected (amplifier, trigger box, remote control etc.): 1. Plug the device into the USB slot 2. Go to the Windows device manager and locate your device there" 3. Select the properties of the device. 4. Select the "Power Management" tab and deactivate the entry "Allow the computer to turn off this device to save power" 5. Repeat steps 2 - 4 for every USB slot and every device.
- The use of many different USB devices on the recording PC may result in an instable system configuration. In that case, please disconnect devices, and/or remove the associated driver through the Windows "Device Manager" control.

### 1.7.3 OPTICAL TRACKING SYSTEM

- Do not use the system in the presence of flammable materials such as anesthetics, cleaning agents, and endogenous gases. Flammable materials may ignite, causing personal injury or death.
- Radio frequency communications equipment, including portable and mobile devices, may adversely affect the performance of the system.
- Do not use the system in an MRI environment.
- Do not transport or store the position sensor outside the recommended storage temperature range, as this may cause the system to go out of calibration. Reliance on data provided by an out of calibration position sensor may lead to inaccurate conclusions and may cause personal injury. A calibration procedure must be performed before using the position sensor after it has been transported or stored outside the recommended storage temperature range.
- **Do not use cables or accessories other than those listed in this guide.**
- All user maintenance must be done by appropriately trained personnel. Maintenance by untrained personnel may present an electric shock hazard.
- Do not attempt to bypass the grounding prong on the power cord by using a three-prong to two-prong adapter. The system must be properly grounded to ensure safe operation. Failure to do so presents an electric shock hazard.
- Do not immerse any part of the system or allow fluid to enter the equipment. If fluids enter any part of the system, they may damage it and present a risk of personal injury.
- Do not sterilize the system as this may cause irreversible damage to its components. Reliance on data provided by a damaged position sensor may lead to inaccurate conclusions.
- Do not use the position sensor without inspecting it for cleanliness and damage both before and during a procedure. Reliance on data provided by an unclean or damaged position sensor may lead to inaccurate conclusions.
- Do not use markers without inspecting them for cleanliness and damage both before and during a procedure. Reliance on data produced by unclean or damaged markers may lead to inaccurate conclusions.
- Do not obstruct the normal flow of air around the position sensor (for example, draping or bagging the position sensor). Doing so will affect the position sensor's operational environment, possibly beyond its recommended thresholds. Reliance on data provided by a position sensor that is outside of recommended thresholds may lead to inaccurate conclusions.

- Take precautions to ensure that people with restricted movement or reflexes (for example, patients undergoing medical procedures) do not look directly into the emitting aperture. Patients undergoing medical procedures may be restricted in the availability of adverse-effects reflexes (turning away eyes and/or head, closing eyes) due to pharmaceutical influences and/or mechanical restraints. The IR emitter on the position sensor emits radiation that is visible and may be harmful to the human eye. Direct viewing of the camera diode emission at close range may cause eye damage.
- Use of camera controls or adjustments other than those specified may result in hazardous radiation exposure.
- Before using the system, make sure the power LED on the position sensor has stopped flashing. The power LED will flash while the position sensor warms up.

## 1.8 OPERATING ENVIRONMENT REQUIREMENTS

- There must be a clear line of sight between the position sensor and the tools to be tracked.
- The tools must be inside the characterized measurement volume.
- Make sure that sources of background IR light in the 800 nm to 1100 nm range (e.g. sunlight, some operating room lights) are minimized. The position sensor is sensitive to IR light. Since the position sensor functions by detecting IR light reflected from, or emitted by, markers, other sources of IR light can interfere with the optical tracking system.
- Make sure that there are no large reflective surfaces within the field of view.
- The environmental conditions must be as follows: the atmospheric pressure must be 70 kPa to 106 kPa, the relative humidity must be 30% to 70%, and the temperature must be 10°C to 40°C.

## 1.9 DISPOSING OF WASTE EQUIPMENT

In Russia after the end of the amplifier service life it must be disposed according to the Sanitary Regulations and Norms 2.1.7.2790-10 class A (epidemiologically non-hazardous wastes that are approximate to urban ore).

The disposal of the system outside the Russian Federation is done in accordance with the current legislation of your region.

## 1.10 CONTACT INFORMATION

For more information and for reporting problems with visor2, please contact your local distributor, Neurosoft LLC using the e-mail: [help@neurosoft.com](mailto:help@neurosoft.com).

## 2 GENERAL INFORMATION FOR VISOR2 SOFTWARE

### 2.1 INTENDED USE

**vIsor2 software** is used to aid 3D navigation for positioning of Transcranial Magnetic Stimulation (TMS) coils over preselected brain regions in scientific research or clinical procedures. **vIsor2 software** is to be used with a compatible TMS stimulator, a compatible optical tracking system and further compatible accessories. Optionally, a compatible amplifier for measuring electrophysiological signals can be connected. The **vIsor2 software** allows to record and view electrophysiological signals from the patient. It visualizes the activated motor cortex area over the anatomical images of the patient's brain based on the analysis of the evoked muscle response after the application of a TMS pulse.

The **vIsor2 software** must be operated by a trained physician, and must not be used during surgical operations.

### 2.2 CONTRAINDICATIONS

See the contraindications for the vIsor2 system in 1.2.

### 2.3 CLASSIFICATION

The **vIsor2 software** is an active device to monitor the position of TMS coils and to help locating motor areas on the brain's surface as assistance in diagnosis. **vIsor2 software** is not intended to be used as therapeutic device, nor for controlling/monitoring the performance of such devices, and is intended to be used only with TMS devices in Class IIa and compatible EEG/EMG amplifiers in Class IIa. The **vIsor2 software** is not intended for direct diagnosis or monitoring of vital physiological processes/parameters.

### 2.4 FEATURES

As a core component of the **vIsor2** system, **vIsor2 software** allows for

- Import of anatomical and functional 3D image data (MRI, CT, fMRI)
- Semi-automated image segmentation of skin and brain
- Coil calibration for 3D navigation using a specially designed calibration board
- Registering the patient's head to the segmented skin using a combined landmark and surface matching algorithm
- Real-time display of the position of one or two TMS coils relative to the patient's 3D image data and segmentation
- Calculation and display of the induced electric field for TMS stimuli
- Display of stimulation sites, i.e. coil positions at the moment the TMS stimulus was invoked
- Navigated targeting at pre-defined targets or existing stimulation sites
- EMG recording and the display of the EMG stimulation responses

## 2.5 LIMITATIONS

- The accuracy **of** the localization of the coil and the other tracked tools is typically within 2mm (RMS), see also Appendix 7: System components for further accuracy information.
- The accuracy **of** the localization depends not only on the overall accuracy of the system, but also on your diligence when using it. Be aware **of** the following:
  - During the whole navigated TMS session, the position **of** the head tracker must not change. If the position **of** the head tracker changes with respect to the position of the patient's head, the patient registration is no more valid.
  - From the accuracy point of view, the patient registration is the most critical step. Be precise when digitizing the positions of the nasion and the left and right preauricular points and perform always the head shape based registration.
  - The segmentation of the scalp in the MRI directly influences the patient registration quality. Check the segmentation result for agreement with the MRI before using it for the neuronavigation session.
  - During coil calibration, be sure that the coil center matches the point of origin on the board. Any displacement will be reflected in the precision of the coil localization.
  - During each TMS session, be sure to validate your coil calibration. This ensures that the coil tracker has not shifted with respect to the coil by more than 2mm. If validation fails or in doubt, recalibrate your coil.
  - A geometrically distorted MR scan will decrease considerably the overall accuracy.
  - The calculation of the electric field induced by the TMS coil is based on an approximation of the head shape. The calculated values are only indicative.
  - The strength of the electric field is calculated on the vertexes of the curvilinear surface. The map is then obtained by interpolation. Be sure to generate a fine curvilinear mesh if you need a precise localization of the field maximum.
  - The probability of excitation of the underlying cells depends not only on the direction and the intensity of the applied field but also on the shape and orientation of the neurons.
  - Depending on the user settings, the extent of the area under a non-vanishing electric field may vary. Be aware that cells outside the area indicated by vIsor2 might be stimulated, and that cells within it might not.

## 2.6 SAFETY WARNINGS

- vIsor2 is to be used only by medical experts, trained on the use of the system, and aware of the vIsor2 indications of use and limitations.
- After loading an MRI scan, always check the correct left-right orientation. After loading an fMRI always, check that the left-right orientation is correct and that the functional scan is correctly co-registered with the anatomical one.
- The use of Analyze as import data format for MRI scans is discouraged, as there is ambiguity about the image orientation. If you need to work with Analyze Images, it is most important to check the left-right orientation for every patient and data set.
- After performing the registration, always verify its goodness. Use the pointer tool to indicate the positions of several anatomical landmarks, and verify if they are correctly displayed on the 3D view.
- If high accuracy of the motor mapping results is required, e.g. as input into pre-surgical procedures, verify that the head reference has not shifted during the session. Use the verification step at the end of the Neuronavigation workflow.
- Any EMG traces recorded in the motor mapping workflow should be reviewed carefully and EMG responses with artifacts should be excluded from further use.
- Motor mapping is only allowed to be performed by trained medical personnel.

- Surgical procedures must not be based solely on data acquired or recorded with visor2 including mapping.
- Perform a screen calibration **if** the absolute EMG / EEG amplitude is to be determined by the height of the EMG / EEG amplitude displayed on screen. Amplitudes on an uncalibrated screen do not match the true amplitudes of the EMG / EEG.

## 2.7 SAFETY MEASURES

### 2.7.1 SPECIAL RECOMMENDATIONS FOR RECORDING

- The visor2 system does not support a 'continuous' recording mode, because of a size limit of CNT files of 2GB. The user should be aware of practical restrictions for the duration of one recording file, for example:  
→ record maximally at most 4 hours with 32 channels at 2048Hz, equally so for 64 channels at 1024Hz; After this time, you need to start a new session.

Notes:

- This is considering 'normal' recordings. As the EEProbe \*.cnt format compresses the data, these values depend on the information content of the EEG and might vary between patients and sessions;
- If you require longer recordings, you can enable the option "Extended mode" in the CNT Export Settings (Tools → Options → CNT Export Settings). This will store the data in the extended EEProbe format, allowing for 17TB (terra bytes, 1TB = 1000GB) of recording data and virtually endless recordings. However, please be aware that this format is not yet supported by most third-party software. For more information, please contact [help@neurosoft.com](mailto:help@neurosoft.com).

## 2.8 PATIENT DATA MANAGEMENT

The first time the **Visor 2 software** is run, it is important to select the root directory where all patient-related data will be stored. The root directory can be selected by changing the Default MRI File Folder under Tools → Options → MRI Import Menu. **It is extremely important to set the root directory only once.** If this directory is changed, access to existing patient data will be lost (these may be retrieved by copying all patient directories into the new root directory).


FIGURE 4: PATIENT DATA STRUCTURE

For each patient, a new patient directory will be created in the root directory, according to the structure shown in Figure 4.

It is important to understand that all navigation data refer to the underlying MRI that was used during data acquisition. That is, Session and Target data are only valid when analyzed using the current MRI. Therefore, if new MRI data is used for a later session, the current Target and Session data is archived together with their respective old MRI in a separate directory called "Archive", see section Load MRI.

## 2.9 APPLICATION WORKFLOWS

The visor2 software is structured in five main workflows, which are accessible through the corresponding toolbar, as displayed in Figure 5.

Setup Patient   TMS Neuronavigation   Offline Review   TMS Settings   Recording Settings

FIGURE 5: THE WORKFLOWS TOOLBAR

When the software is used for the first time, the TMS Settings workflow needs to be run first. It allows the selection of the coil type to be used in the stimulation session and its calibration. The workflow allows calibrating all coils that might be used for stimulation. It is necessary to run this workflow only if a new coil needs to be calibrated: the selection of the active coil is done within the main TMS Neuronavigation workflow.

The Recording Settings workflow needs to be run only **if** the visor2 system is provided with an amplifier. It allows defining the amplifier type, the number and the labels of the recording channels and the sampling rate. It is necessary to run this workflow the first time visor2 is used and each time the configuration of the recording channels needs to be changed.

Online Neuronavigation can only be started if all necessary patient data have been prepared. This is achieved within the Setup Patient workflow that allows to user to prepare all necessary 3D patient data for the navigation, starting from a patient's MR scan or using standard data. This workflow also allows defining stimulation targets for a subsequent Neuronavigation session.

The main application workflow, the TMS Neuronavigation can be started either through the corresponding button in the *Workflows* toolbar, or directly from the Setup Patient workflow, after it has been successfully completed. The TMS Neuronavigation workflow will guide the user to perform all steps that are necessary to start the navigation and to perform e.g., motor mapping in dedicated workflows. The Offline Review workflow can be started either through the corresponding button in the *Workflows* toolbar, or from the TMS Neuronavigation workflow, after this has been finished. Through optimized views this workflow allows the user to easily review and analyze all data that have been generated during a navigated TMS session.

## 3 SETUP PATIENT

Before starting a TMS session with visor2, an MR scan of the patient needs to be imported and segmented, in order to obtain a realistic head model of the scalp and of the brain. The brain surface will be then used to create the curvilinear surface, a mapping surface on which the 3D image texture is projected, which shows the cortical structures at different depths. In order to be able to co-register the MR scan with the head models, and later to register the patient with the corresponding MR scan, it is necessary to define the coordinate markers. All data will be copied to a patient-specific folder for later use in the Neuronavigation session.

Alternatively, if no patient MRI is available, standard data can be used. In this case, an MRI of a healthy adult together with a pre-calculated segmentation is copied to the patient data folder. *Note:* The QT version of visor2 supports only standard MRIs.

The image displays two screenshots of the 'Setup patient | Select MRI' dialog box in the visor2 software.

**Left Screenshot (Create New Patient):**

- Tab:** Create New Patient
- Patient Name:** John Doe
- MRI options:**
  - ☐ Use standard MRI
  - ☒ Import patient MRI
- MRI File:** C:\example\John Doe.nii
- Format:** Nifti
- Next Step:** Specify MRI Fiducials
- Footer:** Save Current Layout, Quit Workflow

**Right Screenshot (Modify Existing Patient):**

- Tab:** Modify Existing Patient
- Patient Name:** test2
- Patient has standard MRI:** (checked)
- MRI options:**
  - ☐ Use standard MRI
  - ☐ Import patient MRI
  - ☒ Do not change MRI
- Next Step:** Define Targets
- Footer:** Save Current Layout, Quit Workflow

In any case, the final step of the Setup Patient workflow allows you to pre-define targets you want to stimulate in a later Neuronavigation session. Targets can be chosen directly on the MRI or curvilinear surface or by referring to atlas coordinates.

### 3.1 CREATING A NEW PATIENT ENTRY

To prepare data of a new patient for visor2, select the option *Create New Patient*. Enter the patient's name in the text field provided. All patient related data including MRI, 3D models, TMS sessions, EMG recordings will be saved in a directory named according to what you type in here. If a patient with the same name already exists, a warning will be issued and you will be asked to enter a different name. In such a case, please add a unique identifier to the name of the patient, e.g. the middle initial or append the date of birth.

Next, you need to select if you want to work with patient specific MRI data or use standard MRI data, see section Load MRI. When all entries have been provided, press the [Specify MRI Fiducials] button to proceed in the workflow.

## 3.2 EDIT EXISTING PATIENT DATA

If you want to use another MRI for an existing patient or if you want to define new targets, choose the *Modify Existing Patient* option. Use the drop-down box to select the patient.

### 3.2.1 DEFINING NEW TARGETS

New targets can be defined based on the existing MRI by selecting the *Do not change MRI* option. Press the [Define Targets] button to arrive directly on the *Define Targets* page.

### 3.2.2 PROVIDE A NEW MRI

Select the *Import patient MRI* option to associate a new MRI data set with the patient. This allows you to import MRI data from file, see section Load MRI and to create a new head model, see section Create Model.

Please be reminded that targets and session data relate to the MRI that has been loaded when performing target definition and running a Neuronavigation session. Therefore, importing a new MRI renders these data invalid. Every time you import a new MRI, visor2 copies the old MRI together with the sessions performed using this MRI and the corresponding targets to a separate Archive folder, see the chapter on patient data management.

## 3.3 LOAD MRI

You can choose whether to associate patient-specific MRI data or standard MRI data with the patient. The third option *Do not change MRI* is only available for existing patients and allows you to define a new set of targets on previously imported MRI data.

### 3.3.1 STANDARD MRI

If you choose this option, an MRI of a healthy adult together with a pre-calculated segmentation is copied to the patient data folder. There is no need to define MRI Fiducials or run a custom segmentation; it has all been done for you already. When pressing the [Define Targets] button, the workflow proceeds directly to the *Define Targets* page (3.6).

We recommend using patient-specific MRI data whenever available. If you choose to work with standard data, please be reminded that

- Navigation accuracy might suffer from imperfect matching of the patient's head shape with the provided standard head
- Neuronal structures differ widely between patients; any anatomical data shown on the standard MRI will not match the true anatomy of the patient
- Mappings performed with the standard data will show the correct hot-spot relative to the registered scalp, but not relative to the standard brain's anatomical structure.

### 3.3.2 PATIENT-SPECIFIC MRI

To import a new MRI data set for the patient, you need to follow these steps:

- Select the type of your MR scan. You can choose between DICOM, Analyze, Nifti or the ASA format, if you already saved your original MR as ASA MRI file type.
- Click the *Browse* (  ) button. Indicate the MR scan you want to import. In case you selected DICOM as the image type, a special DICOM import dialog will appear. For details, see section 7
- *Review Exported DICOM*.
- Click the [Specify MRI Fiducials] button to proceed.

### 3.3.3 MRI REQUIREMENTS

The following restrictions apply to the MRI import:

|                          |                                                                                                                                                                                                                                                                                                                                             |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Supported formats</b> | DICOM; DICOM Secondary Capture Image; Nifti; Analyze; ASA                                                                                                                                                                                                                                                                                   |
| <b>Field of View</b>     | <ul style="list-style-type: none"> <li>• The entire head must be included: the nose (the tip should not be cropped), the ears, the forehead and the chin;</li> <li>• The headrest or the scanning table must not be included in the scan.</li> </ul>                                                                                        |
| <b>Image Properties</b>  | <ul style="list-style-type: none"> <li>• T1 weighted scans;</li> <li>• Maximally 512<sup>3</sup> voxel. If your MR scan contains more than 256 slices in any direction, this direction is automatically down-sampled to 256 slices.</li> <li>• Uniform spacing in each direction is recommended;</li> <li>• Uncompressed format.</li> </ul> |

Please note that any MRI data exported with visor2 in the DICOM Secondary Capture Image format (see chapter 7.1.3 Export to DICOM) is interpreted as monochrome data when imported again in visor2.

The use of Analyze as import data format for MRI scans is discouraged, as there is ambiguity about the image orientation. Please ask your radiology department to provide you with Nifti or DICOM images instead. If you need to work with Analyze images, it is most important to check the left-right orientation for every patient and data set.

On loading a MRI, the grey value data representation will be internally converted to unsigned integer, mapping the grey value minimum to zero. The bit resolution of the existing image is hereby preserved.  
Please note that when exporting the image to DICOM, the converted grey values of the pixel will be used and the format of the DICOM image will be unsigned integer.

### 3.4 SPECIFY MRI FIDUCIALS


FIGURE 6: THE MRI- BASED (LEFT) AND THE NASION-EAR BASED (RIGHT) COORDINATE SYSTEMS

Before setting the coordinate markers, all positions are expressed in the MRI-based coordinate system. The coordinates **of** the crosshair are visible either on the *Status Bar*, or on a tooltip, which becomes visible when dragging the crosshair in the *MRI View*.

The three coordinate markers for nasion and **the** left and right preauricular points are required to define a head-based coordinate frame for the MRI. Based on this coordinate frame, the MRI can be overlaid with the head model and other 3D boundary objects. When the nasion-ear coordinate markers are set, all positions are expressed in the Nasion-ear coordinate system.


Place the crosshair to the position where you want your nasion marker to be placed. This can be done by either dragging the crosshair directly in the *MRI View* or through the slice slider in the *MRI Toolbar*.

- Press the *Nasion Marker* button (  ). A green border indicates that the point has been defined.
- Repeat the procedure for the left and right preauricular points.
- Click the [Create Head Model] button to proceed.

FIGURE 7: CROSSHAIR INDICATING THE NASION

The MRI with the nasion-ear markers will be copied from its original location and stored in ASA file format in the patient directory (please see the "Patient Data Management" section).


- Only trained medical personnel are allowed to identify the positions of the coordinate markers.
- Please set the position of the coordinate markers very carefully. Remind that, during neuronavigation, the same points need to be identified on the patient, through the pointer tool.
- If different persons are responsible for the preparation of the 3D data and for the registration of the patient, the positions of the preauricular points need to be unambiguously defined


### 3.4.1 TALAIRACH MARKERS

The "Specify MRI Fiducials" page also allows you to set Talairach markers if you need to convert positions from and into this standardized head coordinate system during a later TMS session. For more information about the Talairach system and how to set the accompanying image markers, see Appendix 2: MRI coordinate systems.


1. Set the Nasion-Ear-markers as described above. A section that allows you to set the ACPC markers becomes visible (anterior and posterior commissure).

- Use the buttons to set the anterior commissure, posterior commissure and inter-hemispheric point markers: Place the crosshair in the 3d MRI view on the correct location and mark the spot by pressing the corresponding button. After all three points have been defined, the MRI view switches to the ACPC view automatically.

- Use the Talairach Markers buttons to set the boundary points of the brain, defining the limits of the Talairach coordinate space. AP: anterior point, PP: posterior point, SP: superior point, IP: inferior point, LP: left point, RP: right point. To set the two points for each direction, switch ACPC view orientation according to this direction by clicking on the orientation icons of the MRI slice toolbar.

Make sure you set the boundary points according to the ACPC view in the lower right of the 3D image view.

After having defined all markers, press the [Create Head Model] button. The Talairach markers are automatically saved with the MRI and are loaded when using the MRI for a TMS session.

### 3.5 CREATE MODEL

Segmentation is used to divide the MRI into the scalp and brain compartments. The segmented regions will serve as basis for generating the realistic head models of the different head compartments. During a TMS session with visor2, you will need a head model of the scalp and of the brain. The scalp head model is necessary to improve the registration accuracy by performing a so-called shape-based registration, where points digitized on the patient's head by means of a pointer are matched to the patient scalp. The head model of the brain serves as basis for extracting the curvilinear surfaces. These are concentric surfaces, the most outer one lying on the brain head model, which reflect the morphological structure of the brain. The estimated electric field produced by the coil is projected onto the curvilinear surfaces, allowing therefore viewing the distribution of the field at different cortical depths.

In the Segmentation workflow, you have the possibility to choose between manual brain segmentation and automatic segmentation of the brain based on the straightforward scalp segmentation. If you need an accurately segmented brain, we encourage you to perform manual brain segmentation. If the quality of the segmented brain is not crucial, the complete segmentation procedure can be completed very quickly by performing only the scalp segmentation and automatically creating the brain segmentation based on the segmented scalp. Figure 8 shows the result of a manual and automatic segmentation of the brain on the same MR scan.



FIGURE 8: LEFT - AUTOMATIC SEGMENTATION OF THE BRAIN FROM THE SCALP, RIGHT - MANUAL SEGMENTATION OF THE BRAIN

The scalp can be found relatively easy based on the contrast in grey value and so-called 3D region growing algorithms. For volumetric T1-weighted scans, the range of grey values for the scalp is between about 25 and 255 (although very much depending on the scan parameters. The threshold values used for the brain are usually in the range of 30 to 80, depending on the MRI. If the threshold values are inappropriate, the boundary may also contain other parts of the head, such as the scalp. This can in most cases be overcome by adjusting the threshold values. The segmentation of the brain is based on the morphological operation called opening, which breaks up thin connections between regions selected by the threshold criterion. In this way, the brain is separated from the rest of the head. **Since regions not belonging to the brain will be removed, it is very important to define the "seed point" within the brain properly before starting the segmentation**

### 3.5.1 MANUAL BRAIN SEGMENTATION

- Select the lower intensity threshold for the scalp segmentation. Voxels that will be classified as belonging to the scalp, according to the selected threshold values, are rendered in yellow. If the lower voxel intensity threshold is too low, small artifacts connected to the head and in the proximity of the head will appear (Figure 9).
- Try to increase the value of the lower threshold until the artifacts connected to the head disappear. Remind that the artifacts not connected to the head will be removed by performing the segmentation feature.

- Select the lower and the higher Intensity thresholds for the brain segmentation. Voxels that will be classified as belonging **to** the brain, according to the selected threshold values, are rendered in red. The red-rendered brain structure should be disconnected from the rest **of** the head. Increase the value **of** the lower threshold until a sufficient degree of separation between the brain and the remaining head structures is reached (Figure 10).


FIGURE 9: SEGMENTATION OF THE SCALP. LEFT: THE THRESHOLD VALUES ARE NOT CORRECT, AND SMALL ARTIFACTS ARE VISIBLE IN THE PROXIMITY OF THE HEAD. RIGHT: CORRECTLY DETERMINED THRESHOLD VALUES.


FIGURE 10: SEGMENTATION OF THE BRAIN. LEFT: THE LOWER THRESHOLD IS TOO LOW AND THE BRAIN STRUCTURE IS CONNECTED TO THE SCALP. RIGHT: CORRECTLY DETERMINED THRESHOLD VALUES, THE BRAIN IS SUFFICIENTLY SEPARATED FROM THE SCALP

- Select the number of curvilinear surfaces to be generated. The curvilinear surfaces are closed concentric surfaces generated from the brain head model by moving inwards by a 1mm step.
- Place the crosshair within the brain.

- Start the segmentation by clicking the [Create Segments] button. The segmented scalp and brain compartments will be displayed in the *MRI View*, the scalp on the left side, the brain on the right as indicated in Figure 11 .


FIGURE 11: *MRI VIEW* AFTER SUCCESSFULL SCALP AND BRAIN SEGMENTATION

- If the result of the segmentation is not satisfactory, readjust the thresholds, position the crosshair within the brain and perform the segmentation again.
- If necessary, the segmentation of the brain can be readjusted by performing a “dilation operation” or an “erosion operation” of the segmented compartment.
- Create and export the head models of the brain and the scalp by clicking the [Create Model] button. The head models will be displayed in the *MRI View*, as indicated in Figure 12.




FIGURE 12: SCALP AND BRAIN HEAD MODELS, AFTER MANUAL SEGMENTATION

## 3.5.2 AUTOMATIC BRAIN SEGMENTATION

- Select the lower intensity threshold for the scalp segmentation. Voxels that will be classified as belonging to the scalp, according to the selected threshold values, are rendered in yellow. If the lower voxel intensity threshold is too low, small artifacts connected to the head and in the proximity of the head will appear. Try to increase the value of the lower threshold until the artifacts connected to the head disappear. Remind that the artifacts not connected to the head will be removed by performing the segmentation feature.
- Select the number of curvilinear surfaces to be generated. The curvilinear surfaces are closed concentric surfaces generated from the brain head model by moving inwards by a 1mm step.
- Start the segmentation by clicking the [Create Segments] button. The segmented scalp and brain compartments will be displayed in the *MRI View*, the scalp on the left side, the brain on the right as indicated in Figure 11. Please remind that the rendered brain is only the starting scalp surface from which the brain will be automatically generated.

- Create and export the head models of the brain and the scalp by clicking the [Create Model] button. The head models will be displayed in the *MRI View*, as indicated in Figure 13.

After completing the segmentation, the MR image of type ASA (.mri), which includes the positions of the coordinate markers and the head models of the scalp and the brain will be exported in the patient directory.

Click the [Define Targets] button to proceed.

FIGURE 13: SCALP AND BRAIN HEAD MODELS, AFTER AUTOMATIC SEGMENTATION


- Only trained medical personnel is allowed perform the segmentation.
- The quality of the scalp segmentation directly influences the accuracy of the TMS navigation.
- Always check if the segmented scalp surface agrees with the scalp- air-boundary visible in the MRI.

### 3.6 DEFINE TARGETS

#### 3.6.1 SINGLE TARGETS

The Setup Patient workflow gives you the possibility to define targets on the brain anatomy. The defined targets will be saved and exported in the patient directory, and can be loaded within the TMS Neuronavigation workflow. During navigation, the visor2 software will help you targeting the identified markers, as explained in *Stimulate at Targets* section.

Place the crosshair at the desired target position.

This can be done in several ways:

- Directly drag the crosshair in the *MRI View* to the desired position.

- Use the slice sliders in the *MRI Toolbar*.
- Directly enter the coordinates through the *Go to MRI Coordinate* dialog under **Edit / Go to MRI Coordinates....** Depending on the coordinate markers that have been set, you will be able to specify the coordinates in different coordinate systems.


FIGURE 14: THE GO TO MRI COORDINATES DIALOG. NOTE THAT THE TALAIRACH SYSTEM IS ONLY AVAILABLE IF TALAIRACH MARKERS HAVE BEEN SET.

Once you are sure the crosshair is at the right position, press the [Add marker] button. A dialog will appear prompting you to enter the marker label. Enter the target name and Click OK. Repeat to add more targets, if necessary.

If you wish to start a neuronavigation session with the data you just generated, you should directly proceed to the Neuronavigation workflow, by clicking the [Proceed to Neuronavigation] button. Otherwise, terminate the workflow by pressing the [End Workflow] button.

### 3.6.2 TARGET GRID

As an advanced feature, visor2 allows you to generate multiple targets at once, arranged in a regular grid fashion. The target points from the grid are merged with the single targets created as described in 3.6.1 and stored together with them in a single target file.

To create a target grid, switch to the "Target Grid" tab. Set the cross hair at the desired mid-point of the grid. Adjust distance between neighboring grid points ("Spacing") and the extend of the grid ("Grid dimension"). Press [Show preview] to get a preliminary version of the grid presented in the render window. Note that the yellow disks do not yet represent target points, but only indicate where the target points will be created.

The grid is originally aligned with the x-y-plane of the nasion-ear-coordinate system. Use the "Alpha" and "Theta" angles to orient the grid tangential to the head and rotate it the way you prefer. Hereby, "Theta" represents the rotation of the grid around the x-axis of the nasion-ear-system and "Alpha" a rotation around the z-axis of the nasion-ear-system. The theta-rotation is performed before the alpha rotation.

After you have adjusted the grid preview to your needs, press the [Grid -> Target Markers] button. The list of existing targets will now be extended by the points of the grid. If required, you can define a second grid by running through the procedure again: Place the crosshair – press [Show preview] – adjust parameters – press [Grid -> Target Markers].


FIGURE 15 DEFINITION OF A TARGET GRID. ADJUST SIZE AND EXTEND WITH THE SPACING AND DIMENSION PARAMETERS. ADJUST ORIENTATION WITH THE ALPHA AND THETA ANGLES. PRESS GRID -> TARGET MARKERS TO CONVERT THE PRELIMINARY GRID POINTS TO REAL TARGET POINTS.

### 3.7 FILE STRUCTURE

The Setup Patient workflow allows you to load your MR scan as well as to generate and to export the brain and scalp head models. If all the steps are performed correctly, after completing the workflow, the following files will be generated:

- ASA MRI. Your original MR scan will be exported into the ASA .mri type. The ASA MRI file defines also the positions of the coordinate markers.
- The scalp and brain head models and the curvilinear surfaces, organized as follows:
  - A boundary definition file of the scalp and the brain (meg-x.bnd)
  - A boundary definition file, for each of the curvilinear surfaces (clrf-x.bnd).
  - A curvilinear surface definition file, pointing to curvilinear boundaries files (clrf.vol).
  - A head model definition file, pointing to the scalp and brain boundaries files (meg.vol).
- Target markers. A file containing the target coordinates in nasion-ear-coordinates and the corresponding target labels (.asc). The file name is concatenated of the patient name and the date and time when the targets were defined.


## 4 TMS SETTINGS

The TMS Settings workflow needs to be executed each time you need to calibrate a new coil or you connect the position sensor to a different USB port. It should also be the first step when starting visor2 for the first time. The TMS Settings workflow can be started by clicking the corresponding button in the workflows toolbar.

### 4.1 COIL SELECTION

The page allows you to define the coil you need to calibrate or recalibrate:

- If you have already defined this coil for use with visor2, select the entry “Modify Existing Coil” and choose your coil from the drop-down box.
- If you want to define a new coil for the use with visor2, select “Create New coil” and choose:
  - The manufacturer
  - The coil model
  - The calibration method: pointer based (Mag&More only) or using a calibration board
  - The calibration board that will be used later for calibrating the coil, if calibration method “Calibration board” is selected
  - An identifier for the coil. It is not allowed to use the same identifier for two different coils, even if their manufacturer or model differs.
  - The tracking tool. Choose NT-112 for the first coil in a TMS session and NT-143 if you want to use the coil as second coil during the simultaneous navigation of two TMS coils. If you operate visor2 with one coil only, choose NT-112 always.

Press the [Calibrate Tool Faces] button or operate the Play button on your remote control to proceed. You will now be prompted to calibrate the coil. Once you have successfully calibrated the coil, you can either terminate this workflow or you can go back to calibrate another coil.

## 4.2 BEFORE THE FIRST CALIBRATION OF A NEW COIL

Each coil that is used with the visor2 system for the first time has to be provided with a permanent marker on its upper face. Place one of the stickers that came with your visor2 system on the upper right part of the coil's upper face. The sticker serves as a control point to validate the accuracy of a calibration, see *Calibration procedure*. Do not worry about the exact location of the sticker; it is only important that it is easily accessible during the calibration procedure, see Figure 20.

For Mag&More coils, you do not need to attach the sticker. Instead, Point 'B' on the coil can be used as validation point.

## 4.3 CALIBRATION PROCEDURE

In order to be able to track the position of your coil precisely, vIsor2 needs to find out how the tracking tool is attached to the coil. During the calibration procedure, visor2 will calculate the position of the origin of the coordinate system attached to the coil with respect to the coil tracker. Moreover, visor2 will calculate the distance between the control point on the coil and the coil tracker. This distance will be used by vIsor2 as a discriminator of the validity of the coil calibration. During time in fact, although rigidly attached, it might happen that the coil tracker slides along the coil handle. The performed calibration is then no longer valid and the coil needs to be recalibrated. Summarizing, the calibration procedure is necessary:

- Each time a coil tracker is attached to the coil.
- Each time the coil validation procedure fails. During a TMS session with vIsor2, after having selected the stimulating coil, you will be asked to validate the calibration. Even though this step is not mandatory, we strongly encourage you doing so.

Depending on the calibration method selected in the "Coil Selection" step ("calibration points" or "calibration board"), the calibration procedure differs. Both procedures start with calibrating the tool faces of the wing tool, see section 4.4. The "calibration points" method requires then to record the position of the calibration points A, B, C, D on the coil using the pointer tool, see section 4.5. The "calibration board" method requires two additional steps: Recording the position of the validation point, see section 4.6, and calibration of the coil with the calibration board, see section 4.7.

## 4.4 STEP 1: CALIBRATE TOOL FACES

Before the actual calibration of the coil, visor2 checks and refines the geometry of the wing tool. This ensures that the tracking camera can always determine the position of the coil with maximum precision, independently of which wing tool face is presented to the camera.

Please attach the wing tool to your coil now, make sure it cannot slip on the coil handle and present the wing tool frontally to the camera. The indicator lights on the right show the visibility status of the left and right face of the wing tool, respectively. A green light indicates that the tool face is visible to the camera. A red light indicates that the camera cannot determine the position of the face: Either the face is not visible to the

camera or the camera's operation is impaired, e.g. by excess infrared ambient light. In this case, the source of the problem is indicated next to the tool's visibility indicator.

When both faces are visible to the camera, press [Calibrate] or operate the Play button on your remote control. If the calibration succeeds, you will be informed by a green text box and the workflow will automatically proceed to the second step. If the calibration fails, a red error text box will inform you about the source of the problem. Correct the problem and re-try the calibration.


FIGURE 16: CALIBRATION DIALOG, STEP 1 "CALIBRATE TOOL FACES"

#### 4.5 STEP 2 (POINTER): RECORD CALIBRATION POINTS

This is the second and last step of the calibration procedure using the Pointer to mark the calibration points A-B-C-D on a Mag&More coil. Perform the following actions to calibrate:

1. Position the coil such that the wing tool can be seen by the camera. The indicator on the right hand side of the program helps you to verify the visibility.
2. Put the pointer tip into the small hole at calibration point 'A' of the coil. Make sure that pointer and coil are visible and press the [Acquire] button. Alternatively, you can press the Play button on the remote control.
3. Acquire points 'B' to 'D' the same way as you recorded point 'A'.
4. When all points are acquired successfully, the [Calibrate] button will be enabled. Press it or operate the Play button on your remote control to calculate and store the calibration.

If you want to record a point again, press the [Undo Last] button or the Stop button on the remote control until the state of the point goes back from "recorded" to "unknown".

Make sure that you follow the order A-B-C-D as prompted on screen, otherwise a false calibration

will result!


FIGURE 17 COIL CALIBRATION USING THE CALIBRATION POINTS A-B-C-D

MESSAGE AFTER SUCCESSFUL CALIBRATION

#### 4.6 STEP 2 (CALIBRATION BOARD): CALIBRATE CONTROL POINT

This step will only appear if you calibrate using a calibration board. If you calibrate using the pointer, the software will automatically use point 'B' as the control point, see section 4.5.

In the second step, a reference point is defined which allows to check if the wing tool changed its position relative to the coil. Make sure that the control point sticker with the black cross on yellow ground is affixed to the coil. Hold the coil such that the camera can see the wing tool (it does not matter which face of the tool you use – after step 1 both tool faces are calibrated to each other). Hold the pointer such that the tip is on the center of the sticker's black cross and the reflective spheres are visible to the camera.

When ready, press [Calibrate] or operate the Play button of the remote control. If the calibration succeeds, you will be informed by a green text box and the workflow will automatically proceed to the third step. If the calibration fails, a red error text box will inform you about the source of the problem. Correct the problem and re-try the calibration.


FIGURE 18 CALIBRATION DIALOG, STEP 2 "CALIBRATE CONTROL POINT"

#### 4.7 STEP 3 (CALIBRATION BOARD): CALIBRATE COIL CENTER

This step is specific for calibration using a calibration board. For calibration using a pointer, see section 4.5.

In this last step, the position of the coil's center point ("hot-spot") relative to the coil's wing tool is calibrated. Keep in mind that the tracking camera can only measure the position of the wing tool; the position of the coil's hot-spot needs to be calculated from there using the coil center calibration.

To calibrate you will need the calibration board selected for this coil on the Coil Selection page. The calibration board's identifier and image are also presented on the Calibrate Coil Center page again. Please make sure that you use the correct calibration board; otherwise, the calibration results will be inaccurate.

Attach the coil on your calibration board. Place the coil in such a way that the center of the coil matches the point of origin on the board. The coil tracker attached to the coil might prevent from seeing whether the coil is properly centered. It is advisable to control the correct coil placement on the board by viewing it from the backside of the calibration board, as shown in Figure 19. Once you have placed the coil correctly, fix its position through the provided screws. *Note:* If your coil supports a tight fit special calibration board, you only need to place the coil in it and make sure it is properly fitted.


FIGURE 19: BACK VIEW OF THE CALIBRATION BOARD WINDOW WHEN CALIBRATING A COIL

Place the coil and the calibration board such that the camera can see the coil's wing tool and the calibration board, as showed in Figure 20. When ready, press [Calibrate] or operate the Play button of the remote control. If the calibration succeeds, you will be informed by a green text box. If the calibration fails, a red error text box will inform you about the source of the problem. Correct the problem and re-try the calibration.


FIGURE 20: INDICATE THE COIL CONTROL POINT WITH THE POINTER

After successful calibration, you can either return to the Coil Selection page to define and calibrate another coil – e.g. a coil with wing tool NT-143 for a neuronavigation session with two coils – or directly proceed to the Neuronavigation workflow.

## 5 EMG / EEG RECORDING

Note: This chapter applies to visor2 with EMG / EEG license only.

Note: Recording of EEG/EMG with visor2 is an optional feature.

### 5.1 RECORDING SETTINGS

Before you start a Neuronavigation workflow with EMG / EEG recording, you need to select the default montage folder for your amplifier. This needs to be done only once at the first start of the system and whenever you choose to use a different amplifier type. Select Tools → Options ... and then Montages → User Montages File Name. Browse for the montage file of your amplifier. All montages are stored in subfolders of C:\Program Files\Visor2\Data\Montages\visor\. The respective montages reside in:

|                                                          |                       |
|----------------------------------------------------------|-----------------------|
| eego ES-225 amplifier: up to 64 EEG + 6 bipolar channels | EE-225\EE-225.mtg     |
| eego ES-410 amplifier: up to 8 bipolar channels          | EE-410\EE-410.mtg     |
| Refa 72: up to 64 EEG + 4 bipolar + 4 aux channels       | Refa64\ Refa_64.mtg   |
| Refa 32: up to EEG + 4 bipolar + 4 aux channels          | Refa32\Refa_32.mtg    |
| Mobi6: up to 6 bipolar channels                          | MobiMini6\ Mobi_6.mtg |
| MobiMini: up to 2 bipolar channels                       | MobiMini2\ Mobi_2.mtg |

The Recording Settings workflow is dedicated to all amplifier related settings. The recording settings need to be configured before starting the first recording with the amplifier provided within the **visor2 system** and each time you wish to change the recording configuration i.e. the number and labels of the recording channels, the sampling rate and the mapping channels. The Recording Settings workflow can be started by clicking the corresponding button in the workflows toolbar.

The [Edit Montages] button allows to edit the existing montages or to add new ones. The montage defines the configuration of the recording channels i.e. their number and labels. The visor2 software already comes with a predefined set of montages that fit to your amplifier type. If you wish to edit them or add new ones, please see Appendix 1: Montages.


1. Select the amplifier type. The amplifier type determines the range of available sampling rates.
2. Select the sampling rate, depending on your application.
3. Select the recording montage. The recording montage indicates the configuration of the recording channels.
4. Select the mapping montage. The mapping montage indicates those channels that will be used to retrieve the mapping values: only these channels will be processed online, the peak-to-peak values of the corresponding responses will be calculated and projected onto the brain anatomy, as explained in *Finding the hot spot of more muscles*. Please remind that the mapping montage must be either equal to the recording montage or its subset.
5. Select whether to run the Impedance Check or not. This is only relevant in the case EEG is being recorded as well. If this option is checked, you will have the possibility to run the Impedance Check within the main TMS Neuronavigation workflow.
6. If you are using your amplifier for the first time on this PC, you need to calibrate its internal clock against the PC clock. Connect the amplifier, switch it on (if it is an amplifier with power switch) and proceed by pressing the [Calibrate Amplifier] button.

If you have used your amplifier before, you can finish the workflow by pressing the [Save and End Workflow] button.


## 5.2 AMPLIFIER CALIBRATION

This step needs to be performed whenever a new amplifier should be used for working with visor2. It synchronizes the clocks of the PC and the amplifier and allows thus to correlate precisely the recorded EMG signal with the recorded TMS coil position. The amplifier calibration is available as second step in the Recording Settings workflow. Note that you need to connect the amplifier and switch it on (if it is an amplifier with power switch).


After entering the “Calibrate amplifier” step, wait until streaming data are shown. Then press the [Start] button to start the calibration procedure.


We recommend letting the procedure at least for the time of an average TMS-EMG-session, e.g. 90 minutes. The suggested end time is shown on the workflow page. When the time is over, press the [Stop] button.

You receive a message that the calibration was successful. The amplifier can now be used for visor2 neuronavigation. Press the [End Workflow] button to leave the workflow.

If you have already performed a calibration, but you would like to repeat it, e.g. because you plan longer sessions now, just re-run the Calibrate Amplifier workflow step. After pressing [Start], you will be asked if you want to delete the existing calibration. Press the *Yes* button to proceed, press *No* to cancel.

After deleting the old calibration, you need to close and restart the Recording settings workflow. You can now re-calibrate your amplifier as described above.

## 6 NEURONAVIGATION

This section describes how to perform navigated TMS with visor2. It assumes that a valid head model of your patient/patient is already generated (chapter 3 “Setup Patient”) and coils for tracking are calibrated (chapter 4 “TMS Settings”). If you want to run a session with EMG recording (e.g. finding the motor threshold), you also need correct recording settings (chapter 5 “EMG / EEG Recording”).

### 6.1 CHOOSE COIL AND STIMULATOR

In this first step, you need to select the coil(s) that you would like to use for the neuronavigation session. Please select first the stimulator type from the drop down box on the left. Afterwards, select the calibrated coil that is connected to this stimulator and which you would like to use for the session. Please note that only calibrated coils will be available in the dialog.


FIGURE 21. SELECT STIMULATOR TYPE AND COIL FOR THE SESSION


FIGURE 22. OPTIONAL: ADD A SECOND COIL FOR THE NEURONAVIGATION SESSION

If you want to run a session with only one TMS coil, press the button [Validate Single Coil] to proceed. If you want to navigate two coils simultaneously, press the [Add Second Coil] button. In the latter case, the page refreshes and a dialog allows you to select the coil that will be used as a second coil in the session.

*Please note that only coils calibrated with the NT-112 or NT-123 coil tracker can be used as first coil. Coils calibrated with the NT-143 coil tracker can only be used as second coil in a session with two simultaneous coils.*

## 6.2 VALIDATE COIL

You are now prompted to validate your coil calibration. During the calibration procedure, visor2 calculates the distance between the control point indicated on the coil and the center of tracking tool attached to the coil. To ensure that the position of the coil tracker did not shift with respect to the coil, during the validation procedure visor2 measures the distance between the tracker and the control point. If this distance differs from the distance value saved during the calibration by more than 2mm, the validation will fail.

The validation of the coil calibration can be skipped. However, the accuracy of the position of the coil as displayed by visor2, and therefore the estimated position of the induced electric field maximum, depends entirely on the accuracy of the registration and the validity of the coil calibration. We therefore recommend you to validate the calibration of your coil in each TMS session.

*Note:* If you have calibrated a Mag&More coil using the calibration points A-B-C-D, you now need to select point B as control point.


FIGURE 23. COIL VALIDATION DIALOG. PLACE THE POINTER ON THE CONTROL POINT AND PRESS [VALIDATE]

- Place the pointer's tip in the center of the control point on the coil.
- While indicating the coil's control point, press the *Validate* button. Remind that during the validation procedure, the coil tracker and the pointer need to be seen by the position sensor. If any of these tools is not visible to the position sensor while pressing the *Validate* button, you will need to repeat this step.
- If the validation succeeds, the workflow automatically proceeds to the next step. Otherwise, you will see an error message. In this case, you can repeat the validation and try to place the pointer more precisely on the control point.

- If the validation fails again, the relative position of the coil tracker with respect to the coil changed and you need to calibrate your coil again. Quit the workflow and re-calibrate the coil as described in chapter 4 “TMS Settings”.
- If you run a session with two coils, you will first be prompted to calibrate the first coil. After a successful validation, you will be prompted to validate the second coil.

### 6.3 SELECT PATIENT

After you have ensured that the coil calibration is still valid, you are now ready to prepare the patient for the session. In this step, you should load the patient data you have previously prepared in the Setup Patient workflow.


- Select the patient. The drop down box offers you all patients from the patient root directory as explained in the section on **Patient Data Management**. If you run the Neuronavigation workflow right after the Setup Patient workflow, this field will already refer to the patient whose MRI has been processed.
- If patient setup was run with standard data, a reminder will be presented that standard MRI data is used for the Neuronavigation.

After you have selected the patient, choose the options to read stimuli from a previous session, to read target markers or to run a session with EMG / EEG recording (see below). Finally, press the [Define Coord. System] button to proceed to the next workflow step.

#### 6.3.1 READ STIMULI FROM PREVIOUS SESSION

The dialog allows you to optionally read data from a previous session for comparison and use as target information. The stimulus information can be used to exactly reproduce a previous stimulation, see section *Reproduce Stimuli*.

- Select the checkbox *Read stimuli from previous session*.
- Click the *Browse* button (...). A file browser opens.
- Browse for the session you would like to import and select the session's .pos.cnt file in the session's folder.
- The text displayed in the *read stimuli* text field refers to the date and time of the selected study. In the example below, 20140422\_0610 refers to a study started on April 22<sup>nd</sup>, 2014 at 06:10 o'clock.


*Hint:* By default, the file browser opens in the base directory of the selected user. To import a session associated with the current MRI, go to the respective Sessions\SessionDateTime folder. To import an archived session, please go to the Archive\ArchiveX\Sessions\SessionDateTime folder. Please be reminded that when importing a study associated with an archived MRI, the spatial agreement between old data and current MRI depends on the accuracy with which the nasion-ear-markers were set in both MRI.

### 6.3.2 READ TARGET MARKERS

This option allows you to read previously defined targets markers. These markers can be used during the Neuronavigation to exactly stimulate a pre-selected point, see section *Stimulate at Targets*.

- Select the checkbox *Read target markers*.
- Click the *Browse* button (...). A file browser opens.
- Browse for the target file you would like to import.
- The text displayed in the *read target markers* field refers to the patient name and the target definition time. In the example as shown:

20140422\_0610 refers to targets defined on April 22<sup>nd</sup>, 2014 at 06:10 o'clock.

*Hint:* By default, the file browser opens in the base directory of the selected user. To import targets associated with the current MRI, go to the respective Targets folder. To import archived targets, please go to the Archive\ArchiveX\Targets folder. Please be reminded that when importing targets associated with an archived MRI, the spatial agreement between old data and current MRI depends on the accuracy with which the nasion-ear-markers were set in both MRI.

### 6.3.3 SESSION WITH MOTOR MAPPING (EMG)

Options:

☒ Motor mapping (EMG)

FIGURE 24: SELECT OPTION FOR MOTOR MAPPING

If you want to run a session with EMG / EEG recording for motor mapping, tick the box "Motor mapping (EMG)". Depending on the amplifier type and recording settings (see chapter 5.1 "Recording Settings"), the next step will be either Impedance Check or Define Coordinate System.

Please note that using motor mapping may require the purchase of additional software and hardware components. Please contact your local visor2 dealer for more information.

### 6.4 DEFINE COORDINATE SYSTEM

The patient and his 3D representation need now to be registered. This means that the geometrical transformation mapping the physical space to the image space needs to be estimated. This can be done if the coordinates of three points, in this case the nasion and the left and right preauricular points, can be identified both on the MR scan and on the real patient.

In the Setup Patient Workflow, you have already set the nasion-ear coordinate markers on the MRI. The positions of these anatomical landmarks have to be determined now on the real patient, with the pointer tool.


FIGURE 25 RECORD FIDUCIAL POINTS TO DEFINE THE NASION-EAR-COORDINATE SYSTEM

- Identify the position of the nasion on the patient and place the pointer tip on it. Press the *Acquire* button or the *Play* button of the response device. If you accidentally digitize the wrong point, you can always remove the acquired position by pressing the *Remove last* button or pressing the *Stop* button of the response device.
- Repeat the procedure with both left and right preauricular points.
- If the pointer or the head reference tool could not be seen by the camera, an error message is shown. Make sure, both are visible and repeat the sampling.
- After having measured successfully the three positions of the nasion-ear-system, click the [Head Shape Points] button to proceed with refining the registration using sampled scalp points or press [Neuronavigation] or [Prepare mapping] to skip this step and work with the registration defined by the three landmarks.

During the digitization of the points, both the pointer and the head reference tool have to be seen by the tracking camera. During this step, the *Measurement View*, indicating the visibility status of all tracking tools, is displayed next to the Specify Head Coordinate System dialog. If a tracking tool is seen and recognized by the position sensor, it will be indicated in green color, otherwise, in red.

FIGURE 26: TOOL VISIBILITY STATUS IN THE MEASUREMENT VIEW

## 6.5 HEAD SHAPE POINTS GENERATION

When using only the three nasion-ear-markers, performing the transformation between the physical coordinate system to the MRI coordinate system can lack necessary accuracy. If the shape of the patient's head is known both in the physical and the image space, more sophisticated and accurate shape-based algorithms can be used to estimate the mapping transformation.

If MRI data for the patient was available, you have generated an exact head model of the patient in the Setup Patient workflow. The shape of the patient's scalp in the image space is therefore known. During this step, you should generate an approximation of the patient's head shape by **digitizing** as many points on the patient's head as possible. It is important not only to acquire a sufficient number of head shape points (a minimum of 20 points is required) but the area covered by the acquired points should be sufficiently wide. It is therefore important to acquire points not only over the top of the head, but also from the face, in the proximity of the eyes and the nose.


FIGURE 27 HEAD SHAPE POINT ACQUISITION. NOTE THAT THE FIT IS ONLY PERFORMED AFTER THE SAMPLING, SO IT IS NORMAL HERE TO SEE A DISCREPANCY BETWEEN SAMPLED POINTS AND DISPLAYED SCALP.

- Adjust the minimum distance between the points to be acquired (optional).
- If you check the checkbox "Show statistics", the shape-based estimate of the transformation will be displayed when proceeding to the next step.
- Take the pointer tool in the hand and place the tip on the patient's head.
- Press the [Start] button for recording head shape points while moving **continuously** with the tip of the pointer over the patient's head. The digitized head shape points are visualized online in the *HSP Acquisition View*, as shown in Figure 27.
- Note that no sampling will be performed if the pointer or the head reference are not visible to the camera.
- If you want to make a break in order to move the pointer away from the head's surface and resume the digitization in a different region of the head, you can press the [Pause] and resume the digitization by pressing the [Start] button again, once you have placed the pointer correctly.

When you have completed the digitization, press the [Fit] button. The alignment between head model and recorded head shape points will be calculated and applied to all navigation data and you will be forwarded to the next step.

Please note that while moving the pointer, its tip needs always to be in contact with the patient/patient's head.

Hint: To check the outcome of the registration, enable the view of the pointer ("show") in the 3D view properties and align the pointer with a number of feature points on the patient's face such as the nasion. The display should show the pointer at exactly the corresponding position on the head model.

## 6.6 NEURONAVIGATION

After performing the coil validation, you have completed all necessary steps to start the neuronavigation session. Depending whether or not you are performing motor mapping (EMG) and/or other mapping, you will be displayed different views and options.

Any surgical procedures should not be based solely on data acquired or recorded with visor2.

### 6.6.1 NEURONAVIGATION MAIN WORKFLOW

On the top left, you can choose the "active coil". The tasks you select in the central panel are performed for the active coil:

- Stimulate: The electric field is displayed only for the active coil. Note that stimuli are *recorded* for both coils all the time; independent of the active coil and the currently selected task.
- Reproduce stimuli: Navigates the active coil to the position where a previous stimulus was given.
- Stimulate at targets: Navigates the active coil to a position to best stimulate a predefined target.
- Replace coil: Allows you to "hot-swap" the active coil for another previously calibrated coil with the same coil tracking tool type.
- Validate: Check if the coil calibration for the active coil is still valid.

The other tasks in the central panel are independent of the active coil:

- Markers: Allows you to define additional targets.
- rTMS filter: activates grouping of incoming rTMS stimuli into one rTMS stimulus entity.
- EMG settings: Control the time window for the peak-to-peak extraction in the EMG recording after a stimulus.


FIGURE 28. NEURONAVIGATION STEP. VIEW FOR A SESSION WITH TWO TMS COILS AND EMG RECORDING.

## 6.6.2 NEURONAVIGATION WITH EMG RECORDING

The main view consists of the following sub views:

- The *MRI Online View*, displaying the real-time position of the coil respect to the patient/patient's brain.
- The *Average View*, displaying the last evoked muscle responses.
- The *Online Data View*, displaying the free running EMG signal.
- The *Tool Status View*, indicating the visibility status of the tracking tools.


FIGURE 29: NEURONAVIGATION WITH EMG RECORDING

The properties of each sub view can be changed through the View Properties dialog, accessible by double clicking the sub view or through the sub view context menu option "View Properties".

In the *MRI Online View*, a TMS coil is represented with a simplified sketch of its shape, as indicated in Figure 30, and with the projection along the coil axis of its center on the brain, more precisely, on the current curvilinear slice, called Coil target. The loops of the TMS coil are hereby visualized at the downside of the coil casing. The Coil target reflects the orientation of the coil, and it allows therefore determining whether the coil is held tangentially to the brain at the stimulation point.


FIGURE 30: 3D COIL REPRESENTATION. LEFT: SINGLE COIL SESSION. RIGHT: NAVIGATING TWO COILS SIMULTANEOUSLY.

When performing a TMS stimulus with your stimulator (usually via the foot switch or push button at the coil, or through a programmed interface) the visor2 system receives a trigger via the TTL connection cable and trigger unit and sets a stimulus marker at the coil position at this time point. The TMS Markers are shown as 3D lenses, as displayed in Figure 31. Their position corresponds to the projection of the coil's normal line to the curvilinear. Please note that if you select a different curvilinear slice, also the position of the TMS Markers change, as they are only projections of the stimulus position.


FIGURE 31: TMS MARKERS

- Always check EMG signal quality before starting to stimulate in a motor mapping session using the EMG online view.
- In case of bad signal quality, check proper electrode placement and the connection of the reference electrode.
- Any EMG traces recorded in the motor mapping workflow should be reviewed carefully and EMG responses with artifacts should be excluded from further use.
- Any surgical procedures should not be based solely on data acquired or recorded with visor2 including any mapping modality.

### 6.6.3 THE INDUCED ELECTRIC FIELD

During neuronavigation, while moving the coil over the patient's head, you will see the estimated electric field projected onto the cortical structures. The electric field induced by a TMS coil depends on:

- The coil position and orientation
- The coil geometrical structure and its physical properties
- The spatial distribution and electromagnetic properties of the head's tissue.

In visor2, a simplified model of the true physics of electric field generation is used. To allow for real-time display, the head is modeled as a spherical conductor with the center being at the origin of the nasion-ear-system. The coil is approximated by a superposition of magnetic dipoles, located inside the loops representing the coil.

The electric field is displayed on the currently activated curvilinear surface. The field intensity is color-coded by the field strength relative to the field maximum on the current curvilinear slice.


FIGURE 32: THE ELECTRIC FIELD INDUCED BY THE TMS COIL

Through the View Properties, you can adjust different properties of the electric field visualization, such as the threshold value, the color palette etc. You can also select to display always the position of the maximum estimated electric field and the corresponding field strength or to always snap the current slice position to the field maximum.

#### 6.6.4 REPRODUCE STIMULI AND STIMULATE AT TARGETS

The vIsor2 system assists you when you want to stimulate at predefined target positions or if you want to reproduce a previous stimulus. In the *Stimulation* panel, select “Reproduce stimuli” for navigating the coil to the exact position where a previous stimulus was given or “Stimulate at targets” to guide the coil to an optimal place to stimulate the marked spot, respectively.

After the selection of the task, the view layout will change. For “Reproduce stimuli” select the previous stimulus that you would like to reproduce in the TMS Stimulus List at the bottom right of the screen. For “Stimulate at targets”, select the respective target marker in the Image marker list at the bottom right screen.

##### Stimulate at Targets

1. In the “Stimulation” group box, select Stimulate at targets. The view layout will be automatically adjusted. A list view, displaying all available target markers, will be shown.
2. Adjust the Intensity slider to the stimulation intensity of your stimulator.
3. In the target list, select the target you want to use for stimulation. Markers that can be targeted are labeled with the type **Sensor**. The targeting helper will be displayed.
4. When the coil rings match the targeting crosshair, stimulate.


FIGURE 33. STIMULATE AT TARGETS. SELECT THE TARGET IN THE IMAGE MARKER LIST AT THE BOTTOM RIGHT.

### Reproduce Stimuli

1. In the "Stimulation" group box, select Reproduce stimuli. The view layout will be automatically adjusted. A list view, displaying all loaded and previously applied stimuli, will be shown.
2. Adjust the intensity slider to the stimulation intensity of your stimulator.
3. In the stimuli list, select the stimulus you want to reproduce. The targeting helper will be displayed.
4. When the coil rings match the targeting crosshair, stimulate.


FIGURE 34. REPRODUCE STIMULI. SELECT THE STIMULUS IN THE TMS STILULUS LIST AT THE BOTTOM.

For both tasks, target helper rings will guide you to position the coil accurately. As illustrated in Figure 35, the targeting rings indicate the position and the orientation of the target, whereas the crosshair of the targeting helper indicates the position and orientation of the coil, as displayed in Figure 36. When targeting, you should try to position the coil in such a way that the targeting helper and the rings match.


FIGURE 35: THE TARGETING RINGS

Your point of view is attached to the coil: while you move the coil, you will see the brain surface moving with respect to the coil, see Figure 37. As opposed to the *MRI Online View*, you cannot rotate the view, since the viewing perspective is fixed. For zooming in the view, you can still use the mouse wheel.


FIGURE 36: THE COIL RINGS INDICATE THE COIL POSITION AND ORIENTATION


FIGURE 37: DISPLACEMENT TARGET-COIL

The targeting bars help you understand in which direction you need to move the coil to reach the target. While the displayed distance refers to the absolute value, the length of each of three targeting bars represents the displacement target-coil in each of the three directions, as illustrated in Figure 37. The navigation helper rings and bars react immediately to your position changes and guide you to the right spot.

Visor2 displays two figures of merit for the targeting: a distance and an angle value. If you are *reproducing previous stimuli*, visor2 will calculate

- the distance as the difference in position
- the angle as the difference in orientation

of the coil between the previous stimulus orientation and the current one.

If you selected *Stimulate at targets*, visor2 will calculate

- the distance as the distance of the target marker to the normal line of the coil
- the angle as the angle between a line from the center of the head through the target marker and the normal line of the coil

In the View Properties, you can set the threshold values for the distance and the angle. Once the distance and the angle between the coil and the target fall below the defined thresholds, the values will be rendered in green, indicating you that the desired targeting precision is achieved, as displayed in Figure 34. This needs some exercise as the coil has to have the right position *and* orientation and both change while moving it, so do not be frustrated if the results are not perfect in the first session.

### 6.6.5 rTMS

If you are performing rTMS, you should apply the online rTMS filter and set the Inter-train interval according to the settings of your TMS Stimulator. Through this value, visor2 distinguishes between a train of stimuli in rTMS and single pulses applied one after the other. Every train of TMS pulses will be marked as an rTMS event, characterized by its duration, frequency and number of pulses. Use "Pause between pulse trains > 3.0s" if your Inter-train Interval is longer than 3 seconds or if you apply only 1 train. Use "Pause between pulse trains > 0.3s" otherwise. Note that Inter-train Intervals of less than 0.3 seconds are not supported. To switch the rTMS filter off, i.e. to treat every pulse as separate stimulus, select "No rTMS Filter".


FIGURE 38 RTMS FILTER SETTINGS

visor2 automatically displays a statistic about the rTMS pulses received. The TMS Stimulus list below the 3D render view displays the number of stimuli received in the rTMS train and calculates the rTMS frequency out of it.

| Stimulus Id | Stim intensity [%] | TMS Coil | Coil position [x,y,z in mm] | Coil orientation [qx,qy,qz,qw] | TMS stimulus time | rTMS stimulus | rTMS count | rTMS frequency [Hz] |
|-------------|--------------------|----------|-----------------------------|--------------------------------|-------------------|---------------|------------|---------------------|
| 1           | 94                 | 8        | 13, 67, 107                 | -0.382, -0.171, -0.389, 0.820  | 10:42:53.128      | Yes           | 17         | 6.4                 |
| 2           | 94                 | 8        | 47, 76, 96                  | -0.382, -0.171, -0.389, 0.820  | 10:43:42.707      | Yes           | 21         | 3.3                 |

FIGURE 39 RTMS STATISTICS IN THE TMS STILULUS LIST. THE TWO COLUMNS AT THE RIGHT GIVE THE NUMBER OF STIMULI IN A TRAIN AND THE PULSE FREQUENCY, RESPECTIVELY.

Further, the 3D render view can be adjusted to show how the coil moved during the rTMS train. In the 3D image view properties, go to the TMS Marker tab and activate the TMS coil positions -> Shading checkbox. The 3D visualization will now show small, connected dots that indicate the trace of the coil during the rTMS train. Note that this visualization is available during the session as well as in the offline review.


FIGURE 40 THE TMS COIL POSITION -> SHADING OPTION ALLOWS TO DISPLAY THE POSITIONS OF THE TMS COIL DURING A SINGLE RTMS TRAIN. THE POSITIONS ARE INDICATED AS RED DOTS CONNECTED WITH A BLACK LINE.

#### 6.6.6 ONLINE EMG SETTINGS

visor2 allows you to control the visualization of the evoked motor responses and the settings for the extraction of the response peak.

- The number of events is the number of evoked motor responses that will be displayed in the Average View in a waterfall plot. If you select "1", only the last evoked response will be displayed. If you select more, the single sweeps will be displayed in a waterfall plot, the most recent one being on the top.
- The display window is the width of the motor response that will be displayed in the Average View, relative to the time of the incoming TMS pulse.
- The peak search is the time interval where the peak-to-peak of the evoked motor response is searched for, relative to the time of the incoming TMS pulse.

Please remind to click the [Apply Settings] button after changing the online EMG settings.


FIGURE 41: THE ONLINE EMG SETTINGS

To stop the EMG recording, e.g. because you are going to stimulate away from the EMG hotspot, press the [Stop EMG Recording] button. This will close the EMG recording file and disable the amplifier. Please note that for technical reasons it is not possible to restart the recording afterwards!

#### 6.6.7 REPLACE COIL

If you have calibrated several coils, you can switch them during the TMS session. This can be useful if e.g. one of the coils overheats and you want to replace it with a fresh coil. Select the option "Replace Coil". In the drop-down box below, you can now select a coil from the list of all calibrated coils.


FIGURE 42 SWITCH COIL DURING A TMS SESSION

The coil you select in the drop-down box immediately becomes active. Make sure to validate your new coil to ensure correct calibration, see chapter *Validate coil*.

#### 6.6.8 VALIDATE


This method allows you to validate whether the displayed 3D information is still correct. This method helps you check whether:

- The head reference tool that the patient is wearing has shifted during the session and therefore decreased the accuracy of the displayed information.
- The coil tracker has shifted respect to the initial position on the coil and therefore decreased the accuracy of the displayed information.

If possible, validate the accuracy of the session multiple times during the session and every time there is a suspect that the patient moved the head reference tool or that the coil tracker shifted.


#### 6.6.8.1 VALIDATE NASION

- Go to the "Validate" tab, select "Validate Nasion".
- Hold the pointer tool on the nasion of the patient.
- Press the Play button on the remote control or the *Validate* button on the GUI.
- If the distance between the nasion position at the beginning of the session and the current nasion position is smaller than 5 millimeters, you will see the message "Check passed". Otherwise, you will get the message "Check failed".

#### 6.6.8.2 VALIDATE COIL

- Go to the "Validate" tab, select "Validate Active Coil".
- Hold the pointer tool on the coil's control point, as explained in *Validate coil*.
- Press the Play button on the remote control or the *Validate* button on the GUI.
- If the coil tracker has not shifted (threshold distance 5mm), you will see the message "Check passed". Otherwise, you will get the message "Check failed".

In case of a failing accuracy check, we strongly recommend to stop and restart the session. All results obtained so far should be treated with utmost care!

## 6.6.9 FINDING THE MOTOR HOT SPOT


FIGURE 43: TMS MARKERS COLOR CODED WITH A BLUE-WHITE-RED PALETTE

- Select the task "Stimulate".
- Adjust the Online EMG settings according to your preferences.
- When you stimulate, the TMS Marker will be color coded according to the peak-to-peak value of the evoked motor response. Every time you stimulate, the relative value of each TMS Marker will be recalculated, and the new color applied.
- If you want to change any of the properties of the TMS Markers like the mapping color palette, the display threshold or the size, you can do that through the TMS Marker Tab in the MRI Online View Properties dialog.

Please remind that the TMS Markers indicate only the projections of the coil center on the brain, color coded according to the intensity of the evoked response. Please be aware that the real cortical representation of the muscle being mapped might deviate from the location of the TMS Markers associated with the highest response amplitude.

## 6.6.10 FINDING THE HOT SPOT OF MORE MUSCLES

The visor2 system gives you the possibility to map more muscles at the same time, without having to run several neuronavigation sessions. If you are recording the EMG at several muscles, a TMS marker will be placed for each recording channel, and color-coded relatively to the maximum response of that muscle.

TMS Markers associated to different muscles are visually easily distinguishable, having a different pattern printed on the top face of the marker, as shown in Figure 44. visor2 can map simultaneously up to eight muscles at the same time.


FIGURE 44: TMS MARKERS FOR THE FIRST SIX EMG CHANNELS

When mapping several muscles, you can record continuously the EMG at all channels, or you can switch EMG channels ON and OFF during a session. A TMS Marker will be placed only for enabled channels.

If you record the EMG for all channels continuously, a TMS marker will be placed for each TMS pulse and for each EMG channel, as shown in Figure 45. Through the TMS Marker Tab in the MRI Online View Properties dialog, you can select the channels to be displayed and the channel specific display threshold.


FIGURE 45: MAPPING OF THREE MUSCLES. LEFT: SINGLE CHANNEL VIEW, FOR EACH CHANNEL. RIGHT: ALL CHANNELS, WITH A DISPLAY THRESHOLD SET TO ALLOW THE VISUALIZATION ONLY OF THE TMS MARKERS ASSOCIATED WITH THE LARGEST RESPONSE

To be able to switch the channels ON and OFF, please select the “Neuronavigation MRI Online Multichannel” view layout, through the *MRI View* context menu. The Channel List view will be displayed in the right bottom corner, as displayed in Figure 46.


FIGURE 46: MULTICHANNEL MAPPING VIEW LAYOUT

If you change the enabled EMG channels during a neuronavigation session, in the *MRI View*, per default, only the channel/s enabled during the last TMS pulse will be displayed. This can be easily changed through the *TMS Marker Tab* in the *MRI Online View Properties* dialog.

## 6.6.11 STIMULATING AT THE FOUND HOT SPOT

- Select the task “Reproduce stimuli”. The “NeuroNavigation Targeting Stimulus” view will be displayed, as in Figure 47. Per default, you will see only those TMS Markers with a response exceeding the 75% of the maximal evoked response. The TMS Markers will be displayed together with a label, indicating their identification number.
- Select the TMS Marker you want to reproduce from the TMS Marker list, displayed in the bottom, as shown in Figure 47. The identification numbers of the TMS Markers are also displayed in the TMS Marker list. The selected TMS Marker will be highlighted, as shown in the enlargement of Figure 47. The targeting helper will assist you in reproducing the selected stimulus, as explained in Reproduce stimuli section.
- Stimulate. The previously selected TMS Marker will remain the active target until you don't select a new one


FIGURE 47: THE VIEW WHEN TARGETING AT THE HOT SPOT

Please remind that, if you want to apply rTMS at the found hot spot, you need to stop the EMG recording, by clicking the **Stop EMG Recording** button.

## 6.7 TERMINATING THE SESSION

To terminate the session, click the *Finalize* button. You are now invited once again to check if the head reference of the patient has moved during the session. Therefore, hold the pointing device at the nasion of the patient and press the *Validate Nasion* button. If the distance between the nasion position at the beginning of the session and the current nasion position is smaller than 5 millimeters, you will see the message "Check passed". Otherwise, you will get the message "Check failed" together with an estimate of the shift of the head reference tool, informing you that since all positions are measured relative to the head reference tool, the estimated positions of the stimulated sites might not be accurate. Please consider this when reviewing this session. Please remind that while on this step, position information (and EMG if applicable) are still being recorded.


FIGURE 48. FINALIZING THE NEURONAVIGATION SESSION.

In case of a failing accuracy check we strongly recommend to treat all results obtained with utmost care. If in doubt, repeat the session or discard the results!

### 6.7.1 EXPORT TO DICOM

visor2 allows you to define image markers that can be exported to MR images of type DICOM. The image markers can be set after having analyzed the positions of the TMS Markers and the corresponding motor response peak values. visor2 allows you to define and export multiple image markers, which will look like white or colored squares (side length 4mm) in the exported DICOM MR image, as shown in Figure 49. Check the "Export to DICOM" checkbox. The "NeuroNavigation DICOM Export" view will be displayed, as in Figure 48. Per default, you will see only those TMS Markers with a response exceeding the 75% of the maximal evoked response. The TMS Markers will be displayed together with a label, indicating their identification number. Please check chapter 7.1.3 *Export to DICOM* for more details on the colored DICOM Export.

- Position the crosshair at the location where you wish your image marker to be placed. If you want to place the image marker over a TMS Marker, select the TMS Marker in the list. The position of the crosshair will be adjusted according to your selection, as displayed in Figure 50.
- Click the *Add Marker* button. A dialog will appear. Enter the label for the marker and click OK.
- The defined image marker should be visible in the *MRI View*.
- If you want to define and export more image markers before exporting to DICOM, repeat from step 2.
- Export to DICOM by clicking the *Export DICOM* button. The DICOM will be exported in the patient directory, as explained in the "Patient Data Management" section.


FIGURE 49: MARKERS EXPORTED IN DICOM

Depending on your requirements, you can define several image markers and export them all at once into one DICOM image, or you can export several DICOM images, each having a unique image marker.

Either you may now proceed to the Offline Review by pressing the *Proceed to Review* button or you can terminate the session by pressing the *Complete Session* button.


FIGURE 50: POSITIONING THE CROSSHAIR OVER A TMS MARKER

## 7 OFFLINE REVIEW

When performing a navigated TMS session, visor2 records the positions and orientations of all tracking tools. Each time a magnetic pulse is applied; visor2 records the position and the orientation of the coil in the moment of stimulation, the time, the duration and the set stimulation intensity. If the EMG (e.g. motor mapping) has also been recorded during the Neuronavigation session, visor2 records also the latencies and the amplitudes of the evoked motor responses.

You can choose between reviewing the recorded sessions or the exported DICOM images. When reviewing the sessions, you can explore the 3D spatial distribution of the stimuli, the induced electric fields and the evoked responses, if the EMG has been recorded when a motor mapping session was performed. When reviewing the exported DICOM images, you can explore the sites on the patient's original MRI that have been identified and marked.

### 7.1 REVIEW SESSION

#### 7.1.1 LOAD SESSION DATA

- Select the patient, as explained in section Select Patient.
- Select "Review TMS session" button.


FIGURE 51: OFFLINE REVIEW FIRST PAGE FOR SELECTED MOTOR MAPPING SESSION

- Choose the session to review
  - o The entry in the text box refers to the session's date and start time. In the figure above, "20141407\_2023" refers to a session from the 7<sup>th</sup> of April 2014 that was started at 08:23pm.
  - o If you have started the Offline Review workflow directly from a Neuronavigation session, the session field will already refer to the last session and you don't need to browse for it.
  - o Otherwise, click the *Browse* button (⋮). Sessions performed with the current MRI of the patient are stored in the folder called "Sessions" under their start date and time.
  - o A message below the text box shows if EMG ("Contains EMG File") recordings were performed during the selected session.

- Sessions performed with an earlier MRI are stored under the folder “Archive/ArchiveX/Sessions” in the directory of the selected patient. For more information on data storage please refer to the section on Patient Data Management.
- Click the [Review and Export] button to proceed.

The program will automatically load the MRI and head model associated with this session and present it together with the session data.

## 7.1.2 REVIEW SESSION

### 7.1.2.1 NEURONAVIGATION SESSION WITH MOTOR MAPPING (EMG)


FIGURE 52 OFFLINE REVIEW FOR NEURONAVIGATION SESSION WITH EMG

Per default, you will see only those TMS Markers with a response exceeding the 75% of the maximal evoked response. The TMS Markers will be displayed together with a label, indicating their identification number. Centrally at the bottom, the list of all events will be displayed, as shown in Figure 52. TMS Responses are labeled EMG\_X@channel where X is the stimulus number and channel is the label of the EMG channel. Therefore, EMG\_7@EMG1 refers to the 7<sup>th</sup> stimulus on channel EMG1.

The TMS Response view gives you detailed information about the stimulus response:

- Stimulus ID
- Channel label
- Peak to peak response (voltage difference between the global maximum and global minimum) in the search window defined during the Neuronavigation session
- Latency of first peak (can be maximum or minimum, depending on the montage applied), measured from the trigger recorded in the amplifier
- Latency of the second peak
- Time difference between the two latencies
- Stimulation intensities, i.e. the percentage of maximum stimulator output selected for the stimulus
- Label of the coil with which the stimulus was given
- Coil position and orientation in nasion-ear-coordinates
- Time of the stimulus

When a TMS Marker is selected in the list:

- The TMS Marker will be highlighted in the MRI View
- The crosshair will jump to the position of the selected TMS Marker;
- The electric field induced in the moment the TMS pulse was applied will be calculated and displayed;
- The coil, as positioned and oriented in the moment the pulse was applied, will be displayed;
- The free running EMG will jump to the time the pulse was applied.

#### 7.1.2.2 NEURONAVIGATION SESSION WITHOUT EMG


FIGURE 53: OFFLINE REVIEW FOR NEURONAVIGATION SESSION WITHOUT EMG

If a TMS session without EMG is reviewed, the view simplifies to display the stimuli with the stimulator output intensities, see Figure 53. The TMS Stimulus list presents:

- Stimulus ID
- Stimulation Intensities, i.e. the percentage of maximum stimulator output selected for the stimulus
- Label of the coil with which the stimulus was given
- Coil position and orientation in nasion-ear-coordinates
- Time of the stimulus

### 7.1.3 EXPORT TO DICOM

Before any export, please open the "Tools" -> "Options" -> "MRI Export" menu and choose if you want to use a monochrome or colored MRI (RGB) export by default. In this menu, you can also configure the default color for Coordinate markers and Sensors. The default color can later be changed during the export workflow. If the monochrome export option is selected, any coordinate marker or sensor (image markers) will be exported in white color.

**Please note that any TMS stimuli, TMS response (EMG) or TMS classification marker you want to export must be converted into a sensor first.**

To export a marker to DICOM, the following steps need to be performed:

1. Position the crosshair at the location where you wish your image marker to be placed. If you want to place the image marker over a TMS Marker, select the TMS Marker in the sensor, TMS stimuli or TMS response list view. The position of the crosshair will be adjusted according to your selection, as displayed in Figure 53.
2. Click the *Add Marker to MRI* button. A dialog will appear. Enter the label for the marker and select a color (only available with selected RGB MRI export option). Otherwise, the default color configured in the MRI export options will be used.
3. The defined image marker should be visible in the *MRI View*. If you want to define and export more image markers before exporting to DICOM, repeat from step 1.
4. Enable the image markers that you would like to export, disable the ones that you did not want to export to the DICOM file.
5. If you want to change the color of an image marker now, select the marker in the image marker list view and double-click on the shown color in the DICOM export color row. A new export color can now be selected.
6. Export to DICOM by clicking the *Export DICOM* button. The DICOM will be exported in the patient directory, as explained in the "Patient Data Management" section. The markers appear as a bright white cube of voxel size 2 in the exported MRI.
7. If a colored export was performed, the DICOM will be exported in the DICOM Secondary Capture Image format. Please refer to your manufacturer of your DICOM-compatible device if this format is supported.


FIGURE 54: ADD MARKER TO MRI MENU

## 7.2 REVIEW EXPORTED DICOM


FIGURE 55: SELECT PATIENT AND SESSION IN THE OFFLINE REVIEW

- Select the Patient for whom you want to review the exported DICOM data.
- Select the "Review exported DICOM" option.
- Click the [Review and Export] button to proceed. The "Open DICOM Images" dialog will be presented, and all available DICOM images of the selected patient will be shown, as displayed in Figure 56.
- Select the DICOM you want to review and click the [OK] button.
- The DICOM data will be displayed.

Remark: The "Open DICOM Images" dialog allows you also to browse for DICOM data anywhere on your PC:

- If you want to change patient, press the Browse button and select the proper patient directory. Make sure that the "Include Subdirectories In Search" checkbox is checked.
- Press the Refresh button (🔄). All available DICOM images of the selected patient will be displayed.
- Select the DICOM you want to review and click the OK button.


FIGURE 56: OPEN DICOM IMAGES DIALOG

FIGURE 57 THE EXPORTED IMAGE MARKER IS DISPLAYED AS WHITE SQUARE IN THE DICOM

## APPENDIX 1: MONTAGES

## THE MONTAGE EDITOR


FIGURE 58: MONTAGE EDITOR

A montage is an arrangement of EEG or EMG channels, together with their properties. The montage defines channel configuration, their labels and properties like sensitivity, trace color, and filter parameters. The Montage Editor (Figure 58) allows you to edit the already existing montages, to define new montage from scratch or by copying and editing an existing one. The Montage Editor can be accessed through the Recording Settings workflow, by clicking the corresponding button.

The Montage Editor presents four groups of controls, as shown in Figure 59.

To create a new montage, go to the upper left part of the Montage Editor dialog and click *New*. The *New Montage Dialog* appears. In this dialog, the new montage can be labeled, e.g. "EMG 3 Channels" and its storage location defined. If you want to create a montage for data recording, you should check the "Recording Montage" option.

Click OK. The montage information is now displayed in the upper right part of the Montage Editor.

Go to the lower left corner and unfold the *Reference Menu* to select the reference electrode. The average reference is chosen with the option "AVG".

Select the channel or channels of interest in the *Input Sensors* list and then click *Add to output*. In the lower right part of the *Montage Editor*, the selected channels and their reference will be listed. We suggest you select all the channels.


FIGURE 59: MONTAGE EDITOR GROUPS

Click OK to close the dialog and apply the new montage. Once the Montage Editor has been closed, the defined montage will appear in the Recording and Mapping montages list.

If you wish to create a montage similar to an existing one, it is more convenient to create a copy of the existing montage and edit it, rather than defining a new one. This can be done selecting the montage to be copied, clicking the *Copy* button and giving a proper name to the new montage. The newly defined montage will have the same properties and channels as the original one.

The montage allows you to define directly filtering parameters, which will be applied online, once the recording is starting, without need to set them manually. **Please remind that these are only display filters, meaning that the data will be recorded raw, without any filters apply.** Filtering parameters can be set at different levels:

- Montage level: applies to all channels and all channel groups. The montage filtering parameters can be set by clicking the *Edit* button in the left upper part of the Montage Editor after having selected the proper montage.
- Channel level: applies only to the channel. The channel filter parameters overwrite the montage filtering parameters, if any are present. The channel filtering parameters are set in the Montage output space, by selecting the desired channel and unfolding the desired filtering parameter, as shown in Figure 60.
- You have to keep in mind that all the options specified in a montage, including filtering options, will be applied automatically every time the montage is selected overriding any specification done through the toolbars and menus. In other words, there is a hierarchy of view properties:
  - Channel properties
  - Settings in the EEG toolbar.



FIGURE 60: CHANNEL FILTERING PROPERTIES

If a property has been defined in the EEG toolbar (lowest level) and **if** you then apply a montage that specifies this property as well, you will override the specifications made through the toolbar. Due to this property **of** montages, you can review recorded data with different view properties, e.g. filter settings, or you can review a data file with two different montages synchronously.


## APPENDIX 2: MRI COORDINATE SYSTEMS

### 1. Original MRI system

This system reflects the way the MRI data are stored. There are two variants: voxel based and distance based, which distinguish according to the unit used (voxels or millimeters). The origin is always at the corner of the MRI block that represents the first voxel of the first slice. The x axis is defined by the direction the slices are stacked, the y axis corresponds to the slowest changing index within a slice. The visor2 slice coordinates relate to the voxel based coordinate system.

voxel based original MRI system


distance based original MRI system (unit millimeters)

### 2. Fiducial system (Nasion-Ear-System)

This system is defined by so-called fiducial points on the surface of the head. The fiducials used are the nasion and left and right pre-auricular points. The y direction is determined by the connecting line between the two pre-auricular points (pointing left). The origin is then found by projecting the nasion orthogonal onto this line. The x axis is pointing from the origin to the nasion and the z axis is standing perpendicular on the plane spanned by x and y axis.

**Note:** The orientation of the coordinate axes can be changed in the application options. Select the menu option Tools / Options... and choose the tab 3D Coordinates.



Nasion-Ear-system


### **3. AC-PC system**

The AC-PC system is an intermediate step towards the Talairach system.

The origin is determined by the anterior commissure (AC), a small fiber bundle connecting the two hemispheres of the brain. A similar fiber bundle, the posterior commissure (PC) determines the orientation of the x axis, which points into anterior direction. The z axis is oriented upwards and perpendicular to the AC-PC line between the two hemispheres. It lies in a plane defined by the AC-PC line and the interhemispherical point (IHP), which can be chosen at an arbitrary spot in the interhemispherical fissure, provided it does not lie on the AC-PC line itself. The y axis points to the left side, orthogonal to the other two axes.

The AC and PC points are shown in an axial plane (left) and in a sagittal plane (right). The AC is the anterior, more frontal, the PC, the posterior point, more to the back of the head.

Definition **of** the AC-PC coordinate system.

## 5. Talairach system

The Talairach system is created by a piecewise linear transformation **of** the AC-PC system. It transforms the brain to a standardized size. As reference points one uses **the** maximum extension of the cortex along the directions defined by the AC-PC system. The following points are used<sup>1</sup>:


<sup>1</sup> Definitions according to the coordinate axis definition, where x points to the front and y to the **left**.

| Name of point   | Abbreviation | Description                                            |
|-----------------|--------------|--------------------------------------------------------|
| Anterior point  | AP           | point of the cortex with maximum x coordinate          |
| Posterior point | PP           | point of the cortex with maximum negative x coordinate |
| Superior point  | SP           | point of the cortex with maximum z coordinate          |
| Inferior point  | IP           | point of the cortex with maximum negative z coordinate |
| Right point     | RP           | point of the cortex with maximum negative y coordinate |
| Left point      | LP           | point of the cortex with maximum y coordinate          |


These points, as well as the posterior commissure (PC) get new coordinates assigned, according to which the brain is shrunk or stretched. These are:

|    | x / mm | y / mm | z / mm |
|----|--------|--------|--------|
| PC | -23    | 0      | 0      |
| AP | 70     | -      | -      |
| PP | -102   | -      | -      |
| IP | -      | -      | -42    |
| SP | -      | -      | 74     |
| RP | -      | -68    | -      |

In the following figures it is shown how the MRI looks like in vIsor2 after setting each of the Talairach markers (note that the AP marker has been shifted to make the difference to the Nasion-Ear-system more obvious):

MRI after having set the **AP Talairach** marker


MRI after having set the **PP Talairach** marker


MRI after having set the **SP Talairach** marker

MRI after having set the **IP Talairach** marker

MRI after having set **the RP Talairach marker**

MRI after having set the **LP Talairach marker**

Now, all Talairach have been set. The MRI coordinates are now given also in Talairach coordinates. To verify, whether you have set the Talairach markers correctly, scroll through the slices in all three directions. The brain should lie within the red box.

## APPENDIX 3: LICENSING

**visor2 software** is protected from unauthorized use by the implementation of a software licensing technology. A valid license is required to run the software. A license may consist of several options, which are enabled based on the license products that were purchased. The options control access to specific features of the software.

### LICENSE TYPES

The following license types are available

The **single license** unlocks the software on a single computer. The license is provided as a USB dongle (USB key-lock, see Figure 61).

The **single software license** unlocks the software on a single computer. The license is provided as a registration key

The **trial license** provides 30 days of access to all functionality of the software, and is intended for evaluation and demonstration purposes. The trial license is included with each copy of the software.

FIGURE 61: USB DONGLE

### INSTALLATION OF A SINGLE LICENSE

The single license is provided as a USB dongle. The dongle requires a hardware driver for the Windows operating system. The driver is provided by the CodeMeter software, which is installed automatically during the installation.

Installation steps:

1. Keep the USB dongle ready, but do not plug it into the computer yet.
2. Install **visor2 software** on the target computer, this will install the CodeMeter driver automatically. Do not run the **visor2 software** after the installation.
3. The CodeMeter icon now appears in gray color in the Windows System Tray (bottom-right of the monitor screen, embedded in the Windows Taskbar). The gray color indicates that no dongle has been connected to the computer.
4. Plug the USB dongle into the computer. Windows recognizes the dongle as new hardware. The CodeMeter icon now appears in green color. The green color indicates that one dongle has been connected to the computer.
5. The license is now correctly installed.

### TRIAL LICENSE

A trial license is included on every new installation of the **visor2 software**. The trial license grants access to all license options for an evaluation period of 30 days.

Installation steps:

1. Run the software.
2. The trial license window appears. Press the OK button to start the trial license.
3. After 30 days, the license will expire. Contact your local distributor, or the manufacturer for purchasing an official license.

## LICENSING SUPPORT

Contact your local distributor for further support regarding the software licensing.

For specific questions about a dongle, please remember to communicate the serial number, which is located on the side of the dongle.

## APPENDIX 4: SUPPORTED DATA FORMATS

The following table presents supported data formats for EEG/EMG, EEG/EMG data export (no additional manipulations are needed to work with these formats).

The following abbreviations will be used:

- PAT --- Patient information
- EVT --- Events
- MNT --- Montages

**Formats for EEG/EMG data exporting**

| Data format                            | EEG/EMG | PAT | EVT | MNT |
|----------------------------------------|---------|-----|-----|-----|
| ASA (.msr)                             | X       |     |     | X   |
| European Data Format EDF+ (.edf, .rec) | X       | X   | X   |     |
| ASCII                                  | X       |     | X   |     |
| EEProbe (.cnt)                         | X       |     | X   |     |

## APPENDIX 5: WARRANTY, SUPPORT AND RECLAMATIONS

### WARRANTY

The authorized representative **of** the manufacturer in Russia Neurosoft LLC guarantees the systems quality conformance to the technical documentation requirements of the system if the rules of operation, storage, transportation and mounting prescribed in the user documentation are observed and the system is maintained and regularly serviced.

Warranty generally covers 24 months from the day of delivery by Neurosoft to the customer. The delivery date is the date of waybill or other shipment document for the system.

The warranty period **of** components exposed to wear is 30 days.

The warranty does not cover the following:

- Lack of identification labels of the system parts or damage to them.
- Nonobservance of operation rules and conditions, storage, transportation, arrangement and maintenance.
- Expiration of a guarantee period.
- Seal integrity damage without permission of Neurosoft LLC.

Only the authorized representative of the manufacturer can provide repair service.

Neurosoft LLC is obliged to repair the system free of charge in case of failure within the warranty period. Repair service is provided in the Neurosoft Service Center (5, Voronin str., Ivanovo, 153032, Russia) according to the procedure stated in the "Reclamations" section.

The warranty period is extended for a period from filing of reclamation to the repair completion (see the "Reclamation" section).

### SUPPORT

For support please contact the Neurosoft company ([info@neurosoft.com](mailto:info@neurosoft.com)).

## RECLAMATIONS

Please send written notification that includes the following information:

- the customer's name, address and email
- the serial number of the system
- the detailed description of failures. If possible indicate the reasons and circumstances preceding the fault detection (in addition it is recommended to add the test report, the exam data, photos and other materials allowing to solve the problem as soon as possible).

In case of system return to the service center for repair or replacement, the following rules should be observed:

- it must be cleaned and disinfected before sending to service center. Read rules on cleaning and disinfection.
- it must be packed so to exclude the possibility of its damage during the transportation.
- the notice and this manual must be added to the system being returned.

## APPENDIX 6: SHORTCUTS AND MOUSE WHEEL FUNCTIONS

| Shortcut               | Meaning                                                            |
|------------------------|--------------------------------------------------------------------|
| <b>Right</b>           | Go to the next page (segment, event, selected object)              |
| <b>Left</b>            | Go to the previous page (segment, event, selected object)          |
| <b>Up</b>              | Increase sensitivity                                               |
| <b>Down</b>            | Decrease sensitivity                                               |
| <b>End</b>             | Go to the end of the EEG/EMG                                       |
| <b>ENTER</b>           | Start / Stop automatic scrolling                                   |
| <b>+ / -</b>           | Increase / Decrease automatic scrolling velocity                   |
| <b>Control + Right</b> | Go to next second (next sample if page <b>width</b> is too small)  |
| <b>Control + Left</b>  | Go to previous second (previous sample if page width is too small) |
| <b>Control + Up</b>    | Increase page width (less mm / sec)                                |
| <b>Control + Down</b>  | Decrease page width (more mm / sec)                                |
| <b>Control + C</b>     | Copy the view as bitmap to clipboard                               |
| <b>Control + A</b>     | Select all channels and events                                     |
| <b>Control + P</b>     | Print the view                                                     |
| <b>Mouse wheel</b>     | Zoom in 3d view                                                    |

## APPENDIX 7: SYSTEM COMPONENTS AND GUIDANCE RELATED TO EMC

The compliance in the tables below is valid for visor2 systems consisting of:

| ID                                        | Component Name                                                                   | Compliance                                                                                                    |
|-------------------------------------------|----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| M00150                                    | 3D tracking camera<br>(position sensor)                                          | IEC 60601-1<br>IEC 60601-1-2                                                                                  |
| M00262                                    | power supply camera                                                              |                                                                                                               |
| M00261                                    | USB converter camera                                                             |                                                                                                               |
| EE-410                                    | eego amplifier (USB)                                                             | Class IIa, MDD 93/42/EEC                                                                                      |
| 666601                                    | visor2 software                                                                  | Class IIa, MDD 93/42/EEC                                                                                      |
| M01646                                    | isolation transformer                                                            | IEC 60601-1<br>IEC 60601-1-2                                                                                  |
| NH-134                                    | trigger box                                                                      | IEC 60601-1-2 (EMC)                                                                                           |
| M01592                                    | visor2 system compartment                                                        |                                                                                                               |
| M01596                                    | visor2 cart                                                                      |                                                                                                               |
| NT-116                                    | head tracker tool                                                                | Non-medical, conforms to Essential<br>Requirements of Annex I of the<br>Medical Device Directive<br>93/42/EEC |
| NT-115                                    | pointer tool                                                                     | Non-medical, conforms to Essential<br>Requirements of Annex I of the<br>Medical Device Directive<br>93/42/EEC |
| NT-112                                    | coil tracker tool                                                                | Non-medical, conforms to Essential<br>Requirements of Annex I of the<br>Medical Device Directive<br>93/42/EEC |
| M01623                                    | power cord AC OUT2 Display                                                       |                                                                                                               |
| M01593                                    | power cord AC OUT3 Computer                                                      |                                                                                                               |
| M00582.1                                  | USB 2 cable, DC IN 5V, 500mA, camera converter                                   |                                                                                                               |
| M01160                                    | USB 2 cable, DC IN 5V, 500mA, amplifier                                          |                                                                                                               |
| M00582.2                                  | USB 2 cable, DC IN 5V, 500mA, trigger box                                        |                                                                                                               |
| XS-121                                    | BNC cable, Trigger IN1, SIP, 5V 20mA                                             |                                                                                                               |
| XS-294                                    | BNC adapter to double MCX cable, Trigger, SIP, 5V,<br>20mA                       |                                                                                                               |
| XE-292, XE-293, XE-294,<br>XE-295, XE-296 | Bipolar adapters for EE-410 eego amplifier (differ in<br>the number of channels) | Class I, MDD 93/42/EEC                                                                                        |

with following accessories:

| ID     | Accessory Name    | Requirements / Compliances                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      | Manufacturer | Model           |
|--------|-------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|-----------------|
| M00325 | Computer          | <ul style="list-style-type: none"> <li>- EN 60950, EN 62311</li> <li>- Power consumption: &lt;100W</li> <li>- Intel Core Duo 3 GHz or better</li> <li>- 2GB RAM working memory</li> <li>- Video card: DirectX 9.0 Pixel/Vertex shader 2.0 or better (Nvidia or ATI recommended)</li> <li>- VGA interface port</li> <li>- 6 x USB 2.0 ports.</li> <li>- MS Windows 7 or MS Windows 8.1 or Windows 10 operating system in "Pro" variant in 32 bit or 64 mode</li> <li>- Internet Explorer 11 or higher</li> </ul> | Bluechip     | S3300 PC        |
| M00345 | Monitor / Display | <ul style="list-style-type: none"> <li>- EN 60950, EN 62311</li> <li>- Power consumption: &lt;80W</li> <li>- VGA Interface</li> <li>- minimum screen resolution: 1920x1200</li> <li>- minimum size: 24inch</li> </ul>                                                                                                                                                                                                                                                                                           | Ilyama       | E2483HS         |
| M00634 | keyboard          | <ul style="list-style-type: none"> <li>- CE (MDD not required)</li> <li>- English layout</li> <li>- integrated touch pad.</li> <li>- IP22 or higher</li> </ul>                                                                                                                                                                                                                                                                                                                                                  | Active Key   | AK-C4400F-GUS-W |

and port definitions for power:

| ID     | Description                  | Type                                                        | Cable length | Maximum current [A]  | Nominal voltage [V] | Minimal voltage [V] | Maximal voltage [V] |
|--------|------------------------------|-------------------------------------------------------------|--------------|----------------------|---------------------|---------------------|---------------------|
| XG-117 | AC IN, 230 V, 330W           | input, single phase AC-mains, 50/60Hz, not detachable cord  | <2.9m        | 1.3<br>(1.6 at 115V) | NA                  | NA                  | 230                 |
| M01623 | AC OUT, 230 V, 80W Display   | output, single phase AC-mains, 50/60Hz, not detachable cord | <2.9m        | 1.3<br>(1.6 at 115V) | NA                  | NA                  | 230                 |
| M01593 | AC OUT, 230 V, 100W Computer | output, single phase AC-mains, 50/60Hz, not detachable cord | <2.9m        | 1.3<br>(1.6 at 115V) | NA                  | NA                  | 230                 |

and signal ports:

| ID                                     | Description                           | Type                                                 | Cable length | shielded | wireless/<br>RF | number of ports     | service only |
|----------------------------------------|---------------------------------------|------------------------------------------------------|--------------|----------|-----------------|---------------------|--------------|
| M00582.1                               | camera position sensor<br>DC 5V 500mA | USB2, not detachable cord, type A USB connector      | <1m          | yes      | NA              | 1                   | yes          |
| M01160                                 | eego amplifier<br>DC 5V 500mA         | USB2, not detachable cord, type A USB connector      | <1m          | yes      | NA              | 1                   | yes          |
| M00582.2                               | trigger box<br>DC 5V 500mA            | USB2, not detachable cord, type A USB connector      | <1m          | yes      | NA              | 1                   | yes          |
| XS-121                                 | Trigger IN1, TTL, 5V, 20mA            | coaxial, not detachable cord, BNC connector          | <2.9m        | yes      | NA              | 1                   | yes          |
| XE-292, XE-293, XE-294, XE-295, XE-296 | Electrode clip-on lead wires          | Patient lead, not detachable cord, snap on connector | <2.9m        | yes      | NA              | 1-8                 | no           |
| ES600-SIOP-0010                        | tracking volume (field of view)       | infrared, free space                                 | NA           | NA       | infrared        | tracking of 6 tools | no           |

The following guidance and declarations are given with regard to EMC:

| Guidance and manufacturer's declaration – electromagnetic emissions                                                                                                                                  |            |                                                                                                                                                                                                                                                                                                                                                                                                               |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>The visor2 system is intended for use in the electromagnetic environment specified below. The customer and/or user of the visor2 system should ensure that it is used in such an environment.</i> |            |                                                                                                                                                                                                                                                                                                                                                                                                               |
| Emission Test Compliance                                                                                                                                                                             | Compliance | Electromagnetic Environment - guidance                                                                                                                                                                                                                                                                                                                                                                        |
| RF Emissions by CISPR 11 (GOST R 51318.11-2006)                                                                                                                                                      | Group 1    | The visor2 system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.<br><br>The visor2 system is suitable for use in all establishments, including domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF Emissions by CISPR 11 (GOST R 51318.11-2006)                                                                                                                                                      | Class B    |                                                                                                                                                                                                                                                                                                                                                                                                               |
| Harmonic Emissions by IEC 61000-3-2 (GOST 30804.3.2-2013)                                                                                                                                            | Class A    |                                                                                                                                                                                                                                                                                                                                                                                                               |
| Voltage Fluctuations/ Flicker Emissions by IEC 61000-3-3 (GOST 30804.3.3-2013)                                                                                                                       | NA         |                                                                                                                                                                                                                                                                                                                                                                                                               |

| Guidance and manufacturer's declaration – electromagnetic immunity                                                            |                                                                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <i>The visor2 system has been tested as compliant according to standards as indicated below.</i>                              |                                                                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Immunity Test Standard                                                                                                        | Compliance test level                                                    | Electromagnetic Environment - guidance                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2 (GOST 30804.4.2-2013)                                                          | ± 2,4,8,15 kV air<br>± 8 kV contact                                      | Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.                                                                                                                                                                                                                                                                                                                                                                                   |
| Electrical fast transients<br>IEC 61000-4-4 (GOST 30804.4.4-2013)                                                             | ± 2 kV<br>100 kHz repetition frequency                                   | IEC 61000-4-4 The quality of electrical energy should be ensured in accordance with a typical commercial or hospital environment.                                                                                                                                                                                                                                                                                                                                                                                                 |
| Surge<br>IEC 61000-4-5 (GOST 51317.4.5-99)                                                                                    | line to line ± 0,5 kV, ± 1 kV<br>line to ground ± 0,5 kV, ± 1 kV, ± 2 kV | IEC 61000-4-5 The quality of electrical energy should be ensured in accordance with a typical commercial or hospital environment.                                                                                                                                                                                                                                                                                                                                                                                                 |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 (GOST 30804.4.11-2013) | At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°                         | IEC 61000-4-11 The quality of electrical energy should be ensured in accordance with a typical commercial or hospital environment.<br><br>It is recommended to perform system power supply from a battery or an uninterruptible power supply unit if it is necessary to provide continuous running during possible voltage interruption.                                                                                                                                                                                          |
| Power frequency (50/60 Hz) magnetic field<br>IEC 61000-4-8 (GOST 50648-94)                                                    | NA                                                                       | no magnetic sensors or related components included                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Conducted RF<br>IEC 61000-4-6 (GOST 51317.4.6-99)                                                                             | 150 kHz to 80 MHz / 3 Vrms outside ISM band, 6 Vrms in the ISM band      | cable length below 3m                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Radiated RF electromagnetic fields<br>IEC 61000-4-3 (GOST 30804.4.3-2013)                                                     | 3 V/m<br>80 MHz to 2.7 GHz                                               | $d = 1.17\sqrt{P}$ (from 80 to 800 MHz)<br>$d = 2.33\sqrt{P}$ (From 800 MHz to 2.5 GHz)<br>d – recommended space diversity, m;<br>P – rated maximum output of the transmitter set by the manufacturer, W.<br>The field strength in the radio-wave propagation from fixed radio transmitters according to the results of observations over the electromagnetic environment <sup>1)</sup> should be lower than the level of compliance in each frequency band <sup>2)</sup> .<br>There is noise influence near the equipment marked |
| Radiated RF electromagnetic disturbances<br>IEC 61000-4-3 (GOST 30804.4.3-2013)                                               | Band (MHz)                                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                               | 380-390                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                               | 430-470                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                               | 704-787                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                               | 800-960                                                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                               | 1700-1990                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                               | 2400-2570                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                               | 5100-5800                                                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |

|  |  |  |                                                                                              |
|--|--|--|----------------------------------------------------------------------------------------------|
|  |  |  | the  label. |
|--|--|--|----------------------------------------------------------------------------------------------|

<sup>1)</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the system is used exceeds the applicable RF compliance level above, the system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the system.

<sup>2)</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Notes:

1. At 80 MHz and 800 MHz, the higher frequency range applies.
2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

#### Recommended separation distances between portable and mobile RF communications equipment and system

The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the digital system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the digital system as recommended below, according to the maximum output power of the communications equipment.

| Rated maximum output power of transmitter, W | Separation distance according to frequency of transmitter, m    |                                         |                                          |
|----------------------------------------------|-----------------------------------------------------------------|-----------------------------------------|------------------------------------------|
|                                              | 150 kHz up to 80 MHz<br>outside ISM bands<br>$d = 1.17\sqrt{P}$ | 80 MHz to 800 MHz<br>$d = 1.17\sqrt{P}$ | 800 MHz to 2.5 GHz<br>$d = 2.23\sqrt{P}$ |
| 0.01                                         | 0.117                                                           | 0.117                                   | 0.233                                    |
| 0.1                                          | 0.369                                                           | 0.369                                   | 0.738                                    |
| 1                                            | 1.167                                                           | 1.167                                   | 2.333                                    |
| 10                                           | 3.689                                                           | 3.689                                   | 7.379                                    |
| 100                                          | 11.667                                                          | 11.667                                  | 23.333                                   |

Notes:

1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.
3. For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

## APPENDIX 8: MEASURED, DISPLAYED AND ESSENTIAL PARAMETERS

The following general and environmental specifications and parameters hold for the visor2 system.

| Parameter Name                             | Description                       | Unit           | Nominal | Min | Max |
|--------------------------------------------|-----------------------------------|----------------|---------|-----|-----|
| Mode of operation                          | Continuous, no standby            | NA             | NA      | NA  | NA  |
| Maximum power consumption                  |                                   | W              | NA      | NA  | 330 |
| Standby power consumption                  | no standby modus                  | W              | NA      | NA  | NA  |
| Highest Data rates / Frequencies           |                                   | GHz            | NA      | NA  | 1   |
| Installation type                          | transportable (on cart)           | NA             | NA      | NA  | NA  |
| Maximum weight                             |                                   | kg             | NA      | NA  | 100 |
| Maximum volume                             |                                   | m <sup>3</sup> | NA      | NA  | 2   |
| IP class                                   | IP20                              | NA             | NA      | NA  | NA  |
| Operating temperature                      |                                   | °C             | 25      | 10  | 40  |
| Storage temperature                        |                                   | °C             | NA      | 10  | 40  |
| Transport temperature                      |                                   | °C             | NA      | 0   | 50  |
| Humidity                                   |                                   | %              | NA      | 15  | 70  |
| Atmospheric pressure                       |                                   | kPa            | NA      | 70  | 106 |
| Electrical Safety class (system)           | I                                 | NA             | NA      | NA  | NA  |
| Laser class                                | only low power infrared, no laser | NA             | NA      | NA  | NA  |
| Radio/Wireless technologies used in System | no radio/wireless included        | NA             | NA      | NA  | NA  |
| Expected Service Life                      |                                   | years          | 5       | NA  | NA  |

The measured and displayed parameters are:

| Parameter                                               | Description                                           | Unit     | Nominal value                | Minimal Value | Maximal Value |
|---------------------------------------------------------|-------------------------------------------------------|----------|------------------------------|---------------|---------------|
| Maximum mapping stimuli count                           | Maximum number of mapping stimuli                     | NA       | NA                           | NA            | 1000          |
| supported rTMS frequencies                              | any value up to the specified maximum                 | Hz       | NA                           | NA            | 20            |
| Maximum rTMS stimulus count at 1 Hz                     | Maximum number of stimuli                             | NA       | NA                           | NA            | 4000          |
| Maximum rTMS stimulus count at 5 Hz                     | Maximum number of stimuli                             | NA       | NA                           | NA            | 3000          |
| Maximum rTMS stimulus count at 20 Hz                    | Maximum number of stimuli                             | NA       | NA                           | NA            | 2000          |
| EMG potential range, adjustable (eego amplifier EE-410) | refer to user documentation of eego amplifier EE-4xx  | mVpp     | adjustable (150-1000)        | 150           | 1000          |
| EMG potential resolution                                | refer to user documentation of eego amplifier EE-4xx  | bit      | 24                           | NA            | NA            |
| EMG potential error                                     | refer to user documentation of eego amplifier EE-4xx  | %        | NA                           | NA            | 10%           |
| EMG acquisition sampling rate                           | refer to user documentation of eego amplifier EE-4xx  | Hz       | Adjustable (512, 1024, 2048) | 512           | 2048          |
| Position accuracy (Tracking camera)                     | refer to user documentation of position sensor camera | mm (rms) | 0.25                         | NA            | NA            |
| Image acquisition rate                                  | refer to user documentation of position sensor camera | Hz       | NA                           | NA            | 20            |
| Coil validation error                                   | depends on operator                                   | mm       | NA                           | 1mm           | NA            |
| Nasion validation error                                 | depends on operator                                   | mm       | NA                           | 1mm           | NA            |

| Parameter                                               | Description         | Unit  | Nominal value | Minimal Value   | Maximal Value |
|---------------------------------------------------------|---------------------|-------|---------------|-----------------|---------------|
| Coil distance to target /<br>Reproduce stimuli accuracy | depends on operator | mm, ° | NA            | 0.25mm,<br>0.6° | NA            |

The essential performance of the system is defined in the following table:

| Essential Performance Parameter                                                              | Applicable during EMC testing | Description                                                                                                                                                              | Unit  | Nominal value | Minimal Value | Maximal Value |
|----------------------------------------------------------------------------------------------|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------|---------------|---------------|---------------|
| Over all accuracy of coil position relative to the head (absolute value)                     | yes                           | Worst case upper limit of the position mismatch in positioning the coil relative to the head.                                                                            | mm    | NA            | NA            | 10.0          |
| Maximum timing jitter between TMS on set trigger and camera position stream (absolute value) | no                            | Upper limit of the time mismatch between time of stimulus (as TTL pulse) and time when the position of the tracking tools is recorded and associated with that stimulus. | ms    | NA            | NA            | 500.00        |
| Maximum timing jitter between TMS on set trigger and EMG stream (absolute value)             | no                            | Upper limit of the time mismatch between time of stimulus (as TTL pulse) and EMG data recorded.                                                                          | ms    | NA            | NA            | 5             |
| EMG waveform accuracy and dynamic range                                                      | no                            | Input triangular signal at 100Hz with 100 mVpp and +/- 300 mVpp offset is correctly reproduces with 10% accuracy (in time and amplitude).                                | %     | NA            | NA            | 10            |
| EMG input noise                                                                              | yes                           | < 5µVrms for a bandwidth from 40Hz to 200Hz. Notch filter for 50Hz and 60Hz can be applied. Disturbance clearly different from EMG is acceptable.                        | µVrms | NA            | NA            | 5             |
| EMG frequency response                                                                       | no                            | Low pass bandwidth up to 500Hz. Tested with sinusoidal signal of identical amplitude at 40Hz 200Hz and 500Hz shall be reproduced with 10% accuracy.                      | %     | NA            | NA            | 10            |

## APPENDIX 9: EXPLANATION OF SYMBOLS AS USED IN THIS MANUAL AND/OR ON SYSTEM PARTS

| Symbol                                                                              | Meaning                                      | System Components                                                                                                   |
|-------------------------------------------------------------------------------------|----------------------------------------------|---------------------------------------------------------------------------------------------------------------------|
|    | Manual contains important safety information | Position sensor, Host USB converter, tracking tools, amplifier, amplifier power supply, USB stick                   |
|    | Type BF                                      | Amplifier, amplifier power supply                                                                                   |
|    | Consult Instructions for use                 | Tracking tools, amplifier, remote control, trigger box, USB stick                                                   |
|    | Manufacturer                                 | Amplifier, USB stick                                                                                                |
|    | Manufactured date                            | Amplifier, USB stick                                                                                                |
|    | CE marking                                   | Amplifier, USB stick                                                                                                |
|    | Recycle where possible or return             | Position sensor, amplifier, tracking tools, amplifier power supply, amplifier USB unit, remote control, trigger box |
|    | Do not push                                  | vIsor2 cart                                                                                                         |
|    | On (power: connection to mains supply)       | Position sensor, host USB converter                                                                                 |
|    | Status                                       | Position sensor                                                                                                     |
|   | Error                                        | Position sensor, host USB converter                                                                                 |
|  | Connection port                              | Position sensor                                                                                                     |
|  | Direct current                               | Host USB converter                                                                                                  |
|  | USB port                                     | Host USB converter                                                                                                  |
|  | Bluetooth                                    | Mobi amplifier                                                                                                      |

## APPENDIX 10: MAINTENANCE OF THE SYSTEM AND DISINFECTION

The safety measures during maintenance conform to the ones described in sections 1.6 "Safety Warnings" and 1.7 "Safety Measures".

The system assembly and setup should be performed by the authorized representative of the Neurosoft Company. It should be understood that assembling accuracy of the device specifies safety and work quality.

The servicing of bought articles included in the system is performed according to user and technical manuals or standard rules.

If any trouble is detected, can't be eliminated using the control units or by restart, it should be switched off and checked by the specialist.

Maintenance category, amount and frequency are not specified intentionally except mentioned in this section.

The delivery set should match the packing list of the medical equipment.

The user maintenance of the system includes the external examination, check of connectors and cables, decontamination from the surface before each use according to recommendations stated in section 4.2.1 "Disinfection".

### Disinfection

The system (including reusable accessories: a pointer, a reference tracker consisting of head reference tool, clip, mounting tool for head reference as well as coils) is to be disinfected after using according to the methodological guidelines MG-287-113 of disinfection, presterilizing clearing and medical device sterilization (approved by the Department for State Sanitary and Epidemiological Supervision of the Ministry of Health of the Russian Federation on 30 December, 1998).

To avoid setting of organic pollutants by disinfectants on the surface and in the channels of the medical device it is necessary to pre-clean it in compliance with anti-epidemic measures.

The following disinfectants are recommended to use:

- Triosept
- Avansept-spray
- Ultratside
- 0,5% water alcohol solution of chlorhexidine bigluconate

Disinfection with these solutions is recommended to perform according to the methodological guidelines of a specific solution.

Clearing and disinfection of the device is carried out after it is disconnected from the power supply.

Cleaning and disinfection of the device by immersion (in a washing solution, organic solvents, etc.) is not allowed.

Avoid the ingress of liquids into connectors and manufacturing holes of the electronic units.

After cleaning or disinfection wipe the surface of the medical device dry with a cotton cloth.

Disinfection of coils after using is made in accordance with the recommendations set out in the operational documentation for coils.

## APPENDIX 11: PACKING, TRANSPORTATION AND STORAGE

The package should conform to the accepted one when manufacturing and delivering. In case the factory package is damaged, but the long-term system storage and transportation is expected, follow these recommendations:

- The system **with** operational documentation should be packed in plastic sachets and cardboard boxes.
- The cardboard boxes should be covered by the paper tape or pressure sensitive adhesive.

### «VISOR2» TMS Neuronavigation system

**Manufacturer:**  
eemagine Medical Imaging Solutions GmbH  
Gubener Str. 47  
D-10243 Berlin, Germany

RC №.  
Lifetime:

**The authorized representative of the manufacturer: Neurosoft LLC**  
5, Voronin str., Ivanovo, 153032, Russia  
Phone : +7 4932 24-04-80, Fax: +7 4932 24-04-35

**Year of manufacture: 20XX**

The label sample on the package

The system can be transported by all kinds of covered carries, besides unheated plane spaces, according to the rules of goods transportation for each mode of transport.

Transportation by sea should be carried out in accordance with the "Safety Requirements of General Cargos Shipping". Dispatch type: containers and small shipment.

### Transportation:

The system should be transported in a manufacturer package under the following environmental conditions: ambient temperature is from 0 to +50°C and relative humidity is from 17 to 70%. The air should be free from any ad-mixtures which can cause the corrosion.

### Storage:

The system should be stored in a manufacturer package in an enclosed space under the following environmental conditions: ambient temperature is from +10 to +40°C and relative humidity is from 15 to 70% without condensate water. The air should be free from any ad-mixtures which can cause the corrosion.

If long term storage of the system is expected, perform its conservation.

### Conservation:

The system components including accessories and operational documentation should be packed in separate plastic sachets and then placed in a manufacturer package.

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