Instruction for use

T-SEAL MOBILE polymer tube sealer with accessories

Reference:

Date: April 25th, 2019

Version: 2.0

Signature: *DJS Johansson*

Position: *Quality & RA Manager*

Name: *DEBORAH JOHANSSON*

Box 5031
SE-194 05 Upplands Väsby
Sweden

Signature: *DJS Johansson*

I, the undersigned, Jens Ericsson Rindö
Ass. Notary Public of Solna, Sweden,
Hereby certify that *Deborah Johansson*
personally has signed this document
Solna, Sweden. *7/5 2019*
Ex officio.

Jens Ericsson Rindö

Jens Ericsson Rindö

37742
21.06.2019

Contents

1. NAME OF MEDICAL DEVICE.....	3
2. MANUFACTURER'S NAME AND ADDRESS	3
3. DESIGNATION	3
4. INDICATIONS FOR USE, DOMAIN OF USAGE.....	4
5. CONTRAINDICATIONS	4
6. ADVERSE EFFECTS	4
7. WARNINGS AND SAFETY PRECAUTIONS.....	4
8. DESCRIPTION OF DEVICE AND ACCESSORIES.....	6
9. BASIC SPECIFICATIONS	7
+/- 10%.....	8
10. INFORMATION ABOUT THE PHARMACEUTICALS/ANIMAL OR HUMAN MATERIALS IN THE MEDICAL DEVICE	13
11. METHODS OF CLEANING AND DISINFECTION.....	13
12. PRELIMINARY STARTING PROCEDURE.....	14
13. WORKING SEQUENCE	16
14. ELECTROMAGNETIC COMPATIBILITY ACCORDING TO IEC 60601-1-2.....	23
15. STORAGE AND TRANSPORT CONDITIONS.....	23
16. APPLICATION CONDITIONS	23
17. THE LIST OF APPLICABLE STANDARDS	24
18. LABEL: DESCRIPTION OF SYMBOLS.....	24
19. INFORMATION ABOUT RECYCLING OR DISPOSAL. HAZARD CATEGORY OF MEDICAL WASTE	25
20. LIFETIME	26
21. WARRANTY STATEMENT.....	26

1. NAME OF MEDICAL DEVICE

T-SEAL MOBILE polymer tube sealer with accessories (hereinafter – the T-SEAL MOBILE)

I. Composition:

1. Manual sealer gun.
2. Battery (not more than 10 pcs. – If required).
3. Charger (not more than 10 pcs. – If required).
4. Five adapter kit for charger (not more than 10 pcs. – If required).

II. Accessories:

1. Manual gun mount (not more than 10 pcs.).
2. Extension cord (not more than 10 pcs.).

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Authorized representative in Russia

Terumo Russia LLC
Testovskaya Str. 10, 123112, Moscow, Russia, 13-th floor, premise I, room 5
Tel: +7-495-988-4740
Fax: +7-495-988-4739
E-mail: anna.sabirova@terumobct.com
E-mail: quality.EMEA@terumobct.com
www.terumobct.com

3. DESIGNATION

The product is designed to seal polymeric blood bag tubes and the polymeric tubes included in the Blood Packs or the Apheresis Disposable Sets.

4. INDICATIONS FOR USE, DOMAIN OF USAGE

The product is designed to seal polymeric blood bag tubes and the polymeric tubes included in the Blood Packs or the Apheresis Disposable Sets.

The T-SEAL MOBILE is an automatic sealing system for polymeric (PVC- and EVA-) tubes connected to blood bags included either in Blood Packs or in Apheresis Disposable Sets. The sealing can be performed when the donor is still connected to the Blood Pack or the Apheresis Disposable Set.

The T-SEAL MOBILE works with radio frequency (RF) energy to generate heat for sealing.

The T-SEAL MOBILE is used in the same environment as medical equipment (hospitals and blood banks).

Potential consumers of medical product:

The device is to be used by highly qualified personnel only.

5. CONTRAINDICATIONS

Not applicable, there are no known contraindications.

6. ADVERSE EFFECTS

Not applicable, there are no known possible side effects.

7. WARNINGS AND SAFETY PRECAUTIONS

The general safety information in this manual is for the personnel who operate the product. Specific notes, cautions and warnings are found throughout the manual where applicable. Please, read the Instructions For Use carefully before use.

NOTE

It is used to identify the conditions an attention should be drawn to.

CAUTION

It is used to identify the conditions that could result in damage to the equipment.

WARNING

It is used to identify the conditions that could result in personal injury or loss of life.

WARNINGS

The T-SEAL MOBILE must be used in full compliance with all the specifications and operational procedures listed in this manual.

Only trained personnel should be allowed to use the T-SEAL MOBILE Cables and accessories, other than those specified, may result in increased emission or decreased immunity of the equipment or System. Only accessories designed for use with T-SEAL MOBILE T-SEAL MOBILE.

The equipment or system should not be used adjacent to, or stacked with, other equipment. Otherwise, the equipment or system should be monitored constantly to verify normal operation in the configuration in which it will be used. If any of the components of the T-SEAL MOBILE are exposed to blood, they must be cleaned with an appropriate disinfectant solution immediately.

The T-SEAL MOBILE generates RF energy during the sealing procedure and during movement of the electrode. The T-SEAL MOBILE is not to be used in an oxygen-saturated environment.

The T-SEAL MOBILE is not designed to be used with flammable anesthetics or other flammable materials. When using the T-SEAL MOBILE, medical personnel should wear protective clothing and gloves.

CAUTION

ELECTROMAGNETIC INTERFERENCE REGULATIONS

This equipment fulfils EN 60601-1-2:2007 Standards for Electromagnetic

Compatibility. Nevertheless, it should be remembered that this equipment uses radio-frequency (RF) energy to generate tube sealing heat, hence it can affect other Medical Electrical Equipment. See table 1 for guidance. If installation and use are not in accordance with the IFU, there could be interference with radio/TV equipment and instrumentation.

WARNINGS

An operating principle of the T-SEAL MOBILE is based on the use of radio frequency (RF) energy to generate sealing heat. Users should be cautious of potential electrical shocks or hazards while handling this sealer. Always keep your fingers away from the electrodes in the slot. Never place any object other than the tube between the electrodes.

This machine emits a low level of electromagnetic (non-ionizing) radiation while sealing. It should not be used near high frequency sensitive electronic equipment.

In cases of usage of the medical device with the donor still connected to the system:

Do not allow the Hand unit to come directly in contact with the donor.

Do not perform a seal within 8 cm (3 in) of needle to preclude an RF burn at the needle entry point. burn at the needle entry point.

Inspect all parts of the instrument for defects before use. Check sealing pattern if Hand unit is dropped.

CAUTIONS

Always check the seals when the battery capacity is low (20% or less).

8. DESCRIPTION OF DEVICE AND ACCESSORIES

Description of the medical device

The T-SEAL MOBILE is supplied with the Sealing Hand Gun, the Battery, the Charger, and the Five Adapter Kit for the Charger. The Instruction for Use is supplied.

The Hand Gun consists of an RF-generator with intelligent control system, ergonomic handgrip with trigger, moving electrode (which can be removed for cleaning), and cable with connector. The Hand Gun consumes no power between sealing cycles.

The Battery consists of environmental friendly LiIonMn cells, built-in optic charge level indicator, and combined charger/gun cable connector. The charge in the Battery unit is sufficient for approximately 1500 seals.

Description of Sealing

The tube to be sealed is placed in the interelectrode slot of the Hand Gun. When the user presses the trigger on the Hand unit, the tube is pressed together and the sealing procedure will commence automatically. The radio frequency (RF) generator starts and the energy is transferred from the fixed electrode to the tube, which melts to a sterilized welding pattern.

During the whole sealing process, the orange light is lit on the front of the Hand Gun. When this light goes out, the sealing is complete and the trigger may be released. The Hand Gun's intelligent control system is used to detect, control, and adjust the sealing parameters to ensure optimal seal quality, depending on a tube type used. See section 7.4 for a description of the different LED indications on the Hand Unit in case of problems.

Composition

T-SEAL MOBILE comprises the parts listed below:

Description	Hand Unit	Battery	Battery charger without plug	Kit primary adaptors
Part Number	546900000	547900000	544900000	54410120

List of all cables and maximum lengths of cables, transducers and other electrical accessories with which the manufacturer claims compliance (cables, accessories used other than those listed may impact emission/immunity)

Part	REF	Specification
Battery charger without plug	8METSM003	Input: 100-240 VAC/50-60 Hz/700 mA Output: 18VDC/1.66A
Kit primary adaptors	8METSM037	Euro mains connector UK mains connector Australia mains connector USA/Japan mains connector IEC 320 mains connector

Accessories

Description	Holster	Extension cord 1.2 m 5-pol
REF	8METSM004	8METSM005

9. BASIC SPECIFICATIONS

Type polymeric tubes (PVC (EVA) tube)	Different tube types and sizes of up to 6.2 mm (PVC) and 5.4 mm (EVA) in outer diameter are sealed under the control of sophisticated system that adapts sealing time automatically.
Sealing capacity	Approximately 1500 seals/charge with PVC tubes up to 5 mm outer diameter at 20°C (68°F)
Continuous seal capacity	50 seals
Mode of operation	Operation: 25%, Intermittent: 75%
Sealing time	0.4 up to 3 sec. depending on tubing
Input Power	29.6 V DC 1.5 Ah LilonMn Battery
Temperature	Operating: 15 to 35°C (60 to 95°F) Storage and transportation: -20 to 70 °C (-4 to 158 °F)
Humidity	10 to 90% RH (non condensing)

	Storage and transportation: 10 to 90% (non-condensing)
Altitude	Operating: max. 3000 m (9842 feet)
Classification (IEC 60601-1)	Class II ME EQUIPMENT, externally powered
In compliance with	EN 60601-1: 2006 General Requirements for basic safety and essential performance. EN 60601-1-2: 2007 Collateral standards for Electromagnetic Compatibility
Ingress protection	IPX0 (Sealing Hand Gun), IP40 (Charger)
Electrical safety	During charging Class II. When in use: internal power source The T-SEAL MOBILE is used in the same environment as medical equipment (hospitals and blood banks). This instrument must be used by highly qualified personnel.

Battery

Manufacturer (legal and physical)	EVVATECHNOLOGY CO., LIMITED, 4F Bldg A3, New Material Industrial Park of Silicon Valley, Dafu Industrial Zone, Guanian, Longhua, 518110 Shenzhen, PEOPLE'S Republic of CHINA
Size (W x H x D) Battery unit	165 x 106 x 32 mm (6.5 x 4.2 x 1.2 in) +/- 10%
Weight, kg (lb)	0.68 kg (1.50 lb) +/- 10%
Nominal voltage	29.6 VDC
Rated capacity	1500mAh
Fuses	2 A (input) / 10 A (output) (internal)
Std	IEC 62133:2012

Seal Hand Gun

Size (W x H x D)	142 x 253 x 33 (mm) +/- 10%
Weight, kg	0.39 kg (0.86 lb) +/- 10%
RF Output	80 W max. / 50 Ω / 40.680 MHz
Output (Consumption)	Max. 150 W
Force of pressure on key	40-50N

Battery charger

Manufacturer	Friwo Geratebau GmbH Von-Liebig Strasse 11, D-48346 Ostbevern, Germany
Size (W x H x D)	102 x 65 x 30 mm +/- 10%
Weight, kg	0.275 kg (0.606 lb) +/- 10%
Input characteristics	100-240 VAC / 50/60 Hz, 700-350 mA
Output characteristics	18V DC / 1.66 A
Fuses	T2 A 250 V (internal)

Five Adapter Kit for the Charger

Five Adapter Kit for the Charger

Adapter (Euro standard):

Dimensions (H x W x D):	45 mm x 25 mm x 50 mm
Weight:	250g

Allowable deviation: +/-10%

Adapter (USA, Japan standard)

Dimensions (H x W x D):	45 mm x 25 mm x 40 mm
Weight:	250g

Allowable deviation: +/-10%

Adapter (UK Standard)

Dimensions (H x W x D):	50 mm x 40 mm x 45 mm
Weight:	250g

Allowable deviation: +/-10%

Adapter (Australia standard)

Dimensions (H x W x D):	55 mm x 40 mm x 40 mm
Weight:	250g

Allowable deviation: +/-10%

Universal adapter "IEC 320 mains connector"

Dimensions (H x W x D):	45 mm x 55 mm x 25 mm
Weight:	250g

Allowable deviation: +/-10%

Hand Gun Mount

Dimensions (H x W x D):	180 mm x 120 mm x 80 mm
Weight	296 g
	
The Hand Gun Mount allows for belt fixing the gun for more convenient professional work.	

Allowable deviation: +/-10%

Extension cord

Dimensions (H x W x D):	177 mm x 108 mm x 37 mm
Weights:	185 g
	Extension cord 1.2 m 5-pol allows to lengthen the space between hand unit and the battery, allows more freedom and access during the procedure. One end of the cord is plugs to the battery, another to the plug on the hand unit cord.

Allowable deviation: +/-10%

EMS tables

Table 1:

Guidance and manufacturer's declaration – electromagnetic emissions		
The T5460 T-SEAL MOBILE is designed to be used in the electromagnetic environment specified below. The customer or the user of the T5460 T-SEAL MOBILE should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment - guidance
RF emission CISPR 11	Group2	The T5460 T-SEAL MOBILE emits the electromagnetic energy required to perform its intended function. Nearby electronic equipment may be affected. The T5460 T-SEAL MOBILE is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emission CISPR 11	Class A	
Harmonic emission IEC 61000-3-2	Not Applicable	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Not Applicable	

Table 2:

Guidance and manufacturer's declaration – electromagnetic immunity			
The T5460 T-SEAL MOBILE is designed to be used in the electromagnetic environment specified below. The customer or the user of the T5460 T-SEAL MOBILE should assure that it is used in such an environment.			
Interference immunity tests	IEC 60601 Test level	Compliance	Electromagnetic environment - guidance
Electrostatic discharge (ESD)	±6 kV contact	6kV	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at exceed 30 mm.
IEC 61000-4-2	±8 kV air	8kV	
Electrical fast transient/burst	± 2kV for power supply lines	2kV	Mains power quality should be that of a typical commercial or hospital environment.

IEC 61000-4-4	± 1kV for input/ output lines	Not Applicable	
Surge	± 1kV line(s) to line(s)	1kV	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4-5	± 2kV line(s) to earth	2kV	
Voltage dips, short interruptions and voltage fluctuations on power supply input lines	< 5% UT (> 95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) For 5 cycles	< 5% UT (> 95% dip in UT) for 0,5 cycle 40% UT (60% dip in UT) For 5 cycles	Mains power quality should be that of a typical commercial or hospital environment.
IEC 61000-4- 11	70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 s	70% UT (30% dip in UT) for 25 cycles < 5% UT (> 95% dip in UT) for 5 s	
Power frequency (50/60 Hz) magnetic field	3 A/m	3 A/m	Power frequency magnetic fields should not exceed the levels typical for standard commercial or hospital environments.
IEC 61000-4-8			
 NOTE	UT is the a.c. mains voltage prior to application of the test level.		

Table 4

Guidance and manufacturer's declaration – electromagnetic immunity			
The T5460 T-SEAL MOBILE is designed to be used in the electromagnetic environment specified below. The customer or the user of the T5460 T-SEAL MOBILE should assure that it is used in such an environment.			
Interference immunity tests	IEC 60601 Test Level	Compliance Level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the T5460 TSEAL MOBILE including cables than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.17 \sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = 1.17 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.33 \sqrt{P}$ 800 MHz to 2.5 GHz

		Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey^a should be less than the compliance level in each frequency range.^b Interference may occur in the vicinity of equipment marked with the following symbol: 
 NOTES	At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. a. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) tele-phones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the T5460 T-SEAL MOBILE is used exceeds the applicable RF compliance level above, it should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the T5460 TSEAL MOBILE. b. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m	

Table 5

Recommended separation distance between portable and mobile RF communications equipment and the T5460 T-SEAL MOBILE			
The T5460 T-SEAL MOBILE is intended for use in the electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the product can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the T5460 T-SEAL MOBILE as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = 1.17 \sqrt{P}$	80 MHz to 800 MHz $d = 1.17 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.33 \sqrt{P}$
0.01	0.12	0.12	0.23
0.01	0.37	0.37	0.74
1	1.2	1.2	2.3
10	3.7	3.7	7.4
100	11.7	11.7	23.3
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer			
 NOTES	At 80 MHz and 800 MHz, the higher frequency range applies. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

10. INFORMATION ABOUT THE PHARMACEUTICALS/ANIMAL OR HUMAN MATERIALS IN THE MEDICAL DEVICE

The T-Seal Mobile contains no pharmaceutical products as well as materials of animal or human origin.

11. METHODS OF CLEANING AND DISINFECTION

The T-Seal Mobile is delivered unsterile and requires no sterilization.

11.1.1. CLEANING AND DISINFECTION METHODS

This chapter gives cleaning (procedure, frequency) information on the Battery and Hand units. The Sealer requires minimal maintenance for efficient operation. Perform the cleaning procedure as specified below.

WARNINGS

Always disconnect the Hand Gun and the Charger from the Battery for your own safety.

Blood and blood products must be treated as potentially infectious at all times. In the event of blood spillages, appropriate protective clothing and gloves should be worn during cleanup procedures.

After removing residual biological material, surfaces which have been in contact with blood or components must be disinfected using a chemical agent considered to be “sterilizing” (isopropyl Alcohol 70%...). Alternatively, a freshly prepared sodium hypochlorite (household bleach) solution may be used to disinfect surfaces. Diluted solutions of 1 part bleach to 10 parts water may be used.

Regardless of the “sterilant” or disinfection solution used, remember to remove any residue to ensure that surfaces of the equipment are not subject to corrosion or discoloration. Discard all materials in contact with blood according to institutional policies regarding disposal of biohazardous materials.

11.1.2. BATTERY

Cleaning may be required as a result of spilled drops of blood or once per month. If blood spillage occurs, the Battery must immediately be removed from service and cleaned completely before resuming use.

Use a soft lint-free tissue, moistened with a mild detergent to clean the outside of the Battery.

CAUTION

Do not, under any circumstances, submerge the Battery unit in any kind of liquid. This will damage the battery and void the warranty.

11.1.3. SEAL HAND GUN

Clean once per week or more frequently, as required, if spillage of blood occurs.

For cleaning the electrodes remove the moving electrode.

Removal of moving electrode:

1. Pull the blue trigger to “click” to release the moving electrode. Push spring lock forward.
2. First push inward, then pull out the electrode.

Clean the handle and both electrodes with a soft lint free cloth moistened with mild detergent. Dry carefully and ensure that the electrodes are completely dry to prevent sparks. To clean the electrode shaft, use a dry cotton swab. After cleaning, inspect the electrodes for any mechanical damage or wear out. Do not use damaged part.

For proper function, add a drip of light machine oil on the electrode shaft after each cleaning. Assemble the parts in reverse order. Make sure the trigger and spring lock is in full forward position before electrode insertion. Insert the moving electrode ensuring that it is positioned so as to be parallel with the fixed electrode. Press the trigger twice and ensure that the moving electrode moves smoothly, is in position and not loose.

CAUTION

The hand unit is not waterproof. Do not, under any circumstances, submerge the hand unit in any kind of liquid. Doing so will generate electric discharge, cause the device to malfunction, and void the warranty.

NOTES

Some sealing tests are recommended before resuming use.

Do not push the trigger in after removing the electrode since this may block the locking mechanism.

If electrode cannot be inserted after cleaning, the trigger is probably pressed during cleaning.

12. PRELIMINARY STARTING PROCEDURE

This chapter describes the unpacking procedure, temperature requirements and installation of the device.

Unpacking and inspection

1. Visually inspect the cardboard box for damage. Report any damage immediately.
2. Lift the instrument out of the cardboard box and place it on a flat surface.

The device is shipped in one cardboard box that includes:

- One (1) Seal Hand Gun
- One (1) Battery
- One (1) Charger (without adapter)
- One (1) Instruction for Use
- One (1) Five Adapter Kit for the Charger.

3. The above list is subject to change. Refer to the packing list for accurate description of the contents. If any parts are missing or damaged, report it immediately to competent persons. If necessary, individual parts and accessories can be delivered separately and be packed into separate carton boxes:
 1. Manual sealer gun.
 2. Battery (not more than 10 pcs. – If required).
 3. Charger (not more than 10 pcs. – If required).
 4. Five adapter kit for charger (not more than 10 pcs. – If required).
4. Keep all of the packing materials since they may become necessary for repeated transportation at least for the warranty period.

Installation procedure

Device preparation for sealing

1. Select the respective adapter to connect the Charger to mains.
2. The changeable mains adapters for European continent (standard), UK, USA/Canada and Australia/New Zealand, and IEC standard 320 C8 connectors are supplied.
3. Connect the charger (with grey connector) to the combined Hand unit/Battery charger cable connector.
4. Plug the charger into the mains. The charger can be connected to mains voltages between 100 and 240 VAC.
5. It takes about three hours to charge the Battery unit completely. When the charging level indication on the unit indicates 100% and the recharge indication flashes, the battery is fully charged.
6. Charge the battery to 100% before the first use.
7. Connect the Hand unit cable (with grey connector) to the Battery unit connector. Make sure the connector snaps in secured position.
8. If it is necessary to extend the distance between the handgun and the battery, an extension wire can be connected. One end of the wire is connected to the connector on the battery, the other - to the connector on the cable of the handgun
9. Perform a test seal on an empty or water filled tube to ensure proper operation.

CAUTION

Only use the enclosed Battery charger for charging battery in the Battery unit.

NOTE

Ensure that cable to handle is connected before use.

Check of the battery

To check the battery, just press the “TEST” button **③**, and the battery charge level indicator **①** lights up. The capacity indication is in steps of twenty percent. The battery can be recharged at any time or at least when zero (red) remains. The zero charge level is also indicated on the gun with the red light and when seal is attempted, the seal function will be blocked. Recharging the battery:

Connect the Battery charger to the Battery **④** and to the mains. The recharge indicator **②** lights up.

When the recharge indicator starts flashing, the battery is fully charged and the charge current is automatically turned off. The charge can be interrupted at any time without damage to the cells, but a full recharge will take less than three hours and is recommended.

CAUTIONS

Always check the seals when the battery capacity is low (20% or less). See section 7.2 for sealing pattern.

Only use the enclosed Battery charger for charging battery in the Battery unit.

NOTES

To guarantee a long lifetime for the rechargeable battery, please observe the following requirements: Charge the battery at temperature between 5 and 35°C (40-95°F).

If the indicator **2** does not light up, check the connection to the charger and mains.

Preparation before use

1. Place the instrument on a flat surface near the working place or in the holster.
2. Check that there is sufficient battery capacity by pushing the "TEST" button.
3. Connect the Hand Gun cable to the combined hand gun/charger connector.

NOTE


Do not connect to the T-SEAL MOBILE Battery other sealers in addition to T-SEAL MOBILE hand unit manufactured by Conroy Medical AB.

13. WORKING SEQUENCE

This chapter describes how the T-SEAL MOBILE works, where the connectors and the indicators are located, and their functions.

Battery


Front view

Rear View


Appearance of the Sealing Hand Gun

Appearance of the Charger


The device ready for use (hand unit attached to the cord and the battery):

During the sealing process with T-Seal Mobile medical personal should wear protective wear and gloves.

How to seal the tubes


1. Put the tube to be sealed down to the Hand gun slot.
2. Check that the tube is placed in the interelectrode area.


Tube in the hand unit before the sealing

3. Press the trigger to bring the two electrodes closer together until the light on top of the Hand unit comes on. The sealing time is normally 0.4 to 0.9 seconds; after a maximum of 3 seconds the RF is turned off.

Sealing indicator

The indicator is on during the sealing process.

The picture of the device in the process of sealing

4. When the light goes out, sealing is finished. Release the trigger and remove the tubing.
5. The center of the sealed pattern is very thin and pulling both sides will divide the tube into two pieces.
6. Check the tube for leakage.

 **CAUTION**

 **CAUTIONS**

If you should make two or more seals, they should not be within 1 cm (0.4 in) of each other, otherwise the resulting pressure in the tube may cause microscopic cracks and holes in the seals.

In case of sealer malfunction (intermittent operation, poor seal quality, too long or too short sealing time) contact your local authorized representative of the manufacturer for assistance. Periodically check whether the shape of the seal is symmetrical and there is a clear line in the middle of the seal (see picture below).

Good sealing pattern

Bad sealing pattern

NOTES

The tube must be dry on the outside.

Do not pull the tube and keep the trigger fully pressed down until the light goes out. If the trigger is released the sealing process stops.

The

If the sealer does not start

The sealer has several safety functions. Before sealing, these functions control and identify whether it is possible to seal the tube. The table below covers the common probable causes for problems and suggests some recommended actions.

Probable cause	Recommended action
Leaky or wet tube	Dry the tube and try again
Tiny arcs between the electrodes	Dry the electrodes and try again
Wet or dirty electrodes	Clean and dry the electrodes
Low/no capacity in battery	Check/change the battery
Overheated (Hand unit LED red blinking)	Let the Hand unit cool down
No seals	Check Battery unit and connection

If the sealer still doesn't start, refer to the chapter **Troubleshooting** to get acquainted with extra information.

LED indications

Battery

- “Colored LEDs bar indicating battery status.”
- “TEST”: button to check the battery. When pressed, the battery charge level indicator lights up.

Hand gun

- Indicator lights orange when RF-energy is being applied to the fixed electrode and turns off when sealing is completed.
- Indicator flashes red at overheating, sealing function is blocked.
- Indicator lights red one second in case of short circuit, sealing function is blocked.
- Indicator lights red in case of battery low level, sealing function is blocked.

Troubleshooting

Maintenance carried out by the user is limited to changing Battery unit and Hand unit. The following information covers common problems and offers suggested solutions.

Battery:

Problem	Probable cause	Recommended action
Charge level LEDs doesn't light.	Battery too low voltage	Charge battery for 3 hours.
	Blown fuse.	Contact with the manufacturer's authorized representative in Russia.
	Other cause	Contact with the manufacturer's authorized representative in Russia.

Charge Level LEDs do not light up during charge	No voltage	Charge battery for 48 hours
	Blown fuse.	Contact with the manufacturer's authorized representative in Russia.
	Other cause	Contact with the manufacturer's authorized representative in Russia.
Charge Level LEDs do not light up after 48 hours charge	Battery defect	Replace the battery or contact with the manufacturer's authorized representative in Russia.

Hand gun:

Problem	Probable cause	Recommended action
The sealing doesn't start.	Wet electrodes.	Dry the electrodes
	Dirty electrodes.	Clean the electrodes
	Wet tube.	Dry the tube and electrodes, try again.
	Tiny arcs between the electrodes.	Dry the electrodes and try again.
	The trigger does not "click" when it is fully squeezed.	Replace the hand gun or contact with the manufacturer's authorized representative in Russia.
	The trigger cannot fully engage when it is squeezed	
	Defect in the hand gun	
	Battery not charged/low charge.	Charge battery for 3 hours.
	Cable broken	Contact with the manufacturer's authorized representative in Russia.
	Other cause	
Bad seal	Low battery charge	Check the battery charge level. Charge battery if needed.
	Mobile electrode doesn't move smoothly in hand unit.	Clean and oil the electrode shaft.
	Wet electrodes.	Dry the electrodes and try again.
	Dirty electrodes.	Clean the electrodes
	Moving electrode out of position.	Reassemble the electrode.
	Defect in the hand gun	Replace the hand unit or contact with the manufacturer's authorized representative in Russia.
	Other cause	Contact with the manufacturer's authorized representative in Russia.
Hard to divide tube after seal	Mobile electrode doesn't move smoothly in hand unit.	Clean and oil the electrode shaft.
	Wet electrodes.	Dry the electrodes and try again.
	Dirty electrodes.	Clean the electrodes
	Damaged edge of the electrode	Contact with the manufacturer's authorized representative in Russia.
	Other cause	
Intermittent sealing	Defect in the hand gun	Change the Hand unit or contact with the manufacturer's local representative.
	Cable broken	Contact with the manufacturer's authorized representative in Russia.
	Battery defect	

The Hand unit's LED still alight after seal trigger has been released.	Defect in the hand gun	Replace the hand unit or contact with the manufacturer's authorized representative in Russia.
The Hand unit LED doesn't light	Defect in the hand gun	Contact with the manufacturer's authorized representative in Russia.
The Hand unit LED flickers.	Cable broken	Contact with the manufacturer's authorized representative in Russia.
	Other cause	
Electrode cannot be inserted after cleaning	Trigger has been pressed during cleaning which means that internal spring can come out of position.	The spring can be reset to its correct position by pulling out the trigger to the outer locked position. Thereafter carefully insert a pin (Ø3- 4mm, 100mm long) in the front hole where the electrode is inserted and slightly push it downwards. Be careful; do not scratch the surface with the pin.

Preventive Inspection and Maintenance

At least every second year, the T-SEAL MOBILE must be fully inspected:

- Either by the qualified service organization of the original purchaser,
- Or by the manufacturer's authorized representative in Russia.

Interchangeable Parts Replacement (Maintenance performed by user) Maintenance is limited to changing the Hand unit, Battery unit or charger. The information in chapter ***Troubleshooting*** covers common problems and offers suggested solutions.

14.ELECTROMAGNETIC COMPATIBILITY ACCORDING TO IEC 60601-1-2

Class II ME EQUIPMENT, externally powered

15.STORAGE AND TRANSPORT CONDITIONS

The device can be transported by any roofed transport when following conditions of transportation. When transporting or storing T-Seal Mobile or its accessories, place them in the dedicated transport case or original shipping carton. Make sure that no liquid can fall on the case or shipping carton and that the storage conditions are met. Do not load or transport package if damaged.

Storage and transport conditions:

Temperature: -20 to +70°C (-4 to 158°F)

Humidity 10 to 90% (non condensing)

16.APPLICATION CONDITIONS

To keep the instrument operating at its best, please observe the following:

- The instrument should be placed on a flat surface free from dust, solvents and acidic vapor.

- Use the instrument in an area free from vibration and with a room temperature of 15 to 35°C (60 to 95°F), and relative humidity 10% to 90%.
- Handle the components with care in a clean environment.

17. THE LIST OF APPLICABLE STANDARDS

Standards	
EN ISO 13485:2012 + AC:2012	Medical devices – Quality Management Systems - Requirements for regulatory purposes
EN ISO 14971:2012	Medical devices -- Application of risk management to medical devices
EN/(IEC) 61000-4-2:2009	Immunity to electrostatic discharges
EN/(IEC) 61000-4-3:2006+A1+A2	Immunity to radio frequency electromagnetic fields
EN/(IEC) 61000-4-8:2010	Immunity to power frequency magnetic field
CISPR 16-2-3:2010+A1	Measurement of radio frequency electromagnetic field
EN/(IEC) 61000-4-28:2000+A1+A2	Immunity to variations of power frequency
IEC 62133:2012	Safety Requirements for Rechargeable Cells & Batteries used in Portable Devices
UN38.3	Testing-Safety Precautions for the Transportation of Lithium Batteries
IEC 60601-1:2005	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance

18. LABEL: DESCRIPTION OF SYMBOLS

This marking reflects compliance with the Council Directive 93/42/EEC on Medical Devices.

This marking reflects compliance with EN 60601-1 and national standards for USA (ANSI/AAMI ES60601-1:2005) and Canada (CSA C22.2 No. 60601-1) Markets

Class II Battery charger: double insulation during charging.

Battery charger UL recognized component for Canada and U.S.A.

Symbol for "CATALOGUE NUMBER".

Symbol for "WARNING".

Follow instructions.

Symbol for "SERIAL NUMBER".

Symbol for "DATE OF MANUFACTURE".

Symbol for "Manufacturer".

Storage conditions - Symbol for "TEMPERATURE".

Storage conditions - Symbol for "RELATIVE HUMIDITY".

Symbol for "Non-ionizing radiation".

Symbol (WEEE 2002/96/EC) - Do not dispose Product as municipal waste. Collect Product separately. Use collection and return systems available to you. Product brought to EU market after August 13th, 2005.

Indicates information about the battery

19. INFORMATION ABOUT RECYCLING OR DISPOSAL. HAZARD CATEGORY OF MEDICAL WASTE

According to SanPin 2.1.7.2790-10 "Sanitaru-epidemiological requirements for medical waste recycling", this device refers to Class B.

Product Disposal

For disposal of the T-SEAL MOBILE, or parts thereof, (including accessories) at the end of the calculated life cycle of 7 years, please ensure the following:

- Do not dispose T-SEAL MOBILE as unsorted municipal waste.
- Collect the T-SEAL MOBILE separately.
- Use the collection and return systems available to you.

For more information on return, recovery or recycling of the T-SEAL MOBILE, please contact with the manufacturer's local representative in Russia.

20.LIFETIME

The durability or life span of the equipment is seven (7) years from manufacturing year provided that:

- The equipment is handled and stored according to the conditions prescribed on the packaging and in the user manual.
- Periodical preventive maintenance is performed as instructed by the instructions for use and using original spare parts.

Storage life is 7 years.

21.WARRANTY STATEMENT

21.1. Guarantee

Conroy Medical AB guarantees that the equipment shall be free from defects in material and workmanship when delivered to the original purchaser. Conroy Medical AB's sole obligation shall be limited to free-of-charge repair or replacement, at Conroy Medical AB's option, of the defective part or unit for a period of two (2) years following the date of initial delivery to the original purchaser.

The warranty extends only to the original purchaser, is not assignable or transferable, and shall not apply to auxiliary equipment or disposable accessories.

Conroy Medical AB guarantees that the equipment is fit for the purposes and indications described in the labeling when used in accordance with the instructions for use. Unless the equipment is used in accordance with such instructions, this warranty is void and of no effect. Conroy Medical AB shall grant no other warranty, expressed or implied, including any warranty of merchantability or appropriateness for a particular purpose. Conroy Medical AB's sole obligation and original purchaser's exclusive remedy for breach of warranty shall be limited to repair or replacement at Conroy Medical AB's option. Conroy Medical AB shall not be liable for any incidental or consequential damages. Modifications, alterations, recalibrations or abuse, and service by other than a Conroy Medical AB authorized representative will void the warranty.

21.2. Service

Service under warranty period

While under Conroy Medical AB warranty, the instrument must not be opened by unauthorized personnel. Contact your local Sales Office for service and repair information for all T-SEAL MOBILE instruments.

Official representative in Russian Federation for complaints acceptance:

Terumo Russia LLC

Testovskaya Str. 10, 123112, Moscow, Russia, 13-th floor, premise I, room 5

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