The State of Maryland

Office of the Secretary of State

This Apostille is not valid for use anywhere within the United States of America, its territories or possessions.

This Apostille does not certify the content of the document for which it is issued.

Apostille

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1. Country: United States of America
This public document
2. has been signed by **Scott A. Poyer**
3. acting in the capacity of **Clerk of the Circuit Court**
4. bears the seal/stamp of the **Clerk of Circuit Court of Anne Arundel County**

Certified

5. at Annapolis, Maryland
6. the **4th day of March, 2019**
7. by The Secretary of State of Maryland
8. **No. 476876**
9. Seal

10. Signature

Secretary of State

State of Maryland, Anne Arundel County, Sct.

I, **SCOTT A. POYER**, Clerk of the Circuit Court for Anne Arundel County, Maryland, a court of record, do hereby certify that **Linda Davis** was a commissioned/appointed and qualified Notary Public commencing on the 1st day of May, 2017.

In Testimony Whereof, I have hereunto set my hand and affixed the seal of the court March 4, 2019

Scott A. Poyer

Scott A. Poyer, Clerk
Circuit Court for Anne Arundel County, Maryland

Signed *[Signature]*

Subscribed and sworn to before me this

4 day of March, 2019.

[Signature]
Notary Public

LINDA DAVIS
Notary Public-Maryland
Anne Arundel County
My Commission Expires
May 14, 2021

Affidavit

State of Maryland)
)SS:
County of Anne Arundel)

Before me, the undersigned notary public, this day, personally,
appeared to me known or proven, who being duly sworn according to law,
deposes the following:

Yoshiko Yano

I, _____, hereby certify that the attached
document is

- ☒ an original document
☐ a true copy of an original document

to the best of my knowledge and belief.

Signed _____

Subscribed and sworn to before me this

4 day of March, 2019.

Notary Public

LINDA DAVIS
Notary Public-Maryland
Anne Arundel County
My Commission Expires
May 14, 2021

Bayer Medical Care Inc
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.
(412) 767-2400
www.radiology.bayer.com

To whom it may concern,

22 February 2019

To be submitted to the Competent Health Authorities in Federal Republic of Russia

We, Bayer Medical Care Inc., USA (1 Bayer Drive Indianola Pennsylvania 15051-0780, USA) hereby state that attached document (Instruction for use) (Single patient tube 250 cm (SPD 250) for the injection system MEDRAD Stellant) is true.

Best regards,

William Hohn
William Hohn 22-2-19

Regulatory Affairs International Strategist – Europe
22 February 2019

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For use with injector software versions 105.6 or higher and 901.4 or higher
and with Certegra™ Workstation software version 101.05 or higher

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©2017 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer.

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Bayer, the Bayer Cross, MEDRAD, Stellant, MEDRAD Stellant, Certegra, MEDRAD Fluidots and Fluidots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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For use with injector software versions 105.6 or higher and 901.4 or higher and with Certegra™ Workstation software version 101.05 or higher

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©2017-2018 Bayer. This material may not be reproduced, displayed, modified, or distributed without the express prior written consent of Bayer

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Bayer, the Bayer Cross, MEDRAD, Stellant, MEDRAD Stellant, Certegra, MEDRAD FluiDots, and FluiDots are trademarks owned by and/or registered to Bayer in the U.S. and/or other countries.

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UNITED STATES

Bayer Medical Care Inc.
1 Bayer Drive
Indianola, PA 15051-0780
U.S.A.

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ENGLISH

Important Safety Notice

This device is intended to be used by individuals with adequate training and experience in diagnostic imaging studies.

PLEASE FOLLOW YOUR INSTITUTION'S POLICY FOR MEDIA ADMINISTRATION AND ASEPTIC TECHNIQUE.

NOTE: If cleaning of the multi-patient kit is desired, use only ethanol-based disinfectants. Exposure to other cleaning agents or disinfectants may lead to component failure.

NOTE: Refer to the injector system operation manual for additional operating instructions.

Intended Use

The contents of this package are intended to be used in the delivery of contrast media or saline. These devices are indicated for use only with MEDRAD® Stellant CT Injection System.

Warning

Follow aseptic technique when setting up and using disposables. Do not touch the end of the spike or connectors. Do not remove protective covers until ready to make connections. Properly discard the patient tube after each patient injection and the syringes and injector tubing within 12 hours. Contamination can result from failure to follow aseptic technique.

Do not use the multi-patient kit (syringes and injector tubing) for longer than 12 hours. These disposables have been designed and validated for only 12 hours. Extended use of the multi-patient kit poses risks to the patient including contamination and injury due to device malfunction.

Do not exceed the media manufacturer's recommended use time. Refer to the media manufacturer's "Instructions for Use" for specific indications and use time.

For devices labeled for single use, please note: This product is intended for single use only. Do not re-sterilize, reprocess or reuse. The disposable devices have been designed and validated for single use only. Re-use of the single use disposable devices pose risks of device failure and risks to the patient. Potential device failure includes significant component deterioration with extended use, component malfunction, and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be cleaned or re-sterilized.

For devices labeled for multi-use, please note: Do not re-sterilize or reprocess. Potential device failure includes significant component deterioration and system failure. Potential risks to the patient include injury due to device malfunction or infection as the device has not been validated to be reprocessed or re-sterilized.

Do not fill or inject unless the syringe is properly engaged. Patient injury can result if syringe is not properly engaged.


Warning

Do not re-install syringes onto an injector. Syringe sterility will be compromised, and patient infection can result. Use loaded syringes immediately. Discard unused loaded syringes. Contamination can occur if syringes are used to store contrast media.

Visually inspect contents and package before each use. Do not use contents if sterile package is opened or damaged. Patient or operator injury can result.

Do not exceed pressures identified on the front of the package; ensure that the fluid path is open. Use of greater pressures or occlusions in the fluid path can result in fluid leaks or tubing ruptures and patient or operator injury.

1 / 2

Configure Settings:

1. Press **Setup/Fluid Delivery Setup** button.
2. Set **LPCT** type to P8100 LPCT.
3. Set **Auto Load Purge Volume** to **OFF** to reduce fluid waste during syringe filling and refilling.
4. Set **Fluid Control** to **OFF** to decrease resistance felt when turning manual knobs and reduce leaking of fluid.
5. Press **OK** to save settings.

1 / 3

Install Syringes:

1. Insert and push syringes until they snap into place.

MEDRAD® FluidDots Round FluidDots:
Syringe filled with Fluid

Ellipsoidal FluidDots:
Syringe filled with Air

Caution

Ensure all connections are secure; do not over-tighten. This will help minimize leaks, disconnection, and component damage.

Install Fluid Transfer Set:

1. Remove white caps from both syringes.
2. Remove orange caps from fluid transfer sets.
3. Install fluid transfer sets onto syringe tips.

Install Connector Tube:

1. Remove white caps from fluid transfer set.
2. Remove small orange caps from connector tube.
3. Install the end of the connector tube that looks like a "T" onto Syringe A (Contrast).
4. Install opposite tube onto Syringe B (Saline).

Install Use-Time Label:

1. Write the syringe and injector tubing installation and expiration time on the label.
2. Slide the label onto the transfer set or connector tube, so it is visible.

1 / 4

Install Single Patient Tube:

1. Remove a new patient tube from packaging.
2. Remove orange caps from connector tubing and patient tube.
3. Install patient tube onto the injector tubing.

NOTE: Leave white cap on patient tube until ready to connect to patient catheter.

Dust cap on: new patient tube

Dust cap off: used patient tube

Use care in handling and inserting spike into contrast bottle and saline bag. The spike is sharp and can cause personal injury.

Spike Fluids (AS NEEDED):

1. Remove cap from spike and fluid source.
2. Without touching the spike end, insert spike into fluid source.
3. Hang the spiked fluid source.
4. Repeat steps 1-3 for other fluid source.

FILL Syringes (AS NEEDED):

1. Fill syringes either automatically or manually (See Below).
2. Install heat maintainers.
3. Ensure air is expelled from syringes.

AUTOMATIC FILL:


Program a new or recall an existing protocol and press **Auto Load** on injector head. Press **FILL A** or **FILL B** as needed.

NOTE: For first-time auto-fill, air remains in the syringe.

MANUAL FILL:

Press **Move Piston** on injector head and use forward and reverse piston controls to fill syringes to desired amount.

NOTE: Ensure final piston movement is in the forward direction. If the Syringe Undocked error appears, manually advance piston to re-dock syringe.

Prime Tubing:

NOTE: A whistling sound during priming is normal.

NOTE: While priming, tap the check valves on the patient tube to remove air bubbles.

1. Prime tubing.
2. Confirm that air is expelled from tubing.

AUTOMATIC PRIME:

Press **Prime** to prime injector and patient tubing.

MANUAL PRIME:

Either press **Move Piston** and use forward controls or turn manual knobs on injector head.

Warning

Clear all air from the syringes and tubing before connecting a patient to the injector. Air embolization can cause death or serious injury to the patient.

Connect to the Patient:

1. Rotate injector head down.
2. Confirm there is no air in syringes or tubing.
3. Remove cap from patient tube and connect to the patient.
4. Press **I Checked for Air** button on the injector head.

Prepare & Start Injection:

1. Ensure protocol is correct for current patient.
2. **Lock** the protocol at the control room monitor.
3. Press **Arm**.
4. Press yellow **Start** button.

AFTER EACH PATIENT Discard Patient Tube:

1. Disconnect patient tube from the patient catheter.
2. Disconnect patient tube from injector tubing.
3. Discard patient tube in accordance with facility policy.
4. Immediately install new patient tube.
5. Repeat 1-per-patient steps as needed within the 12-hour use life.

WITHIN 12-HOUR USE TIME Discard multi-patient kit:

1. Remove syringe heat maintainers.
2. Remove syringes. Rotate syringes 1/4 turn counter clockwise and pull out.
3. Remove used fluid sources.
4. Discard syringes, tubing, and fluid sources in accordance with facility policy.

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