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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY: The seal of WUXI APEX MEDICAL CO., LTD on the annexed INSTRUCTION FOR USE is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature : SHEN XIAOHUA

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(Date: MAR. 24, 2015)

Instruction for Use

Sterilized Infusion Set for Single Use
Model:FAST SET-II

Wuxi Apex Medical Co., Ltd.


Apex medical Fast Set

Sterilized Infusion Set for Single Use

Model: FAST SET-II

Apex Medical

Refer to instruction for use of pump for installation

Use

Device

The Fast Set infusion set is an infusion system, the proximal end of which is connecting to the reservoir with the medicinal product. Located at the distal end of the system of permanent catheter introduced subcutaneously. The motor drive to push the reservoir, so that insulin will flow into the infusion set. The insulin will be delivered to the human body subcutaneously through the infusion set.

The Fast Set infusion set is an infusion system, the proximal end of which is connecting to the reservoir with the medicinal product. Located at the distal end of the system of permanent catheter introduced subcutaneously. The motor drive to push the reservoir, so that insulin will flow into the infusion set. The insulin will be delivered to the human body subcutaneously through the infusion set.

Indications for use

This product is indicated for the subcutaneous infusion of medication, including insulin, from an infusion pump.

Contraindications

Sterilized infusion sets for single use are indicated for subcutaneous use only and not for intravenous infusion or the infusion of blood or blood products.

Complications

Infection, allergy to stainless steel or teflon, bleeding, needle occlusion, subcutaneous infiltration, needle breakage

Physical Properties

Infusion set: soft cannula-1pc, multiple tubing-1pc, insertion device-1pc
Package: plastic casing, dialyzing paper.

Sterilization Method: ethylene oxide

Content of the set: 10pcs/box


Compatible Pumps

The infusion set should be used with TruCare insulin pump (produced by Apex Medical), Minimed Paradigm MMT-712, Minimed Paradigm MT-722, Minimed Paradigm REAL-TIME MT-722, Minimed Paradigm VEO MMT-754 (produced by Medtronic).

Specification

Reference Number	Catheter Length	Total Tubing Length
FS206051	6mm	50cm
FS206101	6mm	100cm
FS209051	9mm	50cm
FS209101	9mm	100cm

Drawings

1-Connector; 2-Tubing; 3-Guide Needle; 4-Hub; 5-Protective Cap; 6-Separator Cover; 7-Protective Cover; 8-Soft Cannula

Shelf Life: 2 years

Handling: Keep the product away from direct sunlight

Usage time: Max. 72 Hours

Disposal

Dispose in accordance with local regulation

Warnings

Inaccurate medication delivery, infection and/or site irritation may result from improper insertion and maintenance of the infusion site.

Before insertion, clean the insertion site with isopropyl alcohol.

Remove the needle guard before inserting the infusion set.

If using this infusion set for the first time, do the first set-up in the presence of your healthcare professional.

Do not leave air in the infusion set. Prime completely.

Check frequently to make sure the soft cannula remains firmly in place as you may not feel pain if it pulls out. The soft cannula must always be completely inserted to receive the full amount of medication.

If the infusion site becomes inflamed, replace the set, and use a new site until the first site has healed. Replace the infusion set if the tape

, or if the soft cannula
or partially dislodged from the

ie infusion set every 48-72 hours
g to Centers for Disease Control
nes, or per your healthcare
essional's instructions.

Do not reinsert the introducer needle into the
infusion set. Re-insertion may cause tearing
of the soft cannula, which may result in
unpredictable medication flow.

Never prime the set or attempt to free a
clogged line while the set is inserted. You may
accidentally inject too much medication.

Do not put disinfectants, perfumes, or
deodorants on the infusion set as these may
affect the integrity of the set.

Dispose of the infusion set and introducer
needle safely, in a sharps container, after a
single use. Do not clean or re-sterilize.

This device is sterile and non-pyrogenic
unless the package has been opened or
damaged.

Do not use if the package has been opened
or damaged. (For TruCare users, do not use
the infusion set if the tubing connector needle
has been damaged.)

Store infusion sets in a cool, dry place. Do not
leave infusion sets in direct sunlight or inside
a vehicle.

Use aseptic techniques when temporarily
disconnecting the set and consult your
healthcare provider on how to compensate
for missed medication when disconnected.

If infusing insulin, and your blood glucose
level becomes unexplainably high, or an
occlusion alarm occurs, check for clogs
and/or leaks.

If in doubt, change the infusion set because
the soft cannula may be dislodged, crimped
and/or partially clogged. Should any of these
problems arise, make a plan with your
healthcare professional for rapidly replacing
insulin. Test your blood glucose level to make
sure the problem is corrected.

If infusing insulin, do not change the infusion
set just before bedtime unless you can check
your blood glucose 1-3 hours after insertion.
Reuse of the infusion set may cause damage

to the cannula/needle and lead to infection, site
irritation, and/or inaccurate medication delivery.

Transport

Transportation should be conducted in
accordance with the contract which two
sides agreed with.

Storage

The infusion set should be stored in the room
with less than 80% relative humidity,
non-corrosive gases and well-ventilated,
and the products should be fully protected.

Operating Condition

Temperature: 5°C~40°C; Humidity: 20%~80%;
Barometric Pressure: 50kPa~106

Validity

On the condition that comply with the
requirements of the storage and
transportation, sterilized reservoir packing is
valid for two years.

Warranty

Wuxi Apex Medical Co., Ltd. warrants its
infusion sets and reservoir against defects in
materials and workmanship for the period up to
3 days after the individual reservoir/infusion
set sterile packaging has been opened, not to
exceed 6 months from date of shipment of the
infusion set or reservoir to the end user (the
"Warranty Period").

During the Warranty Period, Wuxi Apex
Medical Co., Ltd. will replace any defective
reservoir or infusion set, subject to the
conditions and exclusions stated herein.

This warranty is valid only if the Wuxi Apex
Medical Co., Ltd. Infusion sets and reservoir
are used in accordance with the Instruction for
Use. This warranty will not apply:

If the infusion set or reservoir has been used
for more than a single-time use by a single
end-user;

If damage results during the improper opening
of the sterile package not in conformance with
the procedures outlined in the Instruction for
Use;

If damage results from changes or
modifications made to the infusion set or
reservoir by the user or third persons after the
date of manufacture;

If damage results from service or repairs
performed to any part of the infusion sets or

part of the infusion sets or person or entity other than the

results from a Force Majeure or not beyond the control of the manufacturer; or

age results from negligence or improper including but not limited to improper storage physical abuse such as dropping or otherwise. This warranty shall be personal to the original purchaser. Any sale, rental or other transfer or use of the product covered by this warranty to or for a user other than the original purchaser shall cause this warranty to immediately terminate. The remedies provided for in this warranty are the exclusive remedies available for any breach thereof. Neither Apex Medical nor its supplier or distributors shall be liable for any incidental, consequential, or special damage of any nature or kind caused by or arising out of defect in the product. All other warranties express or implied are excluded, including the warranties of merchantability and fitness for a particular purpose.

Apex Medical Co., Ltd.

Address: Suite 601, 5Bld, No. 7 Changjiang Road, High Tech Area, Wuxi, Jiangsu, China 214028

Phone: +86-510-80100606
Fax: +86-510-80100505
Website: www.apexdiabetes.com

Shanghai International Holding Corp. GmbH (Europe)
Address: Eiffelstrasse 80, 20537 Hamburg, Germany
Phone: +49-40-2513175

Representative in Russia

Diamama, LLC, address: 121471, Moscow, Yabinovaya Str., 43, bldg. 2
Phone: +7 (495) 978-78-17
Email: info@dia-ma.ru


Symbols

Conformité européenne (European Conformity). This symbol means that the device fully complies with European Directive 93/42/EEC.

Caution: Consult Accompanying Documents

Sterilization: Ethylene-Oxide Gas

Use by

Do Not Reuse

Reference Number

Lot Number

Caution: U. S. federal law restricts this device to sale by or on the order of a physician.

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
China Council for the Promotion of International Trade is China Chamber of International Commerce

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


Model:EASY SET-I



Apex Medical

Refer to pump and reservoir user guides for installation and priming.

Disconnecting

Reconnecting

Use

Device

The Easy Set infusion set is used in combination with the insulin pump for the continual injection of insulin with the established and changing speed for the treatment of diabetes mellitus for insulin requiring patients.

The Easy Set infusion set is an infusion system, the proximal end of which is connecting to the reservoir with the medicinal product. Located at the distal end of the system of permanent catheter introduced subcutaneously. The motor drive to push the reservoir, so that insulin will flow into the infusion set. The insulin will be delivered to the human body subcutaneously through the infusion set.

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Complications

Infection, allergy to stainless steel or teflon, bleeding, needle occlusion, subcutaneous irritation, needle breakage

Physical Properties

Infusion set: stainless steel needle-1pc, multiple tubing-1pc
Package: plastic casing, dialyzing paper.

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Content of the set: 10pcs/box


Compatible Pumps

The infusion set should be used with TruCare insulin pump (produced by Apex Medical), Roche Accu-Check D-TRON Plus, Roche Accu-Check Spirit and Roche Accu-Check Combo (produced by Roche).

Specification

Reference Number	Catheter Length	Total Tubing Length
ES106050	6mm	50cm
ES106061	6mm	60cm
ES106100	6mm	100cm
ES106111	6mm	110cm
ES109050	9mm	50cm
ES109061	9mm	60cm
ES109100	9mm	100cm
ES109111	9mm	110cm

Drawings

1-Connector; 2-Tubing; 3-Separator; 4-Hub; 5-Stainless Steel Needle

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if the soft cannula
partially dislodged from the

infusion set every 48-72 hours
to Centers for Disease Control
as, or per your healthcare
professional's instructions.

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Do not prime the set or attempt to free a
clogged line while the set is inserted. You may
accidentally inject too much medication.

Do not put disinfectants, perfumes, or
odorants on the infusion set as these may
affect the integrity of the set.

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The device is sterile and non-pyrogenic
unless the package has been opened or
damaged.

Do not use if the package has been opened
or damaged. (For TruCare users, do not use
infusion set if the tubing connector needle
has been damaged.)

Store infusion sets in a cool, dry place. Do not
store infusion sets in direct sunlight or inside
a vehicle.

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When infusing insulin, and your blood glucose
level becomes unexplainably high, or an
occlusion alarm occurs, check for clogs
and/or leaks.

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When infusing insulin, do not change the infusion
set just before bedtime unless you can check
your blood glucose 1-3 hours after insertion.

Reuse of the infusion set may cause damage to
the cannula/needle and lead to infection, site
irritation, and/or inaccurate medication delivery.

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If damage results from service or repairs

to any part of the infusion sets or by any person or entity other than the manufacturer;
 damage results from a Force Majeure or event beyond the control of the manufacturer; or
 damage results from negligence or improper use, including but not limited to improper storage or physical abuse such as dropping or otherwise. This warranty shall be personal to the original purchaser. Any sale, rental or other transfer of the product covered by this warranty to or by a person other than the original purchaser shall void this warranty to immediately terminate. The remedies provided for in this warranty are the exclusive remedies available for any breach of. Neither Apex Medical nor its supplier or distributors shall be liable for any incidental, consequential, or special damage of any nature and kind caused by or arising out of defect in the product. All other warranties express or implied are excluded, including the warranties of merchantability and fitness for a particular purpose.

Apex Medical Co., Ltd.

Suite 601, 5Bld, No. 7 Changjiang
 Road, High Tech Area, Wuxi, Jiangsu,
 China 214028
 Tel: +86-510-80100606
 Fax: +86-510-80100505
 Email: info@apexdiabetes.com

EC REP

Shanghai International Holding
 GmbH (Europe)
 Eiffeistrasse 80, 20537
 Hamburg, Germany
 Tel: +49-40-2513175

Representative in Russia

Diamama, LLC, address: 121471, Moscow,
 Rabinovaya Str., 43, bldg. 2
 Tel: +7 (495) 978-78-17
 Email: info@dia-ma.ru

Symbols

Conformité européenne (European Conformity).
 This symbol means that the device fully complies
 with European Directive 93/42/EEC.

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Lot Number

Rx only

Caution: U. S. federal law restricts this device
 to sale by or on the order of a physician.

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








or Use

ie Device

Device Fast Set infusion set is using in combination with the insulin continual injection of insulin with the polished and changing speed for the treatment of diabetes mellitus for insulin requiring patients.

The Fast Set infusion set is an infusion system, the proximal end of which is connecting to the reservoir with the medicinal product. Located at the distal end of the system of permanent catheter is introduced subcutaneously. The motor drive to push the reservoir, so that insulin will flow into the infusion set. The insulin will be delivered to the human body subcutaneously through the infusion set.

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Contraindications

Sterilized infusion sets for single use are indicated for subcutaneous use only and not for intravenous (IV) infusion or the infusion of blood or blood products.

Contraindications

Sterilized infusion sets for single use are indicated for subcutaneous use only and not for intravenous (IV) infusion or the infusion of blood or blood products.

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Content of the set: 10pcs/box


Compatible Pumps

The infusion set should be used with TruCare insulin pump (produced by Apex Medical), Minimed Paradigm MMT-712, Minimed Paradigm MMT-722, Minimed Paradigm REAL-TIME MMT-722, Minimed Paradigm VEO MMT-754 (produced by Medtronic)

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or if the soft cannula
partially dislodged from the

infusion set every 48-72 hours
to Centers for Disease Control
es, or per your healthcare
ssional's instructions.

not reinsert the introducer needle into the
fusion set. Re-insertion may cause tearing
of the soft cannula, which may result in
unpredictable medication flow.

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If infusing insulin, do not change the infusion
set just before bedtime unless you can check
your blood glucose 1-3 hours after insertion.

Reuse of the infusion set may cause damage

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erson or entity other than the

results from a Force Majeure or beyond the control of the

ar; or

results from negligence or improper including but not limited to improper use or physical abuse such as dropping or otherwise.

This warranty shall be personal to the original purchaser. Any sale, rental or other transfer or use of the product covered by this warranty to or by a user other than the original purchaser shall cause this warranty to immediately terminate.

This remedies provided for in this warranty are the exclusive remedies available for any breach hereof. Neither Apex Medical nor its supplier or distributors shall be liable for any incidental, consequential, or special damage of any nature or kind caused by or arising out of defect in the product. All other warranties express or implied are excluded, including the warranties of merchantability and fitness for a particular purpose.

Wuxi Apex Medical Co., Ltd.

Add: Suite 601, 5Bld, No. 7 Changjiang Road, High Tech Area, Wuxi, Jiangsu, China 214028

Tel: +86-510-80100606

Fax: +86-510-80100505

www.apexdiabetes.com

Shanghai International Holding

Corp. GmbH (Europe)

Add: Eiffestrasse 80, 20537

Hamburg, Germany

Tel: +49-40-2513175

Representative in Russia

Diamama, LLC, address: 121471, Moscow, Ryabinovaya Str., 43, bldg. 2

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



4



6

Fill empty space in cannula with following
prime of U-100 Insulin:
6mm:0.2u(0.002mL) 9mm:0.3u(0.003mL).




8


10

Disconnecting

Reconnecting

Use

Device

The Easy Set infusion set is a device Easy Set infusion set is used for using in combination with the insulin pump for continual injection of insulin with the established and changing speed for the treatment of diabetes mellitus for insulin requiring patients.

The Easy Set infusion set is an infusion system, the proximal end of which is connecting to the reservoir with the medicinal product. Located at the distal end of the system of permanent catheter introduced subcutaneously. The motor drive to push the reservoir, so that insulin will flow into the infusion set. The insulin will be delivered to the human body subcutaneously through the infusion

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Contraindications

Infused infusion sets for single use are indicated for subcutaneous use only and not for intravenous infusion or the infusion of blood or blood products.

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Physical Properties

Infusion set: soft cannula-1pc, multiple
Insertion device-1pc
Packaging: plastic casing, dialyzing paper.

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Content of the set: 10pcs/box


Incompatible Pumps

The infusion set should be used with TruCare insulin pump (produced by Apex Medical), the Accu-Check D-TRON Plus, Roche Accu-Check Spirit and Roche Accu-Check Combo (produced by Roche).

Specification

Reference Number	Catheter Length	Total Tubing Length
ES206051	6mm	50cm
ES206101	6mm	100cm
ES209051	9mm	50cm
ES209101	9mm	100cm

Drawings

1-Connector; 2-Tubing; 3-Guide Needle; 4-Hub;
5-Protective Cap; 6-Separator Cover;
7-Protective Cover; 8-Soft Cannula

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If the infusion site becomes inflamed, replace the set, and use a new site until the first site has healed. Replace the infusion set if the tape becomes loose, or if the soft cannula becomes

dislodged from the skin.

Infusion set every 48-72 hours

Centers for Disease Control

or per your healthcare

provider's instructions.

Do not reinsert the introducer needle into the

infusion set. Re-insertion may cause tearing

of the soft cannula, which may result in

unpredictable medication flow.

Never prime the set or attempt to free a

clogged line while the set is inserted. You may

accidentally inject too much medication.

Do not put disinfectants, perfumes, or

deodorants on the infusion set as these may

affect the integrity of the set.

Dispose of the infusion set and introducer

needle safely, in a sharps container, after a

single use. Do not clean or re-sterilize.

This device is sterile and non-pyrogenic

unless the package has been opened or

damaged.

Do not use if the package has been opened

or damaged. (For TruCare users, do not use

the infusion set if the tubing connector needle

has been damaged.)

Store infusion sets in a cool, dry place. Do not

leave infusion sets in direct sunlight or inside

a vehicle.

Use aseptic techniques when temporarily

disconnecting the set and consult your

healthcare provider on how to compensate

for missed medication when disconnected.

When infusing insulin, and your blood glucose

level becomes unexplainably high, or an

occlusion alarm occurs, check for clogs

and/or leaks.

In doubt, change the infusion set because

the soft cannula may be dislodged, crimped

and/or partially clogged. Should any of these

problems arise, make a plan with your

healthcare professional for rapidly replacing

insulin. Test your blood glucose level to make

sure the problem is corrected.

When infusing insulin, do not change the infusion

set just before bedtime unless you can check

your blood glucose 1-3 hours after insertion.

Use of the infusion set may cause damage

to the cannula/needle and lead to infection, site irritation, and/or inaccurate medication delivery.

Transport

Transportation should be conducted in accordance with the contract which two sides agreed with.

Storage

The infusion set should be stored in the room with less than 80% relative humidity, non-corrosive gases and well-ventilated, and the products should be fully protected.

Operating Condition

Temperature: 5°C~40°C; Humidity: 20%~80%;
Barometric Pressure: 50kPa~106kPa

Validity

On the condition that comply with the requirements of the storage and transportation, sterilized reservoir packing is valid for two years.

Warranty

Wuxi Apex Medical Co., Ltd. warrants its infusion sets and reservoir against defects in materials and workmanship for the period up to 3 days after the individual reservoir/infusion set sterile packaging has been opened, not to exceed 6 months from date of shipment of the infusion set or reservoir to the end user (the "Warranty Period").

During the Warranty Period, Wuxi Apex Medical Co., Ltd. will replace any defective reservoir or infusion set, subject to the conditions and exclusions stated herein.

This warranty is valid only if the Wuxi Apex Medical Co., Ltd. Infusion sets and reservoir are used in accordance with the Instruction for Use. This warranty will not apply:

If the infusion set or reservoir has been used for more than a single-time use by a single end-user;

If damage results during the improper opening of the sterile package not in conformance with the procedures outlined in the Instruction for Use;

If damage results from changes or modifications made to the infusion set or reservoir by the user or third persons after the date of manufacture;

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...t of the infusion sets or
...son or entity other than the

...results from a Force Majeure or
... beyond the control of the
...rer; or

...ge results from negligence or improper
...including but not limited to improper
...age or physical abuse such as dropping or
...erwise.

...is warranty shall be personal to the original
...rchaser. Any sale, rental or other transfer or
...se of the product covered by this warranty to
...by a user other than the original purchaser
...all cause this warranty to immediately
...minate.

...is remedies provided for in this warranty are
...e exclusive remedies available for any
...each hereof. Neither Apex Medical nor its
...plier or distributors shall be liable for any
...cidental, consequential, or special damage of
...y nature or kind caused by or arising out of
...ect in the product. All other warranties
...press or implied are excluded, including the
...arranties of merchantability and fitness for a
...rticular purpose.

Wuxi Apex Medical Co., Ltd.

Address: Suite 601, 5Bld, No. 7 Changjiang
Road, High Tech Area, Wuxi, Jiangsu,
China 214028
Tel: +86-510-80100606
Fax: +86-510-80100505
www.apexdiabetes.com

EC REP

Shanghai International Holding
Corp. GmbH (Europe)
Address: Eiffestrasse 80, 20537
Hamburg, Germany
Tel: +49-40-2513175

Representative in Russia

Diamama, LLC, address: 121471, Moscow,
Kryabinovaya Str., 43, bldg. 2
Tel.: +7 (495) 978-78-17
E-mail: info@dia-ma.ru

Symbols

Conformité européenne (European Conformity).
This symbol means that the device fully complies
with European Directive 93/42/EEC.

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Reference Number

Lot Number

Caution: U. S. federal law restricts this device
to sale by or on the order of a physician.

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THIS IS TO CERTIFY: The seal of WUXI APEX MEDICAL CO., LTD on the annexed ANNEX TO INSTRUCTION FOR USE is genuine.

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Authorized
Signature : LI JUN

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(Date: SEP. 28, 2015)

ANNEX TO INSTRUCTION FOR USE(FOR RUSSIA)

Fast Set-I, Easy Set I Sterilized Infusion Set for Single Use

Fast Set-II, Easy Set-II Sterilized Infusion Set for Single Use

Fast Set-I Sterilized Reservoir for Single Use

Wuxi Apex Medical Co., Ltd.

ANNEX TO INSTRUCTION FOR USE (FOR RUSSIA)
Fast Set-I Sterilized Reservoir for Single Use

No.	Components		Material Construction
1	Needle		Stainless Steel 304
2	Barrel		Polypropylene
3	Separator	Plastic Part	Poly Carbonate
		Needle	Stainless Steel 304
		Plug	Silicon (Medical-Grade)
4	Injection Part	Plastic Part	Poly Carbonate
		Soft Cannula	Polytetrafluoroethylene
		Metal Needle	Stainless Steel 304
		Protection Cap	Polypropylene
5	Packaging	<u>Dialyzing Paper</u>	Tyvek
		Blister Box	PET
		Size	10.0cmx3.8cmx2.16cm

Additional specification	
Normal/standart injection speed of insulin	0.167 μL/min
The pre-defined min and max injection speed of insulin	0~15 μL/min
Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10)	Class B
Operating Condition:	Temperature: +5℃ ~ +40℃ Humidity:20%~80%; Barometric Pressure: 50kPa~106kPa

ANNEX TO INSTRUCTION FOR USE (FOR RUSSIA)
Fast Set-I, Easy Set-I Infusion Set

No.	Components		Material Construction
1	Tubing		Inner Layer: Poly Ethylene (HDPE, Medical-Grade) Outer Layer: Poly Urethane (Medical-Grade)
2	Connector		Poly Carbonate
3	Separator	Plastic Part	Poly Carbonate
		Needle	Stainless Steel 304
		Plug	Silicon (Medical-Grade)
4	Injection Part	Plastic Part	Poly Carbonate
		Metal Needle	Stainless Steel 304
		Skin Patch	Non-woven Fabrics
		Adhesive layer	Acrylic pressure-sensitive adhesive: Acrylic PSA 50% Ethyl Acetate 15% Isopropyl Alcohol 15.75% Heptane 16.5% Methylbenzene 2.5% 2,4-pentanedione 0.25%
		Protection Cap	Polypropylene
6	Packaging	Dialyzing Paper	Tyvek
		Blister Box	PET
		Size	10.5cmx7.35cmx1.96cm

Additional specification	
Normal/standart injection speed of insulin	0.167 μL/min
The pre-defined min and max injection speed of insulin	0~15 μL/min
Skin Patch	
Degree of adhesion for Skin Patch	≤30%
Linear dimensions for Skin Patch	φ=40mm
Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10)	Class B
Operating Condition:	Temperature: +5℃ ~ +40℃ Humidity:20%~80%; Barometric Pressure: 50kPa~106kPa

ANNEX TO INSTRUCTION FOR USE (FOR RUSSIA)
Fast Set-II, Easy Set-II Sterilized Infusion Set for Single Use

No.	Components		Material Construction
1	Tubing		Inner Layer: Poly Ethylene (HDPE, Medical-Grade) Outer Layer: Poly Urethane (Medical-Grade)
2	Connector		Poly Carbonate
3	Separator	Plastic Part	Poly Carbonate
		Needle	Stainless Steel 304
		Plug	Silicon (Medical-Grade)
4	Injection Part	Plastic Part	Poly Carbonate
		Soft Cannula	Polytetrafluoroethylene
		Metal Needle	Stainless Steel 304
		Skin Patch	Non-woven Fabrics
		Adhesive layer	Acrylic pressure-sensitive adhesive: Acrylic PSA 50% Ethyl Acetate 15% Isopropyl Alcohol 15.75% Heptane 16.5% Methylbenzene 2.5% 2,4-pentanedione 0.25%
		Protection Cap	Polypropylene
5	Packaging	Dialyzing Paper	Tyvek
		Blister Box	PET
		Size	10.5cmx7.35cmx3.76cm

Additional specification	
Normal/standart injection speed of insulin	0.167 μ L/min
The pre-defined min and max injection speed of insulin	0~15 μ L/min
Skin Patch	
Degree of adhesion for Skin Patch	$\leq 30\%$
Linear dimensions for Skin Patch	$\varnothing=40\text{mm}$
Class of hazard for medical wastes by reference to specification of morphological composition (in accordance with Russian Sanitary Regulations and Standards 2.1.7.2790-10)	Class B
Operating Condition:	Temperature: $+5^{\circ}\text{C} \sim +40^{\circ}\text{C}$ Humidity: $20\% \sim 80\%$; Barometric Pressure: $50\text{kPa} \sim 106\text{kPa}$

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The translation of this text was performed by me, the translator Mamedov Timur Djavanshirovich. The requirements regarding the text of the translation (maximum correctness, literate statement) have been explained to me.

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The city of Moscow.

05.10.2015

I, Akimov Gleb Borisovich, Notary of the city of Moscow, certify the authenticity of signature, made by the translator Timur Djavanshirovich in my presence. His identity is established.

It is registered in the register under No. 10 - 110889
Fee collected: 100 rubles

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[Handwritten signature]

Stitched up, numbered and sealed (13) pages

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