

"APPROVE"
Poly Medicure Limited
S. S. Rawat
(name)
Sr. GM-Quality
(position)


(signature)
Stamp

Instruction for medical device use

Infusion sets in different design variants

24 FEB 2020

ATTESTED

NOTARY PUBLIC
(INDIA)

2020

CONTENTS

1	Medical device name	3
2	Medical device manufacturer information	12
3	Medical device use	12
4	Indications for use	12
5	Contraindications	12
6	Precautions	12
7	Possible side effects	13
8	Technical information	13
9	Package contents	25
10	Safety requirements	33
11	Materials	34
12	Preparation for operation and operation procedure	38
13	Service life	38
14	Package, transportation and storage information	43
15	Maintenance	43
16	Cleaning and disinfection	43
17	Warranty statement	44

24 FEB 2020

ATTESTED

(Signature)

NOTARY PUBLIC
(INDIA)

This instruction for use of infusion sets is intended to familiarize you with them and their operation, and contains information about the application, operation and data on the sets.

1. Medical device name

Infusion systems in different design variants (hereafter – "set" or device)

Design variants:

1. System for infusion of intravenous fluids with Polyfusion air duct (polyfusion) as part of:

1.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

1.2 Dropper with an air duct, with a liquid filter with a cell of 15 microns and a length of 57 mm.

1.3 connecting Pipe, in the versions

- inner diameter 2.9 mm, length 150cm;

-inner diameter 2.9 mm, length 180cm;

- inner diameter 2.9 mm, length 200cm.

1.4 flow rate Roller clamp.

1.5 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);

- with Luer slip attachment complete with 21G injection needle (0. 8x40 mm);

· with Luer lock attachment.

1.6 Y-port (or no port).

1.7. Consumer packaging.

1.8 instructions for use.

2. System for the infusion of intravenous fluids without an air duct Polyfusion (polyfusion) in the composition:

2.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

2.2 Dropper without air duct, with a liquid filter with a cell of 15 microns and a length of 57 mm.

2.3 connecting Tube, in versions

- inner diameter 2.9 mm, length 150cm;

-inner diameter 2.9 mm, length 180cm;

- inner diameter 2.9 mm, length 200cm.

2.4 flow rate Roller clamp.

2.5 Injection unit, in versions:

with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);

- with Luer slip attachment complete with 21G injection needle (0. 8x40 mm);

- with Luer lock attachment.

2.6 Y-port (or no port).

2.7. Consumer packaging.

2.8 instructions for use.

3. System for infusion of intravenous fluids with Microfusion microfusion (Microfuge) in the composition:

3.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

3.2 Air filter with a cell of 0.3-3 microns.

3.3 micro-Droplet with a liquid filter with a cell of 15 microns and a length of 57 mm.

3.4 connecting Tube, inner diameter 2.9 mm, length 150cm.

3.5 flow rate Roller clip.

3.6 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);

- with Luer slip attachment complete with 21G injection needle (0. 8x40 mm);

- with Luer lock attachment.

- 3.7 Y-port (or no port).
- 3.8. Consumer packaging.
- 3.9 instructions for use.

4. System for intravenous fluid infusion with polytrol flow controller (Polytrol) in the composition:

- 4.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.
- 4.2 Dropper with a liquid filter with a cell of 15 microns and a length of 63 mm.
- 4.3 connecting Tube, inner diameter 2.9 mm, length 150cm.
- 4.4 flow rate Controller.
- 4.5 Injection unit, in versions:
 - with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
 - with Luer slip attachment complete with 21G injection needle (0. 8x40 mm);
 - with Luer lock attachment.
- 4.6 Y-port (or no port).
- 4.7. Consumer packaging.
- 4.8 instructions for use.

5. System for infusion of intravenous fluids with Novofusion air duct (Novofyuzhn) as part of:

- 5.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.
- 5.2 Dropper with an air duct, with a liquid filter with a cell of 15 microns and a length of 68 mm.
- 5.3 connecting Tube, inner diameter 2.9 mm, length 150cm;
- 5.4 flow rate Roller clamp.
- 5.5 Injection unit, in versions:
 - with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
 - with Luer slip attachment complete with 21G injection needle (0. 8x40 mm);
 - with Luer lock attachment.
- 5.6 Y-port (or no port).
- 5.7. Consumer packaging.
- 5.8 instructions for use.

6. System for infusion of intravenous fluids without an air duct Novofusion (Novofyuzhn) as part of:

- 6.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.
- 6.2 Dropper without air duct, with a liquid filter with a cell of 15 microns and a length of 68 mm.
- 6.3 connecting Tube, inner diameter 2.9 mm, length 150cm.
- 6.4 flow rate Roller clamp.
- 6.5 Injection unit, in versions:
 - with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
 - with Luer slip attachment complete with 21G injection needle (0. 8x40 mm);
 - with Luer lock attachment.
- 6.6 Y-port (or no port).
- 6.7. Consumer packaging.
- 6.8 instructions for use.

7. System for infusion of intravenous fluids with filter and air duct Oncofusion (oncofusion) as part of:

- 7.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.
- 7.2 Dropper with a liquid filter with a cell of 15 microns and a length of 67 mm.
- 7.3 connecting Tube, inner diameter 2.9 mm, length 150cm;

7.4 flow rate Roller clamp.

7.5 Filter with a cell size of 0.2 microns.

7.6 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

7.7 Y-port (or no port).

7.8. Consumer packaging.

7.9 instructions for use.

8. System for infusion of intravenous fluids with a three-way tap and polyfusion air duct (Polyfusion) as part of:

8.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

8.2 Dropper with an air duct with a liquid filter with a cell of 15 microns and a length of 57 mm.

8.3 connecting Pipe, in the versions

- inner diameter 2.9 mm, length 150cm;
- inner diameter 2.9 mm, length 180cm;
- inner diameter 2.9 mm, length 200cm.

8.4 flow rate Roller clamp.

8.5 three-Way crane.

8.6 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

8.7 Y-port (or no port).

8.8. Consumer packaging.

8.9 instructions for use.

9. System for infusion of intravenous fluids with a three-way tap without an air duct Polyfusion (polyfusion) in the composition:

9.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

9.2 Air filter with a cell of 0.3-3 microns.

9.3 Dropper without an air duct with a liquid filter with a cell of 15 microns and a length of 57 mm.

9.4 connecting Pipe, in the versions

- inner diameter 2.9 mm, length 150cm;
- inner diameter 2.9 mm, length 180cm;
- inner diameter 2.9 mm, length 200cm.

9.5 flow rate Roller clamp.

9.6 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

9.7 Y-port (or no port).

9.8. Consumer packaging.

9.9 instructions for use.

10. System for infusion of intravenous fluids with a double Neofusion dropper (Neofyuzhn) in the composition:

10.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

10.2 Double dropper with a liquid filter with a cell of 15 microns and a length of 52 mm.

10.3 Air filter with a cell of 0.3-3 microns.

10.4 connecting Tube, inner diameter 2.9 mm, length 150cm.

10.5 flow rate Roller clamp.

10.6 Filter with a cell size of 0.2 microns.

10.7 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

10.8 Y-port (or no port).

10.9. Consumer packaging.

10.10 instructions for use.

11. Spiked infusion system with Polyfusion-W (Polyfusion-B) wings without air duct in the composition:

11.1 Plastic spike diameter 0.5 mm, length 35 mm with a protective cap, equipped with small wings for additional convenience when piercing the cork.

11.2 Dropper without air duct, with a liquid filter with a cell of 15 microns and a length of 52 mm.

11.3 Air filter with a cell of 0.3-3 microns.

11.4 connecting Tube, inner diameter 2.9 mm, length 150cm.

11.5 flow rate Roller clamp.

11.6 Filter with a cell size of 0.2 microns.

11.7 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

11.8 Y-port (or no port).

11.9. Consumer packaging.

11.10 instructions for use.

12. Spiked infusion system with wings Polyfusion-W (Polyfusion-B) with an air duct in the composition:

12.1 Plastic spike diameter 0.5 mm, length 35 mm with a protective cap, equipped with small wings for additional convenience when piercing the cork.

12.2 Dropper with air duct, with a liquid filter with a cell of 15 microns and a length of 52 mm.

12.3 Air filter with a cell of 0.3-3 microns.

12.4 connecting Tube, inner diameter 2.9 mm, length 150cm.

12.5 Roller clip for flow rate.

12.6 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

12.7 Y-port (or no port).

12.8. Consumer packaging.

12.9 instructions for use.

13. Infusion system for photosensitive preparations Photofusion (Photofusions) in the composition:

13.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

13.2 Dropper with a liquid filter with a cell of 15 microns and a length of 57 mm.

13.3 Air filter with a cell of 0.3-3 microns.

13.4 connecting Tube, inner diameter 2.9 mm, length 150cm.

13.5 flow rate Roller clamp.

13.6 Adapter for the covered Luer tip with MLL cap.

13.7 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

13.8 Y-port (or no port).

13.9. Consumer packaging.

13.10 instructions for use.

14. Infusion system with a micro-dropper and a polytrol-Micro flow controller (Polytrol-Micro) in the composition:

14.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

14.2 Micro-dropper with a liquid filter with a cell of 15 microns and a length of 57 mm.

14.3 Air filter with a cell of 0.3-3 microns.

14.4 connecting Tube, inner diameter 2.9 mm, length 150cm.

14.5 flow rate Controller.

14.6 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

14.7 Y-port (or no port).

14.8. Consumer packaging.

14.9 instructions for use.

15. System for infusion with measurement of the volume of the injected solution of Polyvol (Polyvol) in the composition:

15.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

15.2 graduated Dropper for measuring the volume of the injected solution with a liquid filter with a cell of 15 microns and a length of 160 mm.

15.3 Air filter with a cell of 0.3-3 microns.

15.4 connecting Tube, inner diameter 2.9 mm, length 150cm.

15.5 flow rate Roller clamp.

15.6 Injection unit, in versions:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

15.7 Y-port (or no port).

15.8. Consumer packaging.

15.9 instructions for use.

16. Infusion system with an Autofusion air duct (Autofusion) in the composition:

16.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap.

16.2 Double dropper with a special filter that prevents air from entering the line when emptying a bottle with an air duct with a liquid filter of 15 microns and a length of 68 mm.

16.3 connecting Tube, inner diameter 2.9 mm, length 180cm.

16.4 flow rate Roller clamp.

16.5 the Injection site, in embodiments,:

- with Luer slip attachment complete with 18g injection needle (1, 2x40 mm);
- with Luer slip attachment complete with 21G injection needle (0.8x40 mm);
- with Luer lock attachment.

16.6 Y-port (or no port).

16.7. Consumer packaging.

16.8 instructions for use.

2. Medical device manufacturer information

Poly Medicure Limited

Manufacturing site: UNIT I: Plot No. 104 – 105, Sector – 59, HSIIDC Industrial Area, Ballabgarh – 121 004, Faridabad, Haryana, India.

3. Scope of application of the medical device

In medical centers and health care facilities.

3. Indications for use of the medical device

- Infusion of fluids, such as parenteral nutrition, as well as administration of other medications.
- Support for hydration and/or proper dehydration in patients who are unable to orally take sufficient volume of fluid.

Infusion sets (hereinafter referred to as set of device) **are designed for** gravitational injection of solutions for infusion to a patient from glass and polymer containers.

5. Contraindications to the use of the medical device

Administration of highly viscous fluids.

- Not to be used in patients with known hypersensitivity to any used material.

6. Precautions for use of the medical device

- Do not use for transfusion of blood or blood components.
- Do not use if protective caps are missing or loose.
- Do not use if the package is opened or damaged.
- Use the rubber part for the administration of medications only.
- Dispose of the set after a single use. For injection by gravity feed only.
- If the device contains phthalates, do not use it in children, as well as pregnant and lactating mothers.

Safety instructions

- The procedure is performed according to the attending physician's prescription.
- The sets should not be connected to other infusion or transfusion systems.
- The sets are intended for single use only. Re-use is prohibited.
- Rules of aseptic should be observed when preparing the sets for work and filling.
- The sets are delivered sterile. Resterilization of the sets is prohibited.
- The procedure should be performed only by qualified medical personnel.

7. Possible side effects during the medical device use


Air embolism resulting from: suction of air into open veins with low central venous pressure; faults in airtightness or inadequate filling of the system due to a violation of the infusion technique.

Prevention of side effects

- 1) Accurate compliance with infusion techniques is necessary.
- 2) Monitoring of the absence of air bubbles inside the system and the level of infusion medium in the container during transfusion/infusion, and especially at the end to eliminate air embolism.

8. Technical information

All variants of the Infusion Set contain a plastic spike, a "device" for perforating the stopper (or cap) of the container with a solution for infusion, with a protective cap (cap for a plastic spike). (hereafter – plastic spike or plastic needle).

*Winged spike- "spike with wings"


Infusion sets should be designed in a way to minimize the risk of contamination by particles. The surface of all parts of the set, contacting with the fluid, should be smooth and clean. The number of particles should not exceed the limit value of the pollution index of 90. The infusion system should not cause haemolysis. Infusion sets, with the exception of protective caps, should withstand a static tensile load of at least 20 N.

Infusion set design variant	Length, cm ($\pm 10\%$)	Weight, g ($\pm 10\%$)
Polyfusion IV infusion set with airvent	170.7	25.0
	200.7	26.0
	220.7	27.0
Polyfusion IV infusion set without airvent	170.2	24.0
	200.2	25.0
	220.2	26.0
Microfusion IV infusion set with a micro dropper	170.7	25.0
Polytrol IV infusion set with flow regulator	170.8	47.0
Novofusion IV infusion set with airvent	170.3	25.0
Novofusion IV infusion set without airvent	170.3	24.5
Oncofusion IV infusion set with filter and airvent	184.2	50.0
Polyfusion IV infusion set with three-way stop cock and airvent	173.7	28.0
Polyfusion IV infusion set with three-way stop cock without airvent	173.7	27.5
Neofusion IV infusion set with double dropper	170.7	24.0
Polyfusion-W infusion set with winged spike without airvent	170.7	25.0
	200.7	26.0
	220.7	27.0
Polyfusion-W infusion set with winged spike with airvent	170.7	25.7
	200.7	26.7

	220.7	27.7
Photofusion infusion set for photosensitive drugs	170.7	24.0
Polytrol-Micro infusion set with micro dropper and flow regulator	170.2	45.0
Polyvol measured volume fluid administration set	180.5	40.0
Autofusion infusion set with airvent	175.0	35.0

Protective cap (protective cap for plastic spike)

- Sufficient length to accommodate a plastic spike.
- The material should be transparent or translucent.
- The engagement should not be tight so that the cover can be removed manually, but reliable enough to prevent accidental disconnection during handling and transport.
- The surface should be clean, free from burrs, black particles, mechanical damage, permanent inclusions, etc.
- The cap of the plastic "winged needle" has the wings (which are called "spike with wings", see the picture above) for easy piercing of the stopper of the container with solution for injection.
- Diameter 13 ± 0.1 mm
- The strength of the connection with a plastic spike is not less than 50N.
- The cap of the plastic spike (plastic needle) can be equipped with an airvent with an air filter and a plug, depending on the design variant of the Set.

Airvent


- The airvent and the plastic spike should have a leak-proof connection. The tightness of the connection should be checked at a pressure of 0.3 MPa (3 bar).
- The airvent should be free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- The air filter should provide a bolus feed from the bottle and the fluid into the bottle.
- When connecting the airvent to a rigid container for infusion, air from the container must not enter the flow.

Plastic spike

- The plastic spike should meet the size requirements.
- The tip of the spike should be sharp and be able to easily pierce the cap of the glass and polymer containers of solutions for infusion. Sharp and easily piercing needle should be emptied without residue. It has two holes located one after the other on different sides of the needle at a distance of 9 mm from the tip of the needle.
- The adhesive bond between the spike and the drip chamber should be leak-proof.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- The length of plastic needle with protective cap is $35 \text{ mm} \pm 0.3 \text{ mm}$.
- Diameter of the spike is $0.5 \text{ mm} \pm 0.05 \text{ mm}$.
- The strength of the connection of the plastic spike with the drip chamber should be at least 20N.


Tube (with DEHP/DEHP free, as well as PVC free, as appropriate).

- The tube should be transparent or translucent in order to see the movement of air bubbles.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- There should be no bends or stuck together parts of the tubes that could prevent the normal flow of fluid.
- The place of connection of the tube with the drip chamber (dropper) should be leak-proof and withstand a static tensile force of at least 20 N.
- Amber PVC tubes with UV protection with a wavelength of 290-450 nm for photosensitive medications.
- Inner diameter is $2.9 \text{ mm} \pm 0.1 \text{ mm}$; outer - $4.1 \text{ mm} \pm 0.1 \text{ mm}$.
- The hardness of the tube should be from 75 A to 77 A according to Shore.
- The length of the tube is 100, 110, 150, 180 and 200 cm $\pm 0.05 \text{ cm}$.

Roller clamp for flow rate adjustment (roller and roller body)

- The roller clamp for flow rate adjustment should be able to regulate the infusion fluid flow rate from zero to maximum (from full closure to bolus injection).
- The roller clamp for flow rate adjustment should not cause damage to the tube during prolonged use.
- During storage, there should be no chemical interaction between the roller clamp for flow rate adjustment and the tube.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- The length of the poller clamp is $50.0 \pm 5 \text{ mm}$.

Gate location

Surface should be clean

There should be no burrs on the surface

ation
ation

Corrugation should not be damaged.

Corrugation should not be damaged

Flow regulator

- The flow regulator should be able to regulate the infusion fluid flow rate from full closure to bolus injection.
- The flow regulator should not cause damage to the tube during prolonged use.
- During storage, there should be no chemical interaction between the flow regulator and the tube.
- There should be no burrs, black particles, mechanical damage, etc.
- The flow regulator should be able to regulate the solution flow rate from 5 to 250 ml (maximum) per hour.
- There is no need for frequent flow rate adjustment.
- protection from short interruption of the infusion.
- Control of infusion regardless of the properties of the line.
- Accurate reproduction of the parameters when resuming the infusion.
- Length 35 ± 0.1 mm
- Diameter 2.5 ± 0.1 mm

"Luer lock" and connection "Luer slip"

- The taper and thread profile must meet the requirements of EN 20594-1:1993/ AC-1996/ISO 594-1:1986 and ISO 80369-7:2016.
- The male and female components should have a 6% taper and tight fit with other devices according to EN 20594-1:1993/ AC-1996/ISO 594-1:1986 and ISO 80369-7:2016.
- There should be no burrs, black particles, mechanical damage, etc.
- A closed system should be created, which should be leak-proof when the infusion stops.
- Prevention of the risk of hospital infection.
- Prevention of reverse movement of fluid into the infusion set.
- Suitable for injection needles or catheters.
- "Luer lock" is additionally equipped with a female Luer tip, allowing it to be more tightly fixed on a catheter or needle. Fixation is achieved by rotating the female Luer tip.
- Length 15 ± 0.1 mm (for the "Luer lock" and connection "Luer slip"), 20 ± 0.1 mm (for rotating "Luer lock")

Luer
lock

Luer
slip

Rotating
Luer lock

*Luer lock- "Luer lock"

Luer slip-"Luer slip"

Rotating Luer lock- Rotating "Luer lock"

Injection site

The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

- Tight connection with the head with Luer lock taper and connecting tube.
- Allows bolus injection.
- length of the injection site 42 ± 0.1 mm
- thickness of the injection site 2.2 ± 0.01 mm
- the injection site should be located near the head with a "Luer lock" taper. The injection site should provide self-tightening upon piercing, at the site of the puncture should not be more than one drop of fluid.

Y-shaped injection port


- One end of the injection port is connected to the tube and the other to the male taper. Both connections should withstand the specified static tensile force.
- After perforation (with a needle with an outer diameter of 0.8 mm), the injection port should remain leak-proof, one drop of water is allowed at an air pressure of 50 kPa.
- The material should be transparent or translucent.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- Length 35 ± 0.1 mm

Adapter for a Luer Lock tip

- It should be transparent.
- The male and female components should have a 6% taper and leak-proof fit with other devices according to EN 20594-1:1993/AC-1996 and ISO 80369-7:2016.

Injection needle

- The needle should be made of 304 stainless steel.
- The needle hub should have a 6% female taper to ensure a leak-proof connection with a 6% male end according to BS EN 20594-1:1994 and 594-2:1998.
- The needle is bonded to the needle hub.
- During a visual inspection with 2.5x magnification the tip of the needle looks sharp, without molding edge, burrs and hooks.
- During a visual inspection with 10x magnification the rear cut looks uniform and fully polished on both sides, forming an angle p .

*Section


- The maximum values of acidity and alkalinity, hardness, resistance to destruction and corrosion meet the requirements of the ISO 9626. These tests are carried out at the place of operation and the protocol should be affixed by the manufacturer.
- Needle extract should contain not more than 5 mg/l of total lead, tin, zinc and iron in a microanalytic test, for example, atomic absorption. The cadmium content in the needle extract should not exceed 0.1 mg/l.
- The angle of the lancet bevel has to be clearly defined.
- The needle is sandblasted on all sides and the grit corresponds to the provided example.
- The outer/inner surface should be smooth, elevations, burrs, cracks and other surface damage are not allowed.
- The needle tube should be manufactured according to ISO 9626. If the needle tube is coated with lubricant, the amount of lubricant should be minimal and should not form liquid droplets on the outer and inner surface of the needle during visual inspection. During a visual inspection with 2.5x magnification the tip of the needle should be sharp, without burrs.
- the length of the tube should be 40 mm.

Note — The needle sharpening angle α usually has a cut at an angle of $11 \pm 2^\circ$, but a "short" cut at an angle of $17 \pm 2^\circ$ can also be used.

- The lumen patency should be such that:

- a) a stainless steel mandrel with a diameter according to Table 3 of ISO 7864 passed through the needle, or that
- b) the water flow rate through the needle was not less than 80% of the water flow rate through the tube of the needle with a minimum diameter of the canal according to ISO 9626 and sizes

according to Table 3 of ISO 7864. Water pressure should not exceed 1×10^5 Pa.

1 - needle head; 2 - connective part; 3 - needle tube; 4 - protective cap; l - length of the needle


Injection needle	
Size	Nominal outer diameter, mm
18G	1.2
21G	0.8

Notes:

- 1 The size of the injection needle supplied with the set should correspond to one of the values provided in the table above,
- 2 The length of the cap of the injection needle is $X \pm X$ mm; the inner/outer diameter of the cap of the injection needle is $(X/X) \pm x$ mm.
- 3 The sharpening angle of the needle tip is $XX \setminus X$ mm

Needle cap

- The needle cap should be conical or cylindrical.
- The length should be sufficient to accommodate the male end and the needle.
- It should be translucent or transparent.
- The needle cap should be easily removed, but not disconnected during handling, transportation and storage.

Air filter

- The air filter should be designed so that all incoming air would pass through it.
- It should not significantly reduce the fluid flow rate.
- The pore size should be 0.3-3 micron to prevent the entry of microorganisms.
- The air filter should be equipped with an air valve, which should prevent microbial contamination.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- The design of the air filter should ensure that the air enters the rigid container in such a way that the volumetric flow rate of the liquid is reduced by no more than 20% compared to the volumetric flow rate from the freely ventilated container.

Fluid filter

- It should be located as close as possible to the injection port.
- The porosity of the filter shall be not less than 15 micron.
- Engagement of the lower part of the cap with the spike.
- The engagement should not be too rigid, so it would be possible to remove the cap manually, but sufficiently reliable so that there would be no accidental disconnection during handling, transportation and storage.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.
- Filter for fluids with automatic infusion stop after emptying the drip chamber to prevent air embolism.

0.2 micron fine filter

- Designed for filtration of solutions for injections in order to prevent contamination, bacteria and solid impurities.
- Pore size is 0.2 micron with high throughput.

- Pressure resistance is not less than 0.2 MPa (2 bar).
- Made of polyethylene and polyvinylidene fluoride with nylon filter.
- It should be transparent or translucent.
- Compatible with chemotherapy drugs.
- Control of the absence of air regardless of the location of the filter.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

Adapter for male Luer tip with MLL cap (hereafter – adapter)

- The connector has an injection port and both ends are tightly connected to the feed tube.
- The connector should be transparent and free from any damages or foreign particles.
- The connector should be free from burrs, black particles, etc., The connection of the connector with the lower cap should be leak-proof.


Return valve with a polyethylene cap with filter

- It prevents reverse flow during parallel diffusion.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.

The filter should conform to the characteristics of the fluid filters. It should not leak in the closed position. The connection should be leak-proof.


Locking clamp (hereafter – "clamp")

- The clamp is designed to shut off the supply of liquid from the reservoir to the graduated container. In the closed position the clamp should not allow liquid to pass through.
- The surface should be clean, free from burrs, black particles, mechanical damage, foreign inclusions, etc.


1. The divergence should not exceed 0.2 mm.

2. The divergence should not exceed 0.2 mm.

Lumen should be free

Dropper (drip chamber)

- Clean and flexible.
- Drip chamber should provide the possibility of continuous control of the drop formation process. The liquid should enter the drip chamber through the drop-forming element. The distance

between the end of the drop-forming element and the outlet of the drip chamber should be not less than 40 mm, the distance between the end of the drop-forming element and the fluid filter should be not less than 20 mm. The distance between the wall of the drip chamber and the wall of the drop-forming element should be not less than 5 mm.

- The design of the drop-forming element should ensure the formation of 20 or 60 drops of distilled water from (1 ± 0.1) ml or (1 ± 0.1) g of distilled water at ambient temperature $(23 \pm 2)^\circ\text{C}$ at a volumetric flow rate (50 ± 10) drops per minute. The drip chamber should provide and facilitate the process of primary filling of the infusion line with fluid. (according to GOST 57629-2017 for all design variants of the set, except for the Polyvol measured volume fluid administration set)

- The design of the drop-forming element should provide the formation of 20 drops of (1.0 ± 0.1) g of distilled water at the temperature $(20 \pm 2)^\circ\text{C}$ and the flow rate (50 ± 5) drops per minute (for the design variant: Polyvol measured volume fluid administration set according to GOST 25047-87).

- The fluid filter at the lower part of the drip chamber should meet the specifications described above.


- It has a leak-proof screw connection with a connecting tube.

- Depending on the Set, drip chambers of different sizes are used: drip chamber with a butterfly needle (with a winged spike), "Nova" drip chamber, double drip chamber, drip chamber for measuring the volume of the injected solution, drip chamber with a fluid filter with automatic infusion stop.

- The length of the dropper should be: 57 mm, 68 mm or 63 mm (for the "Novo" dropper), 57 mm (for the micro dropper), 52 mm (for the dropper with butterfly needle), 52 mm and 68 mm (for the dropper with fluid filter and with special filter, which prevents the entry of air in the line when the IV bottle is empty), 160 mm (for the dropper for measured volume fluid administration).

- The dropper for measured volume fluid administration has a clear indelible scale and a silicone float for accurate indication of the volume of solution for infusion. The increment is 1 ml, graduation from 0 to 150 ml.

- The double drip chamber has an elastic lower part (3) and a large separated distributor (1), which provides the formation of droplets connected to each other by a separator (2).

- For dropper with autostop (with a special filter that prevents the entry of air in the line when the IV bottle is empty): in the lower part (3) of the double dropper there is a fluid filter(4) with automatic stop after completion of the infusion. The floating filter (4) automatically stops the infusion after emptying the measuring vessel. Autostop is a unique airtight filter membrane, which not only retains solid impurities, but also prevents the penetration of air. When the IV bottle is emptying, autostop maintains a constant fluid level. The risk of air embolism during infusion is reliably foreseen.

*AIR- air


- Micro dropper has a unit for the formation of microdroplets for a longer flow of fluid for injection.

Tree-way stop cock

- The stop cock should meet the requirements of ISO 594-1:1986 and be free from damage, black particles, burrs or stains.
- The heparin cap profile – the taper and thread should meet the requirements of EN 20594-1:1993/AC-1996/ISO 594-1:1986 and ISO 80369-7:2016.
- "Heparin" cap (heparin is not part of the cap, the name "heparin cap" is used only because heparin solution can be injected through the cap)
- The cap cover and the injection plug connection should be leak-proof. The tightness of the connection is checked at a pressure of 0.3 MPa (3 bar).
- After piercing the injection plug with a needle 0.6 mm in diameter and staying in this position for 15 seconds, there should be no bubbles on it when tested using an air pressure of 0.02 MPa (0.2 bar) under water (after pulling out the needle).

The opening should be centered

- 360° rotation in axial and radial direction.
- Colored rotating blade with arrows indicating the flow direction.
- Prevention of cold deformation in order to avoid pressure on the injection cap.
- The presence of white spots, cracks and other defects on the deforming area is not allowed.
- The surface should be smooth, presence of black particles, signs of lubricant, dirt or impurities is not allowed.
- The component should be transparent and bright. The yellowness, or rust are not allowed.
- Resistance to the dangerous effects of the medications.
- No leakage even after prolonged use.
- Prevention of pollution and air embolism.
 - Reliable access for injection and aspiration.
 - Outer diameter 2.0 ± 0.2 mm
 - Inner diameter 1.0 ± 0.2 mm
 - The strength of the connection is not less than 20N.

9. Package contents

The package contains:

1. Polyfusion IV infusion set with airvent, consisting of:

1.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

1.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

1.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

1.4 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

1.5 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

1.6 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

1.7 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

1.8 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; injection site with Luer attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

1.9 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube inner diameter 2.9 mm, length 200 cm; roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

Polyfusion IV infusion set without airvent, consisting of:

2.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

2.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

2.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

2.4 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with

injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

2.5 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

2.6 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

2.7 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

2.8 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

2.9 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

3. Microfusion IV infusion set with a micro dropper, consisting of:

3.1 Plastic spike 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; micro dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

3.2 Plastic spike 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; micro dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

3.3 Plastic spike 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; micro dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

4. Polytrol IV infusion set with flow regulator, consisting of:

4.1 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (with 15 micron fluid filter 63 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); flow regulator; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

4.2 Plastic spike 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper (with 15 micron fluid filter 63 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); flow regulator; injection site with Luer slip attachment completed with injection needle 21G (0.8x40

mm); Y-site (or without site); consumer packaging; instructions for use.

4.3 Plastic spike 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper (with 15 micron fluid filter 63 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

5. Novofusion IV infusion set with airvent, consisting of:

5.1 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment (or Luer lock or rotating Luer lock) completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

5.2 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment (or Luer lock or rotating Luer lock) completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

5.3 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

6. Novofusion IV infusion set without airvent, consisting of:

6.1 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

6.2 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

6.3 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

7. Oncofusion IV infusion set with filter and airvent, consisting of:

7.1 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (with 15 micron fluid filter 67 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

7.2 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (with 15 micron fluid filter 67 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

7.3 Plastic spike 0.5 mm, length 35 mm with protective cap; dropper (with 15 micron fluid filter 67 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

8. Polyfusion IV infusion set with three-way stop cock and airvent, consisting of:

8.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

8.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

8.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

8.4 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

8.5 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

8.6 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

8.7 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

8.8 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

8.9 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (with airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

9. Polyfusion IV infusion set with three-way stop cock without airvent, consisting of:

9.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

9.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with

15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

9.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

9.4 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

9.5 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

9.6 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

9.7 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

9.8 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

9.9 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; dropper (without airvent, with 5 micron fluid filter 57 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 200 cm); roller clamp for flow rate adjustment; three-way stop cock; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

10. Neofusion IV infusion set with double dropper, consisting of:

10.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; double dropper (with 15 micron fluid filter 52 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer slip attachment with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

10.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; double dropper (with 15 micron fluid filter 52 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer slip attachment with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

10.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; double dropper (with 15 micron fluid filter 52 mm in length); 0.3-3 micron air filter; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with

Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

11. Polyfusion-W infusion set with winged spike without airvent, consisting of:

11.1 Plastic spike 0.5 mm, length 35 mm with protective cap equipped with small wings for additional convenience in piercing the stopper; 0.3-3 micron air filter; dropper (without airvent, with 15 micron fluid filter 52 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

11.2 Plastic spike 0.5 mm, length 35 mm with protective cap equipped with small wings for additional convenience in piercing the stopper; 0.3-3 micron air filter; dropper (without airvent, with 15 micron fluid filter 52 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

11.3 Plastic spike 0.5 mm, length 35 mm with protective cap equipped with small wings for additional convenience in piercing the stopper; 0.3-3 micron air filter; dropper (without airvent, with 15 micron fluid filter 52 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; 0.2 micron filter; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

12. Polyfusion-W infusion set with winged spike with airvent, consisting of:

12.1 Plastic spike 0.5 mm, length 35 mm with protective cap equipped with small wings for additional convenience in piercing the stopper; dropper (with airvent, with 15 micron fluid filter 52 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

12.2 Plastic spike 0.5 mm, length 35 mm with protective cap equipped with small wings for additional convenience in piercing the stopper; dropper (with airvent, with 15 micron fluid filter 52 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

2.3 Plastic spike 0.5 mm, length 35 mm with protective cap equipped with small wings for additional convenience in piercing the stopper; dropper (with airvent, with 15 micron fluid filter 52 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

13. Photofusion infusion set for photosensitive drugs, consisting of:

13.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); adapter for Luer tip with MLL cap; consumer packaging; instructions for use.

13.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); adapter for Luer tip with MLL cap; consumer packaging; instructions for use.

13.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); adapter for Luer tip with MLL cap; consumer packaging; instructions for use.

14. Polytrol-Micro infusion set with micro dropper and flow regulator, consisting of:

14.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; micro dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); flow regulator; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

14.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; micro dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); flow regulator; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

14.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; micro dropper (with 15 micron fluid filter 57 mm in length); connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

15. Polyvol measured volume fluid administration set, consisting of:

15.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper with marking scale for measuring the volume of administered solution (with 15 micron fluid filter 160 mm in length); silicon float; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); locking clamp; consumer packaging; instructions for use.

15.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper with marking scale for measuring the volume of administered solution (with 15 micron fluid filter 160 mm in length); silicon float; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); locking clamp; consumer packaging; instructions for use.

15.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; 0.3-3 micron air filter; dropper with marking scale for measuring the volume of administered solution (with 15 micron fluid filter 160 mm in length); silicon float; connecting tube (inner diameter 2.9 mm, length 150 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); locking clamp; consumer packaging; instructions for use.

16. Autofusion infusion set with airvent, consisting of:

16.1 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; double dropper with special filter, which prevents the entry of air in the line when the IV bottle is empty (with airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 21G (0.8x40 mm); Y-site (or without site); consumer packaging; instructions for use.

16.2 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; double dropper with special filter, which prevents the entry of air in the line when the IV bottle is empty (with airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer slip attachment completed with injection needle 18G (1.2x40 mm); Y-site (or without site); consumer packaging; instructions for use.

16.3 Plastic spike diameter 0.5 mm, length 35 mm with protective cap; double dropper with special filter, which prevents the entry of air in the line when the IV bottle is empty (with airvent, with 15 micron fluid filter 68 mm in length); connecting tube (inner diameter 2.9 mm, length 180 cm); roller clamp for flow rate adjustment; injection site with Luer lock or rotating Luer lock; Y-site (or without site); consumer packaging; instructions for use.

10. Safety requirements

- Requirements for the sanitary protection of air, soil and water bodies should comply with the current sanitary and epidemiological rules and hygienic standards established in the country.
- The possibility of environmental pollution by production waste, generated in the process of production of the sets, should be excluded.
- The technology of production of sets does not involve the formation of liquid waste.
- Disposal of technological losses arising in the process of production of sets, as well as finished products, if it does not meet the requirements, should be carried out in accordance with applicable standards and requirements.
- The operation of the sets does not require special measures to protect the environment from harmful effects.

Disposal:

- ☐ The used devices are classified as class B (epidemiologically unsafe),
- ☐ Sets that have not been used for their intended purpose due to expiration and other reasons are classified as class A epidemiologically safe waste and disposed of as household waste.

11. Materials used in production

All materials from which the components are made meet the specifications and are allowed for use. The device and packaging do not contain components of animal origin.

No.	Item for manufacture of which the material is used	Basic material	Grade	Manufacturer
1	Protective cap air inlet plug	Polypropylene	SM498	HYOSUNG Corporation
2	Airvent filter (Air inlet filter)	Fiberglass	-----	Performance Fibers, Inc.
3	Latex injection site	Latex Natural rubber		TOP GLOVE SDN. BHD.
4	Screw cap	Polycarbonate	1201-15	LG Polycarbonate
5	Plastic spike	Acrylonitrile-butadiene-styrene	IFB700/HF400	TOYOLAC General Purpose
6	Small airtight cap	Low density polyethylene (LDPE) + High density polyethylene (HDPE)	M200056 / LDF 201 FG / HI1600+24FS040	SABIC Indothene-HD
7	Protective cap "winged spike" (large)	LDPE+HDPE	M200056/LDF201 FG / HI1600+24FS040	SABIC Indothene-HD
8	Protective cap "winged spike"	PVC (DEHP free)	IP87ND	Indel polymer
9	Roller clamp	Polypropylene	H 110MA / SM 498	HYOSUNG

No.	Item for manufacture of which the material is used	Basic material	Grade	Manufacturer
	housing			Corporation
10	Flow regulator	Acrylonitrile-butadiene-styrene	IFB920/HF400	TOYOLAC General Purpose
11	Roller clamp wheel	Polypropylene	H 110MA / SM 498	HYOSUNG Corporation
12	Nova drip chamber	Polyvinylchloride compound	PDC 337 K	Indel polymer
13	Dropper for measuring the volume of injected solution (DEHP free)	Polyvinylchloride compound	IP 92 ND	Indel polymer
14	Double drip chamber	Polyvinylchloride compound	PDC 337 K	Indel polymer
15	Micro dropper	Polyvinylchloride compound	PDC 337 K	Indel polymer
16	Male Luer tip	Granules of unplasticized polyvinylchloride (UPVC)	MMP20/B	SABIC
17	MLL cap	LDPE	LDF201FG	SABIC
18	0.2 micron fine filter	Polypropylene	-----	HYOSUNG Corporation
19	Wash Filter	Polyethylene (PE) + Polyvinylidene fluoride (PVDF)	-----	Indothene-HD
20	Y-site (small)	Acrylonitrile-butadiene-styrene (ABS) and Polyisoprene	-----	Styrolution America LLC
21	Y-site (large)	ABS and Polyisoprene	-----	Styrolution America LLC
22	Rotating Luer lock (Housing)	Acrylonitrile-butadiene-styrene	920 555 / TR558A	Styrolution America LLC
23	Luer tip (fixed)	ABS	920 555 / TR558A	Styrolution America LLC
24	Adapter	Polypropylene and UPVC	SM498/MMP 20B	HYOSUNG Corporation Indel polymer
25	Plastic spike	Natural acrylonitrile-butadiene-styrene (ABS)	IFB 920 / HF400	Styrolution America LLC
26	Air inlet stopper / valve	70% HDPE+30% LDPE	50MA180 + 24FS040/ M200056	Indothene-HD SABIC
27	Drip chamber (PVC free)	TPE	-----	SABIC

No.	Item for manufacture of which the material is used	Basic material	Grade	Manufacturer
28	Tube (DEHP free)	Polyvinylchloride compound	IP 31 ND/PDC 320 ND	Indel polymer
29	Tube	Polyvinylchloride compound	IP 31 /BT 31	Indel polymer
30	Tube (DEHP free)	TPE	Not applicable	
31	Roller clamp	Polypropylene	H 110MA/SM498	HYOSUNG Corporation Indel polymer
32	Adapter for Luer tip	Granules of unplasticized polyvinylchloride (UPVC)	MMP 20 B	HYOSUNG Corporation Indel polymer
33	Connecting taper (Luer slip connection)	Polypropylene	SM 498 / 575P	HYOSUNG Corporation Indel polymer
34	Needle hub	Polypropylene	SM 498 / 575 P	HYOSUNG Corporation Indel polymer
35	Protective needle cap	Polypropylene	SM 498 / 575 P	HYOSUNG Corporation Indel polymer
36	Needle	Stainless steel silicone coating	SS 304 dimeticone NF	Henan hzz iron and steel Co., Ltd
37	Fluid filter	Nylon	Nylon 6	Geo-Tech Polymers
40	Y-shaped connector	Clear ABS	PA - 758	Styrolution America LLC
39	Head with Luer lock taper (housing)	Polycarbonate	PC 1201-15/RX-1805	LG Polycarbonate
40	Handle	HDPE	50MA180/M200056	Indothene-HD
41	Luer lock tip	HDPE	M200056	Indothene-HD
42	Injection site	latex free	Polyisoprene	TOP GLOVE SDN. BHD.
43	Pigment	blue pigment light-blue pigment	Pigment No. Blue 15:0: CAS No,147-14-8 (Copper phthalocyanine) Titanium oxide(TiO2): CAS NO,13463-67-7 Pigment No. Blue 15:5: CAS No,147-14-8 (Copper phthalocyanine) Titanium oxide(TiO2)	Focus Pigment Sdn. Bhd.

No.	Item for manufacture of which the material is used	Basic material	Grade	Manufacturer
		yellow pigment	:CAS NO,13463-67-7 Pigment No. Yellow 14: CAS NO.5468-75-7 (Benzidine yellow) Titanium oxide(TiO2) :CAS NO,13463-67-7	
44	Paper for individual packaging of the devices	Packaging paper for medical devices	OCPSxxxxlOx/ OEHPxxxxlOx	AMCOR Flexibles
45	Film for individual packaging of the devices	PE + PP	For medical devices. Non-toxic	AMCOR Flexibles
46	Two-layer box	Cardboard	-----	AMCOR Flexibles
47	Transport package	Semi-craft paper	-----	AMCOR Flexibles
	EtO	20% EtO + 80% CO2	-----	

12. Preparation for operation and operation procedure

ATTENTION: Strictly adhere to the rules of aseptics and follow the general precautions and procedures established by the protocols of medical institutions

12.1 The procedure is carried out in accordance with the guidelines for medical professionals approved by the medical institution, and should include the following:

12.1.1 Preparation of hands for the procedure. Wash your hands according to the algorithm specified in the guideline. This prevents infection.

12.1.2 Prepare sterile tray, cloths, a device for opening metal packaging containers with the solution for infusion, the infusion set, stand, and medications, allowing to perform the procedure.

12.1.3 Check the applicability of the solutions, transparency, shelf life, color, name and concentration.

12.1.4 Prepare your hands for work.

12.1.5 Wipe the metal packaging of the container of medication for infusion to the patient with a cotton ball soaked with alcohol to prevent microbial contamination.

12.1.6 Open the central part of the metal packaging of the container with solution for infusion with a special device.

12.1.7 Wipe the stopper with a cotton ball soaked with alcohol.

12.1.8 Take the infusion set, check the integrity of the package and shelf life. Wipe the place for opening the package with a cotton ball soaked with alcohol and open the package with a special device (scissors).

12.1.9 Place a sterile disposable set on the sterile surface of the tray with sterile forceps.

12.1.10 Take the airvent (if included in the set) with a sterile tweezers, remove the cap from the plastic spike and fully insert the spike into the container plug, secure the airvent with an elastic band along the vial.

12.1.11 Remove the cap from the spike, located closer to the dropper and insert the spike into the container plug.

12.1.12 Close the clamp on the dropper.

12.1.13 Turn the container upside down and secure it on the rack.

12.1.14 Remove the injection needle with sterile forceps together with the cap and place it on the sterile surface of the work tray.

12.1.15 Position the dropper horizontally, open the clamp and fill the dropper halfway.

12.1.16 Close the clamp, position the dropper vertically.

12.1.17 Open the clamp, fill the system until the air is completely removed by directing a flow of the solution to the medical waste tray on the working table.

12.1.18 Close the clamp, connect the injection needle with the cap and secure the system on the rack.

12.1.19 Check the system for air bubbles and the needle for patency.

1.2.2 The attending physician selects the design variant for infusion depending on the infusion therapy (the amount of injected fluid, the therapy program, the conditions of administration, control of administration, etc.).

12. Service life (shelf life) and an indication of the prohibition of use of medical devices after the expiration date (shelf life)

The guaranteed service life of the devices: 5 years from the sterilization date. Do not use the system after the expiration date.

13. Storage and transportation conditions. Packaging.

The devices can be transported by all types of metered vehicles in accordance with the transportation rules,

applicable to this type of transport, while maintaining the temperature.

The transport temperature is between 0 and 40 °C.

Store in a dry place protected from direct sunlight at a temperature from 5 to 40 °C.

Each consumer package should have a label with indication of:

- ☐ manufacturer's name and its location (legal address, including country, and address of manufacture, if it differs from the legal address);
- ☐ manufacturer's trademark (if available);
- ☐ name (designation) of the device;
- ☐ package contents;
- ☐ inscription "For single use";
- ☐ inscription "Sterile", highlighted or printed in a different color;
- ☐ inscriptions "Non-pyrogenic", "Non-toxic";
- ☐ inscription "Do not use if the consumer packaging is damaged";
- ☐ sterilization date (month and year);
- ☐ lot number;
- ☐ inscription "Use before..." (month and year);
- ☐ warning about the need to check the integrity of the package before use;
- ☐ sterilization method;
- ☐ storage conditions;
- ☐ disposal conditions;
- ☐ symbol "presence of instructions";
- ☐ symbol about "the need to read the instructions";
- ☐ mark of conformity.
- ☐ information on the number of drops (20 or 60 drops) of distilled water equivalent to the volume (1±0,1) ml or weight (1±0,1) g, in the form of an inscription or a corresponding graphic symbol

according to GOST R ISO 15223-1;

- For needles, the marking should include: needle size as well as needle tube length(mm) and diameter(mm) (e.g. 18G (1.2x40mm), for needle 21G (0.8x40mm))

Note - the device label should be affixed to the external side of the consumer packaging (the packaged devices should not come into contact with the ink).


Transport marking:

The marking should be made typographically, lithographically or by printing on a paper label firmly glued to the package, indicating:

- product name;
- manufacturer's name and location (legal address, including country, and, in case it differs from the legal address, address(es) of production(s));
- production date (month, year);
- guaranteed validity period;
- lot number or serial number (if available);
- bar code (if available);
- storage and transportation conditions;
- ┘ marketing authorization number and date of issue;
- mark of conformity (if available)

The package should contain the transport label and manipulation signs with the meaning "Protect from moisture", "Fragile.Carefully.»

Description of symbols on the medical device

	See Instructions for use
	Catalogue number
	For single use only
LOT /Ch.B.	Lot number
	Date of manufacture
	Expiration date
	Sterilized by ethylene oxide

Sterilization and packaging

Medical devices in different design variants are packaged in individual sterile packaging that meets the requirements of ISO 11607-1, and are packed in a two-layer secondary packaging (individual packaging paper and individual packaging film).

Packaging is free from holes, cracks, tears, wrinkles, thickenings and attenuations.

Paper for individual package:

Density of the packing material: 65(±2) g/m²

package thickness: 80(±10) micron.

package size: 410 mm x 65 g/m²

Air permeability – 200-400 cm³/m²•min.

The density of glue application – 90(±2) g/m².

The structure of the adhesive layer is continuous.

Adhesion strength – not less than 0.12 MPa.

Adhesive layer width – not less than 8 mm.

Tensile strength – not less than 12 MPa

Sterilization validity period – 5 years.

Film for individual package:

Density of the packing material: 1.32–1.36 g/cm³,

package size: 420 x 0.070 mm

Tensile strength – not less than 0.25 MPa.

Tensile at break – not less than 20 N.

The structure of the adhesive layer is continuous.

The density of glue application – 90(±2) g/m².

Adhesion strength – not less than 0.12 MPa.

Adhesive layer width – not less than 8 mm.

Sterilization validity period – 5 years.

-The device is packed in a heat-sealable package made of paper and polymer film/Soft blister: The device is sterilized by ethylene oxide (EtO) according to a standardized and validated sterilization cycle in accordance with EN ISO 11135-1:2014. The package, in which the device is packed, is made from polypropylene and polyethylene film, and packaging paper for medical devices. Packaging paper for medical devices has a structure that ensures maximum circulation of ethylene oxide, air and moisture. Description of the standard sterilization cycle is provided below and is based on the results of validation:

No.	Detailed information	Maximum values
1	Preliminary treatment (i) Duration (ii) Chamber temperature (iii) Humidity	60 minutes ± 2 minutes 45°C (±5°C) 30% - 90%
2	Conditioning (i) Degree of vacuum and vacuum rate (ii) Vacuum retention time (iii) Chamber temperature (iv) Chamber humidity	0.75 kg/cm ² for 30 min ± 15 min 10 minutes ± 1 minute 45°C (±5°C) 30%-90%
3	Sterilization (i) EtO temperature at the input (ii) EtO concentration (iii) Chamber temperature (iv) Duration of exposure	Not less than 20°C 550 mg/l ± 25 mg/l 45°C (±5°C) 280 minutes
4	Ventilation (i) Degree of vacuum and vacuum	-0.75 kg/cm ² for 40 min ±

rate	20 min
Total number of ventilations procedures	

- The device is packaged in high-molecular (HM) film bags: The device is sterilized with ethylene oxide (EtO) according to a standardized and validated sterilization cycle in accordance with EN ISO 11135-1:2014. The device is packed in high-molecular (HM) bags made of high-density polyethylene, which provide maximum circulation of ethylene oxide, air and moisture. Description of the standard sterilization cycle is provided below and is based on the results of validation:

No.	Detailed information	Maximum values
1	Preliminary treatment	
	(i) Duration	60 minutes $\pm$ 2 minutes
	(ii) Chamber temperature	45°C ($\pm$ 5°C)
	(iii) Humidity	30% ~ 90%
2	Conditioning	
	(i) Degree of vacuum and vacuum rate	-0.70 kg/cm ² for 10 $\pm$ 1 min
	(ii) Vacuum retention time	10 minutes $\pm$ 1 minute
	(iii) Chamber temperature	45°C ($\pm$ 5°C)
	(iv) Chamber humidity	30%-90%
3	Sterilization	
	(i) EtO temperature at the input	Not less than 20°C
	(ii) EtO concentration	550 mg/l $\pm$ 25 mg/l
	(iii) Chamber temperature	45°C ($\pm$ 5°C)
	(iv) Duration of exposure	280 minutes
4	Ventilation	
	(i) Degree of vacuum and vacuum rate	-0.75 kg/cm ² for 40 min $\pm$ 20 min
	Total number of ventilations procedures	2

The appearance and performance of the device does not deteriorate even after three cycles of ethylene oxide sterilization, since the device has been sterilized three times and the integrity of the package has not been violated. The residual amount of ethylene oxide was also within the normal range.

Re-validation is carried out once a year. Calibration of instruments, such as a pressure sensor, a temperature sensor, temperature indicator and regulator, PT 100, etc., is carried out periodically according to the calibration schedule.

Package size:

Name of the infusion set	Package size LxWxH, mm	Quantity, pcs
Polyfusion IV infusion set without airvent	760x480x325	500
Polyfusion IV infusion set with airvent	760x480x325	500
Microfusion IV infusion set with a micro dropper	760x480x325	500
Polytrol IV infusion set with flow regulator	760x480x325	300
Novofusion IV infusion set with airvent	760x480x325	500
Novofusion IV infusion set without airvent	760x480x325	500

Oncofusion IV infusion set with filter and airvent	415x400x255	100
Polyfusion IV infusion set with three-way stop cock and airvent	760x480x325	500
Polyfusion IV infusion set with three-way stop cock without airvent	760x480x325	500
Neofusion IV infusion set with double dropper	760x480x325	400
Polyfusion-W infusion set with winged spike without airvent	760x480x325	500
Polyfusion-W infusion set with winged spike with airvent	760x480x325	500
Photofusion infusion set for photosensitive drugs	415x400x255	150
Polytrol-Micro infusion set with micro dropper and flow regulator	760x480x325	500
Polyvol measured volume fluid administration set	770x370x410	100
utofusion infusion set with airvent	860x480x270	300

14. Maintenance

The sets are ready for use and do not require any maintenance or repair during use.

15. Cleaning and disinfection

Not required. The device is disposable.

16. Manufacturer's warranty

Poly Medicure Limited GUARANTEES THAT THIS DEVICE IS MANUFACTURED ACCORDING TO APPLICABLE STANDARDS AND SPECIFICATIONS. PATIENT'S CONDITION, CLINICAL TREATMENT AND MAINTENANCE OF THE DEVICE MAY AFFECT THIS DEVICE'S PERFORMANCE. THIS DEVICE SHOULD BE USED ACCORDING TO RECOMMENDATIONS OF THE TREATING PHYSICIAN.

The finished sets meet the technical and functional requirements during five years from the sterilization date.

The manufacturer guarantees that the sets meet the requirements of technical documentation in case of the user's compliance with the storage and transportation conditions, as well as the rules of operation.

17. The medical device does not contain any medications.

24 FEB 2020

18. Authorized representative of the manufacturer in the Russian Federation (name, address, telephone, mail)

Authorized representative of the manufacturer: Advamedics LLC

Legal address: Office 15, 2nd floor, block 2, building 10, Korovinskoye shosse, Moscow, 127486, Russia
+7(495) 744-34-79

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