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

page 3

BEFORE USE

Dear User of SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter)!

Read this User's Manual carefully before using the product.

Pay special attention to all warnings and precautions that are marked with sign.

All elements of the System are designed, tested and checked for compatibility in such a way as to ensure accurate results of measurement of blood glucose level.

Use only SelfyCheck Prim test strips and SelfyCheck Prim control solution for SelfyCheck Prim blood glucose meter.

Regular glycemic control allows you and your doctor to monitor the diabetes in a better way. Due to the small size and simplicity of operation the System provides a possibility of self-monitoring of the blood glucose level anywhere and at any time.

Many people, performing the blood analysis themselves, find it useful to start with practising on the control solution (see Section "CHECK OF SYSTEM OPERATION USING THE CONTROL SOLUTION").

Please contact the support service or place of purchase regarding the issues related to the use of this device.

MEDICAL DEVICE NAME

SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter) (the "System", "Device").

MEDICAL DEVICE COMPOSITION

The medical device "SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter)" includes the following versions:

Version 1:

1. SelfyCheck Prim Portable Blood Glucose Meter for Measurement of the Blood Glucose Level – 1 pce.
2. Lancing Device (with AST tip) – 1 pce.
3. SelfyCheck Lancets (10 pcs. in a pack) – 1 pce.
4. User's Manual – 1 pcs.
5. Case – 1 pcs.

Version 2:

1. SelfyCheck Prim Portable Blood Glucose Meter for Measurement of the Blood Glucose Level – 1 pce.
2. SelfyCheck Prim Test Strips (50 pcs. in a pack)
3. Lancing Device (with AST tip) – 1 pce.
4. SelfyCheck Lancets (10 pcs. in a pack) – 1 pce.
5. User's Manual – 1 pcs.
6. Case – 1 pcs.

Version 3:

1. SelfyCheck Prim Portable Blood Glucose Meter for Measurement of the Blood Glucose Level – 1 pce.
2. SelfyCheck Prim Test Strips (50 pcs. in a pack)
3. SelfyCheck Prim (Norm) Control Solution – 1 pce.
4. Lancing Device (with AST tip) – 1 pce.
5. SelfyCheck Lancets (10 pcs. in a pack) – 1 pce.
6. User's Manual – 1 pcs.
7. Case – 1 pcs.

3

page 4

Version 4:

1. SelfyCheck Prim Portable Blood Glucose Meter for Measurement of the Blood Glucose Level – 1 pce.
2. User's Manual – 1 pcs.
3. Case – 1 pcs.

Version 5:

1. SelfyCheck Prim Portable Blood Glucose Meter for Measurement of the Blood Glucose Level – 1 pce.
2. SelfyCheck Prim Test Strips (50 pcs. in a pack) – 1 pce.
3. SelfyCheck Prim (Low) Control Solution – 1 pce.
4. SelfyCheck Prim (Norm) Control Solution – 1 pce.
5. SelfyCheck Prim (High) Control Solution – 1 pce.
6. Lancing Device (with AST tip) – 1 pce.
7. SelfyCheck Lancets (10 pcs. in a pack) – 1 pce.
8. User's Manual – 1 pcs.
9. Case – 1 pcs.

PURPOSE

The medical device "SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter)" is intended to quantify the level of glucose for use as a support for monitoring glucose levels in fresh capillary whole blood (taken from the finger, palm, forearm, upper arm) for self-testing by people with diabetes at home and at stay in clinical settings. The system is not intended to be used for treatment in relation to any symptoms or diseases or testing newborns. The measurement data is for reference only.

SYSTEM COMPONENTS


Blood glucose meter
(with a battery)

Control solutions

Lancing device
(with AST tip)

Sterile lancets – starting set
in the bag

Test strips – starting set in
the container

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page 5

TESTING PRINCIPLE

The device is intended for individual use.

SelfyCheck Prim Blood Glucose Meter as a part of the device is used together with SelfyCheck Prim Test Strips from the starting set included in device composition, and with SelfyCheck Prim purchased separately in the future.

For check of system operation it is necessary to use SelfyCheck Prim Control Solutions with low, normal and high concentration from the starting set included in device components.

Lancing device (with additional AST tip) (included in the device components or similar device by other manufacturer is purchased separately) with installed sterile single-use SelfyCheck Lancet (included in the device components or similar device by other manufacturer is purchased separately) is used to puncture the skin in the point of blood sampling.

For easy storage the device components are placed in the special case.

The device components include User's Manual to ensure correct use of the device, self-study and learning of skill of self-monitoring of blood glucose level.

Determination of glucose in the blood sample is based on the measurement of electric current induced by the reaction of glucose and enzyme on SelfyCheck Prim Test Strip.

After inserting the test strip into the blood glucose meter test port the System is brought to the blood sample with the free end of the test strip. The blood is automatically drawn into the absorption channel where the reaction takes place. The transient electric current is generated during reaction. The system reads the current strength value obtained during the reaction between the enzyme and glucose in the blood sample in the reaction cell of the test strip and calculates the concentration of blood glucose level using the pre-set equation of the software built into the blood glucose meter. The current strength value depends on the blood glucose concentration, batch of test strips, operating voltage and temperature. The result is shown on display. The blood sample for measurement of blood glucose level can be obtained not only from the fingertip, but also from other body parts.

INDICATIONS

Main indications for usage of the medical device – diabetes mellitus type 1 and 2.

CONTRAINDICATIONS

There are no absolute contraindications. There are use limitations.

POSSIBLE SIDE EFFECTS

Not identified.

SCOPE

The medical device "SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter)" is used in the endocrinology.

POTENTIAL USERS

The device can be used for self-monitoring by people at home as well as at stay in clinical settings. The device can be used by lay users with various physical and mental abilities not having special medical, scientific, and technical knowledge with regard to blood glucose level monitoring, having read the User's Manual in advance.

DEVICE USE LIMITATIONS

⚠ The System is designed for external use (for in vitro diagnostics).

⚠ It is not recommended to use the System in dry environment, especially with synthetic materials (synthetic clothes, carpets, etc.).

5

page 6

It can lead to the damage by static discharge and, therefore, to the erroneous results.

⚠ It is prohibited to use the System in close proximity to the sources of strong electromagnetic radiation.

⚠ If the System is used at home, the measurement is performed only on the samples of fresh capillary whole blood; in the medical institutions the healthcare professionals can use the samples of fresh capillary and venous whole blood.

⚠ The System is not intended for use by persons with special needs on their own (children under 18, disabled persons, the elderly). Their usage of the device requires assistance and help.

⚠ The System is not designed for use in combination with other medical devices.

⚠ The System shall not be applied to patients with severe infections, malignant neoplasms, massive oedema.

⚠ The System is not intended for diagnosis of diabetes mellitus.

⚠ The System is not designed for measurement of blood glucose level in the new-born infants.

Possible factors of impact on the results of analysis:

⚠ The checked blood was squeezed out of finger by strong press and therefore contains the interstitial fluid.

⚠ Shelf life of test strips has expired (shelf life is 21 months from the date of manufacture or 6 months from the date of container opening).

⚠ Test strips were stored at improper temperature (proper storage temperature: from +2°C to +30°C, not more than 85% of relative humidity).

⚠ Insufficient moisture protection of test strips.

⚠ Sanitizer, used in the preliminary disinfection of the point of blood sampling, not completely vaporized from the skin.

⚠ The measurement was performed with unwashed hands.

⚠ After washing the hands were not dry for further measurement of blood sugar level.

⚠ The excessively high blood pressure can lead to the distortion and erroneously low measurement results.

⚠ Very high and very low hematocrit (lower than 35% or higher than 60%) can lead to the incorrect results. Hematocrit value shall be within the range from 35% to 60%.

⚠ The results can be inaccurate in case of the patients with very low arterial blood pressure or those who are in shock.

⚠ Endogenous substances, naturally present in the human body (bilirubin, uric acid, glutathione, triglycerides, other medicament chemical compounds), as well as medicinal products do not have a significant impact on the results of blood glucose level measurement using a blood glucose meter if normal concentrations are not exceeded or in case of standard therapeutic doses. High concentration in blood can affect the analysis results. It is recommended to pay special attention to the interpretation of these results.

⚠ In case of patients with impaired peripheral blood circulation, it is not recommended to collect the blood from alternative body parts, because the obtained readings can differ from actual blood glucose level. Such non-conformities can be observed in the following cases: In case of severe dehydration as a result of diabetic ketoacidosis, hyperglycaemia due to the stress, hyperosmolar nonketotic com, shock, decompensated heart failure of class IV according to NYHA scale or peripheral arterial occlusive disease. Measurement errors can lead to the incorrect interpretations and thus to the serious health consequences.

⚠ Some (clinical) states, such as arterial hypotension, state of shock or ketoacidosis, can become the reason of wrong results of measurement.

APPLICATION CONDITIONS

The device is intended for application in hospitals, outpatient departments, medical institutions, ambulance cars as well as at home and in the field.

APPLICATION PRECAUTIONS

⚠ It is necessary to follow the instructions of the User's Manual when using the device.

⚠ Read the manual and perform the dummy test before using the device for measurement of blood

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Black

page 7

glucose level. Perform all quality checks according to the manual.

⚠ The test strips shall be supplied and stored in the plastic container with tightly closed cover. The container cover shall be tightly closed immediately after taking out one test strip. Do not touch the strips with wet or dirty hands.

⚠ The lancets and test strips are designed for single-use, their re-use is prohibited and may cause cross infection. If the previously used test strips are applied, the blood glucose meter displays the error message on the screen.

⚠ Antiseptic rules shall be followed when using the device. It is necessary to wash hands with soap in warm water or use wet disinfecting wipe and let the skin dry before use.

⚠ It is necessary to follow standard precautions when handling blood. All samples and materials that came in contact with the patient are considered to be biohazardous, they shall be handled as potential infection carriers.

⚠ The device contains small details that may cause asphyxia and injuries in young children in case of inhalation and swallowing. The device shall be kept out of reach of children!

⚠ The glucose meter shall not be exposed to the devices which generate strong magnetic fields or ionizing radiation.

⚠ Do not remove the test strip from the test slot of the blood glucose meter during the analysis.

⚠ It is prohibited to disassemble the body of the blood glucose meter on your own.

⚠ Do not apply high force on the body or drop the blood glucose meter.

⚠ Do not use the blood glucose meter located on the cold or hot surface.

⚠ The user shall always have the spare battery in case of replacement.

⚠ Do not use the components which are not specified by the manufacturer.

⚠ Do not use the blood glucose meter in case of improper operation or damage.

⚠ In order to extend the life time of the device it is necessary to ensure the proper service and periodic testing using the control solution. If you are not sure that the measurements are accurate, contact the support service or seller for help.

⚠ Use of this device in dry environment, especially in contact with the synthetic materials (synthetic clothes, carpets, etc.), can lead to the electrostatic discharge which may damage the device, cause the erroneous results and effect the accuracy of device operation.

⚠ It is prohibited to use the blood glucose meter with damaged body or damaged test strips.

⚠ Self-measurement of blood glucose level does not replace the diagnosis by the doctor or regular appointments. Self-testing of blood glucose level shall be used only after receiving the detailed instructions for monitoring and treatment of diabetes from the qualified medical specialist (endocrinologist);

⚠ Not to swallow of control solution;

⚠ Not to allow hit of control solution on skin or eyes in order to avoid irritation.

7


page 8

DESCRIPTION AND CHARACTERISTICS OF SYSTEM COMPONENTS

SelfyCheck Prim Blood Glucose Meter


Appearance

1. Test Slot/ Data Port
2. Front Panel
3. Display
4. Buttons: top/ bottom/ main (M)
5. Battery Compartment Cover
6. Back Panel
7. Ejector Pin of a Test Strip
8. Pad

GLUCOSE METER DISPLAY

1. Test Result
2. Ketonic Warning
3. Error Message
4. Storage Mode
5. Quality Control Mode
6. Measurement Mode
7. Date
8. Time
9. Mean Value per Day

10. Warning Signal Symbol
11. Unit
12. Indicator of High/Low Value
13. Low-Battery Warning
14. Test Strip Symbol
15. Blood Drop Symbol

Specifications

Power source:	1.5V AAA alkaline battery
Battery life:	at least 500 measurement cycles
Display:	LCD, white illumination
Memory:	1,000 measurement results with date and time
External output:	mini-USB cable
Enzyme type:	GDH-FAD (glucose dehydrogenase-flavin-adenine dinucleotide)
Analysis time:	5 s
Measurement units:	mg/dl or mmol/l
Measurement range:	20-650 mg/dl (1.1-36.1 mmol/l)
Resolution (measurement interval):	0.1 mmol/l (1 mg/dl)
Sample size:	0.5 µl
Hematocrit level:	35-60%
Consumed power:	157.08 mW
Dimensions:	52.3 x 90.3 x 18.0 (±0.3) mm
Weight (with battery):	57.4 g ± 15%
Operating mode setting time:	90 sec
Additional functions:	Automatic detection of sample Automatic countdown start Automatic shutoff after 3 minutes with no activity Warning of violation of the temperature conditions

⚠ Always use a special case for storage and transportation of the System.

⚠ Avoid falling and heavy shocks of the blood glucose meter.

⚠ Blood glucose meter shall not be exposed to the direct sunlight and high humidity.

⚠ Consult your doctor regarding the measurement unit used in your therapy.

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page 9

Software Description

The control software used in the blood glucose meter is installed by the manufacturer in plant conditions and constitutes its integral part.

The role of software consists in monitoring of all device functions and provision of the diagnostics information. The main software controlled device functions include the following:

- blood glucose sensor for analysis of the blood glucose level value;
- temperature sensor for temperature measurement;
- keypress event for reading of the user input;
- LCD drive for display of information on LCD panel;
- communications interface for data transfer;
- memory recall for display of previous test results;
- the buzzer driver for switching on or off;
- the battery charge controller for charge determination and control;
- interface for data transfer.

The user controls the software by pressing the buttons or by insertion / removal of the test strip. Thus, the user can change the parameters which control the buzzer, charge controller and LCD.

Interface language:	English
Program seek time:	not more than 0.5 s
Program response time:	not more than 0.1 s
Software hazard class:	B
Software version:	2011_0309 PM 02:08

Operating Environment

Microcontroller:	TAIDOC MD437
ROM capacity:	32 Kbyte
Software name:	TD-4287B
Software Version:	2011_0309 PM 02:08
Programming language:	The C Programming Language

System Requirements

CPU:	TAIDOC MD437
Memory capacity:	RAM of 1,024 bytes; ROM of 32 Kbytes
Battery:	1.5V AAA (1 pce.)

SelfyCheck Prim Test Strips

Specifications	Model: TD-4311
Chemical components:	Enzyme – Glucose Dehydrogenase 8% Potassium phosphate 8% Ferricyanide 55% Carboxymethyl cellulose 29%
Dimensions of test strip:	7.5 x 36.06 x 0.5 (±0.2) mm
Dimensions of container:	53.0 (±0.5) mm, Ø 30.46 (±0.25) mm.
Weight of one test strip:	0.17 (±0.003)g.
Gross weight of container:	- 50 pcs. in container: 24.8 g ± 10%
Sample size:	0.5 µl
Measurement units:	mg/dl or mmol/l
Reaction time:	5 s
Measurement range:	from 1.1 to 36.1 mmol/l (from 20 to 650 mg/dl)
Analyte:	Dextroglucose
Enzyme:	GDH-FAD (glucose dehydrogenase-flavin-adenine dinucleotide)
Shelf life:	21 months from the date of manufacture or 6 months from the date of container opening.
Hematocrit range:	35–60%
Sample type:	fresh capillary whole blood

page 10

The test strip consists of the following parts:

1. Absorbent Hole

Area of contact of test strip and blood sample. In case of contact the blood is absorbed automatically.

2. Confirmation Window

The filling of the window confirms that there is enough blood for analysis.

3. Test Strip Handle

Area for holding the test strip when inserting into the test slot.

4. Contact Bars

Area of contact of the test slot and test strip. When inserting the test strip into the test slot, the blood glucose meter will switch on automatically.

⚠ The test strip shall be inserted firmly until it will go no further with front side facing up.

⚠ If contact bars are not fully inserted into the blood glucose meter, the test result can be incorrect.

⚠ The test strips shall be used within 6 months from the date of first opening of the vial (but not later than 21 months from the date of manufacture). Write the date of first opening of the vial on the label.

⚠ The container cover shall be tightly closed immediately after taking out of one test strip from the container.

⚠ Touch the test strip only with clean and dry hands. Use each test strip immediately after taking it out of the container.

⚠ Do not use the test strip after expiration of the shelf life. It can lead to the incorrect results.

⚠ Do not bend, cut or alter the test strip.

Test strips are produced in two forms of release:

- 1) SelfyCheck Prim Test Strips - with opening control;
- 2) SelfyCheck Prim Test Strips - without opening control.

SelfyCheck Prim Control Solutions

SelfyCheck Prim Control Solution is intended for usage in combination with SelfyCheck Prim Test Strips. SelfyCheck Prim Control Solution is used for quality control of the SelfyCheck Prim Portable Blood Glucose Monitoring System (Blood Glucose Meter). The device is intended for in vitro diagnostics only.

Specifications

Model	SelfyCheck Prim (Low) Control Solution	TD-4907
	SelfyCheck Prim (Norm) Control Solution	TD-4908
	SelfyCheck Prim (High) Control Solution	TD-4909
Nominal vial volume:	4 ml	
Vial dimensions:	58.2(±0.6) mm Ø 17.5 (±0.35) mm	
CV, %:	SD < 5 mg/dL (0.278 mmol/L) at < 100mg/dl (5.55 mmol/l); CV < 5% at 100mg/dl (5.55 mmol/l)	
Weight:	6.6 g ± 10%	
Composition:	Glucose: 0.05–0.2% Inert ingredients: > 99.7%	

⚠ Do not use control solution after the expiration date or 3 months after the initial opening of the vial. Note the date of opening on the vial with the control solution and dispose of the remaining solution after 3 months.

⚠ It is recommended to perform the control solution test at the temperature from 20°C to 25°C and relative humidity from 30% to 85%. Make sure that the control solution, blood glucose meter and test strips are within the specified temperature range before testing.

⚠ Shake the vial before use. Remove the first drop of control solution and wipe the tip of the dispensing hole to ensure sample cleanliness and result accuracy.

Lancing Device (with AST Tip)


The lancing device is used in combination with the lancet to reduce the sense of discomfort due to the puncture of the skin when obtaining the blood sample.

The lancing device is designed so that the fixed lancet is fully hidden under the cap. Due to the compression of the extended spring (activated using the push-button) the lancet needle punctures the skin, then the spring moves back to the initial position with the lancet needle coming back under the lancing device cap.

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
Black

page 11

The device provides six options in terms of depth of skin puncture:

The more the indicator line the deeper the puncture.

Specifications Lancing Device (with AST Tip)	Model: TD-5010
Cocking control spring force:	6.0 H ± 10%
Impact energy (impact on the lancet):	45 mJ ± 10%
Weight:	16.5 g ± 15%
Dimensions:	118.75 (±0.2) mm x 17.45 (±0.1) mm

SelfyCheck Lancets

The lancets are designed for skin puncture for the purpose of obtaining the blood sample. Usually they are used in combination with the lancing device to reduce the sense of discomfort during skin puncture.

SelfyCheck lancet is a sharpened needle of high quality stainless steel SUS304 enclosed in the plastic body. The protective cap of the lancet is a sterile barrier and shall be removed from the needle point only after fixation of the lancet in the lancing device holder slot.

SelfyCheck Lancets	Model: TD-5084
Dimensions of one lancet, mm:	Length of the body with protective cap put on 32.0(±1.5) mm Length of the exposed part of the needle 3.2(±0.3) mm
Needle gauge:	30G/0.3 mm
Net weight, g:	0.4 g ± 5%
Gross weight (10 pcs), g:	4.7±0.5
Roughness:	Ra not more than 0.8 µm
Hardness:	at least 450 HV _{0.05}
Sterilisation:	cobalt-60γ-radiation
Sterilisation dose:	18.7 kGy
Shelf life of the lancets:	5 years after radiation sterilisation

Case

The case is the rectangular box with a zip-fastener, sewed in the sides. Dimensions of the case zipped up, mm: 160 x 100 x 30 (±5).

The case protects the device from the damages, scratches, dirt, dust, ingress of moisture, severe shocks. The zip-fastener runs along three sides and allows for the full opening of the case. The case ensures the space effective storage of the System components.

11

page 12


OPERATION

COMMISSIONING


Before the first use of the blood glucose meter or after replacement of the battery the settings specified below shall be checked or replaced.
Simultaneously press "▲" and "▼" buttons on switched off device (without inserted test strip).

Setting the Date


Procedure of date setting: YEAR → MONTH → DAY.

When the year is flashing press "▲" and "▼" buttons until the required digits appear. Press "M".

When the month is flashing press "▲" and "▼" buttons until the required digits appear. Press "M".

When the day is flashing press "▲" and "▼" buttons until the required digits appear. Press "M".

Setting the Time Format

Press "▲" or "▼" button to select the required time format – 12 h or 24 h. Press "M".

Setting the Time

When the hour is flashing press "▲" or "▼" button until the required digits appear. Press "M".
When the minutes are flashing press "▲" or "▼" button until the required digits appear. Press "M".

Setting the Measurement Units

When the unit of measurement is flashing press "▲" or "▼" button to choose between mg/dl and mmol/l. Press "M".


Setting of the target range of high and low values for the measurement mode

Sequence of setting of the target range of high and low values: Gen low → Gen high → AC low → AC high → PC low → PC high.
Taking in account settings described above, press "▲" or "▼" button until the required target value appears when the sequence is flashing. Press "M".

On/Off sound signal setting

When the vibration indicator is flashing, press "▲" or "▼" button to select between "On" (turn on) and "Off" (turn off) modes. Press "M".

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page 13

Memory Clearance

When "Del" and "M" symbols are shown on a display, press "▲" button, select "No" to save the results to the device memory, then press "M".
To delete all results from the device memory press "▲" button, select "Yes" and press "M"

Setting the Reminder Alarms

The device has four sound reminder alarms. "On" or "Off" and "🔊" symbols will be displayed on the screen. If there is no need to set the alarm, press "▲" or "▼" button to select "OFF" symbol, then press "M" to skip this step. Or select "On" symbol and press "M" button to proceed with setting. When clock digits flash one by one, press the main button to select the correct time. Press SET to finish.

When hour/minute setting indicator is flashing, press "▲" or "▼" button to set the required time values. Press "M" to go to the next alarm setting.

⚠ When alarm is activated, press "M" to deactivate it. Otherwise, the device will issue a sound signal within 2 minutes.

⚠ These parameters can be changed ONLY in the setting mode.

⚠ SelfyCheck Prim Blood Glucose Meter shall be used only with SelfyCheck Prim Test Strips. The use of other test strips can lead to the incorrect results.

CHECK OF SYSTEM OPERATION USING THE CONTROL SOLUTION

It is recommended to perform the independent check of the joint operation of the blood glucose meter and test strip using SelfyCheck Prim Control Solution (supplied in the stating set). The control solution is put into the vial with the seal cap.

The vial with "Low" marking is required for check and display of the low glucose level; the vial with "Norm" marking is required for check and display of normal glucose level; the vial with "High" marking is required for check and display of high glucose level.

The allowable range of results in the control solution analysis is indicated on the container with test strips. To obtain accurate results, it is necessary to conduct such diagnostics regularly.

The control solution test shall be performed:

- if you suspect that the blood glucose meter or test strips work improperly;
- if your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate;
- for practising the testing process;
- if you have dropped or think you may have damaged the device.

PERFORMING A CONTROL SOLUTION TEST

For performance of a control solution test you will need:


- blood glucose meter
- test strips
- control solution

13

page 14

1. Insert a test strip

Insert the test strip into the slot with the contact strips forwards and upwards. The Blood Glucose Meter will turn on automatically with "🔊" and "🔋" symbols appearing on the display.

2. Mark the test with the control solution

Press the "M" button to mark the test with the control solution. When the "QC" symbol is displayed on the device, the test results are stored in memory. If you decide not to carry out the test with the control solution, press the "M" button again and the "QC" symbol will disappear.

⚠ When conducting a test with a control solution, it is necessary to switch to the "QC" mode so that the result is not stored in the device memory along with the TEST RESULTS with blood. Otherwise, the results of the blood glucose test will be mixed in the device memory with the results of the control solution test.

3. Bring the test strip to the drop of control solution

Turn the closed vial over several times (do not shake) before opening. The drop of the control solution shall be applied to the clean vial cap; then the test strip, inserted into the blood glucose meter, shall be brought to the drop. Once the absorbent hole of the test strip touches the drop of the control solution, the drop will be drawn in automatically.

⚠ To avoid contamination of the control solution, do not apply the solution directly to the test strip.

4. Read and verify the result


The confirmation window of the test strip will be filled in case of sufficient sample size; after that the blood glucose meter will produce a sound signal and start the test. After the end of the countdown from 5 to 0, the test result with the control solution will appear on the display. Compare the result with the range indicated on the test strip container. The result must be within this range.

If the result is outside the specified range

If the repeated test results continue to fall outside the range specified on the container, the device can work or be used improperly, or the blood glucose meter and test strips can be incompatible.

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page 15

DO NOT USE this kit for testing of a blood sample. Contact the support service or place of purchase for help.

⚠ The range of values for the test with the control solution, indicated on the container of test strips, is intended only for evaluating the results of the test with the control solution. This range is not a recommendation for blood glucose.

Performance of Blood Sample Test

For performance of the blood sample test you will need:

- SelfyCheck Prim Blood Glucose Meter
- SelfyCheck Prim Test Strips
- SelfyCheck Prim Control Solutions (Low, Norm, High) - For Control Solutions Test
- Lancing device (with additional AST tip)
- SelfyCheck Lancets.

⚠ Complete set of the System is specified on the package.


Components that can be purchased additionally:

- ⚠ Never share a lancet or the lancing device with other people;
- ⚠ Always use new sterile lancet. Lancets are for single use only;
- ⚠ Avoid getting lotions, oils, dirt and debris inside or on the lancets and the lancing device.

1. Setting up the Lancing Device for Blood Testing

1. Turn the lancing device cap and remove it by pulling it up.
2. Insert the lancet into the lancet holder and press it in order to fix it cap from the lancet firmly.
3. Remove the protective holder and press it in order to fix it cap from the lancet firmly.


4. Put the device cap back and turn it until clicking sound to fix it (see figure). The cap is fixed if the fixing indicator is in the extreme right position.
5. Hold the cap and choose the puncture depth turning the adjustable tip in any direction so that the arrow on the cap points at the required puncture depth.
6. Pull by the slope adjustment area until click. When the device is ready for use the release button colour will be changed.

⚠ If you did not hear the click possibly the device jammed before the lancet fixation.


15

page 16

2. Insert the Test Strip into the Blood Glucose Meter to Switch It On

Wait for " " and " " symbols to appear on the display.

3. Select the Suitable Measurement Mode

The device has four measurement modes: Gen, AC, PC and QC.

MODE	CONDITIONS FOR USE
Gen	any time of day, without regard to timing of food ingestion
AC	8 hours before meals minimum
PC	2 hours after a meal
QC	diagnostics using the control solution

The switching between the modes is performed in the following way:


1. " " symbol or "Gen" will appear on the display after switching on

2. Press " " button to switch the modes Gen, AC and PC.

3. Press "M" button to switch to QC mode.

Obtaining the Blood Sample

Preparation of the Puncture Site

Stimulating blood circulation has a significant impact on the quality of the obtained result of blood glucose level measurement; to do so rub the puncture site before blood draw. The value of the glucose level in the blood collected from the site that has not been rubbed preliminarily differs essentially from the value of glucose level in the blood collected from the finger. Rubbing the puncture site before the blood draw allows for significant reduction of this difference.

Please follow the instructions below before collecting a drop of blood:

- ⚠ Wash and dry your hands before starting the procedure;
- ⚠ Select the puncture site at fingertip or on other part of the body;
- ⚠ Rub the puncture site for about 20 seconds before penetration;
- ⚠ Clean the puncture site using disinfecting wipe and let the skin dry.

Blood Drawing from the Finger

1) Puncture the finger

Press the lancing device tip to the bottom part of the finger pad tightly. In order to make puncture press the release button. After the puncture completion the click will sound.

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page 17

2) Remove the device from the puncture site.
Use the lancing device to prick the skin. After the puncture completion remove the first blood drop using clean cotton. It is necessary to obtain not less than 0.5 µl of blood sample. Press to the puncture site carefully to collect another blood drop. Be careful in order NOT to smear the blood sample.

- ⚠ Choose a different spot each time you perform a test. Repeated punctures at the same spot may cause pain feelings and calluses.
- ⚠ Consult the doctor before collecting the blood from the alternative sites.

5. Apply the Blood Sample

It is necessary to bring the test strip inserted in the glucose meter to the drop of blood at an angle.
Once the absorbent hole of the test strip touches the drop of blood, the blood will be drawn in automatically. The confirmation window of the test strip will be filled in case of sufficient blood sample size; after that the blood glucose meter will produce a sound signal and start the test.

- ⚠ DO NOT REMOVE the finger until you hear the sound signal.
- ⚠ Do not press the puncture site to the test strip and do not smear the blood.
- ⚠ If you do not apply the blood sample to the test strip within 3 minutes, the device will switch off automatically. To start a new measurement remove and replace the test strip.
- ⚠ The blood shall completely fill the confirmation window to start the count down. DO NOT try to add more blood to the test strip after the first drop is absorbed. Remove the used test strip and repeat the measurement using new test strip.
- ⚠ If you have trouble filling the confirmation window, consult your doctor or support service.

6. Reading of Result

The result of blood glucose test will appear after the countdown from 5 to 0 and will be saved to the device memory automatically (100 mg/dl = 5.5 mmol/l).

- ⚠ Severe dehydration and excessive water loss can lead to the readings lower than the actual values. Immediately contact your doctor in case of suspicion of severe dehydration.
- ⚠ If the blood glucose level readings are lower or higher than the normal value, and the corresponding symptoms of the disease are absent, repeat the test first. Follow the recommendations of your doctor in case of symptoms or if the repeated tests continue to show the results lower or higher than the normal value.
- ⚠ To measure the blood glucose level use only fresh whole blood samples. In case of use of other samples the results will be incorrect.
- ⚠ If you have the symptoms that do not correspond to the test results and subject to adherence to the instructions provided in this Manual.
- ⚠ Check, if the result of analysis is within the acceptable range of results. Refer to your doctor, if the value exceeds 14 mmol/l.

17

page 18

The System allows obtaining a result equivalent to the value for blood plasma.

Time of measurement	Range of normal values of blood plasma glucose level for people without diabetes (mg/dl)
Fasting and before meal	< 100 mg/dl (5.6 mmol/l)
2 hours after a meal	< 140 mg/dl (7.8 mmol/l)

⚠ Consult your doctor to define the target range of values appropriate for you.
Blood or plasma glucose concentration level can be expressed in units of weight (mg/dl) or molar concentration (mmol/l). Principle of approximate conversion of mg/dl value to mmol/l value is as follows:

mg/dl	divide by 18	= mmol/l
mmol/l	multiply by 18	= mg/dl

For example, 1. 120 mg/dl ÷ 18 = 6.6 mmol/l;
2. 7.2 mmol/l × 18 = 129 mg/dl
(conversion results are approximate).


7. Disposal of Used Materials

Removal of the used test strip

Remove the test strip by moving the glucose meter ejector pin slider to the end position; the glucose meter will switch off automatically.

Removal of the used lancet

1) Twist off the device cap and remove the lancet. Dispose the used lancet following precautions: put the lancet protective cap on a hard surface and insert the used needle in it.

2. In order to remove the used lancet move the ejector forward. Dispose the used lancet according to the regulations of local legislation.

⚠ Used lancets and test-strips are considered biohazardous waste. Dispose of the used lancet according to the regulations of the local legislation.

⚠ Use a transparent tip for blood draw from alternative sites (not from the finger pad).

BLOOD DRAW FROM THE ALTERNATIVE SITE

You can perform the analysis of blood sample collected from the alternative sites.

⚠ There are limitations in terms of possibility of blood draw from the alternative sites. Please consult your physician before using this method.

Replace the cap of the delivery device with a transparent cap for blood draw from the alternative sites. Pull the cocking control back until it engages. When making injections to the forearm, upper arm or hand, avoid sites with well defines veins, to avoid excessive bleeding.


When should the method of blood draw from the alternative sites be used?

Dietary habits, treatment, physical activity level or stress and diseases may influence blood glucose level. Capillary blood in finger pads reflects such changes faster than blood in other parts of the body. Therefore **only the blood samples from fingers** should be used when measuring blood glucose level during or immediately after a meal, exercise or any other actions or events.
It is strongly recommended to perform blood draw from the alternative sites ONLY under the following conditions:

- Before a meal or on an empty stomach (not earlier than 2 hours after the last meal).
- Not earlier than 2 hours after the insulin administration.
- Not earlier than 2 hours after physical exercise.

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
page 19

DO NOT USE blood draw from the alternative sites in case:

- You assume low level of blood glucose.
- You are not sure whether you have hypoglycaemia.
- You are performing test for hyperglycaemia.
- The result of blood glucose level measurement from the alternative sites does not correspond to your state of health.
- The results of your ordinary checks are subject to fluctuations.

Alternative sites for the blood draw:

1. finger
2. palm
3. forearm
4. upper arm

GLUCOSE METER MEMORY

The device stores last 1,000 results of blood glucose level measurements with indication of date and time. To view the results saved in the device memory **start with the switched off device.**

1. Activation of Memory Access Mode

Press and release "M" or "▲" button on a switched off device. "M" symbol and first result displayed on the screen will appear – last measurement result with a date, time and measurement mode. Press "▲" or "▼" buttons to show the measurement results stored in the device memory. To switch the glucose meter off, press and hold "M" button.

2. Viewing Mean Measurement Values of Blood Glucose Level

- Press and release "▼" button to enter the device storage mode. Mean values will be displayed when "M" and "DAY AVG" symbols appear on the screen. Release "M" button. Mean measurement results for 7 days will be displayed on the screen in General mode.
- Press "▲" button to view mean measurement results stored in each of the measurement modes in the order of Gen, AC, PC for 14, 21, 28, 60 and 90 days.

19

page 20

3. Exiting the Memory Access Mode

After display of the last measurement result, hold "M" button to switch off the glucose meter.

△ In order to exit the memory access mode at any time, press and hold the "M" button for 3 seconds or do not press any buttons for 3 minutes. The blood glucose meter will switch off automatically.

△ If the device memory is filled, the very first measurement result will be automatically deleted and replaced by new one.

△ Measurement results using control solution are NOT included into mean measurement results.

PERFORMANCE CHARACTERISTICS

Precision

Precision shall correspond to the requirements of ISO 15197:2015 standard:

Glucose Concentration	≥ 5.55 mmol/l (100 mg/dl)	< 5.55 mmol/l (100 mg/dl)
Result	± 5%	± 0.28 mmol/l (5 mg/dl)

The System has acceptable intermediate precision and reproducibility with the value of coefficient of variation lower than 5% for each concentration level lower than 100 mg/dl and the standard deviation value within 5 mg/dl for each concentration level lower than 100 mg/dl.

Linearity

The linearity of the System corresponds to the acceptance criteria within the measured range of concentration from 20 to 650 mg/dl.

Accuracy

The System accuracy shall correspond to the requirements of ISO 15197:2015 standard:

- for glucose concentration < 5.55 mmol/l (< 100 mg/dl):
95% of results shall be within ± 0.83 mmol/l (± 15 mg/dl).
- for glucose concentration of ≥ 5.55 mmol/l (≥ 100 mg/dl):
95% of results shall be within ± 15%.
- 99% of measured values of glucose shall be within A and B zones of Consensus Error Grid.

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page 21

Interference

Substances and concentrations exceeding the permissible limit and causing the interference:

Substance	Maximum permissible concentration (mg/dl)	Therapeutic / normal concentration range (or upper limit) (mg/dl)
Acetaminophen (Paracetamol)	> 6.25	0.45-3
Ascorbic acid (Vitamin C)	> 5	2
Bilirubin (unconjugated)	> 20	0-2
Dopamine	> 1.25	0.03
Levodopa	> 0.7	0.02-0.28
Methyldopa	> 1.875	0.1-0.5
Reduced glutathione	> 30	24.25-32.2
Tolazamide	> 12.5	2-2.5
Uric acid	> 10	2-8
Galactose	> 250	< 5
Mannitol	> 5000	0.0128
Mannose	> 125	1.15
Xylose	> 3.125	Not applicable
Lipemic samples (triglycerides)	> 3000	30-300

Substances and concentrations causing the interference:

Substance	Maximum permissible concentration (mg/dl)	Concentration in the blood draw tube (mg/dl)
Sodium fluoride/potassium oxalate	< 250/200	250/200

Compounds that can lead to increase of results of blood glucose level measurement within the therapeutic concentration range:

Substance	Tested concentration (mg/dl)	Therapeutic range (or upper limit) (mg/dl)
Pralidoxime	> 5	~ 10 (the dose for intravenous administration is 500 mg)

21

page 22

PRODUCTION MATERIALS OF DEVICE CONTACTING SURFACES

When producing the device the components of animal origin are not used. The product does not contain substances that can be considered as drugs or pharmaceuticals. Materials of biological origin intended for recycling and re-use are not used in the device.

The blood glucose meter, test strips, control solutions, lancing device, lancets, storage case are in contact with the skin surface (fingers, palm, wrist) for a short time when in use.

The test strips, lancing device and lancets are in contact with a damaged skin (fingers).

The lancets are in transient intrusive contact with internal environment of the body. The period of continuous contact is less than 1 sec. The lancets are sterile, apyrogenic, non-toxic.

LABELING

Labels of the consumer packaging contain the following information: name, trademark, completeness, name and address of the manufacturer, address and contact information of the authorized representative, "In vitro diagnostics medical device" symbol, "Prevent exposure to sunlight" symbol, "Protect from moisture" symbol, "Temperature range" symbol, "Warning! Refer to the User's Manual" symbol; device bar code; summarized information on the possibilities and technical parameters of the device; unique serial number of device; shelf life; MA number and date.

Marking of the body of blood glucose meter contains the following information: name, model, trademark, name and address of the manufacturer, serial number, QR code; internal label code; "Warning! Refer to the User's Manual" symbol; "Refer to the User's Manual" symbol, "In vitro diagnostics medical device" symbol; "Direct current" symbol, CE conformance mark, battery parameters.


Label on the container body with test strips contains the following information: name, model, trademark, name and address of the manufacturer, name and address of the authorised representative, hot line, serial number, place for indication of the date of opening, number of test strips in the pack, "Warning! Refer to the instructions for use" inscription; "Best before" inscription (in YYYY-MM format); "In vitro diagnostics medical device" inscription; "Temperature range" inscription; "Humidity range" inscription; "Re-use is forbidden" inscription; range for qualitative interpretation of measured values in mmol/l and mg/dl units (low, normal, high levels) for control solution (in tabular form); "Use only with SelfyCheck Prim System and SelfyCheck Prim Control Solutions" inscription; "Use within 6 months after the initial opening" inscription; "Do not freeze" inscription; internal label code; date and number of the Marketing Authorisation; opening control (strip, glued to container for strips cover) depending on the release form.

Label on the control solution vial body contains the following information: name, model, trademark, name and address of the manufacturer, serial number, place for indication of the date of opening, number of the test strips in the package, "Warning! Best before" inscription (in YYYY-MM format); "In vitro diagnostics medical device" inscription; "Temperature range" inscription; "Humidity range" inscription; "Re-use is forbidden" inscription; range for qualitative interpretation of measured values in mmol/l and mg/dl units (low, normal, high levels) for control solution (in tabular form); "Use only with SelfyCheck Prim System and SelfyCheck PrimTest Strips" inscription; "Use within 6 months after the initial opening" inscription; "Do not freeze" inscription; internal label code; date and number of the Marketing Authorisation.

Marking of the lancet package contains the following information Marking of the lancet package (bag) contains the following information: device name and designation; model, name, trademark, address of the manufacturer and manufacturing facility; serial number; number of devices per package; indication of the needle gauge 30G/0.3 mm; "Warning! Refer to the User's Manual" symbol; "Radiation sterilisation" symbol; "Use by" symbol (in YYYY-MM format); "Re-use is forbidden" symbol; CE quality mark; internal label code; MA date and number.

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page 23

Symbols Used in the Device Marking					
Symbol	Value	Symbol	Value		
	In vitro diagnostics medical device		Manufacturer		
	Refer to the User's Manual		Serial number		
	Temperature range		Batch code		
	Best before		Humidity range		
	Warning! Refer to the User's Manual		Disposal with electrical and electronic equipment		
	"Radiation sterilisation" symbol (lancet use only)		CE marking		

PRINCIPLE OF CODING OF SYSTEM COMPONENTS LOTS

Blood glucose meter

Code	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16
Description	Model No.				Model version 1=A, 2=B, etc.	Year of manufacture (last 2 digits 00-99)			Month of manufacture	Serial number				Model version		
Example	4	1	8	3	1	1	8	2	0	0	0	0	0	1	8	B

Test strips

Code	1	2	3	4	5	6	7	8	9	10	11
Description	Production year (2 last digits in 00-99)				Month A-L	Lot	Date of manufacture 01-31		-	Q-ty A: 10 pcs. B: 25 pcs. C: 50 pcs.	Category code
Example	T	D	1	7	B	1	1	5	-	B	

1: T – Taoyuan (Place of manufacture)
2: D – test strip with carbon electrode
6: Lot number (daily) (No. 1,2)
10-11 Category code (A-Z, 0-9), in the absence of code – 0

Control solutions

Code	1	2	3	4	5	6	7	8
Description	Colour code	Identification code	Model version	Production year (2 last digits in 00-99)		Month A-L	Serial number	
Example	R	A	A	1	2	A	0	1

Code 1: Colour code: Y (yellow) – SelfyCheck Prim Low (low level)
G (green) – SelfyCheck Prim Norm (normal level)
R (red) – SelfyCheck Prim High (high level)

23

page 24

Lancets

Code	1	2	3	4	5	6	7
Description	Production year (2 last digits in 00-99)		Month 01-12		Date 01-31		T – for identification of lancets with twist-off cap
Example	1	7	0	5	2	3	T

MAINTENANCE

Battery


The blood glucose meter set includes one 1.5 V AAA alkaline battery.

Low Battery Level


The analyser will display one of the messages listed below as an alert on the low level of electric charge of the battery:

The symbol appears on the display together with other messages: The device is functional and the results remain accurate, but the battery shall be changed.

The symbol is displayed with "E-b", "A" and "Low" or message: The power is not sufficient to perform a diagnostics. Please, replace the battery.

Battery Changing

- To change the battery, make sure that the blood glucose meter is switched off.
- Pull the edge of the battery compartment cover and open it.
- Remove the old battery and replace it with new 1.5 V AAA alkaline battery.
- Close the battery compartment. In case of correct installation of the battery, you will hear the sound signal.

Battery changing does not affect the measurement results stored in the device memory.

As with all small parts, the battery shall be kept away from children. If swallowed, immediately seek medical assistance.

Batteries may leak chemicals if unused for a long time. Remove the battery for the hold time if the use of the blood glucose meter is not required for a long period (for example, 3 months or more).

Dispose the battery according to local legislation.

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Black

page 25

SYMBOLS USED BY THE SYSTEM FOR ERROR MESSAGE, TROUBLESHOOTING

The symbols used by the system for error message are provided below. If the error symbol appears on the blood glucose meter display, follow the instructions specified in the table below. If it is impossible to fix the problem, please contact the customer support service for assistance.

Test Result

SYMBOL	MEANING		
Lo	< 20 mg/dl (1.1 mmol/l)		
LOW ☹️	20–69 mg/dl (1.1–3.8 mmol/l)		
😊	AC 🍏	PC 👤	Gen
	70–129 mg/dl (3.8–7.1 mmol/l)	70–179 mg/dl (3.8–9.9 mmol/l)	70–119 mg/dl (3.8–6.6 mmol/l)
HIGH ☹️	AC 🍏	PC 👤	Gen
	130–239 mg/dl (7.2–13.2 mmol/l)	180–239 mg/dl (10–13.2 mmol/l)	120–239 mg/dl (6.6–13.2 mmol/l)
KETONE?	≥ 240 mg/dl (13.3 mmol/l)		
Hi	> 650 mg/dl (36.1 mmol/l)		

Error Messages

SYMBOL	MEANING	WHAT TO DO
E-b ⚠️ Low	Insufficient capacity of the battery	Replace the battery immediately
E-U ⚠️	Used test strip is inserted	Repeat the measurement with new test strip
E-t ⚠️	Air temperature is above/below the system operation range	The system operates at the temperature from 10°C to 40°C. Repeat the measurement when the temperature of blood glucose meter and test strip is within the specified range

25

page 26

E-E ⚠️	Device malfunction	Repeat the measurement using new test strip. If the device still does not work, please contact the customer support service for assistance
E-D ⚠️		
E-R ⚠️		
E-C ⚠️		
E-F ⚠️	Test strip is removed during the count down, or the blood sample volume is insufficient	Refer to the manual and repeat the measurement using new test strip. If the device still does not work, please contact the customer support service for assistance

In case of doubts regarding the obtained results, please refer to the instructions for use of test strips.

TROUBLESHOOTING

1. If the meter does not display the message after inserting a test strip

POSSIBLE CAUSE	WHAT TO DO
Low battery	Replace the battery
Test strip is inserted incorrectly or incompletely	Insert the test strip into the device with contact bars end first and facing up
Defective meter or test strips	Please contact the support service

2. The test does not start after inserting a test strip with blood sample

POSSIBLE CAUSE	WHAT TO DO
Insufficient blood sample	Repeat the measurement using new test strip and larger volume of blood sample
Defective test strip	Repeat the measurement using new test strip
The blood sample is applied after automatic switch-off (3 minutes after last action of the user)	Repeat the measurement using new test strip. Apply the sample only when 🍏 symbol appears on the screen
Defective blood glucose meter	Please contact the support service

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Black

page 27

3. Results of control solution test are out of the specified range

POSSIBLE CAUSE	WHAT TO DO
Error when performing the test	Read instructions thoroughly and repeat the test
Control solution vial was poorly shaken	Shake the control solution vial properly and repeat the test
Expired or contaminated control solution	Check the expiry date of the control solution.
Control solution is too warm or too cold	Control solution, blood glucose meter and test strips shall be stored at the room temperature from 20°C to 25°C before testing
Defective test strip	Repeat the test using new test strip
Blood glucose meter malfunction	Please contact the support service
Malfunction of blood glucose meter with test strip	Please contact the support service

27

page 28

SAFETY INFORMATION

The manufacturer does not guarantee the correct operation of SelfyCheck Prim Blood Glucose Meter in case of use of any test strips other than SelfyCheck Prim Test Strips. Use of other test strips with SelfyCheck Prim Blood Glucose Meter can lead to the erroneous results.
The product cannot be used until it has been tested with the SelfyCheck Prim Control Solution and the results obtained with the test are within the required range.
To avoid cross-contamination, it is not recommended to use the same lancing device by several people. It is prohibited to use the lancet that has been used by another person. To avoid infection, a new sterile lancet is to be used every time.
Any items that come in contact with blood, even after treatment, can be sources of vector-borne viral infections.

REQUIREMENTS FOR ENVIRONMENT PROTECTION

The medical device is safe and does not affect the environment.
Potential biohazard: doctors or people who use this product for analysis for several patients should be aware that any items that come in contact with blood, even after treatment, can be sources of vector-borne viral infections.
The product does not generate its own energy and does not absorb the user's energy. The product is not toxic at room temperature, does not release toxic substances into the environment, does not have any harmful effect on the human body through direct contact.
The product does not contain liquid and gaseous components in its composition.

QUALITY STANDARDS

The device corresponds to the international standards: IEC 62366:2007, BS EN 62366:2008, IEC 60601-1-6: 2010, IEC 60601-1-11:2015, ISO 14971, IEC 62366-1:2015, EN ISO 15197:2015, ISO 23640: 2011, EN ISO 17511, EN 13532:2002, EN ISO 15223-1:2014, GOST ISO 11137-2:2011, EN ISO 18113-4:15, EN ISO 18113-5, EN ISO 13485, CE MDD (93/42/EEC), GOST IEC 61010-1:2014, GOST IEC 61010-2:101-2013, IEC 61326-1.

STERILITY

The device components, except for lancets, are supplied non-sterile. The lancets are supplied sterile. Sterilisation method: radiation (cobalt-60γ-radiation). Sterilisation dose is 18,7 kGy.

METHODS AND MEANS OF CLEANING

It is expected that daily cleaning of the external surfaces of the blood glucose meter, lancing device, storage case will be performed in case of contaminations.

Blood glucose meter

Cleaning: The external surface of the blood glucose meter shall be cleaned using a cloth wetted with main water or mild detergent and wiped with a soft dry cloth; the device shall not be rinsed under water or cleaned using the organic solvents.

Disinfection: in case of individual use at home disinfection is not required;

Lancing Device


Cleaning: the lancing device is cleaned using mild soap and water; the device shall not be rinsed under water or washed in a dish washer or using detergents.

Disinfection: in case of individual use at home disinfection is not required;

Test strips and lancets are for single use only, therefore their disinfection and cleaning is not provided. For control solutions cleaning or disinfection is not required.

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page 29

STORAGE CONDITIONS

Name	Storage conditions
SelfyCheck Prim Blood Glucose Meter	From -20°C to 60°C, not more than 95% of relative humidity.
SelfyCheck Prim Test Strips	From 2°C to 30°C, not more than 85% of relative humidity. Store the device in original container only, do not transfer to another package. In a dry cool place. Avoid direct sunlight and heat. Do not freeze. Keep the container with test strips away from children.
SelfyCheck Prim Control Solutions (Low, Norm, High)	From 2°C to 30°C, relative humidity from 30% to 85%. Store the control solution in the tightly closed vial. Do not freeze.
Lancing device (with AST tip)	From -10°C to 50°C, not more than 95% of relative humidity.
SelfyCheck lancets	From -20°C to 60°C, not more than 95% of relative humidity.

Storage Conditions for Kits (Complete Sets)

Set	Storage conditions
Version 1	From -10°C to 50°C, not more than 95% of relative humidity
Version 2	From 2°C to 30°C, not more than 85% of relative humidity
Version 3	From 2°C to 30°C, not more than 85% of relative humidity
Version 4	From -20°C to 60°C, not more than 95% of relative humidity
Version 5	From 2°C to 30°C, from 30% to 85% of relative humidity

Warranty Shelf Life

Guaranteed storage life ("Use by") shall be specified on the package of the sets and accessories for each element of the set.

Name	Warranty Shelf Life
SelfyCheck Prim Blood Glucose Meter	5 years
SelfyCheck Prim Test Strips	21 months
SelfyCheck Prim Control Solutions (Low, Norm, High)	24 months
Lancing device (with AST tip)	2 years
SelfyCheck lancets	5 years after radiation sterilisation

29

page 30

OPERATING CONDITIONS

Name	Operating conditions
SelfyCheck Prim Blood Glucose Meter	From 10°C to 40°C, not more than 85% of relative humidity. Avoid falling and heavy shocks. Prevent exposure to the direct sunlight and high humidity. Use of this device in dry environment, especially with synthetic materials (synthetic clothes, carpets, etc.), can lead to the damage by static discharge and, therefore, to the erroneous results. It is prohibited to use the blood glucose meter in close proximity to the sources of strong electromagnetic radiation, as it can affect the accuracy of device operation. The device shall not be used in case of improper operation or damage.
SelfyCheck Prim Test Strips	From 10°C to 40°C, not more than 85% of relative humidity. The container cover shall be tightly closed immediately after taking out of one test strip. Touch the test strip only with clean and dry hands. Use each test strip immediately after taking it out of the container. Never bend, cut or damage the test strip. Do not use the test strip after expiration of the shelf life. It may affect accuracy of the results.
SelfyCheck Prim Control Solutions (Low, Norm, High)	From 20°C to 25°C, from 30% to 85% of relative humidity. Before the testing, make sure that control solution, glucose meter and the test strips are within the set temperature range.
Lancing device (with AST tip)	From -10°C to 50°C, not more than 85% of relative humidity.
SelfyCheck lancets	From 10°C to 40°C, not more than 85% of relative humidity. Do not use if the protective cap is removed.

In case of violation of operation conditions specified by the manufacturer, protection used in the device may be damaged.

TRANSPORTATION CONDITIONS


Transportation can be carried out using all types of transport in covered vehicles or general-purpose containers protected from sunlight in accordance with cargo transportation rules applicable to them. The device is transported with all storage conditions observed.

PROCEDURE FOR DISPOSAL AND DESTRUCTION OF THE MEDICAL DEVICE

The device is disposed and destroyed in accordance with the classification, rules of collection, usage, decontamination, placement, storage, transportation, record keeping and disposal of medical waste established by the competent authority of the importing country. Precautions and local regulations shall be observed in case of disposal of the used batteries. It is prohibited to throw batteries into a fire. Used test strips and lancets can be infectious, if the user has contagious infectious diseases. It is crucial that the used materials are disposed of as infectious and potentially biohazardous wastes of Class B.

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page 31

The lancet shall be put into the sharps container.
The blood glucose meter is subject to mechanical corrosion with disposal as the production or household waste.
Upon usage, expiration of shelf life and in case of packaging damage and impossibility of usage of the devices as intended, such devices should be discarded under SanPin 2.1.7.2790 as waste of Class B. Unutilized devices are referred to class A waste, they do not contain harmful substances and can be discarded with household waste in accordance with SanPin 2.1.7.2790-10.

WARRANTY OBLIGATIONS AND SERVICE LIFE

Compliance of SelfyCheck Portable Blood Glucose Monitoring System (Blood Glucose Meter) to the requirements established in the specification is ensured for the entire period of guarantee, provided that the device is used as intended, and its use, storage and transportation is carried out according to the indicated conditions.

Product	Period of guarantee
SelfyCheck Prim Portable Rapid-Response Analyser (Blood Glucose Meter) for measurement of blood glucose level	For the period of the guaranteed storage life.
SelfyCheck Prim Test Strips	For the period of the guaranteed storage life.
SelfyCheck Prim Control Solutions (Low, Norm, High)	For the period of the guaranteed storage life.
Lancing device (with AST tip)	For the period of the guaranteed storage life.
SelfyCheck lancets	For the period of the guaranteed storage life.

Product	Guaranteed storage life
SelfyCheck Prim Blood Glucose Meter	5 years from first usage
SelfyCheck Prim Test Strips	6 months after first opening
SelfyCheck Prim Control Solutions (Low, Norm, High)	90 days after first opening
Lancing device (with AST tip)	2 years from first usage
SelfyCheck lancets	5 years after radiation sterilisation

31

page 32

Appendix A. Electromagnetic Compatibility

Electromagnetic compatibility is ensured by compliance with the requirement of GOST R IEC 60601-1-2-2014.

ATTENTION!

DEVICES require special measures to ensure EMC and shall be installed and put into operation according to the information on EMC provided in the operating documentation.

ATTENTION!

Application of mobile radio-frequency communication means can affect the **DEVICE**.

ATTENTION!

Use of accessories, converters and cables, not specified in the list, except for converters and cables supplied by the manufacturer of **DEVICE** as replacement parts for internal components, can lead to the increase of electromagnetic emission or reduction of noise immunity of the **DEVICE**.

Table B.1

Manual and Declaration of Manufacturer. Electromagnetic Emission		
DEVICE is designed for use in electromagnetic environment defined below. The buyer and user of the DEVICE shall ensure the application in the specified environment		
Test for electromagnetic emission	Compliance	Electromagnetic environment – instructions
Radio frequency interference as per CISPR 11	Group 1	DEVICE uses radio-frequency energy only for performance of internal functions. The level of radio frequency interference emission is low and probably will not lead to the malfunction of the electronic equipment located nearby.
Radio frequency interference as per CISPR 11	Class B	DEVICE is suitable for use in any locations, including domestic houses and buildings directly connected to the power distribution network ensuring supply for the domestic houses
Voltage fluctuations and flicker as per IEC 610003-3	Fail	

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page 33

Table B.2

Manual and Declaration of Manufacturer. Immunity			
DEVICE is designed for use in electromagnetic environment defined below. The buyer and user of the DEVICE shall ensure the application in the specified environment			
Immunity test	Test level according to IEC 60601	Compliance level	Electromagnetic environment – instructions
Electrostatic discharge (ESD) as per IEC 61000-4-2	±6 kV – contact discharge	±6 kV – contact discharge	The floor shall be made of wood, concrete or ceramic tile. If the floor is covered with synthetic material, the relative humidity shall be minimum 30%
	±8 kV – air discharge	±8 kV – air discharge	
Fast transient burst as per IEC 61000-4-4	±2 kV for power lines	±2 kV for power lines	Quality of electrical energy in the electrical network of the building shall correspond to the typical conditions of the commercial or hospital environment
	±1 kV for input/output lines	±1 kV for input/output lines	
Microsecond high energy pulse disturbance as per IEC 61000-4-5	±1 kV – in case of generation of disturbances via "wire-wire" circuit	±1 kV – in case of generation of disturbances via "wire-wire" circuit	Quality of electrical energy in the electrical network of the building shall correspond to the typical conditions of the commercial or hospital environment
	±2 kV – in case of generation of disturbances via "wire-ground" circuit	±2 kV – in case of generation of disturbances via "wire-ground" circuit	
Dips, interruptions and voltage variations in the input power lines as per IEC 61000-4-11	< 5% U_N (voltage dip > 95%) during 0,5 of the period	< 5% U_t (voltage dip > 95%) during 0,5 of the period	Quality of electrical energy in the electrical network of the building shall correspond to the typical conditions of the commercial or hospital environment. If the user of Device requires the continuous operation under the conditions of the circuit voltage interruptions, it is recommended to ensure the power supply of the Device from the battery or uninterruptible power supply
	40% U_N (voltage dip of 60%) during five periods	40% U_t (voltage dip of 60%) during five periods	
	70% U_N (voltage dip of 30%) during 25 periods	70% U_t (voltage dip of 30%) during 25 periods	
	< 5% U_N (voltage dip > 95%) for 5 s	< 5% U_t (voltage dip > 95%) for 5 s	
Power frequency magnetic field as per 61000-4-8	3 A/m	3 A/m	The levels of the power frequency magnetic field shall correspond to the typical conditions of the commercial or hospital environment.

Note: U_N — AC network voltage prior to application of the test level.

33

page 34

Table B.3

Manual and Declaration of Manufacturer. Immunity			
DEVICE is designed for use in electromagnetic environment defined below. The user shall ensure that the DEVICE is used in such environment			
Immunity test	Test level according to IEC 60601	Compliance level	Electromagnetic environment – instructions
Conducted disturbances induced by RF fields as per IEC 61000-4-6	3 V (mean square value) from 150 kHz to 80 MHz out of frequency range for industrial, scientific and medical high-frequency devices ^{a)}	3 V	The distance between used portable/mobile radio facility and any part of DEVICE , including the cables, shall not be less than recommended separation distance calculated according to the presented formula in respect to the transmitter frequency. Recommended separation distance is: $d=1,2 \cdot \sqrt{P}$
Radio-frequency electromagnetic field as per IEC 61000-4-3	3 V/m from 80 MHz to 2.5 GHz	3 V/m	$d=1,2 \cdot \sqrt{P}$ (from 80 MHz to 800 MHz) $d=2,3 \cdot \sqrt{P}$ from (from 800 MHz to 2.5 GHz) where d – recommended separation distance, m ^{b)} ; P – nominal value of the maximum output in watts according to the value, set by the manufacturer. Intensity of the field in case of radiowave propagation from the fixed radio transmitters according to the results of observations over the electromagnetic environment ^{a)} shall be lower than the level of compliance in each frequency band ^{b)} . The disturbances can be observed near the equipment, marked with sign


- a) Intensity of the field in case of radiowave propagation from the fixed radio transmitters, such as base stations of radio telephone networks (cellular/wireless) and mobile ground radio stations, amateur radio stations, AM and FM broadcast transmitters, TV transmitters cannot be determined by calculation with adequate accuracy. To do this, it is necessary to perform the practical measurements of field intensity. If the measured values in the **DEVICE** location exceed the applicable levels of compliance, the observations over the **DEVICE** operation shall be carried out in order to check the proper functioning. If the deviation from the normal operation is detected in the course of observation, it is necessary to take additional measures, such as reorientation or movement of **DEVICE**.
- b) The field intensity of less than 3 V/m shall be provided outside the frequency band from 150 kHz to 80 MHz.

Notes

- On frequencies of 80 and 800 MHz, the greater value of the field intensity shall be used.
- The presented formulas are not always applicable. The electromagnetic wave propagation is impacted by absorption or reflection from the structures, objects and people.

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page 35


Table B.4

Recommended values of the separation distance between portable and mobile radio-frequency communication means and Device			
DEVICE is designed for use in the electromagnetic environment, in case of which the radiated disturbance level control is ensured. The buyer or user of the DEVICE can avoid the electromagnetic noise effect by means of ensuring the minimum separation distance between the portable and mobile radio-frequency communication means (transmitters) and DEVICE, as recommended below, taking into account the power output of the communication means			
Nominal maximum output of the transmitter, W	Separation distance, m, depending on the transmitter frequency		
	$d=1.2 \cdot \sqrt{P}$ in the band from 150 kHz to 80 MHz	$d=1.2 \cdot \sqrt{P}$ in the band from 80 kHz to 800 MHz	$d=2.3 \cdot \sqrt{P}$ in the band from 800 MHz to 2.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
Notes 1. On frequencies of 80 and 800 MHz, the greater value of the field intensity shall be used. 2. When determining the recommended values of the separation distance for transmitters with nominal maximum power output, not specified in the table, the nominal maximum output in watts, indicated in the documentation of transmitter manufacturer, is inserted into the presented formulas. 3. The presented formulas are not always applicable. The electromagnetic wave propagation is impacted by absorption or reflection from the structures, objects and people.			

TaiDoc Technology Corporation
B1-7F, No.127, Wugong 2nd Rd., Wugu Dist.,
24888 New Taipei City, Taiwan
www.taidoc.com

For self-testing

Address of the authorized representative:

NEXT Bio LLC
5 Akademika Pavlova Street, lit. B, room 46-H, 197022 Saint Petersburg, Russia

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