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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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IS TO CERTIFY THAT: the seal of HANGZHOU BIOTEST
TECH CO., LTD. on the annexed OPERATION MANUAL is
genuine.

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Operation manual

Reagent kit "Rapid test for qualitative detection of D-dimer in human plasma, serum or whole blood using immunochromatographic methods" (D-dimer rapid test cassette)

INTENDED USE

A kit "Rapid test for qualitative detection of D-dimer in human plasma, serum or whole blood using immunochromatographic methods" (D-dimer rapid test cassette) is used for qualitative detection of D-dimer in human plasma, serum or whole blood in order to facilitate diagnosis of deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC) by immunochromatographic method.

For a single use of the kit for the intended purpose.

For the in vitro clinical laboratory diagnosis.

Potential users of the product

Laboratory units of healthcare centers.

Professional level of potential users

The doctor of clinical laboratory diagnostics, medical laboratory technician.

Indications

At the diagnosis of deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC).

Contraindications

When using by specially trained personnel and taking into account the intended use contraindications are not identified.

Possible adverse effects

When using by specially trained personnel and taking into account the intended use adverse effects are not identified.

Description of the target analyte, information on its scientific validity

D-dimer is a biomarker that generally indicates the activation of hemostasis and fibrinolysis. It is the product of the disruption of fibrin, which is formed during the disintegration of cross-linked fibrin under the influence of plasmin-induced fibrinolytic activity. Fibrin D-dimer being the marker of fibrin forming

ruption is widely used in clinical assessment procedures to detect activation of the blood coagulation system. Since plasma levels of D-dimer are elevated after the clot formation, D-dimer measurement is usually used in combination with clinical parameters during the initial evaluation in case of suspected acute venous thromboembolism (VTE). Elevated levels of D-dimer are prognostic criteria in various conditions, including venous thromboembolism, disseminated intravascular coagulation, cardiovascular disease, infectious diseases and cancer. The concentration of D-dimer normally ranges from 0 to 500 ng/ml. Rapid test for qualitative detection of D-dimer in human plasma, serum or whole blood using immunochromatographic methods (D-dimer rapid test cassette) is used for qualitative determination of D-dimer in specimens of human whole blood, serum or plasma. The test uses monoclonal antibodies for selective determination of elevated D-dimer levels in human whole blood, serum, or plasma.

Kit composition and complete set

Test cassette	A rectangular plate 73x25x5,5 made of white plastic with two windows, inside of which goat anti-mouse IgG and D-dimer-McAT were applied to the membrane and dried on the control line; and monoclonal antibodies against D-dimer conjugated with gold were applied on one of the filters; the plate is packaged hermetically in an individual laminated foil package with desiccant ¹ .	10 units
Buffer solution for sample dilution	Solution containing sodium chloride 0.02%/sodium phosphate basic/casein, sodium salts/ kanamycin sulfate 0.025%; sodium azide 0.02% transparent colorless liquid, pH 7.4+/-0.1	1 unit (3,0 ml)
Disposable pipettes	Pipette of colorless transparent plastic with balloon-like extension at one end	10 units

Note:

- One package of laminated foil includes 1 test cassette and 1 disposable pipette.
- There are no materials in the device that come into direct or indirect contact with the patient's body and personnel using the device, while meeting the requirements of the operational documentation (instructions for use).
- The device does not contain medicines and pharmaceutical substances.

¹ including disposable pipette

The kit components are packed in a box, the operation manual is enclosed in the box.

The package includes: a kit of reagents, operation manual, certificate.

GENERAL CONSUMER PROPERTIES

The basic version of the kit allows the study of 10 samples.

Sterilization of the device

The device does not require the sterilization.

Software for device work

Absent.

Maintenance and repair of the device

ice does not require maintenance and is not repairable.

PRINCIPLE OF OPERATION

ner rapid test cassette is an express test for carrying out a chromatographic immunoassay that makes an anti-D-dimer antibody conjugated with colloidal gold particles, an antibody to the D-dimer on test line and an IgG-class antibody on the control line. During testing, samples of whole blood, serum or plasma react with particles coated with an antibody to the D-dimer coated on the test strip. The mixture then migrates upward on the membrane by capillary action and reacts with the antibody to the D-dimer applied in the region of the test line. If the sample contains a D-dimer, a single purple line appears in the test line region. The presence of this purple line indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a purple line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

SPECIMEN COLLECTION AND PREPARATION

25 mcl of human serum or plasma (K_2 EDTA, sodium citrate 3.2%, sodium heparin, lithium heparin) or 50 mcl of whole blood.

❖ For collection fingertip whole blood specimens:

- Wash the patient's hand with soap and warm water or clean with an alcohol swab. Allow to dry.
- Massage the hand without touching the puncture site by rubbing down the hand towards the fingertip of the middle or ring finger.
- Puncture the skin with a sterile lancet. Wipe away the first sign of blood.
- Gently rub the hand from wrist to palm to finger to form a rounded drop of blood over the puncture site.
- Place a sample of fingertip whole blood to the test cassette using a capillary tube:

Touch the tip of the capillary tube to a drop of blood to fill the tube with approximately 50 mcl. Avoid the air bubbles formation.

Place the rubber bag in the top of the capillary tube, then squeeze the bag to expel the whole blood on the specimen well of test-cassette.

Separate serum or plasma from blood immediately to avoid hemolysis. Use only clear, non-hemolyzed specimens.

Testing should be performed immediately after specimen collection. Do not leave the specimens at room temperature for a long time. Serum and plasma specimens may be stored at 2-8°C for up to 3 days. For long-term storage, specimens should be kept below -20°C. Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 2 days of collection. Do not freeze whole blood specimens. Whole blood collected from fingertip should be tested immediately.

Bring specimens to room temperature prior to testing. Frozen specimens must be completely thawed and mixed well. Specimens should not be frozen and thawed repeatedly.

If specimens are to be shipped, they should be packed in compliance with local regulations covering the transportation of etiologic agents.

Collection of blood samples should be performed in accordance with good practice by structure (ГОСТ Р [Russia State Standard] 53079.4-2008 "Ensuring the quality of clinical laboratory tests").

ANALYTICAL AND DIAGNOSTIC CHARACTERISTICS

Analytical sensitivity (the minimum detectable concentration) is 380 ng / ml.

Sensitivity of the kit is determined from the reference positive biological samples containing the D-dimer, as a percentage of the samples determined by the kit as positive, and is 100%.

Specificity of the kit is determined from control negative biological samples that do not contain D-dimer, as a percentage of samples identified by the kit as negative, and is 100%.

PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after expiration date.
- The test device must remain in the sealed pouch until use.
- Do not eat, drink or smoke in the area where the specimens or kits are handled.
- Handle all specimens as if they contain infectious agents. During all procedures observe established precautions against microbiological hazards, and follow the standard procedures for proper disposal of specimens.
- Wear protective clothing such as laboratory coats, disposable gloves or eye protection when specimens are being tested.
- The used test, specimens and contaminated materials must be discarded according to local regulations.
- Humidity and temperature can adversely affect results.
- Do not mix components from different lots.
- The kit is biologically safe; however, it is necessary to deal with the samples as a potentially infectious material.
- Reagent kit, in accordance with Sanitary Regulations and Standards 2.1.7.2790-10 "Sanitary and epidemiological requirements for the handling with medical waste" refers to Class A - waste that does not have contact with the biological fluids of patients, infectious patients. Class A - epidemiologically safe waste, approximated in composition to solid domestic waste (SDW).
- Recovery or destruction, and disinfection of reagent kits should be carried out in accordance with Sanitary Regulations and Standards 2.1.7.2790-10 "Sanitary and epidemiological requirements for the handling with medical waste" and Methodical Guidelines 287-113 "Guidelines for disinfection, pre-sterilization cleaning and sterilization of medical devices".
- The list of risks identified during the risk analysis process, and methods for risk management are presented in a separate document (see Risk Analysis).

DIRECTIONS FOR USE

Materials required but not provided

Containers for collection of plasma, serum or venous blood specimens

ary centrifuge

lancets (for fingertip whole blood only)

rinized capillary tubes and dosing rubber bag (for fingertip whole blood only).

ow the test device, specimen, buffer solution and/or controls to reach room temperature (15-30°C) or to testing.


Bring the pouch to room temperature before opening it. Remove the test device from the sealed pouch and use it for an hour.

2. Place the test device on a clean and level surface.

For serum or plasma specimens: Hold the dropper vertically and transfer 1 drop of serum or plasma (approximately 25 mcl) to the specimen well then add a drop of buffer (approximately 40 mcl) and start the timer. See figure below.

For venipuncture whole blood specimens: Hold the dropper vertically and transfer 2 drops of venipuncture whole blood (approximately 50 mcl) to the specimen well and a drop of buffer (approximately 40 mcl), then start the timer. See figure below.

For fingertip whole blood specimens: Using capillary tube: Fill the capillary tube and transfer approximately 50 mcl of fingertip whole blood to the specimen well in the test device, add 1 drop of buffer (approximately 40 mcl) and start the timer. See figure below.

3. Wait for the purple line(s) to appear. Read results at 15 minutes. Do not interpret results after 20 minutes.

INTERPRETATION OF RESULTS

(Refer to the figure above)

POSITIVE RESULT:* Two distinct lines appear. One purple line should be in the control line region (C) and another apparent purple line should be in the test line region (T).

The intensity of the color in the test line region (T) will vary depending on the concentration of in the specimen. Therefore, any shade of color in the test line region (T) should be considered.

INVALID RESULT: The appearance of one purple line in the control line region (C). No line appears in the test line region (T).

INVALID RESULT: Control line (C) fails to appear. Insufficient specimen volume or incorrect procedural techniques are the most likely reasons for control line failure. Review the procedure and repeat the test with a new test cassette. If the problem persists, discontinue using the test kit immediately and contact your local distributor.

QUALITY CONTROL

An internal procedural control is included in the test. A purple line appearing in the control line region (C) is considered as an internal tool for monitoring the test procedure. It confirms adequate membrane wicking.

Control standards are not supplied with this kit; however it is recommended that positive and negative controls be tested as a good laboratory practice to confirm the test procedure and to verify proper test performance. Samples for external quality control should be tested in accordance with the quality control standards established in your laboratory.

LIMITATIONS

1. D-dimer rapid test cassette is for the in vitro diagnostic use only. This test should be used for the detection of D-dimer in whole blood, serum or plasma specimens only. Neither the quantitative value nor the rate of increase in D-dimer level can be determined by this qualitative test.
2. D-dimer rapid test cassette will only indicate the qualitative level of D-dimer in the specimen and should not be used as the sole criteria for the diagnosis of deep vein thrombosis (DVT), pulmonary embolism (PE) and disseminated intravascular coagulation (DIC).
3. To confirm the obtained results, additional specimen analyzes should be performed in accordance with the guidelines of local health authorities, for example, enzyme immunoassay.
4. As with other diagnostic tests, all results must be interpreted together with other clinical data available to the physician.
5. This test is for screening purposes only.
6. A negative test result does not exclude the possibility of thrombosis. Therefore, in order to establish an accurate diagnosis, the result obtained with the help of D-dimer rapid test cassette should be used in conjunction with clinical data.
7. False positive result is possible in cases of liver disease, high rheumatoid factor, pregnancy, injury.
8. False negative results are found in cases where the analysis is taken too quickly after the thrombus formation, in the presence of anticoagulants in the blood.

Attention!

This device will give the real result only in the case of strict compliance with the instructions for use.

EXPECTED VALUES

Results of D-dimer rapid test cassette using have been compared to the results of leading specially available rapid tests (products that are trusted in the market), and during the comparison, overall test accuracy was confirmed as 97%.

PERFORMANCE CHARACTERISTICS

Sensitivity and specificity of samples obtained from apparently healthy blood donors, as well as of samples obtained from patients

D-dimer rapid test cassette was validated using samples obtained from apparently healthy blood donors, as well as by means of samples obtained from patients. The results show that the sensitivity of the D-dimer rapid test cassette is 97.4%, and the specificity is 98.1% in compare to other rapid tests.

Method		Other rapid test		Total
D-dimer rapid test cassette	Results	Positive result	Negative result	results
	Positive result	76	2	78
	Negative result	2	101	103
Total results		78	103	181

Precision

For the results of the same assay

Intra-serial precision was determined for 15 series of the following three samples: levels of D-dimer in samples were 0 mcg / ml, 0.5 mcg / ml and 5 mcg / ml. The specimens were correctly identified in more than 99% of cases.

For the results of different assays.

Inter-serial precision was determined for 3 independent assays using the same three samples: 0 mcg / ml, 0.5 mcg / ml and 5 mcg / ml of D-dimer. Three different batches of D-dimer rapid test cassette (whole blood/serum/plasma) were tested using these samples. Samples were correctly identified in more than 99% of cases.

Cross reactivity

D-dimer rapid test cassette was tested with the following positive samples: fibrinogen, fibrin monomer, 5 types of hepatitis B virus (HBsAg, HBsAb, HBeAg, HBeAb, HBcAb), syphilis, antibodies to hepatitis C virus, antibodies to HIV and antibodies to H.pylori. According to the results of the testing cross-reactivity was not revealed.

Interfering substances

The following potentially interfering substances were added to the negative and positive samples for the D-dimer.

Acetaminophen: 20 mg/dL	caffein: 20 mg/dL	bilirubin: 1000 mg/dL
Acetylsalicylic acid: 20 mg/dL	creatine: 200 mg/dL	Oxalic acid: 600 mg/dL
Gentisinic acid: 20 mg/dL	Ascorbic acid: 20 mg/dL	cholesterol: 800 mg/dL

hemoglobin: 1000 mg/dL

Albumin: 10,500 mg/dL

triglycerides: 1600 mg/dL

Interfering to the results of the testing no one substance in said concentrations is interfering.

EXPIRE DATE

The shelf life of the kit is 2 years.

A kit with expired shelf life is not a subject to application.

STORAGE AND STABILITY

Storage (stability)

Storage of the kit in the manufacturer's packaging must be at room temperature or in the refrigerator (at a temperature 2 - 30 ° C) ensuring a regulated temperature regime with a daily temperature record throughout the shelf life. DO NOT FREEZE.

The test device is stable through the expiration date printed on the sealed pouch.

Storage (stability) after opening

Opened buffer solution must be stored at room temperature or in the refrigerator (at a temperature 2 - 30 ° C) during the entire shelf life indicated on the package. DO NOT FREEZE. The contents of one individual package for single use. The kit is stable until the expiration date.

Storage (stability) imitation of shipment conditions

It should be transported in the original packaging of the manufacturer, which protects the product from external influences, on all modes of transport in covered vehicles in accordance with the rules of transportation operating in this type of transport at a temperature of 2 to 30 ° C for the entire shelf life. When carrying out loading and unloading operations and transportation, the requirements of handling marks on packaging and shipping containers should be strictly observed. If the transportation and storage conditions are not observed, the medical device cannot be used!

SALE CONDITIONS

For health care institutions.

MANUFACTURER WARRANTY

The production is certified for compliance with the requirements of international standards Directive 98/79 / EC, ISO 13485 quality management system.

For questions related to the quality and handling of the reagent kit "Rapid test for the qualitative determination of D-dimer in human plasma, serum and whole blood by immunochromatography" (D-dimer rapid test cassette), please contact: Limited Liability Company "ROTANA", OOO "ROTANA", 117342, Russia, Moscow, General Antonov str., 3A, tel. / Fax: +7 (495) 646-28-61, med@rotana-rf.ru

Manufacturer:

Biotech Co., Ltd, China

Hangzhou, 311121-P.R. China

+86-571-88630899

+86-571-88633388

Authorized representative of the manufacturer in the territory of the Russian Federation

United Liability Company "ROTANA", OOO "ROTANA"

17342, Russia, Moscow, General Antonov str., 3A, tel. / Fax : +7 (495) 646-28-61, med@rotana-rf.ru

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