The «03» of the August, 2020г.

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(signed at)

Printed Name: Alan Zettin

Title with Diaco Biofarmaceutici S.r.l.: Chief executive officer

DIACO BIOFARMACEUTICI S.R.L.
Via Flavia, 124 - 34147 Trieste
C.F. e P.IVA 01221530320

INSTRUCTION FOR USE

Device:

Injectable hyaluronic acid-based implants XELA REDERM, approaches

1. Injectable hyaluronic acid-based implant XELA REDERM 1,1% – in the 1 ml syringe - 1 pce, package contents with needles 30G ½" (0,3 x 13 mm) - 2 pcs, tracking labels - 2 pcs;
2. Injectable hyaluronic acid-based implant XELA REDERM 1,1% – in the 2 ml syringe. - 1 pce, package contents with needles 30G ½" (0,3 x 13 mm) - 2 pcs., tracking labels - 2 pcs;
3. Injectable hyaluronic acid-based implant XELA REDERM 1,8% - in the 1 ml syringe. - 1 pce, package contents with needles 30G ½" (0,3 x 13 mm) - 2 pcs., tracking labels - 2 pcs;
4. Injectable hyaluronic acid-based implant XELA REDERM 1,8% - in the 2 ml syringe. - 1 pce, package contents with needles 30G ½" (0,3 x 13 mm) - 2 pcs., tracking labels - 2 pcs;
5. Injectable hyaluronic acid-based implant XELA REDERM 2,2% - in the 1 ml syringe. - 1 pce, package contents with needles 30G ½" (0,3 x 13 mm) - 2 pcs, tracking labels - 2 pcs;
6. Injectable hyaluronic acid-based implant XELA REDERM 2,2% – in the 2 ml syringe. – 1 pce, package contents with needles 30G ½" (0,3 x 13 mm) – 2 pcs, tracking labels – 2 pcs;

Composition:

Active substance: sodium hyaluronate.

Other components: water for injection, sodium succinate, succinic acid, sodium chloride

The contents of the syringes are sterilized by steam.

The needles are sterilized with ethylene oxide.

Package contents:

- A syringe containing a sterile injectable implant XELA REDERM 1 ml or 2 ml - 1 piece;
- Needles for hypodermic injections (0.3 x 13 mm), 30G ½" manufactured by Becton, Dickinson and Company Limited – 2 pcs;
- Tracking labels – 2 psc. for tracing of injections to the patient. The tracking labels contain device information (lot number –(LOT) and expiration date – (EXP).

Device description:

Hyaluronic acid is a natural polysaccharide and an important building block for the skin and connective tissues.

Injectable implant XELA REDERM -is a gel, colorless, transparent, consisting of a sterile, apyretic hyaluronic acid of non-animal origin with optimal pH.

Indications:

Implants are intended to restore skin moisture balance, improve skin structure and the elasticity of the skin, and also for correction of facial wrinkles and folds, improving skin color. Implants should be injected into the dermal layer of the skin, preferably into the deeper parts of dermis.

Contraindications:

XELA REDERM should not be used:

- for other routes of administration, except for intradermal;
- in combination with peeling, laser or ultrasound procedures;
- if the patient has skin diseases, inflammation, infected areas in the procedure area or areas adjacent to it;
- with a known hypersensitivity of the patient to hyaluronic acid, if there is a history of allergy or anaphylactic shock;
- during pregnancy and breastfeeding;
- aged under 18 years;
- do not inject into blood vessels;
- when taking high doses of vitamin E, aspirin, anti-inflammatory products and anti-coagulants, less than a week before the procedure;
- in the case of autoimmune diseases of the patient;
- for skin areas where the procedures of with implants (fillers) have been carried out, and biodegradation of this implants (filler) has not yet occurred.

Precautions for Use:

1. Do not inject intravenously.
2. Do not use if the color and clarity of the solution deviate from those declared by the manufacturer.
3. Do not use if package integrity is compromised.
4. Use immediately after opening the package; Dispose the remaining unused implant.
5. Do not reuse. Only for single use in one patient during one procedure. After the first injection, the implant remaining in the syringe is no longer sterile. The remnants of the implant cannot be reused after the first use, due to a violation of sterility.
6. Do not re-sterilize. Repeated sterilization can lead to cross-contamination from the patient and the doctor.
7. Do not allow air to enter the syringe.
8. Do not freeze.
9. It is not recommended to warm up in a microwave oven.
10. Do not take orally.
11. To keep from taking high doses of vitamin E, aspirin, anti-inflammatory products and anticoagulants within one week prior to the procedure.
12. Do not use makeup for at least 12 hours after the procedure.
13. Avoid sudden changes in temperature (strong cold, sauna, etc.) within 72 hours after the procedure.

Special instructions:

1. When using an injection implant, you must follow the aseptic techniques.
2. Local anesthesia in the injection areas is allowed.
3. Do not heat the injection site.
4. Avoid exposure to laser, chemical peelings or other similar means that can cause an inflammatory reaction at the injection site, until the skin is completely normal after the procedure.
5. Since interaction with other known injectable products has not been investigated, XELA REDERM should not be used together with other injections in the perimeter or in the immediate area of the procedure.

Side effects:

Application of an injectable implant can be unpleasant for patients, as well as any other applications of injections.

After the procedure, rashes on the skin (redness, swelling, bruising, itching) and minor soreness at the injection sites may occur. These symptoms usually disappear within 72 hours after the injection.

In rare cases, there may be an injection site mass with a loss of skin sensitivity.

In the literature, extremely rare cases of necrosis, papilliform swelling of the tissue, formation of granulomas, hypersensitivity reactions and the formation of abscess after the injection of products containing hyaluronic acid are described.

If any of these undesirable effects occur, you should immediately contact your doctor and notify the authorized representative in the Russian Federation or the manufacturer of the injectable implant.

Method of administration and dosage:

XELA REDERM injectable implants should be used only by medical specialists who have been special trained in the safe and effective use of injectable dermal preparations. Before the procedure, the patient's history should be carefully read and the patient should be warned about the expected outcome of the procedure and the possible undesirable effects of the product.

The choice of concentration depends on the degree of skin photoaging and clinical picture of aging.

Degree of skin photoaging according to R. Glogau	Age	Morphotypes of aging according to I.I. Kolgunenko		
		Tired	Deformative	Fine-wrinkled
Degree I (minor damage)	20–40 yo	XELA REDERM 1.1% XELA REDERM 1.8%	XELA REDERM 1.1%	XELA REDERM 1.8%
Degree II (moderate damage)	40–50 yo	XELA REDERM 1.1% XELA REDERM 1.8%	XELA REDERM 1.1% XELA REDERM 1.8%	XELA REDERM 1.8%
Degree III (pronounced damage)	50–60 yo	XELA REDERM 1.1% XELA REDERM 1.8% XELA REDERM 2.2%	XELA REDERM 1.1% XELA REDERM 1.8%	XELA REDERM 1.8% XELA REDERM 2.2%
Degree IV (very pronounced damage)	Over 60 yo	XELA REDERM 1.1% XELA REDERM 1.8% XELA REDERM 2.2%	XELA REDERM 1.1% XELA REDERM 1.8% XELA REDERM 2.2%	XELA REDERM 1.8% XELA REDERM 2.2%

In oily and combined skin, as well as in very dry and thin skin, it is recommended to use XELA REDERM 1.1%, and XELA REDERM 1.8 and 2.2% in dry skin.

Method of administration of XELA REDERM: intradermal.

Zones of administration: any areas of face, arms, neck, and décolleté.


The following methods (injection techniques) are used to inject the implants for all zones:

- **Papule technique:** the product is administered by microinjections until small papules with a diameter of 1 - 2 mm occur at a distance of 0.5 - 1 cm from each other. The needle is introduced into the upper dermal layer at a depth of 2 - 3 mm, with the needle cut facing upwards;

- **Line technique (tunnelling and tracing)** the whole length of the needle is inserted parallel to the skin at a depth of 2 - 3 mm. The needle is moved along a straight line (tunnelling technique) or from side to side (tracing technique). The needle is inserted into the derma then slowly retracted, injecting the product gradually into the skin.

As a rule, it is necessary to perform 3-5 procedures with an interval of 1 - 2 weeks between procedures. Regular doses of additional injections will maintain the desired long-term effect.

Attention! The injection technique must be selected only by a medical specialist performing similar procedures.

Recommendations for use and procedure:

1. Check the expiration date on the package. Do not use the product after the expiration date.
2. Check the integrity of the package. Do not use if the package is damaged.
3. Attach the tracking label (with the product data: batch code and expiration date) to the patient's medical record, to whom the procedure is carried out.
4. Open the sterile package.
5. Remove the plug from the tip of the syringe.
6. Put a hypodermic needle on the syringe. Make sure that the needle fits tightly on the syringe.
7. Perform the injections and follow the procedure according to the plan.
8. At the end of the procedure, dispose of the unused product and the device for its administration in accordance with the requirements for disposal of medical devices.

Shelf life:

The product shelf life is two years provided in the original packaging, provided the storage conditions are observed.

Storage and transport conditions:

Store at temperatures between $+5^{\circ}\text{C}$ and $+30^{\circ}\text{C}$. Do not expose to direct rays of light and high temperatures. Store in closed original packaging.

The expiry date is valid if the storage conditions in the original packaging are observed.

Implants are transported by all types of covered transport under the conditions of storage, in accordance with the rules of transportation of goods, applicable for each mode of transport.


Fig. 1 Operating procedures with syringe and needle after the package is opened

Ways and measures for disposal:

An unused medical product is classified as Class A waste used in Class B and must be disposed of by local public law organizations in accordance with SanPin 2.1.7.2790-10 "Sanitary and Epidemiological Requirements for the Management of Medical Waste".

A syringe and any unused materials should be disposed immediately after the procedure. Any residues of materials after the procedure cannot be used, stored or re-sterilized for use in several patients.

Needles are to be used once, and must be disposed immediately after use. Needles after application carry a biological hazard.

After use, needles should be disposed of in a purpose made puncture-proof and moisture-proof container with lid:

1. The puncture proof container should be collected and disposed by qualified professional sub-contractors.
2. Do not handle or touch needles once they have been disposed.
3. Puncture proof needle containers should not be reused.

It is prohibited:

- To destroy and cut medical devices manually once they have been used;
- To remove the needle from the syringe manually after it has been used, recapping the needle after the injection;
- To pour (overload) unpacked the needle from one container to another;
- Ram;
- To carry out any operation without gloves or necessary personal protective equipment and clothing;
- To use a soft disposable packaging for the collection of sharp medical instruments and other sharp objects;
- To establish disposable and reusable containers for waste collection at a distance of less than 1 m from heating appliances.

ATTENTION! ANY ADVERSE EVENTS SHOULD BE REPORTED TO THE AUTHORIZED REPRESENTATIVE ON THE TERRITORY OF THE RUSSIAN FEDERATION OR DIRECTLY TO THE MANUFACTURER.

WORLDWIDE DISTRIBUTED BY

INSTITUTE HYALUAL GMBH, RISCH, Birkenstrasse 47, 6343, Rotkreuz, Switzerland.

E-mail: info@hyalual.com, www.hyalual.com

MANUFACTURER: Diaco Biofarmaceutici S.r.l, Italy Via Flavia, 124, 34147 Trieste (TS), Italy

Manufacturing sites addresses:

- Diaco Biofarmaceutici S.r.l, Via Flavia, 124, 34147 Trieste (TS), Italy
- "Yuria-Pharm" Ltd, 18030, Cherkassy, str. Kobzarskaya, 108, Ukraine

Authorized Representative over Russian Federation territory:



Limited Liability Company "RedermaSwiss" (LLC "RedermaSwiss»), 105094, Russian Federation, Moscow, Zolotaya str., 11 build. 1, suits 8A17; 8A22 phone (499) 653-71-12, e-mail:


info@hyalual.ru

Registration certificate No.

Table 1. Explanation of designation symbols

Symbol	Title
Symbols on the secondary package	
	Prohibition of reuse Do not reuse - only for single use in one patient during one procedure
	Protect from moisture
	Fragile, handle with care

	Do not expose to sunlight
	Non-pyrogenic
	Sterile liquid inside the product. Steamed or sterilized by dry heat .
	Do not use if the package is damaged
	Do not re-sterilize
	Please read the instructions for use carefully
	Maximum permissible temperature ranges for storage
	Humidity limits
	Mark of circulation within the EU market and number of the authorized certification body
	Producer (Manufacturer)
	Production date (Manufacturing date)
	Shelf life, best before
	Batch code
Symbols on the sticker-grading of the syringe	
	Sterilized by steam or dry heat
	Prohibition of reuse Do not reuse - only for single use in one patient during one procedure

	Batch code
	Maximum permissible temperature ranges for storage
	Mark of circulation within the EU market and number of the authorized certification body
	Production date (Manufacturing date)
	Producer (Manufacturer)
	Shelf life, best before
The symbols on the back of the needle blister	
	Sterilization by ethylene oxide
	Mark of circulation within the EU market and number of the authorized certification body
	Caution! Please refer to the instructions
	Prohibition of reuse Do not reuse - only for single use in one patient during one procedure
	Producer (Manufacturer)
	Shelf life, best before
	Batch code
	Catalog reference
	Needle dimensions

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