Instructions for Use of medical device Nasal spray Quixx Barrier

Medical Device Name:

Nasal Spray Quixx Barrier 20 ml

Description and composition of the medical device:

Nasal spray QUIXX barrier is a sterile, clear, colourless, viscous solution without odour.

1 ml solution contains:

iota-carrageenan (Carragelose ®) (1,2 mg), kappa-carrageenan (0.4 mg), sodium chloride (5.0 mg) and purified water (up to 1 ml).

Does not contain preservatives.

Intended use

The nose is the gateway to our breathing passages and one of the protective factors against upper respiratory tract infections (URTI, flu). The nasal mucosa, which is covered by a mucin layer acts as a physiological barrier. It is lined with a special epithelial tissue, which provides the moisturizing of breathing air and its filtration from mechanical particles (like dust or pollen) and infectious particles (like for instance the cold virus).

Cold viruses can attach to the nasal mucosa where they multiply and spread. This results in URTI or flu diseases having the typical symptoms like a blocked or runny nose, a sore throat and cough etc.

Nasal Spray QUIXX barrier is used against the virus infections of upper respiratory tract at an early stage of disease. **Nasal Spray QUIXX barrier** helps to protect against colds and infections caused by respiratory viruses, while visiting public places, including presence in public transport and in contacts with people who may be infected.

Nasal Spray QUIXX barrier helps to

- shorten the duration of common colds
- relieve symptoms of common colds
- accelerate recovery
- reduce the risk of resurgence of virus-related symptoms of a cold by reducing viral load.

Nasal Spray QUIXX barrier also serves to moisten the nasal mucous membrane when it is dry or irritated. A moistened nasal mucosa is essential for the prophylaxis of infections caused by respiratory viruses and the support of the natural physiological functions of the nose.

Properties and mode of action:

Nasal Spray Quixx Barrier is a sterile aqueous saline solution free from preservatives. The main component is iota-carrageenan (Carragelose ®), a natural polysaccharide

extracted from marine red algae that keeps the nasal mucosa moist and forms a film around the cold virus. Thus the cold virus is immobilized. Hence, Carragelose ® prevents cold viruses from penetration into the body, thereby decreasing the viral burden on the epithelial tissue of the mucosa nasal cavities (Figure 1). The particles trapped by the protective film are further removed from the nasal cavity. **Nasal spray QUIXX barrier** due to its saline content helps to maintain the natural moisture of the nasal mucosa and thus supports the body's own defence against viruses and bacteria.

Fig.1

Keeping the
nasal mucosa
moist

Trapping and
immobilising/
blocking of
viruses

Reducing symptoms of
infections caused by
respiratory viruses and
removing the viruses

Protective properties of **Nasal spray QUIXX barrier**, which forms the barrier against the most common viruses of acute respiratory diseases, are proved by results of clinical studies.

Nasal spray QUIXX barrier contains none of drug preparations for medical usage, materials of animal and/or human origin within the composition and causes no pharmacological, immunological, or metabolic effects on the human body. Its action is based on the physical process, i.e., formation of a film (barrier) protecting the nasal mucosa from the penetration of viruses.

Nasal spray QUIXX barrier is a nasal spray of topical action and is not absorbed by the body. Spray may be used for a long time, if necessary.

Indications

- Prevention of infections caused by respiratory viruses while visiting public places, including presence in public transport and in contacts with people who may be infected.
- Moistening the nasal mucous membrane when it is dry or irritated.

Dosage and method of administration

How often and how much **Nasal spray QUIXX barrier** should be used?

At the first signs of a cold (e.g. a tingling or runny nose, etc.): spray 1 puff into each nostril at least 3 up to 6 times a day depending on your individual requirements.

Usage of **Nasal spray Quixx Barrier** should begin as early as possible after the first symptoms of a cold occur and be continued until the symptoms have abated. If the symptoms persist, you should consult a doctor. Nasal passages should be cleaned before use of **Nasal spray QUIXX barrier**.

To moisten a dry nasal mucous membrane, it is recommended to spray 1 puff into each nostril 3 times per day.

How should Nasal spray QUIXX barrier be used?

Remove the protective cap from the bottle. Before the first usage pump a few times until a uniform atomised spray is discharged. **Nasal spray QUIXX barrier** is now ready for use. Hold the bottle upright when spraying. In order to avoid injury, insert the spray head carefully and not too deep into the nose and press the spray nozzle. The solution should be sprayed alternately into each nasal passage. After use, cover the spray bottle with a protective cap.

Nasal spray Quixx Barrier is intended for use at home. **Nasal spray Quixx Barrier** is an individual use product.

Who can use Nasal spray QUIXX barrier?

Adults and children from 1 year old, and also women during pregnancy and in breast-feeding period.

Contraindications

Known hypersensitivity to any of the ingredients of the product.

Special precautions for use

Do not use **Nasal spray Quixx Barrier** if the package has any visible signs of damage. For hygienic reasons and to avoid the transmission of pathogenic agents, one spray bottle should always be used **by one and the same person**.

Possible side effects

No side effects of **Nasal spray QUIXX barrier** usage have been reported till now.

Interaction with other products

No interactions with other products have been reported till now.

Presentation form

20 ml glass bottle (142 puffs) with spray nozzle and protective cap.

One bottle with instruction for use in a carton box.

Basic technical and performance specifications of the medical device:

Nasal spray Quixx Barrier is manufactured in compliance with the global standards: Directive of EC 93/42/EEC, EN ISO 13485:2012 +AC:2012, EN ISO 14971.

The solution is sterile. «A» class of sterility.

Sterilization of the product is performed by the filtration sterilization method through the filter having porosity of 0.2 µm; sterilization of primary packaging: by dry heat or γ-radiation. No re-sterilization is needed. This medical device is sealed.

Designation and marking

	Carefully read the instruction before usage.
	Medical device was sterilized using the methods of aseptic processing.
	Temperature range, within which the medical device is safely stored.
	Medical device shall not be exposed to the sunlight.

Additional advertising information on the packaging is allowed.

Safety measures requirements

Safety of the product is ensured during the shelf life.

Environmental conditions requirements

Materials, which are used in manufacturing of this medical device, are not life-threatening, dangerous for life and environment, either during the maintenance of device or after its termination.

Storage conditions and terms of use

Keep out of the reach and sight of children! Store protected from light at temperature from 5°C to 30°C.

Do not use **Nasal spray Quixx Barrier** if the package has any visible signs of damage. For hygienic reasons and to avoid the transmission of pathogens, one spray bottle should always be used **by one and the same person**.

Warranty

The product is designed, manufactured, tested and packaged in compliance with all applicable requirements. The manufacturer guarantees the quality of the medical device until the expiration of its shelf life, provided the packaging is not damaged and all storage conditions of the product are observed.

Transportation conditions

Transportation of the product is possible by all types of the roofed transport in accordance with the rules of the cargo transportation for each type of transport.

Transportation of the medical device shall be performed at the temperature of 5 to 30°C. Do not freeze.

The procedure for disposal and destruction

Disposal of the product shall be performed in accordance with the requirements / rules of the buyer's country. On the territory of the Russian Federation: in terms of the hazard, the product belongs to class A (epidemiologically safe waste) and is subject to disposal together with household waste.

Shelf life

36 months.

The solution can be used for 12 months after first opening.

Do not use after the expiry date.

Sold without a prescription.

Manufacturer

MoNo chem-pharm Produkte GmbH, LeystraÙe 129, 1200 Vienna, Austria

Tel.: +43-1-330 06 71

Fax: +43-1-330 06 71 - 38

mail@mono.co.at

Produced:

1. MoNo chem-pharm Produkte GmbH, LeystraÙe 129, 1200 Vienna, Austria.
2. Laboratoires Pharmaster, Z.I. de Krafft, 67150 Erstein, France.

Organization authorized to accept claims:

OOO «Berlin-Chemie/A. Menarini», 123112, Moscow, Presnenskaya quay 10,
Business Center «Tower on Quay», Block B

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
Certification Register Number: 448/19

I herewith certify the authenticity and genuineness of the
foregoing signature which is acknowledged today before me
of **DI Dr. Bernhard WITTMANN**, born on 14.02.1973 (fourteenth
day of February nineteen hundred and seventy-three), A-1200
Vienna, Leystraße 129, acting as Managing Director with
individual signing power for, in the name and on behalf of
MoNo chem-pharm Produkte GmbH, with principal place of
business at A-1200 Wien, Leystraße 129. -----

According to today's inspection of the Companies' Register
kept by electronic data processing I herewith certify
pursuant to Section 89a (eighty-nine paragraph a) Notary
Act that the above mentioned Company is duly registered
with the Companies' Register of the Vienna Commercial Court
under Court File Number FN 189118 s and that **DI Dr.
Bernhard WITTMANN** is today duly authorized to sign
individually for, in the name and on behalf of the said
Company in his capacity as above mentioned. -----

Furthermore, I certify that the party has verified to have knowledge of the content of the document and to sign it free from compulsion. -----

Given at Vienna this 04.03.2019 (fourth day of March two thousand and nineteen) -----

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