
APPROVED BY
NEW LIFE RADIOLOGY S.R.L.
July, 14, 2017

DIGISENS

Intraoral Sensor

simplify your work

User Manual

NEW LIFE RADIOLOGY s.r.l.
Via Latina 17 - Telef. 011/781.994
10095 GRUGLIASCO (TO)
Partita IVA 05989400014

User Manual

DIGISENS Intraoral Sensor

MADE IN ITALY

Version: 5 January 2016 (Rev. 0)

TABLE OF CONTENTS

1	INTRODUCTION	4
2	GENERAL SYSTEM DESCRIPTION	6
3	SYSTEM USABILITY	8
4	SAFETY ASPECTS	9
	4.1.1 Electromagnetic environment	10
	4.1.2 Disposal hazards	14
5	CLEANING AND DISINFECTING	14
	5.2.1 First level disinfection	16
	5.2.2 Second level disinfection	16
6	IDENTIFICATION AND DESCRIPTION	16
7	TECHNICAL SPECIFICATIONS	17
8	PRECAUTIONS FOR USING THE SENSOR	20
9	SOFTWARE INSTALLATION	22
10	INSTALLATION AND PERIODIC CHECKS	31
	10.2.1 Frequency of monitor quality tests	33
	10.3.1 Frequency of spatial resolution test	34
	10.4.1 Frequency of low contrast resolution test	34
11	DOCUMENT STATUS	35

The original language of the DIGISENS user manual is Italian.

PRODUCTION

We thank you for trusting our company and choosing DIGISENS as your intraoral sensor. We appreciate your support and hope that DIGISENS serves you well. Our continued commitment lies in the complete customer satisfaction for each and every product we produce. Our DIGISENS intraoral sensor is designed to provide many years of reliable service during its use and handle. This user's manual will assist you on the installation and operation of your DIGISENS. Please read carefully the warnings and instructions, keep them for future reference.

Symbols used in this manual

Following icons are used in this manual:

SYMBOL	DESCRIPTION
	Indicates a "NOTE"; all text marked with this icon is very important and should be read carefully.
	Indicates a "WARNING"; all text marked with this icon relates to safety issues for the patient and/or operator.

Symbols used

Following symbols are used in this manual and on the DIGISENS packaging:

SYMBOL	DESCRIPTION
	Device with Type BF applied parts
	This symbol indicates that the DIGISENS sensor contains electrostatic-sensitive electronic parts susceptible to damage by electrostatic discharge. Refer to the section on Precautions for Use.
	The device contains solid materials which, at the end of its life cycle, must be disposed of at authorised recovery centres according to local regulations in order to prevent human health and environmental damages caused by improper disposal.
	NON-STERILE. DIGISENS is a non-sterile product and cannot be sterilised.

DIGISENS

	Temperature limits
	Product identification code
	Serial number
	Date of manufacture (year and month)
	Name and address of manufacturer
	Consult accompanying documents
	Conforms with EC Directive 93/42 and its amendments and supplements

1.3 General Introduction

The purpose of this manual is to provide the user with the instructions necessary for safe, secure, and efficient operation of this device.

SYMBOL	DESCRIPTION
	NOTE This manual is updated to correspond to the product status it is sold with, to ensure that the user has appropriate reference on operating the device and all safety-related aspects. This manual may not reflect any product variation that has no impact on operating procedures and safe use.

DIGISENS

	DIGISENS is an electro medical device and, as such, it may only be used by dentists, radiologists or other legally qualified professionals. Proper operation requires that it be used in combination with radiographic equipment and by persons with the necessary knowledge on radiation (x-rays) protection.
	The device must be used in compliance with the procedures included in this manual and never be used for purposes other than those provided for herein.
	The user is responsible for fulfilling all legal requirements related to the installation and use of the device.
	Read this manual carefully before using the device. Keep this manual near the device for future reference.

GENERAL SYSTEM DESCRIPTION

DIGISENS, manufactured by New Life Radiology S.r.l. is an easy-to-use digital intraoral system has the following unique features:

- Beveled edges
- Smooth casing
- Outstanding cable flexibility,
- Direct connection to personal computer for image acquisition via USB direct connector.

DIGISENS is based on CMOS technology, allowing the user to obtain a small pixel size (20 µm) thereby ensuring excellent spatial resolution; the use of a fiber optic plate (FOP) also helps to ensure a high quality image and, at the same time, increases the life of the sensor protecting it against incident x-rays.

The scintillator layer can be supplied in either Cesium Iodide (CsI) or Gadolinium Oxy-Sulfide (GOS) technology. The sensitive (active) area is available in two formats: standard size 1 and size 2 for intraoral film. All of these features make DIGISENS the perfect tool for any type of positioning, while ensuring maximum patient comfort. DIGISENS uses ARCHIMED SUITE image capturing and management software which saves images directly in DICOM format.

Intended use
The system is designed for taking intraoral x-rays of human teeth; the device can be used on patients, regardless of their type, gender or race.

DIGISENS NEW LIFE RADIOLOGY

The device is designed for both the dental and radiology market. It can be installed both in dedicated facilities (hospitals or clinics) and in residential structures equipped with appropriate shielding systems.

2.2 Classification of the device in accordance with EC 93/42
DIGISENS, in all of its configurations, is an active medical device, invasive through natural orifices, for temporary use, and intended for diagnostic purposes. This device falls into Class IIa according to the classification rules set out in Annex IX of EC 93/42, amended by Directive 2007/47/EC.

2.3 Applicable Standards
The standards applicable to the device mainly concern rules on general safety (for the patient and operator) and electromagnetic compatibility. The following standards apply:

Reference standard	Description
CEI EN 60601-1-2005	Medical electrical equipment Part 1: General requirements for basic safety and essential performance
EN 60601-1-2:2007	Medical electrical equipment Part 1: General requirements for basic safety and essential performance – Collateral standard: Electromagnetic compatibility – Requirements and tests

2.4 Type of Installation
DIGISENS is used as a temporary device and is not connected to the mains power source; power is supplied directly from the computer via the USB connection.

2.5 Address of manufacturer
The address of manufacturer is as follows:

	New Life Radiology Srl Via Latina 17 10095 Grugliasco (TO) - Italy website: www.newliferadiology.it
---	---

DIGISENS NEW LIFE RADIOLOGY

Optical principles of operation

DIGISENS sensor works like a regular digital sensor, i.e. it transforms the measured dose strikes each element of the sensor (pixel) in an electric signal that can be processed by an analog digital converter.

Conversion process includes the following steps:

Conversion of incident x-rays into visible light; this conversion takes place in the CsI or GOS sensitive layer.

Visible light is transferred, through the Optical Fibre, onto the sensitive layer of the CMOS.

The CMOS sensor converts the light rays into electric charges which are stored in special structures until reading.

In this way, each picture element (pixel) accumulates a number of charges proportional to the quantity of incident light beams and to the exposure time.

SYSTEM USABILITY

User profile and knowledge requirements

DIGISENS Medical System is designed to be used in radiography and dental facilities.

In most cases, the primary user is a professional who has the knowledge required to properly assess the risks and benefits associated with their radiological imaging technologies.

Users must have basic knowledge about:

Use of ionising radiation emissions

Harmful biological effects related to excessive use of ionising radiation

Methods to reduce the risk of excessive exposure to radiation as a patient (use of lead shields, etc.)

The operator should be familiar with using personal computers (PC) and the related programs, in order to be able to use the functions on the PC easily.

2. Training

After the system has been installed the operator will receive training on using the system and the image acquisition and viewing programs.

Training does not involve the use of special tools, only the digital sensor and related acquisition program.

3. Patient profile

The machine is suitable for use on any type of patient.

The different procedures for carrying out each exam based on the type of patient depend on the x-ray system used and are not included in this user manual.

DIGISENS

3.4 Usability

All documentation supplied with the DIGISENS system has been designed to help the operator in performing the operations.

The information contained in this manual call upon the knowledge requirements described for the user profile.

Information on using the acquisition, storage and processing system for images captured using the DIGISENS sensor is available in a specific manual, which should be read for further details.

Read this manual carefully before using the device. Keep this manual near the device for future reference.

4 SAFETY ASPECTS

WARNING

This chapter contains very important information concerning system, operator and patient safety. Read this chapter very carefully.

New Life Radiology S.r.l. designs and manufactures these devices in compliance with all relevant safety requirements. It also provides all necessary information for appropriate use and warnings on the risks associated with using X-rays for diagnostic purposes.

New Life Radiology S.r.l. shall not be held liable for:

1. Use of the DIGISENS device for any purpose other those for which it has been designed.
2. Damage to the device, injuries to the operator or patient caused by either incorrect installation or maintenance that does not follow the procedures contained in the User and Service Manuals provided with the device, as well as incorrect operating techniques.
3. Mechanical and/or electrical changes, made during or after installation, that differ from those listed in the Service Manual.
4. Only personnel authorised by New Life Radiology S.r.l. may install and perform technical work on the device.

DIGISENS

Warnings

	WARNING DIGISENS must be used in radiography, dental or hospital facilities.
	WARNING DIGISENS and its accessories are supplied non-sterile and cannot be heat sterilized. Follow the instructions below to ensure proper protection of the patient and operator.
	WARNING Do not sterilize the DIGISENS sensor in an autoclave or using dry heat as this could cause serious damage to the sensor. Do not sterilize with UV units.
	WARNING Do not immerse the USB connector of the sensor in cleaning fluids.
	WARNING For DIGISENS to work properly, it must be connected to a personal computer designed for image acquisition and image processing. The dedicated software must be installed on the personal computer.

4.1 Electromagnetic environment

	WARNING DIGISENS meets the requirements of IEC 60601-1-2 concerning electromagnetic emissions; therefore it may be used in electromagnetic environments that comply with the conditions described below. The customer or user of the system should assure that it is used in such an environment
---	---

DIGISENS

4.1.1.1 Electromagnetic emissions

	WARNING This device is intended for use by healthcare professionals only. The device may cause radio interference and may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orientation or relocating DIGISENS or shielding the location.
	WARNING The system is supplied with all cables required for its operation and connection to the computer. The use of other cables not supplied by the manufacture or additional to those supplied by the manufacturer may have a significant effect on the electromagnetic behaviour of DIGISENS, increasing its emissive capacity or decreasing its susceptibility.

Electromagnetic Emissions

The DIGISENS device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.

EMISSIONS TEST	COMPLIANCE	EMC ENVIRONMENT - GUIDANCE
RF Emissions CISPR 11	Group 1	The DIGISENS device uses RF energy only for its internal functions. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	DIGISENS is suitable for use in domestic establishments and in those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations /flicker emissions IEC 60000-3-3	Complies	

60601-1-2 Table

DIGISENS

Immunity

The DIGISENS device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601-1-2 TEST LEVEL	COMPLIANCE LEVEL	EMC ENVIRONMENT - GUIDANCE
Static electric (ESD) IEC 61000-4-2	± 6 kV contact ± 8 kV air	± 6 kV contact ± 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Radio frequency electromagnetic field (RF) IEC 61000-4-3	2 V/m for power supply lines 1 V/m for input/output lines > 3 mm	2 V/m for power supply lines 1 V/m for input/output lines > 3 mm	Mains power quality should be that of a typical commercial or hospital environment.
RF IEC 61000-4-3	± 1 V/m line to line ± 2 V/m line to earth	± 1 V/m line to line ± 2 V/m line to earth	Mains power quality should be that of a typical commercial or hospital environment.
Surge (dip, short duration) IEC 61000-4-11	<5% U _T (> 95% dip in U _T) for 0.5 cycles 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (> 95% dip in U _T) for 5 s	<5% U _T (> 95% dip in U _T) for 0.5 cycles 40% U _T (60% dip in U _T) for 5 cycles 70% U _T (30% dip in U _T) for 25 cycles <5% U _T (> 95% dip in U _T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the phototherapy device be powered from an uninterruptible power supply or battery.
Power frequency magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

60601-1-2 Table

Electromagnetic Immunity

The DIGISENS device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.

IMMUNITY TEST	IEC 60601-1-2 TEST LEVEL	COMPLIANCE LEVEL	EMC ENVIRONMENT - GUIDANCE RECOMMENDED SEPARATION DISTANCE
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 Vrms 150 kHz to 80 MHz	$d = 1.2 \cdot \sqrt{P}$
Radiated emissions IEC 61000-4-6	3 V/m 80 MHz to 2.5 GHz	3 V/m	from 80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P}$ from 800 MHz to 2.5 GHz

Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and "d" is the recommended separation distance in meters (m).
Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.
Interference may occur in the vicinity of equipment marked with the symbol shown on the right.

NOTE 1) At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2) These guidelines may not apply in all situations.
Electromagnetic protection is provided by absorption and reflection from structures, objects and people.

60601-1-2 Table

Recommended separation distances between portable and mobile RF communications equipment and DIGISENS devices

The DIGISENS device is intended to use in an electromagnetic environment in which radiated electromagnetic disturbances are controlled.

To prevent electromagnetic interference, maintain a minimum distance between the portable and mobile RF communications equipment (transmitters) and all parts of DIGISENS, including the power cords, as recommended below, according to the maximum output power of the transmitter.

Rated maximum output power of transmitter

Rated maximum output power of transmitter	150 kHz to 80 MHz $d = 1.2 \cdot \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \cdot \sqrt{P}$	800 MHz to 2.5 GHz $d = 1.2 \cdot \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

60601-1-2 Table

DIGISENS Intelligent Sensing

DIGISENS Intelligent Sensing

Disposal hazards

Parts of this device contain solid materials which, at the end of its life cycle, must be disposed of at authorised recovery centres according to local regulations.

Typically, the device contains the following materials and/or components:
Lead, electronic boards, electronic components.

Information on proper disposal of the system or parts thereof

Available in the European Union and in other European countries with separate waste collection systems
In accordance with Article 13 of Legislative Decree no. 151 of 25 July 2005 "Implementation of Directives 2002/95/EC, 2002/96/EC and 2003/108/EC, on reducing the use of hazardous substances in electrical and electronic equipment, and on the disposal of waste"

The crossed-out wheeled bin symbol on the device and/or its packaging indicates that it was put on the market after 13 August 2005, it complies with European Directives on the disposal of electrical and electronic equipment; at the end of its useful life, the product must be collected separately from other waste and must not be handled as domestic waste but disposed of separately and delivered to an appropriate collection point for recycling electrical and electronic equipment. Scrap it according to local waste disposal regulations.

Separate collection of this device when it reaches the end of its life is organised and managed by the manufacturer. The user should contact the manufacturer or its representative when it is time to dispose of the device and follow the system adopted by the manufacturer for the separate collection of end of life equipment.

5 CLEANING AND DISINFECTING

5.1 Disinfection products compatible with DIGISENS

Use 70% isopropyl alcohol to properly clean and disinfect the DIGISENS sensor.

Tests carried out have shown that the DIGISENS sensor can be immersed in the disinfection solution mentioned above without suffering any damage.

DIGISENS

	WARNING Do not use liquid disinfectant containing aldehydes (formaldehyde, glutaraldehyde, etc.).
	NOTE New Life Radiology Dental recommends that you only use disinfectants that are in compliance with EC Directive 93/42 on Medical Devices and that they bear the CE marking.

Scrupulously following the procedures provided below in order to ensure thorough cleaning and hygiene.

5.2 Disinfection procedures

	WARNING Cleaning operations must be performed with the device disconnected from the computer.
	WARNING DIGISENS and its accessories are supplied non-sterile and cannot be heat sterilized. To protect the health and safety of patients and prevent possible risks of contamination and/or cross-infection, please read and carefully follow the general guidelines provided below.
	Before each use, it is essential that you apply a protective device on the sensor. This protective device usually consists of a disposable sheath to be slipped on the sensor before use.
	Use a new protective sheath for each patient. For optimum performance, use protective sheaths specifically designed for the size of your sensor.
	Do not remove the disposable sheath by pulling on the sensor cable. Remove the protective sheath by carefully cutting it or using the tear-off strip (if provided on the sheath).
	The sensor, cable (sensor side only), and any accessories used must be carefully disinfected before each use.

DIGISENS

There are two levels of disinfection provided for the DIGISENS sensor. Disinfect the sensor using one of the two procedures described below depending on the conditions observed.

6.1 First level disinfection

The disinfection procedure must be carried out the first time the sensor is used and any time the protective sheath shows signs of damage.
Remove the protective sheath from the sensor and accessories and make sure that there are no residues of blood, saliva, tissue or secretions on it.
Carefully check the sensor and accessories used to make sure that there are no traces of organic matter on them.
Prepare the disinfecting solution according to the manufacturer's instructions.
Carefully disinfect the sensor, following the instruction provided by the manufacturer of the disinfectant.

6.2 Second level disinfection

This procedure is required when one or more of the following issues are noted during a visual inspection:

- the protective sheath is torn
- residues of organic matter (blood, saliva, other types of secretions, residual tissues) are found on the sensor and/or on the sensor-PC cable connection.

If any of these issues appear, follow the steps below for the second level disinfection procedure:

Wash the sensor thoroughly with soap and water to remove all organic matter. Be careful not to immerse the end part of the cable with the USB connection.

Prepare the disinfecting solution according to the manufacturer's instructions.

Carefully disinfect the sensor, following the instruction provided by the manufacturer of the disinfectant.

6 IDENTIFICATION AND DESCRIPTION

DIGISENS is a device designed for taking intraoral x-rays of the human tooth structure.

Normal use and operation of this device does not involve:

- the administration of biological substances
- the sterilisation of parts of the product, since only regular cleaning is required
- the interpretation of the final results
- the updating and modification of the control software.

6.1 Nameplates

Nameplates cannot be applied to the sensor due to its size and also due to hygienic reasons.

DIGISENS

The sensor is sold in a package that shows all the data necessary for its correct identification, as shown in the image below.

Additionally, a nameplate is applied near the USB connector and bears the serial number of the sensor.

6.2 Versions

DIGISENS is available in the following configurations:

- DIGISENS 1, which identifies a sensor size 1
- DIGISENS 2, which identifies a sensor size 2

Both of these versions can be equipped with either a CsI or GOS sensitive layer.

7 TECHNICAL SPECIFICATIONS

7.1 General characteristics

CHARACTERISTIC	VALUE
Type of device	DIGISENS
	New Life Radiology Srl 10095 Grugliasco (TO) - Italy
Class (according to IEC 60601-1 classification)	Class I with Type BF applied parts 
Degree of Protection	Equivalent to IP67
Sensitive (active) surface	30 x 20 mm (size 1) 33 x 25 mm (size 2)
Pixel size	20 x 20 µm
Number of Pixels	1500 x 1000 (size 1) 1650 x 1250 (size 2)
Spatial resolution	20 lp/mm typical (theoretical 25 lp/mm)
Dynamic range	58 dB
Scintillator type:	CsI or GOS
Photodiode detectors	Positioned outside of sensitive area
Maximum absorb. dose	50 Gy

7.2 Mechanical characteristics

CHARACTERISTIC	VALUE
Mechanical dimensions	25 x 39 x 12.5 mm (size 1) (L x W x H)
USB cable length	2 m

DIGISENS

10 Electrical characteristics

CHARACTERISTIC	VALUE
Power voltage	5 V DC
Power Supply	Directly via USB port
Maximum absorbed current	275 mA
Frame rate	0.7 fps

11 Electrical and optical characteristic

CHARACTERISTIC	VALUE
Local dark current (@23°C)	350 LSB/s
Integration level pixel (@70 kV)	340 µGy
Dynamic Range	57 dB
Grey levels	65535 (16 bit)
Sensitivity	15 LSM/µGy
RAW response non uniformity (XRNU)	± 30 %
Total dose irradiation	50 Gy

12 Environmental characteristics

CHARACTERISTIC	VALUE
Operating temperature	0 ± 35°C
Storage temperature	-20 ± +70°C
Sensor protection rating	Equivalent to IP67
Storage humidity	less than 75% non-condensing

13 Connection characteristics

DIGISENS sensor connects directly to the type A USB port of the computer.

The sensor is compatible with the standard USB 2 and the system is able to recognise up to three different sensors connected at the same time to the PC.

DIGISENS Intraoral Sensor

7.7 Computer specifications

The minimum system requirements for the computer used to acquire and store images are provided below.

COMPONENT	REQUIREMENTS
Operating System	Windows XP® SP3, Windows® 7 (32/64 bit) SP1, Windows® 8 64 bit Not compatible with previous version of O.S.
CPU	Intel i5-Z520M 2.5 GHz or higher
RAM Memory	2 GB or higher
Hard drive	10 GB minimum
USB Port	2.0
Video controller	1024x768 resolution in 65,000 colours (Ideally 1280x1024 - 16 million colours, 32 bit)

The hard drive size must be proportionate to the size of the archive that has to be used, with particular reference to the images that will be stored in it.

It is suggested to have a free space on the hard disk, after installing the program, considering the following data:

- Each radiographic image occupies disk space that varies between 500KB and 4MB.
- An image from the camera is set to 500 KB to 2 MB.

7.8 Monitor specifications

A high contrast and high definition monitor that has at least a 17" screen is preferable. Check the quality of the monitor using specific images (see paragraph 10.2)

Monitor with a good quality is essential for a correct diagnosis of the exam by analysing the image. A monitor that does not pass the required tests should not be considered a diagnostic-grade monitor. Using a non-diagnostic monitor will render the DIGISENS product non-diagnostic.

7.9 Hardware Installation

After installation, the system must not pose any risk to the patient or the operator.

Therefore, please follow the safety instructions provided below.

DIGISENS Intraoral Sensor

	The computer where our imaging systems are installed, along with all equipment connected to it, must bear the CE marking (IEC 950)
	The computer and all other associated equipment must be placed outside of the patient environment (about 1.5 m away from the chair).
	Only the sensor and x-ray generator are allowed to be placed in the so-called patient environment when using the sensor.
	Do not plug the PC used for the DIGISENS sensor in a power strip.
	No additional ground receptacle is required because DIGISENS is in compliance with safety standard EN 60601-1 for applied parts classified as type BF.

PRECAUTIONS FOR USING THE SENSOR

	For correct use of the sensor, consult the user manual.
	Use a disposable protective device (not supplied) to prevent patients from exposure to infections. Failure to use a disposable protective device may cause serious danger to the health of the patient.
	Do not sterilize the product using dry heat or autoclaving or UV devices.
	Before using the sensor, make sure it is in good condition (no cracks in the protective part of the sensor, cord ripped, etc.).
	Clean the sensor and the cable (at the output of the sensor) with a cloth moistened with 70% isopropyl alcohol for disinfection. Do not use other liquids or disinfectants and do not use too much rubbing.

DIGISENS

	The DIGISENS sensor has some electro-static sensitive parts: make sure to observe the precautions for use. Do not touch the sensor and computer screen at the same time. Do not touch the USB connector.
	When the sensor is not in use, store it away from static electricity.
	In case of any problems, do not use the product. In case of any fault during use, discontinue use of the sensor.
	Do not use USB connectors/ports if they are dusty or damp.
	Handle this product with care. If the sensor is not connected to PC, be careful not to damage the USB connector.
	Do not use a wet cloth or spray on the USB connector because it will deteriorate with moisture and can cause harm to the patient and / or operator.
	Do not attach or hang anything on the sensor, especially on the sensitive part.
	Do not exert pressure on the head of the sensor (tight sensor holders, clamps, etc.)
	Do not forcibly twist, bend, pull or pinch the cable. Do not pull the cable to remove the disposable protective sheath.
	When connecting/disconnecting the sensor, grasp the connector itself—never pull on the cable.
	The temperature of the sensor will rise considerably (even by 10°C) if it remains in operation for an extended period of time. Take care to use it only when the temperature is below 35° C.
	Although the sensor has been designed and engineered to be resistant to the entrance of liquids and powders, do not let the sensor immersed in liquid disinfectant, water or other chemicals for a long time.

DIGISENS

SOFTWARE INSTALLATION


DIGISENS sensor comes with ARCHIMED SUITE, an image acquisition and management software program specially designed for DIGISENS. ARCHIMED SUITE software is supplied on CD/DVD ROM.

Installation process for this system entails:

- Installation of main program
- Installation of the protection key drivers
- Installation of sensor drivers

Insert the CD/DVD in the CD drive of the computer where the program is to be installed; wait a few seconds until the start-up window shown below automatically appears.

In the main screen of the start-up procedure, it is possible to select the language that will be used during the installation. Please select your preferred language.

1.1 Installation of program

Select button "1" to install the ARCHIMED SUITE program on your computer; the installation procedure involves a number of different steps, all of which are noted via special windows and messages.

Press the "Next" button in each window to confirm the task to perform.

NOTE

During installation you will be asked if you would like to create a desktop shortcut (to add a program icon on your desktop). We recommend that you confirm the creation of this desktop shortcut for easy access to the ARCHIMED SUITE program.

DIGISENS

When the program is loading a progress bar will appear showing the percentage of completion of installation; when installation is completed, a window will appear confirming that the installation is complete; in this window you must select "End" to return to the main installation screen.

1.2 Installation of the protection key drivers

The ARCHIMED SUITE program is protected by a hardware protection key; without this key, the program will not work or remain in demo mode and it will not be possible to capture images.

Install the drivers for the protection key by pressing button "2"; the installation procedure is similar to that of the main program.

1.3 End of installation process

At the end of the installation process press the button "4" to exit the program and return to the main screen of the operating system in use.

1.4 Installation of sensor drivers

The sensor requires a special program (driver) for it to work properly; follow the following procedure:

1. Open the install disk with the command explorer
2. Navigate to the folder "Drivers_I-View_H"
3. Double click on the file "Driver Setup_32bit.exe" or "Driver Setup_64bit.exe" depending on the type of PC available.

NOTE

The selection of the sensor type is irrespective of the size of the sensor (size 1 or size 2) because it is automatically detected by the driver. Always select "I-View H" for proper operation of the acquisition program, regardless of the size of the sensor (or sensors) used.

Also here the driver installation process is similar to those above.

1.5 Using the ARCHIMED SUITE program

Instructions for using the ARCHIMED SUITE are available on the software CD/DVD. The instructions are provided in PDF format and can be accessed using Adobe Reader. If this software is not installed on your computer, you can download it for free at <http://get.adobe.com/it/reader/>

Read the documentation on using the ARCHIMED SUITE software carefully before you start to use it.

DIGISENS

Software configuration

Installation procedure has to be performed on the first use of the program, because ArchiMED Suite needs to be properly configured.

Insert USB dongle protection key on the PC and start the program

ArchiMED Suite is the main screen

Application is protected by user identification with different permission levels.

At the first run of the software it will appear this window

Click on **OK** and the program will automatically use the software starting drive as default database

The application will ask to you the preferred language.

Select your preferred language and click on **OK**

DIGISENS Medical Imaging

Then, if you don't have installed the demo database, this window will open


Click on **Yes** and ArchiMED Suite® will automatically create the environment for begin

NOTE


We recommend you to perform the functional verification of the monitor to verify that the monitor is proper for the visualization of the radiographic images. This procedure automatically started the first time that is accessed the images module.

At the first access at the program, you'll see this window

If you don't have registered the application yet, the follow window will be displayed.

DIGISENS Medical Imaging


NOTE

If you are using the software on the PC that is not the final one, please do not fill the above data.

You can enter now your cabinet data. If you don't want to write the data now, click on **Close**. The window will be displayed every time you open the application, until you insert the data. Please remember that you cannot change the registration data after you store them in the software protection key, so do not register any information before it is agreed with the customer.

If you are running for the first time ArchiMED Suite®, the password is **12345678**. Right after the program will ask you to change it, using this window.


Insert the old password one time, then the new for two times. Click on **OK** then insert the new password in the above window. Click on **Access** to enter in the program.

DIGISENS Individual Solution

Accessing ArchiMED Suite its main window will appear.


To configure properly the software, you must create a new dummy patient; click on **New** as showed on the next figure.

DIGISENS Universal Service

Following windows will be displayed


To create a new patient you must insert some important information. The needed information are highlighted with a darkest colour: name, last name, birth date. It is possible to complete these data in another moment too.

Before save a new patient you have to select the destination database (only if you have more than one active database). Click on the destination database and an arrow near the database icon will appear.

We suggest you to create a dummy patient, named Dummy Dummy on the following example.

Click on **Save** to save it and the following window will be displayed, where you have to click on **Images** to access the image acquisition module.


DIGISENS Intraoral Scanner

On the Image main module screen, click on 'Options'


On the 'Options' windows, click on Intraoral button, as showed.


DIGISENS Intraoral Scanner

Following image shows that no device has been selected.


Click on the down arrow "▼" of "Device Type" to select a list of available sensor types. Scroll down the list up to "I-View H" as displayed on the following figure, and save the selection by clicking to "OK" and return to main acquisition window.


DIGISENS Intraoral Sensor

Now, selecting the Intraoral 1 button, the acquisition screen will be displayed.

Once performed the emission, the sensor status windows will become green, as follow:

10 INSTALLATION AND PERIODIC CHECKS

After installing the DIGISENS intraoral sensor, perform the following steps to ensure that everything is correctly installed and that the sensor functions properly.

1. Installation check
2. Monitor quality test
3. Spatial resolution test

10.1 Installation check

This check is performed by launching the ARCHIMED SUITE software and verifying that it works properly and recognizes the DIGISENS sensor when it is connected to the computer.

Read the documentation on using the ARCHIMED SUITE program carefully before you start to use it.

NOTE

The selection of the language to be used by ARCHIMED SUITE is possible in the program's main screen and selecting 'Utilities'. In the new window, you can choose the desired language.


DIGISENS Intraoral Sensor

Monitor quality test


A good quality monitor is essential for a correct diagnosis of the exam by analysing the image. A monitor that does not pass the required checks should not be considered diagnostic. Using a non-diagnostic quality monitor will render the DIGISENS product non-diagnostic.

ARCHIMED SUITE program has a guided procedure to test the quality of the monitor

We highly recommend that you perform the functional test for the monitor to make sure it is suitable for viewing radiographic images. This procedure is automatically launched anytime you start the images module, as shown in the following window

Guided procedure. After confirmation the following screen will appear which contains descriptions of the tasks to perform.

DIGISENS Personal Monitor

Each task is associated to an appropriate image, accompanied by the related operating instructions.

Before continuing, the following general conditions must be met:

- monitors, especially CRT monitors, must be turned on for at least 30 minutes;
- no light sources are reflected on the screen;
- monitor surface is clean.

10.2.1 Frequency of monitor quality tests

We highly recommend that you perform the monitor quality test:

1. The first time using the program
2. Once every six months
3. Any time the monitor is replaced

10.3 Spatial resolution test

A spatial resolution phantom is required to perform this test. The captured images must not be adjusted in contrast or brightness.

1. Launch the ARCHIMED SUITE program to select images and create a new patient (e.g. "Periodic Check").
2. Place the phantom, with an additional filter made of 6 mm aluminium, touching the sensor, so as to cover all of the sensitive area. The sensor and phantom must be in contact with the collimator, so that they are in the centre of the beam.
3. On the x-ray generator, select the lowest exposure time possible.
4. If necessary, enable the sensor and then take the exposure
5. Check that the image obtained is properly exposed; if it is not, adjust the exposure time and repeat acquisition. In case of overexposure, increase the distance between the sensor/phantom and collimator
6. Repeat the correct exposure test; if necessary repeat steps 4 and 5 above.
7. Once a correct image is obtained you should cancel the previous images
8. Zoom in on the image and check that the line pairs are clearly visible. If necessary, adjust the contrast and brightness for optimal viewing.
9. Note the parameters used for exposure (date of execution, exposure time, focus-sensor distance, any additional filters and their characteristics, spatial resolution detected) of the image.

The test is considered passed if the spatial resolution is greater than 10 lp/mm1.

This test must be carried out on all x-ray generators that use the digital sensor, identifying the images for each sensor.

DIGISENS Personal Monitor

If two or more digital sensors are used, the test must be carried out on each digital sensor, identifying the images for each sensor.

Frequency of spatial resolution test

Spatial resolution test must be performed once every 6 months.

When following the procedure described in the paragraph above, the following items must be taken into consideration during the periodic tests:

The results must be the same as those described for the initial system check.

If the visual performance appears to be degraded:

1. Perform the monitor test, checking that it corresponds to the initial quality

2. If the results of the monitor test are positive, call the Technical Service Centre so they can check the performance of the sensor and/or x-ray generator.

Low contrast resolution test

The program used to test low contrast resolution also generates the image necessary to test low-resolution contrast.

Follow the same procedure described above, noting also the minimum value of the step when you are able to distinguish it from the background.

The test is passed if you are able to distinguish the step from 1.5 mm

Frequency of low contrast resolution test

The test must be performed with the same frequency as the spatial resolution test; the results must be the same as those of the acceptance test.

If the periodic low contrast resolution test is not passed, the same considerations for the spatial resolution test apply.

Periodic maintenance

The sensor DIGISENS does not require special maintenance except for the above controls.

It is recommended, however, perform the following checks on a monthly basis:

Check for presence and legibility of the S/N on the plate located near the USB connector;

DIGISENS digital sensor

- if it appears damaged or not readable, request a copy from the manufacturer.
2. Check the integrity of the USB cable and the connector.

11 DOCUMENT STATUS


Version	Date	Page N°	Parts / Modified Pages
0	05/01/2016	----	First Version

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