Processing instructions KaVo ARCTICA C-Temp-Blanks


Approved by

«Kaltenbach & Voigt GmbH»

Date : 16.09.2013

Stefan Wankel (Signature)
Name of Responsible Person

Post
Stamp

Always be on the safe side.

KaVo. Dental Excellence.

Sold by:
Dental GmbH
Bismarckring 39
88400 Biberach
+49 7351 56-0
+49 7351 56-1488

Manufacturer:
Kaltenbach & Voigt GmbH
Bismarckring 39
D-88400 Biberach
www.kavo.com

Table of contents

1 User instructions 2

1.1 User guide 2

1.1.1 Abbreviations 2

1.1.2 Symbols 2

1.1.3 Target group 2

2 Safety 3

2.1 Description of safety instructions 3

2.1.1 Warning symbol 3

2.1.2 Structure 3

2.1.3 Description of hazard levels 3

3 Information on materials 4

3.1 Technical specifications of KaVo ARCTICA C-Temp Blanks 4

3.2 Material numbers of KaVo ARCTICA C-Temp Blanks 4

3.3 Intended use 4

3.4 Indications for the use of KaVo ARCTICA C-Temp blanks 5

3.5 Contraindications 5

3.6 Customised products 5

3.7 Preparation basics 5

3.8 Transportation and storage 5

3.9 Disposal 6

4 Processing 7

4.1 C-Temp-specific design parameters 7

4.2 Manufacturing 7

4.2.1 Preparing C-Temp blanks for manufacturing 7

4.2.2 Setting the connectors 7

4.2.3 Setting pins 8

4.2.4 Separating C-Temp workpieces from the blank 8

4.3 Final tasks 8

4.3.1 Finishing the C-Temp piece 8

4.3.2 Veneering C-Temp frameworks 8

4.3.3 Shipping C-Temp frameworks 9

4.3.4 Incorporating the C-temp framework 9

5 Identification plates 10

1 User instructions

1.1 User guide

1.1.1 Abbreviations

IFU	Instructions for Use
CI	Care instructions
AI	Assembly instructions
TI	Technician's instructions
SC	Safety checks
IEC	International Electrotechnical Commission
RI	Repair instructions
EMC	Electromagnetic compatibility
VA	Processing instructions

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1.1.3 Target group

This document is for dental technicians, dentists, and trained professionals.

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3 Information on materials

Fibreglass-reinforced polymer from a materials science perspective

The new generation of dental framework materials are fibreglass-reinforced polymers. However, there are crucial differences in composition, strength, and aging. KaVo ARCTICA C-Temp blanks are industrially-manufactured polymer blanks with ideal and homogeneous structure. They consist of continuous molecular chains that increase the strength of the polymer to over 500 MPa. This strength is sufficient to support large spans and is suitable for long-term temporary restorations (for durations between 3 and 12 months).

The solubility and water absorption of the polymer is low that no difficulties arise while worn in the mouth.

The large amount of glass in this material can ensure a very strong adhesive bond between the framework material and conventional veneer plastics.

C-temp workpieces made using CAD/CAM have the following additional advantages:

- high comfort of wear owing to the low weight
- low reactivity with physiological solutions
- reproducible results
- no change in the material properties
- manufacturing with optimised properties (micro-texture)
- can be veneered with commercial composite veneers
- translucent to x-rays

3.1 Technical specifications of KaVo ARCTICA C-Temp Blanks

Flexural strength	500 MPa
chemical solubility	0.55 µg/cm³
Water absorption	9.6 µg/cm²
Bending modulus	20,000 MPa
Adhesive bond (GC/Gradia)	> 20 MPa
Adhesive bond (Ivoclar/Adoro)	> 18 MPa

3.2 Material numbers of KaVo ARCTICA C-Temp Blanks

2 x C-Temp Blank (70/40/20) Mat. no. 1.008.7465

3.3 Intended use

KaVo ARCTICA C-Temp Blanks are medical devices. Dental restorations (fixed dental prosthesis) are manufactured from these medical products in the KaVo ARCTICA Engine milling and grinding unit. The holders screwed to KaVo ARCTICA C-Temp blanks can be reused.

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- Do not prepare the product again or re-use.

- Provisional veneered crowns for the anterior and posterior dentition
- Provisional veneered bridges up to a span of 60 mm are possible.
- Long-term provisional of crowns and bridges wearable up to 12 months
- Provisional veneered bridges up to a span of two pontics.

3.5 Contraindications

- Duration of use > 12 months.
- Non-veneered, non-painted crowns and bridges.
- Bar cross section of less than 7 mm² in the anterior area.
- Bar cross section of less than 9 mm² in the posterior area.
- Temporary abutments.

3.6 Customised products

3.7 Preparation basics

All types of preparations can be used for a C-Temp restoration.
All transitions from the axial to the occlusal or incisal surfaces are to be rounded off with a radius of 1 mm or blocked out on the model using Scan wax.

See also: KaVo preparation basics (Mat. no. 1.004.9297)

3.8 Transportation and storage

The C-Temp material is nonhazardous as defined by national and international regulations.

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- Store KaVo ARCTICA C-Temp blanks in a dry and clean place.

3.9 Disposal

- ▶ KaVo ARCTICA C-Temp blanks can be disposed with household garbage in compliance with local regulations (waste code: 200301).
- ▶ If permitted by regulations, the packaging can be recycled (used paper, cardboard) or disposed with household garbage.

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- Thickness of edges: at least 0.2 mm
- Thickness of framework: at least 0.5 mm
- Cross-section of the bars in the anterior area: at least 7 mm²
- Cross-section of the bars in the lateral tooth area: at least 9 mm²
- Transition area: no more than 1.2 mm

4.2 Manufacturing

Please use only the following ZS dental burs for processing:

- Mat. no. 1.008.7095ZS Milling Pin D2 long shaft
- Mat. no. 1.008.7094ZS Milling Pin D1 long shaft
- Mat. no. 1.008.7093ZS Milling Pin D1
- Mat. no. 1.008.7092ZS Milling Pin D 0.5
- Mat. no. 1.008.7091ZS Milling Pin D2

4.2.1 Preparing C-Temp blanks for manufacturing

- ▶ Screw the C-Temp blank to the ARCTICA exchangeable holder using the ARCTICA universal tool.
- ▶ Tighten the internal hexagon screws at the proper torque.
- ▶ Insert the C-Temp blank into the tool holder of the KaVo ARCTICA Engine until it reaches the stop.
- ▶ Tighten the clamping bolt of the tool holder at the proper torque using the universal tool.
- ▶ Use the "KaVo CSS" nesting software.

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4.2.3 Setting pins

- ▶ Pins do not have to be set for C-Temp jobs.

4.2.4 Separating C-Temp workpieces from the blank

- ▶ Separate workpieces carefully from the blank using a suitable carbide cutter. Direct your special attention to the edges of the restorations.

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⚠ CAUTION

Irritant dust during processing of the polymers.
Irritation of the eyes and mucous membranes.

- ▶ Process polymers in the presence of local dust extraction only.
- ▶ Make sure the product does not contact the eyes or mucous membranes during processing.
- ▶ Do not inhale dust from grinding.
- ▶ Wear a dust mask and protective goggles.
- ▶ Do not eat or drink during work.

4.3.2 Veneering C-Temp frameworks

Note

The C-Temp frameworks must be veneered. The surfaces to be veneered do not have to be blasted.

- ▶ The remaining veneer procedure can be found in the processing instructions of the veneer ceramic manufacturer.

The outstanding bond of GC Gradia and IVOCLAR Adoro was demonstrated according to ISO standard 10 477.

4.3.3 Shipping C-Temp frameworks

Note
Only send frameworks that have been prepared in accordance with the currently relevant hygiene guidelines for shipping.

4.3.4 Incorporating the C-temp framework

	CAUTION
	The material is a medical product that contacts body surfaces (see DIN EN ISO 7405,4.1.2.). Infections ▶ The material should not be allowed to contact pulp and/or bone.

- ▶ Insert the framework using the customary fastening materials for temporary restorations.

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10 Kg max.

Sample

Type	ARCTICA C-Temp Blank (X) Material, (dimension in mm)
	Quantity
REF	Reference number or material number
	CE mark according to Medical Devices Directive EC 93/42
LOT	Batch ID
	Manufacturing date Year - month
	Comply with electronic processing instructions

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
Sefan Trampel (Signature)

Name of Responsible Person

Head of Quality Management
& Regulatory Affairs

Post

Stamp

Always be on the safe side.

KaVo. Dental Excellence.

Table of contents

User instructions	2
1.1 User guide	2
1.1.1 Target group	2
1.1.2 Abbreviations	2
1.1.3 Symbols	2
Safety	3
2.1 Description of safety instructions	3
2.1.1 Warning symbol	3
2.1.2 Structure	3
2.1.3 Description of hazard levels	3
Information on materials	4
Processing	7
4.1 Planning the model	7
4.2 Manufacturing	7
4.3 Final tasks	9
4.3.1 Reworking of frameworks produced by KaVo ARCTICA ZS	9
4.3.2 Veneering of frameworks produced by KaVo ARCTICA ZS	10
4.3.3 Sending frameworks produced by KaVo ARCTICA ZS	10
4.3.4 Integration of frameworks produced by KaVo ARCTICA ZS	10
Identification plates	12
Transportation and storage	14

1 User instructions

1.1 User guide

1.1.1 Target group

This document is for dental technicians and laboratory personnel.

1.1.2 Abbreviations

Abbreviation	Explanation
IfU	Instructions for Use
CI	Care instructions
AI	Assembly instructions
TI	Technician's instructions
SC	Safety checks
IEC	International Electrotechnical Commission
RI	Repair instructions
EMC	Electromagnetic compatibility
VA	Processing instructions

1.1.3 Symbols

	See the Safety/Warning Symbols section
	Important information for users and technicians
	CE mark (European Community). A product bearing this mark meets the requirements of the applicable EU directive.
	Action required

2 Safety

2.1 Description of safety instructions

2.1.1 Warning symbol

	Warning symbol
---	----------------

2.1.2 Structure

	 DANGER The introduction describes the type and source of the hazard. This section describes potential consequences of non-compliance. ► The optional step includes necessary measures for hazard prevention.
--	---

2.1.3 Description of hazard levels

Safety instructions with three hazard levels are used in this document to prevent personal and property damage.

	 CAUTION CAUTION indicates a hazardous situation that can cause damage to property or mild to moderate injuries.
--	---

	 WARNING WARNING indicates a hazardous situation that can lead to death or fatal injury.
--	---

	 DANGER DANGER indicates a maximal hazard due to a situation that can directly cause death or fatal injury.
--	--

3 Information on materials

Zirconium dioxide as defined by materials science

Zirconium dioxide is a high-performance ceramic and is an oxide ceramic.

Zirconium dioxide possesses an excellent biocompatibility and has therefore been used for implants for years in addition to dentistry. In addition, its elasticity modulus is high, and resists breaking. This quality is achieved by doping with yttrium oxide. When external energy is applied, for example in a crack, the volume of the material increases at this site at room temperature in a stabile phase.

The arising pressure counteracts the crack and prevents it from spreading and the ceramic from failing.

The transformation is the basis of the high strength, the impressive toughness and the excellent long-term stability.

CAD/CAM fabricated zirconium oxide dental work additionally offers the following benefits:

- Corrosion resistance
- Very fine particle size
- No glass phase as the particle binder
- Extremely high density
- Porosity-free
- Conventional incorporation
- Adhesive or non-adhesive

KaVo ARCTICA ZS-Blanks consist of pre-sintered, yttrium-stabilised zirconium oxide. They are easy to process in this soft state, and are fabricated larger on the KaVo ARCTICA Engine and then sintered to achieve full density in a special high-temperature oven. They shrink to form a very strong, precisely fitting framework.

Technical data on KaVo KaVo ARCTICA ZS-Blanks:

CTE (25°C to 500°C)	10.0 * 10 ⁶ K ⁻¹
Flexural strength	1200 MPa
Elasticity modulus	210 GPa
Fracture toughness	8 MN/m ^{3/2}
Chemical solubility	10 µg/cm ³
Thermal conductivity	2.5 W/ mK

Material numbers of KaVo ARCTICA ZS-Blanks

- ARCTICA ZS-Blank (B70/12) (Mat. no. 1.008.7461)
- ARCTICA ZS-Blank (B70/16) (Mat. no. 1.008.7462)
- ARCTICA ZS-Blank (B70/20) (Mat. no. 1.008.7463)
- ARCTICA ZS-Blank (B70/25) (Mat. no. 1.008.7464)

Intended use

KaVo ZS-Blanks are medical devices. Dental restorations (fixed dental prosthesis) are manufactured from these medical products in the KaVo ARCTICA Engine milling and grinding unit. The holders attached to the KaVo ZS-Blanks are disposed of after use.

Indications for the KaVo ARCTICA ZS-Blanks

⚠ CAUTION

The product is not re-usable.
Loss of biocompatibility, danger of infection.

- ▶ Do not prepare the product again or re-use.

- Crown frameworks in the anterior and posterior areas
- Bridge frameworks in the anterior and posterior areas
- Primary conical crowns and telescopic crowns
- Inlay frameworks
- Implant abutments based on titanium bases in the anterior and premolar regions

Contraindications for the KaVo ARCTICA ZS-Blanks

- insufficient oral hygiene
- bruxism
- tangential preparations
- unsatisfactory preparation results
- post and cores
- unveneered crowns, inlays, onlays and bridges
- when there are more than 2 pontics
- with insufficient bar cross-sections
 - 9 mm² with 1 pontic
 - 12 mm² with 2 pontics
- Implant abutments without a titanium base
- Implant abutments in the molar region
- Cantilevered bridges
- Attachments
- with provisional integration

It is the responsibility of the dentist to determine whether to use customized products such as long lower jaw bridges.

The crown and bridge frameworks made of ZS-Blanks can be individually veneered with the veneer ceramic developed especially for zirconium dioxide.

Preparation basics

For an appropriate zirconium oxide treatment, the following preparation foundations must be observed. The preparation can be made either with a 90° hollow neck, or a shoulder with rounded inner edge. Try to achieve a circular preparation depth of 1 mm.

All transitions from the axial to occlusal to incisal surfaces should be rounded with a radius of 1 mm.

Transportation and storage

The ZS material is nonhazardous as defined by national and international regulations.

Transportation/storage temperature	-20 °C to +55 °C
Humidity	5 % to 95 %

- ▶ Store ZS-Blanks in a clean and dry location.
- ▶ Store the blanks in the original packaging at room temperature (+20 °C).

Disposal

- ▶ ZS-blanks can be disposed with household waste while observing local regulations (waste code no.: 200301).
- ▶ If permitted by regulations, the packaging can be recycled (used paper, cardboard) or disposed of with household waste.

4 Processing

4.1 Planning the model

The following CAD parameters must be adhered to guarantee a safe design:

- Margin thickness: at least 0.3 mm
- Framework thickness: at least 0.55 mm
- The cross-section of the connecting bars is dependent on the span width of the restoration:
 - 3-member: at least $9 \text{ mm}^2 = 1 \text{ ZWG}$
 - 4-member: at least $12 \text{ mm}^2 = 2 \text{ ZWG}$
- Transition area: max. 0.5 mm

The shape of the connecting bars has a direct influence on the stability of the work. A round or oval cross-section provides better static reinforcement than a square one.

The height of the bar is more important than its width.

4.2 Manufacturing

Note


Check the concentration, purity and level of the cooling lubricant weekly, see instructions for use combination and base test (**Mat. no. 1.003.0825**).

For shaping, use the following ZS burrs:

- ZS Milling Pin D2 long shaft (**Mat. no. 1.008.7095**)
- ZS Milling Pin D2 (**Mat. no. 1.008.7091**)
- ZS Milling Pin D1 long shaft (**Mat. no. 1.008.7094**)
- ZS Milling Pin D1 (**Mat. no. 1.008.7093**)
- ZS Milling Pin D0.5 (**Mat. no. 1.008.7092**)


Nest restorations

- ▶ Position larger restorations at the margin of the blank in such a way that a socket with an even installation surface ① remains.

Restoration to be manufactured in the nesting software

- ▶ Set three connectors ②, ③ and ④ for the connection of the restoration with the socket ①.
- ▶ Set further connectors to stabilise the restoration for the production process.


Restoration to be manufactured in the nesting software

- ▶ Set the cutting depth of the slit cuts ① at 100 %.
- ▶ After the production, separate the restoration with the cutting disk from the blank at position ②.


Sinter the restorations

- ▶ Place the restoration on a socket in the oven.

- Sinter the restoration in the KaVo Everest Therm with the „Sinter ZSII” program.

Restoration after the sintering process

Restoration of a posterior bridge prepared for the sintering process

4.3 Final tasks

4.3.1 Reworking of frameworks produced by KaVo ARCTICA ZS

Note

The entire framework does not have to be reworked as is the case with metal frameworks.

- Observe the processing instructions for the respective veneer ceramic.

⚠ CAUTION

Hot preheating oven
Burning hazard!

- Be particularly cautious when working with the preheating oven.

- After sintering, the ceramic bowl can be removed from the oven at 300 °C.
- Remove the piece from the bowl only when it has cooled down and is safe to touch.
- Do not cool the restoration with compressed air. This can cause cracks.

- ▶ Avoid inhaling the zirconium sand.
- ▶ If necessary, the cavity can be blasted with aluminium oxide 50 µm.
- ▶ Observe the additional information given by the manufacturer of the fixing material.
- ▶ If you need to refine a sintered ZS piece, use water cooling and a turbine. Do not apply pressure, and use tools with new diamonds.

Note

Avoid the formation of sparks!

Note

Do not blast the areas to be veneered! Blasting is not required.

- ▶ Preheating before veneering as instructed by Vita is recommended to relieve stress in the material and degrease the surface.

t. [°C]	Hold [min]	Heat up [min]	Heat up [°C/min]	Temperature [°C]	Hold [min]	VAC [min]
1000	-	5.00	as low as possible	approx. 1000	15.00	-

(source: Vita firing table for Vita VM9).

4.3.2 Veneering of frameworks produced by KaVo ARCTICA ZS

Note

ZS frameworks must be veneered. The surfaces to be veneered do not have to be blasted.

- ▶ The additional steps necessary to prepare for veneering can be found in the processing instructions of your veneer plastic manufacturer.

4.3.3 Sending frameworks produced by KaVo ARCTICA ZS

- ▶ Only send frameworks that have been prepared in accordance with the currently relevant hygiene guidelines for dispatch.

4.3.4 Integration of frameworks produced by KaVo ARCTICA ZS

Note

Provisional cementing is highly inadvisable.

- ▶ Clean and disinfect the framework according to the relevant hygiene guidelines before inserting into the mouth of the patient.
- ▶ The framework is permanently inserted using the usual fastening materials for adhesive and non-adhesive fastening.

- ▶ Additional information can be obtained from the preparation guide (Mat. no. 1.004.9297).

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**2012 - 07**

Shrinkage
20%

see operating instructions
www.kav-o.com

55°C

Transport

Storage

-20°C

95%

5%

700 hPa

1060 hPa

10 Kg max.

2015 - 07

Caution: Federal law restricts this device to sale by or on the order of a health care professional/ dentist. For dental use only.

Made in Germany (European Community)

**ARCTICA ZS-Blank (B70/12)**
REF: 1.008.7461

55°C

Transport

Storage

-20°C

95%

5%

700 hPa

1060 hPa

10 Kg max.

2015 - 07

Sample

Material number	Abbreviation of the dimension
1.008.7461	B70/12
1.008.7462	B70/16
1.008.7463	B70/20
1.008.7464	B70/25


Type	KaVo ARCTICA ZS-Blank (X) Material, (dimension in mm)
	Quantity
REF	Reference number or material number
	CE mark according to EC Directive 93/42 for medical devices
LOT	Batch ID

12/14

 ООО ЦРМИ  INFO@CIRMI.RU  WWW.GOSZDRAVNADZOR.RU

	Manufacturing date Year - month
	Reference to maximum period of use (Same as "Can be used until...") Year - month
	Observe the electronic processing instructions

6 Transportation and storage

	Air pressure range 700 hPa to 1060 hPa
	Humidity range 5 % to 95 %
	Temperature limit -20 °C to +55 °C
	Permissible stacking load 10 kg
	Fragile - protect against impact!
	Protect from moisture!
	Protect from radiation and sunlight.

Manufacturer : "Kaltenbach&Voigt GmbH", Bismarckring, 39, D-88400, Biberach, Germany.

Representative : "KaVo Dental Russland GmbH", nab.Fontanki, 130 A, 190005 St.Peteraburg. Russia.

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Approved by
«Kaltenbach & Voigt GmbH»



Date : 16.09.2013

Stefan Trampler (Signature)
Name of Responsible Person

Head of Quality Management
Regulatory Affairs

Post
Stamp

Always be on the safe side.

KaVo. Dental Excellence.

Contents

1	User instructions	2
1.1	User guide	2
1.1.1	Abbreviations	2
1.1.2	Symbols	2
1.1.3	Target group	2
2	Safety	3
2.1	Description of safety instructions	3
2.1.1	Warning symbol	3
2.1.2	Structure	3
2.1.3	Description of hazard levels	3
3	Information on materials	4
4	Processing	7
4.1	Titanium-specific design parameters	7
4.2	Manufacturing	7
4.2.1	Preparing ARCTICA T-Blanks for manufacturing	7
4.2.2	Setting the connectors	7
4.2.3	Setting pins	8
4.2.4	Separating the titanium work from the blank	8
4.3	Final tasks	8
4.3.1	Refining the titanium restoration	8
4.3.2	Veneering titanium pieces	9
4.3.3	Shipping frameworks	9
4.3.4	Incorporating the titanium restoration	9
5	Rating plates	10

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2.1.2 Structure

	 DANGER The introduction describes the type and source of the hazard. This section describes potential consequences of non-compliance. ▶ The optional step includes necessary measures for hazard prevention.
--	---

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3 Information on materials

Titanium in the view of materials science

Used in medicine and dentistry for many years, titanium, a transition metal, is a pure metal with a degree of purity of 2. Pure titanium has proven to be valuable in medicine for decades as an implant material and is also time-proven in dental applications as an implant and structural material.

Titanium shows excellent corrosion properties which are based on the passivation layer that is produced immediately through a reaction with oxygen. The layer thus formed affords protection and electrical insulation. This means that other metals being incorporated into the oral environment along with titanium do not react with titanium.

The substance shows excellent biocompatibility. Bone and tissue cells show perfectly natural and unperturbed behaviour in the immediate vicinity of titanium producing, e.g., a normal calcified bone.

There are no undesirable biological reactions at crown edges, e.g. as gingival inflammation or discolouration.

Titanium parts made by CAD/CAM show the following additional advantages:

- high comfort of wear owing to the low weight
- X-ray translucence
- very low heat conductivity
- fairly inert to solutions and acids
- no alpha-case layer
- lasting high quality
- reproducible results
- no change in the material properties
- Manufactured with optimised properties (micro-texture)

KaVo ARCTICA T blanks are industrially-manufactured titanium blanks with an ideal and homogeneous structure.

Technical specifications of KaVo ARCTICA T-Blanks

CTE (25 °C to 500 °C)	9.2 * 10 ⁶ K ⁻¹
Chemical solubility	0.10 µg/cm ³
Modulus of elasticity	91 GPa
Flexural strength	322 MPa
Tensile strength	461 MPa
Vickers hardness	148

Material numbers of KaVo ARCTICA T-Blanks

2 x T-Blank (12) Mat. no. 1.008.7467

2 x T-Blank (16) Mat. no. 1.008.7468

Intended use

KaVo ARCTICA T-Blanks are medical devices. These medical devices are used in the milling and grinding unit, KaVo ARCTICA engine, to produce dental restorations (permanent and temporary restorations). The holders screwed to the KaVo ARCTICA T blanks can be reused.

Indications for the use of KaVo ARCTICA T blanks

	 CAUTION
	The product is not re-usable. Loss of biocompatibility, danger of infection. ▸ Do not prepare the product again or re-use.

- Crowns in the area of the anterior and lateral teeth
- Primary conical crowns and telescopic crowns
- Cantilevered bridges
- Maryland bridges
- Inlay bridges
- Inlays / onlays / full crowns
- Superstructures and tertiary frameworks
- Implant abutments

Note
It is the responsibility of the dental technician or dentist to determine whether to use customised products such as long low jaw bridges, cantilever bridges, attachments.

Frameworks made from KaVo ARCTICA T blanks can be veneered in custom-manner using veneer ceramics developed especially for titanium.

Contraindications

- More than 3 pontics per unit

Preparation basics

All types of preparations can be used for a titanium restoration. All transitions from the axial to the occlusal or incisal surfaces are to be rounded off with 1 mm or blocked out on the model using Scan wax.

See also: CAD/CAM Preparation basics (Mat. no. 1.004.9297)

Transportation and storage

T material is non-hazardous as defined by national and international regulations.

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- ▶ Store KaVo ARCTICA T blanks in a dry and clean place.

Disposal

- ▶ KaVo ARCTICA T blanks can be disposed with household garbage in compliance with local regulations (waste code: 200301).
- ▶ If permitted by regulations, the packaging can be recycled (used paper, cardboard) or disposed with household garbage.

4 Processing

4.1 Titanium-specific design parameters

In order to ensure the safety of the design, please comply with the following parameters:

- Thickness of edges: at least 0.2 mm
- Thickness of framework: at least 0.5 mm
- Cross-section of the bars in the anterior area: at least 5 mm²
- Cross-section of the bars in the lateral tooth area: at least 7 mm²
- Transition area: no more than 1.2 mm

4.2 Manufacturing

Note

Check the concentration, purity and level of the cooling lubricant weekly, see instructions for use combination and base test (Mat. no. 1.003.0825).

Please use only the following dental burs for processing:

- Titanium Milling Pin D0.5 (Mat. no. 1.008.7098)
- Titanium Milling Pin D1 (Mat. no. 1.008.7096)
- Titanium Milling Pin D1 long stem (Mat. no. 1.008.7099)
- Titanium Milling Pin D2 short (Mat. no. 1.008.7100)
- Titanium Milling Pin D2 (Mat. no. 1.008.7097)
- ARCTICA Exchangeable Holder (Mat. no. 1.008.7469)

Perform the following preparations before each use:

- ▶ Clean the chip drawer and drain area of the cooling lubricant in the machine.
- ▶ Check the cooling lubricant level.

4.2.1 Preparing ARCTICA T-Blanks for manufacturing

- ▶ Screw the T-Blank to the ARCTICA exchangeable holder for titanium blanks using the ARCTICA universal tool.
Tighten the internal hexagon screws at the proper torque.
- ▶ Push the ARCTICA exchangeable holder to the stop in the blank holder of the KaVo ARCTICA engine.
- ▶ Screw the ARCTICA exchangeable holder in place in the blank holder using the universal tool.
Tighten the clamping bolt at the proper torque.
- ▶ Nesting work with the nesting software KaVo CSS.

See also: Instructions for Use KaVo ARCTICA Engine (Mat. no. 1.009.4851)

4.2.2 Setting the connectors

The number of connectors and the maximal depth of the split cuts assuring the safety of the fabrication is always a function of the size and geometry of the restoration. For this reason, the minimal and maximal values specified in the following paragraph are for purposes of orientation only.

Coping:

- At least 2 connectors, depth of split cut ≤ 50 %.

Bridges:

- Anterior teeth: At least 1 connector on every other unit, depth of split cut ≤ 100 %
- Premolars/molars: At least 1 connector per unit, depth of split cut ≤ 100 %
- Terminal units: At least 2 connectors, depth of split cut < 100 %.
- Bridges from 6 units: additional connectors with depth of split cut < 100 %

4.2.3 Setting pins

- ▶ Do not insert pins in titanium work.

4.2.4 Separating the titanium work from the blank

CAUTION

Irritating dust when cutting out the titanium piece
Irritation of the eyes and mucosa!

- ▶ Do not inhale dust from grinding.
- ▶ A dust mask and protective glasses are required.

- ▶ Carefully separate the work from the blank (e.g. using a suitable hard-metal carbide cutter). Direct your special attention to the edges of the restorations.

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- ▶ Observe the processing instructions for the respective veneer ceramic.

CAUTION

Irritating dust from working with titanium
Irritation of the eyes and mucosa!

- ▶ Titanium must be worked on using a local dust extraction system.
- ▶ Make sure the milling dust does not contact the eyes or mucous membranes.
- ▶ Grinding dust must not be inhaled!
- ▶ A dust mask and protective glasses are required!
- ▶ Do not eat or drink while working!

Note

Local overheating during processing leads to a darkly discoloured, slightly rough oxide layer.
No sufficient bonding can be attained at this surface later on.

▶ Exclude local overheating.

- Additional requirements for trimming:
- Comply with the maximal speed of the handpiece you use.
 - Work at low contact pressure.
 - Work with clean milling handpieces only.
 - The production of sparks must be prevented.
 - Rework milled pieces with a special cross-cut titanium bur.

4.3.2 Veneering titanium pieces

- ▶ Consult the documentation of the respective veneer manufacturer in regard to preparing a framework for veneering.

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4.3.4 Incorporating the titanium restoration

	⚠ CAUTION The material is a medical device contacting body surfaces (see DIN EN ISO 7405, 4.1.2). Infection hazard. ▶ Do not contact the material directly to pulp tissue and/or bone.
--	---

- ▶ Make sure to clean and disinfect the framework sufficiently according to the applicable hygiene regulations before inserting the framework in the patient's mouth.
- ▶ Insert the framework using the customary fastening materials for titanium work.

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**2012 - 07**

**CE**

see operating instructions
www.kavo.com

-20°C

55°C

5%

95%

700 hPa

1060 hPa

Transport

Storage

10 Kg max.

Caution: Federal law restricts this device to sale by or on the order of a health care professional/ dentist.
For dental use only.

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**2012 - 07**

**CE**

see operating instructions
www.kavo.com

-20°C

55°C

5%

95%

700 hPa

1060 hPa

Transport

Storage

10 Kg max.

Sample

Type	ARCTICA T-Blank (X) Material, (dimension in mm)
	Quantity
REF	Reference number or material number
	CE mark according to Medical Devices Directive EC 93/42
LOT	Batch ID
	Manufacturing date Year - month
	Comply with electronic processing instructions

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 Kaltenbach & Voigt GmbH Bismarckring 39 D-88400 Biberach		
 ARCTICA Titan-Blank (B70/12) REF: 1.008.7467		
LOT	123456789	
 2012 - 07	 0124	 see operating instructions www.kavo.com
Transport -20°C — 55°C 5% — 95% 700 hPa Storage 1060 hPa 10 Kg max. 10 Kg max.   		
Caution: Federal law restricts this device to sale by or on the order of a health care professional/ dentist. For dental use only.		
Made in Germany (European Community)		
 ARCTICA Titan-Blank (B70/12) REF: 1.008.7467		
Transport -20°C — 55°C 5% — 95% 700 hPa Storage 1060 hPa 10 Kg max. 10 Kg max.   		

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