NC0J9708HB

NOTARIAL CERTIFICATE

TO ALL TO WHOM these presents shall come

I, Lee Meng Mew, NOTARY PUBLIC duly admitted, authorised to practise in the Republic of Singapore, DO HEREBY CERTIFY

that I was present on the 27th day of November 2019 and I did see the person, **PAULINE LIOW**, sign and execute the **attached Operational documentation for MEDICAL DEVICE as a certified true copy** and that the signature thereto subscribed is the true signature of **PAULINE LIOW** who signed the document in her capacity as Senior Regulatory Affairs Specialist, Merit Medical Singapore Pte. Ltd.

IN FAITH AND TESTIMONY whereof I the said notary have subscribed my name and set and affixed my seal of office at Singapore, this 27th day of November 2019.

NOTARY PUBLIC
SINGAPORE

By virtue of Rule 8(3)(c) of the Notaries Public Rules, a Notarial Certificate must be authenticated by the Singapore Academy of Law in order to be valid.

This **Authentication Certificate** only certifies the authenticity of the signature and the capacity of the person who signed the Notarial Certificate.

This **Authentication Certificate** is not valid if the seal of the Singapore Academy of Law is removed or altered in any way whatsoever. This Certificate does not authenticate or confirm the content of the Document attached to the annexed Notarial Certificate.

Verification code:
25215367

Authentication

- | | |
|---|---------------|
| 1. Country: | Singapore |
| 2. This public document has been signed by: | Lee Meng Mew |
| 3. Acting in the capacity of: | Notary Public |
| 4. Bears the seal/stamp: | Notary Public |

Certified


- | | |
|-----------------------------|-----------------------------------|
| 5. Authentication Cert No.: | AC0J970EEG |
| 6. At: | Singapore Academy of Law |
| 7. The: | 27th November 2019 |
| 8. By: | Melissa Goh, Deputy Director, SAL |
| 9. SAL Certification Seal: | 10. Signature: |

Melissa

CERTIFIED TRUE COPY

Pauline Liow
Senior Regulatory Affairs Specialist
Merit Medical Singapore Pte., Ltd.
27 November 2019

**Operational documentation for
MEDICAL DEVICE**

«Sets of vascular catheters Careflow™, Hydrocath Assure™»

Manufactured by

Merit Medical Singapore Pte., Ltd., Singapore

Manufacturer:

"Merit Medical Singapore Pte., Ltd. ", Singapore
10 Collyer Quay, # 10-01 Ocean Financial Centre, Singapore, 049315

The place of production of the medical device:

"Merit Medical Singapore Pte., Ltd. ", Singapore
Merit Medical Singapore Pte Ltd., 198 Yishun Avenue 7, Singapore 768926

The authorized representative of the manufacturer on the territory of the Russian Federation:

"Merit technologies"

Of the Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

phone: +7 (495) 221-89-02

1. DEVICE DESCRIPTION

1.1 The name and structure of the medical device

Sets of vascular catheters Careflow™, Hydrocath Assure™:

1. The catheter is a central venous Careflow™ in the set:

1.1 Careflow™ catheter with one lumen, 2.5 Fr (0.90 mm) in diameter, 60 mm, 100 mm, 150 mm, 200 mm, 300 mm long; diameter 3.0 Fr (1.10 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm; diameter 4.0 Fr (1.45 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm; diameter 5.0 Fr (1.75 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm; diameter 7.0 Fr (2.40 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm - 1 pc.;

1.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0,6 mm diameter, 450 mm long; diameter 0.9 mm, length 450 mm, 500 mm, 600 mm, 700 mm - 1 pc.;

1.3 Device for fixing a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);

1.4 Introducer needle, diameter 1.4 mm, length 70 mm - 1 pc.;

1.5 Coaxial introducer, diameter 1.70 mm, length 70 mm; diameter 1.27 mm, length 60 mm; diameter 1.0 mm, length 45 mm; with a diameter of 0.75 mm, a length of 38 mm - 1 pc. (if necessary);

1.6 Introducer FloSwitch™, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);

1.7 Dilator, 4Fr diameter (1.33 mm), 88 mm long; 5Fr diameter (1.7 mm), 88 mm long; 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);

1.8 Cap - 1 pc. (if necessary);

1.9 Floswitch™ rotating connector - 1 pc. (if necessary);

1.10 Clamp for catheter, 18 mm long, 25 mm - not more than 2 pcs. (if necessary);

1.11 Syringe, volume 5 ml - 1 pc. (if necessary);

1.12 Guiding Syringe, volume 5 ml - 1 pc. (if necessary);

1.13 Scalpel - 1 pc. (if necessary);

1.14 Surgical sheet with a dedicated surgical field for the patient - 1 pc. (if necessary);

1.15 Instructions for use - 1 pc

2. The catheter is a central venous Careflow™ in the set:

2.1 A Careflow™ catheter with a two-lumen diameter of 4.0 Fr (1.45 mm), 60 mm long, 100 mm, 150 mm, 200 mm, 300 mm; Diameter 5.0 Fr (1.75 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm; Diameter 7.0 Fr (2.40 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm - 1 pc.;

2.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0,6 mm diameter, 450 mm long; diameter 0.9 mm, length 450 mm, 500 mm, 600 mm, 700 mm - 1 pc.;

2.3 Device for fixing a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);

2.4 Introducer needle, diameter 1.4 mm, length 70 mm - 1 pc.;

2.5 Coaxial introducer, diameter 1.70 mm, length 70 mm; diameter 1.27 mm, length 60 mm; diameter 1.0 mm, length 45 mm; with a diameter of 0.75 mm, a length of 38 mm - 1 pc. (if necessary);

2.6 Introducer FloSwitch™, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);

2.7 Dilator, 4Fr diameter (1.33 mm), 88 mm long; 5Fr diameter (1.7 mm), 88 mm long; 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);

2.8 Cap - not more than 2 pcs. (if necessary);

2.9 Floswitch™ rotating connector - not more than 2 pcs. (if necessary);

2.10 Clamp for catheter, 18 mm long, 25 mm - not more than 2 pcs. (if necessary);

2.11 Syringe, volume 5 ml - 1 pc. (if necessary);

2.12 Guiding Syringe, volume 5 ml - 1 pc. (if necessary);

2.13 Scalpel - 1 pc. (if necessary);

2.14 Surgical sheet with a dedicated surgical field for the patient - 1 pc. (if necessary);

2.15 Instructions for use - 1 pc

3. The catheter is a central venous Careflow™ in the set:

3.1 Careflow™ catheter with three lumen, diameter, 5 Fr (1.75 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm, diameter 7 Fr (2.40 mm), length 60 mm, 100 mm, 150 mm, 200 mm, 300 mm - 1 pc.;

3.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0,6 mm diameter, 450 mm. 650 mm long; diameter 0.9 mm, length 450 mm, 500 mm, 600 mm, 700 mm - 1 pc.;

3.3 Device for fixing a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);

3.4 Introducer needle, diameter 1.4 mm, length 70 mm - 1 pc.;

3.5 Coaxial introducer, diameter 1.70 mm, length 70 mm; diameter 1.27 mm, length 60 mm; diameter 1.0 mm, length 45 mm; with a diameter of 0.75 mm, a length of 38 mm - 1 pc. (if necessary);

3.6 Introducer FloSwitch™, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);

3.7 Dilator, 4Fr diameter (1.33 mm), 88 mm long; 5Fr diameter (1.7 mm), 88 mm long; 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);

3.8 Cap - not more than 3 pcs (if necessary);

3.9 Floswitch™ rotating connector - not more than 3 pcs. (if necessary);

3.10 Clamp for catheter, 18 mm long, 25 mm - not more than 3 pcs. (if necessary);

3.11 Syringe, volume 5 ml - 1 pc. (if necessary);

3.12 Guiding Syringe, volume 5 ml - 1 pc. (if necessary);

3.13 Scalpel - 1 pc. (if necessary);

3.14 Surgical sheet with a dedicated surgical field for the patient - 1 pc. (if necessary);

3.15 Instructions for use - 1 pc

4. The catheter is a central venous Careflow™ in the set:

4.1 Careflow™ catheter with four lumen, diameter 8.5 Fr (2.80 mm), length 150 mm, 200 mm - 1 pc.;

4.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0.9 mm, length 450 mm, 500 mm, 700 mm - 1 pc.;

4.3 Device for fixing a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);

4.4 Introducer needle, diameter 1.4 mm, length 70 mm - 1 pc.;

- 4.5 Coaxial introducer, diameter 1.70 mm, length 70 mm; diameter 1.27 mm, length 60 mm; diameter 1.0 mm, length 45 mm; with a diameter of 0.75 mm, a length of 38 mm - 1 pc. (if necessary);
 - 4.6 Introducer FloSwitchTM, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);
 - 4.7 Dilator, 4Fr diameter (1.33 mm), 88 mm long; 5Fr diameter (1.7 mm), 88 mm long; 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);
 - 4.8 Cap - not more than 4 pcs (if necessary);
 - 4.9 FloswitchTM rotating connector - not more than 4 pcs. (if necessary);
 - 4.10 Clamp for catheter, 18 mm long, 25 mm - not more than 4 pcs. (if necessary);
 - 4.11 Syringe, volume 5 ml - 1 pc. (if necessary);
 - 4.12 Guiding Syringe, volume 5 ml - 1 pc. (if necessary);
 - 4.13 Scalpel - 1 pc. (if necessary);
 - 4.14 Surgical sheet with a dedicated surgical field for the patient - 1 pc. (if necessary);
 - 4.15 Instructions for use - 1 pc.
5. The catheter is a central venous CareflowTM in the set:
- 5.1 CareflowTM catheter with five lumen, diameter 9.5 Fr (3.22 mm), length 150 mm, 200 mm - 1 pc.;
 - 5.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0.9 mm, length 450 mm, 500 mm, 600 mm, 700 mm - 1 pc.;
 - 5.3 Device for fixing a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);
 - 5.4 Introducer needle, diameter 1.4 mm, length 70 mm - 1 pc.;
 - 5.5 Coaxial introducer, diameter 1.70 mm, length 70 mm; diameter 1.27 mm, length 60 mm; diameter 1.0 mm, length 45 mm; with a diameter of 0.75 mm, a length of 38 mm - 1 pc. (if necessary);
 - 5.6 Introducer FloSwitchTM, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);
 - 5.7 Dilator, 4Fr diameter (1.33 mm), 88 mm long; 5Fr diameter (1.7 mm), 88 mm long; 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);
 - 5.8 Cap - not more than 5 pcs (if necessary);
 - 5.9 FloswitchTM rotating connector - not more than 5 pcs. (if necessary);
 - 5.10 Clamp for catheter, 18 mm long, 25 mm - not more than 5 pcs. (if necessary);
 - 5.11 Syringe, volume 5 ml - 1 pc. (if necessary);
 - 5.12 Guiding Syringe, volume 5 ml - 1 pc. (if necessary);
 - 5.13 Scalpel - 1 pc. (if necessary);
 - 5.14 Surgical sheet with a dedicated surgical field for the patient - 1 pc. (if necessary);
 - 5.15 Instructions for use - 1 pc
6. The catheter is a central venous Hydrocath AssureTM in the set:
- 6.1 Hydrocath AssureTM catheter with one lumen, diameter 5 Fr (1.75 mm), length 200 mm - 1 pc.;
 - 6.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0.9 mm, length 450 mm, 500 mm, 600 mm, 700 mm - 1 pc.;
 - 6.3 Device for fixation of a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);

- 6.4 Introducer needle, diameter 1.4 mm, length 70 mm; - 1 pc.;
- 6.5 Coaxial introducer, diameter 1.70 mm, length 70 mm - 1 pc. (if necessary);
- 6.6 Introducer FloSwitch™, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);
- 6.7 Dilator, 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);
- 6.8 Cap - 1 pc. (if necessary);
- 6.9 Floswitch™ rotating connector - 1 pc. (if necessary);
- 6.10 Clamp for catheter, 18 mm long, 25 mm - not more than 2 pcs. (if necessary);
- 6.11 Syringe, volume 5 ml - 1 pc. (if necessary);
- 6.12 Surgical sheet with a dedicated surgical field for the patient - 1 pc. (if necessary);
- 6.13 Instructions for use - 1 pc

7. The catheter is a central venous Hydrocath Assure™ in the set:

- 7.1 Hydrocath Assure™ catheter with two lumen, diameter 7 Fr (2.40 mm), length 150 mm, 200 mm, 300 mm - 1 pc.
- 7.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0.9 mm, length 450 mm, 500 mm, 600 mm, 700 mm - 1 pc.;
- 7.3 Device for fixation of a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);
- 7.4 Introducer needle, diameter 1.4 mm, length 70 mm; - 1 pc.;
- 7.5 Coaxial introducer, diameter 1.70 mm, length 70 mm - 1 pc. (if necessary);
- 7.6 Introducer FloSwitch™, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);
- 7.7 Dilator, 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);
- 7.8 Cap - not more than 2 pcs (if necessary);
- 7.9 Floswitch™ rotating connector - not more than 2 pcs (if necessary);
- 7.10 Clamp for catheter, 18 mm long, 25 mm - not more than 2 pcs. (if necessary);
- 7.11 Syringe, volume 5 ml - 1 pc. (if necessary);
- 7.12 Instructions for use - 1 pc

8. The catheter is a central venous Hydrocath Assure™ in the set:

- 8.1 Hydrocath catheter Assure™ with three lumen, diameter 7 Fr (2.40 mm), length 150 mm, 200 mm, 300 mm - 1 pc.;
- 8.2 Guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger, diameter 0.9 mm, length 450 mm, 500 mm, 600 mm, 700 mm - 1 pc.;
- 8.3 Device for fixation of a catheter, 25 mm wide, 30 mm - 1 pc. (if necessary);
- 8.4 Introducer needle, diameter 1.4 mm, length 70 mm; - 1 pc.;
- 8.5 Coaxial introducer, diameter 1.70 mm, length 70 mm - 1 pc. (if necessary);
- 8.6 Introducer FloSwitch™, 2.0 mm diameter, 70 mm long; diameter 1.7 mm, length 70 mm; diameter 1.20 mm, length 60 mm; diameter 1.0 mm, long 45 mm - 1 pc. (if necessary);
- 8.7 Dilator, 6Fr diameter (2.0 mm), 88 mm long; diameter 8Fr (2.69 mm), length 110 mm; diameter 9Fr (3.0 mm), length 110 mm; diameter 10Fr (3.33 mm), length 110 mm - 1 pc. (if necessary);
- 8.8 Cap - not more than 3 pcs. (if necessary);
- 8.9 Floswitch™ rotating connector - not more than 3 pcs (if necessary);
- 8.10 Clamp for catheter, 18 mm long, 25 mm - not more than 3 pcs. (if necessary);

8.11 Syringe, volume 5 ml - 1 pc. (if necessary);

8.12 Instructions for use - 1 pc

1.2 Purpose of the medical device

Sets of vascular catheters Careflow™, Hydrocath Assure™ consists of separate sets of catheters: a catheter of a central venous Careflow™, a catheter of a central venous Hydrocath Assure™ is an intravenous catheter that is intended for short term use (no longer than 30 days) to access the human circulatory system via insertion through either the internal Jugular or Sub-Clavian Vein using Seldinger technique, where the catheter tip resides in the superior vena cava.

The conditions for using the medical device are in medical institutions.

1.3 Description of the principles of medical devices

Sets of vascular catheters Careflow™, Hydrocath Assure™ for short term use (30 days) in order to access to the circulatory system of man by Seldinger method.

Products used in hospitals.

Sets of vascular catheters Careflow™, Hydrocath Assure™ are intravascular catheters designed in single-lumen or multi-lumens with a radiopaque polyurethane catheter tubing available in gauge sizes from 2.5 to 9.5 Fr and in lengths from 60 mm to 300mm. The catheter tubing for the Careflow™ Central Venous Catheter is uncoated while the catheter tubing for Hydrocath Assure™ Central Venous Catheter has a hydrophilic antimicrobial coating of Benzalkonium Chloride (BAC).

The Hydrocath Assure™ CVC with the hydrophilic coating becomes lubricious when wet. This eases the insertion of the catheter into the patients and provides a thromboresistant surface. The Benzalkonium Chloride (BAC) is impregnated into the Hydromer layer, which leaches out from the Hydromer layer during use to prevent catheter related sepsis. The catheter has a tapered atraumatic tip which reduces the risk of venous damage and improves insertion. The catheter hub also includes a luer lock connection for accessory attachment. The over-needle introducer is also siliconised to ensure a smooth insertion.

The catheter is supplied sterile in a catheterization kit containing below listed components to facilitate catheter placement. These components are packed into the catheterization tray in various configurations:

- | | |
|--------------------------------|----------------------------------|
| • Introducer needle | • Shortstop Floswitch™ plug |
| • introducer catheter | • Secondary Fixation Device |
| • Guide wire with J-shaped tip | • Fenestrated drape |
| • Straight guide wire | • Scalpel |
| • vessel dilator | • Rotatable connector Floswitch™ |
| • Guiding Syringe | • Introducer Floswitch™ |
| • Syringe | |

The needle introducer supplied as part of kits is necessary for percutaneous puncture of the vein or arteries in order to facilitate the insertion of the catheter. The puncture of a vein is made by means of a needle with the attached syringe. When the needle is properly positioned in a vein or artery, blood flows freely into the syringe.

Coaxial introducer is intended for access to the vascular system.

The guide wire with the J-shaped tip is inserted through the needle for further insertion of the catheter. A straight conductor reduces the risk of injury and subsequent thrombosis, and allows for a number of manipulations to correct the position of the catheter. The dilator is used to extend the puncture of the skin. Also, if necessary, further expansion of the puncture site used a scalpel.

The Floswitch™ Introducer is used to enter conductors used with central or peripheral venous catheters.

Secondary Fixation Device the catheter is designed to securely hold the catheter and attach it to the skin to significantly reduce the risk of mechanical damage to the inner wall of the vessel.

The Guiding Syringe is included in some kits and is intended for insertion of conductors (allows insertion of a conductor into the vessel lumen through the piston), this method reduces the number of manipulations when installing a catheter, eliminates the likelihood of needle dislocation, effectively prevents air embolism, minimizes blood loss and physician contact with the patient's blood.


The connector is required to create a tightness in the catheter.

A surgical sheet with a dedicated surgical field for the patient, which is part of some sets is necessary to close the surface of the patient's body, has a special opening to delimit the puncture site.

The rotating connector and the Floswitch™ flow interrupter are designed to be mounted on the skin and short-term overlap of the lumen of the outer portion of the catheter prevents air embolism or bleeding when the catheter is manipulated.

The Central Venous Catheter and kit are individually packed into a rigid thermoformed blister in which the components are securely held. The Catheterization kits are supplied sterile and non-pyrogenic for single use only.

catheters scheme shown below

Picture 1 - The Central Venous Catheter Careflow™ one lumen (without extender)


Figure 2 - The Central Venous Catheter Careflow™ / Hydrocath Assure™ one lumen (with extension)


Figure 3- The Central Venous Catheter Careflow™ / Hydrocath Assure™ with two lumens


Figure 4- The Central Venous Catheter Careflow™ / Hydrocath Assure™ three lumens


Figure 4- The Central Venous Catheter Careflow™ with four lumens


Figure 5- The Central Venous Catheter Careflow™ with five lumens

Kits vascular catheters Careflow™, Hydrocath Assure™ have different length, diameter, number of lumens. The following types of lumens: distal, proximal 1, 2 proximal, middle 1, midline.

The number of gaps introduced depends on the compatibility of fluids, as well as the functionality (e.g., blood collection / SCN / PM infusion, pressure monitoring, large volume of infusion).

2. The method of sterilization of medical products

The company "Merit Medical Singapore" sterilizes the product described in this technical documentation file, ethylene oxide (EO). EO sterilization cycle validation is performed in accordance with the international standard ISO 11135. Annual re-validation evaluates the worst case loading on the basis of reference density packaging configuration that allows for sterilization standard load pallet mixed. The products of "Merit Medical Singapore" meet the requirements of ISO 10993-7 standard for residual content of ethylene oxide after sterilization. Minutes and reports on the validation of the sterilization process are stored in the company's document control system "Merit Medical Singapore."

This procedure involved contractor, certified according to international harmonized standards for sterilization of medical devices.

The new products are evaluated and, when justified, shall be approved and implemented in the existing EO sterilization cycles in accordance with the applicable standard AAMI TIR28 «Approval of the product and the process equivalent to ethylene oxide sterilization." These reports are stored in the department of sterilization / laboratory services.

Sets of vascular catheters CareflowTM, Hydrocath AssureTM products are less complex than the standard process control device (SAR), which is the worst-case configuration, and easily incorporated into a validated sterilization cycle of "Merit Medical Singapore." Forsets vascular catheters CareflowTM, Hydrocath AssureTM It does not provide for re-sterilization user forces.

The sterilization cycle consists of the steps incorporated in one cycle: exposure within the chamber, sterilization and venting with nitrogen gas and air. These steps are conducted in a special chamber. Sterility of the product is verified using *Bacillus atrophaeus* bacteria (previously called *Bacillus subtilis*) as biological indicators. Regular monitoring of the level of bioburden, endotoxins, and bacterial content of residual ethylene oxide / ethylene chlorohydrin is conducted in accordance with the requirements of US Pharmacopoeia or methods manuals AAMI, ANSI and ISO.

The sterilization process was validated to match the guaranteed level of sterility 10^{-6} .

3. PRODUCT SPECIFICATIONS

3.1 central venous catheter Careflow TM in the set.

CareflowTM catheter with one lumen (Without extender)	
The outer diameter of catheter	2,5 F (0,90 mm); 3,0 F (1,10 mm); 4,0 F (1,45 mm); 5,0 F (1,75 mm); 7,0 F (2,40 mm)
The inner diameter of the lumen of the tube	22 G (0,70 mm); 20G (0,90 mm); 18G (1,25 mm); 16G (1,65 mm); 14G (2,10 mm)
The working length of the catheter	100 mm 150 mm 200 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy

Flow rate, ml / min, at least	For catheter 2,5 F, 100 mm long - 10 ml / min For catheter 3 F, 150 mm long - 16 ml / min For catheter 3 F, 200 mm long - 15 ml / min For catheter 4 F, 200 mm long - 36 ml / min For the catheter 5 F, 150 mm long - 71 ml / min For the catheter 5 F, 200 mm long - 63 ml / min For catheter 7 F, length 150 mm - 190 ml / min For catheter 7 F, length 200 mm - 136 ml / min
The configuration of the tip	tonchik distal end of the catheter is smooth, rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter	Not more than 15 N
CareflowTM catheter with one lumen (With extension)	
The outer diameter of catheter	2,5 F (0,90 mm); 3,0 F (1,10 mm); 4,0 F (1,45 mm); 5,0 F (1,75 mm); 7,0 F (2,40 mm)
Inner catheter diameter / lumen diameter tube	22 G (0,70 mm); 20G (0,90 mm); 18G (1,25 mm); 16G (1,65 mm); 14G (2,10 mm)
The working length of the catheter	60 mm 100 mm 150 mm 200 mm 300 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
The flow rate of at least	For catheter 2,5 F, length 200 mm - 7 ml / min For catheter 2,5 F, 60 mm - 12 ml / min For catheter 2,5 F, 100 mm long - 10 ml / min For catheter 3 F, 150 mm long - 18 ml / min For catheter 3 F, 200 mm long - 15 ml / min For catheter 4 F, 60 mm - 47 ml / min For catheter 4 F, 200 mm long - 35 ml / min For the catheter 5 F, 150 mm long - 65 ml / min For the catheter 5 F, 200 mm long - 60 ml / min For the catheter 5 F, 300 mm long - 46 ml / min For catheter 7 F, length 150 mm - 159 ml / min For catheter 7 F, length 200 mm - 149 ml / min For catheter 7 F, 300 mm - 136 ml / min
The configuration of the tip	Tip of distal end of the catheter should be smooth rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage

Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter	Not more than 15 N
catheter Careflow™ with two lumens	
The outer diameter of catheter	4,0 F (1,45 mm); 5,0 F (1,75 mm); 7,0 F (2,40 mm)
The inner diameter of the lumen of the tube	For catheter diameter 4,0 F: 22 G (0,70 mm); 20G (0,90 mm); For catheter diameter 5,0 F: 18G (1,25 mm); 20G (0,90 mm); For catheter diameter 7,0 F: 16G (1,65 mm); 18G (1,25 mm)
The working length of the catheter	60 mm 100 mm 150 mm 200 mm 300 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
The flow rate of at least	For catheter 4 F, 60 mm long distal lumen - 20 ml / min, proximal lumen 1 - 8 ml / min For catheter 4 F, length 100 mm, lumen distal - 19 ml / min, proximal lumen 1 - 6 ml / min For catheter 4 F, length 150 mm, lumen distal - 16 ml / min, proximal lumen 1 - 4 ml / min For catheter 4 F, length 200 mm, lumen distal - 17 ml / min, proximal lumen 1 - 4 ml / min For the catheter 5 F, 60 mm long distal lumen - 20 ml / min, proximal lumen 1 - 8 ml / min For the catheter 5 F, length 150 mm, lumen distal - 29 ml / min, proximal lumen 1 - 6 ml / min For the catheter 5 F, length 200 mm, lumen distal - 30 ml / min, proximal lumen 1 - 6 ml / min For catheter 7 F, 150 mm in length, distal lumen - 74 ml / min, proximal lumen 1 - 24 ml / min For catheter 7 F, 200 mm in length, distal lumen - 73 ml / min, proximal lumen 1 - 21 ml / min For catheter 7 F, 300 mm in length, distal lumen - 61 ml / min, proximal lumen 1 - 19 ml / min
The configuration of the tip	Tip of distal end of the catheter should be smooth, rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter	Not more than 15 N

catheter Careflow™ three lumens	
The outer diameter of catheter	5,0 F (1,75 mm); 7,0 F (2,40 mm)
The inner diameter of the lumen of the tube	For catheter diameter 5,0 F: 18G (1,25 mm); 20G (0,90 mm); 20G (0,90 mm); For catheter diameter 7,0 F: 16G (1,65 mm); 18G (1,25 mm); 18G (1,25 mm)
The working length of the catheter	60 mm 100 mm 150 mm 200 mm 300 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
The flow rate of at least	For the catheter 5 F, 60 mm length: Distal lumen - 40 ml / min, the median clearance 1 - 10 ml / min, proximal lumen 1 - 12 ml / min For the catheter 5 F, 100 mm in length: The distal lumen - 32 ml / min, the median clearance 1 - 7 ml / min, proximal lumen 1 - 8 ml / min For the catheter 5 F, 150 mm in length: The distal lumen - 29 ml / min, the median clearance 1 - 5 ml / min, proximal lumen 1 - 6 ml / min For the catheter 5 F, 200 mm in length: The distal lumen - 26 ml / min, the median clearance 1 - 5 ml / min, proximal lumen 1 - 6 ml / min For catheter 7 F, 150 mm in length: The distal lumen - 66 ml / min, the median clearance 1 - 22 ml / min, proximal lumen 1 - 23 ml / min For catheter 7 F, 200 mm in length: The distal lumen - 59 ml / min, the median clearance 1 - 16 ml / min, proximal lumen 1 - 17 ml / min For catheter 7 F, 300 mm in length: The distal lumen - 53 ml / min, the median clearance 1 - 11 ml / min, proximal lumen 1 - 9 ml / min
The configuration of the tip	Tip of distal end of the catheter should be smooth rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter	Not more than 15 N
catheter Careflow™ with four lumens	
The outer diameter of catheter	8,5 F (2,80 mm)
The inner diameter of the lumen of the tube	14G (2,10 mm) 16G (1,65 mm); 16G (1,65 mm);

	18G (1,25 mm)
The working length of the catheter	150 mm 200 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
The flow rate of at least	For catheter 8,5 F, 150 mm in length: The distal lumen - 64 ml / min, the median clearance 1 - 82 ml / min, proximal lumen 1 - 20 ml / min, median lumen 2 - 57 ml / min For catheter 8,5 F, 200 mm in length: The distal lumen - 58 ml / min, the median clearance 1 - 76 ml / min, proximal lumen 1 - 19 ml / min Median lumen 2 - 45 ml / min
The configuration of the tip	Tip of distal end of the catheter should be smooth, rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter catheter Careflow™ with five lumens	Not more than 15 N
The outer diameter of catheter	9,5 F (3,22 mm)
The inner diameter of the lumen of the tube	14G (2,10 mm) 14G (2,10 mm) 16G (1,65 mm); 16G (1,65 mm); 18G (1,25 mm)
The working length of the catheter	150 mm 200 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
The flow rate of at least	For catheter 9,5 F, 150 mm in length: The distal lumen - 73 ml / min, the median clearance 1 - 61 ml / min, proximal lumen 1 - 19 ml / min, median lumen 2 - 65 ml / min, proximal lumen 2 - 21 ml / min For catheter 9,5 F, 200 mm in length: The distal lumen - 68 ml / min, the median clearance 1 - 58 ml / min, proximal lumen 1 - 13 ml / min, median lumen 2 - 60 ml / min, proximal lumen 2 - 12 ml / min
The configuration of the tip	Tip of distal end of the catheter should be smooth, rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction

The strength of the extraction catheter	Not more than 15 N
kit components	
guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger	Diameter of the guide wire 0,6 mm, length 450 mm; 650 mm Diameter of the guide wire 0,9 mm, length 450 mm, 500mm, 600 mm, 700mm
straight guide wire	The diameter of the guide wire: 0,5 mm guide wire length: 350 mm 450 mm
Introducer needle	Outer Diameter: 1,40 mm Full length of needle tube: 70 mm On the outer surface should be free of surface defects and process.
coaxial introducer	1) Outer sheath diameter of introducer sheath: 1.70 mm Introducer sheath effective length: 70 mm Outer diameter of the introducer needle (included in the introducer sheath): 1.24 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). 2) Outer sheath diameter of introducer sheath: 1.27 mm Introducer sheath effective length: 60 mm Outer diameter of the introducer needle (included in the introducer sheath): 0.85 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). 3) Outer sheath diameter of introducer sheath: 1.00 mm Introducer sheath effective length: 45 mm Outer diameter of the introducer needle (included in the introducer sheath): 0.5 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). 4) Outer sheath diameter of introducer sheath: 0.75 mm Introducer sheath effective length: 38 mm Outer diameter of the introducer needle (included in the introducer sheath): 0.40 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer).
Introducer FloSwitch™	1) Outer sheath diameter of introducer sheath: 2.00 mm Introducer sheath effective length: 70 mm Outer diameter of the introducer needle (included in the introducer sheath): 1.40 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N

n;
n,

defects

70 mm

ed in the

ffects on
h is part

1.27mm

ed in the

defects on
ich is part

: 1.00 mm

uded in the

defects on
hich is part

h: 0.75 mm

luded in the

e defects on
hich is part

th: 2.00mm

cluded in the

ce defects on
(which is part

nan 10 N

	2) Outer sheath diameter of introducer sheath: 1.70mm Introducer sheath effective length: 70 mm Outer diameter of the introducer needle (included in the introducer sheath): 1.20 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N 3) Outer sheath diameter of introducer sheath: 1.20mm Introducer sheath effective length: 60 mm Outer diameter of the introducer needle (included in the introducer sheath): 0.86 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N 4) Introducer sheath outer diameter: 1.00 mm Introducer sheath effective length: 45 mm Outer diameter of the introducer needle (part of the introducer sheath): 0.50 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N																					
vessel dilator	<table><tr><th>External diameter</th><th>Length products</th></tr><tr><td>4F (1,33 mm)</td><td>88 mm</td></tr><tr><td>5F (1,7 mm)</td><td>88 mm</td></tr><tr><td>6F (2,0 mm)</td><td>88 mm</td></tr><tr><td>8F (2,69 mm)</td><td>110 mm</td></tr><tr><td>9F (3,0 mm)</td><td>110 mm</td></tr><tr><td>10F (3,33 mm)</td><td>110 mm</td></tr><tr><td>The inner diameter of the distal (narrowed) end of dilator</td><td>For dilator 4F: 0,66 mm for dilator 5F: 0,66 mm for dilator 6F: 0,66 mm for dilator 8F: 0,96 mm for dilator 9F: 0,96 mm for dilator 10F: 0,96 mm</td></tr><tr><td>The geometric shape of the distal end of the dilator</td><td>The presence of tapering, polished edge. Edge should be sharp. Without cracks, scratches.</td></tr><tr><td>radio-opacity</td><td>radiopaque</td></tr></table>		External diameter	Length products	4F (1,33 mm)	88 mm	5F (1,7 mm)	88 mm	6F (2,0 mm)	88 mm	8F (2,69 mm)	110 mm	9F (3,0 mm)	110 mm	10F (3,33 mm)	110 mm	The inner diameter of the distal (narrowed) end of dilator	For dilator 4F: 0,66 mm for dilator 5F: 0,66 mm for dilator 6F: 0,66 mm for dilator 8F: 0,96 mm for dilator 9F: 0,96 mm for dilator 10F: 0,96 mm	The geometric shape of the distal end of the dilator	The presence of tapering, polished edge. Edge should be sharp. Without cracks, scratches.	radio-opacity	radiopaque
External diameter	Length products																					
4F (1,33 mm)	88 mm																					
5F (1,7 mm)	88 mm																					
6F (2,0 mm)	88 mm																					
8F (2,69 mm)	110 mm																					
9F (3,0 mm)	110 mm																					
10F (3,33 mm)	110 mm																					
The inner diameter of the distal (narrowed) end of dilator	For dilator 4F: 0,66 mm for dilator 5F: 0,66 mm for dilator 6F: 0,66 mm for dilator 8F: 0,96 mm for dilator 9F: 0,96 mm for dilator 10F: 0,96 mm																					
The geometric shape of the distal end of the dilator	The presence of tapering, polished edge. Edge should be sharp. Without cracks, scratches.																					
radio-opacity	radiopaque																					
The rotating connector Floswitch™	Length: 8 mm Width 11 mm Height 13 mm																					
Clamp the catheter	1 embodiment clip: Length 25 mm The width of 9 mm 2 embodiment clip: Length 18 mm width 6 mm																					

	height 12 mm
A syringe volume of 5 ml	- nominal volume of the syringe: 5 ml; - purity of the syringe surfaces: a syringe the surface must be clean and free of foreign particles by visual inspection; - scale must be calibrated graduations. the calibration line must have digital values corresponding to the capacity values (1 ml, 2ml, 3ml, 4 ml, 5 ml) - division of the scale of the syringe – 0,2 ml - scale length 45 mm
Guiding Syringe	- nominal volume of the syringe: 5 ml; - purity of the syringe surfaces: a syringe the surface must be clean and free of foreign particles by visual inspection; - the speed of aspiration of water: ≥ 4 g / min - scale must be calibrated graduations. the calibration line must have digital values corresponding to the capacity values (1 ml, 2ml, 3ml, 4 ml, 5 ml) - division of the scale of the syringe – 0,2 ml - scale length 37 mm
Scalpel	Working length scalpel 24,0 mm (cutting portion) The total length of the scalpel 147,0 mm The width of the cutting edge - less than 3 microns The connection strength of the polymer handle and the blade is not less than 20 N
Surgical sheet with a dedicated surgical field for the patient	Length 914 mm Width 609 mm Diameter of the surgical field: 102 mm Material Density: 0.70 g / cm^3 Material strength: not less than 5 N Adhesive Layer: None White colour

1.2 The catheter is a central venous Hydrocath Assure™ in the set:

Catheter Hydrocath Assure™ with a clearance	
The outer diameter of catheter	5,0 F (1,75 mm);
The inner diameter of the catheter tube lumen	16G (1,65 mm);
The working length of the catheter	200 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
Flow rate, ml / min	56 ml / min
The configuration of the tip	Tip of distal end of the catheter is smooth, rounded, without bends.

Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter	Not more than 15 N
Catheter Hydrocath Assure™ with two lumens	
The outer diameter of catheter	7,0 F (2,40 mm)
The inner diameter of the catheter tube lumen	16G (1,65 mm); 18G (1,25 mm)
The working length of the catheter	150 mm 200 mm 300 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type
The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
flow rate	For catheter 7 F, 200 mm in length, distal lumen - 69 ml / min, proximal 1 - 22 ml / min For catheter 7 F, 300 mm in length, distal lumen - 61 ml / min, proximal 1 - 17 ml / min For catheter 7 F, 150 mm in length, distal lumen - 71 ml / min, proximal 1 - 23 ml / min For catheter 7 F, 200 mm in length, distal lumen - 69 ml / min, proximal 1 - 22 ml / min
The configuration of the tip	Tip of distal distal end of the catheter should be smooth, rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter	Not more than 15 N
Catheter Hydrocath Assure™ three lumens	
The outer diameter of catheter	7,0 F (2,40 mm)
The inner diameter of the catheter tube lumen	16G (1,65 mm); 18G (1,25 mm); 18G (1,25 mm)
The working length of the catheter	150 mm 200 mm 300 mm
Compatible connection Luer-type	The hub of the catheter should be compatible with rotary adapters and fixed parts of the male type

The distance between the marks depth	10 mm
radio-opacity	The catheter should be viewed under fluoroscopy
The flow rate of at least	For catheter 7 F, 150 mm in length, distal lumen - 68 ml / min, proximal 1 - 24 ml / min, median 1 - 22 ml / min For catheter 7 F, 200 mm in length, distal lumen - 58 ml / min, proximal 1 - 19 ml / min, median 1 - 18 ml / min For catheter 7 F, 300 mm in length, distal lumen - 50 ml / min, proximal 1 - 13 ml / min median 1 - 13 ml / min
The configuration of the tip	Tip of distal end of the catheter should be smooth, rounded, without bends.
Surface	On the outer surface should be free of surface defects and process.
Resistance to liquid leakage	Cannula, fitting joint assembly and other parts of the catheter shall not allow fluid leakage
Air Leak Resistance	Cannula, fitting assembly and other parts of the catheter shall prevent air leakage during suction
The strength of the extraction catheter	Not more than 15 N
kit components	
guide wire with J-shaped tip of the device for the passage of the catheter Venaguide Seldinger	The diameter of the guide wire: 0,9mm guide wire length: 450 mm, 500 mm, 600 mm, 700 mm
Device for fixation of the catheter	option 1: Length: 25 mm Maximum width: 25 mm Pitch: 5 mm Option 2: Length: 25 mm Maximum width: 30 mm Height 6 mm
Introducer needle	Outer Diameter: 1,40 mm Full length of needle tube: 70 mm On the outer surface should be free of surface defects and process.
Introducer FloSwitch™	1) Outer sheath diameter of introducer sheath: 2.00mm Introducer sheath effective length: 70 mm Outer diameter of the introducer needle (included in the introducer sheath) 1.40 mm There should be no technological or surface defects on the outer surface the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N 2) Outer sheath diameter of introducer sheath: 1.70mm Introducer sheath effective length: 70 mm Outer diameter of the introducer needle (included in the introducer sheath) 1.20 mm There should be no technological or surface defects on the outer surface the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N

	3) Outer sheath diameter of introducer sheath: 1.20mm Introducer sheath effective length: 60 mm Outer diameter of the introducer needle (included in the introducer sheath): 0.86 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N 4) Introducer sheath outer diameter: 1.00 mm Introducer sheath effective length: 45 mm Outer diameter of the introducer needle (part of the introducer sheath): 0.50 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer). Effort of switching of a stream: no more than 10 N															
coaxial introducer	Introducer sheath outer diameter: 1.70 mm Introducer sheath effective length: 70 mm Outer diameter of the introducer needle (included in the introducer sheath): 1.24 mm There should be no technological or surface defects on the outer surface of the introducer needle (which is part of the introducer).															
The rotating connector Floswitch TM	Length: 8 mm Width 11 mm Height 13 mm															
vessel dilator	<table><tr><th>External diameter</th><th>Length products</th></tr><tr><td>6F (2,0 mm)</td><td>88 mm</td></tr><tr><td>8F (2,69 mm)</td><td>110 mm</td></tr><tr><td>10F (3,33 mm)</td><td>110 mm</td></tr><tr><td>The inner diameter of the distal (narrowed) end of dilator</td><td>For dilator 6F: 0,66 mm For dilator 8F: 0,96 mm For dilator 10F: 0,96 mm</td></tr><tr><td>The geometric shape of the distal end of the dilator</td><td>The presence of tapering, polished edge. Edge should be sharp. Without cracks, scratches.</td></tr><tr><td>radio-opacity</td><td>radiopaque</td></tr></table>		External diameter	Length products	6F (2,0 mm)	88 mm	8F (2,69 mm)	110 mm	10F (3,33 mm)	110 mm	The inner diameter of the distal (narrowed) end of dilator	For dilator 6F: 0,66 mm For dilator 8F: 0,96 mm For dilator 10F: 0,96 mm	The geometric shape of the distal end of the dilator	The presence of tapering, polished edge. Edge should be sharp. Without cracks, scratches.	radio-opacity	radiopaque
External diameter	Length products															
6F (2,0 mm)	88 mm															
8F (2,69 mm)	110 mm															
10F (3,33 mm)	110 mm															
The inner diameter of the distal (narrowed) end of dilator	For dilator 6F: 0,66 mm For dilator 8F: 0,96 mm For dilator 10F: 0,96 mm															
The geometric shape of the distal end of the dilator	The presence of tapering, polished edge. Edge should be sharp. Without cracks, scratches.															
radio-opacity	radiopaque															
connector	Connector size: Length 8 mm The outer diameter of 10 mm															
Clamp the catheter	1 embodiment clip: Length 25 mm The width of 9 mm 2 embodiment clip: Length 18 mm width 6 mm height 12 mm															
A syringe volume of 5 ml	- nominal volume of the syringe: 5 ml; - purity of the syringe surfaces: a syringe the surface must be clean and free of foreign particles by visual inspection; - scale must be calibrated graduations. the calibration line must have digital values corresponding to the capacity values (1 ml, 2ml, 3ml, 4 ml, 5 ml) - division of the scale of the syringe – 0,2 ml - scale length 45 mm															
Guiding Syringe	- nominal volume of the syringe: 5 ml;															

	- purity of the syringe surfaces: a syringe the surface must be clean and free of foreign particles by visual inspection; - scale must be calibrated graduations. the calibration line must have digital values corresponding to the capacity values (1 ml, 2ml, 3ml, 4 ml, 5 ml) - division of the scale of the syringe – 0,2 ml - scale length 37 mm
--	---

4. DESCRIPTION OF PACKAGING ARTICLES WITH THE AMOUNT IN PACKAGING

Standard warehouse storage is acceptable for Merit products unless specific storage conditions are required as noted in labeling and as driven by the clinical risk assessment and/or FMEA. The products are packaged and labelled as sterile, single use devices. Sterile barrier systems are specified per ISO11607 to prevent ingress of microorganisms and allow aseptic presentation of products at point of use. Package integrity testing is performed as part of the production process. Device Expiration of Careflow™ Central Venous Catheter is Five (5) years and Hydrocath Assure™ Central Venous Catheter is Two (2) years.

Packaging Materials and Configuration

The Hydrocath Assure™ Central Venous Catheters and kit components are packaged in a Tyvek pouch. The Careflow Central Venous Catheters and kit components are packaged in a rigid thermoformed bottom web tray with Tyvek Lid. Multiple units of the sealed package are packed into a shelf-box, which is subsequently packed into a corrugated fiber board shipping carton for Ethylene Oxide sterilization.

The packaging materials are shown as below:

Packaging Material Specifications

Components	Material
Packaging Pouch	Tyvek , High-density Polyethylene
Tray	Polyethylene Terephthalate Glycol (PETG)

b) Packaging Qualification

Packaging qualification of Merit Medical products is conducted to the requirements of the following standards:

ISO 11607: Packaging for terminally sterilized medical devices

ISTA-2A: Partial simulation performance test

The standard Careflow™ and Hydrocath Assure™ CVC kit models are packed in thermoformed bottom web tray with Tyvek lid using Form Fill Seal (FFS3) packaging equipment and ancillary systems. Packaging for this product was qualified according to approved protocols and reports.

Note that the PVC inner/outer trays have been phased out in year 2008.

The labels of the standard packaging are printed and adhere to the Tyvek Top.

c) Shipping Qualification

A model of Central Venous Catheter (Hydrocath Assure™) that has longest length, highest complexity and highest weight was selected as the representative test article in the shipping qualification study to represent the models for the Hydrocath Assure™ and Careflow™ Central Venous Catheter. The representative test articles of Hydrocath Assure™ CVC were sealed in their unit packaging, packed in the shelf-box and then in the corrugated fiber board shipping carton per the assembly procedure. Samples were inspected pre-sterilized and post-sterilized (3x sterilization) after the simulated shipment tests.




5. INFORMATION ON LABELING OF PRODUCTS



Label Formats products have been developed and approved for the individual articles, sterile packaging and the transport packaging. The approved labeling regulations is in the process of the product.

The company "Merit Medical Singapore Pte. Ltd." Uses the following standards for the maintenance and formation of marking symbols in accordance with the current requirements of the Medical Devices Directive:

- EN 980: symbols used in the labeling of medical devices.
- EN 1041: Information supplied by the manufacturer.
- ISO 15223-1: medical products. Symbols used on labels for labeling of medical devices and in the information provided - Part 1: General requirements.
- ISO 15223-2: medical products. The symbols used on medical devices, labeling and the information provided - Part 2: Design, selection and validation of characters.
- FDA guidance on the labeling of medical devices for patients; final guidance for industry and FDA reviews (19 April 2001)

Description of the marking symbol is represented in the table.

	Symbol	Decoding the symbol
1.		Authorized Representative in the European Community
2.		Manufacturer
3.		It sold only to doctors and hospitals
4.		A ban on re-use.
5.		REF number
6.		use up
7.		lot number
8.		sterilized using ethylene oxide
9.		Attention! Refer to the instructions for use
10.		pyrogen-free
11.		Do not re-sterilize
12.		CE mark

13.		Protect from moisture
14.		Hooks do not take Do not turn over
15.		Hooks do not take Do not turn over
16.		Fragile, handle with care
17.		Do not use if packaging is damaged
18.		Hooks do not take
19.		Do not turn over
20.		Number of lumen (lumen 3)
21.		flow rate
22.		Refer to the instructions for use.
23.	The inscription "non-toxic"	
24.	Date of manufacture / sterilization	

6. CONDITIONS, storage and transportation.

Storage conditions:

Temperature + 5 ° C to +27 ° C (environment)

Relative humidity 10% to 90% (non-condensing)

Atmospheric pressure 84.0 kPa - 106.7 kPa (630 - 800 mm Hg)

Transportation Conditions:

Temperature -35 ° C to +55 ° C (environment)

Relative humidity 10% to 90% (non-condensing)

Atmospheric pressure 84.0 kPa - 106.7 kPa (630 - 800 mm Hg)

Ambient temperature not exceeding 25 ° C.

Temperature range for use - the body temperature of 32 to 42 ° C

Relative humidity: 100% (blood contact)

Atmospheric pressure 86-106 kPa.

7. SHELF LIFE.

Shelf life catheters venous Careflow™ in sets of five (5) years, and the shelf life of catheters central venous Hydrocath Assure™ - 2 (two) years.

Based on the European regulirovanii symbol "Use By" or "date of manufacture" can be used for sterile products. In the case of "use by" date is indicated, in which the product can be used (expiration date), after that date, it is impossible to use the product. In this case, to determine the date of production is necessary to take the general term storage of the product, manufacturer's declaration of "Use By" date.

8. RECYCLING

Used sets disposed of as hazardous waste epidemiologically class B (the presence of contact with patient's blood and mucous tissues).

At the end of the shelf life of unused sets, not in contact with blood and mucous membranes should be disposed of as hazardous waste epidemiologically class B close in composition to the solid waste. Dispose of used products shall be in accordance with the requirements of laws and regulations on the disposal of products made in the country of use of the medical device.

9. WARRANTY

Limited Warranty and Disclaimer

The company "Merit Medical Singapore Pte. Ltd." Guarantees the conformity of the product specifications at the time of shipment to the Buyer, and be free from defects in materials and workmanship for the period specified in the written limited warranty supplied separately to each product. Warranty period of storage catheters, venous CareflowTM in sets of five (5) years, and the catheter central venous Hydrocath AssureTM - 2 (two) years. This warranty does not apply to products which were: (i) modified, altered, repaired, recovery, recycling or other changes in the company, in addition to "Merit Medical Singapore Pte. Ltd."; (ii) misused, subjected to improper treatment, reuse or resterilization, suffered an accident, violation of the rules, negligence or tampering; (iii) is damaged by physical impact or electrical nature, as well as altered, damaged or removed serial number; (iv) used in conjunction with any other article; (v) used in accordance with the instructions for use, is installed an authorized organization, used in accordance with the field of application specified in the specification, or under conditions not foreseen for this kind of products. Any technical advice on behalf of the company "Merit Medical Singapore Pte. Ltd." Provided without liability or compensation, the company "Merit Medical Singapore Pte. Ltd." Does not assume any obligation or responsibility, all these recommendations are given and accepted under the responsibility of the buyer. Save as aforesaid, the company "Merit Medical Singapore Pte. Ltd." MAKES REPRESENTATIONS OR WARRANTIES OF ANY KIND WITH RESPECT TO PRODUCT, WHETHER EXPRESS OR IMPLIED, buyers and end users to third parties. "Merit Medical Singapore Pte. Ltd." DISCLAIMS AND EXCLUDES ANY AND ALL EXPRESS, IMPLIED, OR LEGALLY UNDER WARRANTY, INCLUDING THE WARRANTIES AGAINST INFRINGEMENT OR IMPLIED WARRANTIES OF MERCHANTABILITY OR FITNESS

10. CLAIM

An organization that adopts the territory of the Russian Federation claims from consumers on product quality:

Limited Liability Company "Merit Technologies" (LLC "Merit Technologies"),

BIN 1117746679533, INN 7705960130

Legal address: Of the Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

Tel / fax: +7 (495) 221-89-02

Instructions for use
Sets of vascular catheters Careflow™

Description

Sets of vascular catheters Careflow™ includes Careflow™ Catheter is a radio-opaque polyurethane catheter.

Intended Use

The Central Venous Catheter Careflow™ is an intravenous catheter that is intended for short term use (no longer than 30 days) to access the human circulatory system via insertion through either the internal Jugular or Sub-Clavian Vein using Seldinger technique, where the catheter tip resides in the superior vena cava. The Central Venous Catheter Careflow™ is intended for infusion of drugs, Total Parental Nutrition (TPN) fluids, large volume infusions, repeated blood sampling and to monitor the central venous pressure.

Contraindications

Percutaneous puncture of veins/ arteries may be contraindicated in patients with pulmonary hypertension.

Patients with suspected hypersensitivity to nickel must pass the test to assess the sensitivity of the skin before using the conductors at the location of central venous catheters.

Note:

Use only as directed by a physician.

A full list of components is shown on the package lid.

General Warnings

1. Physicians must be familiar with the complications associated with central venous catheterisation, i.e. vessel perforation, air embolism, catheter embolism, pleural and mediastinal damage, septicaemia, thrombosis and cardiac tamponade secondary to vessel wall or atrial perforation.
2. Complications are associated with right atrial and inadvertent right ventricular catheterisation. Physicians must be aware of these complications before advancing the catheter beyond the depth required for normal vena cava placement. Do not advance the catheter past this depth unless procedure requires right atrial placement. If catheter is advanced beyond normal vena cava placement depth, monitor electrocardiogram during insertion and confirm final position by chest X-Ray.
3. It is recommended that patients be placed in a slight Trendelenburg position during insertion procedure to reduce the risk of air embolism.
4. We recommend that Careflow™ Catheter lumens (Central Venous, Peripheral Venous Catheters and Radial, Femoral Artery Catheter) are flushed with sterile saline solution prior to catheter insertion.
5. All catheter placements must be inspected for flow rates, security of dressing and security of luer connections.
6. To reduce or eliminate the potential for catheter migrations, we advise that every catheter placement be secured by suturing at the eyelets of the junction boot housing/hub and that, where the use of the secondary fixation device is necessary, it should be used as additional support and not the only means of fixation. Additionally, the security of catheter fixation and position of the catheter tip should be checked throughout use.
7. When removing dressings at or close to catheter sites, care must be taken to avoid severing the catheter.

8. Acetone must not come into contact with the catheter as the material may weaken and this may result in leakage or aspiration.
9. Exposure of product componentry to topicals containing alcohol is not recommended.
10. Do not attempt to re-insert a partially or completely withdrawn introducer cannula.
11. Percutaneous puncture of a central vein may be contraindicated in patients with pulmonary hypertension.
12. Use of a syringe smaller than 5 mL to irrigate or de-clot an occluded catheter may cause intraluminal or catheter rupture.
13. Syringes are supplied for blood aspiration only.
14. Luer connections: as standard practice the security of luer connections must be checked routinely.
15. Physicians should be aware that central venous catheters are intended for use up to 30 days.
16. Patients with suspected hypersensitivity to nickel should undergo skin test to assess hypersensitivity prior to use of Guidewires in the placement of central venous catheters.

Warnings - Seldinger Technique

1. Do not withdraw the guidewire against needle bevel, as this increases the risk of severing the guidewire.
2. During insertion do not reinsert a partially or completely withdrawn needle into the cannula.
3. Ensure the flexible end of the guidewire is advanced into the vein.
4. Ensure the guidewire moves freely in the needle introducer.
5. A firm grip must be maintained on the guidewire at all times.
6. When using the 'J' wire straightener maintain a firm grip on the plastic sleeve.
7. Ensure the dilator is removed prior to catheter advancement.
8. The moveable suture devices are designed as additional support and must not be used as the only means of fixation.
9. Potential for guidewire breakage. Although the incidence of guidewire breakage is extremely uncommon, physicians must be aware of the potential of guidewire breakage if undue force is applied to the wire. If resistance is met when attempting to remove the guidewire after central venous placement, the wire may be kinked within the area of the catheter tip and the vessel. Undue force may cause the wire to break. If resistance is encountered, withdraw the catheter relative to the guidewire (2- 3 cm) and attempt to remove the wire. If resistance is still apparent, remove the wire and the catheter simultaneously.
10. Physicians should be aware that the guidewire can pick up material from the vein. This may prevent the guidewire from being withdrawn through the catheter.
11. Do not force the guidewire. If resistance is met, carefully withdraw the guidewire and attempt insertion.
12. If femoral approach is used, place patient in supine position for insertion procedure.

Suggested Procedure Seldinger Technique (Venous Access)

1. Prepare the insertion site using the full aseptic technique required for venous access.
2. The guidewire may be introduced either by using a thin walled steel needle or introducer cannula.
3. Feed the desired flexible end of the guidewire through the introducer into the vein. If a 'J' wire is to be used, the 'J' can be straightened prior to insertion by using the plastic inserter sleeve.
4. Remove the introducer.
5. A vessel dilator may be used to enlarge the cutaneous puncture site. If further enlargement of the puncture site is necessary, use a scalpel.
6. When using multi lumen catheters, lumens other than distal should be flushed and attached to a desired fluid administration set or, alternatively, heparin locked as standard hospital practice.

1. Release
2. Straighten
3. Insert to depth.

WARNING
To prevent

7. Pass the distal tip of the catheter over the guidewire. (Sufficient guidewire length must remain exposed at the hub end of the catheter to maintain a firm grip on the guidewire.) Grasping the catheter near the skin, push the catheter into the vein with a slight twisting motion and advance into the final position.

8. Hold the catheter in position, withdraw the guidewire and aspirate with a syringe to ensure correct placement.


9. The distal hub should be connected to the appropriate fluid administration set. If a hub with Floswitch™ is being used this can now be switched off.

10. The catheter may now be secured by suturing eyelets of the junction boot housing/hub to the skin.

11. Apply sterile dressing as appropriate.

Secondary Fixation Device (7Fr and 9.5Fr catheters)

Note: This device is found in certain models only.

1. Ensure the Secondary Fixation Device is in place at the point of suture.


2. To close the device, snap the top plate in place.

3. To ensure fixation, suture at each eyelet on the wings of the Secondary Fixation Device.

4. To open the device, hold the lock device down and pull the plate up, using the lever.

Venaguide

Note: This device is found in certain models only.

1. Release guidewire by removing the guidewire cap.

2. Straighten guidewire 'J' by retracting into introducer system with thumb.

3. Insert into hub of introducer needle and advance guidewire into vein. Advance to required depth.

WARNING:

To prevent over-insertion stop when last depth mark reaches the venepuncture site.

Guiding Syringe

Note: This accessory is found in certain models only. This device is used for introducing guidewires for use with Merit central or peripheral venous catheters.

1. Insert introducer needle attached to Guiding Syringe into vessel and aspirate.
2. Feed the desired flexible end of the guidewire through the hole in rear of Guiding Syringe plunger into the vein. If a 'J' wire is to be used, straighten guidewire 'J' by retracting into introducer system with thumb.
3. Hold the guidewire in place and remove introducer needle and Guiding Syringe.

WARNING: Do not aspirate with guidewire in place or air may enter the syringe

Caution:

Do not reinfuse blood to minimize the risk of blood leakage from rear of syringe.

Floswitch™ Introducer

This device is used for introducing guidewires for use with central or peripheral venous catheters.

Using a Floswitch™ introducer cannula, puncture the vein, aspirate and advance the cannula into the vein and remove the needle. At this stage the Floswitch™ may be switched off. The switch must be opened again to facilitate guidewire or catheter passage.

Rotatable connector Floswitch™

This device may be attached to a luer hub of a central or peripheral venous catheter.

1. As standard practice the catheter should be secured to the skin.
2. To attach the Floswitch™, insert the Floswitch™ luer into the catheter hub. Tighten rotating collar ensuring that 'hand tight' connection has been made and that the catheter hub is fully advanced into the Floswitch™ collar (see figure 1.1).
3. The Floswitch™ should be fixed in position by either suturing or taping the wings. The recommended procedure for taping is shown in figure 2.
4. The black markings indicate flow status. When visible the catheter is open and there is flow. When covered, the catheter is closed.

Figure 1. Floswitch™ Inspection Procedure



Figure 2. Floswitch™ Taping Procedure

Warnings - Floswitch™ Products

1. The Floswitch™ should not be switched off (black marks covered) before the needle has been fully withdrawn. This applies to Floswitch™ introducers only.
2. The Floswitch™ must not be switched off (black marks covered) until the guidewire has been fully withdrawn.
3. When the catheter is not in use for infusion or aspiration, the Floswitch™ must be switched off and a suitable luer cap locked into the hub.
4. When using an intermittent injection bung attached to a Floswitch™, only use short needles. Ensure Floswitch™ is in 'ON' position prior to injecting or aspirating. Do not insert the needle more than 8 mm into the Floswitch™.
5. As standard practice, the security of the luer connection must be checked routinely. This is essential when lubricious substances such as Intralipids are being used. This applies to Floswitch™ luer lock attachment only.

Storage Conditions Store in cool dry place away from direct sunlight.

STERILE and non-pyrogenic in unopened, undamaged package. For single use only. Check integrity of the individual package before use. Dispose of product after use. Do not resterilize. Re-use may lead to infection or other illness/ injury. For reordering information or assistance please contact local representative.

Instructions for use
Sets of vascular catheters Hydrocath Assure™

In
The
sh
the
ca
inf
sar
Ca
few
ven
col
sign
cath
cath
insu
cath

Con

1. No
Hydr
paedi
than 2
2. On
3. Thi
benza
4. NO
been c
shown
are not
5. No c
Hydroc
through
be avoi

Description

Sets of vascular catheters Hydrocath Assure™ includes Hydrophilic Anti- Microbial Central Venous Catheter.

Intended Use

The Hydrocath Assure™ Central Venous Catheter is an intravenous catheter that is intended for short term use (no longer than 30 days) to access the human circulatory system via insertion through either the internal Jugular or Sub-Clavian Vein using Seldinger technique, where the catheter tip resides in the superior vena cava.

The Central Venous Catheter is intended for infusion of drugs, Total Parental Nutrition (TPN) fluids, large volume infusions, repeated blood sampling and to monitor the central venous pressure.

The Hydrocath Assure™ central venous catheter has a hydrophilic (Poly Vinyl Pyrrolidone), anti-microbial surface which is intended to help prevent colonisation of the catheter.

Anti-Microbial Activity: The antiseptic, benzalkonium chloride is integrated into the catheter's hydrogel coating system, internal and external to the catheter lumens.

Indications:

The Central Venous Catheter Hydrocath Assure™ is an intravenous catheter that is intended for short term use (no longer than 30 days) to access the human circulatory system via insertion through either the internal Jugular or Sub-Clavian Vein using Seldinger technique, where the catheter tip resides in the superior vena cava. The Central Venous Catheter is intended for infusion of drugs, Total Parental Nutrition (TPN) fluids, large volume infusions, repeated blood sampling and to monitor the central venous pressure. The Hydrocath Assure™ Central Venous Catheters are intended to assist in the prevention of colonisation of the catheter during the first few days after insertion. In a clinical trial which compared Hydrocath Assure™ with a central venous catheter not coated with hydromer or benzalkonium chloride the incidence of colonisation of the tip and subcutaneous segments by coagulase-negative staphylococci was significantly lower for Hydrocath Assure™. This significant benefit was observed only for catheters which had been in situ for less than 5 days. Hydrocath Assure™ central venous catheters have not been shown to reduce the rate of infection at the insertion site. There are insufficient data to conclude that Hydrocath Assure™ central venous catheters reduce the rate of catheter-related systemic infections.

Contraindication for use:

1. NOT FOR PAEDIATRIC USE Due to a risk of immediate-type hypersensitivity reaction, Hydrocath Assure™ catheters are not for use in paediatric patients. Based on reports received, paediatric patients are most likely to experience such a reaction, with lower body weight (less than 20 kg) and younger (less than 6 years) patients appearing to be at highest risk.
2. Only one Hydrocath Assure™ should be used per patient at any one time.
3. This product is contraindicated for use in patients with a known hypersensitivity (allergy) to benzalkonium chloride or any benzalkonium chloride-containing product.
4. NOT RECOMMENDED FOR USE FOR PREGNANT WOMEN No controlled studies have been conducted in pregnant women with Hydrocath Assure™ catheters. Animal studies have shown that Benzalkonium Chloride is embryo lethal and teratogenic, therefore, these catheters are not recommended during pregnancy.
5. No controlled studies have been conducted to investigate potential incompatibilities between Hydrocath Assure™ and intravenous fluids. Literature searches suggest that administration through the catheter of Bicarbonate Solutions, Nitrates, organic acids and Sulphonamides should be avoided.

6. Percutaneous puncture of a central vein may be contraindicated in patients with pulmonary hypertension. 7. Patients with suspected hypersensitivity to nickel must pass the test to assess the sensitivity of the skin before using the conductors at the location of central venous catheters.

Note:

Use only as directed by a physician.

A full list of components is shown on the package lid.

General Warnings

1. Rare cases of immediate-type hypersensitivity reactions have been reported when Hydrocath Assure™ catheters are inserted into patients. As a precautionary measure, the possibility of such reactions must be considered when using these catheters.
2. Physicians must be familiar with the complications associated with central venous catheterisation, i.e. vessel perforation, air embolism, catheter embolism, pleural and mediastinal damage, septicaemia, thrombosis and cardiac tamponade secondary to vessel wall or atrial perforation.
3. Complications are associated with right atrial and inadvertent right ventricular catheterisation. Physicians must be aware of these complications before advancing the catheter beyond the depth required for normal vena cava placement. Do not advance the catheter past this depth unless the procedure requires right atrial placement. If catheter is advanced beyond normal vena cava placement depth, monitor electrocardiogram during insertion and confirm final position by chest X-Ray.
4. It is recommended that patients be placed in a slight Trendelenburg position during insertion procedure to reduce the risk of air embolism.
5. We recommend that Merit Central Venous Catheter lumens are flushed with sterile saline solution prior to catheter insertion.
6. All catheter placements must be inspected for flow rates, security of dressing and security of luer connections.
7. To reduce or eliminate the potential for catheter migration, we advise that every catheter placement be secured by suturing at the eyelets of the junction boot housing/hub and that, while the use of the secondary fixation device is necessary, it should be used as additional support and not the only means of fixation. Additionally, the security of catheter fixation and position of catheter tip should be checked throughout use.
8. When removing dressings at or close to catheter sites, care must be taken to avoid severing catheter.
9. Acetone must not come into contact with the catheter as the material may weaken and may result in leakage or aspiration.
10. Exposure of product componentry to topicals containing alcohol is not recommended.
11. Do not attempt to re-insert a partially or completely withdrawn introducer cannula.
12. Percutaneous puncture of a central vein may be contraindicated in patients with pulmonary hypertension.
13. Use of a syringe smaller than 5 mL to irrigate or de-clot an occluded catheter may cause intraluminal or catheter rupture.
14. Syringes are supplied for blood aspiration only.
15. Physicians should be aware that central venous catheters are intended for use for up to 14 days. The maximum duration of residence of Hydrocath Assure™ catheters in clinical trials is 14 days. It should be noted that the anti-microbial efficacy diminishes with time.
16. Patients with suspected hypersensitivity to nickel should undergo skin test to assess hypersensitivity prior to use of Merit Guidewires in the placement of central venous catheters.

Warnings- Seldinger technique

1. Do not withdraw the guidewire against needle bevel to avoid the possibility of severing of the guidewire.
2. During insertion do not reinsert a partially or completely withdrawn needle into the cannula.
3. Ensure the flexible end of the guidewire is advanced into the vein.
4. Ensure the guidewire moves freely in the needle introducer.
5. A firm grip must be maintained on the guidewire at all times.
6. When using the 'J' wire straightener maintain a firm grip on the plastic sleeve.
7. Ensure the dilator is removed prior to catheter advancement.
8. The moveable suture devices must be used as additional support and not the only means of fixation.
9. Potential for guidewire breakage. Although the incidence of guidewire breakage is extremely uncommon, physicians must be aware of the potential of guidewire breakage if undue force is applied to the wire. If resistance is met when attempting to remove the guidewire after central venous placement, the wire may be kinked within the area of the catheter tip and the vessel. Undue force may cause the wire to break. If resistance is encountered, withdraw the catheter relative to the guidewire (2-3 cm) and attempt to remove the wire. If resistance is still apparent remove the wire and the catheter simultaneously.
10. Physicians should be aware that the guidewire can pick up material from the vein. This may prevent the guidewire from being withdrawn through the needle.
11. Do not force the guidewire. If resistance is met, carefully withdraw the guidewire and re-attempt insertion.
12. If femoral approach is used, place patient in supine position for insertion procedure.
13. HighFlow Central Venous Catheter Kit Seldinger Technique - This device is used for high volume infusions and caution should be applied.


Suggested Procedure Seldinger Technique

1. Prepare the insertion site using the full aseptic technique required for central venous access.
2. The guidewire may be introduced either by using a thin walled steel needle or introducer cannula.
3. Feed the desired flexible end of the guidewire through the introducer into the vein. If a 'J' wire is to be used, the 'J' can be straightened prior to insertion by using the plastic insertion sleeve.
4. Remove the introducer.
5. A vessel dilator may be used to enlarge the cutaneous puncture site. If further enlargement of puncture site is necessary, use a scalpel.
6. When using multi lumen catheters, lumens other than distal should be flushed and attached to desired fluid administration set or, alternatively, heparin locked as standard hospital practice.
7. Pass the distal tip of the catheter over the guidewire. (Sufficient guidewire length must remain exposed at the hub end of the catheter to maintain a firm grip on the guidewire.) Grasping the catheter near the skin, push the catheter into the vein with a slight twisting motion and advance into the final position.
8. Hold the catheter in position, withdraw the guidewire and aspirate with a syringe to ensure correct placement.
9. The distal hub should be connected to the appropriate fluid administration set. If a hub with Floswitch™ is being used this can now be switched off.
10. The catheter may now be secured by suturing eyelets of the junction boot housing to the skin.
11. Apply sterile dressing as appropriate.

Secondary Fixation Device (5Fr -8 Fr catheters)

Note:


This device is found in certain models only.

1. Position the secondary fixation device on catheter. To close the device, press down on each wing. A 'click' sound confirms that the device is secured and cannot be moved easily.
2. To ensure fixation, suture at each eyelet on the wings of the Secondary Fixation Device.
3. To open the device, press down on the 2nd hinge, pull up on each wing.

Secondary Fixation Device (7Fr and 9.5Fr catheters)


Note: This device is found in certain models only.

1. Ensure the Secondary Fixation Device is in place at the point of suture.
2. To close the device, snap the top plate in place.
3. To ensure fixation, suture at each eyelet on the wings of the Secondary Fixation Device.
4. To open the device, hold the lock device down and pull the plate up, using the lever.

Venaguide

Note: This device is found in certain models only. This device is not for sale in the USA.

1. Release guidewire by removing the guidewire cap.
2. Straighten guidewire 'J' by retracting into introducer system with thumb.
3. Insert into hub of introducer needle and advance guidewire into vein. Advance to required depth.

WARNING: To prevent over insertion stop when last depth mark reaches the venepuncture site.

Guiding Syringe

Note: This accessory is found in certain models only. This device is used for introducing guidewires for use with Merit central or peripheral venous catheters.

1. Insert introducer needle attached to Guiding Syringe into vessel and aspirate.
2. Feed the desired flexible end of the guidewire through the hole in rear of Guiding Syringe plunger into the vein. If a 'J' wire is to be used, straighten guidewire 'J' by retracting into introducer system with thumb.
3. Hold the guidewire in place and remove introducer needle and Guiding Syringe.

WARNING: Do not aspirate with guidewire in place or air may enter the syringe.

Caution: Do not reinfuse blood to minimize the risk of blood leakage from rear of syringe.

Floswitch™ Introducer

This device is used for introducing seldinger guidewires or with Merit peripheral and central venous catheters. Using a Floswitch™ introducer cannula, puncture the vein, aspirate and advance the cannula into the vein and remove the needle. At this stage the Floswitch™ may be switched off. The switch must be opened again to facilitate guidewire or catheter passage.

Floswitch™ Luer Lock Attachment

This device may be attached to a luer hub of a central venous catheter.

1. As standard practice the catheter should be secured to the skin.
2. To attach the Floswitch™, insert the Floswitch™ luer into the catheter hub. Tighten the rotating collar ensuring that 'hand tight' connection has been made and that the catheter hub is fully advanced into the Floswitch™ collar (see figure 1.1).
3. The Floswitch™ should be fixed in position by either suturing or taping the wings. The recommended procedure for taping is shown in figure 2.
4. The black markings indicate flow status. When visible the catheter is open and there is free flow. When covered, the catheter is closed.

Figure 1. Floswitch™ Inspection Procedure



Figure 2. Floswitch™ Taping Procedure

Warnings - Floswitch™ Products

1. The Floswitch™ should not be switched off (black marks covered) before the needle has been fully withdrawn. This applies to Floswitch™ introducers only.
2. The Floswitch™ must not be switched 'OFF' (black marks covered) until the guidewire has been fully withdrawn.
3. When the catheter is not in use for infusion or aspiration, the Floswitch™ must be switched off and a suitable luer cap locked into the hub.
4. When using an intermittent injection bung attached to a Floswitch™, only use short needles. Ensure Floswitch™ is in 'ON' position prior to injecting or aspirating. Do not insert the needle more than 8 mm into the Floswitch™.
5. As standard practice, the security of the luer connection must be checked routinely. This is essential when lubricious substances such as Intralipids are being used. This applies to Floswitch™ luer lock attachment only.

Storage Conditions

Store in cool dry place away from direct sunlight.

STERILE and non-pyrogenic in unopened, undamaged package. For single use only. Check integrity of the individual package before use. Dispose of product after use. Do not resterilize. Re-use may lead to infection or other illness/ injury. For reordering information or assistance please contact local representative.

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