21 November 2018

DECLARATION

We, SYSMEX CORPORATION, having its principal place of business at 1-5-1, Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan, hereby declare that the following document is truly prepared, and the content of the document is certified true and correct.

- Instruction for Use for CELLCLEAN AUTO

For and on behalf of
SYSMEX CORPORATION,

Hitoshi Nishiumi
Vice President
Regulatory Affairs
Regulatory Affairs and Quality Assurance

Instruction for Use

**CELLCLEAN AUTO cleaning solution for Sysmex
hematology instruments for in vitro diagnostics**

**manufactured by
SYSMEX CORPORATION, Japan**

Version 1 dated 01. 03.2018

Only for Russian Federation

1. Name of the medical device

CELLCLEAN AUTO cleaning solution for Sysmex hematology instruments for in vitro diagnostics (hereinafter referred to as solution, reagent, CELLCLEAN AUTO).

2. Intended use

For in vitro diagnostic use only

CELLCLEAN AUTO is to be used as a strong alkaline detergent to remove SYSMEX lysing reagent, cellular residuals and blood proteins remaining in the hydraulics of automated hematology analyzers XN series, XN-L-series and automatic system for blood smear preparation and staining SP-10.

3. Field of application

CELLCLEAN AUTO is used in clinical laboratory where diagnostic tests are performed. The device is intended for use only by professionals having appropriate qualification. Prior to use, the operators should carefully read the Instruction for Use.

4. Indication

CELLCLEAN AUTO is necessary for hematology instrument cleaning after performing an analysis.

5. Contraindications

Contradictions to the use of CELLCLEAN AUTO are absent.

6. Possible side effects

Known side effects related to the use of the CELLCLEAN AUTO are absent.

7. Sterility.

CELLCLEAN AUTO cleaning solution is supplied nonsterile.

8. Principles of use

One vial is used for one cleaning sequence.

CELLCLEAN AUTO is a detergent to prevent protein precipitation in the Reaction Chambers, the Waste Chambers, the Detectors, the Piercers, the Aspiration tubes and the HGB Flow Cell.

If this product is used on any other Sysmex instrument, it is unrecognized as detergent and the content solution may leak into the instrument. When using CELLCLEAN AUTO, do not peel off the upper film.

9. Medical devices for use with the cleaning solution

CELLCLEAN AUTO is intended for use only with analyzers Automated Hematology Analyzer XN series (model: XN-10 and XN-20), with Automated Hematology Analyzer XN-L series (model: XN-350, XN-450, XN-550) and Automated Hematology Slide Preparation Unit SP-10 manufactured by SYSMEX CORPORATION, Japan.

Prior to use, please read the analyzer, system and reagent Instructions for Use.

10. Package contents

CELLCLEAN AUTO cleaning solution is supplied as the following set:

- 1) CELLCLEAN AUTO cleaning solutions, 20 vials 4 ml each in package.
- 2) Instruction for Use – 1 pcs.

11. Warnings and precautions

Avoid contact with skin and eyes. In case of skin contact, flush the area with water. In case of contact with eyes, rinse immediately with plenty of water and seek medical advice immediately. Small white dots may be visible in the seal. This is normal and does not affect the quality of the product or integrity of the seal.

Danger

H314 Causes severe skin burns and eye damage.

P264 Wash skin thoroughly after handling.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P303+P361+P353 IF ON SKIN (or hair): Take off immediately all contaminated clothing.

Rinse skin with water/shower.

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes.

Remove contact lenses, if present and easy to do. Continue rinsing.

P310 Immediately call a POISON CENTER/doctor.

P501 Dispose of contents/container in accordance with local and national regulations.

EUH031 Contact with acids liberates toxic gas.

Limitations

Refer to devices Instruction for Use.

CELLCLEAN AUTO should be only used on Automated Hematology Analyzer XN series (model: XN-10, XN-20), Automated Hematology Analyzer XN-L series (model: XN-550, XN-450, XN-350), Automated Hematology Slide Preparation Unit SP-10

When using CELLCLEAN AUTO, do not peel off the upper film. One vial is used for one cleaning sequence.

12. Examination procedure

Scheduled cleaning should be performed as per corresponding analyzer instruction for Use. Scheduled cleaning procedure includes CELLCLEAN AUTO aspiration and hydraulics flushing.

Cleaning procedure is specified in instrument's Instruction for Use.

13. Device characteristics

Device components

Sodium hypochlorite (available chlorine concentration 5.0 %).

Physical properties of the device

Physical state/color	Liquid/yellowish green
----------------------	------------------------

14. Conditions of use, storage and transportation. Reagent stability.

Conditions of use

The reagent is used in clinical laboratories where diagnostic tests are performed within temperature range of 1 °C to 30 °C

Storage conditions

If stored at 1-30 °C the CELLCLEAN AUTO solution in unopened vial is stable within 12 months from the date of manufacture. CELLCLEAN AUTO should be stored in a dark place. Avoid exposing to direct sunlight, or chlorine component will deteriorate and the effectiveness of this detergent will be lost depending upon the time period of exposure.

Storage of the medical device in the original package should be performed in appropriate rooms providing the necessary temperature range (less than 30 °C) with daily recording of the air temperature.

If the storage conditions are not met the medical device should not be used!

If the storage conditions are met the CELLCLEAN AUTO solution in unopened vial is stable within 12 months from the date of manufacture.

Transportation conditions

Transportation is allowed within the shelf life at 1-30 °C in all types of closed vehicles (trucks, planes, trains) according to the transportation regulations valid for particular vehicle. Freezing is not allowed.

Transportation at temperature above 30 °C (but less than 45 °C) is allowed within 42 days.

Other conditions:

- (1) Do not freeze the device
- (2) The devices should not contact the cold package (do not freeze)
- (3) Protect from sunlight/ small particles (dust)

Other conditions:

- (1) Do not freeze the device
- (2) The devices should not contact the cold package (do not freeze)
- (3) Protect from sunlight/ small particles (dust)

Reagent stability

Shelf life after first opening (open vial stability)

Not applicable (this is a single-use device).

Long-term stability

CELLCLEAN AUTO is stable within 12 months at 1-30 °C.

15. Labeling symbols

Symbol	Designation
	Manufacturer
	Manufacturer logo
	<i>in vitro</i> diagnostics medical device
	Warning! Causes corrosion
	Temperature limitation
	Keep away from sunlight
	Consult instructions for use
	Compliance with the essential requirements of the European Directives: 98/79/EC on <i>in vitro</i> diagnostic medical devices
	Use by
	Batch code
	Authorized representative in European Community
	Catalogue number

16. Cleaning, disinfection, sterilization

Not applicable.

17. Contact with human body

There is no direct contact between the product and the human body. When working with the reagent the analyst should wear medical gloves.

18. Disposal

Disposal procedure:

CELLCLEAN AUTO should not be disposed in general waste but should be disposed with infectious medical waste.

Disposal by incineration is recommended. Requirements of applicable local regulations must be considered.

In order to prevent possible contamination of the environment the expired CELLCLEAN AUTO as well as used reagents should be disposed according to the applicable local regulations and hospital regulations.

Appropriate disposal and waste recycling procedures are necessary for natural resources maintenance and protection of human health.

Disposal of waste packaging materials:

CELLCLEAN AUTO package is made to minimize the environment contamination, while ensuring product integrity during transportation and storage. Used packaging materials should be disposed in local places for collection of the packages into special containers situated in municipalities.

Product recycling, disposal of used and expired reagents:

Upon product and used/expired reagents disposal it is necessary to follow hygienic regulations and formal recommendations on recycling and disposal of potentially infectious and chemical materials applied on the territory of the Russian Federation and in hospitals.

Key points of the waste managing strategy include optimization of the management and approaches for the finished products as well as minimization of the environment pollution.

Detailed information is available at manufacturer/authorized representative.

Disposal should be performed according to the local regulatory standards as well as healthcare facility protocols where the medical devices were used for data collection. CELLCLEAN AUTO is disposed on the territory of the Russian Federation according to the SanPiN 2.1.7.2790-10 as potentially infectious medical waste (class B).

Do not dispose with general waste.

19. Requirements for environmental protection

The reagents are developed, manufactured and tested according to the national and international safety and quality standards. The manufacturing process is also performed according to the national laws on environmental protection.

The manufacturer guarantees high level of safety for the environment if the medical device is used according to the intended use and is properly disposed.

20. Warranty

Manufacturer (Sysmex) guarantees that CELLCLEAN AUTO conforms to the specifications indicated on the label. Manufacturer is solely responsible for device performance only if it used according to the instructions and the storage and transportation of the device is performed according to the conditions and regulations specified. Manufacturer guarantees safety of the

medical devices, absence of the unacceptable risk regarding human health as well as environment if the reagent set is used according to the conditions specified by manufacturer. The products should be used up to the expiration date specified on the package. This limited warranty replaces any other explicit or implicit warranties including any implicit warranties of applicability for commercial or any other use. Sysmex takes no responsibility for any accidental, indirect or intentional material losses.

On the topics regarding the quality and performance characteristics of the product, please contact the authorized representative of the company on the territory of Russian Federation:

Sysmex RUS LLC

Registered address: 11 1 Magistral'niy tupik, building 10, floor 2, space V, rooms 1-23, 25-39 Moscow, 123290

Tel: +7(495)7816772

E-mail: support.rus@sysmex-europe.com

Quality claims (complaints):

Complaints mechanisms and corresponding replies are regulated by the law. Only the complaints not related to the product acceptance, manufactured according to the agreement conditions, can be submitted.

On the issues related to the device use one should contact the authorized representative of the manufacturer.

In case of complaint, contact manufacturer and/or the corresponding authorized representative.

21. International standards

CELLCLEAN AUTO cleaning solution is developed and manufactured according to the following standards:

Reference	Title
EN ISO 14971:2012	Medical devices. Application of risk management to medical devices
EN ISO 18113-1:2009	In vitro diagnostic medical devices — Information supplied by the manufacturer (labeling) — Part 1: Terms, definitions and general requirements
EN ISO 18113-2:2009	In vitro diagnostic medical devices -- Information supplied by the manufacturer (labelling) -- Part 2: In vitro diagnostic reagents for professional use
EN ISO 13485: 2003	Medical devices - Quality management systems - Requirements for regulatory purposes
EN ISO 9001: 2008	Quality management systems - Requirements
EN 23640: 2013	In vitro diagnostic medical devices. Evaluation of stability of in vitro diagnostic reagents
EN 980: 2008	Symbols for use in the labeling of the medical devices

22. Contact information

Manufacturer and developer:

SYSMEX CORPORATION

1-5-1 Wakinohama-Kaigandori, Chuo-Ku, Kobe, 651-0073, Japan

Tel: +81 78-265-0500

Fax: +81 78-265-0524

Authorized representative of the manufacturer in Russia:

Sysmex RUS LLC

Registered address: 11 1 Magistral'niy tupik, building 10, floor 2, space V, rooms 1-23, 25-39

Moscow, 123290

Tel: +7(495)781-67-72

Authorized representative in Europe

SYSMEX EUROPE GmbH

Bornbarch 1, 22848 Norderstedt, Germany

Register Number : 737

NOTARIAL CERTIFICATE

This is to certify that Ms. Risa Mochizuki an agent of Mr. Hitoshi Nishiumi, Vice President, Regulatory Affairs, Regulatory Affairs and Quality Assurance of Sysmex Corporation has stated in my very presence that the said Mr. Hitoshi Nishiumi acknowledged himself to have signed to the foregoing document.

Dated this 22nd day of November, 2018

Makoto Okada 

NOTARY Makoto Okada
Kobe District Legal Affairs Bureau
44, Akashimachi, Chuo-ku,
Kobe, Japan

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Registration № 1331

CERTIFICATE

This is to certify that the annexed Notarial Certificate has been executed by Notary, duly authorized and practising in Kobe, Japan, and that the Official seal appearing on the same is genuine.

Date November 22, 2018

Director of the Kobe District Legal Affairs Bureau

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Lighting the way with diagnostics

7 →

GHS Requirements
Danger
H314 Causes severe skin burns and eye damage.
P264 Wash skin thoroughly after handling.
P280 Wear protective gloves/protective clothing/eye protection/face protection.
P303+P361+P353 IF ON SKIN: If on skin: Take off immediately all contaminated clothing. Rinse skin with water/shower.
P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.
P310 Immediately call a POISON CENTER or doctor.
P501 Dispose of contents/container in accordance with local, national regulations.

18 →

Detergent for Fully Automated Hematology Analyzers
Détérgeant pour analyseurs d'hématologie entièrement automatisés
Detergente para analizadores hematológicos totalmente automatizados

19 →

Distributed in the U.S.A. by
Sysmex America, Inc.
577 Apt. 430 P.O. Box 110009, U.S.A.
Tel: +1 224 541-9500

Sysmex is Registered in U.S. Pat. & Tm. Off.

Components
Sodium Hypochlorite
Concentration: 5.0% available chlorine concentration 5.0%.

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
Supplementary Instructions for Using the Medical Product
CELLCLEAN AUTO Cleaning Solution for Sysmex Hematology
Equipment for In Vitro Diagnostics [IVD]

Manufactured by the Company
SYSMEX CORPORATION, Japan

For and on behalf of
SYSMEX CORPORATION

Hitoshi Nishiumi
Vice President
Regulatory Affairs
Regulatory Affairs & Quality Assurance

Sysmex Corporation

www.sysmex.co.jp

1-5-1 Wakinohama-Kaigandori, Chuo-ku, Kobe 651-0073, Japan
Tel 81-78-265-0500 Fax 81-78-265-0524

1. Medical product category code

Category code of the medical product CELLCLEAN AUTO Cleaning Solution for Sysmex Hematology Equipment for in Vitro Diagnostics according to the nomenclature classification of medical products by category:

Category code - 109770.

Name - Washing/cleaning solution for IVD, used in automated/ semi-automated systems.

Description - Concentrated washing or cleaning solution intended for use as a consumable on automatic or semi-automatic devices that prepare, stain and/or analyze clinical laboratory samples.

2. Risk class

Risk class of the medical product CELLCLEAN AUTO Cleaning Solution for Sysmex Hematology Equipment for In Vitro Diagnostics according to the nomenclature classification of medical products that separates them into classes based their risk potential during use - 2a.

3. pH value

The pH value of CELLCLEAN AUTO cleaning solution: ≥ 12.0

4. The medical product is labeled as follows:

4.1 The vial labeling contains (see Figure 1):

- manufacturer's logo (see № 1 in Figure 1)
- name of manufacturer, address, telephone (see № 2 in Figure 1)
- special labeling (EC) EUH031 (see № 3 in Figure 1)
- abbreviated name of the product (see № 4 in Figure 1)
- alphanumeric packaging code (see № 5 in Figure 1)
- vial volume (see № 6 in Figure 1)

- lot number (see № 7 in Figure 1)
- expiration date (see № 8 in Figure 1)
- symbol “Only for In Vitro Diagnostics” (see № 9 in Figure 1)
- warning “Danger: causes severe skin burns and eye damage” (see № 10 in Figure 1)
- codification of hazard characteristics P264, P280, P303+P361+P353, P305+P351+P338, P310, P501, EUH031 (see № 11 in Figure 1)
- symbol “Danger: causes corrosion and burns” (№ 12 in Figure 1)
- symbol “Temperature range” (see № 13 in Figure 1)
- symbol “Consult instructions for use” (see № 14 in Figure 1)
- symbol “Do not expose to sunlight” (see № 15 in Figure 1)
- symbol of compliance with EU directives and harmonized standards (see № 16 in Figure 1)
- bar code (see № 17 in Figure 1)
- components: contain sodium hypochlorite (see № 18 in Figure 1)


Figure 1. Layout of vial label

4.2 Labeling on the outer box displays the following (see Figures 2-7):

- manufacturer's logo (see № 1 in Figures 2, 5, 6)
- name of manufacturer, address, (see № 2 in Figure 2)
- name of company representative in EU, address (see № 3 in Figure 2)
- abbreviated name of the product (see № 4 in Figures 2, 5, 6)
- alphanumeric packaging code (see № 5 in Figures 4,7)
- symbol "Consult instructions for use" (see № 6 in Figure 2)
- codification of hazard characteristics H314, P264, P280, P303+P361+P353, P305+P351+P338, P310, P501, EUH031c (see № 7 in Figures 2,3)
- symbol "Danger. Causes corrosion and burns" (see № 8 in Figure 2)
- components of the product (see № 9 in Figure 3)
- catalog number (see № 10 in Figure 5)
- symbols "Only for In Vitro Diagnostics" (see № 11 in Figure 5)
- symbol "Do not expose to sunlight" (see № 12 in Figure 5)
- symbol of compliance with EU directives and harmonized standards (see № 13 in Figure 5)
- bar code (see № 14 in Figure 5,7)
- symbol "Temperature range" (see № 15 in Figure 5)
- expiration date (see № 16 in Figure 7)
- lot number (see № 17 in Figure 7)
- description of the product (see № 18 in Figure 3)
- information about USA representative (see № 19 in Figure 3)
- warnings and precautionary measures (see № 20 in Figure 3)
- package contents (see № 21 in Figure 5)
- information in Japanese: description of the product, warnings and precautionary measures, components of the product, manufacturer's name, address, telephone (see № 22 in Figure 4)
- information in Chinese: description of the product, warnings and precautionary measures, components of the product, manufacturer's name, address, telephone (see № 23 in Figure 4)

Figure 2. Layout of the manufacturer's labeling on the top of the outer box, in which the medical product is packed

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Figure 3. Layout of the manufacturer's labeling on the side of the outer box, in which the medical product is packed


Figure 4. Layout of the manufacturer's labeling on the side of the outer box, in which the medical product is packed

Figure 5. Layout of the manufacturer's labeling on the side of the outer box, in which the medical product is packed


Figure 6. Layout of the manufacturer's labeling on the side of the outer box, in which the medical product is packed


Figure 7. Layout of the manufacturer's labeling on the side of the outer box, in which the medical product is packed

On every outer box, there will also be handling instructions that apply to the territory of the Russian Federation.

The layout of an additional label is given below:

CELLCLEAN AUTO Cleaning Solution for Sysmex Hematology Equipment for In Vitro Diagnostics
Contents: 1) Cleansing solution CELLCLEAN AUTO, 20 vials of 4 ml. each/pkg. 2) Instructions for use.
Registration certificate: Roszdravnadzor № XXXX/XXX
Series number, expiration date, storage requirements, contents and graphic symbols are indicated on the original packaging
Manufacturer Sysmex Corporation Address: 1-5-1 Wakinohama-Kaigandori, Chuo-Ku Kobe, Japan 651-0073 Telephone: +81 78-265-0500 Fax: +81 78-265-0525 Country of manufacture: Sysmex International Reagents Co., Ltd, Ono Factory, 17 Takumidai, Ono, Hyogo 675-1322, Japan
Manufacturer's authorized representative: Sysmex RUS, LLC, Russia Address: 123290 Moscow, 1 Magistral Cul-de-sac, 11 Building 10, 2nd floor, Suite V, Rooms 1-23, 25-39 Telephone: +7(495)781-67-72.
Danger: can seriously burn the skin and damage the eyes Components: sodium hypochlorite

4.3 Labeling of a multi-unit package (overpack) indicates:

- the name of the medical product
- the manufacturer's name and address
- the product's expiration date (year, month inclusively);

- the date when product was manufactured (year, month);
- symbol "Danger. Causes corrosion and burns"
- symbol "Do not expose to sunlight"
- conditions for storing and transporting the product
- symbol of compliance with EU directives and harmonized standards
- symbol "Temperature range"
- catalog number
- lot
- weight
- Number and issuance date of registration certificate

Labeling of a multi-unit package (overpack):

CELLCLEAN AUTO Cleaning Solution for Sysmex Hematology Equipment for In Vitro Diagnostics	
Date of manufacture: XX./XXXX	 use until XX.XX./XXXX
REF XX-XXX-XXX	LOT XXXXXXXX
  	 Temperature range from XX to XX
Weight: XXXX	
Registration certificate: Roszdravnadzor № XXXX/XXX	
Manufacturer: Sysmex Corporation Address: 1-5-1 Wakinohama-Kaigandori, Chuo-Ku Kobe, Japan 651-0073 Telephone +81 78-265-0500 Fax +81 78-265-0525 Country of manufacture: Sysmex International Reagents Co., Ltd, Ono Factory, 17 Takumidai, Ono, Hyogo 675-1322, Japan	

1


Register Number : 618

NOTARIAL CERTIFICATE

This is to certify that Kaho Uetani an agent of Hitoshi Nishiumi, Vice President, Regulatory Affairs, Regulatory Affairs & Quality Assurance of Sysmex Corporation has stated in my very presence that the said Hitoshi Nishiumi acknowledged himself to have signed to the foregoing document.

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