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
Abbott

Abbott Ireland Diagnostics
Division,
Lisnamuck,
Longford
Co. Longford
Ireland

Tel. +353 43 3331000
Fax +353 43 3331001

ABBOTT DIAGNOSTICS DIVISION LONGFORD

I, Nora Fitzgerald, Regulatory Affairs Specialist, Abbott Ireland Diagnostics Division, Longford, Ireland, hereby certify that the attached documents that keep the number of pages indicated below are the Instruction for Use. These documents are in a full accordance with actual manuals of manufacturer.

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Signed: Nora Fitzgerald

Date: 14/04/2019

Nora Fitzgerald
Regulatory Affairs Specialist
Abbott Ireland Diagnostics Division
Lisnamuck, Longford
Co Longford
Ireland

ABBOTT DIAGNOSTICS

SLIGO CHAMBER OF COMMERCE (INC)
CERTIFIED CORRECT IN ALL GOOD FAITH

Mary Hart 08/04/2019

Authorized Official

Instru
canni
instru

■ N

Alinity

■ IF

The /
micro
deter
plasr

■ S

3,5,3'
weigh
circul
boun
preall
bindir
TBG.
these
in sol
physi
It has
role i
meat
pane
cond
of se
hyper
disen
T₃ m
Reco
asso-
(T₄),
the li
defic
desp
serui
Horn
with
iodin
eleva
indiv
T₃ le
cont
prec
occi
ther.
salic
abni
the
abni
Ser.
or s
wha
norr
assa

Alinity i

Total T₃ Reagent Kit

Read Highlighted Changes: Revised February 2018.

en

Total T₃

07P94

G71277R04

B7P940

07P9420

07P9430

Instructions must be carefully followed. Reliability of assay results cannot be guaranteed if there are any deviations from these instructions.

■ NAME

Alinity i Total T₃ Reagent Kit

■ INTENDED USE

The Alinity i Total T₃ (TT3) assay is a chemiluminescent microparticle immunoassay (CMIA) used for the quantitative determination of total triiodothyronine (Total T₃) in human serum and plasma on the Alinity i analyzer.

■ SUMMARY AND EXPLANATION OF THE TEST

3,5,3' Triiodothyronine (T₃) is a thyroid hormone with a molecular weight of 651 daltons¹ and a half-life in serum of 1.5 days.² T₃ circulates in the blood as an equilibrium mixture of free and protein bound hormone.³ T₃ is bound to thyroxine binding globulin (TBG), prealbumin, and albumin. The actual distribution of T₃ among these binding proteins is controversial as estimates range from 38-80% for TBG, 9-27% for prealbumin, and 11-35% for albumin.⁴ The binding of these proteins is such that only 0.2-0.4% of the total T₃ is present in solution as unbound or free T₃.⁵ This free fraction represents the physiologically active thyroid hormone.³

It has become apparent in recent years that T₃ plays an important role in the maintenance of the euthyroid state. Serum T₃ measurements can be a valuable component of a thyroid screening panel in diagnosing certain disorders of thyroid function as well as conditions caused by iodine deficiency. Clinically, measurements of serum T₃ concentration are especially valuable in diagnosing hyperthyroidism and in following the course of therapy for this disorder.^{2, 6, 7} Under conditions of strong thyroid stimulation, the T₃ measurement provides a good estimation of thyroid reserve.² Recognition of a thyroid dysfunction called T₃-thyrotoxicosis, associated with an increased serum T₃ level but normal thyroxine (T₄), free T₄, and *in vitro* Uptake results have further highlighted the importance of serum T₃ measurements.^{2, 8-11} Dietary iodine deficiency results in inadequate production of thyroid hormones despite the presence of normal thyroid tissue. In these cases, the serum T₄ concentration is often low while the Thyroid Stimulating Hormone (TSH) concentration is elevated. Elevated TSH associated with low T₄ is normally indicative of hypothyroidism. However, in iodine deficiency, these results together with normal or slightly elevated serum T₃ are indicative of euthyroid status in most individuals.¹²

T₃ levels are also affected by conditions which affect TBG concentration.¹³⁻¹⁵ Slightly elevated T₃ levels may occur in pregnancy or during estrogen therapy, while depressed levels may occur during severe illness, malnutrition, in renal failure and during therapy with anti-thyroid drugs, propranolol and propylthiouracil and salicylates.^{2, 16, 17} In patients with severe or chronic illnesses, many abnormalities of thyroid hormone balance occur. T₄ production and the extent of serum thyroid hormone binding may be independently abnormal, resulting in a low, normal or high free T₄ estimate. Serum T₃ concentrations are often low, TSH levels may be normal or slightly elevated. Total T₃ measurements may be valuable when hyperthyroidism is suspected and the free T₄ estimate is normal.¹³ The Alinity i Total T₃ assay is to be used as an aid in the assessment of thyroid status.

■ BIOLOGICAL PRINCIPLES OF THE PROCEDURE

This assay is a two-step immunoassay for the quantitative determination of total triiodothyronine (Total T₃) in human serum and plasma using chemiluminescent microparticle immunoassay (CMIA) technology.

Sample and anti-T₃ coated paramagnetic microparticles are combined and incubated. The T₃ present in the sample binds to the anti-T₃ coated microparticles. The mixture is washed. T₃ acridinium-labeled conjugate is added to create a reaction mixture and incubated. Following a wash cycle, Pre-Trigger and Trigger Solutions are added.

The resulting chemiluminescent reaction is measured as relative light units (RLUs). There is an inverse relationship between the amount of Total T₃ in the sample and the RLUs detected by the system optics.

For additional information on system and assay technology, refer to the Alinity ci-series Operations Manual, Section 3.

■ REAGENTS

Kit Contents

Alinity i Total T₃ Reagent Kit 07P94

NOTE: Some kit sizes are not available in all countries. Please contact your local distributor.

Volumes (mL) listed in the table below indicate the volume per cartridge.

	07P9420	07P9430
Tests per cartridge	100	600
Number of cartridges per kit	2	2
Tests per kit	200	1200
MICROPARTICLES	6.6 mL	32.1 mL
CONJUGATE	4.2 mL	16.3 mL

MICROPARTICLES Anti-T₃ (sheep) coated microparticles in MES buffer with sheep IgG stabilizers. Minimum concentration: 0.05% solids. Preservative: ProClin 300.

CONJUGATE T₃ acridinium-labeled conjugate in citrate buffer with NaCl and Triton X-100 stabilizers. Minimum concentration: 0.33 ng/mL. Preservative: ProClin 300.

Warnings and Precautions

- **IVD**
- For *In Vitro* Diagnostic Use
-

Safety Precautions

CAUTION: This product requires the handling of human specimens. It is recommended that all human-sourced materials be considered potentially infectious and handled in accordance with the OSHA Standard on Bloodborne Pathogens. Biosafety Level 2 or other appropriate biosafety practices should be used for materials that contain or are suspected of containing infectious agents.¹⁸⁻²¹

Abbott

The following warnings and precautions apply to: **MICROPARTICLES / CONJUGATE**

WARNING	Contains methylisothiazolones.
H317	May cause an allergic skin reaction.
Prevention	
P261	Avoid breathing mist / vapors / spray.
P272	Contaminated work clothing should not be allowed out of the workplace.
P280	Wear protective gloves / protective clothing / eye protection.
Response	
P302+P352	IF ON SKIN: Wash with plenty of water.
P333+P313	If skin irritation or rash occurs: Get medical advice / attention.
P362+P364	Take off contaminated clothing and wash it before reuse.
Disposal	
P501	Dispose of contents / container in accordance with local regulations.

Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

For a detailed discussion of safety precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 8.

Reagent Handling

- Upon receipt, gently invert the unopened reagent kit by rotating it over and back for a full 180 degrees, 5 times with green label stripe facing up and then 5 times with green label stripe facing down. This ensures that liquid covers all sides of the bottles within the cartridges. During reagent shipment, microparticles can settle on the reagent septum.
 - Place a check in the square on the reagent kit to indicate to others that the inversions have been completed.
- After mixing, place reagent cartridges in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- If a reagent cartridge is dropped, place in an upright position for 1 hour before use to allow bubbles that may have formed to dissipate.
- Reagents are susceptible to the formation of foam and bubbles. Bubbles may interfere with the detection of the reagent level in the cartridge and cause insufficient reagent aspiration that may adversely affect results.

For a detailed discussion of reagent handling precautions during system operation, refer to the Alinity ci-series Operations Manual, Section 7.

Reagent Storage

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Unopened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright, gently invert the cartridge 10 times and place in an upright position for 1 hour before use.
Onboard	System Temperature	30 days	

	Storage Temperature	Maximum Storage Time	Additional Storage Instructions
Opened	2 to 8°C	Until expiration date	Store in upright position. If cartridge does not remain upright during storage, discard the cartridge. Do not reuse original reagent caps or replacement caps due to the risk of contamination and the potential to compromise reagent performance.

Reagents may be stored on or off the system. If removed from the system, store reagents with new replacement caps in an upright position at 2 to 8°C. For reagents stored off the system, it is recommended that they be stored in their original trays or boxes to ensure they remain upright.

For information on unloading reagents, refer to the Alinity ci-series Operations Manual, Section 5.

Indications of Reagent Deterioration

Deterioration of the reagents may be indicated when a calibration error occurs or a control value is out of the specified range. Associated test results are invalid, and samples must be retested. Assay recalibration may be necessary.

For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.

INSTRUMENT PROCEDURE

The Alinity i Total T₃ assay file must be installed on the Alinity i analyzer prior to performing the assay.

For detailed information on assay file installation and viewing and editing assay parameters, refer to the Alinity ci-series Operations Manual, Section 2.

For information on printing assay parameters, refer to the Alinity ci-series Operations Manual, Section 5.

For a detailed description of system procedures, refer to the Alinity ci-series Operations Manual.

Alternate Result Units

Edit assay parameter "Result Units" to select an alternate unit.

Conversion formula:

(Concentration in Default result unit) x (Conversion factor) = (Concentration in Alternate result unit)

Default Result Unit	Conversion Factor	Alternate Result Unit
ng/mL	1.536	nmol/L
	100.0	ng/dL

SPECIMEN COLLECTION AND PREPARATION FOR ANALYSIS

Specimen Types

The specimen types listed below were verified for use with this assay.

Other specimen types and collection tube types have not been verified with this assay.

Specimen Types	Collection Tubes
Serum	Serum Serum separator
Plasma	Sodium heparin Lithium heparin Potassium EDTA

- When serial specimens are being evaluated, the same type of specimen should be used throughout the study.

- Performance has not been established for the use of neonatal specimens.
- The instrument does not provide the capability to verify specimen types. It is the responsibility of the operator to verify that the correct specimen types are used in the assay.

Specimen Conditions

- Do not use:
 - heat-inactivated specimens
- For accurate results, serum and plasma specimens should be free of fibrin, red blood cells, and other particulate matter. Serum specimens from patients receiving anticoagulant or thrombolytic therapy may contain fibrin due to incomplete clot formation.
- Ensure that complete clot formation in serum specimens has taken place prior to centrifugation. Some specimens, especially those from patients receiving anticoagulant or thrombolytic therapy may exhibit increased clotting time. If the specimen is centrifuged before a complete clot forms, the presence of fibrin may cause erroneous results.
- To prevent cross contamination, use of disposable pipettes or pipette tips is recommended.

Preparation for Analysis

- Follow the tube manufacturer's processing instructions for collection tubes. Gravity separation is not sufficient for specimen preparation.
- Specimens should be free of bubbles. Remove bubbles with an applicator stick before analysis. Use a new applicator stick for each specimen to prevent cross-contamination.

To ensure consistency in results, recentrifuge specimens prior to testing if

- they contain fibrin, red blood cells, or other particulate matter.
- NOTE: If fibrin, red blood cells, or other particulate matter are observed, mix by low speed vortex or by inverting 10 times prior to recentrifugation.

Prepare frozen specimens as follows:

- Frozen specimens must be completely thawed before mixing.
- Mix thawed specimens thoroughly by low speed vortex or by inverting 10 times.
- Visually inspect the specimens. If layering or stratification is observed, mix until specimens are visibly homogeneous.
- If specimens are not mixed thoroughly, inconsistent results may be obtained.
- Recentrifuge specimens.

Recentrifugation of Specimens

- Transfer specimens to a centrifuge tube and centrifuge.
- Transfer clarified specimen to a sample cup or secondary tube for testing. For centrifuged specimens with a lipid layer, transfer only the clarified specimen and not the lipemic material.

Specimen Storage

Specimen Type	Temperature	Maximum Storage Time	Special Instructions
Serum/ Plasma	2 to 8°C	6 days	If testing will be delayed more than 6 days, specimens should be frozen at -10°C or colder.

If testing will be delayed more than 24 hours, remove serum or plasma from the clot, serum separator or red blood cells.

Specimens stored frozen at -10°C or colder for 6 days showed no performance difference.

Avoid multiple freeze/thaw cycles.

Specimen Shipping

Package and label specimens in compliance with applicable state, federal, and international regulations covering the transport of clinical specimens and infectious substances.

PROCEDURE

Materials Provided

07P94 Alinity i Total T₃ Reagent Kit

Materials Required but not Provided

- Alinity i Total T₃ assay file
- 07P9401 Alinity i Total T₃ Calibrators
- Commercially available controls containing Total T₃
- 07P9440 Alinity i Total T₃ Manual Diluent
- Alinity Trigger Solution
- Alinity Pre-Trigger Solution
- Alinity i-series Concentrated Wash Buffer

For information on materials required for operation of the instrument, refer to the Alinity ci-series Operations Manual, Section 1.

For information on materials required for maintenance procedures, refer to the Alinity ci-series Operations Manual, Section 9.

Assay Procedure

For a detailed description of how to run an assay, refer to the Alinity ci-series Operations Manual, Section 5.

- If using primary or aliquot tubes, refer to the Alinity ci-series Operations Manual, Section 4 to ensure sufficient specimen is present.
- To minimize the effects of evaporation, verify adequate sample cup volume is present prior to running the test.
- Maximum number of replicates sampled from the same sample cup: 10
 - Priority:
 - Sample volume for first test: 70 µL
 - Sample volume for each additional test from same sample cup: 20 µL
 - ≤ 3 hours on the reagent and sample manager:
 - Sample volume for first test: 150 µL
 - Sample volume for each additional test from same sample cup: 20 µL
 - > 3 hours on the reagent and sample manager:
 - Replace with a fresh aliquot of sample.
- Refer to the Alinity i Total T₃ calibrator package insert and the commercially available control package insert for preparation and usage.
- For general operating procedures, refer to the Alinity ci-series Operations Manual, Section 5.
- For optimal performance, it is important to perform routine maintenance as described in the Alinity ci-series Operations Manual, Section 9. Perform maintenance more frequently when required by laboratory procedures.

Sample Dilution Procedures

Samples with a Total T₃ value exceeding 6.00 ng/mL (9.22 nmol/L) are flagged with the code "> 6.00 ng/mL" (>9.22 nmol/L) and may be diluted with the Manual Dilution Procedure.

Manual Dilution Procedure

Suggested dilution: 1:2

It is recommended that dilutions not exceed 1:2.

Add a minimum of 75 µL of the sample to 75 µL of Alinity i Total T₃ Manual Diluent.

To avoid contamination of Alinity i Total T₃ Manual Diluent, dispense several drops of the diluent into a clean test tube prior to pipetting.

The operator must enter the dilution factor (2) in the Specimen or Control tab of the Create Order screen. The system will use this dilution factor to automatically calculate the concentration of the sample and report the result. The dilution should be performed so that the reported result reads greater than 1.0 ng/mL (1.5 nmol/L).

If the operator does not enter the dilution factor, the result must be manually multiplied by the appropriate dilution factor before reporting the result. If a diluted sample result is ≤ 0.5 ng/mL (0.8 nmol/L), do not report the result. Rerun using an appropriate dilution.

For detailed information on ordering dilutions, refer to the Alinity ci-series Operations Manual, Section 5.

Calibration

For instructions on performing a calibration, refer to the Alinity ci-series Operations Manual, Section 5.

Each assay control must be tested to evaluate the assay calibration. Once a calibration is accepted and stored, all subsequent samples may be tested without further calibration unless:

- A reagent kit with a new lot number is used.
- Daily quality control results are outside of statistically-based quality control limits used to monitor and control system performance, as described in the Quality Control Procedures section of this package insert.
 - If statistically-based quality control limits are not available, then the calibration should not exceed a 30-day limit for recalibration frequency.

This assay may require recalibration after maintenance to critical parts or subsystems or after service procedures have been performed.

Quality Control Procedures

The recommended control requirement for the Alinity i Total T₃ assay is that a single sample of each control level be tested once every 24 hours each day of use.

Additional controls may be tested in accordance with local, state, and/or federal regulations or accreditation requirements and your laboratory's quality control policy.

To establish statistically-based control limits, each laboratory should establish its own concentration target and ranges for new control lots at each clinically relevant control level. This can be accomplished by assaying a minimum of 20 replicates over several (3-5) days and using the reported results to establish the expected average (target) and variability about this average (range) for the laboratory. Sources of variation that should be included in this study in order to be representative of future system performance include:

- Multiple stored calibrations
- Multiple reagent lots
- Multiple calibrator lots
- Multiple processing modules (if applicable)
- Data points collected at different times of the day

Refer to published guidelines for information or general control recommendation, for example Clinical and Laboratory Standards Institute (CLSI) Document C24-A3 or other published guidelines, for general quality control recommendations.²²

- If more frequent control monitoring is required, follow the established quality control procedures for your laboratory.
- If quality control results do not meet the acceptance criteria defined by your laboratory, sample results may be suspect. Follow the established quality control procedures for your laboratory. Recalibration may be necessary. For troubleshooting information, refer to the Alinity ci-series Operations Manual, Section 10.
- Review quality control results and acceptance criteria following a change of reagent or calibrator lot.

Commercial controls should be used according to the guidelines and recommendations of the control manufacturer. Concentration ranges provided in the control package insert should be used only for guidance.

For any control material in use, the laboratory should ensure that the matrix of the control material is suitable for use in the assay per the assay package insert.

Quality Control Guidance

Refer to "Basic QC Practices" by James O Westgard, Ph.D. for guidance on laboratory quality control practices.²³

Verification of Assay Claims

For protocols to verify package insert claims, refer to Verification of Assay Claims in the Alinity ci-series Operations Manual.

RESULTS

Calculation

The Alinity i Total T₃ assay utilizes a 4 Parameter Logistic Curve fit data reduction method (4PLC, Y-weighted) to generate a calibration and results.

For information on alternate result units, refer to the INSTRUMENT PROCEDURE, Alternate Result Units section of this package insert.

Flags

Some results may contain information in the Flags field. For a description of the flags that may appear in this field, refer to the Alinity ci-series Operations Manual, Section 5.

Measuring Interval

Measuring interval is defined as the range of values in ng/mL (nmol/L) which meets the limits of acceptable performance for linearity, imprecision, and bias.

The measuring interval of the Alinity i Total T₃ assay is 0.40 to 6.00 ng/mL (0.61 to 9.22 nmol/L).

LIMITATIONS OF THE PROCEDURE

- For diagnostic purposes, results should be used in conjunction with other data; e.g., symptoms, results of other thyroid tests, clinical impressions, etc.
- If the Total T₃ results are inconsistent with clinical evidence, additional testing is suggested to confirm the result.

EXPECTED VALUES

This study was performed on the ARCHITECT i System.

Representative performance data are provided in this section.

Results obtained in individual laboratories may vary.

It is recommended that each laboratory determine its own reference range based upon its particular locale and population characteristics.

A normal range of 0.64 - 1.52 ng/mL (Central 95% interval) was obtained by testing serum specimens from 260 individuals determined as normal by ARCHITECT Anti-Tg, Anti-TPO and TSH assays.

SPECIFIC PERFORMANCE CHARACTERISTICS

Representative performance data are provided in this section.

Results obtained in individual laboratories may vary.

The Alinity i analyzer and the ARCHITECT i System utilize the same reagents and sample/reagent ratios.

Unless otherwise specified, all studies were performed on the Alinity i analyzer.

Precision

Within-Laboratory Precision

A study was performed based on guidance from CLSI EP05-A2. Testing was conducted using 1 lot of the Alinity i Total T₃ Reagent Kit, 1 lot of the Alinity i Total T₃ Calibrator Kit, and 1 lot of commercially available controls and 1 instrument. Three serum panels were assayed in a minimum of 2 replicates at 2 separate times per day on 20 different days.²⁴

Sample	n	Mean (ng/mL)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
			SD	%CV	SD	%CV
Panel 1	120	0.36	0.018	5.0	0.018	5.1
Panel 2	120	1.01	0.023	2.2	0.023	2.3
Panel 3	120	5.67	0.500	8.8	0.550	9.7

Sample	n	Mean (nmol/L)	Within-Run (Repeatability)		Within-Laboratory (Total) ^a	
			SD	%CV	SD	%CV
Panel 1	120	0.55	0.027	4.9	0.028	5.0
Panel 2	120	1.56	0.035	2.2	0.035	2.3
Panel 3	120	8.72	0.768	8.8	0.845	9.7

^a Includes within-run, between-run, and between-day variability.

Lower Limits of Measurement

A study was performed based on guidance from CLSI EP17-A2. Testing was conducted using 3 lots of the Alinity i Total T₃ Reagent Kit on each of 2 instruments over a minimum of 3 days. The Limit of Blank (LoB), Limit of Detection (LoD), and Limit of Quantitation (LoQ) values are summarized below. These representative data support the lower limit of the measuring interval.²⁵

	ng/mL	nmol/L
LoB ^a	0.00	0.00
LoD ^b	0.05	0.08
LoQ ^{c,d}	0.30	0.50

^a The LoB represents the 95th percentile from $n \geq 60$ replicates of zero-analyte samples.

^b The LoD represents the lowest concentration at which the analyte can be detected with 95% probability based on $n \geq 60$ replicates of low-analyte level samples.

^c The LoQ is defined as the lowest concentration at which a maximum allowable precision of 20 %CV was met.

^d This value represents the observed LoQ on the ARCHITECT System. The LoQ observed on the Alinity i analyzer supports this LoQ.

Linearity

A study was performed based on guidance from CLSI EP06-A.²⁶ This assay is linear across the measuring interval of 0.40 to 6.00 ng/mL (0.61 to 9.22 nmol/L).

Analytical Specificity

This study was performed on the ARCHITECT i System. Analytical specificity was evaluated with the ARCHITECT Total T₃ assay and determined to have a mean analytical specificity of $\leq 0.1\%$ cross reactivity with thyroxine (T₄) at a concentration of 1100 ng/mL.

Interference

This study was performed on the ARCHITECT i System.

Potentially Interfering Endogenous Substances

The testing demonstrated $\leq 10\%$ mean interference at the levels indicated below.

Potentially Interfering Substance	Interferent Level
Hemoglobin	≤ 500 mg/dL
Bilirubin	≤ 20 mg/dL
Triglycerides	≤ 2000 mg/dL
Protein	≤ 12 g/dL

Method Comparison

A study was performed based on guidance from CLSI EP09-A3 using the Passing-Bablok regression method.²⁷

Assay	Sample Type	Units	n	Correlation Coefficient	Intercept	Slope	Concentration Range
Alinity i	Serum	ng/mL	123	1.00	0.06	0.90	0.46-5.44
Total T ₃ vs ARCHITECT	Serum	nmol/L	123	1.00	0.09	0.90	0.70-8.36
Total T ₃							

BIBLIOGRAPHY

1. Budavari S, editor. *Merck Index* (11th Ed.). Rahway, NJ: Merck and Co., Inc.; 1989:868.
2. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972;21:1073-1092.

3. Ekins RP, editor. *Methods for the Measurement of Free Thyroid Hormones*. Amsterdam: Excerpta Medica Foundation; 1979:72-92.
4. Robbins J, Rail JE. The Iodine-Containing Hormones. In: *Hormones in Blood* (3rd Ed.). London: Academic Press; 1979:1:632-667.
5. DeGroot LJ, Larsen PR, Refetoff S, Stanbury JB. Transport of Thyroid Hormone and Cell Uptake. In: *The Thyroid and Its Diseases*. New York: Wiley and Sons; 1984:62-66.
6. Wahner HW, Gorman CA. Interpretation of Serum Triiodothyronine Levels Measured by the Sterling Technic. *N Engl J Med* 1971;284:225-230.
7. Marsden P, McKerron CG. Serum Triiodothyronine Concentration in the Diagnosis of Hyperthyroidism. *Clin Endocrinol* 1975;4:183-189.
8. Ivy HK, Washner HW, Gorman CA. Triiodothyronine (T₃) Toxicosis: Its Role in Graves' Disease. *Arch Intern Med* 1971;128:529-534.
9. Hollander CS, Mitsuma T, Nihel N, Shenkman L, Burday SZ, Blum M. Clinical and Laboratory Observations in Cases of Triiodothyronine Toxicosis Confirmed by Radioimmunoassay. *Lancet* 1972;1:609-611.
10. Sterling K, Refetoff S, Selenkow HA. T₃ Thyrotoxicosis: Thyrotoxicosis Due to Elevated Serum Triiodothyronine Levels. *JAMA* 1970;213: 571-575.
11. Hollander CS, Mitsuma T, Shenkman L, Stevenson C, Pinada G, Silva E. T₃ Toxicosis in an Iodine-Deficient Area. *Lancet* 1972;2:1276-1278.
12. Ermans AM. Disorders of Iodine Deficiency. In: Ingbar SH, Braverman LE, editors. *The Thyroid* (5th Ed.). Philadelphia: JB Lippincott Co.; 1986:705-721.
13. Kaplan MM, Larsen PR, Crantz FR, Dzau VJ, Rossing TH, Haddow JE. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. *Am J Med* 1982;72:9-16.
14. Bermudez F, Surks MI, Oppenheimer JH. High Incidence of Decreased Serum Triiodothyronine Concentration in Patients with Nonthyroid Disease. *J Clin Endocrinol Metab*. 1975;41:27-40.
15. Oppenheimer JH. Thyroid Function Tests in Nonthyroidal Disease. *J Chronic Dis* 1982;35:697-701.
16. Abuid J, Larsen PR. Triiodothyronine and Thyroxine in Hyperthyroidism: Comparison of the Acute Changes During Therapy with Antithyroid Agents. *J Clin Invest* 1974;54:201-208.
17. Felig P, Baxter JD, Broadus AE, Frohman LA, editors. *Endocrinology and Metabolism* (2nd Ed.). New York: McGraw-Hill Book Co.; 1987:408-416.
18. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens
19. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
20. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
21. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
22. Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures: Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
23. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI: Westgard Quality Corporation; 2010.
24. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
25. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
26. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
27. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

Note for number formatting:

- A space is used as thousands separator (example: 10 000 specimens).
- A period is used to separate the integer part from the fractional part of a number written in decimal form (example: 3.12%).

■ Key to Symbols

ISO 15223 Symbols	
	Consult instructions for use
	Manufacturer
	Sufficient for
	Temperature limitation
	Use by/Expiration date
	In Vitro Diagnostic Medical Device
	Lot Number
	List Number
	Serial number

Other Symbols	
	Conjugate
	Distributed in the USA by
	Information needed for United States of America only
	Inversions Performed
	Microparticles
	Product of Ireland
	For use by or on the order of a physician only (applicable to USA classification only).

Alinity and ARCHITECT are trademarks of Abbott Laboratories in various jurisdictions. All other trademarks are property of their respective owners.

 Abbott Ireland
Diagnostics Division
Lisnamuck, Longford
Co. Longford
Ireland
+353-43-3331000

Abbott Laboratories
Abbott Park, IL 60064 USA

Customer Service: Contact your local representative
or find country-specific contact information on
www.abbottdiagnostics.com

Revised February 2018.
©2017, 2018 Abbott Laboratories

ADC
to the
SHIP
Recon
contain
mainle
MAB
for use
vehicle
carrier
WAF
Manuf
stored
month
in mar
MAII
Not re
REP
Not re
DES
Destru
SanPir
Requir
Transr
with lo
LIST
98/79/1
80/181
CLAS
ISO 90
EN ISO
1, EN
POTI
The Al
laboral
person
INFC
COM
Does r
STEI
Medic
INFC
APP
Medic
MAT
Materi
Data
repre
INDI
Not a
STC
Stora
room
tempe

ADDENDUM

to the Instruction for Use Alinity i Total T3 Reagent Kit (for Russia)

SHIPPING

Recommended shipping conditions for device is 2-8 °C in single-use thermo containers containing coolant elements, or in reusable thermo containers with automatically maintained temperature or in autorefriigerators using thermal indicators. Device is stable until the expiration date indicated on the package. Device is no longer suitable for use if shipping temperature limits have been violated. Ship in an enclosed carrier vehicle of any kind in accordance with shipping rules, applicable to the respective carrier.

WARRANTY

Manufacturer guarantees that device remains stable until the expiration date when stored and handled as directed. The maximum shelf life of the medical device is 11 months from date of manufacturing. Device must be discarded if stored with deviations in manufacturer's requirements.

MAINTENANCE

Not required.

REPAIRS

Not required.

DESTRUCTION AND DISPOSAL OF A MEDICAL DEVICE

Destruct and dispose of a medical device in accordance with the requirements of SanPin 2.1.7.2780-10 dated December 12, 2010 "Sanitary-epidemiological Requirements for the Management of Medical Wastes".

Transmit and dispose of a product and its components and packaging in accordance with local regulations.

LIST OF APPLICABLE STANDARDS

98/79/EC, 1272/2008/EC (CLP) as amended, 1907/2006/EC (REACH) as amended, 80/181/EEC, 2009/3/EC, Regulation (EC) 1069/2009, Regulation (EC) 142/2011, CLASS 2013 (as amended), GHS (as amended), WHSR Part 7.1 (as amended), EN ISO 9001, ISO 13485, EN ISO 14971, EN 13641, EN 13612, EN ISO 15223-1, EN 980, EN ISO 18113-1, EN ISO 18113-2, EN ISO 18113-3, EN ISO 23840, EN ISO 20776-1, EN 13975, EN ISO 17511, ISO 80000-1, EN ISO 15193, EN ISO 15194

POTENTIAL USER

The Alinity i Total T3 Reagent Kit is intended to be used in a variety of clinical and laboratory conditions, including, without limitation, hospitals and laboratories, whose personnel have been trained in the use of the Alinity i analyzer.

INFORMATION ABOUT PRESENCE OF DRUGS IN COMPOSITION OF MEDICAL DEVICE

Does not contain any drugs.

STERILIZATION

Medical device is not required to be sterile.

INFORMATION ABOUT POSSIBILITY OF MULTIPLE APPLICATION

Medical device cannot be re-applied.

MATERIALS OF ANIMAL ORIGIN

Materials of animal origin are tested and risks of using MD are indicated in SDS. Safety Data Sheets are available at www.abbottdiagnostics.com or contact your local representative.

INDICATIONS AND CONTRA-INDICATIONS FOR USE

Not applicable.

STORAGE AND USE CONDITION

Storage of Medical device at a temperature of 2 °C to 8 °C should be carried out in cold rooms or in refrigerators providing a regulated temperature regime with a daily temperature record.

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Date Issued September 2018

Product

Product Name	List Number	Lot Numbers
Alinity i Total T ₃ Reagent Kit	07P9420	All
	07P9430	

Reason

When Alinity i Total T₃ (LN 07P94) Controls and Patient samples are run on board with the Alinity i 25-OH Vitamin D assay (LN 08P45) depressed results may be observed as a result of reagent carryover. This kit stuffer is being issued to provide steps to ensure the performance integrity of Alinity i Total T₃ test results is maintained.

Current Status

Samples tested with Alinity I Total T₃ list numbers 07P9420 and 07P9430 may show falsely depressed control and patient sample results when run in combination with the Alinity I 25-OH Vitamin D assay (LN 08P45).

Action

Actions for use of reagents:

- Separate the Alinity i Total T₃ assay from the Alinity i 25-OH Vitamin D assay by running these tests on different systems.
Or
- Perform daily maintenance procedure 2500 Daily Maintenance (i-series) found in section 10 of the Alinity ci-series Operations Manual prior to performing batch testing for all Total T₃ samples.

Contact

If you or any of the health care providers you serve have any questions regarding this information, please contact your local area Customer Service.

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1. Budavari S, editor. *Merck Index* (11th Ed.). Rahway, NJ: Merck and Co., Inc., 1989:668.
2. Larsen PR. Triiodothyronine: Review of Recent Studies of Its Physiology and Pathophysiology in Man. *Metabolism* 1972;21:1073-1092.
3. Ekins RP, editor. *Methods for the Measurement of Thyroid Hormones*. Amsterdam: Excerpta Medica Foundation, 1979:72-92.
4. Robbins J, Rall JE. The Iodine-Containing Hormones. In: *Hormones in Blood* (3rd Ed.). London: Academic Press, 1979:1:632-667.
5. DeGroot LJ, Larsen PR, Refetoff S, Stanbury JB. Transport of Thyroid Hormone and Cell Uptake. In: *The Thyroid and Its Diseases*. New York: Wiley and Sons, 1984:62-66.
6. Wahner HW, Gorman CA. Interpretation of Serum Tri-iodothyronine Levels Measured by the Sterling Technic. *N Engl J Med* 1971;284:225-230.
7. Marsden P, McKernon CG. Serum Triiodothyronine Concentration in the Diagnosis of Hyperthyroidism. *Clin Endocrinol* 1975;4:183-189.
8. Ivy HK, Washner HW, Gorman CA. Triiodothyronine (Т₃) Toxicosis: Its Role in Graves' Disease. *Arch Intern Med* 1971;128:529-534.
9. Hollander CS, Mitsuma T, Nhei N, Shenkman L, Burday SZ, Blum M. Clinical and Laboratory Observations in Cases of Triiodothyronine Toxicosis Confirmed by Radioimmunoassay. *Lancet* 1972;1:609-611.
10. Sterling K, Refetoff S, Selenkow HA. Т₃ Thyrotoxicosis: Thyrotoxicosis Due to Elevated Serum Triiodothyronine Levels. *JAMA* 1970;213: 571-575.
11. Hollander CS, Mitsuma T, Shenkman L, Stevenson C, Pineda G, Silva E. Т₃ Toxicosis in an Iodine-Deficient Area. *Lancet* 1972;2:1276-1278.
12. Ermans AM. Disorders of Iodine Deficiency. In: Ingbar SH, Braverman LE, editors. *The Thyroid* (5th Ed.). Philadelphia: JB Lippincott Co., 1986:705-721.
13. Kaplan MM, Larsen PR, Crantz FR, Dzau VJ, Rossing TH, Haddow JE. Prevalence of Abnormal Thyroid Function Test Results in Patients with Acute Medical Illnesses. *Am J Med* 1982;72:9-16.
14. Bermudez F, Surks MI, Oppenheimer JH. High Incidence of Decreased Serum Triiodothyronine Concentration in Patients with Nonthyroid Disease. *J Clin Endocrinol Metab*. 1975;41:27-40.
15. Oppenheimer JH. Thyroid Function Tests in Nonthyroidal Disease. *J Chronic Dis* 1982;35:697-701.
16. Abuld J, Larsen PR. Triiodothyronine and Thyroxine in Hyperthyroidism: Comparison of the Acute Changes During Therapy with Antithyroid Agents. *J Clin Invest* 1974;54:201-208.
17. Felig P, Baxter JD, Broadus AE, Frohman LA, editors. *Endocrinology and Metabolism* (2nd Ed.). New York: McGraw-Hill Book Co., 1987:408-416.
18. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
19. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
20. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
21. Clinical and Laboratory Standards Institute (CLSI). *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. CLSI Document M29-A4. Wayne, PA: CLSI; 2014.
22. Clinical and Laboratory Standards Institute (CLSI). *Statistical Quality Control for Quantitative Measurement Procedures; Principles and Definitions; Approved Guideline—Third Edition*. CLSI Document C24-A3. Wayne, PA: CLSI; 2006.
23. Westgard JO. *Basic QC Practices*. 3rd ed. Madison, WI. Westgard Quality Corporation: 2010.

24. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Precision Performance of Quantitative Measurement Methods: Approved Guideline—Second Edition*. CLSI Document EP05-A2. Wayne, PA: CLSI; 2004.
25. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures: Approved Guideline—Second Edition*. CLSI Document EP17-A2. Wayne, PA: CLSI; 2012.
26. Clinical and Laboratory Standards Institute (CLSI). *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach: Approved Guideline*. CLSI Document EP06-A. Wayne, PA: CLSI; 2003.
27. Clinical and Laboratory Standards Institute (CLSI). *Measurement Procedure Comparison and Bias Estimation Using Patient Samples: Approved Guideline—Third Edition*. CLSI Document EP09-A3. Wayne, PA: CLSI; 2013.

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 **DISTRIBUTED IN THE USA BY**
Abbott Laboratories
Abbott Park, IL 60064 USA

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