
Instructions for Use
ReSolve Non-Locking Catheter
Dated September 9, 2019
Consisting of 19 pages

ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

On this 9th day of September, 2019, before me, Jennifer Fordham, a notary public, personally appeared Robert Wolfarth, Senior Manager, Regulatory Affairs proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.

Notary Public

My Commission Expires: June 27, 2021

Robert Wolfarth
Senior Manager, Regulatory Affairs
Merit Medical Systems, Inc.
09 September 2019

**MEDICAL DEVICE
INSTRUCTIONS FOR USE**

«Kits of The ReSolve Non-Locking Drainage Catheter (RNL)»


**Manufactured by
Merit Medical Systems, Inc., the USA**

PRODUCT NAME

Kits of The ReSolve Non-Locking Drainage Catheter (RNL)

DESCRIPTION OF THE PRODUCT

The ReSolve Non-Locking Drainage Catheter (RNL) with hydrophilic coating is a radiopaque catheter with multiple side holes used for percutaneous drainage. The components of the catheter allow for introduction and placement using a trocar stylette or over-the-wire method. Kits of the ReSolve Non-Locking Drainage Catheter (RNL) may be packaged in a pouch with the following components: One (1) The ReSolve Non-Locking Drainage Catheter (RNL) with hydrophilic coating (A) One (1) Metal stiffening cannula (B) One (1) Trocar stylette (C)

INDICATIONS FOR USE

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) are designed for percutaneous drainage of fluids:

- when removing the contents of the abscess in different locations;
- abdominal cavity with the formation of serous fluid due to hemodynamic disorders, such as right ventricular failure, portal hypertension (liver disease);
- abdominal cavity with irritations due to infection;
- abdominal cavity due to metastases of some tumors;
- with the termination of outflow of urine in the kidney or the proximal part of the ureter;
- to gain access to the pelvis, ureter, or bladder for a diagnostic or therapeutic procedure;
- with the formation of stoma with damage to the ureter.

CONTRAINDICATIONS

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) is contraindicated for use where percutaneous drainage catheterization is unacceptable.

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) is contraindicated for intravascular use.

RESTRICTED DEVICE

Rx only: Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of this device.

PRECAUTIONS

- Read manufacturer's instructions prior to use.
- Contents are sterile (via ethylene oxide) and non-pyrogenic.

- Do not use if packaging is opened, damaged or broken.
- For single patient use only. Do not reuse, reprocess, or resterilize. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device and/ or lead to device failure, which in turn may result in patient injury, illness, or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness, or death of the patient.
- Follow universal precautions when inserting and maintaining this device. Due to the risk of bloodborne pathogens, healthcare professionals should always use standard blood and body fluid precautions in the care of all patients. Sterile technique should always be followed.
- Do not use after expiration date.
- Store in a cool, dry place.

WARNING:

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) is not to be used to deliver nutritional supplements.

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) is not to be used in the billing system.

ADVERSE REACTIONS

- Septic shock
- Bacteremia
- Hemorrhage
- Superinfection
- Bowel transgression
- Pleural transgression
- Vascular injury
- Pneumothorax
- Skin Infection
- Catheter Occlusion
- Catheter Dislodgement

INSTRUCTIONS FOR USE PLACEMENT TECHNIQUES OPTION 1: DIRECT PUNCTURE USING TROCAR STYLET

1. Remove the stiffening cannula and trocar stylette assembly from the Kits of The ReSolve Non-Locking Drainage Catheter (RNL).
2. Ensure that the distal portion of the catheter is wet prior to placement.

CAUTION: To maximize the advantages of the hydrophilic coating on the surface of the distal portion of the catheter, wet the catheter prior to use with sterile water or saline. Keep catheter wet during placement.

WARNING: DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating.

3. Flush catheter prior to use.
4. Place the metal stiffening cannula into the catheter and tighten the Luer lock fittings.

See Figure 1. Figure 1

Figure 1

5. Remove the paper spacer from the trocar stylette. Advance the trocar stylette through the metal stiffening cannula and tighten the Luer lock fittings.

See Figure 2. Figure 2

Figure 2

6. Place the catheter/cannula/trocar assembly into the fluid collection site using standard insertion technique. Placement should be confirmed with diagnostic imaging.

7. After placement is confirmed, remove the trocar stylette and stiffening cannula.

8. Kits of The ReSolve Non-Locking Drainage Catheter (RNL) is now ready to be connected to appropriate drainage bag or tubing

9. A flush regimen should be designed for the circumstances of each patient and the protocol of the physician.

Note: Instruct patient or other healthcare personnel in appropriate device function and/ or maintenance.

INSTRUCTIONS FOR USE PLACEMENT TECHNIQUES OPTION 2:

SELDINGER ENTRY TECHNIQUE OR GUIDE WIRE EXCHANGE

1. Remove the stiffening cannula and trocar stylette assembly from catheter.

2. Ensure that the distal portion of the catheter is wet prior to placement.

CAUTION:

To maximize the advantages of the hydrophilic coating on the surface of the distal portion of the catheter, wet the catheter prior to use with sterile water or saline. Keep catheter wet during placement.

WARNING: DO NOT wipe catheter with dry gauze or any solvents because it may damage the catheter coating.

3. Flush catheter prior to use.

4. Place the metal stiffening cannula into the catheter and tighten the Luer lock fittings.

See Figure 3. Figure 3

Figure 3

5. Place catheter over appropriate guide wire into the fluid collection site. The catheter accommodates an 0.038" (0.97 mm) wire.

See Figure 4. Figure 4

Figure 4

6. Advance the catheter into the site while holding the guide wire in place. Placement should be confirmed with diagnostic imaging. After placement is confirmed, remove cannula and guide wire.

7. The catheter is now ready to be connected to appropriate drainage bag or tubing.

8. A flush regimen should be designed for the circumstances of each patient and the protocol of the physician. Note: Instruct patient or other healthcare personnel in appropriate device function and/ or maintenance.

CATHETER EXCHANGE OR REMOVAL

1. Disconnect catheter from drainage tubing or bag.

2. For catheter exchange or if access is to be maintained, advance appropriate guide wire through catheter; use diagnostic imaging to confirm wire placement. Guide wire will maintain access to drainage site, if necessary.

3. Carefully remove the catheter. Proceed with either catheter exchange or skin closure.

The FastBreak connector.

The FastBreak Breakaway Connector is a two-piece coupling with a patient side connecting to the hub of the drainage catheter, and a second side connecting to the drainage bag during placement. In use, the two (2) connected components allow fluids to pass freely. Upon application of variable strength against the catheter tubing, the coupling will detach. The FastBreak connector lessens the chance of the drainage catheter being inadvertently detached. The FastBreak connector should not be used with biliary catheters.

FASTBREAK INSTALLATION

1. Connect the male luer end of the FastBreak connector to the open end of the allpurpose drainage catheter and twist in a clockwise direction. (See Step 2)

2. Connect the open female end of the FastBreak connector to the male luer end of the external fluid reservoir and twist in a clockwise direction. (See Step 3)

3. Connect the FastBreak male to female connectors together until you hear a "click" noise. (See Steps 4 and 5)


INDICATIONS FOR USE

The FastBreak connector is intended to serve as a conduit between a general drainage catheter and an external fluid reservoir.

PRECAUTIONS

- Follow instructions for use and institutional guidelines for disposal.
- Do not store at extreme conditions.
- Do not use in billiary application REUSE PRECAUTION STATEMENT For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

	Do not use if package is damaged				
	Do not re-use. Do not resterilize				
	Consult instructions for use				
	Not for individual sale.				
Rx ONLY	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.				
	Orange (light)		Green (medium)		Yellow (strong)

**ATTACHMENT INSTRUCTIONS FOR USE
MEDICAL DEVICE**

«Kits of The ReSolve Non-Locking Drainage Catheter (RNL)»

**Manufactured by
Merit Medical Systems, Inc., the USA**

1. DEVICE DESCRIPTION

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) is packaged with the following components:

1) Kit of the ReSolve Non-Locking Drainage Catheter (RNL), catheter diameter 6,5 F (2,1mm); 7,5 F (2,4 mm); 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); length 20 cm, 25 cm consisting of:

- ReSolve Non-Locking Drainage Catheter (RNL) 'J' tip with hydrophilic coating - 1 pc;
- metal stiffening cannula – 1 pc;
- trocar - 1 pc.

2) Kit of the ReSolve Non-Locking Drainage Catheter (RNL), catheter diameter 6,5 F (2,1mm); 7,5 F (2,4 mm); 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); length 20 cm, 25 cm consisting of:

- ReSolve Non-Locking Drainage Catheter (RNL) 'J' tip with hydrophilic coating - 1 pc;
- metal stiffening cannula – 1 pc;
- trocar - 1 pc.
- FastBreak connector - not more than 6 pcs. / packaging

3) Kit of the ReSolve Non-Locking Drainage Catheter (RNL), catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); length 20 cm, 25 cm consisting of:

- ReSolve Non-Locking Drainage Catheter (RNL) straight tip with hydrophilic coating - 1 pc;
- metal stiffening cannula – 1 pc;
- trocar - 1 pc.

4) Kit of the ReSolve Non-Locking Drainage Catheter (RNL), catheter diameter 6,5 F (2,1 mm); 7,5 F (2,4 mm); 8,5 F (2,8 mm); 10F (3,3 mm); 12F (4,0 mm); length 20 cm, 25 cm consisting of:

- ReSolve Non-Locking Drainage Catheter (RNL) straight tip with hydrophilic coating - 1 pc;
- metal stiffening cannula – 1 pc;
- trocar - 1 pc.
- FastBreak connector - not more than 6 pcs. / packaging

Kits of the ReSolve Non-Locking Drainage Catheter (RNL) are designed for percutaneous drainage of fluids from body cavities with abscesses, abdominal ascites, nephrostomy. The products are intended for use in clinics, hospitals and other medical institutions.

Kits of the ReSolve Non-Locking Drainage Catheter (RNL) is designed for long-term use (up to 30 days).


The ReSolve Non-Locking Drainage Catheter (RNL) is intended for percutaneous drainage of fluids from body cavities. The RNL is offered in 6.5, 7.5, 8.5, 10 and 12 French sizes with a catheter length of 20 or 25 cm. The RNL is a single lumen catheter with a non-locking 'J' or straight tip configuration and 4 to 5 drainage holes (depending on French Size) in the distal end. The 6.5F, 7.5F, 8.5F and 10F catheters have 4 drainage holes while the 12F catheter has 5 drainage holes. The catheter shaft is radiopaque and visible under X-ray to help in catheter placement. The components of the catheter allow for introduction and placement using a trocar or the Seldinger's technique.

The 'J' shape of the catheter tip may help to keep the drainage site from collapsing as fluid is being removed.

A paper spacer is attached to the trocar to prevent catheter damage.

The connector metal cannula shows the serial number of the product within one batch

Figure 1 - RESOLVE NON-LOCKING CATHETER (J-tip Configuration)

The The ReSolve Non-Locking Drainage Catheter (RNL) is packaged with the following components:

- One (1) straight or 'J' tip drainage catheter with hydrophilic coating
- One (1) metal stiffening cannula
- One (1) trocar

All sets can be equipped with a FastBreak connector (see Figure 2). The FastBreak Breakaway Connector is a two-piece coupling with a patient side connecting to the hub of the drainage catheter, and a second side connecting to the drainage bag during placement. In use, the two (2) connected components allow fluids to pass freely. Upon application of variable strength against the catheter tubing, the coupling will detach. The FastBreak connector lessens the chance of the drainage catheter being inadvertently detached. The FastBreak connector should not be used with biliary catheters. Connectors are available in the following colors - yellow, orange, green, depending on the degree of fixation.


Figure 2

2. STERILIZATION METHOD

Merit Medical Systems, Inc. utilizes Ethylene Oxide (EO) sterilization for kit of the ReSolve Non-Locking Drainage Catheter (RNL). The EO sterilization cycle is validated in accordance to the requirements of The European International Standard ISO 11135. Annual re-validation occurs assessing a worst-case reference load based on density and boxing configuration, which allows routine sterilization of a mixed-pallet load. Merit's products meet ISO 10993-7 requirements for ethylene oxide sterilization residuals. The residual ethylene oxide content does not exceed 10^{-6} .

The protocols and reports validating the sterilization process are maintained in Merit's Document Control system and are referenced on Merit's Sterilization Summary.

This device is not intended to be re-sterilized by the user.

Kits of the ReSolve Non-Locking Drainage Catheter (RNL) is less complicated than a standard process challenge device (PCD) that represents a worst case configuration and is readily adopted into Merit's validated sterilization cycle.

The dosage test and the bioburden check according to EN ISO 11137-2 were routinely performed to confirm that the dose for sterilization continues to reach the required SAL (10^{-6}).

3. TECHNICAL CHARACTERISTICS

3.1 ReSolve Non-Locking Drainage Catheter (RNL)

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F
Catheter OD	0,085" (2,16 mm)	0,098" (2,49 mm)	0,112" (2,84 mm)	0,131" (3,33 mm)	0,158" (4,01 mm)
Catheter ID	0,070" - 0,073" (1,78mm - 1,85 mm)	0,070" - 0,073" (1,78mm - 1,85 mm)	0,070" - 0,073" (1,78mm - 1,85 mm)	0,089" ± 0,002" (2,26 mm ±0,051 mm)	0,106" ± 0,002" (2,69 mm ±0,051 mm)
Catheter working Length (cm)	20	25	25	25	25
Number of Drainage Holes	4	4	4	4	5
Drain hole width	1,20 mm - 1,90 mm	1,20 mm - 1,90 mm	1,20 mm - 1,90 mm	1,20 mm - 1,90 mm	1,20 mm - 1,90 mm
Drain hole length	2,00 mm - 4,45 mm	2,00 mm - 4,45 mm	2,00 mm - 4,45 mm	2,00 mm - 4,45 mm	2,00 mm - 4,45 mm
Hole spacing	4,45- 8,65 mm	4,45- 8,65 mm	4,45- 8,65 mm	4,45- 8,65 mm	4,45- 8,65 mm
The distance from the distal end of the catheter	9,60 mm- 15,80 mm	9,60 mm- 15,80 mm	9,60 mm- 15,80 mm	9,60 mm- 15,80 mm	9,60 mm- 15,80 mm
Liquid flow rate, ml / min	At least 30 ml / min	At least 50 ml / min	At least 50 ml / min	At least 60 ml / min	At least 80 ml / min
Radiographic visibility	The catheter should be viewed under fluoroscopy				
Compatibility guidewire	guidewire of 0.038 inches (0.97 mm) with an installed cannula needs to pass through the entire structure in the Assembly.				

Compatibility connection, Luer-type	The hub of the catheter must be compatible with swivel and fixed adapters items covered by type (hemostatic valves, Luer lock connection, syringes, catheter valves, etc.).
A cone-shaped tip	At the distal end of the tip should be tapered.
Safety of magnetic resonance imaging (MRI)	The catheter should not interfere with the MRI.
Surface	On the outer surface should not be technology and surface defects. This surface should be minimally traumatic to the vessel during use.
Outer diameter connection, Luer-type	6,60 mm
The inner diameter of the connection, Luer-type	4,0 mm
The taper connection of Luer-type	Not less 6%

3.2 Metal Stiffening Cannula

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F
Cannula Size	18	18	18	15	14
Cannula OD	0,050" (1,27 mm)	0,050" (1,27 mm)	0,050" (1,27 mm)	0,072" (1,83 mm)	0,083" (2,11 mm)
Cannula ID	0,033" (0,84 mm)	0,033" (0,84 mm)	0,033" (0,84 mm)	0,054" (1,372 mm)	0,063" (1,60 mm)
Exposed Length	6,7	8,8	8,8	9,2	11,1
Allowable angle of deflection, without compromising the integrity	15°				
Surface	On the outer surface should not be technology and surface defects.				

3.3 Trocar

Characteristic	6,5 F	7,5 F	8,5 F	10 F	12 F
The diameter of the trocar	0,038" (0,97 mm)	0,038" (0,97 mm)	0,038" (0,97 mm)	0,038" (0,97 mm)	0,038" (0,97 mm)
The length of the trocar	9,279" 23,5 cm	11,744" 29,8 cm	11,252" 28,5 cm	11,258" 28,5 cm	11,262" 28,6 cm
Surface	On the outer surface should not be technology and surface defects.				
Type of sharpening	triangular				

3.4 FastBreak

length of the connector, type "male"	20,82 mm
length of female connector	20,32 mm
outer diameter of Luer coupling	6,85 mm , thread pitch - 1,4mm

inner diameter of Luer's standard plug connection,	7,23 mm thread pitch - 1,4mm
tightness	the connector must be sealed

4. PACKAGING QUALIFICATION

For products "Merit Medical Systems, Inc." Acceptable standard warehousing, unless special conditions of storage, in accordance with the risk management reports are not listed on the packaging. The product is packaged and labeled as a sterile disposable article. The sterile barrier systems specified in the standard ISO 11607, are designed to prevent the ingress of microorganisms and ensures aseptic preservation of the product to the final application. Checking package integrity is carried out within the production process on the package.

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) drainage catheter consist of the following components: one (1) ReSolve Non-Locking Drainage Catheter (RNL) drainage catheter, one (1) metal cannula, one (1) trocar placed in an impact-resistant polystyrene tray with lid PET (polyethylene terephthalate) to support products, which are then placed in a bag of Tyvek® material sealed with a validated thermal welding process, designed to cover and protect the product during sterilization, loading, transportation, manipulation and storage.

Also, kits of The ReSolve Non-Locking Drainage Catheter (RNL) can be supplied with the FastBreak connector. The FastBreak connector is packaged in an individual package in the form of a packet of Tyvek® material, sealed with a validated thermal welding process.

Individual packings of kits of The ReSolve Non-Locking Drainage Catheter (RNL) and FastBreak connectors are connected by an external fixing package with a common label on it, which shows the composition of the kit.

Individual packaging of the drainage catheter and FastBreak connector are connected by an external fixing package with a common label applied.

Instructions for use are included in the group packaging of the sets. Sets of catheters in a group box are placed in a transport box.

Kits of The ReSolve Non-Locking Drainage Catheter (RNL) packaging qualification was conducted to meet the requirements of the following standards:

ISO 11607-1: Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems.

ISO 11607-2: Packaging for terminally sterilized medical devices - Part 2: Validation Requirements for forming, sealing and assembly processes

ASTM D4169: Standard Practice for Performance Testing of Shipping Containers and Systems.

5. LABEL

Product label formats have been developed and approved for unit product, sterile packages and transport packaging. The approved labeling is located in the product Device Master Record.

Merit utilizes the following standards for labeling content and symbol generation in order to meet current requirements of the MDD:

EN 980: Symbols for use in the Labeling of Medical Devices.



EN 1041: Information Supplied by the Manufacturer.

ISO 15223 - 1: Medical Devices - Symbols to be used with medical device labels, labeling and information to be supplied Part 1: General Requirements.

ISO 15223 - 2: Medical Devices - Symbols to be used with medical devices, labeling, and information to be supplied Part 2: Symbol Development, Selection and Validation.

FDA Guidance on Medical Device Patient Labeling; Final Guidance for Industry and FDA Reviews (4/19/2001).

Description of symbols labeling shown in the table.

	Symbol.	Explanation of symbols
1.		EC official
2.		manufacturer
3.		It sold only to doctors and hospitals
4.		For single use only.
5.		catalog number
6.		shelf life * Based on the European Regulation, the symbol for "Use By" or "date of manufacture" can be used to sterile products. In the case of "use by" date is indicated, in which the product can be used (the expiration date), after that date, it is impossible to use the product. In this case, to determine the date of production must be deducted the total period of storage of the product, manufacturer's declaration of "Use By" date.
7.		lot number
8.		sterilized ethyleneoxide
9.		attention
10.		CE Mark
11.		Pyrogen-free
12.	 MAX GUIDE WIRE (1 D38" (0.97 mm))	Compatibility with the guidewire
13.		Not for individual sale

14.		Do not re-sterilize
15.		Protect from moisture
16.		This side up
17.		Fragile, handle with care
18.		Do not use if package is damaged
19.		Refer to the instructions for use.
20.	Надпись «Нетоксично»	The inscription "non-toxic"
21.		Caution! To prevent catheter damage, fully seat cannula before removing this bag!

6. TRANSPORTATION, STORAGE AND HANDLING

Storage conditions:

Products released from the transport packaging should be kept at an ambient temperature of +5°C to +27°C and a relative humidity of 10% to 90% non-condensing.

Transport conditions:

Temperature - 35 °C to +55°C and a relative humidity of 10% to 90% non-condensing.

Terms of Use:

Ambient temperature not exceeding 25 ° C.

Temperature range for use - the body temperature of 32 ° C to 42 ° C

Relative humidity: 100% (blood contact)

Atmospheric pressure 86-106 kPa.

The shelf life of sterile - 36 months from the date of sterilization.

7. UTILIZATION

Use a set of catheters disposed of as hazardous waste epidemiologically class B (presence of contact with blood and mucosal tissues of the patient).

After the expiry date the unused sets of catheters not in contact with blood and mucous membranes should be disposed of as hazardous waste epidemiologically, close in composition to the solid waste. Dispose of used conductors shall be in accordance with the laws and regulations on the disposal of products made in the country of use of the medical device.

8. WARRANTY

Limited Warranty and Disclaimer

Merit represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Guaranteed storage period for sets of drainage catheters is 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

9. COMPLAINT

The organization, which takes in the Russian Federation claims from consumers on the quality of the product:

Limited Liability Company "Merit Technologies" (LLC "Merit Technologies"), BIN 1117746679533, TIN 7705960130

Legal address: Of the Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

Tel / fax: +7 (495) 221-89-02

Component	Description of materials
Drainage Catheter	
Luer	polycarbonate CALIBRE™ Megarad™ 2081, Trinseo LLC, USA
	Dye - polycarbonate granules PC 2081-15 CALIBRE™ Megarad™, Trinseo LLC, USA Color: Yellow, Pantone 109C classification code (for 12,5F diameter)
	Dye - pellets 10050274 TRANS PURPLE PolyOne™, Trinseo LLC, USA Color: lilac, classification code Pantone 15-3817 (for 10,5F diameter)
	Dye - pellets 10200442 TRANS BLUE PolyOne™, POLYONE CORPORATION, USA Color: blue, Pantone 3005U Classification code (for 8,5F diameter)
	acetal composite PolyOne Red 500P Acetal X HC 98354 DuPont™, POLYONE CORPORATION, USA Dye - solids fraction of titanium oxide Xenoy EXN X5600HP, SABIC Innovative Plastics B.V. USA Color: Red, Classification 186S Pantone code (for the diameter of 7,5 F)
	Dye - pellets 10200443 TRANS GREEN PolyOne™, POLYONE CORPORATION, USA Color: Green, Pantone 354U Classification code (for 6,5F diameter)
Cap	Polycarbonate CALIBRE™ Megarad™ 2081, Trinseo LLC, USA
	Pigment - titanium oxide solid fraction A32329C PolyOne™, POLYONE CORPORATION, USA Color: bright white, code classification Pantone 11-0601 TCX
	Red ink INK-PG Pad Printing TYPE UV 2002 IWT Trans Tech, ITW-PM Trans Tech, USA Color: black, code classification Pantone BLACK WITH
Catheter shaft, all sizes	The polyurethane polymer CARBOTHANE™ PC-3595A-B20 MM, THE LUBRIZOL CORPORATION, USA
Catheter shaft tip sleeve (6.5F &	Polyurethane CARBOTHANE™ PC-

8.5F sizes only)	3572D-B20 MM, THE LUBRIZOL CORPORATION, USA
Hydrophilic Coating	The hydrophilic urethane ester-based HydroMed™ D640, AdvanSource Biomaterials Corporation, USA
stylet trocar	
hub	polycarbonate CALIBRE™ Megarad™ 2081, Trinseo LLC, USA
Trocar	stainless steel alloy AISI 304 (composition: C -0.08%, Mn - 2.0%, Si - 1.0%, Cr - 18.0 / 20.0%, Ni - 8.0 / 10.5%, Fe - 70.92-66.42%) / AISI 316 (Composition: C - 0.08%, Mn - 2.0%, Si - 1.0%, Cr - 16.0 / 18.0%, Ni - 10.0 / 14.0%, Mo - 2.0-3.0%, Fe - 68.92-61.92%), Ulbrich Stainless Steels & Special Metals, Inc, USA
The metal cannula	
hub	polycarbonate CALIBRE™ Megarad™ 2081, Trinseo LLC, USA
	Dye - polycarbonate granules PC 2081-15 CALIBRE™ Megarad™ Color: Yellow, Pantone 109C classification code (for 12,5F diameter), Trinseo LLC, USA
	Dye - pellets 10050274 TRANS PURPLE PolyOne™, POLYONE CORPORATION, USA Color: lilac, classification code Pantone 15-3817 (for 10,5F diameter)
	Dye - pellets 10200442 TRANS BLUE PolyOne™, POLYONE CORPORATION, USA Color: blue, Pantone 3005U Classification code (for 8,5F diameter)
	acetal composite PolyOne Red 500P Acetal X HC 98354 DuPont™, POLYONE CORPORATION, USA Dye - solids fraction of titanium oxide Xenoy EXN X5600HP. SABIC Innovative Plastics B.V. USA Color: Red, Classification 186S Pantone code (for the diameter of 7,5 F)
	Dye - pellets 10200443 TRANS GREEN PolyOne™, POLYONE CORPORATION, USA Color: Green, Pantone 354U Classification code (for 6,5F diameter)
The metal cannula	stainless steel alloy

	AISI 304 (composition: C -0.08%, Mn - 2.0%, Si - 1.0%, Cr - 18.0 / 20.0%, Ni - 8.0 / 10.5%, Fe - 70.92-66.42%) / AISI 316 (Composition: C - 0.08%, Mn - 2.0%, Si - 1.0%, Cr - 16.0 / 18.0%, Ni - 10.0 / 14.0%, Mo - 2.0-3.0%, Fe - 68.92-61.92%), Ulbrich Stainless Steels & Special Metals, Inc, USA
lubricant	Medical Silicone DOW CORNING (R) MDX4-4159 + 50% lubricant medical purpose, Dow Corning Corporation, USA
FastBreak	Polycarbonate copolymer CALIBRE™ Megarad™ 2081 and polyester resin, Trinseo LLC, USA. Dye - a solid fraction of titanium oxide Xenoy HX5600HP - 8H9D370, SABIC Innovative Plastics B.V. USA. dye ratio of not more than 4% Color: yellow, Pantone 109C code
	Polycarbonate copolymer CALIBRE™ Megarad™ 2081 and polyester resin, Trinseo LLC, USA Dye - a solid fraction of titanium oxide Xenoy HX5600HP - 8H9D370, SABIC Innovative Plastics B.V. USA. dye ratio of not more than 4% Color: orange, code according to Pantone 1645C
	Polycarbonate copolymer CALIBRE™ Megarad™ 2081 and polyester resin, Trinseo LLC, USA. Dye - a solid fraction of titanium oxide Xenoy HX5600HP - 8H9D370, SABIC Innovative Plastics B.V. USA. dye ratio of not more than 4% Color: green, Pantone 355C code
	Polycarbonate copolymer CALIBRE™ Megarad™ 2081 and polyester resin, Trinseo LLC, USA. Dye - a solid fraction of titanium oxide Xenoy HX5600HP - 8H9D370, SABIC Innovative Plastics B.V. USA. dye ratio of not more than 4% Color: blue, Pantone Blue C code
	Polycarbonate copolymer CALIBRE™ Megarad™ 2081 and polyester resin, Trinseo LLC, USA. Dye - a solid fraction of titanium oxide Xenoy HX5600HP -

	8H9D370, SABIC Innovative Plastics B.V. USA. dye ratio of not more than 4% Color: red, Pantone 711C code
O-Ring	Ethylene Propylene Rubber EPDM, Dow Corning Corporation, USA
	Silicone Dow Corning 360, 20 CST, Dow Corning Corporation, USA

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