IFU

11107

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Registered No. 599

NOTARIAL CERTIFICATE

This is to certify that Satoshi Kabaya , an agent of Masashige Kuwayama , CEO of REPROLIFE Inc. , has stated in my very presence that he acknowledged himself to have signed to the attached document.

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

2. has been signed by

Motohiro Yano

3. acting in the capacity of Notary of the Tokyo Legal Affairs Bureau

4. bears the seal/stamp of

Motohiro Yano, Notary  
Certified

5. at Tokyo

6. FEB. 27. 2019

7. by the Ministry of Foreign Affairs

8. 19-N0017599

9. Seal/stamp:

10. Signature

T. TANAKA

Toshie TANAKA

For the Minister for Foreign Affairs

## DECLARATION

I, Masashige Kuwayama, CEO of REPROLIFE Inc. (2-5-3 9F, Shinjuku, Shinjuku-ku, 160-0022 Tokyo, Japan), the undersigned, do hereby solemnly and sincerely declare:


That the attached documents:

CRYOTEC IFU

is a true copy issued for me.

And I make this solemn declaration conscientiously believing the same to be true and correct.

Dated this 27<sup>th</sup> day of February, 2019

  
\_\_\_\_\_  
Masashige Kuwayama  
CEO

Intention: For RU

The name of the product

CRYOTEC IFU (English) - Russian Market

Medical device manufacturer

REPROLIFE Inc., Japan

1-5-101 Shibuya, Shibuya-ku, Tokyo, 150-8501

Japan

Tel: +81 (0)3 5561 1111

info@reprolife.jp

Address of place of production of medical

products

REPROLIFE Inc., Japan

1-5-101 Shibuya, Shibuya-ku, Tokyo, 150-8501

Japan

Tel: +81 (0)3 5561 1111

REPROLIFE Inc., Japan

# CRYOTEC IFU

For Russian Market only

## Instructions for HCP

### The name of the product

Carrier for vitrification "CRYOTEC", different colors.

#### Medical device manufacturer

REPROLIFE, Inc., Japan  
2-5-3-9F Shinjuku, Shinjuku-ku Tokyo, 160-0022,  
Japan  
Tel. +813 59258931  
info@reprolife.jp

#### Addresses of places of production of medical products

REPROLIFE, Inc., Japan  
2-5-5 9F Shinjuku, Shinjuku-ku Tokyo, 160-0022,  
Japan  
Tel. +813 59258931

### Intended use of the medical device

The carrier is used for vitrification in liquid nitrogen, storage and subsequent defrosting of oocytes and embryos of the patient during the use of assisted reproductive technologies (art), ensuring the survival of the biomaterial.

### Potential consumer

Carriers for vitrification "CRYOTEC" are intended for use in medical institutions accredited to assist in the field of assisted reproductive technologies (art), such as clinics, laboratories and IVF centers.

Work with the products should be carried out by qualified specialists-embryologists with the necessary specialization, with relevant professional knowledge and skills, as well as experience in the field of assisted reproductive technologies (art).

It is recommended that special training in vitrification and thawing, because the method has a lot of special technical nuances that are important for minimising cellular damage.

### Indications for use of medical devices

For placement and storage of patient's oocytes and embryos in liquid nitrogen during the application of assisted reproductive technologies (art)

### Contraindications to the use of medical devices

There are no special contraindications.

### Conditions of use

Carriers for vitrification "CRYOTEC" can be used only in a specialized laboratory with a cryopreservation IVF with conditions for working with liquid nitrogen. Use of the product in other conditions is not allowed.

Carriers for vitrification "Cryotech" are used at a temperature of 10 to 35°C, while the product has a long - term resistance to low temperatures (- 196 ° C, the temperature of liquid nitrogen in the Dewar vessel).

CRYOTEC vitrification media should be stored and used in accordance with the instructions on the labels and in the instructions for use.

### Sterilely

Carrier for vitrification "Cryotech" - radiation sterilization. Re-sterilization is not allowed.

Check the expiration date on the package. The use of products with expired shelf life is not allowed.

Read the product safety data Sheet!

It is strictly forbidden to mix these solutions with solutions for vitrification of oocytes and embryos and defrosting after vitrification of other manufacturers.

### Precautionary measures

The staff of the institution must comply with the generally accepted safety measures and precautions in the laboratories, use personal protective equipment – special clothing, gloves and goggles.

**Do not use the product if the packaging is damaged**

Not subject to reuse and sterilization

**The order and conditions of use of the product**

The product is used only during the process of using assisted reproductive technologies (art). The user must first become familiar with the procedure for the preparation of biomaterial, which is described in the relevant instructions for the use of kits for vitrification and defrosting of oocytes and embryos.

Carefully inspect the package for visible damage. Do not use the product if there are obvious signs of opening and/or damage to the integrity of the package. Inspect the marking of the product and make sure that the expiration date has not expired. Do not use expired products.

Carriers for vitrification "CRYOTEC" can be used only in a specialized laboratory with a cryopreservation IVF with conditions for working with liquid nitrogen. Use of the product in other conditions is not allowed.

**Vitrification**

- \* Warning! The process of preparing the oocyte/embryo for freezing (vitrification) is described in detail in the instructions for use of the vitrification kit.
- \* Select the prepared oocyte/embryo with the tip of the pipette and place it on the tip of the CRYOTEC carrier along with the minimum volume of the VS solution.
- \* Immediately immerse the CRYOTEC in liquid nitrogen.
- \* Close the CRYOTEC cap and store in a Dewar vessel.

Please keep CRYOTEC in liquid nitrogen during the whole period of storage.

**Thawing**

- \* Warning! The process of preparing the oocyte / embryo for defrosting and further actions are described in detail in the instructions for use of the defrosting kit.
- \* Quickly remove from liquid nitrogen and place the CRYOTEC carrier in a well with the appropriate solution.

**Attention!**

1 CRYOTEC carrier may carry more than one oocyte/embryo (the amount recommended by the manufacturer does not exceed 3-4 oocytes or 4 cell embryos or 1 blastocyst).

**Storage conditions**

Storage temperature 10-20 °C, avoiding direct sunlight. Use caution when carrying.  
Do not use after the expiration date.

**Transportation conditions**

temperature range -80°C and +60°C  
-relative humidity 40-60%,  
-atmospheric pressure – (84,0-106,7) kPa (630-800 mm Hg).  
Transportation is possible by all modes of transport in covered vehicles in compliance with the conditions of transportation and in accordance with the requirements and rules of transportation of goods applicable to each mode of transport.

**Terms of disposal of products.**

Used products are subject to disposal as potentially infected waste, in accordance with state standards (SanPiN 2.1.7.728-99 "Rules for collection, storage and disposal of waste of health facilities") and internal rules adopted in the medical institution.

Unused products with expired shelf life are subject to disposal as non-hazardous waste in accordance with the state norms and rules adopted in the medical institution.

Products cannot be disposed of with household waste.

**Maintenance and repair requirements**

The medical device does not contain components subject to maintenance and repair.

All products are designed for single use.

**Recommended reading:**

- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. *Theriogenology* 67, 732007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. *FertilSteril.* J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. *Reproductive BioMedicine Online* 14, 5
- Vajta G, Kuwayama M. Improving cryopreservation systems. 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. *Reproductive BioMedicine Online* 11:300
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured invitro following vitrification. *J Reprod Dev.* 50:685
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexual maturity on vitrification of minke whale oocytes at germinal vesicle stage. *Zygote* 12,333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. *J. Reprod. Develop.* 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro- matured oocytes. *Biology of Reproduction.* 71, 432

This product should only be used in ART procedures.

Authorized representative of the Manufacturer in the Russian Federation: LLC "Intergen", Russia  
Address: 4A Nagornaya str., Moscow, 117186,  
Tel. (495) 724-38-62, t / Fax (499) 969-83-38  
info@intergen.ru

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- Kuwayama M. Highly efficient vitrification for cryopreservation of human oocytes and embryos: The CryoTop method. Theriogenology 67, 732007.
- Cobo A, Kuwayama M. Comparison of concomitant outcome achieved with fresh and cryopreserved donor oocytes vitrified by the Cryotop method. Fertil Steril. J89(6): 1657-64, 2007.
- Antinori M, Licata E, Dani G, Cerusico F, Versaci C, Antinori S. vitrification of human oocytes results in high survival rate and healthy deliveries. Reproductive BioMedicine Online 14, 5
- Vajta G, Kuwayama M. Improving cryopreservation systems. 65(1), 236-44, 2006.
- Kuwayama M. Highly efficient vitrification method for cryopreservation of human oocytes. Reproductive BioMedicine Online 11:300
- Ushijima J, Kuwayama M. High survival rate of bovine oocytes matured in vitro following vitrification. J Reprod Dev. 50:685
- Fukui Y, Kuwayama M. Effect of cryodevice type and donor's sexal maturity on vitrification of minke whale oocytes at germinal vesicle stage. Zygote 12,333-338, 2004.
- Hochi S, Kuwayama M. Improved Survival of Vitrified in vivo embryos. J. Reprod. Develop. 50, 481-486, 2004.
- Esaki R, Kuwayama M. Cryopreservation of porcine embryos derived from invitro- matured oocytes. Biology of Reproduction. 71, 432

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