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1.1 General Information of the Product

Performance and usage

Dental Implant System Oneplant with kits is intended to be surgically placed in the bone of the upper or lower jaw arches to provide support for prosthetic devices, such as artificial teeth, and to restore the patient's chewing function.

Purpose for implants

Inserted to support prosthesis like artificial dental implants. It is a substructure used for patient masticatory function recovery in 2-stage use.

Purpose for components

It is the implant to be inserted to maintain the dental prosthesis such as artificial teeth for the recovery of dental implant patient's chewing function. This is upper structure of screw, Cylinder, abutment, etc.

Purpose for instruments and kits

-Surgical Drill for Dental Implant system.

This product is a surgical drill for dental implant system which is used for inserting the fixture variously depending on the bone quality of the patients, Designed to make safer path and easier to secure the bone chip through low speed drilling technology.

-Surgical Screw Driver for Dental Implant system.

This product is a surgical screw driver for dental implant system which is used for inserting or attaching the implant fixture and abutment.

Its shape is composed of Ratchet type, Square type, and Handpiece type. It is designed to twist for the insertion of fixture, prosthetics, and drill.

-Surgical instrument for dental implants

This product is a surgical instrument for dental implants, which is designed to measure or induce the location, depth, and direction of the fixtures and prosthetics depending on the bone quality of the patients for the successful surgery.

Indications

- 1) Traumatic tooth loss without alveolar bone injury
- 2) Root treated tooth planned for extraction
- 3) Root fracture
- 4) Tooth loss due to severe caries or periodontal disease without pyogenic exudate or acute
- 5) Congenital permanent tooth deficit and remaining milk tooth with resorbed root
- 6) Internal and external root resorption
- 7) Root resorption after orthodontia

- 8) Sufficient and healthy soft tissue able to be sutured
- 9) If there is sufficient bony tissue proximal to the tooth extraction wound

Contraindications

In following conditions you should reconsider implant operation

- 1) Intraoral contraindications A. In cases with insufficient bony tissue so that severe bone resorption is predicted, or if there is insufficient remaining bone for early-fusion in the proximal tooth extraction wound. B. Disorder in mastication or functional relation C. Pathologic condition of the alveolar bone D. Prior radiotherapy on jawbone E. Xerostomia F. Pathologic change of oral mucosa (vitiligo, lichen planus, stomatitis) G. Macroglossia H, if vital anatomical structures are nearby I. Cellulitis in surrounding soft tissues J. If there are not sufficient soft tissues or its condition is poor
- 2) Transient contraindications A. Acute inflammatory disease or infection B, Pregnancy C. Temporary effect of specific **drugs** (anticoagulant, immune-suppressant) D. Mental, physical fatigue
- 3) Psychological contraindications 'A. Poor compliance B. Alcohol or other substance abuse C. Neurosis, psychosis patient D. Troublesome patient
- 4) General medical contraindications A. General/nutritional condition - age (obesity, cachexia, Syr survival rate) B. Current medications (corticosteroid, long-term antibiotic treatment) C. Metabolic disorder (pubertal diabetes, overt hyperglycemia ($>300\text{mg/dl}$)) D. Hematologic disorder (disorder of RBC, WBC, coagulation) E. Cardiovascular diseases (atherosclerosis, overt hypertension ($>300\text{mmHg}$)) F. Metabolic disorder of skeletal system (osteomalacia, Paget's disease, menopausal osteoporosis) G. Connective tissue disease (dermatosclerosis, rheumatoid arthritis) H. Implant as potential infection focus (prosthetic valve, bacterial endocarditis)

Warning

- 1) Implant operation should be performed by skilled dental surgeon because mishandled procedure may damage the implant or recipient bone
- 2) Implant is not to be recycled and it should be used for its original purpose
- 3) Damaged or mishandled implant should be removed
- 4) Inappropriate implant selection and improper implantation site or unstable fixation may shorten the lifespan of the implant
- 5) Defective product should be withdrawn,
- 6) Handle the implant carefully to prevent any damage or deformation
- 7) When using the ID excessive torque may lead to fixture and operation. failure. We recommend using less than 80Ncm of torque.
- 8) Never Reuse, reclean or resterilize the dental implant. These activities can adversely affect implant materials and alter the surface characteristics, which may result in poor function and implant failure.

Cautions for the usage

- 1) The operation must be performed by well-trained qualified dental specialist.
- 2) While drilling, you must follow the procedure written in the catalog and the fixation part should be sufficiently implanted.

- 3) Take care not to allow any soft tissues to come between the implant and prosthesis and disrupt the attachment.
- 4) If worn instruments are used during the operation, other components may be damaged. Therefore, inspect the condition of the instruments before every operation.
- 5) The product is wrapped in aseptically via gamma ray sterilization so you should open the packaging just before use (and use it immediately after the packaging has been opened).
- 6) If the package has been damaged, you must discard the product because aseptic condition cannot be maintained. (unable to use)
- 7) Bone diseases (osteoporosis, osteomalacia), bone metabolic disorders, etc. If there are, consider them carefully before treatment
- 8) The IU Implant System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of TU Leet in the MR environment is unknown. Scanning a patient who has this device may result in patient injury

Complications General complications after intraoral implant surgery include local hemorrhage, edema and hematoma. Transient loss of taste sense and masticatory function may occur. Additionally, following complications may develop: -iatrogenic trauma of surrounding tissues (lower alveolar nerve injury or sensory change, injury or hemorrhage in maxillary sinus or nasal cavity) -Insufficient or failed bony fusion - Wound dehiscence on sutured site -Delayed recovery, edema due to anesthesia - Mucositis around implant due to insufficient adhesive soft tissue -Incomplete implant placement due to insufficient bone removal or overt compression -General hypersensitivity reaction

Instructions for operation and use

1. Preparation before use

- 1) Before clinical use, the dentist must be well acquainted with the surgical procedure and the product, and he has to inform the patient about the limitations of the implant. The patient should be well aware of the functional and aesthetic limits of the implant.
- 2) Because proper selection and fixation of the implant is closely related to the life-span of the implant, you must follow the indications, contraindication, cautions and recommendations strictly.
- 3) When handling the implant, without careful review of the patient's condition and establishment of proper diagnosis and therapeutic plans, you may damage the implant.
- 4) You should choose the appropriate product based on careful review of the patient's X-ray picture.
- 5) Check whether the product's expiration date hasn't expired, packaging is sound and the product isn't damaged.

- 6) As this product is supplied aseptically, you may not use it if the packaging is damaged or torn.
- 7) The contamination of surgical instruments may lead to loss of the implant so check the hygiene and preparatory state of the proper surgical instruments.
- 8) Check for foreign-body before use.

2. Instructions and procedural sequence During the diagnostic and therapeutic planning, you must exclude any patient with local lesions or other contraindications and choose candidates who have proper bony condition to undergo implant surgery. Before proceeding the operation, you must sterilize operation room and patient's oral cavity and perioral area thoroughly. After proper draping, perform local anesthesia and make an incision on the implant site and form a flap. Exposure the implant site sufficiently and proceed to implant surgery.

- 1) Implant site preparation to implant the fixture, an adequate hole must be made in the recipient bone. To place the fixture accurately in the selected site, a hole must be made according to the size of the artificial dental prosthesis, using the respective instruments (drilling, tapping). Rotatory speed during these procedures must be adjusted taking the recipient bone condition and type of equipment used into consideration. The maximal permissible rotatory speed for the drill is generally 1,500~2,000rpm and 20~30 rpm for the tap drill. The procedure should be performed using adequate normal saline to prevent heat necrosis of the alveolar bone.
- 2) Placing the fixture Remove the wrapping and place the fixture into the epiphysis using the placement driver and adapter. Insert the fixture into the prepared implant site to the planned depth, using a handpiece or hand wrench with a clockwise rotation of 20030rpm velocity and 20~45Ncm rotatory force If the maximal torque is below 25Ncm, use equipment one unit larger.
- 3) Inserting the cover screw After the fixture has been placed, attach the cover screw using a 0480 driver below 10Ncm torque. Make sure there are no foreign bodies inside, and suture the operation site.
- 4) Connecting the abutment Bony fusion of the fixture requires 3~4 months for the mandible and 6~8 months for the maxilla. During this period, expose the implant and connect the healing abutment to enhance mucosal healing.
- 5) Prosthesis attachment After a healing period of between 2~4 weeks, connect the impression post to obtain an impression and manufacture a dental mockup. Attach the prosthesis produced in a dental lab to the oral abutment.

5. Complications General complications after intraoral implant surgery include local hemorrhage, edema and hematoma. Transient loss of taste sense and masticatory function may occur. Additionally, following complications may develop: -iatrogenic trauma of surrounding tissues (lower alveolar nerve injury or sensory change, injury or hemorrhage in maxillary sinus or nasal cavity) -Insufficient or failed bony fusion -Wound dehiscence on sutured site -Delayed recovery, edema due to anesthesia -Mucositis around implant due to insufficient adhesive soft tissue -Incomplete implant placement due to insufficient bone removal or overt compression -General hypersensitivity reaction

Precautionary Measures to Take before Getting Dental Implants

1. Talk It Out

Discuss with your dental care provider to **get** all your doubts related to dental implants cleared. Start making a list of the surgery-related questions before the consultation. This will help you **gain** confidence in the procedure and **minimize** your fear.

Also, it is important that your dental care provider is aware of your complete medical history including any allergies, current medical conditions, and any medicines and/or over-the-counter supplements you take.

2. Don't Forget Your Medicines

Your dental care provider will prescribe certain medication for you before the surgery. Please make sure that you follow the instructions as not taking these medicines can have lasting effects on your health post the surgery.

For example, if you have a heart condition, your dental care provider may ask you to take certain antibiotics before the surgery to mitigate the risk of infection. In case of any confusion/queries with respect to any medicine, call your dental clinic.

3. Eat Well before the Procedure

You will not be able to consume any eatables for a couple of hours post the surgery, which is why you need to eat a wholesome meal beforehand. After the surgery, do not try to eat anything until your dental care provider gives you the permission. The surgery site will be numb for quite a while after the procedure.

4. No Heavy Workout

You will not be in a position to exert energy in any manner or perform heavy-duty physical workouts for a few days after the surgery. So, avoid planning such activities in advance.

5. Stop Smoking

If you are a smoker, you will need to abstain from smoking for a week to ten days before the date of your surgery.

6. Bring the Companion with You

You're going to be under sedation during the surgery. Therefore, it is suggested to have someone you trust close by as you will probably not be able to go back to your home by yourself after the procedure.

7. Discuss about Pain Control Options

Discuss the pain control measures with your dental care provider in advance to know exactly what you're getting into. Your options may include local anesthesia, sedation or general anesthesia. Depending on the decision, your dental care provider can advise you about your eating and drinking schedule before the surgery.

8. Make Sure Your Teeth are Bright Clean

It is important that your teeth are cleaned thoroughly before the dental implant's procedure. Teeth need to be free of bacteria, infections, and plaque as only then the surgery will be successful. If you forget this precaution, your implants will fail and instigate you unnecessary pain.

9. Prepare Your Treatment Plan

Apart from dealing with the comprehensive dental exam, x-rays, and models of your teeth and mouth, you need to plan the treatment plan with your dental care provider. The plan will depend on certain conditions such as the number of teeth that need to be replaced and the health of your jawbone.

Various types of dental care specialists, such as periodontists and maxillofacial surgeons, may be involved when finalizing the treatment plan.

10. Choose the Right Care Provider

When getting dental implants, the last thing you want is to travel long distances to and from the dental specialist's clinic. Therefore, it is recommended to pick a locally skilled practitioner and make the whole process that much convenient for yourself. For example, if you reside around Moscow, pick someone who does dental implants in Moscow.

Instruction for Use&Order of Procedure

Dental Implant Fixture

1) Preparation before use

- ① Before clinical use, the doctor must fully understand the order of procedure and the product and must explain the limitation of implant to the patient. patient must recognize the function and esthetical limitation of the implant.
- ② Appropriate selection and fixation is major factor affecting the life of the implant. It is recommended to strictly abide by the indications, contraindications, precautions and warnings,
- ③ It is recommended to analyze the detailed status of patient, pre surgery diagnosis and treatment planning before the procedure. If not properly followed, it may cause the failure of the implant.
- ④ By reading X-ray images, applicable product may be selected in accordance to each patient's characteristics by reading X-ray image.
- ⑤ Packaging condition, expiration date and any other damages should be observed and confirmed.

⑥ This product is supplied in sterilized condition. If packaging is damaged or opened in any way, it is prohibited to use.

⑦ Implant failure may be caused by contaminated or infected surgical instruments. Check for the sanitation and availability of the surgical instruments in accordance to each treatment planning.

⑧ Sterilized medical devices (Fixture) method Gamma sterilization

2) Instruction for Use

① X-Ray image of patient is used to determine proper osteotomy site for implant.

② On selected surgical site, incision is made on the gum using the scalpel.

③ Once bone is exposed, the guide drill is used to position the osteotomy site.

④ Initial drill is used to drill initial hole to match the length of the fixture. The depth of the hole may be measured using depth gauge.

⑤ If more than one osteotomy site is required, use the positioning guide on the first osteotomy site and continue with the drilling of the next site.

⑥ Final drill is selected according to the diameter of the fixture and start drilling by maintaining the appropriate drilling speed.

⑦ Confirm the accurate connection between the Driver connector and fixture located within ampoule. Afterward, fixture is placed by rotating clockwise in the osteotomy site. It is recommended not to exceed 50Ncm of torque.

⑧ If unclear about the depth of the osteotomy site, use Torque Wrench to adjust the depth.

⑨ After placing fixture, suitable cover screw or closing screw is selected to protect the fixture. Connect the selected screw with the torque below 15Ncm and suture the surgical site.

⑩ To prevent contamination and infection in the oral cavity, recommended maintaining the cleanliness and prohibit drinking and smoking.

⑪ After 3~6 months of healing time, verify the osseointegration and proceed with the prosthesis.

Prosthetic Components

1) Preparations before use

- ① Operators must understand the operating procedures and products completely before clinical use. The operators must explain the limitation of implants to patients and the patients must clearly understand the limitations of the implants functions & esthetics.
- ② The packaging must be opened immediately before use in front of the patient to check for abnormality of product. The product cannot be used when impurity or foreign substance is found.
- ③ Check abnormality and foreign substance in the product prior to use.
- ④ As selecting and fixing implants properly have great influence on the life of the implants, it is essential to follow indicants, contraindications, cautions, and warnings precisely.
- ⑤ Detailed conditions of the patient must be inspected and plans for diagnosis and treatment need to be sufficiently devised before the operation.
- ⑥ In treatment of graft, oral health of the patient must be inspected in detail prior to operation to check for biological and somatological factors that may impede the outcome of operation.
- ⑦ Plans for diagnosis and treatment should precede operation. Loss of graft may occur when patient evaluation and plans for diagnosis and treatment are not carried out properly.
- ⑧ Upper structure transplanted by high pressure steam sterilization method (pursuant to ISO 17665-1, 2) is used after sterilizing for 3 minutes at 135 °C using high pressure steam sterilizer.
- ⑨ After sterilizing with a high pressure steam sterilizer, the implants should be dried for 30 minutes.
- ⑩ The treatment of all equipment to be used, must be sterilized in advance.

2) Instruction for Use

- ① After the soft tissues are cured, use the impression materials to take an impression.
 - ② Compare the impression and X-Ray image. Select the type and size of abutment considering the relationship with adjacent teeth.
 - ③ Attach the selected abutment in the mouth.
 - ④ Cut the attached abutment to an appropriate length to modify its shape.
 - ⑤ Take the final impression from the modified abutment and manufacture final prosthesis.
- During manufacture of prosthesis, plastic sleeve is burnt and removed.

⑥ After mounting the final prosthesis on top of graft or abutment in the mouth, screw is connected and occlusion between upper and lower jaws is checked for completion.

KITS – Drill & Instruments

① For operation procedures and detailed instructions and procedures for each product, please see our manual, catalogue, or website. (www.warantec.com/customers/catalog)

② For product code, specification, date of manufacture, expiration date, please see the label on our product.

3. Precaution during use

Dental Implant Fixture

1) Insertion Torque requirements for placement in the oral cavity

① Guide Drill, Initial Drill, Lindemann Drill, Pilot Drill – Recommended Speed: 1200 ~ 1500 rpm, Insertion Torque: 15 Ncm

② Tapered Drill, Straight Drill, Profile Drill – Recommended Speed: 1200 ~ 1500 rpm, Insertion Torque: 15 Ncm

③ Tap Drill - Recommended Speed: 700 ~ 1000 rpm

③ Fixture Driver – Recommended Speed: 30 rpm, Insertion Torque: 25 ~ 30 Ncm

④ Torque Wrench – Recommended Speed: Minimum, Insertion Torque: 25 ~ 30 Ncm

※ Fixture Driver and Torque Wrench are for Fixture insertion.

Prosthetic Components

1) Handling

① Cover screw and healing abutment are products non-sterilized by radiation.

② Before you use non-sterilized products, make sure to sterilize them at 135 °C for 3 minutes with a high pressure steam sterilizer.

③ After sterilizing with a high pressure steam sterilizer, the implants should be dried for 30 minutes.

④ As plastic products may transform during high pressure steam sterilization, more than 160 °C of temperature and 0.45 MPa of pressure should not be added.

⑤ UCLA Gold /CCM Abutment and Mill Abutment are to be fabricated with at least a minimum post height of 4 mm above the transmucosal collar.

- ⑥ Type of method used for steam sterilization is “Pre-Vaccum Displacement”
 - ⑦ Non-sterile device components are to be packaged in a “Wrapped Tray” that is a FDA-cleared sterilization accessory prior to undergoing steam sterilization.
 - ⑧ UCLA Gold /CCM Abutment and Mill Abutment are to be straight abutments only during the customization.
 - ⑨ After use of investment, it should be cooled down gradually to avoid transformation of a fine contraction and expansion when the investment is coagulated.
 - ⑩ They should be treated only with exclusive operating instruments. It is recommended to clamp with recommended torque.
- mouth are disposable, they are not be reused.

KITS – Drill & Instruments

1. Cleaning and sterilization for re-use

Surgery instruments must be cleaned and sterilized before re-use. Please follow the following procedure.

(Cleaning)

- ① Remove pollutants (skin tissue, blood, pieces of bones, etc.) off the surface of the used instrument.

Till the next step is implemented, place the instrument in a humid place.

- ② Cleaning surgery instruments with hands

- a. Put surgery instruments in water mixed with enzyme cleansing solution or similar cleansing solutions.
- b. Scrub the surface of the surgery instrument with a soft brush till all visible pollutants are removed. Wash small parts of the instrument with a cleaning needle.
- c. Remove cleansing solutions off the surface of the instrument with water.

- ③ Automatic cleaning with an ultrasonic cleaner.

- a. Put surgery instruments in a container of an ultrasonic cleaner and fill it with water completely.
- b. Wash it at room temperature for 10 or 15 minutes and drain the water.
- c. Wash the surgery instrument with water again and drain the water.

- ④ Follow the drying process of the cleaner and dry the surface of the instrument.

⑤ Repeat the cleaning process till the instrument is clean. Then, carry out a visual inspection.

⑥ Dry the cleaned surgery instrument with non-abrasive cloth and keep it before sterilization.

(Sterilization)

① Make sure to sterilize a non-sterilized product with a high pressure steam sterilize at 132 °C for 30 minutes before use.

② After high pressure steam sterilization, make sure to dry it for 30 minutes.

③ Make sure to store the sterilized surgery instrument in a clean place at 15-30 °C.

4. Disclaimer of liability

This product is part of the Dental Implant System and may only be used in conjunction with the corresponding original components and instruments according to WARANTEC instructions and recommendations. Use of products and tools made by third parties detracts from the efficient functioning of the Dental Implant System and will void any warranty or other obligation, expressed or implied by WARANTEC. Advice on the use of our products is provided verbally, in writing, by electronic media or in demonstrations. It does not relieve the user of WARANTEC products of the responsibility of determining whether or not any product is suitable for the intended purpose, indications and procedures. Use of this product is not within the control of WARANTEC. It is the responsibility of the user. All liability for loss or damage attributable to the use of this product is excluded. Within the framework of the Conditions at Sale and Delivery of WARANTEC, we warrant the quality of the Dental Implant System

5. Warranty

All WARANTEC products are guaranteed to be free from defects in material and workmanship at the time of shipping. All of our products are designed and manufactured to meet the highest quality standards. We cannot accept any liability for failure of products which have been modified in any way from their originals, for any failure due to misuse or application which is not in accordance with the designer's intentions. Products will lose their right to carry the CE sign in case they will be changed or repaired at another place than a company authorized by WARANTEC to do the change or repair.

6. Availability

Some items in the Dental Implant System are not available in all countries.

Copyright and trademarks

The products and logos from WARANTEC that are mentioned here are trademarks or registered trademarks of WARANTEC or its affiliates.

7. Conditions for use, transport and storage

Storage and shelf life.

Storage of products in the manufacturer's packaging under the following environmental conditions

Temperature: from - 5 ° C to + 30 ° C

Relative humidity: keep dry

Keep away from sunlight

The shelf life of implant is 7 years from the time of sterilization (expiration date shelf life on the package). Store in a dry place in closed packaging.

For components and instruments not shelf life.

8. Disposal

After the use the device has to be disposed as a contaminated healthcare waste. Class A.

9. Maintenance after initial use

For one time use only. Reuse is prohibited

10. prospective patients

patients with missing tooth or patients who need to remove teeth and re-implant

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
23, Dec. 1992 Under Law No.3790.

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