Operational Documentation
CorVocet Biopsy System
Dated April 21, 2020
Consisting of 18 pages


ACKNOWLEDGEMENT

State of Utah)
 §
County of Salt Lake)

On this 21st day of April 2020, before me, Jennifer Fordham, a notary public, personally appeared Robert Wolfarth, Senior Manager, Regulatory Affairs proved on the basis of satisfactory evidence to be the person whose name is subscribed to the attached instrument, and acknowledge he executed the same.

Witness my hand an official seal.

Notary Public
My Commission Expires: June 27, 2021

Robert Wolfarth
Senior Manager, Regulatory Affairs
Merit Medical Systems, Inc.
21 April 2020

Operational documentation for

MEDICAL DEVICE

«CorVocet™ biopsy systems»

Manufactured by

Merit Medical Systems, Inc., the USA

Instruction for use

DEVICE DESCRIPTION

CorVocet™ biopsy systems are a sterile single patient use device comprised of a core needle biopsy device, optional coaxial introducer and stylet. It is available in several needle gauge sizes and lengths. The top and rear actuator buttons on the device are color coded according to the various gauge sizes, Yellow=20 gauge, Pink=18 gauge, Purple=16 gauge and Green=14 gauge

INDICATIONS FOR USE

CorVocet™ biopsy systems are intended for use as a guiding needle in obtaining core biopsy samples from soft tissues such as liver, kidney, prostate, spleen, breast, lung, lymph nodes and various soft tissue tumours.

CONTRAINDICATIONS:

Physician judgment is required when considering biopsy on patients with bleeding disorders or receiving anticoagulant medications.

R_x only Federal (U.S.A.) law restricts this device to sale by or on the order of a physician trained and/or experienced in the use of this device

WARNINGS

- After use, this product may be a potential biohazard. Handle and dispose of in accordance with acceptable medical practice and with applicable local, state, and federal laws and regulations.
- Unusual force applied to the needle or unusual resistance against the needle may cause the needle to deform or bend. A bent or deformed needle may not function properly.
- Post-biopsy care may vary with the biopsy technique utilized and the individual patient's physiological condition. Observation of vital signs and other precautions should be taken to avoid and/or treat potential complications that may be associated with biopsy procedures.
- The collection of multiple core biopsy samples may help to ensure the detection of any cancer tissue. A "negative" biopsy in the presence of suspicious radiographic findings does not preclude the presence of carcinoma.
- Corvocet™ Biopsy System is not intended for use in bone.

PRECAUTIONS

- This product should be used by a physician who is completely familiar with the indications, contraindications, limitations, typical findings and possible side effects of core needle biopsy, in particular, those relating to the specific physiology being biopsied.
- For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- Inspect the product for damage. If damaged, appropriately discard the entire device and prepare a new instrument.
- These instructions for the Corvocet™ are NOT meant to define or suggest any medical or surgical technique. The individual physician is responsible for the proper procedure and techniques to be used with this product.
- The introduction of the Corvocet™ Coaxial Introducer into the body should be carried out under imaging control (ultrasound, X-Ray, CT, etc.). This product has not been tested for MR imaging compatibility.

NOTE. : If collecting multiple samples, inspect the needle for a damaged point, bent shaft or other imperfections after each sample is collected. Do not use the needle if any imperfection is noted.

POTENTIAL COMPLICATIONS/ADVERSE EFFECTS

Potential complications associated with core biopsy procedures are site specific and include, but are not limited to:

- Air Embolism
- Bleeding
- Unintended Organ Injury
- Track Seeding
- Infection
- Perforation
- Hematoma
- Hemorrhage
- Adjacent tissue injury
- Pain
- Hemoptysis
- Hemothorax
- Non-target tissue, organ or vessel perforation
- Pneumothorax

Air embolism is a rare but serious potential complication of lung biopsy procedures. Rapid deterioration of neurological status and/or cardiac arrhythmia may be indicative of air embolism. Prompt diagnosis and treatment must be considered if the patient exhibits signs or symptoms of air embolism.


INSTRUCTIONS FOR USE

BIOPSY DEVICE PREPARATION


1. Using aseptic technique remove the biopsy device from its package. Remove protective needle sheath that is secured to device and pull back slightly on the top trigger to remove packaging clip


2. Adjust device throw length to desired level (usable range is 10 mm – 25 mm). Each “click” is 0.5 mm. The colored line on the adjustment dial indicates the penetration depth.

3. Energize (prime) device by pulling back on the top or rear trigger until you feel a click. When the device is primed, a yellow indicator will show in the indicator window (labeled “ready”).


4. Advance the rear or top trigger until you feel a slight increase in resistance.
NOTE: Advancing the trigger past this point of increased resistance will cause the device to fire.

5. Device is ready to fire (primed).

BIOPSY PROCEDURE (WITHOUT OPTIONAL COAXIAL INTRODUCER)


1. Using aseptic technique, prepare site as required. Adequate anesthesia should be considered prior to procedure.
2. Verify instrument is ready to fire (primed).
3. Advance the needle tip of the biopsy device to a position proximal to the desired biopsy site.
NOTE: Once fired, the biopsy needle will travel to the preset penetration depth beyond the biopsy needle tip.
4. While maintaining the device's position and needle orientation, depress the rear trigger or push the top trigger distally to cause the needle to automatically advance.

NOTE: For devices that include a safety, it must be disengaged before the device will fire. The safety will toggle to red when ready to fire.

5. Remove the device from the patient.

6. Pull proximal on the top or rear trigger until you feel a click. See Figure below. This motion expels the specimen.

NOTE: This will prime the biopsy device.

7. Remove the specimen.

8. If additional biopsies are required, advance the rear or top trigger until you feel a slight increase in resistance.

NOTE: Advancing the triggers past this point of increased resistance will cause the device to fire. 9. Device is ready to fire (primed). Repeat steps 2-8 for all additional biopsies.

NOTE: If collecting multiple samples, inspect the needle for damaged point, bent shaft or other imperfections after each sample is collected. Do not use needle if any imperfection is noted.

10. After your last sample is collected it is recommended that you disengage the device by dry firing it before you discard it.

BIOPSY PROCEDURE (WITH OPTIONAL CORVOCET COAXIAL INTRODUCER)

Select the appropriate Corvocet™ Coaxial Introducer to correspond with the gauge and length of the Corvocet™ Biopsy System in use. Refer to catalog for product numbers and descriptions.

NOTE: An optional blunt tip stylet may be included with 18 and 20 gauge Corvocet™ Biopsy Systems. The blunt tip stylet may be used to manipulate through soft tissue and around vasculature or other organs to minimize the risk of unintentional damage to these areas.

1. Using aseptic technique, prepare site as required. Adequate anesthesia should be considered prior to procedure.

2. Using aseptic technique, remove the Corvocet™ Coaxial Introducer from its package.

3. If preferred, advance the supplied depth stop to the predetermined placement depth.

NOTE: The depth stop should be adjusted so that the coaxial introducer is in proper position when the depth stop is in contact with the skin. This will help stabilize the Corvocet™ Coaxial Introducer. 4. Using image guidance, insert the tip of the Corvocet™ Coaxial Introducer proximal to the lesion to be biopsied using the depth stop as an aid for proper placement and adjust as necessary.


5. Hold the guiding cannula needle hub. Squeeze the stylet hub and pull proximally to remove the stylet from the outer cannula. Leave the cannula in place as a guide for the placement of the Corvocet™ biopsy device.

6. Verify instrument is ready to fire (primed).

7. Place the Corvocet™ biopsy device all the way through the guiding cannula to the lesion to be biopsied.

NOTE: There should not be a gap between the hub of the Corvocet™ Coaxial Introducer and the proximal end of the biopsy device.

8. Ensure device needle tip is at the correct location. NOTE: Once fired the biopsy needle will travel to the preset penetration depth beyond the biopsy needle tip.

9. While maintaining the device's position and needle orientation, depress the rear trigger or push the top trigger distally to cause the needle to automatically advance.

NOTE: For devices that include a safety it must be disengaged before the device will fire.

10. Remove the device from the patient.

11. Pull proximal on the top or rear trigger until you feel a click. This motion expels the sample.

NOTE: This will prime the biopsy device.

12. Remove the specimen.

13. If additional biopsies are required, advance the rear or top trigger until you feel a slight increase in resistance.

NOTE: Advancing the triggers past this point of increased resistance will cause the device to fire. 14. Device is ready to fire (primed). Repeat steps 6-12 for all additional biopsies. NOTE: If collecting multiple samples, inspect the needle for damaged point, bent shaft or other imperfections after each sample is collected. Do not use needle if any imperfection is noted.

15. After your last sample is collected it is recommended that you disengage the device by dry firing it before you discard it. Recommendation: When collecting multiple samples, wipe the core tissue biopsy needle with sterile moist gauze prior to reinsertion into the guiding cannula. This will aid in proper movement of the Corvocet™ biopsy needle within the cannula

	Caution: Consult accompanying documents
	Single use only
	Do not use if package is damaged
	Catalog Number

	Lot Number
	Use By
	Sterilized Using Ethylene Oxide

Manufacturer:

Merit Medical Systems, Inc., 1600 West Merit Parkway
South Jordan, Utah 84095 U.S.A. 1-801-253-1600
U.S.A. Customer Service 1-800-356-3748

Authorized Representative:

Merit Medical Ireland Ltd,
Parkmore Business Park West, Galway, Ireland

www.merit.com

CE 0086

401157003EN_001 2016-02-11

**APPENDIX TO INSTRUCTION FOR USE
FOR MEDICAL DEVICE**

«CorVocet™ biopsy systems»

**Manufactured by
Merit Medical Systems, Inc., the USA**

Manufactured by:

Merit Medical Systems, Inc.
1600 West Merit Parkway South Jordan, Utah, 84095, USA

The place of production of the medical device

1. Merit Medical Systems, Inc., 1600 West Merit Parkway South Jordan, Utah 84095 USA
2. Merit Maquiladora Mexico, S. DE R.L. DE C.V., Avenida Sor Juana Ines De L a Cruz 19970 Interior B, Edificio 2, Parque Industrial Frontera, Tijuana, Baja California, MEXICO 22630.

The authorized representative of the manufacturer on the territory of the Russian Federation:

"Merit technologies"

Of the Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

phone: +7 (495) 221-89-02

1. DEVICE DESCRIPTION

CorVocet™ biopsy systems:

1. CorVocet™ biopsy systems, consisting of:

- biopsy needle, with a diameter of 14 G (2.1 mm), 16 G (1.6 mm), 18 G (1.25 mm), 20 G (0.9 mm), 10 cm long, 15 cm, 20 cm, 25 cm;
- coaxial introducer with a diameter of 13 G (2.45 mm), 15 G (1.8 mm), 17 G (1.5 mm), 19 G (1.1 mm); length: 4.4 cm; 4.5 cm; 4.6 cm; 9.4 cm; 9.5 cm; 9.6 cm; 14.5 cm; 14.6 cm; 19.5 cm - no more than 2 pcs. (if necessary);
- Stylet - no more than 2 pcs. (if necessary);

2. CorVocet™ biopsy systems, consisting of:

- coaxial introducer with a diameter of 13 G (2.45 mm), 15 G (1.8 mm), 17 G (1.5 mm), 19 G (1.1 mm); length: 4.4 cm; 4.5 cm; 4.6 cm; 9.4 cm; 9.5 cm; 9.6 cm; 14.5 cm; 14.6 cm; 19.5 cm;
- Stylet - no more than 2 pcs. (if necessary);

CorVocet™ biopsy systems are intended for use in obtaining core biopsy samples from soft tissues such as liver, kidney, prostate, spleen, breast, lung, lymph nodes and various soft tissue tumors.

CorVocet™ biopsy systems are sterile disposable products consisting of a biopsy needle, an optional coaxial introducer and stylet.

The procedure for sampling tissue from a patient takes no more than 60 minutes. Thus, it is believed that this medical device is mandatory for brief (less than 24 hours).

Biopsy needle is an automatic device that uses a spring coupled to a cutting needle to core soft tissue. Biopsy needle is available in several needle gauge sizes and lengths. The top and rear actuator buttons are color coded according to the various gauge sizes (e.g. yellow = 20G, pink = 18G, purple = 16G, and green = 14G).

The Coaxial Introducer is intended for use as a guiding needle in obtaining core biopsy samples from soft tissues such as liver, kidney, prostate, spleen, breast, lung, lymph nodes and various soft tissue tumors. The coaxial biopsy introducer is a three-part disposable device consisting of an outer (guiding) cannula with an attached female Luer lock hub, an inner stylet with an attached hub, and a flexible slip ring-style depth stop. The Coaxial Introducer is color coded to match the correlating gauge of biopsy device. A blunt-tipped stylet will be included with all 17G (1.5 mm) and 19G (1.1 mm) coaxial introducers to manipulate through soft tissue vasculature and minimize the risk of unintentional damage.

Fig. 1 CorVocet CorVocet™ biopsy systems, consisting of a stylet with a blunt tip (top), a biopsy needle (in the middle) and a coaxial introducer with a depth limiter (bottom)

Manufactured by:

Merit Medical Systems, Inc.

1600 West Merit Parkway South Jordan, Utah, 84095, USA

The place of production of the medical device

Merit Medical Systems, Inc., 1600 West Merit Parkway South Jordan, Utah 84095 USA

Merit Maquiladora Mexico, S. DE R.L. DE C.V., Avenida Sor Juana Ines De L a Cruz

19970 Interior B, Edificio 2, Parque Industrial Frontera, Tijuana, Baja California, MEXICO 22630.

2. PRODUCT SPECIFICATIONS

2.1 biopsy needle.

Table 1

Characteristic	14G	16G	18G	20G
The outer diameter of the cannula	2,1 mm	1,60 mm	1,25 mm	0,90 mm
Working length of the cannula	10 cm 15 cm 20 cm 25 cm	10 cm 15 cm 20 cm 25 cm	10 cm 15 cm 20 cm 25 cm	10 cm 15 cm 20 cm 25 cm
The distance between the marks depth	10mm			
bending resistance	Multiple needle deflection on 15° tube must withstand at least 20 tests			
The shape of the cannula tip	At the distal end of the cannula tip must type Francine			
Emission length (working range) of the needle	10-25 mm Step 0,5 mm			
blade edge	It must be sharpened and free of visible defects, such as burrs, by external			

Characteristic	14G	16G	18G	20G
	examination			
The surface of the cannula	On the outer surface should be free of surface defects and process.			
Maximum arm length	17 0mm			
The maximum outer diameter of handles, where the product is held by hand	17,7 mm			
The outer diameter of the stylet, trocar	1,50 mm	1,05 mm	0,85 mm	0,58 mm
Surface	On the outer surface should be free of surface defects and process.			
The product must withstand a certain number of starts without giving up	12 cycles (12 starts)			

2.2 introducer coaxial 2

Table

Characteristic	13G	15G	17G	19G
The outer diameter of the cannula	2,45mm	1,80 mm	1,45 mm	1,10 mm
Working length of the cannula	4,4 cm 4,5 cm 4,6 cm 9,4 cm 9,5 cm 9,6 cm 14,5 cm 14,6 cm 19,5 cm	4,4 cm 4,5 cm 4,6 cm 9,4 cm 9,5 cm 9,6 cm 14,5 cm 14,6 cm 19,5 cm	4,4 cm 4,5 cm 4,6 cm 9,4 cm 9,5 cm 9,6 cm 14,5 cm 14,6 cm 19,5 cm	4,4 cm 4,5 cm 4,6 cm 9,4 cm 9,5 cm 9,6 cm 14,5 cm 14,6 cm 19,5 cm
The distance between the marks depth	10mm			
The shape of the cannula tip	At the distal end of the cannula tip line			
The surface of the	On the outer surface should be free of surface defects and process.			

Characteristic	13G	15G	17G	19G
cannula				
Luer connector handles	The outer diameter of 7.60 mm The inner diameter of 4.0 mm The taper of 6%			
The outer diameter of the stylet	2,09 mm	1,48 mm	1,24 mm	0,88 mm
Surface	On the outer surface should be free of surface defects and process.			

2.3 Stylet

TABLE 3

Characteristic	Supplied with introducer 17G	Supplied with introducer 19G
The outer diameter of the stylet	1,27 mm	0,92 mm
Tip shape stylet	straight	straight
Surface	On the outer surface should be free of surface defects and process.	

3. STERILIZATION METHOD

Merit Medical Systems, Inc., utilizes Ethylene Oxide (EO) sterilization for the products described in this Technical File. The EO sterilization cycle is validated according to European International Standard ISO 11135. Annual re-validation occurs assessing a worst-case reference load based on density and boxing configuration, which allows routine sterilization of a mixed-pallet load. Merit's products meet ISO 10993-7 requirements for ethylene oxide sterilization residuals, no more than 10^{-6} .

This product does not provide a user re-sterilization.

4. DESCRIPTION OF PACKAGING

For products of the company "Merit Medical Systems, Inc." permissible standard warehousing, unless special conditions of storage, in accordance with the risk management reports are not listed on the packaging. The product is packaged and labeled as a sterile disposable article. Sterile barrier systems specified in the standard ISO 11607, are designed to prevent the ingress of microorganisms and provide an aseptic preservation of products to end use. Checking package integrity is carried out within the production process.

All products (biopsy needle, biopsy introducer, stylet) are packaged in a protective case. The biopsy needle is packed in a protective case, fastened with a packing clip.

The needle of the Corvocet™ Biopsy System is packaged in a protective sheath that is held on with a packaging clip. The entirety of the device is then placed in a thermoformed PETG blister tray sealed with a Tyvek® lid. Some configurations also include the coaxial introducer and

blunt stylet in the tray. Five (5) units are packaged per shelf carton and four (4) shelf cartons per shipper.

The stand-alone Corvocet™ Coaxial Introducer needle is packaged in a protective sheath and placed in a Tyvek. Some configurations also include one (1) blunt stylet. The coaxial introducer is packaged five (5) units per shelf carton and four (4) shelf cartons per shipper.

Instructions for use placed in group packaging.

Certification of packaging products was carried out in accordance with the requirements of the following standards:

ISO 11607-1: packaging of medical devices to be sterilized at the final stage of production -

Part 1: Requirements for materials, sterile protection systems and packaging systems.

ISO 11607-2: packaging of medical devices to be sterilized in the final stages of production -

Part 2: Requirements for validation processes forming, sealing and assembly.

ASTM D4169: Standard method of testing the efficiency of transport of containers and systems.

5. RISK ANALYSIS SYSTEM

Potential risks are evaluated during the design / product development and manufacturing processes according to ISO 14971. The appropriate risk assessment tool (FMEA, Clinical Risk Assessment - CRA.) is applied in order to identify potential failure modes and potential effect on patients or users. Appropriate control and prevention mechanisms are defined and implemented. Risk Analysis reports are revisited and updated as necessary and controlled under Merit's Document Control system. A cross-functional team, including persons with clinical experience, reviewed the Risk Analysis reports to ensure that all potential risks were addressed. *An analysis of the relationship between benefits and risks is included in the risk analysis plan.*

Important Hazards

The conclusions drawn by the cross-functional team are located in the Risk Assessment reports. The most important hazards include air embolism, vascular perforation, systemic infection, moderate blood loss, chemical toxic effects – limited or moderate organ involvement, and tract seeding.

Manufacturer statement regarding safety

Based on the documentation provided and/or referenced in this technical file describing the design, manufacture, packaging and sterilization of the Corvocet™ Biopsy System and the assessment of potential risks compared to the benefits of the use of these products through formal risk analysis, clinical evaluation, and product/process qualifications, Merit Medical Systems has determined that the Essential Requirements of the European Community Council Medical Device Directive 93/42/EEC or the General Safety and Performance Requirements of the Medical Device Regulations (EU) 2017/745 have been satisfied and that the Corvocet™ Biopsy System is safe and effective for its intended use as labeled.

Drug substance that are parts of the medical device:

The manufacturer Merit Medical Systems, Inc., USA, declares with full responsibility that **there are NO** drug substances used during the manufacture of medical device «Corvocet™ Biopsy System».




Materials of animal or human origin with indication information on their biological compatibility:

Materials of animal and human origin **are not used**.

6. INFORMATION ON LABELING OF PRODUCTS

Product label formats have been developed and approved for unit product, sterile packages and transport packaging. The approved labeling is located in the product Device Master Record. For label content and symbol generation, Merit utilizes, EN 1041, ISO 15223-1, EN 15986, and FDA Guidance on Medical Device Patient Labeling; Final Guidance for Industry and FDA Reviews (4/19/2001), in order to meet current requirements of the MDD or the MDR.

Description of the symbols used on the medical device's label are represented in the table.

	Symbol	Explanation of symbols
1.		Authorized representative in the European Community
2.		Manufacturer
3.		Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.
4.		For single use only
5.		Catalog number
6.		Use by * Based on the European Regulation, the symbol for "Use By" or "date of manufacture" can be used to sterile products. In the case of "use by" date is indicated, in which the product can be used (the expiration date), after that date, it is impossible to use the product. In this case, to determine the date of production must be deducted the total period of storage of the product, manufacturer's declaration of "Use By" date.
7.		Batch code
8.		Sterilized using ethylene oxide
9.		Caution! Consult accompanying documents
10.		This symbol means that product complies with the essential requirements of the

		Directive 93/42/EEC concerning medical devices
11.		Protect from moisture
12.		This side up
13.		Fragile, handle with care
14.		Do not use if package is damaged
15.		Refer to the instructions for use.

7. TRANSPORTATION, STORAGE AND HANDLING

Storage conditions:

The device should be stored in a cool, dry and dark place. Do not store catheters where they are directly exposed to organic solvents, ionizing radiation, ultraviolet light and high temperatures. Rotate inventory so that device is used prior to the Use By date on the package label. Products released from the transport packaging should be kept at an ambient condition:

Storage conditions:

Temperature +5°C to +27 °C (environment)

Relative humidity 10% to 90 % (non-condensing)

Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Transport conditions:

Temperature -35°C to +55 °C (environment)

Relative humidity 10% to 90 % (non-condensing)

Atmospheric pressure 84.0 kPa – 106.7 kPa (630 – 800 mmHg)

Usage environment:

Ambient temperature not exceeding 25 °C.

Temperature range for use - the body temperature of 32 °C to 42 °C

Relative humidity: 100% (blood contact)

Atmospheric pressure 86-106 kPa.

The shelf life of sterile devices - 36 months from the date of sterilization.

8. DISPOSAL

To prevent the formation of health risks when disposal of medical devices and also formation of environment risks connected with disposal of medical devices, the personnel have to dispose MD per protocol adopted in particular organization.

Medical Device "Corvocet™ Biopsy System" according to the Sanitary Rules and Regulations 2.1.7.2790-10 is related to Class B contaminated and potentially pathological wastes (epidemiologically hazardous wastes).

9. MAINTENANCE

Not subject to maintenance and repair. There are no parts that are serviced on site. This means that repair in conditions operation is not carried out (non-applicable).

10. WARRANTY

Limited Warranty and Disclaimer

Merit Medical Systems, Inc. represents and warrants that the products shall materially conform to the specifications at the time of shipment to Customer and will remain materially free of defects in materials and workmanship for the period specified on the written limited warranty separately provided with each product. Warranty shelf life of medical device Corvocet™ Biopsy System is 3 years. This warranty shall not apply to products that have been (i) modified, changed, repaired, refurbished, reprocessed, or altered by anyone other than Merit; (ii) subjected to misuse, mishandling, re-use, re-sterilization, accident, abuse, neglect or tampering; (iii) damaged by excessive physical, environmental or electrical stress; or that have had a serial number altered, defaced or removed; (iv) combined with any other product; (v) used outside the approved "Indications for Use" as cleared by the relevant competent authority, used contrary to the use outlined in the product specifications or in an application or environment for which such product was not designed or contemplated. Any technical product advice furnished by Merit is provided without duty or compensation and Merit assumes no obligation or liability, all such advice being given and accepted at Customer's risk. EXCEPT AS EXPRESSLY SET FORTH HEREIN MERIT MAKES NO REPRESENTATION OR WARRANTY OF ANY KIND, EXPRESS OR IMPLIED, TO CUSTOMER, ITS END-USERS OR ANY THIRD PARTIES WITH RESPECT TO THE PRODUCTS. MERIT DISCLAIMS AND EXCLUDES ANY AND ALL OTHER EXPRESS, IMPLIED, AND STATUTORY WARRANTIES, INCLUDING WITHOUT LIMITATION, WARRANTIES AGAINST INFRINGEMENT, AND THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

Limitation of Liability

UNDER NO CIRCUMSTANCES SHALL MERIT BE LIABLE TO CUSTOMER OR ANY OTHER PERSON OR ENTITY FOR ANY PUNITIVE, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, WHETHER BASED UPON LOST GOODWILL, LOST RESALE PROFITS, WORK STOPPAGE, PRODUCT FAILURE, IMPAIRMENT OF OTHER PRODUCTS OR OTHERWISE AND WHETHER IN A CONTRACT, TORT OR OTHER ACTION FOR OR ARISING OUT OF BREACH OF WARRANTY, BREACH OF CONTRACT, DELAY, NEGLIGENCE, STRICT LIABILITY, OR OTHERWISE.

11. COMPLAINT

The organization, which takes claims from consumers about the product quality on the Russian Federation territory:

Limited Liability Company "Merit Technologies"

(LLC "Merit Technologies"), BIN 1117746679533, TIN 7705960130

Legal address: Russian Federation, 119017, Moscow, 1-st Kazachiy pereulok 7, floor 1, room 6

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