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
Name Mary Simchik

Date 07 FEB 2020

**Instruction for use
for Medical Devices (MD)**

Non-absorbable arthroscopic suture anchors (fixation implants)

MULTIFIX® P Knotless Fixation Device

Instructions for Use

DESCRIPTION

The MULTIFIX® P Knotless Fixation Device facilitates the attachment of tissue to bone and consists of an implant and associated instruments. The implant can accommodate 2, 3 or 4 suture strands or a maximum volume of 2 strands of suture and 2 strands of ULTRATAPE®. The implant is manufactured from polyetheretherketone (PEEK) and is preloaded onto a disposable handle designed to facilitate direct insertion into bone. The suture anchor is a knotless fixation device, in other words surgical knots are not necessary for fixation of suture to tissue.

INDICATIONS FOR USE

The MULTIFIX P Knotless Fixation Device is indicated for use in the fixation of soft tissue to bone. Examples of such procedures include:

Shoulder: Bankart Repair, SLAP lesion repair, acromio-clavicular separation, rotator cuff repair, capsule shift/capsule-labral reconstruction, biceps tenodesis and deltoid repair.

Ankle: Lateral instability, medial instability, Achilles tendon repair/reconstruction and midfoot reconstruction

Foot: Hallux valgus reconstruction

Elbow: Tennis elbow repair, biceps tendon reattachment

Knee: Extra-capsular repairs; reattachment of medial collateral ligament, posterior oblique ligament or joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions.

CONTRAINDICATIONS

- Surgical procedures other than those listed in the INDICATIONS FOR USE section.
- Pathologic bone conditions, such as cystic changes or severe osteopenia, which would impair the secure fixation of the implant.
- Pathological changes in the affected soft tissues which would prevent secure suture fixation.
- Comminuted bone surfaces which could compromise secure fixation of the implant.
- Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitation, previous infections, etc.
- Conditions such as senility or mental illness which may impact the patient's ability or willingness to restrict activities or comply with physician instructions during the healing period.
- The implant is not designed for and should never be used to attach artificial ligaments or other implants.
- Consider the potential risk of a second surgical procedure before deciding to remove the device after wound closure. Implant removal should be followed by adequate postoperative management.
- This device has been designed for use with a 3.3mm drill or 3.4mm punch. Use of other drill bits or punches may compromise implant insertion and/or retention.

WARNINGS AND PRECAUTIONS

Warnings

Do not use if package is opened or damaged.

- Do not resterilize or reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases.
- The implant with inserter handle is a single use device to be used only once for a single patient.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Do not implant the anchor in poor quality bone or where bone quantity is limited. Incomplete insertion or poor bone quality may result in implant pullout.
- Failure to distribute proper tension through both suture ends may result in inadequate suture lock and potential failure.
- Do not insert an implant into a bone hole where an existing implant is present. This may cause the previously placed implant to extrude out through the bone.
- Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.
- Detailed instructions on the use and limitations of the device should be given to the patient. The implant with inserter handle is a single use device to be used only once for a single patient.

Precautions

- Implantation of this device requires an appropriate level of surgical skill and experience. Surgeons who plan to use this device should carefully review the instructions for use prior to clinical use. Additionally, completion of a skills laboratory and training by the manufacturer or one of its appointed representatives, or observation of/ assistance with similar implantation procedures, is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device.
It is the responsibility of the surgeon to advise the patient of the known risks and complications associated with the surgical procedure and device.
- Any decision to remove the implant should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- A back-up device should always be available to avoid any delay in surgery should it be required.
- Do not insert an implant into a bone hole where an existing implant is present. This may cause the previously placed implant to extrude out through the bone.

POTENTIAL ADVERSE EVENTS

As with any surgical procedure, there is risk involved. Potential complications accompanying such implant surgery may include, but are not limited to the following: complications with anesthesia; pain at the incision or surgical site; infection, both deep and superficial; bone damage or fracture; injury to surrounding tissues or vasculature; embolus or blood clotting issues (e.g., pulmonary embolus, deep vein thrombosis, etc.); nerve injury or palsy; implant or treatment failure. Secondary surgical interventions may include, but are not limited to: removal of implant; placement of new or additional implant(s).

OPERATOR TRAINING REQUIREMENTS

- Use of this medical device requires an appropriate level of surgical skill and experience.

- Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/assistance with similar surgical procedures is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using this medical device.
- The surgeon should advise the patient of the known risks and complications associated with the surgical procedure and use of this device.

PREPARATION FOR USE

Caution: This device has been qualified for use with #2 MAGNUMWIRE®, ULTRABRAID®, and ULTRATAPE sutures. Use of other suture is not recommended and may compromise suture locking integrity and result in potential failure.

Caution: When required, use only the approved drill, bone punch or bone tap instruments. Use of other ancillary instrumentation may compromise implant insertion and/or retention.

DIRECTIONS FOR USE

- 1) Suture tissue using a suture passer utilizing the appropriate suture.
 - 2) Using the Insertion Guide with Obturator is recommended. When using the insertion guide in the portal, align the slit of the insertion guide in the direction of the stitch in the tissue.
 - 3) Assess the patient's bone quality prior to use. If in the physician's opinion the bone is hard, create a bone hole using a 3.3mm drill bit or a 3.4mm punch prior to anchor insertion. Advance the drill or the punch until the shoulder reaches the bone surface. Alternately, a 4.5mm tapered bone punch advanced to the 15mm depth marker may be used. **Caution:** Care must be taken to avoid creation of a shallow bone hole. A proper hole depth is accomplished when the 4.5mm tapered punch is advanced to the 15mm mark.
 - 4) Thread the suture ends through the snare loop. Ensure about 6 cm of each suture end extends past the snare (see diagram 1). **Note:** The device is intended to be used with 2, 3 or 4 suture strands or a maximum volume of 2 strands of suture and 2 strands of ULTRATAPE.
 - 5) To snare the sutures through the anchor, use a syringe type motion to pull the suture loader cross bar (see diagram 2) away from the shaft and toward of the handle at approximately a 45° angle. Discard the suture loader after the sutures have been passed through the anchor.
 - 6) Holding the free ends of the suture, advance the loaded anchor either percutaneously or through a cannula to the surface of the bone. Take up suture slack by pulling on the suture tails until suture can be clearly seen between the tissue and the implant. Confirm that the suture is not wrapped around the device. If wrapped, pull the device out and repeat this step.
 - 7) Place the tip of the anchor at the desired location lateral to where the tendon is to be secured. Establish axial alignment of the inserter shaft perpendicular to the bone during insertion. Assure the eyelet of the anchor tip orients the sutures toward the tissue (see diagram 3). Do not split the sutures with the implant. If a bone hole was created, establish axial alignment of the inserter handle with the bone hole.
 - 8) Use a mallet to tap the proximal end of the inserter handle to drive the anchor tip until the horizontal depth mark (see diagram 4) is flush or below the surface of the bone.
- Caution:** If the anchor tip does not penetrate the bone on the first strike, create a bone hole

according to step 3. Use a new anchor into the bone hole.

Caution: Use care to properly align the implant and inserter handle with the bone hole during pound-in. Avoid excessive probing. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.

Caution: Incomplete insertion or poor bone quality may result in the implant pulling out.

9) In order to maneuver the tissue to the desired location, pull the free suture tails in line with the axis of the inserter while maintaining downward pressure on the inserter handle. When desired tension is achieved, secure the sutures in the jam cleat (see diagram 5).

Caution: Do not over tension the sutures as breakage or tissue pull through may occur.

10) While maintaining downward pressure on the inserter handle, deploy the suture lock by turning the deployment knob

(see diagram 6) two full turns clockwise. Continue to turn the deployment knob clockwise until it spins freely.

11) Check the location of the horizontal depth mark on the shaft (see diagram 7). If the mark is visible, use the mallet to tap the device further into the bone until the depth mark on the shaft is flush or below the surface of the bone.

12) Disengage sutures from the jam cleats prior to extracting the inserter handle. Trim remaining loose suture ends at bone hole.

Caution: Do not pull on loose suture ends. This may result in implant pullout. If there is a need to check fixation, probe the repaired tissue.

13) If placing more than one implant, assure that the anchors are located at least 7mm apart and repeat steps 1 through 12.

MATERIALS

- Implant: PEEK (polyetheretherketone)
- Suture: Implant grade polyethylene and polypropylene.
- Inserter Handle: Surgical grade stainless steel and medical grade plastic.
- The product is not made with natural rubber latex.
- The product does not contain DEHP.

HOW SUPPLIED

- The device is provided STERILE for SINGLE USE ONLY. The device is only sterile if the packaging is not opened, damaged or broken.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Do not resterilize or reuse the device as this may result in product malfunction, failure or patient injury, and may also expose the patient to the risk of infectious diseases.

MAGNETIC RESONANCE (MR) IMAGING

The implant is MR safe.

STORAGE

Store away from moisture and direct heat.

DISPOSAL

International and US regulations require controlled disposal of used and unused medical devices. Used devices must be disposed of in accordance with applicable government regulations and Medical Waste Management Plans/Acts.

CUSTOMER SERVICE

Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.

Product Complaints & Return Authorization All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized Customer Service representative.


Manufacturer

ArthroCare Corporation

7000 West William Cannon Drive Austin, TX 78735 USA

Phone: (800) 343-5717

www.arthrocare.com

Authorized European Representative

Smith & Nephew Medical Limited 101 Hessle Road

Hull HU3 2BN

United Kingdom

ec.rep@smith-nephew.com

SYMBOLS LEGEND	
Symbol	Definition
	Catalogue number
	Lot number
	Use by
	MR Safe
	Consult instructions for use
	Do not reuse
	Do not resterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Caution
	Sterilization using irradiation
	Manufacturer
	Authorized representative in the European Community
Rx Only	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.
	CE mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).

° Trademark of Smith & Nephew.


These products may be covered by one or more U.S. patents.

See smith-nephew.com/patents to view a listing of patents that cover our products.

© Smith & Nephew 2017.

P/N 62536 Rev.C

February 2017

MULTIFIX® S Knotless Fixation Device

Instructions for Use

DESCRIPTION

The MULTIFIX® S Knotless Fixation System consists of an implant (anchor) and associated instruments which facilitate the attachment of tissue to bone. The implant, which is available in 5.5mm and 6.5mm sizes, can accommodate 2, 3 or 4 suture strands or a maximum volume of 2 strands of suture and 2 strands of ULTRATAPE®. The implant is manufactured from polyetheretherketone (PEEK) and is pre-loaded onto a disposable handle designed to facilitate direct insertion into bone. The suture anchor is a knotless fixation device, in other words surgical knots are not necessary for fixation of suture to tissue.

INDICATIONS FOR USE

The MULTIFIX S Knotless Fixation Device is indicated for use in the fixation of soft tissue to bone. Examples of such procedures include:

Shoulder: Bankart Repair, SLAP lesion repair, acromioclavicular separation, rotator cuff repair, capsule shift/capsule-labral reconstruction, biceps tenodesis and deltoid repair.

Ankle: Lateral instability, medial instability, Achilles tendon repair/reconstruction and midfoot reconstruction

Foot: Hallux valgus reconstruction

Elbow: Tennis elbow repair, biceps tendon reattachment

Knee: Extra-capsular repairs; reattachment of medial collateral ligament, posterior oblique ligament or joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions.

CONTRAINDICATIONS

- Surgical procedures other than those listed in the INDICATIONS FOR USE section.
- Pathologic bone conditions such as cystic changes or severe osteopenia, which would impair secure fixation of the implant.
- Pathological changes in the affected soft tissues, which would prevent secure fixation.
- Inadequate bone stock or comminuted bone surface(s), which would compromise secure fixation of the implant.
 - Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitation, previous infections, etc.
 - Conditions such as senility or mental illness which may impact the patient's ability or willingness to restrict activities or comply with physician instructions during the healing period.
 - The implant is not designed for and should never be used to attach artificial ligaments or other implants.

⚠ WARNINGS AND PRECAUTIONS

Warnings

- Do not use if package is opened or damaged.
- Do not resterilize or reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases.
- The implant with inserter handle is a single use device to be used only once for a single patient.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Do not implant the anchor in poor quality bone or where bone quantity is limited. Incomplete insertion or poor bone quality may result in implant pullout.
- Failure to distribute proper tension through both suture ends may result in inadequate suture lock and potential failure.
- Do not insert an implant into a bone hole where an existing implant is present. This may cause the previously placed implant to extrude out through the bone.
- Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.
- Detailed instructions on the use and limitations of the device should be given to the patient. The implant with inserter handle is a single use device to be used only once for a single patient.

Precautions

- Implantation of this device requires an appropriate level of surgical skill and experience. Surgeons who plan to use this device should carefully review the instructions for use prior to clinical use. Additionally, completion of a skills laboratory and training by the manufacturer or one of its appointed representatives, or observation of/ assistance with similar implantation procedures, is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device.
It is the responsibility of the surgeon to advise the patient of the known risks and complications associated with the surgical procedure and device.
- Any decision to remove the implant should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- A back-up device should always be available to avoid any delay in surgery should it be required.
- Do not insert an implant into a bone hole where an existing implant is present. This may cause the previously placed implant to extrude out through the bone.

POTENTIAL ADVERSE EFFECTS

As with any surgical procedure, there is risk involved. Potential complications accompanying such implant surgery may include, but are not limited to the following: complications with anesthesia; pain at the incision or surgical site; infection, both deep and superficial; bone damage or fracture; injury to surrounding tissues or vasculature; embolus or blood clotting issues (e.g., pulmonary embolus, deep vein thrombosis, etc.); nerve injury or palsy; implant or treatment failure. Secondary surgical interventions may include, but are not limited to: removal of implant; placement of new or additional implant(s).

OPERATOR TRAINING REQUIREMENTS

- Use of this medical device requires an appropriate level of surgical skill and experience.

- Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/assistance with similar surgical procedures is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using this medical device.
- The surgeon should advise the patient of the known risks and complications associated with the surgical procedure and use of this device.

PREPARATION FOR USE

Caution: This device has been qualified for use with #2 MAGNUMWIRE®, ULTRABRAID®, and ULTRATAPE sutures. Use of other suture may compromise suture locking integrity and result in potential failure.

Caution: When required use only the approved bone punch or bone tap instruments. Use of other ancillary instrumentation may compromise implant insertion and/or retention.

DIRECTIONS FOR USE

- Suture tissue using a suture passer utilizing the appropriate suture.
- Using the Insertion Guide with Obturator is recommended. When using the insertion guide in the portal, align the slit of the insertion guide in the direction of the stitch in the tissue.
- Assess the patient's bone quality prior to use. If in the physician's opinion the bone is hard, create a bone hole using a 4.5mm Tapered Bone Punch prior to anchor insertion. A bone hole is also required if the implant is to be placed at an angle other than perpendicular to the bone surface. Advance the punch until the 25mm mark is flush with the bone surface.

Caution: Care must be taken to avoid creation of a shallow bone hole. A proper hole depth is accomplished when the punch is advanced to the 25mm mark.

- Thread the suture ends through the snare loop. Ensure about 3cm of each suture end extends past the snare loop (see diagram 1).

Note: The device is intended to be used with 2, 3 or 4 suture strands or a maximum volume of 2 strands of suture and 2 strands of ULTRATAPE.

- To snare the sutures through the anchor, use a syringe type motion to pull the suture loader cross bar (see diagram 2) away from the shaft and toward the handle at approximately a 45° angle. Discard the suture loader after the sutures have been passed through the anchor.
- Holding the free ends of the suture, advance the loaded anchor either percutaneously or through a cannula to the surface of the bone. Take up suture slack by pulling on the suture tails until suture can be clearly seen between the tissue and the implant. Confirm that the suture is not wrapped around the device. If wrapped, pull the device out and repeat this step.
- Place the pound-in tip of the anchor at the desired location lateral to where the tendon is to be secured. Assure the eyelet of the pound-in tip orients the sutures toward the tissue (see diagram 3). Do not split the sutures with the implant. Establish axial alignment of the inserter shaft perpendicular to the bone. If a bone hole was created, establish axial alignment of the inserter handle with the bone hole.
- Remove desired slack in the sutures by pulling on the suture tails.
- Use a mallet to tap the proximal end of the inserter handle to drive the pound-in tip into the bone. Tap the proximal end of the inserter until the laser mark on the leading thread of the screw body is sub-cortical (see diagram 4).

Caution: If the pound-in tip does not penetrate the bone on the first strike, create a bone hole according to step 3. Use a new anchor in the bone hole.

Caution: Use care to properly align the implant and inserter handle with the bone hole during pound-in. Avoid excessive probing. Do not bend or twist the inserter handle during

and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.

Caution: Incomplete insertion or poor bone quality may result in the implant pulling out.

- Tension will be applied to the tissue during advancement of the pound-in tip. If additional tension is required individually pull the free suture tails in line with the axis of the inserter while maintaining downward pressure on the inserter handle. Suture may be placed in the suture cleats for suture management if desired (see diagram 5).

Caution: Individually tensioning the sutures may cause the device to back out of the bone hole. If this occurs, advance the pound-in tip until the laser mark on the leading edge of the screw is sub-cortical before applying tension to the next suture.

Caution: Do not over tension the sutures as breakage or tissue pull through may occur. Additionally, over tensioning the suture may result in incomplete insertion or implant pull out.

- While continuing to maintain downward pressure on the inserter handle, deploy the screw by turning the deployment knob (see diagram 6) clockwise until the screw is flush or below the surface of the bone (see diagram 7) and the deployment knob has come to a hard stop.
- Disengage the inserter handle by turning the end cap counterclockwise a minimum of six turns (see diagram 8).
- If implant extraction is indicated at this point, turn deployment knob counterclockwise to remove the implant screw and pull on the suture tails to remove the implant. If implant has detached from inserter handle, or for post-operative repair, insert the Anchor Extraction tool and turn counterclockwise to remove the implant screw and then pull on the sutures to remove the implant.
- Disengage sutures from the cleats if needed prior to retracting the inserter handle. Trim remaining loose suture ends at the bone hole.
- If placing more than one implant, assure that the anchors are located at least 7mm apart and repeat steps 1 through 14.

MATERIALS

- Implant: PEEK (polyetheretherketone)
- Suture: Implant grade polyethylene and polypropylene. The suture is USP size 2 (metric size 5) except for diameter.
- Inserter Handle: Surgical grade stainless steel and medical grade plastic.
- The product is not made with natural rubber latex.
- The product does not contain DEHP.

HOW SUPPLIED

- The Device is provided **STERILE** for **SINGLE USE ONLY**. The device is only sterile if the packaging is not opened, damaged or broken.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Do not resterilize or reuse the device as this may result in product malfunction, failure or patient injury, and may also expose the patient to the risk of infectious diseases.

MAGNETIC RESONANCE (MR) IMAGING

The implant is MR safe.

STORAGE

Store away from moisture and direct heat.

DISPOSAL

International and US regulations require controlled disposal of used and unused medical devices. Product must be disposed of in accordance with applicable government regulations and Medical Waste Management Plans/Acts.

CUSTOMER SERVICE

Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. DO NOT REUSE. **THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR**

STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.

Product Complaints & Return Authorization

All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized Customer Service representative.

Manufacturer

ArthroCare Corporation

7000 West William Cannon Drive Austin, TX 78735 USA

Phone: (800) 343-5717

www.arthrocare.com


Authorized European Representative

Smith & Nephew Medical Limited

101 Hessle Road Hull HU3 2BN United Kingdom

ec.rep@smith-nephew.com

SYMBOLS LEGEND	
	Catalogue number
	Lot number
	Use by
	MR Safe
	Consult instructions for use
	Do not reuse
	Do not resterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Caution
	Sterilization using irradiation
	Manufacturer
	Authorized representative in the European Community
Rx Only	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.
	CE mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).


° Trademark of Smith & Nephew.

These products may be covered by one or more U.S. patents.

See smith-nephew.com/patents to view a listing of patents that cover our products.

© Smith & Nephew 2017.

P/N 62538 Rev. C April 2017

MULTIFIX® S ULTRA Suture Anchor

Instructions for Use

DESCRIPTION

The MULTIFIX® S ULTRA Knotless Fixation System consists of an implant (anchor) and associated instruments which facilitate the attachment of tissue to bone. The implant, which is available in 5.5mm and 6.5mm sizes, can accommodate 2 to 4 strands of suture thread or a maximum volume of 4 strands of ULTRATAPE® suture. The implant is manufactured from polyetheretherketone (PEEK) and is preloaded onto a disposable handle designed to facilitate direct insertion into bone. The suture anchor is a knotless fixation device, in other words surgical knots are not necessary for fixation of suture to tissue.

INDICATIONS FOR USE

The MULTIFIX S ULTRA Knotless Fixation Device is indicated for use in the fixation of soft tissue to bone. Examples of such procedures include:

Shoulder: Bankart Repair, SLAP lesion repair, acromio-clavicular separation, rotator cuff repair, capsule shift/capsule-labral reconstruction, biceps tenodesis and deltoid repair

Ankle: Lateral instability, medial instability, Achilles tendon repair/reconstruction and midfoot reconstruction

Foot: Hallux valgus reconstruction

Elbow: Tennis elbow repair, biceps tendon reattachment

Knee: Extra-capsular repairs; reattachment of medial collateral ligament, posterior oblique ligament or joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions

CONTRAINDICATIONS

- Surgical procedures other than those listed in the INDICATIONS FOR USE section.
- Pathologic bone conditions such as cystic changes or severe osteopenia, which would impair secure fixation of the implant.
- Pathological changes in the affected soft tissues, which would prevent secure fixation.
- Inadequate bone stock or comminuted bone surface(s), which would compromise secure fixation of the implant.
- Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitation, previous infections, etc.
- Conditions such as senility, mental illness or alcoholism, which may impact the patient's ability or willingness to restrict activities or comply with physician instructions during the healing period.
- The implant is not designed for and should never be used to attach artificial ligaments or other implants.

WARNINGS AND PRECAUTIONS

Warnings

- Do not use if package is opened or damaged.
- Do not resterilize or reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases.
- The implant with inserter handle is a single use device to be used only once for a single patient.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- Do not implant the anchor in poor quality bone or where bone quantity is limited. Incomplete insertion or poor bone quality may result in implant pullout.
- Failure to distribute proper tension through both suture ends may result in inadequate suture lock and potential failure.
- Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.

Precautions

- Implantation of this device requires an appropriate level of surgical skill and experience. Surgeons who plan to use this device should carefully review the instructions for use prior to clinical use. Additionally, completion of a skills laboratory and training by the manufacturer or one of its appointed representatives, or observation of/assistance with similar implantation procedures, is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. It is the responsibility of the surgeon to advise the patient of the known risks and complications associated with the surgical procedure and device.
- Any decision to remove the implant should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- A back-up device should always be available to avoid any delay in surgery should it be required.
- Do not insert an implant into a bone hole where an existing implant is present. This may cause the previously placed implant to extrude out through the bone.

POTENTIAL ADVERSE EFFECTS

As with any surgical procedure, there is risk involved. Potential complications accompanying such implant surgery may include, but are not limited to the following: complications with anesthesia; pain at the incision or surgical site; infection, both deep and superficial; bone damage or fracture; injury to surrounding tissues or vasculature; embolus or blood clotting issues (e.g., pulmonary embolus, deep vein thrombosis, etc.); nerve injury or palsy; implant or treatment failure. Secondary surgical interventions may include, but are not limited to: removal of implant; placement of new or additional implant(s).

OPERATOR TRAINING REQUIREMENTS

- Use of this medical device requires an appropriate level of surgical skill and experience.
- Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/assistance with similar surgical procedures is

recommended.

- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using this medical device.
- The surgeon should advise the patient of the known risks and complications associated with the surgical procedure and use of this device.

PREPARATION FOR USE

Caution: This device has been qualified for use with #2 MAGNUMWIRE®, ULTRABRAID®, and ULTRATAPE sutures. Use of other suture may compromise suture locking integrity and result in potential failure.

Caution: When required use only the approved bone punch or bone tap instruments. Use of other ancillary instrumentation may compromise implant insertion and/or retention.

DIRECTIONS FOR USE

- Suture tissue using a suture passer utilizing the appropriate suture.
 - Using the Insertion Guide with Obturator is recommended. When using the insertion guide in the portal, align the slit of the insertion guide in the direction of the stitch in the tissue.
 - Assess the patient's bone quality prior to use. If in the physician's opinion the bone is hard, create a bone hole using a 4.5mm Tapered Bone Punch prior to anchor insertion. A bone hole is also required if the implant is to be placed at an angle other than perpendicular to the bone surface. Advance the punch until the 25mm mark is flush with the bone surface.
- Caution:** Care must be taken to avoid creation of a shallow bone hole. A proper hole depth is accomplished when the punch is advanced to the 25mm mark.
- Insert the suture ends through the first snare loop. Ensure about 3cm of each suture extends through the snare loop (see diagram 1). **Note:** The first snare loop can accommodate 2 suture threads, 2 suture tapes or 1 suture tape and 1 suture thread.
 - To snare the sutures through the anchor, use a syringe type motion to pull the suture loader cross bar (see diagram 2) toward the handle.
 - If additional strands of suture are to be used, thread those suture ends through the second snare loop. Ensure about 3cm of each suture thread or tape extends through the second snare loop (see diagram 3). **Note:** The volume of suture the second snare loop can accommodate is dependent on the volume of suture already threaded through the first snare loop. If the first snare loop had 2 or less suture threads, do not place more than 2 suture threads into the second loop. If the first snare loop already snared 2 suture tapes do not place more than 2 suture tapes or 2 suture threads through the second loop. The total volume of suture through both snare loops should not exceed 4 strands of suture thread, 4 strands of suture tape or 2 strands of suture tape and 2 strands of suture thread.
 - To snare the sutures through the anchor, pull the cross bar of the suture loader handle away from the shaft of the anchor inserter handle at approximately a 45° angle and pull the suture loader handle away from the anchor inserter. Discard the suture loader after the sutures have been passed through the anchor.
 - Holding the free ends of the suture, advance the loaded anchor either percutaneously or through a cannula to the surface of the bone. Take up suture slack by pulling on the ends of the suture until suture can be clearly seen between the tissue and the implant. Confirm that the suture is not wrapped around the device. If wrapped, pull the device out and repeat this step.
 - Place the pound-in tip of the anchor at the desired location lateral to where the soft tissue is to

be secured. Assure the eyelet of the pound-in tip orients the sutures toward the tissue (see diagram 5). Do not split the sutures with the implant. Establish axial alignment of the inserter shaft perpendicular to the bone. If a bone hole was created, establish axial alignment of the inserter handle with the bone hole.

„Remove desired slack in the sutures by pulling on the ends of the suture.

„Use a mallet to tap the proximal end of the inserter handle to drive the pound-in tip into the bone. Tap the proximal end of the inserter until the laser mark on the leading thread of the screw body is sub-cortical (see diagram 6).

Caution: If the pound-in tip does not penetrate the bone on the first strike, create a bone hole according to step 3. Use a new anchor in the bone hole.

Caution: Use care to properly align the implant and inserter handle with the bone hole during pound-in. Avoid excessive probing. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.

Caution: Incomplete insertion or poor bone quality may result in the implant pulling out.

„Tension will be applied to the tissue during advancement of the pound-in tip. If additional tension is required individually pull the free ends of the suture in line with the axis of the inserter while maintaining downward pressure on the inserter handle. Suture may be placed in the suture cleats for suture management if desired (see diagram 7).

Caution: Individually tensioning the sutures may cause the device to back out of the bone hole. If this occurs, advance the pound-in tip until the laser mark on the leading edge of the screw is sub-cortical before applying tension to the next suture.

Caution: Do not over tension the sutures as breakage or tissue pull through may occur. Additionally, over tensioning the suture may result in incomplete insertion or implant pull out.

„While continuing to maintain downward pressure on the inserter handle, deploy the screw by turning the deployment knob (see diagram 8) clockwise until the screw is flush or below the surface of the bone (see diagram 9) and the deployment knob has come to a hard stop.

„Disengage sutures from the cleats if needed prior to retracting the inserter handle.

„Disengage the inserter handle by turning the end cap counterclockwise a minimum of six turns (see diagram 10).

„Trim remaining loose suture ends at the bone hole.

„After complete deployment or disengagement of inserter from anchor body, if revision/removal of implant is deemed necessary, then proceed with the following steps: advance inserter handle until hex pattern is reengaged with implant. Turn deployment knob counterclockwise to remove the implant screw and pull on the ends of the suture to remove the remaining implant eyelet. If implant is not able to be reengaged with the inserter handle, or for post-operative repair, insert the Anchor Extraction tool and turn counterclockwise to remove the implant screw and then pull on the sutures to remove the remaining implant eyelet.

„If placing more than one anchor, assure that the anchors are located at least 7mm apart and repeat steps 1 through 16.

MATERIALS

• Anchor: PEEK [polyetheretherketone].

• Suture: Implant grade polyethylene and polypropylene. The suture is USP size 2 (metric size 5) except for diameter.

• Inserter Handle: Surgical grade stainless steel and medical grade plastic.

• The product is not made with natural rubber latex.

- The product does not contain DEHP.

HOW SUPPLIED

- The Device is provided **STERILE** for **SINGLE USE ONLY**. The device is only sterile if the packaging is not opened, damaged or broken.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Do not resterilize or reuse the device as this may result in product malfunction, failure or patient injury, and may also expose the patient to the risk of infectious diseases.

MAGNETIC RESONANCE (MR) IMAGING

The implant is MR safe.

STORAGE

Store away from moisture and direct heat.

DISPOSAL

International and US regulations require controlled disposal of used and unused medical devices. Product must be disposed of in accordance with applicable government regulations and Medical Waste Management Plans/Acts.

CUSTOMER SERVICE

Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. **DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

Product Complaints & Return Authorization

All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized Customer Service representative.



Manufacturer

ArthroCare Corporation
7000 West William Cannon Drive Austin, TX 78735 USA
Phone: (800) 343-5717
www.arthrocure.com

Authorized European Representative

Smith & Nephew Medical Limited 101 Hessle Road
Hull HU3 2BN United Kingdom
ec.rep@smith-nephew.com

SYMBOLS LEGEND	
	Catalogue number
	Lot number
	Use by
	MR Safe
	Consult instructions for use
	Do not reuse
	Do not resterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Caution
	Sterilization using irradiation
	Manufacturer
	Authorized representative in the European Community
Rx Only	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.
	CE mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).




° Trademark of Smith & Nephew.

These products may be covered by one or more U.S. patents.

See smith-nephew.com/patents to view a listing of patents that cover our products.

© Smith & Nephew 2017.

P/N 65101 Rev. B April 2017

Q-FIX® Soft Suture Anchor Instructions for Use

DESCRIPTION

The Q-FIX® Soft Suture Anchor System consists of an implant (anchor) and associated instruments which facilitate the attachment of tissue to bone. The implant, which is available in 1.8mm (single loaded) and 2.8mm (double loaded) sizes, is manufactured from braided polyester with ultra high molecular weight polyethylene suture and is pre-loaded into a disposable tool designed to facilitate direct insertion into a pre-drilled/pre-punched bone hole.

INDICATIONS FOR USE

The Q-FIX Soft Suture Anchor System is intended to be used for soft tissue to bone fixation for:

Shoulder: Bankart lesion repair; SLAP lesion repair; acromio-clavicular repair; capsular shift/capsulolabral reconstruction; deltoid repair; rotator cuff tear repair; biceps tenodesis.

Foot & Ankle: Medial/Lateral repair and reconstruction; midfoot and forefoot repair; Hallux valgus reconstruction; Metatarsal ligament/tendon repair or reconstruction; Achilles tendon repair.

Elbow: Ulnar or radial collateral ligament reconstruction; lateral epicondylitis repair; biceps tendon reattachment.

Knee: Extra-capsular repair: medial collateral ligament (MCL), lateral collateral ligament (LCL) and posterior oblique ligament; iliotibial band tenodesis (IBT); patellar tendon repair; vastus medialis obliquus advancement (VMO); joint capsule closure.

Hip: Acetabular labral repair.

CONTRAINDICATIONS

Surgical procedures other than those listed in the INDICATIONS FOR USE section.

- Pathologic bone conditions such as cystic changes or severe osteopenia, which would impair secure fixation of the implant.
- Pathological changes in the affected soft tissues, which would prevent secure fixation.
- Insufficient quality or quantity of bone, inadequate bone stock or comminuted bone surface(s) which would compromise secure fixation of the implant.
- The use of this device may not be suitable for patients with insufficient or immature bone. The physician should carefully assess bone quality before performing orthopedic surgery on patients who are skeletally immature. The use of this device and the placement of implants must not bridge, disturb or disrupt the growth plate.
- Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitations, active or previous infections, etc.
- Foreign body sensitivity to implant materials. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Conditions such as senility, mental illness or alcoholism, which may impact the patient's ability or willingness to restrict activities or comply with physician instructions during the healing period.
- The implant is not designed for and should never be used to attach artificial ligaments or other implants.
- Use of suture other than ArthroCare® USP #2 MAGNUMWIRE suture.

WARNINGS AND PRECAUTIONS **Warnings**

- Do not use if package is opened or damaged.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- Do not resterilize or reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases. The implant is a single use device to be used only once for a single patient.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Any decision to remove the device should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- Do not implant the anchor in poor quality bone or where bone quantity is limited. Incomplete insertion or poor bone quality may result in implant pullout or suture breakage.
- Do not bend, apply excessive torque, or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.
- Do not attempt to drill an anchor out of the bone hole. Should a placed anchor need to be abandoned, trim the suture ends flush with the bone to avoid potential tissue irritation.
- Do not insert an implant into a bone hole where an existing implant is present. This may cause the previously placed implant to extrude out of the bone.
- Do not use the inserter to probe the tissue or bone as damage may occur to the device resulting in potential patient injury.

Precautions

- A back-up device should always be available to avoid any delay in surgery should it be required.

POTENTIAL ADVERSE EFFECTS

As with any surgical procedure, there is risk involved. Potential complications accompanying such implant surgery may include, but are not limited to the following: complications with anesthesia; pain at the incision or surgical site; infection, both deep and superficial; bone damage or fracture; injury to surrounding tissues or vasculature; embolus or blood clotting issues (e.g., pulmonary embolus, deep vein thrombosis, etc.); nerve injury or palsy; treatment or implant failure; and secondary surgical intervention to address complications associated with surgery or treatment.

OPERATOR TRAINING REQUIREMENTS

- Use of the Q-FIX Soft Suture Anchor System requires an appropriate level of surgical skill and experience.
- Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/assistance with similar surgical procedures is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using the Q-FIX Soft Suture Anchor System.
- The surgeon should advise the patient of the known risks and complications associated with the surgical procedure and use of this device.

PREPARATION FOR USE

Caution: Use only ArthroCare drill bits designed for use with these anchors. Use of other drill bits

will compromise implant retention integrity.

Caution: Use only the ArthroCare Reusable Bone Punch designed for use with these anchors. Use of other punches may compromise implant insertion and/or retention.

DIRECTIONS FOR USE

- Using standard surgical technique, expose the bone at the desired insertion location for the implant.
- Attach the appropriate size drill bit to a standard orthopedic drill. Be sure to chuck the drill no farther than the flats on the proximal end of the drill.
- Position the tip of the Drill Guide with Obturator against the bone surface at the desired position (see figures, Item 1), then remove the Obturator (see figures, Item 2).

Caution: Do not advance the Drill Guide into the bone.

Doing so may compromise the cortical surface of the bone and reduce implant retention capability.

- Insert the drill into the drill guide and drill the bone hole. Advance the drill until the shoulder on the drill bottoms in the drill guide (see figures, Item 3). Alternatively, insert the reusable punch into the drill guide and drive the punch forward until the distal edge of the handle is aligned with the proximal edge of the drill guide (see figures, Item 4).

Caution: Do not use the drill or punch without the drill guide as the bone surface can be damaged.

Caution: Care must be taken to not create too deep or too shallow a bone hole. To accomplish this, be sure to drill the hole until the shoulder on the drill bottoms in the drill guide. If using the punch, insure the handle lines up with the proximal edge of the drill guide. Care must be taken to not advance the Drill Guide into the bone as the implant retention force may be compromised.

Caution: Confirm that bone hole is clear of material by advancing and retracting the drill until the flutes of the drill are visually clear of bone fragments upon withdrawal.

Failure to clear material from the bone hole may result in difficulty inserting the implant which could compromise retention force.

- Remove the drill or punch. A PATHFINDER® may be used to confirm bone hole location and depth. **Caution:** Do not attempt to insert the implant into the bone without first creating a bone hole.
- Insert the Q-FIX implant through the drill guide into the bone hole (see figures, Item 5). Be sure that the implant inserter is fully engaged with the drill guide throughout the deployment of the anchor (see figures, Item 6). This can be accomplished by holding the drill guide and inserter together by hand or by applying downward pressure to maintain the engagement during deployment.

Also, be sure to maintain the alignment of the drill guide to the bone hole during insertion.

Caution: Do not use excessive force or pound the inserter into the bone hole as damage to the inserter/implant may occur or result in injury to the patient. If resistance is met on insertion, turn the inserter handle clockwise while keeping a loose hold on the drill guide. If continued resistance is met, remove the inserter and re-drill the existing bone hole.

Caution: Use care to properly align the implant and inserter handle with the bone hole during insertion. A PATHFINDER may be used to confirm bone hole location and depth. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy bent or damaged inserter.

Caution: Be sure to advance the inserter fully into the bone hole. Incomplete insertion or poor bone quality may result in implant pullout.

- Rotate the activation knob clockwise on the proximal end of the implant inserter to deploy the implant (see figures, Item 7). Rotate the knob until a hard stop is reached and the suture cleat is exposed at the proximal end of the inserter (see figures, Item 8).
- Release the sutures from the cleat (see figures, Item 9) and remove the inserter and drill guide. Pull on the suture tails to ensure that the implant is set in the bone hole.
- If placing more than one implant, assure that the anchors are located at least 7 mm apart and repeat steps 1 through 9.

MATERIALS

Anchor implant: Braided Polyester.

Suture: Ultra High Molecular Weight Polyethylene.

Inserter Handle: Surgical grade stainless steel and medical grade plastic.

The product is not made with natural rubber latex.

The product does not contain DEHP.

HOW SUPPLIED

The Device is provided **STERILE for SINGLE USE ONLY**. The device is only sterile if the packaging is not opened, damaged or broken.

Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.

Do not resterilize or reuse the device as this may result in product malfunction, failure or patient injury, and may also expose the patient to the risk of infectious diseases.

The device is sterilized using ethylene oxide.

MAGNETIC RESONANCE (MR) IMAGING

The implant is MR safe.

STORAGE

Store away from moisture and direct heat.

DISPOSAL

International and US regulations require controlled disposal of used and unused medical devices. Product must be disposed of in accordance with applicable government regulations and Medical Waste Management Plans/Acts.

CUSTOMER SERVICE

Warranty information

This product is guaranteed for materials, function, and workmanship for single patient use only.

DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.

Product Complaints & Return Authorization

All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized Customer Service representative.

**Manufacturer**

ArthroCare Corporation

7000 West William Cannon Drive Austin, TX 78735 USA

Phone: (800) 343-5717

www.arthrocare.com**Authorized European Representative**

Smith & Nephew Medical Limited 101 Hessle Road

Hull HU3 2BN United Kingdom

ec.rep@smith-nephew.com

SYMBOLS LEGEND	
	Catalogue number
	Lot number
	Use-by date
	MR Safe
	Consult instructions for use
	Do not re-use
	Do not resterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Caution
	Sterilization using ethylene oxide
	Manufacturer
	Authorized representative in the European Community
	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.
	CE mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).

°Trademark of Smith & Nephew.

This product may be covered by one or more US Patents.

See smith-nephew.com/patents for details.

© Smith & Nephew 2017.

April 2017 P/N 62526 Rev. B


Diagram 1 (Figure 1-2)


Diagram 2 (Figure 3)


Diagram 3 (Figure 4)


Diagram 4 (Figure 5-7)


Diagram 5 (Figure 8-9)

SPEEDLOCK® KNOTLESS FIXATION DEVICE

INSTRUCTIONS FOR USE

DESCRIPTION

The SPEEDLOCK® Knotless Fixation Device is an implant that facilitates the attachment of tissue to bone. The SPEEDLOCK system consists of a fixation device pre-loaded with a suture snare in an inserter handle. The SPEEDLOCK Implant is a knotless fixation device, In other words surgical knots are not necessary for fixation of suture to tissue.

INDICATIONS FOR USE

The SPEEDLOCK Knotless Fixation Device is indicated for use in the fixation of soft tissue to bone. Examples of such procedures include:

Shoulder: Bankart Repair, SLAP lesion repair, acromio-clavicular separation, rotator cuff repair, capsule shift/ capsule-labral reconstruction, biceps tenodesis and deltoid repair.

Ankle: Lateral instability, medial instability, Achilles tendon repair/reconstruction and midfoot reconstruction.

Foot: Hallux valgus reconstruction.

Elbow: Tennis elbow repair, biceps tendon reattachment.

Knee: Extra-capsular repairs; reattachment of medial collateral ligament, posterior oblique ligament or joint capsule closure to anterior proximal tibia; extra capsular reconstruction, ITB tenodesis; patellar ligament and tendon avulsions.

CONTRAINDICATIONS

- Surgical procedures other than those listed in the INDICATIONS section.
- Pathologic bone conditions such as cystic changes or severe osteopenia, which would impair secure fixation of the implant.
- Pathological changes in the affected soft tissues, which would prevent secure fixation.
- Inadequate bone stock or comminuted bone surface, which would compromise secure fixation of the implant.
- Physical conditions that would eliminate or tend to eliminate adequate implant support or retard healing, i.e., blood supply limitation, previous infections, etc.
- Conditions such as senility, mental illness or alcoholism, which may impact the patient's ability or willingness to comply with physician instructions during the healing period.
- The implant is not designed for and should never be used to attach artificial ligaments or other implants.
- Use of suture other than ArthroCare USP #2 (metric size 5) MAGNUMWIRE suture.

WARNINGS AND PRECAUTIONS

Warnings

- Do not use if package is opened or damaged.
- Do not resterilize or reuse as this may result in product malfunction, failure, or patient injury and may also expose the patient to the risk of infectious diseases.
- The implant with inserter handle is a single use device to be used only once for a single patient.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Detailed instructions on the use and limitations of the device should be given to the patient.
- Do not implant the anchor in poor quality bone or where bone quantity is limited. Incomplete

insertion or poor bone quality may result in implant pullout.

- Failure to distribute proper tension through both suture ends may result in inadequate suture lock and potential failure.
- Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.
- The SPEEDLOCK implant system requires tension to be distributed through both suture ends to achieve suture lock. The use of a locking suture stitch (i.e. modified Mason- Allen) may cause tensioning to be more difficult due to the inability of the suture ends to move freely through the tissue.

Precautions

- Implantation of this device requires an appropriate level of surgical skill and experience. Surgeons who plan to use this device should carefully review the instructions for use prior to clinical use. Additionally, completion of a skills laboratory and training by the manufacturer or one of its appointed representatives, or observation of/assistance with similar implantation procedures, is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques and proper selection and placement of the implant, are important considerations in the successful utilization of this device. It is the responsibility of the surgeon to advise the patient of the known risks and complications associated with the surgical procedure and device.
- Any decision to remove the implant should take into consideration the potential risk to the patient of a second surgical procedure. Implant removal should be followed by adequate postoperative management.
- A back-up device should always be available to avoid any delay in surgery should it be required.
- Do not insert an implant into a bone hole where an existing implant is present. This may cause the previously placed implant to extrude out through the bone.

POTENTIAL ADVERSE EFFECTS

As with any surgical procedure, there is risk involved. Potential complications accompanying such implant surgery may include, but are not limited to the following: complications with anesthesia; pain at the incision or surgical site; infection, both deep and superficial; bone damage or fracture; injury to surrounding tissues or vasculature; embolus or blood clotting issues (e.g., pulmonary embolus, deep vein thrombosis, etc.); nerve injury or palsy; implant or treatment failure. Secondary surgical interventions may include, but are not limited to: removal of implant; placement of new or additional implant(s).

OPERATOR TRAINING REQUIREMENTS

- Use of this medical device requires an appropriate level of surgical skill and experience.
- Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/assistance with similar surgical procedures is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using this medical device.
- The surgeon should advise the patient of the known risks and complications associated with the surgical procedure and use of this device.

PREPARATION FOR USE

Do not use sutures other than ArthroCare #2 MAGNUMWIRE suture.

Use of other suture may compromise suture locking integrity and result in potential failure.

DIRECTIONS FOR USE

- Suture tissue using #2 MAGNUMWIRE suture. Ensure both ends of the suture are of equal length.
- Insert the drill guide and obturator at the desired bone hole site.
- Remove obturator and drill bone site using a compatible 3.0mm or 3.0mm+ drill.
Caution! Care must be taken to not drill the bone hole too deep or too shallow. To accomplish this, be sure to drill the hole until the proximal edge of the depth mark on the drill is fully subcortical.

- Insert PATHFINDER® instrument to maintain alignment with bone hole.
- Thread the two suture ends through the proximal snare wire loop (Item 1). Pull the distal snare (Item 2).
- Ensure about 6 cm of each suture end extends past the proximal snare (Item 3).
- Rotate suture knob (Item 4) to pull suture ends through the SPEEDLOCK implant, until suture is visible around the suture reel (Item 5). If required, secure the suture ends onto a fixation slit on the suture reel (Item 5).

- Align slit on the drill guide in the direction of the stitch in the tissue.
- Remove PATHFINDER instrument prior to inserting SPEEDLOCK implant (Item 6).
- Insert the SPEEDLOCK implant into the drill guide oriented such that the suture leading from the implant is aligned with the slit on the drill guide (Item 7). Advance the implant, continuing to take up the suture slack, until it is at the opening of the drilled hole.

- Confirm alignment of the implant such that the slot on the plug guide (Item 8) orients the suture towards the tissue (Item 9). If not aligned, reposition prior to implant insertion into the drilled hole. Insert/drive the implant into the drilled hole to the laser line (Item 10) on the plug guide.

Note: Check that the entire length of the suture is visible, from the tissue to the implant. Do not split the sutures with the implant.

Caution! Failure to align the suture exiting the implant to the stitch may result in a loose repair or suture breakage.

Caution! Use care to properly align the implant and inserter handle with the bone hole during insertion. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy bent or damaged implant.

Caution! Ensure the suture does not wrap around the inserter or drill guide when inserting the implant into bone.

Caution! Incomplete insertion or poor bone quality may result in implant pullout.

Caution! Inserting the implant beyond the laser line may result in difficulty in removing the inserter.

- Rotate the suture knob to take up the slack in the suture and tension the tissue to the desired position (Item 11). If suture tension needs to be released, push the corresponding ratchet release button.

Caution! Do not over-tension the suture.

- When desired tension is achieved, deploy the suture lock and detach the handle by turning the activation knob (Item 12) in the direction of the arrow approximately 1.5 complete revolutions until it stops.
- Remove handle. Trim remaining loose suture ends at drilled hole.
- If placing more than one implant, assure drilled holes are at least 5 mm apart and repeat steps

1 through 14.

MATERIALS

- Implant: PEEK [polyetheretherketone].
- Suture: Implant grade polyethylene and polypropylene.
- Inserter Handle: Surgical grade stainless steel and medical grade plastic.
- The product is not made with natural rubber latex.
- The product does not contain DEHP.

HOW SUPPLIED

- The Device is provided **STERILE** for **SINGLE USE ONLY**. The device is only sterile if the packaging is not opened, damaged or broken.
- Do not use the product after the "use by" date printed on the label. Device performance and patient safety may be compromised if the product is used after its expiration date.
- Do not resterilize or reuse the device as this may result in product malfunction, failure or patient injury, and may also expose the patient to the risk of infectious diseases.

MAGNETIC RESONANCE (MR) IMAGING

The implant is MR safe.

STORAGE

Store away from moisture and direct heat.

DISPOSAL

International and US regulations require controlled disposal of used and unused medical devices. Product must be disposed of in accordance with applicable government regulations and Medical Waste Management Plans/Acts.

Warranty Information

This product is guaranteed for materials, function, and workmanship for single patient use only. **DO NOT REUSE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESS, IMPLIED, AND/OR STATUTORY, INCLUDING BUT NOT LIMITED TO, WARRANTIES OF MERCHANTABILITY, FITNESS, AND/OR SUITABILITY FOR A PARTICULAR PURPOSE, AND OF ALL OTHER OBLIGATIONS OR LIABILITIES ON THE PART OF ARTHROCARE.**

Product Complaints & Return Authorization

All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized Customer Service representative.

Manufacturer


ArthroCare Corporation

7000 West William Cannon Drive Austin, TX 78735 USA

Phone: (800) 343-5717





www.arthrocare.com

EC REP Authorized European Representative
 Smith & Nephew Medical Limited 101 Hessle Road
 Hull HU3 2BN United Kingdom
ec.rep@smith-nephew.com

SYMBOLS LEGEND	
	Catalogue number
	Lot number
	Use-by date
	MR Safe
	Consult instructions for use
	Do not re-use
	Do not re-sterilize
	Do not use if package is damaged
	Keep away from sunlight
	Keep dry
	Caution
	Sterilization using ethylene oxide
	Manufacturer
	Authorized representative in the European Community
Rx Only	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician.
	CE mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC).

°Trademark of Smith & Nephew.
 This product may be covered by one or more US Patents.
 See smith-nephew.com/patents for details.
 © Smith & Nephew 2017.

March 2017 PN 62513 Rev. B

**Addendum to the instruction for use
for Medical Devices (MD)**
Non-absorbable arthroscopic suture anchors (fixation implants)

When delivering a product to the Russian Federation, the information contained in this addendum to the instruction for use is a priority.

1. Purpose of the medical device

1.1. Name of the medical device

Non-absorbable arthroscopic suture anchors (may be referred to as Suture anchors, medical devices, devices)

Please consider the names of the variants of execution, specified in the application instructions, identical names, presented in the table below:

Name in the IFU	Identical name
The MULTIFIX P Knotless Fixation Device	MultiFIX P 4.5 mm fixation implant
The MULTIFIX S Knotless Fixation Device	MultiFIX S 5.5 mm fixation implant MultiFIX S 6.5 mm fixation implant
The MULTIFIX S Knotless Fixation System	
MULTIFIX S ULTRA Suture Anchor The MULTIFIX S ULTRA Knotless Fixation System	MultiFIX S ULTRA 5.5 mm fixation implant MultiFIX S ULTRA 6.5 mm fixation implant
Q-FIX Soft Suture Anchor The Q-FIX Soft Suture Anchor System	Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE # 2 suture (thread) (white/blue co-braid) Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)
The SPEEDLOCK system	SPEEDLOCK fixation implant

1.2. The purpose of the medical device

Non-absorbable arthroscopic suture anchors are designed for fastening soft tissue to bones during arthroscopic surgery.

1.3. Operating conditions

The product is designed for use in the internal environment of the body and in the environment in the operating system, presented below:

Environmental requirements upon using this device in patient's body	
Temperature	+32°C to + 42°C
Relative humidity	100%
Atmospheric pressure	86 kPa - 106 kPa

Environmental requirements for surgical intervention	
Temperature	+15°C to + 25°C
Relative humidity	30% - 75%
Atmospheric pressure	86 kPa - 106 kPa

1.4. Conditions for medical device use in medical organizations, mobile transport media, field environment or home conditions, public or individual application and potential device users

These medical devices should be used by specially trained personnel at medical and preventive treatment institutions.

This Medical Device (MD) is intended for use by trained medical personnel **only**:

- Use of this medical device requires an appropriate level of surgical skill and experience.
- Completion of a skills laboratory, training by the manufacturer or one of its appointed representatives, and/or observation of/assistance with similar surgical procedures is recommended.
- Preoperative and operative procedures, including proper patient selection, knowledge of surgical techniques, and proper device selection are important considerations when using this medical device.
- The surgeon should advise the patient of the known risks and complications associated with the surgical procedure and use of this device.

This medical device has no antigenic properties and is apyrogenic.

1.5. Indication of possibilities and peculiarities of medical device application for patients with implanted medical devices, pregnant women, women during lactation and breast feeding periods, children and adults with chronic diseases

Use of this implant is allowed upon medical prescription and is forbidden in the event of contraindications. This device may be used for patients of all age, sex and ethnicity.

2. Information on the medical device manufacturer

2.1. Name

ArthroCare Corporation

2.2. Legal address

7000 West William Cannon Drive Austin, TX 78735 USA

2.3. Manufacturing site address

1. ArthroCare Corporation

7000 West William Cannon Drive, Austin, TX 78735 USA

2. ArthroCare Corporation

15285 Alton Parkway, Suite 200, Irvine, CA 92618, USA

3. ArthroCare Corporation

B32.1, St2, Zona Franca Coyol, Coyol, Costa Rica, 20101, Costa Rica

2.4. Information on the authorized representative of the medical device manufacturer
LLC "Smith and Nephew"
Russia, 105120, Moscow, 2nd Syromiatnicheskiy lane, building 1, 9th floor, room1, office 1.

3. Medical device classification

Class under potential risk of medical device application	2B
Electrical class	Not applicable. Inactive medical device
Device class under type and duration of contact	constant contact with internal body environment
Frequency of use	single-use medical device
Invasiveness	Invasive medical device, surgically installed
Sterility	Sterile medical device

4. Description of the principles underlying operation of the medical device and its features

This device is implanted into the bone and attaches soft tissues by means of pre-installed or intraoperatively installed sutures. An aperture is made in the bone prior to this device implanting. Polyetheretherketone (PEEK) implants will be secure after the bone tissue regenerates. Implants do not require further retrieval, except in the case described below:

- Loosening the product in the bone canal in violation of the rules of rehabilitation of the duration of the immobilization of the joint appointed by the doctor, depending on the type of damage;
- Incorrect installation of the product (deep, angle of introduction) in case of non-compliance with the rules of installation of the product.

5. Instructions for the safe removal of MI from the patient's body

Used an MRI scanner for the diagnosis and the localization of the shoulder implant on the humerus or joint trough of the shoulder blade.

Using a surgical scalpel, a skin incision of about 5 cm is performed, then the same instrument makes a layered incision of fascia until the moment of visual detection of the strands of the shoulder implant and the channel with the shoulder retainer installed. The strands installed in the shoulder retainer are cut off by a surgical scalpel. A hollow drill/trepan is used to extract the shoulder fixer, the diameter of which is 2-3 mm larger than the diameter of the installed shoulder retainer. The hollow drill/trepan is centered against the shoulder retainer, the bone block containing the shoulder retainer is drilled. This bone block is then removed from the patient's body.

If necessary, another part of the bone is undergoing an audit operation to refix soft tissues.

A control MRI examination is being carried out.

6. Method of application

MultIFIX P 4.5 mm fixation implant

PREPARATION FOR USE

Caution! This device has been qualified for use with #2 MAGNUMWIRE, ULTRABRAID, and ULTRATAPE sutures. Use of other suture is not recommended and may compromise suture locking integrity and result in potential failure.

Caution! When required, use only the approved drill, bone punch or bone tap instruments. Use of other ancillary instrumentation may compromise Implant insertion and/or retention.


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5


Figure 6



Figure 7

DIRECTIONS FOR USE

- 1) Suture tissue using a suture passer utilizing the appropriate suture.
- 2) Using the Insertion Guide with Obturator is recommended. When using the insertion guide in the portal, align the slit of the insertion guide in the direction of the stitch in the tissue.
- 3) Assess the patient's bone quality prior to use. If in the physician's opinion the bone is hard, create a bone hole using drill bit or punch (OM-9920) prior to anchor insertion. Advance the drill or the punch until the shoulder reaches the bone surface. Alternately, a 4.5mm tapered bone punch advanced to the 15mm depth marker may be used.

Caution!

Care must be taken to avoid creation of a shallow bone hole. A proper hole depth is accomplished when the 4.5mm tapered punch is advanced to the 15mm mark.

- 4) Thread the suture ends through the snare loop. Ensure about 6 cm of each suture end extends past the snare (see figure 1).

Note: The device is intended to be used with 2, 3 or 4 suture strands or a maximum volume of 2 strands of suture and 2 strands of ULTRATAPE.

- 5) To snare the sutures through the anchor, use a syringe type motion to pull the suture loader cross bar (see figure 2) away from the shaft and toward of the handle at approximately a 45° angle. Discard the suture loader after the sutures have been passed through the anchor.

- 6) Holding the free ends of the suture, advance the loaded anchor either percutaneously or through a cannula to the surface of the bone. Take up suture slack by pulling on the suture tails until suture can be clearly seen between the tissue and the implant. Confirm that the suture is not wrapped around the device. If wrapped, pull the device out and repeat this step.

- 7) Place the tip of the anchor at the desired location lateral to where the tendon is to be secured. Establish axial alignment of the inserter shaft perpendicular to the bone during insertion. Assure the eyelet of the anchor tip orients the sutures toward the tissue (see figure 3). Do not split the sutures with the implant. If a bone hole was created, establish axial alignment of the inserter handle with the bone hole.

- 8) If it is required, tighten threads.

- 9) Use a mallet to tap the proximal end of the inserter handle to drive the anchor tip until the horizontal depth mark (see figure 4) is flush or below the surface of the bone.

Caution! If the anchor tip does not penetrate the bone on the first strike, create a bone hole according to step 3.

Use a new anchor into the bone hole.

Caution! Use care to properly align the implant and inserter handle with the bone hole during pound-in. Avoid excessive probing. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.

Caution! Incomplete insertion or poor bone quality may result in the implant pulling out.

- 10) In order to maneuver the tissue to the desired location, pull the free suture tails in line with the axis of the inserter while maintaining downward pressure on the inserter handle. When desired tension is achieved, secure the sutures in the jam cleat (see figure 5).

Caution! Do not over tension the sutures as breakage or tissue pull through may occur.

- 11) While maintaining downward pressure on the inserter handle, deploy the suture lock by turning

the deployment knob (see figure 6) two full turns clockwise. Continue to turn the deployment knob clockwise until it spins freely (see figure 7).

12) Prior to inserter tool handle extraction, release suture material from the clamping brackets. Cut off loose suture ends at near the bone hole.

13) If implant extraction is indicated at this point, turn deployment knob counterclockwise to remove the implant screw and pull on the suture tails to remove the implant. If implant has detached from inserter handle, or for post-operative repair, insert the Anchor Extraction tool and turn counterclockwise to remove the implant screw and then pull on the sutures to remove the Implant.

14) Disengage sutures from the jam cleats prior to extracting the inserter handle.

Trim remaining loose suture ends at bone hole.

Caution! Do not pull on loose suture ends. This may result in implant pullout. If there is a need to check fixation, probe the repaired tissue.

15) In the event of multiple implants installation the anchors should be placed at not less than 7 mm intervals between each other. Repeat steps 1 – 14 of the above procedure.

MultIFIX S 5.5 mm fixation implant, MultIFIX S 6.5 mm fixation implant

PREPARATION FOR USE

Caution! This device has been qualified for use with #2 MAGNUMWIRE, ULTRABRAID, and ULTRATAPE sutures. Use of other suture is not recommended and may compromise suture locking integrity and result in potential failure.

Caution! When required, use only the approved drill, bone punch or bone tap instruments. Use of other ancillary instrumentation may compromise implant insertion and/or retention.


Figure 1


Figure 2


Figure 3


Figure 4

Figure 5


Figure 6


Figure 7


Figure 8

DIRECTIONS FOR USE

- 1) Suture tissue using a suture passer utilizing the appropriate suture.
 - 2) Using the Insertion Guide with Obturator is recommended. When using the insertion guide in the portal, align the slit of the insertion guide in the direction of the stitch in the tissue.
 - 3) Assess the patient's bone quality prior to use. If in the physician's opinion the bone is hard, create a bone hole using a 4.5mm Tapered Bone Punch (OM-9220) prior to anchor insertion. A bone hole is also required if the implant is to be placed at an angle other than perpendicular to the bone surface. Advance the punch until the 25 mm mark is flush with the bone surface.
Caution! Care must be taken to avoid creation of a shallow bone hole. A proper hole depth is accomplished when the punch is advanced to the 25mm mark.
 - 4) Insert the suture ends through the first snare loop. Ensure about 3cm of each suture extends through the snare loop (see figure 1).
- Note:** Note: The device is intended to be used with 2, 3 or 4 suture strands or a maximum volume of 2 strands of suture and 2 strands of ULTRATAPE.

- 5) To snare the sutures through the anchor, use a syringe type motion to pull the suture loader cross bar (see figure 2) away from the shaft and toward of the handle at approximately a 45° angle. Discard the suture loader after the sutures have been passed through the anchor.
- 6) Holding the free ends of the suture, advance the loaded anchor either percutaneously or through a cannula to the surface of the bone. Take up suture slack by pulling on the ends of the suture until suture can be clearly seen between the tissue and the implant. Confirm that the suture is not wrapped around the device. If wrapped, pull the device out and repeat this step.
- 7) Place the pound-in tip of the anchor at the desired location lateral to where the soft tissue is to be secured. Assure the eyelet of the pound-in tip orients the sutures toward the tissue (see figure 3). Do not split the sutures with the implant. Establish axial alignment of the inserter shaft perpendicular to the bone. If a bone hole was created, establish axial alignment of the inserter handle with the bone

hole.

8) Remove desired slack in the sutures by pulling on the ends of the suture.

9) Use a mallet to tap the proximal end of the inserter handle to drive the pound-in tip into the bone. Tap the proximal end of the inserter until the laser mark on the leading thread of the screw body is sub-cortical (see figure 4).

Caution! If the pound-in tip does not penetrate the bone on the first strike, create a bone hole according to step 3. Use a new anchor in the bone hole.

Caution! Use care to properly align the implant and inserter handle with the bone hole during pound-in. Avoid excessive probing. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.

Caution! Incomplete insertion or poor bone quality may result in the implant pulling out.

10) Tension will be applied to the tissue during advancement of the pound-in tip. If additional tension is required individually pull the free ends of the suture in line with the axis of the inserter while maintaining downward pressure on the inserter handle. Suture may be placed in the suture cleats for suture management if desired (see figure 5).

Caution! Individually tensioning the sutures may cause the device to back out of the bone hole. If this occurs, advance the pound-in tip until the laser mark on the leading edge of the screw is sub-cortical before applying tension to the next suture.

Caution! Do not over tension the sutures as breakage or tissue pull through may occur. Additionally, over tensioning the suture may result in incomplete insertion or implant pull out.

11) While continuing to maintain downward pressure on the inserter handle, deploy the screw by turning the deployment knob (see figure 6) clockwise until the screw is flush or below the surface of the bone (see figure 7) and the deployment knob has come to a hard stop.

12) Disengage sutures from the cleats if needed prior to retracting the inserter handle and disengage the inserter handle by turning the end cap counterclockwise a minimum of six turns (see figure 8).

13) After complete deployment or disengagement of inserter from anchor body, if revision/removal of implant is deemed necessary, then proceed with the following steps: advance inserter handle until hex pattern is reengaged with implant. Turn deployment knob counterclockwise to remove the implant screw and pull on the ends of the suture to remove the remaining implant eyelet. If implant is not able to be reengaged with the inserter handle, or for post-operative repair, insert the Anchor Extraction tool and turn counterclockwise to remove the implant screw and then pull on the sutures to remove the remaining implant eyelet.

14) Trim remaining loose suture ends at the bone hole.

15) If placing more than one anchor, assure that the anchors are located at least 7mm apart and repeat steps 1 through 14.

Multifix S ULTRA 5.5 mm fixation implant, Multifix S ULTRA 6.5 mm fixation implant

PREPARATION FOR USE

Caution! This device has been qualified for use with #2 MAGNUMWIRE, ULTRABRAID, and ULTRATAPE sutures. Use of other suture is not recommended and may compromise suture locking integrity and result in potential failure.

Caution! When required, use only the approved drill, bone punch or bone tap instruments. Use of

other ancillary instrumentation may compromise implant insertion and/or retention.


Figure 1

Figure 2


Figure 3


Figure 4

Figure 5


Figure 6


Figure 7


Figure 8

Figure 9


Figure 10

DIRECTIONS FOR USE

- 1) Suture tissue using a suture passer utilizing the appropriate suture.
- 2) Using the Insertion Guide with Obturator is recommended. When using the insertion guide in the portal, align the slit of the insertion guide in the direction of the stitch in the tissue.
- 3) Assess the patient's bone quality prior to use. If in the physician's opinion the bone is hard, create a bone hole using a 4.5mm Tapered Bone Punch (OM-9220) prior to anchor insertion. A bone hole is also required if the implant is to be placed at an angle other than perpendicular to the bone surface.

Advance the punch until the 25mm mark is flush with the bone surface.

Caution! Care must be taken to avoid creation of a shallow bone hole. A proper hole depth is accomplished when the punch is advanced to the 25 mm mark.

4) Insert the suture ends through the first snare loop. Ensure about 3cm of each suture extends through the snare loop (see figure 1).

Note: The first snare loop can accommodate 2 suture threads, 2 suture tapes or 1 suture tape and 1 suture thread.

Note: The volume of suture the second snare loop can accommodate is dependent on the volume of suture already threaded through the first snare loop. If the first snare loop had 2 or less suture threads, do not place more than 2 suture threads into the second loop. If the first snare loop already snared 2 suture tapes do not place more than 2 suture tapes or 2 suture threads through the second loop. The total volume of suture through both snare loops should not exceed 4 strands of suture thread, 4 strands of suture tape or 2 strands of suture tape and 2 strands of suture thread.

5) To snare the sutures through the anchor, use a syringe type motion to pull the suture loader cross bar (see figure 2) away from the shaft and toward of the handle at approximately a 45° angle. Discard the suture loader after the sutures have been passed through the anchor.

6) Holding the free ends of the suture, advance the loaded anchor either percutaneously or through a cannula to the surface of the bone. Take up suture slack by pulling on the ends of the suture until suture can be clearly seen between the tissue and the implant. Confirm that the suture is not wrapped around the device. If wrapped, pull the device out and repeat this step.

7) Place the pound-in tip of the anchor at the desired location lateral to where the soft tissue is to be secured. Assure the eyelet of the pound-in tip orients the sutures toward the tissue (see figure 4). Do not split the sutures with the implant. Establish axial alignment of the inserter shaft perpendicular to the bone. If a bone hole was created, establish axial alignment of the inserter handle with the bone hole.

8) Remove desired slack in the sutures by pulling on the ends of the suture.

9) Use a mallet to tap the proximal end of the inserter handle to drive the pound-in tip into the bone. Tap the proximal end of the inserter until the laser mark on the leading thread of the screw body is sub-cortical (see figure 5).

Caution! If the pound-in tip does not penetrate the bone on the first strike, create a bone hole according to step 3. Use a new anchor in the bone hole.

Caution! Use care to properly align the implant and inserter handle with the bone hole during pound-in. Avoid excessive probing. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy a bent or damaged implant.

Caution! Incomplete insertion or poor bone quality may result in the implant pulling out.

10) Tension will be applied to the tissue during advancement of the pound-in tip. If additional tension is required individually pull the free ends of the suture in line with the axis of the inserter while maintaining downward pressure on the inserter handle. Suture may be placed in the suture cleats for suture management if desired (see figure 7).

Caution! Individually tensioning the sutures may cause the device to back out of the bone hole. If this occurs, advance the pound-in tip until the laser mark on the leading edge of the screw is sub-cortical before applying tension to the next suture.

Caution! Do not over tension the sutures as breakage or tissue pull through may occur. Additionally, over tensioning the suture may result in incomplete insertion or implant pull out.

- 11) While continuing to maintain downward pressure on the inserter handle, deploy the screw by turning the deployment knob (see figure 8) clockwise until the screw is flush or below the surface of the bone (see figure 9) and the deployment knob has come to a hard stop.
- 12) Disengage sutures from the cleats if needed prior to retracting the inserter handle and disengage the inserter handle by turning the end cap counterclockwise a minimum of six turns (see figure 10).
- 13) After complete deployment or disengagement of inserter from anchor body, if revision/removal of implant is deemed necessary, then proceed with the following steps: advance inserter handle until hex pattern is reengaged with implant. Turn deployment knob counterclockwise to remove the implant screw and pull on the ends of the suture to remove the remaining implant eyelet. If implant is not able to be reengaged with the inserter handle, or for post-operative repair, insert the Anchor Extraction tool and turn counterclockwise to remove the implant screw and then pull on the sutures to remove the remaining implant eyelet.
- 14) Trim remaining loose suture ends at the bone hole.
- 15) If placing more than one anchor, assure that the anchors are located at least 7mm apart and repeat steps 1 through 14.

SPEEDLOCK Fixation Implant

DIRECTIONS FOR USE

Caution! Do not use sutures other than ArthroCare #2 MAGNUMWIRE suture.

Use of other suture may compromise suture locking integrity and result in potential failure.


Figure 1 (Item 1-6)


Figure 2 (Item 7-11)


Figure 3 (item 8)


Figure 4 (item 12)

PREPARATION FOR USE

1. Suture tissue using #2 MAGNUMWIRE suture. Ensure both ends of the suture are of equal length.
2. Insert the drill guide and obturator at the desired bone hole site.
3. Remove obturator and drill bone site using a compatible 3.0mm drill.

Caution!

Care must be taken to not drill the bone hole too deep or too shallow. To accomplish this, be sure to drill the hole until the proximal edge of the depth mark on the drill is fully subcortical.

4. Insert PATHFINDER instrument to maintain alignment with bone hole.
5. Thread the two suture ends through the proximal snare wire loop (Figure 1). Pull the distal snare (Item 2).
6. Ensure about 6 cm of each suture end extends past the proximal snare (Item 3).
7. Rotate suture knob (Item 4) to pull suture ends through the SPEEDLOCK implant, until suture is visible around the suture reel (Item 5). If required, secure the suture ends onto a fixation slit on the suture reel (Item 5).
8. Align slit on the drill guide in the direction of the stitch in the tissue.
9. Remove PATHFINDER instrument prior to inserting SPEEDLOCK implant (figure 4).
10. Insert the SPEEDLOCK implant into the drill guide oriented such that the suture leading from the implant is aligned with the slit on the drill guide. Advance the implant, continuing to take up the suture slack, until it is at the opening of the drilled hole.
11. Confirm alignment of the implant such that the slot on the plug guide (Item 8) orients the suture towards the tissue (figure 4). If not aligned, reposition prior to implant insertion into the drilled hole. Insert/drive the implant into the drilled hole to the laser line on the plug guide. Note: Check that the entire length of the suture is visible, from the tissue to the implant. Do not split the sutures with the implant.

Caution! Failure to align the suture exiting the implant to the stitch may result in a loose repair or suture breakage.

Caution! Use care to properly align the implant and inserter handle with the bone hole during insertion. Do not bend or twist the inserter handle during and after insertion as damage to the implant

or incomplete insertion may result. Do not deploy bent or damaged implant.

Caution! Ensure the suture does not wrap around the inserter or drill guide when inserting the implant into bone.

Caution! Incomplete insertion or poor bone quality may result in implant pullout.

Caution! Inserting the implant beyond the laser line may result in difficulty in removing the inserter.

12. Rotate the suture knob to take up the slack in the suture and tension the tissue to the desired position (Item 11). If suture tension needs to be released, push the corresponding ratchet release button.

Caution! Do not over-tension the suture.

13. When desired tension is achieved, deploy the suture lock and detach the handle by turning the activation knob (Item 12) in the direction of the arrow approximately 1.5 complete revolutions until it stops.
14. Remove handle. Trim remaining loose suture ends at drilled hole.
15. If placing more than one implant, assure drilled holes are at least 5 mm apart and repeat steps 1 through 14.

Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE # 2 suture (thread) (white/blue co-braid)

Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)

PREPARATION FOR USE

Caution! Use only ArthroCare drill bits designed for use with these anchors. Use of other drill bits will compromise implant retention integrity.

Caution! Use only the ArthroCare Reusable Bone Punch designed for use with these anchors. Use of other punches may compromise implant insertion and/or retention.


Diagram 1 (Figure 1-2)


Diagram 2 (Figure 3)


Diagram 3 (Figure 4)


Diagram 4 (Figure 5-7)


Diagram 5 (Figure 8-9)

DIRECTIONS FOR USE

- 1) Using standard surgical technique, expose the bone at the desired insertion location for the implant.
 - 2) Attach the appropriate size drill bit to a standard orthopedic drill. Be sure to chuck the drill no farther than the flats on the proximal end of the drill.
 - 3) Position the tip of the Drill Guide with Obturator against the bone surface at the desired position (see diagram, Item 1), then remove the Obturator (see diagram, Item 2).
- Caution!** Do not advance the Drill Guide into the bone. Doing so may compromise the cortical surface of the bone and reduce implant retention capability.
- 4) Insert the drill into the drill guide and drill the bone hole. Advance the drill until the shoulder on the drill bottoms in the drill guide (see diagram, Item 3). Alternatively, insert the reusable punch into the

drill guide and drive the punch forward until the distal edge of the handle is aligned with the proximal edge of the drill guide (see diagram, Item 4).

Caution! Do not use the drill or punch without the drill guide as the bone surface can be damaged.

Caution! Care must be taken to not create too deep or too shallow a bone hole. To accomplish this, be sure to drill the hole until the shoulder on the drill bottoms in the drill guide. If using the punch, insure the handle lines up with the proximal edge of the drill guide. Care must be taken to not advance the Drill Guide into the bone as the implant retention force may be compromised.

Caution! Confirm that bone hole is clear of material by advancing and retracting the drill until the flutes of the drill are visually clear of bone fragments upon withdrawal. Failure to clear material from the bone hole may result in difficulty inserting the implant which could compromise retention force.

5) Remove the drill or punch. A PATHFINDER may be used to confirm bone hole location and depth.

Caution! Do not attempt to insert the implant into the bone without first creating a bone hole.

6) Insert the Q-FIX implant through the drill guide into the bone hole (see figures, Item 5). Be sure that the implant inserter is fully engaged with the drill guide throughout the deployment of the anchor (see figures, Item 6). This can be accomplished by holding the drill guide and inserter together by hand or by applying downward pressure to maintain the engagement during deployment. Also, be sure to maintain the alignment of the drill guide to the bone hole during insertion.

Caution! Do not use excessive force or pound the inserter into the bone hole as damage to the inserter/implant may occur or result in injury to the patient. If resistance is met on insertion, turn the inserter handle clockwise while keeping a loose hold on the drill guide. If continued resistance is met, remove the inserter and re-drill the existing bone hole.

Caution! Use care to properly align the implant and inserter handle with the bone hole during insertion. A PATHFINDER may be used to confirm bone hole location and depth. Do not bend or twist the inserter handle during and after insertion as damage to the implant or incomplete insertion may result. Do not deploy bent or damaged inserter.

Caution! Be sure to advance the inserter fully into the bone hole. Incomplete insertion or poor bone quality may result in implant pullout.

7) Rotate the activation knob clockwise on the proximal end of the implant inserter to deploy the implant (see diagram, Item 7). Rotate the knob until a hard stop is reached and the suture cleat is exposed at the proximal end of the inserter (see diagram, Item 8).

8) Release the sutures from the cleat (see diagram, Item 9) and remove the inserter and drill guide. Pull on the suture tails to ensure that the implant is set in the bone hole.

9) If placing more than one implant, assure that the anchors are located at least 7 mm apart and repeat steps 1 through 9.

7. Configuration options

1. MultiFIX P 4.5 mm Fixation implant
2. MultiFIX S 5.5 mm Fixation implant
3. MultiFIX S 6.5 mm Fixation implant
4. MultiFIX S ULTRA 5.5 mm Fixation implant
5. MultiFIX S ULTRA 6.5 mm Fixation implant

6. Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread) (white/blue co-braid)
7. Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)
8. SPEEDLOCK Fixation implant

8. Means and possibilities of integration with other medical devices

MultIFIX S ULTRA 5.5 mm, 6.5 mm fixation implant, and MultIFIX S 5.5 mm, 6.5 mm fixation implant

These devices are certified for use with № 2 MAGNUMWIRE suture thread and ultra-high molecular polyethylene suture thread. Use of other suture materials may have an adverse effect on anchoring and result in a failed outcome.

MultIFIX P 4.5 mm fixation implant

This fixation implant is intended for use with a 3.0 mm drill or a 3.4 mm awl for MULTIFIX P. Use of other tools may adversely affect implant insertion and/or anchoring (fixation).

Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread) (white/blue co-braid)

Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)

SPEEDLOCK fixation implant

These fixation implants are used together with a PathFinder obturator.

9. Medical device manufacturing materials and its core components

a. List of materials

Device (unit) name	Materials
1. MultIFIX P 4.5 mm fixation implant	RIBBON- SUTURE LOADER : Stainless steel 302 DIE : Stainless steel 316 TIE ROD : Stainless steel 304 ANCHOR, SLEEVE, SUTURE LOCKING : Polyetheretherketone LOADER : Nylon HANDLE (Top Half, bottom half) : PC KNOB, ACTUATION : PC SCREWS : Stainless steel 18-8 BLOCK, DRIVER, PTFE LOADED : PC DOWEL PIN : Stainless Steel SAE 30302/30304 CAM, DRIVER,PTFE LOADED : PC CLIP DIE : Stainless Steel 301
2. MultIFIX S 5.5 mm fixation implant MultIFIX S 6.5 mm fixation implant	RIBBON-SUTURE LOADER : Stainless steel 302 DIE : Stainless steel 17-4 PH ANCHOR, SLEEVE, SUTURE LOCKING, SCREW : Polyetheretherketone SUTURE LOADER : Nylon HANDLE (Top Half, bottom half) : POLYCARBONAT KNOB, ACTUATION : POLYCARBONAT

	SCREWS: Stainless steel 18-8 DIE BLOCK: Nylon END CAP : ABS SCREW DRIVE: ABS DOWEL PIN : Stainless Steel 18-8
3. MultiFIX S ULTRA 5.5 mm fixation implant MultiFIX S ULTRA 6.5 mm fixation implant	RIBBON- SUTURE LOADER: Stainless steel 302 DIE: Stainless steel 17-4 PH TIE ROD: Stainless steel 304 ANCHOR, SLEEVE, SUTURE LOCKING, SREW : Polyetheretherketone SUTURE LOADER: Nylon HANDLE (Top Half, bottom half): Polycarbonate KNOB, ACTUATION: Polycarbonate SCREWS: Stainless steel 18-8 DIE LOCK : Nylon END CAP : ABS SCREW DRIVE: ABS DOWEL PIN : Stainless Steel 18-8
4. Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread) (white/blue co-braid) Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)	Polyester Braided, Precut: Braided suture made of PET Suture, white/blue cobraid (thread 1): ultrahigh molecular weight medical polyethylene (UHMWPE). Suture, white/black co-braid (thread 2): ultrahigh molecular weight medical polyethylene (UHMWPE). HANDLE (Top Half, bottom half): Polycarbonate DRIVER ASSEMBLY : PC TIE ROD: Stainless steel 304 SCREWS: Stainless steel 18-8 CLIP DIE : Stainless Steel 301 or 304 DRIVER BLOCK : PC CAGE, FORCE LIMITING: ABS SPOOL: PC RATCHET PAWL, TWIST KNOB: NYLON KNOB, DEPLOYMENT TWIST: PC
5. SPEEDLOCK fixation implant	IMPLANT/FLAT PIN: polyetheretherketone DIE: Stainless steel 304 HANDLE (Top Half, bottom half): Polycarbonate KNOB, ACTUATION: PC TIE ROD : Stainless steel 302 or Stainless steel 304 WIRE: Stainless steel 304 DISTAL SNARE RING: Stainless steel 302 or Stainless Steel 316 SCREW : Stainless Steel 18-8 EXTENSION, DRIVER PUSH ROD: Stainless Steel 316 DRIVER, LASER CUT: Stainless steel 316 DIE BLOCK : Polyamide BLOCK, DRIVER, PTFE LOADED: PC REEL, SUTURE : PC PLUG: Polyetheretherketone


	RATCHET RELEASE: Polyamide KNOB : Polyamide REEL RATCHET : PC MAGNUM PI SPRING: Stainless Steel 302 CAM, DRIVER: Nylon.
--	--

10. Medical device basic parameters and performance characteristics

Tolerance for all characteristics specified therein constitutes $\pm 5\%$ unless indicated otherwise

MultiFIX P 4.5 mm fixation implant

Design drawing:

Weight	95 g
Implant weight	0,2 g
Total length	$286,766 \pm 0,254$ mm
Handle diameter	$38,1 \pm 0,254$ mm
Handle length	38,1 mm
Shaft length	$145,29 \pm 0,254$ mm
Shaft diameter	$3,175 \pm 0,0127$ mm
Implant length	17,02 mm
Implant diameter	$4,5 \pm 0,051$ mm
Hinge length	$17,018 \pm 0,762$ mm
Hinge hardness	Not less than 40 HRC
Ra hinge roughness	Not more than $0,8 \mu\text{m}$
Ra shaft roughness	Not less than $0,8 \mu\text{m}$
Shaft hardness	Not less than 41 HRC
DIE Length	215,3 mm
Outer diameter DIE	3,18 mm
Ra die roughness	Not more than $0,8 \mu\text{m}$
Die hardness	Not less than 41 HRC
Length of the loader of suture material	66,4 mm
Width of the loader of suture material	0,3 mm
Height of the loader of suture material	4,0 mm
Diameter KNOB, ACTUATION	38,1 mm
Length KNOB, ACTUATION	47,1 mm
Length CUP/DIE	14,7 mm

Width CLIP/DIE	1,2 mm
Height CLIP/DIE	8,9 mm
Ra clip/die roughness	Not more than 0,8 μm
Clip/die hardness	Not less than 36 HRC
Length BLOCK / DRIVER/PTFE LOADED	38,4 mm
Diameter BLOCK/DRIVER/PTFE LOADED	20,2 mm
Length PIN/SPIROL	9,5 mm
Diameter PIN/SPIROL	1,6 mm
Ra pin/spirol roughness	Not more than 0,8 μm
Pin/spirol hardness	Not less than 40 HRC
Dimensions of the CAM, DRIVER,PTFE LOADED	23,0x24,9x12,4

MultIFIX S 5.5 mm fixation implant

Design drawing:

Weight	98 g
Implant weight	0,3 g
Total length	291 mm
Handle length	38,1 mm
Handle diameter	38 mm
Shaft length	137 mm
Shaft diameter	4,1 mm
Implant length	14,99 mm
Implant diameter	4,2 mm
Hinge length	56 mm
Hinge hardness	Not less than 40 HRC
Ra hinge roughness	Not more than 0,8 μm
Ra shaft roughness	Not more than 0,8 μm
Shaft hardness	40 - 45 HRC
Die length	284,8 mm
Die Outer diameter	2,6 mm
Ra die roughness	Not more than 0,8 μm
Die hardness	Not less than 40 HRC
Length of the loader of suture material	66,4 mm
Width of the loader of suture material	0,3 mm

Height of the loader of suture material	4,0 mm
Diameter of the KNOB, ACTUATION	38,1 mm
Length of the KNOB, ACTUATION	47,1 mm
Length of the End Cap	21,6 mm
Diameter of the End Cap	9,5 mm
Diameter of the Die Lock	7,4 mm
Length of the Die Lock	7,9 mm
Length of the Screw Drive	84,8 mm
Inner diameter of the Screw Drive	3,2 mm
Outer diameter of the Screw Drive	11,3 mm
Length of the Dowel Pin	9,5 mm
Diameter of the Dowel Pin	1,6 mm
Ra dowel pin roughness	Not more than 0,8 μ m
Dowel Pin hardness	Not less than 90 HRB

MultiFIX S 6.5 mm fixation implant

Design drawing:


Weight	98 g
Implant weight	0,27 g
Total length	291 mm
Handle length	38,1 mm
Handle diameter	38 mm
Shaft length	137 mm
Shaft diameter	4,1 mm
Implant length	22 mm
Implant diameter	5,5 mm
Hinge length	56 mm
Hinge hardness	Not less than 40 HRC
Ra hinge roughness	Not more than 0,8 μ m
Ra shaft roughness	Not more than 0,8 μ m
Shaft hardness	40 – 45 HRC
Die length	284,8 mm
Die Outer diameter	2,6 mm
Ra die roughness	Not more than 0,8 μ m

Die hardness	Not less than 40 HRC
Length of the loader of suture material	66,4 mm
Width of the loader of suture material	0,3 mm
Height of the loader of suture material	4,0 mm
Diameter of the KNOB, ACTUATION	38,1 mm
Length of the KNOB, ACTUATION	47,1 mm
Length of the End Cap	21,6 mm
Diameter of the End Cap	9,5 mm
Diameter of the Die Lock	7,4 mm
Length of the Die Lock	7,9 mm
Length of the Screw Drive	84,8 mm
Inner diameter of the Screw Drive	3,2 mm
Outer diameter of the Screw Drive	11,3 mm
Length of the Dowel Pin	9,5 mm
Diameter of the Dowel Pin	1,6 mm
Ra dowel pin roughness	Not more than 0,8 μ m
Dowel Pin hardness	Not less than 90 HRB

MultIFIX S ULTRA 5.5 mm fixation implant

Design drawing:

Weight	97 g
Implant weight	0,3 g
Total length	291 mm
Handle length	35,9 mm
Handle diameter	38 mm
Shaft length	137 mm
Shaft diameter	4,1 mm
Implant length	22 mm
Implant diameter	5,5 mm
Hinge length	116 mm
Hinge hardness	Not less than 40 HRC
Ra hinge roughness	Not more than 0,8 μ m
Ra shaft roughness	Not more than 0,8 μ m

Shaft hardness	40 – 45 HRC
Die length	284,8 mm
Die Outer diameter	2,6 mm
Ra die roughness	Not more than 0,8 μ m
Die hardness	Not less than 40 HRC
Length of the loader of suture material	124,6 mm
Width of the loader of suture material	4,2 mm
Height of the loader of suture material	0,3 mm
Diameter of the KNOB, ACTUATION	38,1 mm
Length of the KNOB, ACTUATION	47,1 mm
Diameter of the Die Lock	7,4 mm
Length of the Die Lock	7,9 mm
Length of the End Cap	21,6 mm
Diameter of the End Cap	9,5 mm
Length of the End Cap	84,8 mm
Diameter of the End Cap	3,2 mm
Outer diameter of the Screw Drive	11,3 mm
Length of the Dowel Pin	9,5 mm
Diameter of the Dowel Pin	1,6 mm
Ra dowel pin roughness	Not more than 0,8 μ m
Dowel Pin hardness	Not less than 90 HRB

MultiFIX S ULTRA 6.5 mm fixation implant

Design drawing:


Weight	97 g
Implant weight	0,27 g
Total length	291 mm
Handle length	35,9 mm
Handle diameter	38 mm
Shaft length	137 mm
Shaft diameter	4,1 mm
Implant length	22 mm
Implant diameter	6,5 mm

Hinge length	116 mm
Hinge hardness	Not less than 40 HRC
Ra hinge roughness	Not less than 0,8 μm
Ra shaft roughness	Not less than 0,8 μm
Shaft hardness	40 – 45 HRC
Die length	284,8 mm
Die Outer diameter	2,6 mm
Ra die roughness	Not less than 0,8 μm
Die hardness	Not less than 40 HRC
Length of the loader of suture material	124,6 mm
Width of the loader of suture material	4,2 mm
Height of the loader of suture material	0,3 mm
Diameter of the KNOB, ACTUATION	38,1 mm
Length of the KNOB, ACTUATION	47,1 mm
Diameter of the Die Lock	7,4 mm
Length of the Die Lock	7,9 mm
Length of the End Cap	21,6 mm
Diameter of the End Cap	9,5 mm
Length of the End Cap	84,8 mm
Diameter of the End Cap	3,2 mm
Outer diameter of the Screw Drive	11,3 mm
Length of the Dowel Pin	9,5 mm
Diameter of the Dowel Pin	1,6 mm
Ra dowel pin roughness	Not less than 0,8 μm
Dowel Pin hardness	Not less than 90 HRB


Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread)

Design drawing:

Weight	82 g
Implant weight	0,3 g
Total length	365 mm

Shaft length	220 mm
Shaft diameter	2,11 mm
Handle length	95,1 mm
Handle diameter	38 mm
Thread type	Braided multifilament, Non-absorbable
Threads quantity	2
Thread (suture) color	Thread 1: White/blue Thread 3: white
Thread 1 length	1219,2±12,7mm
Thread 1 diameter	0,500-0,599 mm
Thread 3 length	12,31 ±0,6 mm
Thread 3 diameter	1,8 ± 0,1 mm
The minimum diameter of the thread 3 after tightening (knot diameter)	3,8
Durability of threads 1 in simple knot	>63,5 H
Durability of thread 3 in simple knot	>22 H
Ra shaft roughness	Not more than 0,8 µm
Shaft hardness	Not less than 33 HRC
Length of the Fixing part of the handle	50,8 mm
Diameter of the Fixing part of the handle	25,4 mm
Dimensions of the Clip, Die (L×W×H)	9,0×0,8×6,0
Ra Clip/die roughness	Not more than 0,8 µm
Clip/die hardness	Not less than 26 HRC
Dimensions of the Driver Block (L×W×H), mm	83,2×17,2×16,1
Dimensions Cage, Force Limiting (L×W×H), mm	13,7×10,2×13,0
Length of the Spool, Jam Cleat	46,0 mm
Outer diameter of the Spool, Jam Cleat	14 mm
Inner diameter of the Spool, Jam Cleat	2,2 mm
Dimensions of the Ratchet Pawl, (L×W×H), mm	17,5×9,3×6,3
Outer diameter Knob, Deployment Twist	38,1 mm
Inner diameter Knob, Deployment Twist	14,2 mm
Length of the Knob, Deployment Twist	63,9 mm
Maximum allowable pressure on the suture # 1	Not less than 178 N
The Force required to tighten (fasten) the suture # 1 (Suture, MagnumWire, white/blue co-braid)	Not less than 70 N
The Force required to tighten (fasten) the suture # 3 (Polyester Braided 2.8mm, Precut)	Not less than 140 N
Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads)	
Design drawing:	


Weight	84 g
Implant weight	0,3 g
Total length	368 mm
Shaft length	224 mm
Handle length	95,1 mm
Handle diameter	38 mm
Type of threads (suture)	Braided multifilament, Non-absorbable
Quantity of threads	3
Thread colors:	Thread 1: White/blue Thread 2: White/black Thread 3: White
Thread1 length	1219,2±12,7 mm
Thread 1 diameter	0,500-0,599 mm
Thread 2 length	1219,2±12,7mm
Thread 2 diameter	0,500-0,599 mm
Thread 3 length	14,78±0,8 mm
Thread 3 diameter	2,8 ±0,15 mm
The minimum diameter of the thread 3 after tightening (knot diameter)	4,65
Durability of threads 1 and 2 in simple knot	>63,5 H
Durability of thread 3 in simple knot	>22 H
Ra shaft roughness	Not more than 0,8 µm
Shaft hardness	Not less than 33 HRC
Length of the Fixing part of the handle	50,8 mm
Diameter of the Fixing part of the handle	25,4 mm
Shaft Diameter	3,0 mm
Dimensions of the Clip, Die (L×W×H)	9,0×0,8×6,0
Ra Clip/die roughness	Not more than 0,8 µm
Clip/die hardness	Not less than 26 HRC
Dimensions of the Driver Block (L×W×H), mm	83,2×17,2×16,1
Dimensions Cage, Force Limiting (L×W×H), mm	13,7×10,2×13,0
Length of the Spool, Jam Cleat	46,0 mm
Outer diameter of the Spool, Jam Cleat	14 mm

Inner diameter of the Spool, Jam Cleat	2,2 mm
Dimensions of the Ratchet Pawl, (L×W×H), mm	17,5×9,3×6,3
Outer diameter Knob, Deployment Twist	38,1 mm
Inner diameter Knob, Deployment Twist	14,2 mm
Length of the Knob, Deployment Twist	63,9 mm
Maximum allowable pressure on the suture # 1	Not less than 178 N
Maximum allowable pressure on the suture # 2	Not less than 178 N
The Force required to tighten (fasten) the suture # 1 (Suture, MagnumWire, white/blue co-braid)	Not less than 70 N
The Force required to tighten (fasten) the suture # (2 Suture, MagnumWire, white/black co-braid)	Not less than 70 N
The Force required to tighten (fasten) the suture # 3 (Polyester Braided 2.8mm, Precut)	Not less than 140 N

SPEEDLOCK fixation implant

Design drawing:

Weight	93 g
Total length	345 mm
Handle length	95,1 mm
Handle diameter	38 mm
Shaft length	170 mm

Shaft diameter	3,4 mm
Implant length	12,3 mm
Implant diameter	3,43 mm
Hinge diameter	3,4 mm
Hinge hardness	Not less than 38 HRC
Ra hinge roughness	Not more than 0,8 μm
Ra shaft roughness	Not more than 0,8 μm
Shaft hardness	Not less than 31 HRC
Die length	226, 2 mm
Die diameter	3,4 mm
Ra die roughness	Not more than 0,8 μm
Die hardness	Not less than 30 HRC
Diameter of the Knob, Actuation	38,1 mm
Length of the Knob, Actuation	47,1 mm
Hinge hardness	Not less than 29 HRC
Length of the Extension, Driver Push Rod	29,0 mm
Diameter of the Extension, Driver Push Rod	2,1 mm
Ra extension, driver push rod roughness	Not more than 0,8 μm
Extension, driver push rod hardness	Not less than 40 HRC
Length the Driver, Laser Cut	244,6 mm
Outer diameter of the Driver, Laser Cut	2,1 mm
Ra driver, laser cut roughness	Not more than 0,8 μm
Driver, laser cut hardness	Not less than 30 HRC
Dimensions of the die block (L×W×H), mm	35,0×15,2×12,4
Length of the Block, Driver, PTFE Loaded, mm	38,4 mm
Diameter of the Block, Driver, PTFE Loaded, mm	20,2 mm
Diameter of the Reel, Suture, mm	26,9 mm
Height of the Reel Suture, mm	50,8 mm
Dimensions of the plug (L×W×H), mm	8,3×1,3×1,3
Dimensions of the Ratchet release (L×W), mm	26,6×27,3×15,2
Diameter of the Knob	24,9 mm
Height of the Knob	15,9 mm
Dimensions of the Reel Ratchet (L×W×H)	8,0 × Ø26,9
Inner diameter of the Distal Snare Ring	22,0 mm
Width of the Distal Snare Ring	3,0 mm
Ra distal snare ring roughness	Not more than 0,8 μm
Distal snare ring hardness	Not less than 43 HRC

11. Draft Labeling in Russian language

Label in Russian contains the following information:

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*

Information about "name of the medical product" is provided in accordance with the registration certificate, after it's received.

Information about "size" for configuration options MultiFIX P 4.5 mm Fixation implant, MultiFIX S 5.5 mm Fixation implant, MultiFIX S 6.5 mm Fixation implant, MultiFIX S ULTRA 5.5 mm Fixation implant, MultiFIX S ULTRA 6.5 mm Fixation implant,

Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread) (white/blue co-braid), Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)

"Precautions" for all variants is provided in accordance with the data from the technical documentation for this type of products.

For the MultiFIX P 4.5 mm Fixation implant, MultiFIX S 5.5 mm Fixation implant, MultiFIX S 6.5 mm Fixation implant, MultiFIX S ULTRA 5.5 mm Fixation implant, MultiFIX S ULTRA 6.5 mm Fixation implant, "Sterilization method" is provided like "radiation sterilization".

For the Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread) (white/blue co-braid), Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid), SPEEDLOCK Fixation implant is provided like "ethylene oxide".

The words "Apirogenic", "Sterile", "non-toxic" also presented on the marking for the all variants. Information about the "storage conditions" – according to the data from the technical documentation for this type of product.

Information about material of suture, suture structure, size, length for the Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread) (white/blue co-braid), Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid) is indicated according to the technical documentation.

For each variant of execution the separate label is glued

11.1 Labeling symbols

Symbol	Description
	Manufacturer
	Consult IFU
	Radiation-sterilized
	ETO-sterilized
	Single-use only
	Do not re-sterilized
	Do not use if package is damaged
	Use by (YYYY/MM/DD: year/month/day)
	Catalog number
	Serial number/Lot
	MR Safe
	Caution: Federal (USA) Law restricts this device to sale by or on the order of a physician
	Keep away from sunlight
	Keep dry
	Authorized representative in the EC

Symbol	Description
	CE mark and Identification number of Notified Body. The product meets the essential requirements of Medical Device Directive (93/42/EEC)
	Quantity of devices
MADE IN	Made in

12. Information about medical substances included to composition of the medical device

Non-absorbable arthroscopic suture anchors do not contain drug substances.

13. List of animal and(or) human-derived components included to composition of the medical device

Non-absorbable arthroscopic suture anchors do not contain human or animal-derived components.

14. Package information

Device is placed in a plastic tray and sealed with Tyvek cover. Next, the device (in sterile package) is placed inside consumer packaging along with an IFU. Finally, consumer packages are placed in shipping cartons.

Package properties (for all types of packing):

Transport packaging (for all types of packing)	
Transport packaging dimensions and weight	Dimensions depend on the order, but not exceeding 200x200x200 cm. Transport package weight depends on its dimensions and content but not exceeding 100 pounds (45.36 kg).
Transport package material	Package: cardboard boxes

15. Transportation and storage conditions

Keep in dry places away from heat sources.

	Storage conditions	Transportation conditions
Temperature	From + 5°C till + 30 °C	From -10°C till + 35 °C
Humidity	Up to 75%	Up to 80%
Pressure	86 kPa - 106 kPa	86 kPa - 106 kPa

16. Service life, shelf life

Devices are intended for permanent implantation at normal use conditions. Implant extraction is a surgical operation compatible with implant installation procedure. Observe aseptic requirements and compliant surgical methods. Implant extraction decision should be properly decided with assessment of all possible risks related with additional surgical intrusion.

Medical device shelf life: 3 years.

17. Safety requirements and precaution measures *

MultiFIX P 4.5 mm fixation implant

Warnings

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due to changes in the functioning of the implant.

Precaution measures

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due to changes in the functioning of the implant.

MultiFIX S 5.5 mm fixation implant, MultiFIX S 6.5 mm fixation implant *

MultiFIX S ULTRA 5.5 mm fixation implant, MultiFIX S ULTRA 6.5 mm fixation implant

Warnings

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due to changes in the functioning of the implant.

Precaution measures

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due to changes in the functioning of the implant.

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Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)*

Warnings

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due to changes in the functioning of the implant.

Precaution measures

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due

to changes in the functioning of the implant.

SPEEDLOCK fixation implant *

Warnings

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due to changes in the functioning of the implant.

Precaution measures

- There are no possible warnings in the normal daily activities of the patient after implantation.
- There are no warnings to prevent harmful effects on the patient's health after implantation due to changes in the functioning of the implant.

*Please follow the original IFU.

18. Equipment methods of decontamination and pre-sterilization cleaning

Non-absorbable arthroscopic suture anchors are supplied sterile and do not require sterilization. Suture anchors are single use for one patient only.

19. Sterilization methods and conditions

Device	Sterilization method
MultiFIX P 4.5 mm fixation implant	Radiation
MultiFIX S 5.5 mm fixation implant	Radiation
MultiFIX S 6.5 mm fixation implant	Radiation
MultiFIX S ULTRA 5.5 mm fixation implant	Radiation
MultiFIX S ULTRA 6.5 mm fixation implant	Radiation
Q-Fix 1.8 mm Fixation implant with 1 MAGNUMWIRE #2 suture (thread) (white/blue co-braid)	Ethylene oxide
Q-Fix 2.8 mm Fixation implant with 2 MAGNUMWIRE #2 sutures (threads) (white/blue and white/black co-braid)	Ethylene oxide
SPEEDLOCK fixation implant	Ethylene oxide

Sterility period is at least 3 years old.

20. Maintenance information and validation data

Not applicable to this device.

21. Environment protection

21.1 Requirements to environmental protection during medical equipment application

Non-absorbable arthroscopic suture anchors have no environmental impact in process of use, transportation and storage.

21.2 Disposal and utilization order, indications (if any) for special precaution measures, required in process of non-used medical devices disposal

Non-absorbable arthroscopic suture anchors should be disposed of in accordance with applicable national regulations and plans/manuals for medical waste disposal.

If the device package is damaged this device is disposed of in accordance with Class A medical waste disposal procedures under SanPiN 2.1.7.2790-10.

Used medical device is disposed of according to Class B medical waste under SanPiN 2.1.7.2790-10.

22 Risks management summary

Summary (Overall Residual Risk):

The risk analysis performed provides assurance that the appropriate risks have been identified and either eliminated or reduced as low as possible. For the risks that remain, risk controls appropriate to each risk are provided to minimize the hazard from potential failures during the projected expected service life of the device.

Results of the risk analysis demonstrate that use of the device as intended does not adversely affect the health or safety of the patient, user, or other persons.

The Overall Residual Risk evaluation analysis shows that there are no risks which require further controls.

23 List of national and international regulatory documents/standards applicable to the medical device

Identification №	Title
93/42/EEC last revised and amended by 2007/47/EC	Medical Devices Directive
BS EN ISO 13485:2016	Medical Devices – Quality management systems – Requirements for regulatory purposes
BS EN ISO 14971:2012	Application of Risk Management to Medical Devices
BS EN 1041:2008+A1:2013	Information supplied by the manufacturer with medical devices
BS EN 980:2008	Symbols for use in the labelling of medical devices
BS EN ISO 15233-1:2012 +CORR:2012	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ASTM D4169:2014	Standard Practice for Performance Testing of Shipping Containers and Systems
ASTM F88/F88M:2009	Standard Test Method for Flexible Barrier Materials
ASTM F1839 8 th edition: 2012	Standard Specification for Rigid Polyurethane Foam for use as a Standard Material for Testing Orthopedic Devices and Instruments
BS EN ISO 14644-1:2015	Cleanrooms and associated controlled environments – Part 1: Classification of air cleanliness
BS EN ISO 10993-1:2010 +CORR:2010	Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process
EN ISO 10993-5:2009	Biological Evaluation of Medical Devices – Part 5: Tests for Cytotoxicity

Identification №	Title
ISO 10993-10:2013	Biological Evaluation of Medical Devices – Part 10: Tests for Irritation and Delayed Hypersensitization
EN ISO 10993-6:2009	Biological Evaluation of Medical Devices – Part 6: Tests for local effects of implantation
BS EN ISO 10993-3:2014	Biological Evaluation of Medical Devices – Part 3: Tests for Genotoxicity, Carcinogenicity & Reproductive Toxicity
BS EN ISO 10993-11:2009	Biological Evaluation of Medical Devices – Part 11: Tests for Systemic Toxicity
BS EN 14630:2012	Non-active surgical implants – General Requirements
BS EN 556-1:2001	Sterilization of medical devices – Part 1: Requirements for medical devices to be designated "STERILE"
BS EN ISO 11607-1:2009	Packaging for terminally sterilized medical devices– Part 1: Requirements for materials, sterile barrier systems and packaging systems
BS EN ISO 11607-2:2010	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
BS EN ISO 11737-1:2006 +CORR1:2008, +CORR2:2011	Sterilization of medical devices – Part 1: Microbiological methods. Determination of a population of microorganisms on products
BS EN ISO 11737-2:2009	Sterilization of medical devices. Microbiological methods – Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
ASTM F1980 7 th editions: 2011	Standard for Accelerate Aging of Sterile Barrier Systems for Medical Devices
BS EN ISO 10993-7:2008	Biological evaluation of medical devices—Part 7: Ethylene oxide sterilization residuals
AAMI TIR 28:2009	Product adoption and process equivalency for ethylene oxide sterilization
ISO 11135:2014	Sterilization of health-care products — Ethylene oxide — Requirements for the development, validation and routine control of a sterilization process for medical devices
EN ISO 11137-2:2007	Sterilization of health care products – Radiation. Part 2 – Establishing the sterilization dose
BS EN 62366-1:2015	Medical devices – application of risk management to medical devices
BS EN ISO 14155:2011	Clinical Investigation of medical devices for human subjects – Good Clinical practice

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The manufacturer guarantees conformity of declared characteristics, term of preservation of sterility,

guarantee of safety and integrity of packing at correct transportation, operation and storage of the product during 3 years from the moment of sterilization Products.

25 Complaints

All questions, complaints, or return authorization requests should be directed to Customer Service or an authorized representative.

Address for reclamations:

LLC "Smith & Nephew"

Russia, 105120, Moscow, 2nd Syromiatnicheskiy lane, house 1, 9th floor, room 1, office 1.

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«АртроКэр Корпорейшн»
(ArthroCare Corporation)
(Organization/Organization)

7000 West William Cannon Drive
Austin TX 78735-8531, USA (США)
(Address/Address)

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