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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Suzhou Semorr Medical Technology Co., Ltd. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized

Signature:

Yu Ying

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(Date: Nov. 25, 2019)

证明书查询网址: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

All-Russian product
classification code
OKPD2 26.70.22.150

Approved:
Signature/SEAL
President
06.10.2019

[Handwritten signature]

INSTRUCTION FOR USERS for Medical Device:

Operation Microscope *Just Vision micro JVM*
Configurations: JVMTME, JVMA, JVMB, JVMC, JVMD,
JVME

Manufacturer: *Suzhou Semorr Medical Technology Co., Ltd*

Complies with the requirements of international and
Russian standards

GOST R 50444-92, GOST R IEC 60601-1-2010

Contents

1. Product Designation
2. Supply Package
3. General Use and Operation Information
4. Indications
5. Contraindications
6. Product Description
7. Product Materials
8. Condition and rules of product operation
9. Product Specifications
10. Labeling
11. Packaging
12. Transportation
13. Storage
14. Maintenance
15. Minor Repairs and Warranties
16. Disposal and Destruction Information
17. List of National and International Standards
18. List of Data Concerning International Certificates and Declarations of Conformity
19. The main steps of the engineering and manufacturing processes
20. Electromagnetic compatibility table
21. Manufacturer's guidelines and statement
22. Manufacturer's guidelines and statement
23. Recommended distance between portable and movable equipment

1. PRODUCT DESIGNATION

Operation microscope is designed for image scaling in course of surgeries.

2. SUPPLY PACKAGE

Operation microscope Just Vision micro JVM, configurations: JVMTME, JVMA, JVMB, JVMC, JVMD, JVME

Set of the product:

1. Microscope
2. Binocular loupes - (if necessary)
3. Lens objectives – (if necessary)
4. Magnification knob – (if necessary)
5. Protection shield– (if necessary)
6. Swing arm – (if necessary)
7. Cross arm - (if necessary)
8. Suspension arm – (if necessary)
9. Standard handle – (if necessary)
10. Oblique handle – (if necessary)
11. T-handle – (if necessary)
12. Eyepiece – (if necessary)
13. Optical fiber – (if necessary)
14. Optical double ray divider – (if necessary)
15. Imaging camera – (if necessary)
16. Adapter – (if necessary)
17. Photoflash ring – (if necessary)
18. Swing ring – (if necessary)
19. Extender with a swing ring - (if necessary)
20. Splitter with a vertically adjustable head - (if necessary)
21. Splitter without a double divider/with double divider - (if necessary)
22. Iris diaphragm – (if necessary)
23. Amplifier – (if necessary)
24. Bracket – (if necessary)
25. Lens with a focusing device F175 – (if necessary)
26. Lens with a focusing device F200 – (if necessary)
27. Lens with a focusing device F250 – (if necessary)
28. Lens with a focusing device F300 – (if necessary)
29. Lens with a focusing device F350 – (if necessary)
30. Lens with a focusing device F400 – (if necessary)
31. Lens with changeable focus F200-F400 – (if necessary)
32. Set of silicone plates – (if necessary)
33. Zoom lens (with adjustable working distance) – (if necessary)
34. Assistant binocular unit – (if necessary)
35. Portable case – (if necessary)
36. Mobile rolling column – (if necessary)
37. Instruction for users

(hereinafter – Operation Microscope)

3. GENERAL USE AND OPERATION INFORMATION

Product use conditions: general, plastic and dental surgery units, ENT units of healthcare facilities.

Operation microscope represents a multiple-use product.

Operation microscope represents a non-sterile medical product.

In terms of potential application risk, the operation microscope grades as class 2a (according to the Decree of the Ministry of Health of the Russian Federation No. 4n dated 06.06.2012).

All-Russian product classification code by types of economic activities: OKPD2 26.70.22.150.

Connection of any products not included in the operation microscope supply package is strictly forbidden.

The operation microscope is used in a patient environment.

Patient touching of any parts of the operation microscope under maintenance is strictly forbidden.

In accordance with GOST R 56169, the product refers to the operation microscope group.

The product is not designed for eye surgery!

The product refers to I-class products in re electric shock hazard protection and to B-type working section protection degree in accordance with GOST R IEC 60601-1.

WARNING: "Caution! To avoid electric shock hazard, the product must be connected only to a line supply having a protective grounding."

The microscope is connected to power supply via standard network cable (not supplied as part of package).

4. INDICATIONS

Operation microscope is designed for image scaling in course of surgeries.

5. CONTRAINDICATIONS

Do not use the operation microscope to light the eye area.

6. PRODUCT DESCRIPTION

Operation microscope Just Vision micro JVM, configurations: JVMTME, JVMA, JVMB, JVMC, JVMD, JVME with accessories

Features of the microscope configurations:

JVMTME: 5-steps magnification

JVMA: 6-steps magnification

JVMB: zoom magnification

JVMC: 5-steps magnification

JVMD: 6-steps magnification

JVME: zoom magnification

Figure 1. General view of the microscope

Availability of the accessories is stipulated by a particular filed of application and individual preferences of attending specialists.


The microscope is designed for image scaling.

Accessories:

- Binocular loupes, lens, magnification knob, protection shield, eyepiece, optical double ray divider, optical fiber, lenses with focusing devices, zoom lenses, iris diaphragm, and amplifier are used to scale a particular area of an attending specialist's interest depending on a particular clinical situation.
- Handles and swing ring are used to adjust the position of the microscope in course of work with the medical device.
- Imaging camera, adapter, and splitters are used to record the scaling area.
- Set of silicone plates is used for image stabilization.
- Swing arm, cross arm, suspension arm, bracket, and extender with a swing ring are used to place the microscope conveniently and to fix it in an operation office.
- Assistant binocular unit is an auxiliary unit of the microscope with binocular loupes for joint operation of an attending specialist and an assistant. The image on the assistant unit is duplicates the image transmitted to the main area of work of the attending specialist.
- Photoflash ring is used to fix the photoflash to the lens of the device.

- Portable case is used to transport the microscope among healthcare facility units.
- Mobile rolling column is used to fix the accessories and auxiliary in course of operational researches.

Configuration variant of the assembled microscope with accessories:

1. Wheels to move the microscope
2. Wheelbase
3. Bracket
4. Lock knob to prevent uncontrolled rotation
5. Cross arm
6. Lock knob of the swing arm
7. Swing arm
8. Balance knob
9. Height adjuster
10. Adapter
11. Microscope arm
12. Optical fiber and optical double ray divider
13. Handle
14. Magnification knob
15. Lenses (lens with a focusing device or zoom lens)
16. Lens focusing knob
17. Protection shield
18. Microscope body
19. Lock screw for binocular loupes
20. Binocular loupes
21. Distance adjustment knob

- 22. Eyepiece
- 23. Swing ring
- 24. Splitter
- 25. Brightness adjustment knob
- 26. Light socket
- 27. Imaging camera socket
- 28. Power switch
- 29. HDMI socket
- 30. Power socket

7. CONDITIONS AND RULES OF PRODUCT OPERATION

Installation and preparation for use:

Open the product's protective package. Mount the microscope bracket on the wheelbase so that the slot of the column is aligned with the cylindrical pin on the base. Tighten the screw using a 10-mm wrench. Base installation is completed. Mount the cross arm, the swing arm and the optical fiber on the bracket. Mount the handle in the package using locking screws. Fix the handle to the microscope head using a 3-mm wrench. Mount the protection shield and lens objectives, whereupon connect the cables to the power supply.

Before use:

Connect the microscope to the main power supply. Move the power switch to the ON position, making sure that the microscope is functioning properly. After inspection, shut down the power supply.

Operation of the microscope:

- Loose lock knobs to test flexibility of all mounts.
- Holding the handle, move the microscope to its focal plane (250mm from eyepiece lenses to the operation surface). Then fasten the lock knob to get a stable observing position.
- Fasten the height adjuster of the swing arm to prevent it from moving.
- Loosen the height adjuster; holding the handle, move the microscope up and down. Rotate the balance knob counterclockwise to mount the lock knob on the as shown in Fig. 4. Repeat the step until the up and down movement of the swing arm is flexible.


Fig. 3



Fig. 4

Adjustment during use:

- Move the microscope to the required location. Lock the brakes of the wheels to fix the microscope.
- Connect the power cable to the power supply. Move the power switch to the ON position.
- Using the brightness adjustment knob, set the required brightness. Minimum brightness is set by clockwise rotation, maximum brightness is set by counterclockwise rotation.
- Set the image focus using the magnification knob.
- Diopter regulation: turn the pupil distance adjustment knob under the eyepieces [21] until a clear and sharp stereo image appears on the microscope.

- Filter switch: apply filters in or out of the optical path using the splitter.

After use:

- Fold the microscope arms on the bracket. Tighten lock knobs of all joints to fix the microscope.
- Shut down the power supply.

- Loosen the brakes of the wheels and, holding the cross arm, carefully move microscope, avoiding collision and dump.

8. PRODUCT MATERIALS

Operation Microscope Just Vision micro JVM Configurations: JVMTME, JVMA, JVMB, JVMC, JVMD, JVME with accessories

No.	System, Component Name	Color Marks and Materials	Material/ Trademark	Name	Human Contact Type
1	Operation microscope Just Vision micro JVM	N/A	Stainless steel ABS resin Polycarbonate Aluminum Cast iron	1. SUS304, China 2. 765A, China 3. SC1004A-EM59P, China 4. Al6061, China 5. HT250, China	Brief, with uninjured skin
2	Binocular loupes	N/A	Optical flint Polycarbonate Stainless steel	1. ZK7, China 2. SC1004A-EM59P, China 3. SUS304, China	Brief, with uninjured skin
3	Lenses, Lens with a focusing device F175, Lens with a focusing device F200, Lens with a focusing device F250, Lens with a focusing device F300, Lens with a focusing device F350, Objective with a focusing device F400, Lens with changeable focus F200-F400, Zoom lens, eyepiece, iris diaphragm	N/A	Optical flint Stainless steel	1. ZK7, China 2. SUS304, China	Brief, with uninjured skin
4	Magnification knob	N/A	Polycarbonate	SC1004A-EM59P, China	Brief, with uninjured skin
5	Protective shield, portable case	N/A	ABS resin	765A, China	Brief, with uninjured skin
6	Swing arm, cross arm, suspension arm, standard handle, oblique handle, T-handle	N/A	Stainless steel	SUS304, China	Brief, with uninjured skin
8	Optical fiber, optical double ray divider	N/A	1. Optical flint 2. ABS resin 3. Stainless steel	1. ZK7, China 2. 765A, China 3. SUS304, China	Brief, with uninjured skin
9	Imaging camera	N/A	1. Optical flint 2. Polycarbonate 3. ABS resin	1. ZK7, China 2. SC1004A-EM59P, China 3. 765A, China	Brief, with uninjured skin
10	Adapter, photoflash ring, swing ring, extender with a swing ring, splitter with a vertically adjustable head, splitter without double divider/with double divider, amplifier, bracket, mobile rolling column	N/A	Stainless steel ABS resin	1. SUS304, China 2. 765A, China	Brief, with uninjured skin
11	Set of silicone plates	N/A	Silicone	Sun Group LCM 1234, China	Brief, with uninjured skin
12	Assistant binocular unit	N/A	1. Optical flint 2. Polycarbonate 3. Stainless steel 4. Aluminum	1. ZK7, China 2. SC1004A-EM59P, China 3. SUS304, China 4. Al6061, China	Brief, with uninjured skin

9. PRODUCT SPECIFICATIONS

- 9.1 Dimensions and weight of the assembled operation microscope:
- Length: 9450 ± 50 mm
 - Width: 1390 ± 50 mm
 - Height: 1670 ± 50 mm
 - Weight: 115kg max.
- 9.2 Eyepiece inclination: $0-210^\circ$
- 9.3 Adjustable pupil distance: 50-75mm
- 9.4 Microscope magnification (for all microscope models): 2.7x-24x
- 9.5 Focal distance: up to 250mm
- 9.6 FOV diameter: from 9mm to 83mm
- 9.7 Illumination diameter: 85mm
- 9.8 The operation microscope must provide illumination of at least 60000lx at a working distance of 250mm
- 9.9 Built-in LED-lighting with a maximum power of 20W
- 9.10 Arm gyration max: 460mm
- 9.11 Input power: 130 ± 10 W
- 9.12 Input voltage: 220 ± 22 V
- 9.13 Frequency: $50 \pm 0,5$ Hz
- 9.14 Built-in filter: orange
- 9.15 Microscope illumination sport diameter: no more than 85mm
- 9.16 Cracks, blisters, nicks, scratches, chips, burrs, drosses, grinding and polishing material particles are not admissible.
- 9.17 Products must withstand temperature and air humidity effects in the process of transportation and storage under group-5 conditions (OZH4 – particularly severe) according to GOST 15150.
- 9.18 During operation, products must be resistant to environment climatic aspects effect in accordance with the requirements to products designed for UKHL 4.2 (moderately cool climate) according to GOST 15150.
- 9.19 Products in a transportation package must be resistant to mechanical during transportation and comply with GOST-R 50444.
- 9.20 Requirements to optomechanical characteristics:
- Allowance for compound magnification – $\pm 7.5\%$
 - Difference in magnification of right and left optical systems – not over 1.5%
 - Divergence of right and left optical systems' axes:
 - vertical – not over 15 degrees
 - horizontal: convergence – not over 45 degrees, divergence – not over 10 degrees

- Focusing planes displacement at magnification alteration
 - Difference in focusing of right and left optical systems
 - At maximum magnification, resolution in the center FOV must comprise not over
 - difference in image rotation in right and left optical systems – not over 1 degree
- Eyepiece:

- difference in height of exit pupils of right and left optical systems – not over 0.5mm
- no diopter scale calibration error
- minimum regulation interval:
- overall – $\pm 5D$
- by eye apex – $\pm 2D$

9.21. Rendering camera specifications and ports:

- The rendering camera is connected to the microscope via HDMI cable (not supplied as part of package). The camera also has a USB port for connection to the microscope and powering from it. The USB cable is not supplied as part of package.

Camera specifications:

Video format	JPEG color images, MPEG color videos
Video resolution, px	1920 x 1080
Framerate, rates per second	60
Overall dimensions, mm	$64 \pm 1 \times 66 \pm 1 \times 29 \pm 1$
Weight, g, not over	1000

10. LABELLING

10.1. The product labeling must comply with indications in the user manual and the requirements of GOST R 50444-92, GOST R IEC 60601-1-2010.

10.2. Every operation microscope must have a plate containing:

- name of the manufacturer,
- product name,
- month and year of manufacture,
- serial number,
- voltage (V);
- network frequency, Hz
- input power (VA),

- sign "Throwing out is forbidden. Special disposal is required"
- B-type working section protection sign,
- manufacturer's address and phone number,
- package contents description,
- lot (batch) number.

Additional product information inclusion is admissible.

11. PACKAGING

Individual packaging consists of a carton box with inserted carton dividers and foam fillers.

12. TRANSPORTATION

The product is packed into an appropriate packaging together with operational documentation. The product transportation conditions must correspond to storage conditions 5 according to GOST 15150 (at an ambient temperature ranging from -50 to +50°C and relative humidity of up to 100%)

13. STORAGE

Conditions for storage of the product in a transportation package at the manufacturer's (user's) warehouses must correspond to storage conditions 2 according to GOST 15150 (at an ambient temperature ranging from -40 to +50°C and relative humidity of up to 98%)

14. MAINTENANCE

Cleaning:

Binocular loupes:

- The binocular loupes must be cleaned with a watered or alcohol-wetted soft cloth. Do not use any solvents or detergents.
- Sterilization of the loupes with any other agents is forbidden.

Operation microscope:

- Make sure the device is switched off before cleaning.
- The outer surface of the product may be cleaned with a damp cloth, 3-% penozone with 0.5-% solution of a detergent like *Lothos* or 1-% chloramine solution.
- To disinfect the outer surface, use a soft sterilizing agent in accordance with the manufacturer's instructions.
- Application of abrasives or tough detergents is forbidden.

- Avoid liquid getting into the microscope.
- Sterilization of the product with any other agents is forbidden.

Fuse replacement:

When fuses blow (fuses blow at 250V), replace them with fuses having identical rating (M10AL250V). To do that, press the fuse holder. The holder will move out as shown in Figure 1.

Figure 1. Fuse holder in the open position

Thereupon, take the fuses (2 pieces) out from the holder by pushing them (the fuses will slide out) and removing; then insert new fuses in their place and press the holder, which will return to the initial position, as shown in Figure 2.


Figure 2. Fuse holder in the open position

15. MINOR REPAIRS AND WARRANTIES

The user does not carry out minor repairs.

Warranty for materials and quality of the microscope is provided for one (1) year from the initial trade date. If the product breaks down due to defects in the material or improper quality within the one (1)-year warranty period, the product is repaired or replaced.

In case of warranty event, contact the manufacturer's agent:

Dr. Kim Ltd., ul. Turistskaya 28/1-390, Saint Petersburg 197082

16. DISPOSAL AND DESTRUCTION INFORMATION

Medical devices, packing material and accessories must be disposed in accordance with applicable national laws and standards. In accordance with SanPiN 2.1.7.2790, the product grades as A-class – epidemiologic non-hazardous wastes.

17. LIST OF NATIONAL AND INTERNATIONAL STANDARDS


The list of applicable RF national standards and international standards for safety, efficiency and quality assurance, as well as standards with test methods is as follows: GOST R 50444-92, GOST-R IEC 60601-1-2010, EN ISO 13485:2016, FDA certificate No. 3013154160, CE certificate No. 02290.

The product is manufactured in accordance with the FDA and CE certificates.

18. LIST OF DATA CONCERNING INTERNATIONAL CERTIFICATES AND DECLARATIONS

- Certificate EN ISO 13485:2016 dated 06.11.2017
- Certificate FDA No. 3013154160
- Company registration certificate

19. THE MAIN STEPS OF THE ENGINEERING AND MANUFACTURING PROCESSES

20. ELECTROMAGNETIC COMPATIBILITY TABLE

Electromagnetic Compatibility Tables

Table 1 – MANUFACTURER'S Guidelines and Statement concerning ELECTROMAGNETIC RADIATION for MEDICAL EQUIPMENT and MEDICAL SYSTEMS

Manufacturer's Guidelines and Statement – Electromagnetic Radiation			
Radiation Check	Value	Compatibility	Electromagnetic Environment – Guidelines
RF radiation CISPR11		Group 1	RF energy is used only for internal operation. Therefore, RF radiation is very insignificant and will probably not cause interferences for electronic installations nearby. Suitable for use in all institutions other than household facilities – thus, the device is connected directly to a low-voltage public power network providing power supply to buildings used for household purposes.
RF radiation CISPR11		Class A	
Harmonic IEC 61000-3-2		Class A	
Flicker IEC 61000-3-3		Compatible	

21. MANUFACTURER'S GUIDELINES AND STATEMENT CONCERNING ELECTROMAGNETIC IMMUNITY FOR MEDICAL EQUIPMENT AND MEDICAL SYSTEMS

Table 2 – MANUFACTURER'S Guidelines and Statement Concerning ELECTROMAGNETIC IMMUNITY for MEDICAL EQUIPMENT and MEDICAL SYSTEMS

Manufacturer's Guidelines and Statement – Electromagnetic Immunity			
Immunity Check	EN/IEC 60601 Test Level	Immunity Requirements Compliance Level	Electromagnetic Environment – Guidelines
Static electricity EN/IEC 61000-4-2	±6kV contact ±8kV air	±6kV contact ±8kV air	The floor must be wooden, concrete or tiled. For synthetic floors, relative humidity must be at least 30%.
Microsecond pulse disturbances EN/IEC 61000-4-4	±2kV network ±1kV inputs/outputs	±2kV network ±1kV inputs/outputs	Mains supply quality must correspond to mains supply quality in an industrial and household or hospital environment.
Surge voltage/Surge current EN/IEC 61000-4-5	±1kV differential ±2kV general	±1kV differential ±2kV general	Mains supply quality must correspond to mains supply quality in an industrial and household or hospital environment.
Voltage dips / Failures EN/IEC 61000-4-11	Dip>95% for 0.5 cycle Dip60% for 5 cycles Dip30% for 25 cycles Dip>95% for 5 seconds	Dip>100% for 0.5 cycle Dip60% for 5 cycles Dip30% for 25 cycles Note 1	Mains supply quality must correspond to mains supply quality in an industrial and household or hospital environment. Whenever continuous operation is required during power interruptions, it is recommended to connect to an uninterruptible power supply or an accumulator battery.
Industrial frequency 50/60Hz	3A/m	3A/m	Power frequency magnetic fields must correspond to

Magnetic field EN/IEC 61000-4-8			standard power frequency magnetic fields in an industrial and household or hospital environment.
------------------------------------	--	--	---

Note 1: The equipment under test turns off, and the operator must manually restart it in 5 seconds after the AC power outage.

22. MANUFACTURER'S GUIDELINES AND STATEMENT CONCERNING ELECTROMAGNETIC IMMUNITY FOR MEDICAL EQUIPMENT AND MEDICAL SYSTEMS THAT DO NOT REPRESENT LIFE-SUPPORT EQUIPMENT OR SYSTEMS

Table 4 – MANUFACTURER'S Guidelines and Statement concerning ELECTROMAGNETIC IMMUNITY for MEDICAL EQUIPMENT and MEDICAL SYSTEMS that do not represent LIFE-SUPPORT EQUIPMENT OR SYSTEMS

Manufacturer's Guidelines and Statement – Electromagnetic Immunity

Immunity Check	EN/IEC 60601 Test Level	Immunity Requirements Compliance Level	Electromagnetic Environment – Guidelines
Induces radio waves/interferences EN/IEC 61000-4-6	3V (r.m.s. voltage) from 150kHz to 80MHz	(3)V (r.m.s. voltage)	The below calculated/specified minimal distance must be observed between the portable and movable communications equipment: $D = (3.5/V1)(P \text{ sq. root})$ $D = (3.5/E1)(P \text{ sq. root})$ from 80 to 800 MHz $D = (7/E1)(P \text{ sq. root})$ from 800 MHz to 2.5 GHz
Radiated radio waves EN/IEC 61000-4-3	3V (r.m.s. voltage) from 80MHz to 2.5GHz	(3)V/m	Where P – max. power in Watts, and D – recommended distance in meters. The field intensity created by transmitters and determined via site EM survey, must not exceed the immunity requirements compliance levels (V1 and E1). Interferences may arise near devices equipped with transmitters.

23. Recommended distances between portable and movable communications equipment and equipment that does represent not a life-support equipment or system

Table 6 – Recommended Distances between Portable and Movable Communications Equipment and EQUIPMENT that does represent not a LIFE-SUPPORT EQUIPMENT or SYSTEM

Max output power (W)	Distance (m) from 150 kHz to 80 MHz $D = (3,5/V1)(P \text{ sq. root})$	Distance (m) from 150 kHz to 80 MHz $D = (3,5/E1)(P \text{ sq. root})$	Distance (m) from 150 kHz to 80 MHz $D = (7/E1)(P \text{ sq. root})$
0.01	0.1166	0.1166	0.2333
0.1	0.3689	0.3689	0.7378
1	1.1666	1.1666	2.3333
10	3.6893	3.6893	7.3786
100	11.6666	11.6666	23.333

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