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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of CHONGQING JINSHAN
SCIENCE & TECHNOLOGY (GROUP) CO., LTD on the annexed
DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明真伪网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED BY

General Manager

Wung imshan

“ ” 9.10 2019

Operator's Manual

Hysteroscopy System for Vacuum Aspiration Abortion

Chongqing Jinshan Science & Technology (Group) Co., Ltd., China

1. Overview of the medical device

The medical device name is "Hysteroscopy system for vacuum aspiration abortion" (hereinafter referred to as the device).

Manufacturer: Chongqing Jinshan Science & Technology (Group) Co., Ltd., China
Manufacturer's address: No. 18, Nishang Road, LiangLu Industrial City, 401120 Yubei District, Chongqing, China
Phone: +86-23-86098099
Manufacturing site: Chongqing Jinshan Science & Technology (Group) Co., Ltd., China
No. 18, Nishang Road, LiangLu Industrial City, 401120 Yubei District, Chongqing, China

Development, manufacture, and sale of the device are carried out in compliance with the requirements of ISO 13485. Notified body: TUV Rheinland LGA Products GmbH.

The device is CE marked. Notified body: TUV Rheinland LGA Products GmbH.

Russian Classification of Products by Economic Activities (OK 034-2014): 32.50.50.190

The class of potential risk of using the medical device in accordance with the nomenclature classification of medical devices approved by the Ministry of Health of the Russian Federation of June 06, 2012, No. 4H: 2a

The type of the medical device in accordance with the nomenclature classification of medical devices approved by the Ministry of Health of the Russian Federation of June 06, 2012, No. 4H: 284910

Nature of contact with the body:

- Image processor, AC power cord, USB cable – short-term contact with intact skin.
- Cannula: JSX-T-5 cannula, JSX-T-6 cannula, JSX-T-7 cannula, JSX-T-8 cannula, JSX-T-9 cannula, JSX-T-10 cannula – short-term contact with an intact mucous membrane;

Medical use conditions: in the environment of medical and treatment facilities.

Potential consumers are patients who need vacuum aspiration. The device should be used by physicians having undergone appropriate training and having basic knowledge in the field of endoscopy, physiology, and pathology.

2. Intended use

The device is designed for transcervical aspiration of the uterine cavity used together with the third party vacuum aspiration abortion system.

3. Full name of the device

Hysteroscopy system for vacuum aspiration abortion comprised of:

- Image processor – 1 pc;
- AC power cord – 1 pc;
- USB cable – 1 pc;
- JSX-T-5 cannula (if necessary) – not more than 100 pcs;
- JSX-T-6 cannula (if necessary) – not more than 100 pcs;
- JSX-T-7 cannula (if necessary) – not more than 100 pcs;
- JSX-T-8 cannula (if necessary) – not more than 100 pcs;
- JSX-T-9 cannula (if necessary) – not more than 100 pcs;
- JSX-T-10 cannula (if necessary) – not more than 100 pcs;
- CD with software - 1 pc;
- Operator's Manual – 1 pc.

4. Operating principle of the medical device

Extraction of the fetus is carried out by means of a special cannula which has a built-in camera with a LED allowing controlling the procedure. The image processor, a cannula is connected to, processes the received photos and video recordings. Software installed to the operator's PC allows creating a patient database thus recording the procedure results.

Fetus extraction lasts for 2 – 5 minutes and is carried out under local or general anesthesia.

The device is used together with a vacuum aspiration system, which is purchased separately and is designed to create of uniform negative pressure in the cannula* ranging from 400 to 500 mm Hg.

*The cannula shall mean all the design options including JSX-T-5 cannula, JSX-T-6 cannula, JSX-T-7 cannula, JSX-T-8 cannula, JSX-T-9 cannula, JSX-T-10 cannula;

5. Appearance, device design, and designation of its components

5.1 Appearance of the medical device and its components

The device and its component parts are given in Table 1 below.

Table 1 – Appearance of the medical device and its components


Component	Image
Cannula	

Fig. 1 – General view of the cannula


Fig. 2 – Connector for vacuum aspiration system


Fig. 3 – Cannula components including the suction port, the protective head with a light source and the 8P8C connector

Image processor


Fig. 4 –Front panel of the image processor


Fig. 5 – Image processor, top view

Fig. 6 –Back panel of the image processor

Fig. 7 –Image processor, left-side view

Fig. 8 –Image processor, right-side view


Fig. 9 –Image processor, bottom view

USB cable



Fig. 10 – USB cable

5.2 Device design and designation of its components

All cannula models have the same designation and design.
Overall view of the cannula is shown in Figure 12 below.


Figure 12 – Overall view of the cannula

1 – 8P8C connector; 2 – connector for vacuum aspiration system; 3 – data cable; 4 – cannula handle; 5 – metal tube; 6 – suction port; 7 – protective head with a light source, lenses, and a built-in camera;

Designation of the cannula components is given in Table 2

Table 2 – Designation of cannula components

Cannula component	Designation
1) 8P8C connector	Designed to connect the cannula to the image processor and transmit the signal from cannula camera to the image processor
2) Connector for vacuum aspiration system	Designed to connect the cannula to a compatible vacuum aspiration system
3) Data cable	Designed to transmit the signal from cannula camera to the image processor
4) Cannula handle	Designed to insert the cannula into the uterine cavity in a guided way. The cannula body has marking applied showing the cannula insertion depth.
5) Metal tube	
6) Suction port	Designed to extract the fetus out of the uterine cavity by an action of uniform negative pressure created by the vacuum aspiration system.
7) Protective head with a light source, lenses and a built-in camera	Designed to visualize the procedure of transcervical aspiration of the uterine cavity with the help of inbuilt LED and camera.

Overall view of the image processor and its components is shown in Figure 13


Fig. 13 – Overall view of the image processor

1 – power indicator; 2 – status indicator; 3 – cannula connection port to the image processor; 4 – power button; 5 – AC power cord socket connector; 6 – fuse drawer; 7 – power ON/OFF button; 8 – marking of the image processor; 9 – RS232 port; 10 – USB port; 11 – equipotential connection;

Designation of the image processor components is given in Table 3.

Table 3 – Designation of the image processor components

Component	Designation
1) Power indicator	Shows status of the image processor connection to the power supply network. Green glowing LED means that the image processor is connected;
2) Status indicator	Shows status of the cannula connection to the image processor. Green glowing LED means that the cannula is connected;
3) Cannula connection port to the image processor	Designed to connect cannula to the image processor;
4) Power button	Designed to activate and deactivate the image processor;
5) AC power cord socket connector	Designed to connect AC power cord
6) Fuse drawer	Designed to replace fuses of the image processor
7) Main switch button	Button for manual switching ON and OFF the image processor power
8) Marking of the image processor	Shows the image processor model, the range of operating voltages, frequency, special marking symbols, information about the manufactures and its authorized representative
9) RS232 port 10) USB port	Ports for connection of the image processor to the user PC

11) Equipotential connection	In case several active devices are connected to the patient, the equipotential connection is designed to equalize their potentials.
------------------------------	---

Designation of the AC power cord, USB cable, CD with software, fuses are given in Table 4

Overall view of USB cable is shown in Figure 14


Fig. 14 – Overall view of USB cable

USB 3.0 cable is designed to connect the image processor to PC and carry out data interchange between them. The manufacturer supplies USB 3.0 cable compatible with the USB 2.0 cable connector provided in the image processor.

USB 3.0 connectors are connected to the USB port of the image processor and USB port of the user PC.

Overall view of AC power cord is shown in Figure 15


Fig. 15 – Overall view of AC power cord

AC power cord is designed to be connected to the image processor and for electrical supply.

Fuse view is shown in Figure 16


Fig. 16 – Fuse view

Fuses are designed to protect the image processor from inadmissible load currents or short circuit currents. They are supplied in a fuse drawer of the image processor.

CD with software is designed to install the program for creation of the patient database and visual control of the procedure by means of displaying the photo images and videos of the procedure.
PC software version is V 1.0 dated 01/2017

6 Weight-size and technical parameters of the medical device

6.1 Weight-size and technical parameters of the cannula

Weight-size and technical parameters of the cannula are given in Table 5 in accordance with Figure 17


Figure 17 – Overall dimensions of the cannula

Table 5.1 – Weight-size and technical parameters of the cannulas

Cannula model	L , mm	l , mm	l_1 , mm	a , mm	D_1 , mm	D , mm	Examination direction angle	Weight, g
JSX-T-5 cannula	305 ± 5	60 ± 2	20.0 ± 1	15.0 ± 1	5.0 ± 1	9.5 ± 1	$30^\circ \pm 10^\circ$	$60 \pm 5\%$
JSX-T-6 cannula	305 ± 5	60 ± 2	20.0 ± 1	15.0 ± 1	6.0 ± 1	9.5 ± 1	$30^\circ \pm 10^\circ$	$60 \pm 5\%$
JSX-T-7 cannula	305 ± 5	60 ± 2	20.0 ± 1	15.0 ± 1	7.0 ± 1	9.5 ± 1	$90^\circ \pm 10^\circ$	$60 \pm 5\%$
JSX-T-8 cannula	305 ± 5	60 ± 2	20.0 ± 1	15.0 ± 1	8.0 ± 1	9.5 ± 1	$30^\circ \pm 10^\circ$	$60 \pm 5\%$
JSX-T-9 cannula	305 ± 5	60 ± 2	20.0 ± 1	15.0 ± 1	9.0 ± 1	9.5 ± 1	$30^\circ \pm 10^\circ$	$60 \pm 5\%$
JSX-T-10 cannula	305 ± 5	60 ± 2	20.0 ± 1	15.0 ± 1	10.0 ± 1	9.5 ± 1	$30^\circ \pm 10^\circ$	$60 \pm 5\%$

The internal diameter of the cannula metal tube and the internal diameter of the connector for vacuum aspiration system are given in Table 5.2 in accordance with Figures 16.1 and 16.2.


Figure 18.1 – Showing internal diameter of the cannula metal tube


Fig. 18.2 – Showing internal diameter of the connector for vacuum aspiration system

Table 5.2 – Internal diameter of the cannula and internal diameter of the connector for vacuum aspiration system

Cannula model	D ₂ , mm	D ₃ , mm
JSX-T-5 cannula	4.4 ± 0.05	9.5 ± 0.5
JSX-T-6 cannula	5.4 ± 0.05	9.5 ± 0.5
JSX-T-7 cannula	6.4 ± 0.05	9.5 ± 0.5
JSX-T-8 cannula	7.4 ± 0.05	9.5 ± 0.5
JSX-T-9 cannula	8.4 ± 0.05	9.5 ± 0.5
JSX-T-10 cannula	9.4 ± 0.05	9.5 ± 0.5

Leakage current protection of the applied part: cannula – type BF applied part (without defibrillator discharge protection).

Electromagnetic compatibility: group I, class A.

Ingress protection rating: IPX7;

Operating mode: continuous.

Shelf-life of the cannula: 2 years.

Optical specifications of the device are given in Table 6

Table 6 – Optical specifications of the cannula

Parameter	Value
Light source	2 LEDs
Light intensity (center intensity for 5 mm), mm	≥300
Field-of-view angle (air), °, ±15%	100°
Examination depth, mm	0-3
Examination depth (air), mm	0-2
The best focal length, mm	0
Resolution (for 5 mm)	>6 line pairs/mm
Resolution (0 mm)	>20 line pairs/mm
Color rendering index (CRI) (R _a)	≥80

6.2 Weight-size and technical parameters of the image processor

Weight-size parameters of the image processor are given in Table 7 in accordance with Figure 19.


Fig. 19 – Size parameters of the image processor

Table 7 – Weight-size parameters of the image processor

Image processor			
Height (H)	Length (L)	Width (W)	Weight, kg
50 ± 2 mm	314 ± 5 mm	21 ± 1 mm	1.960 ± 5%

Electrical safety of the device: class I of electrical safety

The device is not suitable for use in the environment of a flammable mixture of anesthetics with air, oxygen or nitrogen oxide

Electromagnetic compatibility: group I, class A.

Ingress protection rating: IPX0 – Image processor.

Operating mode: continuous.

Expected service life: 5 years.

Self-life of the image processor: 5 years.

Inbuilt software – version V 1.00 dated 01/2017

BF type applied part

Supply voltage: 110-240 V ± 5 %

Supply frequency: 50/60 Hz;

Power: 25 VA;

6.3 Weight-size and technical parameters of USB cable

Maximum allowable current for the transmission through USB cable– 900 mA;

Rated voltage is 5 V ± 5 %;

Cable length – 150 ± 1 cm;

Cable external diameter – 5.5 mm ± 5 %;

Weight of USB cable – 320 g ± 5 %;

Connector type – USB 3.0 type B;

Shelf-life of USB cable – 5 years;

6.4 Weight-size and technical parameters of AC power cord

Maximum allowable current for transmission through AC power cord – 10 A;

Rated voltage is 220 V ± 0.5 %;

Power cord length is 150.5 ± 1 cm;

Number of cores, pcs – 3

Core diameter, mm² – 0.75 mm² ± 5 %;

Power cord external diameter – 7 mm ± 5 %;

Weight of AC power cord – 210 g ± 5 %;
Connector type at the image processor side – Schuko CEE 7/7;
Connector type at the power socket side – IEC 60320 C13 connector;
Shelf-life of AC power cord – 5 years;

6.5 Technical specifications of a program for PC

Video record format: AVI,
Video record resolution: 480*480 pixels;
Video streaming rate is 51 KB/sec; (3060 KB/minute);
Photo image format: BMP;
Photo image resolution: 480*480 pixels;
Size of one photo image: 675 KB.
Maximum time for the recording of one video – 12 minutes.
The software version is V 1.0 dated 01/2017.

6.6 Technical parameters of compatible vacuum aspiration systems

The medical device shall be compatible with the vacuum aspiration system providing for a pressure of 400-500 mm Hg. Besides, the vacuum aspiration system shall be securely connected to the cannula with a connector of the diameter specified in Table 5.2.

7. General indications, contraindications for use and adverse events

Indications for use:	• Patient's wish to terminate a pregnancy; • Physical disability or genetic pathology of the previously born baby; • Missed miscarriage; • Anembryonic gestation; • Childbearing is dangerous for the patient's life and health (diabetes mellitus, severe hepatic and renal pathology, malignant tumors); • The inefficiency of the previously carried out medical abortion confirmed by the uterus ultrasound investigation;
-----------------------------	--

Contraindications for use:	• Period of pregnancy of more than 10 weeks; • Infection in acute form; • Uterine malformations; • Unconfirmed pregnancy; • Labor activity within the last 6 months; • Fever (high body temperature); • Inability to tolerate operative intervention;
Adverse events (complications):	• Uterine wall damage during the procedure; • Introduction of infection during the procedure; • Menstrual disorders; • Incomplete aspiration of the gestational sac;

8. Preventive measures

- Prior to the device use, the operator shall be well knowledgeable in the principles of its operation;
- All the provisions of the device operator's manual shall be strictly followed;
- Act with due care when carrying out vacuum aspiration in order not to inflict injury to a patient;

- Negative pressure shall be lower than 500 mm Hg to prevent uterine perforation;
- The device shall be kept away from the magnetic fields;
- The device shall be kept away from the high-frequency or laser equipment.
- The device is designed to be operated only by a medical professional of appropriate competence;
- Prior to the procedure, the device shall be carefully examined for any damages;
- The image processor shall be kept away from liquids. Please consult the manufacturer in case of liquid ingress in the device;
- Prior to the procedure, an operator shall check the cannula for integrity, smoothness of the body and availability of real-time imaging;
- The cannula shall not be deformed, subject to sudden impact action or fall, the force of over 100 N shall not be applied to it, as well;
- The device can be examined and repaired only by the qualified specialists of the manufacturer;
- The cannula is disposable and not intended for re-use;
- Before disconnecting the cannula after abortion procedure first switch off the cannula power supply on the front panel;
- Prevent AC power cord and USB cable bend and stretch.

9 Connection and setting-up the medical device

9.1 Connecting the image processor to PC

1) USB cable shall be connected to appropriate connectors of the image processor and PC, and then AC power cord shall be connected to an appropriate socket connector of the image processor and the power socket as shown in Figure 20.


Fig. 20 – Connection diagram for the medical device components

2) Then power ON/OFF button located on the back panel of the image processor shall be activated and the power indicator on the front panel will start glowing green.


Fig. 21 – Connection diagram for the medical device components

9.2 Connecting the cannula to the image processor and vacuum aspiration system

Diagram for connection of the cannula to the image processor and vacuum aspiration system is shown in Figure 22


Figure 22 – Cannula connection diagram

9.3 General connection diagram for the medical device components

General connection diagram for the medical device components is shown in Figure 23


Fig. 23 – General connection diagram for the medical device components.

Vacuum aspiration system is not included in the scope of delivery and shall be purchased separately.

9.4 Installation of the software for medical device

In order to start operating the medical device, the software shall be installed to the user PC;
Requirements to the user PC for installation of the medical device software are the following:

Operating system	Windows 7, 32 or 64 bit; Windows 10, 32 or 64 bit;
Processor	Intel Core I3 or higher;
RAM	2 GB or more;
Hard disk	250 GB
Input/output ports	USB 2.0 or USB 3.0 connector

The software version is V1.0 dated 01/2017.

Microsoft.NET Framework 3.5 or higher program shall be available on PC to install the software;

The manufacturer supplies the Microsoft.NET Framework 3.5 program together with the software, and this program will be automatically installed in case it is not available on the user PC.

The following shall be done to install the software:

1) Insert a CD with the software in the drive

2) Press icon  setup.exe , then a dialog window for preparing to installation will be displayed as shown in Figure 24.


Fig. 24 Dialog window for preparing to install a program

3) Then read and accept the License Agreement (option "I accept the terms of the license agreement") and then press the Next button;


Fig. 25 Dialog window with the License Agreement

4) In the new dialog window please enter User Name and Company Name and then press the Next button;


Fig. 26 Dialog window with a User Name and Company Name

5) In the next dialog window setup type will be offered to you, please choose a Complete option for program setup and press the Next button;


Fig. 27 Dialog window with the program setup type

6) Then please press the Install button and wait till the program installation completion


Fig. 28 Dialog window for installation of the program

7) Upon installation completion, a window with the information on successful program installation will be displayed. To complete the process, please press the Finish button;


Fig. 29 Dialog window of the program installation completion

9.5 Logging in the program and patient registration

Upon successful installation of the program an icon  will be displayed on the desktop, please double-click it with a left mouse button to start the program.

1) The language shall be set up during first logging in the system as shown in Figure 30. Please choose language and press OK.


Fig. 30 Language selection dialog window

2) Enter User Name and Password and then press the Login button to log in the program. In order to remember the password and not to type it in again when re-launching the program, please choose Remember password option;


Fig. 31 Dialog window to finish the program installation

A window with the patient database is shown in Figure 32. In order to add a new patient please press the Add button.


Fig. 32 – Window with the patient database

Upon clicking Add button a form with the patient data will be displayed. Lines marked with * are obligatory for filling in. To save the data entered please press the Save button.

The screenshot shows a dialog window for adding patient data. It has a vertical sidebar on the left with an 'Add' button. The main area contains the following fields: 'In-patient ID' (text box), 'Serial Number' (text box with value '20150721115234' and an asterisk), 'Last name' (text box with an asterisk), 'First name' (text box with an asterisk), 'Birthday' (calendar icon with value '2015/ 7/21' and an asterisk), 'Profession' (text box), 'Company' (text box), 'Phone' (text box), 'PO Box' (text box), and 'Home Address' (text box). There are also labels for 'Chief Complaint' and 'On-call Physician' on the right. At the bottom right, there are 'Save' and 'Cancel' buttons. An arrow points from the 'Add' button in the sidebar to the dialog window.


Fig. 33 Dialog window with the procedure for an addition of the patient data

In order to save personal data, please press the Save button.

10 Vacuum aspiration procedure

10.1 Preparation to vacuum aspiration procedure

1) Please press the Power button on the front panel of the image processor and the status indicator will start glowing

2) Please select a patient in the program database (see Figure 34).


Fig. 34 – Dialog window for selection of a patient in the database

Then press Start procedure button to switch to a video channel.


Fig. 35 Window with the video channel

Functions of the buttons in the window with a video channel:

Camera function is used to receive an image;

Start Video is used to record a video;

Finish is used to stop the video recording;

Full screen is used to display the image on full PC screen;

Original size is used to return to the original size of the image;

Settings option is used to set up photo image and video record;

1.5 is used to display a screen in one and half of its actual size;

10.2 Methods for carrying out vacuum extraction procedure

When carrying out the procedure, at first external genitalia shall be treated as appropriate. The uterine cervix shall be treated with a disinfectant. Then the uterine cavity length is measured with help of a uterine probe.

A diameter of a cannula (see Table 5) is selected by a doctor depending on the patient's period of pregnancy.

Period of pregnancy	Cannula diameter
From 5 to 7 weeks	5 mm
From 7 to 9 weeks	6 mm
From 9 to 12 weeks	From 7 to 12 mm

Cannula size is also selected depending on the size of the gestational sac, uterus size, and uterine cervix length. In case the procedure course is slow, the doctor may select a cannula of a larger diameter.

The operation lasts from 2 to 5 minutes and is carried out under local or general anesthesia. The cannula shall be inserted slowly and horizontally turning it for 180° in each direction, slightly moving it away from and towards oneself. The suction port of the cannula shall be located behind the internal mouth of the uterus otherwise the vacuum will be lost. The video received with the help of inbuilt camera of the cannula and transmitted to the software allows controlling the procedure.

By means of vacuum aspiration system uniform negative pressure ranging from 400 to 500 mm Hg is created in the cannula, thus allowing easy separation of the gestational sac from uterine walls not depending on its location. And in this case blood vessels are damaged and bleeding occurs, blood together with the gestational sac is sucked out of the uterine cavity into the collection tank.

Pink frothing in the cannula without tissue and the uterine contractions around the cannula evidence of the completion of the procedure. Before removing the cannula out of the uterine cavity, it is necessary to switch off the power supply of the image processor by pressing power ON/OFF button on the front panel of the image processor.

Actions to be done upon vacuum aspiration:

- 1) To disconnect and discard the cannula;
- 2) To close the window with a video;
- 3) To fill in and print out the report form (see section 10.3);

101 Working with the application

11.1 Personal setup

Please press Personal Setup option in the main program interface to edit information about the current user

The screenshot shows the 'Uterus Visualization System' main interface. At the top, there is a header bar with the system name, version (1.1), and navigation buttons: 'Personal Setup' (selected), 'User', and 'Start P'. Below the header, the 'Personal Setup' window is open, displaying two sections: 'User Information Setup' and 'Security Setup'. In the 'User Information Setup' section, the 'Login Name' is set to 'admin', and there is a 'Detail' button with a left arrow. The 'Security Setup' section is divided into 'Change Password' and 'Login Setup'. Under 'Change Password', there are three input fields for 'Original Password', 'New Password', and 'Confirm Password', followed by a 'Save' button. Under 'Login Setup', there is a checked checkbox for 'Remember Password'.

Press the Save button upon completion of the information editing process.

11.2 User information setup

Please press User option in the main interface to add, edit or delete the user.

User Information

Basic Information

Login Name: admin * Password: 100000

Name: admin * Age: 42

Gender: ☐ Male ☒ Female Telephone: 1388888888

Detailed Information

Number: 100001

Department: Gynecology * Procedural Physician *

Profession: Pharmacy Healer

Work Age: 9

Permission: ☒ User Management ☒ Procedural Management ☒ Reporting Management ☒ Hospital Management *

Individual resume

☒ Save

11.3 Report setup

Please select those points in the Report Setup window to be printed out.

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☐ Uterus Position
☐ Cervical Depth (Pre-Pos) / cm

☐ Cervical Size (Pre-Pos) / cm
 ☐ Suction Tube Size

☐ Negative Pressure
☐ Aspirates Content

☐ Villus
☐ Embryonic Sac

☐ Size of Aspirated Embryonic Sac mm
 ☐ Blood Loss ml

☐ Curettage
☐ Treatment

☐ Drug
☐ Menstrual Cycle

☐ Special Notes

☐ Contraceptive Device Placement
 ☐ Size mm
 ☐ Type

☐ Duration of Placement
☐ Tail

☐ Length of the tail mm
 ☐ Model

☐ Specification

☐ Others

Procedure Picture ☒

☐ ☐ ☐ ☐ ☐ ☐

Post-op Notes ☒

Fig. 36 Window with report data

Then press the Save button upon completion of the editing process.

11.4 Hospital information setup.

Please choose Facility setup option in the main program interface to log in the Hospital information interface.

Then press Save to save the information entered.

11.5 Working with the database

1) Creating backup copies of the database.

The user can transfer data with help of data management function.


A required option of the drop-down menu shall be selected and then Backup Data button shall be pressed to create a backup copy of all the patient, user, and hospital information.

Upon completion of the backup copying process, all the existing data will be deleted.


To restore the data, please choose a location for restored data and then press Restore Data button to start the restoration.

2) Database transfer

Database transfer is necessary when the current database has reached its limits. This action will result in a transfer of all the images and videos from the current database to the specified folder and the database will be cleared as a consequence.


The user can add new data by pressing the Open button

11.6 Searching for the patients

Please select Item option for quick search for the patient. For instance, select Name option in the drop-down menu. Content field allows carrying out the search by the words entered upon pressing the Search button.

11.7 Access to the patient information

If double-click the patient option in the main program interface you can access the patient information. General inspection: please record the results of the general inspection the patient has already been subject to.

Dialog window Procedure Status allows recording information on the procedures carried out.

Detailed Patient Information

Basic Information General Inspection **Procedure Status** Procedure Picture

Date of Procedure: Attending Physician:

Uterus Position: Cervical Depth (Pre Pos): / cm Cervical Size (Pre Pos): / cm

Suction Tube Size: Negative Pressure: Aspires Content:

Vitrus: Embryonic Sac: Size of Aspirated Embryonic Sac: mm

Blood Loss: ml Curettage: Treatment:

Drug: Menstrual Cycle: Special Notes:


Contraceptive Device Placement: Choose 'Yes' to unlock the following tests

Size: mm Type: Duration of Placement:

Tail: Length of the tail: mm Model:


Specification: Others:

Dialogue window with the procedure displayed.


11.8 Viewing video and images

In order to view all the images and video done for a definite patient, please double-click the Patient information option in the main window and then press the Folder button.


11.9 Quick view of the patient video

In order to quickly view the video, please double-click the option with the patient name in the main window and then press the View button. Video playback is started with a double-click.


11.10 Deleting the patient information

Please press Delete in the program interface to delete the information about the selected patient.


11.11 Uninstalling the application


The program can be uninstalled from the list of the programs installed on PC.


You may also enter the folder with the program installed, and double-click uninstall icon.

Then you shall confirm the program uninstallation from PC.

Upon completion of all the actions specified a note about successful uninstallation of the program will be displayed

11 Information on the device materials

The device materials are given in Table 8

Table 8 – Materials of the device and its components

Device or component and accessory (under test)	Materials used during manufacture	Nature of contact with the body
Image processor:		Short-term contact with intact skin
- body	Stainless steel Wuxi Daya AISI 304 with applied white powder polymer coating Wuxi Daya WH1008 manufactured by Wuxi Daya Special Steel Co. Ltd., China	
- power button	Stainless steel Wuxi Daya AISI 304 manufactured by Wuxi Daya Special Steel Co. Ltd., China	
- AC power cord connector, fuse drawer cover, a main switch button	ABS plastic Fukuang ABS9010 with a black colorant Fukuang Black10 manufactured by Fukuang Plastic Technology Co.,Ltd., China	
- equipotential connection	Brass Nickel Plated 08 manufactured by Taizhou Jiaxin Metals Co., Ltd., China	
AC power cord	PVC Zhengzhou SPVC 64TS with a black colorant Zhengzhou Black233 manufactured by Zhengzhou San Techchem Co., Ltd., China	
USB cable	PVC Zhengzhou SPVC 60S with a blue colorant Zhengzhou Blue122 manufactured by Zhengzhou San Techchem Co., Ltd., China	

JSX-T-5 cannula, JSX-T-6 cannula, JSX-T-7 cannula, JSX-T-8 cannula, JSX-T-9 cannula, JSX-T-10 cannula:		Short-term contact with an intact mucous membrane
- Protective head with a light source	Polycarbonate Guangzhou Yuemei PC2458, transparent manufactured by Guangzhou Yuemei Technology Materials Co., Ltd., China	
- Cannula body	Steel of Shenzhen Kingship SUS316 grade manufactured by Shenzhen Kingship Machinery & Electronic Co., Ltd., China	
- Cannula handle	ABS plastic Fukuang ABS1104 with a white colorant Fukuang White04 and red colorant Fukuang Red04 manufactured by Fukuang Plastic Technology Co.,Ltd., China;	
Connector for vacuum aspiration system	Polyamide Zhangjiagang Huachuang PA66 with a white colorant Zhangjiagang Huachuang White008 manufactured by Zhangjiagang Huachuang Engineering Plastics Co., Ltd., China	
- Cable for connection to the processor	PVC Zhengzhou SPVC 60TS with a white colorant Zhengzhou White44 and black colorant Zhengzhou Black84 manufactured by Zhengzhou San Techchem Co., Ltd., China	
Individual package for all the models of hysteroscopy cannulas:		Short-term contact with intact skin
Lower part	Non-woven fabric Tyvek 1073B manufactured by DuPont, USA	
Upper part	Polyester /2 mil polyethylene; manufactured by Beacon Converters, Inc., USA	

12. Package and marking

The medical device is supplied in separate carton boxes for the cannulas and the image processor.

Complete set of the transportation package for the image processor is as follows:

- Image processor – 1 pc
- AC power cord – 1 pc
- USB cable – 1 pc
- CD with software – 1 pc
- Operator's Manual – 1 pc

10 multiple packages are included in the transportation package. Multiple packages include 10 cannulas in individual packages of the same design.

Each item of the supply set is individually packed in a polyethylene bag.

12.1 Marking of the image processor

Marking of the image processor contains the following items:

- Device name;
- Device component;
- Supply voltage and frequency;
- Rated power;
- Electrical safety class;

Including the following symbols and conventional signs:

-  Serial number with the specification of the device serial number;
-  Manufacture date with the specification of the manufacture date;
-  Best before (with the specification of the year and month);
-  CE certification mark with the specification of the notified body number;
-  WEEE Directive;
-  BF type applied parts;
-  Caution! Refer to the Operator's Manual;
-  Refer to the Operator's Manual;

Marking of the medical device is shown in Figure 16.


Figure 16 – Marking of the image processor

12.2 Marking of the cannula

Marking of the cannula contains the following items:

- Name of the device component;
- Manufacturer's name;
- Electrical equipment protection class IPX7;

Including the following symbols and conventional signs:

-  WEEE Directive;
-  BF type applied parts;
-  Caution! Refer to the Operator's Manual;
-  For single use only
-  Ethylene oxide sterilized
-  Refer to the Operator's Manual;
-  CE certification mark with the specification of the notified body number;

Photographs of the marking of a cannula for vacuum aspiration are given in Figures 18.1 and 18.2. Marking of all the other cannula models is the same.

Fig. 18.1 – Photograph of JSX-T-7 cannula marking


Fig. 18.2 – Photograph of JSX-T-7 cannula marking


12.3 Marking of the individual cannula package


The cannula is supplied in a branded manufacturer's package protecting it from climatic and mechanical impacts. Each item of the supply set is individually packed in a polyethylene bag.

Marking of the individual cannula package contains the following items:

- Device name;
- Device component;
- Number of device components;
- The diameter of the tube;
- Electrical equipment protection class IPX7;
- Information on the device sterility;
- Warning "Do not use in case of the damaged package"
- Information on non-toxicity and non-pyrogenicity of the medical device;

Including the following symbols and conventional signs:

-  Information on the manufacturer (with the specification of the manufacturer's name and address);
-  Serial number with the specification of the device serial number;
-  Lot code;
-  WEEE Directive;
-  BF type applied parts;
-  Caution! Refer to the Operator's Manual;
-  For single use only
-  Ethylene oxide sterilized

-  Refer to the Operator's Manual;
-  CE certification mark with the specification of the notified body number;
-  Manufacture date with the specification of the manufacture date;
-  Best before (with the specification of the year and month);

Photograph of JSX-T-7 cannula package marking is given in Fig. 17
Marking of all the other cannula models is the same.


Fig. 17 – Photograph of JSX-T-7 cannula individual package marking

12.4 Marking of the multiple cannula packages

Marking of the multiple cannula packages contains the following items:

- Device name;
- Device component;
- Number of device components;
- Electrical equipment protection class IPX7;
- Information on the device sterility;
- Warning "Do not use in case of the damaged package"
- Information on non-toxicity and non-pyrogenicity of the medical device;
- Number of the device components in the package;
- Transportation and storage conditions;

Including the following symbols and conventional signs:

-  Information on the manufacturer (with the specification of the manufacturer's name and address);
-  Serial number with the specification of the device serial number;
-  Lot code;
-  WEEE Directive;
-  BF type applied parts;
-  Caution! Refer to the Operator's Manual;
-  For single use only


-  Ethylene oxide sterilized
-  Refer to the Operator's Manual;
-  CE certification mark with the specification of the notified body number;
-  Manufacture date with the specification of the manufacture date;
-  Best before (with the specification of the year and month);

Photograph of the marking for multiple cannula packages is shown in Figure 19.

Device:


Hysteroscopy system for vacuum aspiration abortion

Device component:

- | | | |
|--|--|---|
| ☐ JSX-T-5 cannula | ☐ JSX-T-6 cannula | ☒ JSX-T-7 cannula |
| ☐ JSX-T-8 cannula | ☐ JSX-T-9 cannula | ☐ JSX-T-10 cannula |

- The device is sterile, do not to use in case of the damaged package;
- Non-pyrogenic;
- Non-toxic

Number: 10 pcs

Chongqing Jinshan Science & Technology (Group) Co., Ltd.
China
No. 18, Nishang Road, LiangLu Industrial City,
401120 Yubei District, Chongqing, China

Transportation and storage conditions:

To keep at air temperature ranging from +10 to +35°C and relative humidity of not more than 80% and atmosphere pressure ranging from 86.6 KPa to 106.7 KPa



Fig. 19 – Marking of the multiple cannula packages.

12.5 Marking of the transportation cannula package


Marking of the transportation cannula package contains the following items:

- Device name;
- Device component;
- Transportation package number;
- Electrical equipment protection class IPX7;
- Information on the device sterility;
- Transportation and storage conditions;

Including the following symbols and conventional signs:

-  Information on the manufacturer (with the specification of the manufacturer's name and address);
-  WEEE Directive;
-  BF type applied parts;
-  Caution! Refer to the Operator's Manual;
-  For single use only
-  Ethylene oxide sterilized
-  Refer to the Operator's Manual;
-  CE certification mark with the specification of the notified body number;
-  Manufacture date with the specification of the manufacture date;
-  Best before (with the specification of the year and month);

Layout for marking of the transportation cannula package:

Device: Hysteroscopy system for vacuum aspiration abortion	
Device component: ☐ JSX-T-5 cannula ☐ JSX-T-6 cannula ☐ JSX-T-7 cannula ☐ JSX-T-8 cannula ☐ JSX-T-9 cannula ☐ JSX-T-10 cannula	
The device is sterile	
Transportation package number: 1 of 2;	
     IPX7   0197	
	Chongqing Jinshan Science & Technology (Group) Co., Ltd., China No. 18, Nishang Road, LiangLu Industrial City, 401120 Yubei District, Chongqing, China http://english.jinshangroup.com Phone: +86-23-86098099
 03.2017	 03.2019
	+35 C° Transportation and storage conditions: To keep at air temperature ranging from +10 to +35°C and relative humidity of not more than 80% and atmosphere pressure ranging from 86.6 KPa to 106.7 KPa

12.6 Marking of the transportation package for image processor

Marking of the transportation package for image processor contains the following items:

- Device name;
- Device component;
- Supply voltage and frequency;
- Rated power;
- Electrical safety class;

Including the following symbols and conventional signs:

-  Serial number with the specification of the device serial number;
-  Manufacture date with the specification of the manufacture date;
-  Best before (with the specification of the year and month);
-  CE certification mark with the specification of the notified body number;
-  WEEE Directive;
-  BF type applied parts;
-  Caution! Refer to the Operator's Manual;
-  Refer to the Operator's Manual;

Layout for marking of the transportation package for image processor:

Device: Hysteroscopy system for vacuum aspiration abortion	
Device component: Image processor	
Transportation package number: 2 of 2;	
Electrical safety – Class I.	
Delivery scope: - Image processor – 1 pc - AC power cord – 1 pc - USB cable – 1 pc - CD with software – 1 pc - Operator's Manual – 1 pc	    
Supply voltage: 110-240 V; Supply frequency: 50/60 Hz; Rated power: 25 VA;	


Chongqing Jinshan Science & Technology (Group) Co., Ltd., China
No. 18, Nishang Road, LiangLu Industrial City,
401120 Yubei District, Chongqing, China
<http://english.jinshangroup.com> Phone: +86-23-86098099

SN

XYG046427029

10.01.2017

+35 C° Transportation conditions:
To keep at air temperature
ranging from +10 to +35°C and
relative humidity of not more
than 80% and atmosphere
pressure ranging from 86.6 KPa
to 106.7 KPa

10.01.2022

13 Compliance with the Russian Federation standards

GOST 31214-2016. "Medical devices. Requirements for samples and documentation presented for toxicological tests, sanitary and chemical analyses, tests for sterility and pyrogenicity".

GOST R 52770-2016. "Medical devices. Safety requirements. Methods of sanitation-chemical and toxicological tests".

GOST ISO 10993-1-2011. "Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing".

GOST R ISO 10993-2-2009. "Medical devices. Biological evaluation of medical devices. Part 2. Animal welfare requirements".

GOST ISO 10993-5-2011. "Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity".

GOST ISO 10993-7-2016. "Medical devices. Evaluation of the biological effects of the medical devices. Part 7. The residual content of ethylene oxide after sterilization".

GOST ISO 10993-10-2011. "Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and sensitization".

GOST ISO 10993-12-2015. "Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials".

The State Pharmacopoeia of the Russian Federation, XIV Edition. General monograph 1.2.4.0003.15. Sterility.

GOST 31209-2003 "Containers for blood and its components. Requirements for chemical and biological safety and methods of testing".

MUK 4.1.742-99 "Inversion voltammetric measurement of zinc, cadmium, lead, and copper ion concentration in the water".

GOST 4011-72. Drinking water. Methods for determination of total iron (with Amendments No. 1, 2)

MUK 4.1.3166-14 "Gas chromatography for determination of hexane, heptane, acetaldehyde, acetone, methyl acetate, ethyl acetate, methanol, isopropanol, acrylonitrile, n-propanol, n-propyl acetate, butyl acetate, isobutanol, n-butanol, benzene, toluene, ethylbenzene, m-, o- and p-xylenes, isopropylbenzene, styrene, α -methyl styrene content in the water and water extracts of materials of various composition".

I-880-71 "Instructions on sanitary and chemical examination of the devices manufactured from polymer and other synthetic materials contacting with food products".

MR 1436-76 "Instructional guidelines on a determination of diphenylol propane and some phenols in its presence in the course of sanitary and chemical tests of the devices manufactured from polymer materials contacting with food products".

GOST 31866-2012 "Drinking water. Determination of the content of elements by stripping voltammetry method".

GOST 31870-2012 "Drinking water. Determination of the element content by atomic spectrometry (as amended)".

GN 2.3.3.972-00 "Maximum allowable quantities of chemicals released from materials in contact with food".

MUK 4.1.3169-14 "Gas chromatography determination of dimethyl phthalate, dimethyl terephthalate, diethyl phthalate, dibutyl phthalate, butyl benzyl phthalate, bis(2-ethylhexyl)phthalate and dioctyl phthalate content in the water and water extracts of various composition".

MR 1503-76 "Instructional guidelines on a determination of hexamethylenediamine content in the water in the course of sanitary and chemical tests of polymer materials used in the food and textile industry".

GOST R 50444-92 "Medical instruments, apparatus and equipment. General specifications".

GOST R IEC 60601-1-2010 "Medical electrical equipment. Part 1. General requirements for basic safety and essential performance".

GOST R ISO 15223-1-2014 "Medical devices. Symbols to be used with medical device labels, labeling, and information to be supplied. Part 1. General requirements".

GOST R IEC 62304-2013 "Medical device software. Software life cycle processes".

GOST R ISO/IEC 12119-2000 "Information technology. Software packages. Quality requirements and testing".

GOST R ISO/IEC 9126-93 "Information technology. Software product evaluation. Quality characteristics and guidelines for their use".

GOST R ISO 9127-94 "Information processing systems. User documentation and cover information for consumer software packages".

GOST 28195-89 "Quality control of software systems. General principles".

GOST R 51188-98 "Information security. Software testing for the presence of computer viruses. The sample manual".

GOST R IEC 60601-2-18-2014 "Medical electrical equipment. Part 2-18. Particular requirements for the basic safety and essential performance of endoscopic equipment".

GOST R 53469-2009 "Optics and optical instruments. Medical endoscopes and endotherapy devices. Part 1. General requirements".

GOST 26332-84 "Medical endoscopes. Marking, packing, transportation and storage".

GOST R 55037-2012 "Optics and optical instruments. Medical endoscopes and endotherapy devices. Particular requirements. Test methods".

GOST R IEC 60601-1-2-2014 "Medical electrical equipment. Part 1-2. General requirements for basic safety and essential performance. Collateral standard. Electromagnetic compatibility. Requirements and tests".

GOST ISO 11607-1-2018 "Packaging for medical devices subject to final sterilization. Part 1. Requirements for materials, sterilization barrier systems and packaging systems".

14 Control measures according to the Russian Federation standards

GOST 31214-2016. "Medical devices. Requirements for samples and documentation presented for toxicological tests, sanitary and chemical analyses, tests for sterility and pyrogenicity".

GOST R 52770-2016. "Medical devices. Safety requirements. Methods of sanitation-chemical and toxicological tests".

GOST ISO 10993-1-2011. "Medical devices. Biological evaluation of medical devices. Part 1. Evaluation and testing".

GOST R ISO 10993-2-2009. "Medical devices. Biological evaluation of medical devices. Part 2. Animal welfare requirements".

GOST ISO 10993-5-2011. "Medical devices. Biological evaluation of medical devices. Part 5. Tests for in vitro cytotoxicity".

GOST ISO 10993-7-2016. "Medical devices. Evaluation of the biological effects of the medical devices. Part 7. The residual content of ethylene oxide after sterilization".

GOST ISO 10993-10-2011. "Medical devices. Biological evaluation of medical devices. Part 10. Tests for irritation and sensitization".

GOST ISO 10993-12-2015. "Medical devices. Biological evaluation of medical devices. Part 12. Sample preparation and reference materials".

The State Pharmacopoeia of the Russian Federation, XIV Edition. General monograph 1.2.4.0003.15. Sterility.

GOST 31209-2003 "Containers for blood and its components. Requirements for chemical and biological safety and methods of testing".

MUK 4.1.742-99 "Inversion voltammetric measurement of zinc, cadmium, lead, and copper ion concentration in the water".

GOST 4011-72. Drinking water. Methods for determination of total iron (with Amendments No. 1, 2)

MUK 4.1.3166-14 "Gas chromatography for determination of hexane, heptane, acetaldehyde, acetone, methyl acetate, ethyl acetate, methanol, isopropanol, acrylonitrile, n-propanol, n-propyl acetate, butyl acetate, isobutanol, n-butanol, benzene, toluene, ethylbenzene, m-, o- and p-xylenes, isopropylbenzene, styrene, α -methyl styrene content in the water and water extracts of materials of various composition".

I-880-71 "Instructions on sanitary and chemical examination of the devices manufactured from polymer and other synthetic materials contacting with food products".

MR 1436-76 "Instructional guidelines on a determination of diphenylol propane and some phenols in its presence in the course of sanitary and chemical tests of the devices manufactured from polymer materials contacting with food products".

GOST 31866-2012 "Drinking water. Determination of the content of elements by stripping voltammetry method".

GOST 31870-2012 "Drinking water. Determination of the element content by atomic spectrometry (as amended)".

GN 2.3.3.972-00 "Maximum allowable quantities of chemicals released from materials in contact with food".

MUK 4.1.3169-14 "Gas chromatography determination of dimethyl phthalate, dimethyl terephthalate, diethyl phthalate, dibutyl phthalate, butyl benzyl phthalate, bis(2-ethylhexyl)phthalate and dioctyl phthalate content in the water and water extracts of various composition".

MR 1503-76 "Instructional guidelines on a determination of hexamethylenediamine content in the water in the course of sanitary and chemical tests of polymer materials used in the food and textile industry".

GOST R 50444-92 "Medical instruments, apparatus and equipment. General specifications".

GOST R IEC 60601-1-2010 "Medical electrical equipment. Part 1. General requirements for basic safety and essential performance".

GOST R ISO 15223-1-2014 "Medical devices. Symbols to be used with medical device labels, labeling, and information to be supplied. Part 1. General requirements"

GOST R IEC 62304-2013 "Medical device software. Software life cycle processes"

GOST R ISO/IEC 12119-2000 "Information technology. Software packages. Quality requirements and testing";

GOST R ISO/IEC 9126-93 "Information technology. Software product evaluation. Quality characteristics and guidelines for their use"

GOST R ISO 9127-94 "Information processing systems. User documentation and cover information for consumer software packages"

GOST 28195-89 "Quality control of software systems. General principles".

GOST R 51188-98 "Information security. Software testing for the presence of computer viruses. The sample manual"

GOST R IEC 60601-2-18-2014 "Medical electrical equipment. Part 2-18. Particular requirements for the basic safety and essential performance of endoscopic equipment"

GOST R 53469-2009 "Optics and optical instruments. Medical endoscopes and endotherapy devices. Part 1. General requirements"

GOST 26332-84 "Medical endoscopes. Marking, packing, transportation and storage"

GOST R 55037-2012 "Optics and optical instruments. Medical endoscopes and endotherapy devices. Particular requirements. Test methods"

GOST R IEC 60601-1-2-2014 "Medical electrical equipment. Part 1-2. General requirements for basic safety and essential performance. Collateral standard. Electromagnetic compatibility. Requirements and tests".

GOST ISO 11607-1-2018 "Packaging for medical devices subject to final sterilization. Part 1. Requirements for materials, sterilization barrier systems and packaging systems".

15 Transportation, storage, operation

The medical device can be transported by all the transport means in covered vehicles in accordance with the shipping rules valid for these transport means under the following conditions:

Transportation conditions: air temperature ranging from +10 to +35°C;
relative humidity of not more than 80%;
atmosphere pressure ranging from 86.6 KPa to 106.7 KPa;

Operating and storage conditions for the medical device:

Storage conditions: air temperature ranging from +10 to +35°C;
relative humidity of not more than 80%;
atmosphere pressure ranging from 86.6 KPa to 106.7 KPa;

Operating conditions for the device (with exception of an optical part of the cannula): air temperature ranging from +10 to +35°C;
relative humidity of not more than 80%;
atmosphere pressure ranging from 86.6 KPa to 106.7 KPa;

Operating conditions for an optical part of the cannula: air temperature ranging from +10 to +35°C;
relative humidity of not more than 80%;
atmosphere pressure ranging from 86.6 KPa to 106.7 KPa;

16 Information on sterility and device cleaning

16.1 Sterility of the cannula

The cannula is supplied sterile. Sterilization with ethylene oxide is carried out in accordance with GOST ISO 11135-2012. Guaranteed sterility level is 10^{-6} . The cannula is for single use only, re-sterilization and re-use is forbidden.

16.2 Cleaning

In order to clean the image processor, AC power cord, and USB cable, do the following:

- use dry gauze or the other cloth to remove slight dirt accumulations.
- use gauze and neutral cleaning agent to remove severe dirt accumulations.

Cleaning of the image processor not complying with the specified recommendations may lead to water ingress in electronic components of the device and its damage.

16.3 Disinfection

The surface of the image processor, AC power cord, and USB cable may be cleaned with some alcohol (ethyl alcohol of not less than 70%).

17 Information on repair and maintenance of the medical device

17.1 Information on the medical device repair

Repair, replacement of the medical device and its components shall be carried out only by the manufacturer's authorized personnel. In case of the device failure, please contact the manufacturer using the following details:

China
Chongqing Jinshan Science & Technology (Group) Co., Ltd., China
Address: No. 18, Nishang Road, LiangLu Industrial City, 401120 Yubei District, Chongqing,
Phone: +86-23-86098099;
E-mail: international@jinshangroup.com;

17.2 Information on the medical device maintenance

17.2.1 Periodic checks

It is recommended to check the medical device function:

- every 3 months;
- every 6 months and in case the device has not been operated for 6 months, it is necessary to check the main functions of the device with a cannula connected to the image processor;
- annually;

The medical device function shall be checked only by the manufacturer's authorized personnel. The annual check shall include the device integrity, grounding, electronic component, and device function inspection.

To carry out periodic check, repair, and replacement of the medical device, please contact the manufacturer:

China
Chongqing Jinshan Science & Technology (Group) Co., Ltd., China
Address: No. 18, Nishang Road, LiangLu Industrial City, 401120 Yubei District, Chongqing,
Phone: +86-23-86098099;
E-mail: international@jinshangroup.com;

17.2.2 Replacement of the fuses

Please contact the manufacturer to have fuses replaced.

18 Information on disposal

Disposal of the medical device shall be carried out in accordance with the Sanitary Rules and Regulations on class B wastes (epidemiologically hazardous wastes) SanPiN 2.1.7.2790-10. After application of the instrumental decontamination methods with the use of physical methods and changing appearance of the wastes excluding the possibility of their re-use, class B wastes can be accumulated, temporarily stored, transported, disposed and dumped together with class A wastes (epidemiologically safe wastes similar to solid household wastes in their composition). Package of the decontaminated class B medical wastes shall have marking proving their decontamination.

19. Guarantees

This guarantee provided by manufacturer Chongqing Jinshan Science & Technology (Group) Co., Ltd. covers the medical device "Hysteroscopy system for vacuum aspiration abortion" as to its repair, replacement or replacement of its components. A faulty device shall be returned to the manufacturer or its authorized distributor within the specified service life of the device. This guarantee shall not cover the device or its components in case they have been repaired by the persons other than the manufacturer's authorized personnel and if such repair could impact the device stability, safety and lead to its improper use or an accident.

This guarantee shall not be transferred to and shall not cover the original purchaser of this medical device. Chongqing Jinshan Science & Technology (Group) Co., Ltd. disclaims all the responsibilities connected with the product sale and any other indirect violations of law or consequential losses.

Self-life of the image processor is 5 years;
Self-life of the cannulas is 2 years;

Manufacturer: Chongqing Jinshan Science & Technology (Group) Co., Ltd.,
China
Address: No. 18, Nishang Road, LiangLu Industrial City, 401120 Yubei
District, Chongqing, China
Phone: +86-23-86098099
E-mail: international@jinshangroup.com

The authorized representative of the
manufacturer at the territory of the
Russian Federation: BW Company, LLC

Address: 129085, Moscow, 101 Mira Avenue, bld. 1, premise 17
Phone: +7 (499) 281-67-68
E-mail: info@beawire.com

20. Information on the electromagnetic compatibility of the hysteroscopy system for vacuum abortion according to GOST R IEC 60601-1-2-2014

Table 1 - Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emission		
The hysteroscopy system for vacuum abortion * is intended for use in the electromagnetic environment defined below. The buyer or user of the hysteroscopy system for vacuum-aspiration abortion should ensure that the working environment complies with these characteristics. * Further in the text - hysteroscopy system		
Radiation test	Compatibility	Electromagnetic environment guidance
Radio interference by CISPR 11	Group 1	The hysteroscopy system uses radio frequency energy only for its internal functions, so its radio emission is small and most likely will not create interference for other equipment. The hysteroscopy system is suitable for use in all locations, except for residential buildings and buildings directly connected to the distribution network that feeds residential buildings.
Radio interference by CISPR 11	Class A	
Harmonic radiation IEC 61000-3-2	Not applicable	
Voltage / Radiation Fluctuations IEC 61000-3-3	Not applicable	

Table 2 - Guidance and Manufacturer's Declaration - Immunity

Guidance and Manufacturer's Declaration - Immunity
The hysteroscopy system is intended for use in the electromagnetic environment defined below. The customer or the user of the hysteroscopy system must ensure that the working environment complies with these

characteristics.			
Noise immunity test	Test level according to IEC 60601	Compliance Level	Electromagnetic environment - instructions
Electrostatic discharge (ESD) according to IEC 61000-4-2	± 6 kV contact discharge ± 8 kV air discharge	± 6 kV contact discharge ± 8 kV air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity of the air should be at least 30%.
Nanosecond pulse interference according to IEC 61000-4-4	± 2 kV - for power lines ± 1 kV - for input / output lines	± 2 kV - for power lines ± 1 kV - for input / output lines	The quality of the electrical network must comply with the requirements for electrical networks of medical institutions.
High energy microsecond pulse interference according to IEC 61000-4-5	± 1 kV - when applying interference under the "wire-to-wire" scheme ± 2 kV - in case of wire-to-ground interference	± 1 kV - when applying interference under the "wire-to-wire" scheme ± 2 kV - in case of wire-to-ground interference	The quality of electrical energy in the electrical network of a building must meet the typical conditions of a commercial or hospital environment.
Dips, interruptions and voltage changes in the input power supply lines according to IEC 61000-4-11	$< 5\% U_T$ (voltage dip $> 95\% U_T$) within 0,5 period	$< 5\% U_T$ (voltage dip $> 95\% U_T$) within 0,5 period	The quality of electrical energy in the electrical network of a building must meet the typical conditions of a commercial or hospital environment. If the user of the hysteroscopy system needs continuous operation in the conditions of mains voltage interruptions, it is recommended to provide power to the hysteroscopy system from the battery or uninterruptible power supply.
	$40\% U_T$ (voltage dip $60\% U_T$) during 5 periods	$40\% U_T$ (voltage dip $60\% U_T$) during 5 periods	
	$70\% U_T$ (voltage dip $30\% U_T$) during 25 periods	$70\% U_T$ (voltage dip $30\% U_T$) during 25 periods	
	$< 5\% U_T$ (voltage dip $> 95\% U_T$) during 5 seconds	$< 5\% U_T$ (voltage dip $> 95\% U_T$) during 5 seconds	
Magnetic field of industrial frequency according to IEC 61000-4-8	3 A/m	3 A/m	The levels of the magnetic field of industrial frequency should be provided in accordance with the typical conditions of a commercial or hospital environment.
Note: U_T voltage level of the electrical network until the moment of application of the test exposure			

Table 4 – Guidance and Manufacturer's Declaration - Immunity

Guidance and manufacturer's declaration - protection against electromagnetic fields			
The hysteroscopy system is intended for use in the electromagnetic environment defined below. The customer or the user of the hysteroscopy system must ensure that the working environment complies with these characteristics.			
Noise immunity test	Test level according to IEC	Compliance Level	Electromagnetic environment - instructions

	60601		
Conductive interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6	3 V (root mean square value) 150 kHz – 80 MHz	3 V (root mean square value)	Portable and mobile radio equipment should be used no closer than a certain distance to any part of the hysteroscopy system, including cables. This distance is calculated using the equation applicable to the transmitter frequency. Recommended distance: $d = 1,2\sqrt{P}$
Radiated radio frequency electromagnetic field according to IEC 61000-4-3	3 V/m in the strip from 80 MHz to 2.5 GHz	3 V/m	$d = 1,2\sqrt{P}$ (from 80 MHz to 800 MHz); $d = 2,3\sqrt{P}$ (from 800 MHz to 2,5 GHz), where P is the maximum transmitter output power in watts (W) according to the manufacturer, and d is the recommended distance in meters (m). The electric field strength from stationary radio transmitters, as determined by the survey of the electromagnetic section, must be less than the level of compliance in each frequency range. There may be interference in the vicinity of equipment marked with the following symbol.  The field strength during propagation of radio waves from stationary radio transmitters, according to the results of observations of the electromagnetic environment a), should be lower than the level of compliance in each frequency band b).
a) The field strength in propagation of radio waves from stationary radio transmitters, such as base stations of radio telephone networks (cellular / wireless) and land mobile radio stations, amateur radio stations, AM and FM broadcast transmitters, television transmitters, cannot be determined by calculation with sufficient accuracy. To do this, practical measurements of field strength must be made. If the measured values at the location of the hysteroscopy system are more than the applicable levels of compliance, then the operation of the hysteroscopy system should be monitored to verify its normal operation. If in the process of observation a deviation from normal functioning is detected, then additional measures should be taken, such as reorienting or moving the hysteroscopy system. b) Outside the frequency range from 150 kHz to 80 MHz, field strengths of less than 1 V / m should be ensured. Notes: 1. At 80 and 800 MHz, a higher field strength value is used. 2. These expressions are not applicable in all cases. Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.			

Table 6 - Recommended spacing values between portable and mobile radio frequency communications and model 006

Recommended values of the separation distance between portable and mobile radio frequency communications and hysteroscopy system

The hysteroscopy system is intended for use in an electromagnetic environment in which emitted interference levels are monitored. The customer or the user of the hysteroscopy system can avoid the effects of electromagnetic interference by providing minimal spatial separation between portable and mobile radio frequency communications (transmitters) and the hysteroscopy system, as recommended below, given the maximum output power of the communications.

Nominal maximum output power, W	Spatial separation, m, depending on the frequency of the transmitter		
	$d = 1.2\sqrt{P}$ in the strip from 150 kHz to 80 MHz	$d = 1.2\sqrt{P}$ in the frequency band from 80 to 800 MHz	$d = 2.3\sqrt{P}$ in the 800 MHz to 2.5 GHz band
0.12	0.12	0.12	0.23
0.38	0.38	0.38	0.73
1.2	1.2	1.2	2.3
3.8	3.8	3.8	7.3
12	12	12	23

Notes:

- 1 At frequencies of 80 and 800 MHz, a higher field strength value is used.
- 2 These expressions are not applicable in all cases. Electromagnetic propagation is affected by absorption or reflection from structures, objects and people.
- 3 When determining the recommended values of the spatial separation d for transmitters with a rated maximum output power not specified in the table, the nominal maximum output power P in watts specified in the documentation of the transmitter manufacturer is substituted into the above expressions.

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

The screenshot shows a software interface with a top menu bar containing 'User', 'Start Procedure', 'Report Setup', 'Facility Setup', and 'Data'. Below the menu is a search bar with a dropdown set to 'All' and a 'Search' button. The 'Data Management' window is open, displaying a 'Data Management' section with a 'Backup Savepath' dropdown and a 'Backup Data' button. An arrow points from the 'Data' menu item to the 'Data Management' window, and another arrow points from the 'Backup Data' button to a 'Notice' dialog box. The 'Notice' dialog box contains the text: 'Are you sure to archive the data to the designated folder? The image data will be cleared from the current database!' and has 'Yes(Y)' and 'No(N)' buttons.

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
Detailed Patient Information

Basic Information General Inspection Procedure Status Procedure Picture

GP Last Pregnancy Terminal Result of Pregnancy Breast Feeding

Menstrual Period Amount of Menses Dysmenorrhea Last Menstruation Contraception History

Past History Allergy Blood Pressure mmHg

Body Temperature Heartbeat Respiration Pulse

Vulva Vagina Cleanliness Cleanliness

Uterus Size Aden Blood Routine

Urine Pregnant Test Trichomonas Candida

Diameter of Embryonic Sac under Ultrasound Diagnosis

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Uterus Position: Cervical Depth (Pre-Pos): cm Cervical Size (Pre-Pos): cm

Suction Tube Size: Negative Pressure: Aspirates Content:

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Contraceptive Device Placement: Choose "Yes" to unlock the following tests

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