	INSTRUCTION FOR USE iVascular luminor	КОПИЯ RUSSIA/Edition 2
	PACLITAXEL ELUTING PTA BALLOON DILATATION CATHETER	DATE: July 2019
		Page 1 of 32

PEDRO-ÁNGEL CASADO MARTÍN
PUBLIC NOTARY
Registered with number 655/2019 in
my Registered Book of Certifications.

INSTRUCTION FOR USE (RUSSIA)

Paclitaxel eluting PTA balloon dilatation catheter iVascular Luminor

APPROVED

M. Eugenia Villanueva

Qualified Person

LIFE VASCULAR DEVICES BIOTECH, S.L.
C.I.F. B-65405169
Cami de Ca n'Ubach, 11
(Pol. Ind. Les Fallúies)
08620 SANT VICENÇ DELS HORTS (Barcelona)

1. DEVICE NAME

Paclitaxel eluting PTA balloon dilatation catheter iVascular Luminor.

Modifications:

- luminor 14m
- luminor 18
- luminor 35

2. MANUFACTURER

LIFE VASCULAR DEVICES BIOTECH, S.L.

Camí de Can Ubach, 11;

Pol. Ind. Les Fallulles

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Barcelona, Spain.

Tel/fax +34936724711.

E-mail: info@ivascular.es

3. AUTHORIZED REPRESENTATIVE

RAYMED TRADING GROUP PVT, LLC

141001, Moscow region, Mytishchinsky district, Mytishi, 1st Vokzalnaya street, bld. 8

Tel: (495) 956-83-16;

Fax: (499) 706-83-16

e-mail: raymedoffice@mail.ru


4. GENERAL DESCRIPTION

Brief description of the device with operating principles

The family of Paclitaxel-eluting balloon catheters is designed to dilate a stenosis or peripheral arterial obstruction with the aim of artery recanalization.

The procedure is conducted by minimally invasive techniques also known as percutaneous transluminal angioplasty (PTA), where a puncture on the external skin allows access to the blood flow and to the lesion.

The medical device is a two lumen balloon catheter (double lumen from the connector to the tip, also known as OTW), designed for percutaneous transluminal angioplasty of peripheral arteries. In the case of ***luminor 35***, the two lumen are parallel and ***luminor 14m&18*** have a coaxial design with an internal and an external lumen along all the catheter shaft.

EX0651248

I, **PEDRO-ANGEL CASADO MARTIN**, Notary of the **Catalonian Association of Notaries**, living in the city of **Barcelona (SPAIN)**,
CERTIFY:-----

That the previous signature belongs to **Mrs. M^a EUGENIA VILLANUEVA PÉREZ**, of legal age, of Spanish nationality and with DNI/NIF. (National Identity Card) number **37323062L.-----**

Mrs. M^a EUGENIA VILLANUEVA PÉREZ acts in name and on behalf of the company **"LIFE VASCULAR DEVICES BIOTECH, S.L."**, a company validly existing and incorporated under the laws of Spain; with registered address at **Sant Vicenç dels Horts, industrial state Les Fullulles, Camí de Can Ubach, number 11, C.P. 08620; with Tax Identification Number ("CIF") B-65405169; and registered in the Commercial Registry of Barcelona Volume 42171, folio 171 and sheet number B-402656.-----**

She acts as attorneys by virtue of the power of attorney granted by the notary public of Barcelona, Mr. **Pedro Angel Casado**

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Martín, on 20 of
November of the year
2013, with number of
official record 2281. --

I, the Public
Notary, attest that I
understand the English
language.-----

Likewise I, the
Notary, certify the
identity of the
signatory party and by
its virtue I state that
the previous signature
belongs to Mrs.
M^aEUGENIA VILLANUEVA
PÉREZ.-----

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
Figure 1. The figure above shows the structure of all models, the figure below shows the internal shaft of luminor 35 on the left and luminor 14m& 18 on the right

At the distal part of the catheter, before the tip, is located the balloon or inflatable segment to different pressures that dilates the artery when it is inflated by means of contrast liquid infusion inside it. The balloon is coated with the drug Paclitaxel included in a physiologically innocuous matrix. Drug dosage form is $3 \mu\text{g} / \text{mm}^2$ of balloon surface and it is intended for avoiding cell proliferation, while re-intervention index subsequently decreases.

Paclitaxel coating technique does not involve chemical adherence or a controlled release mechanism, so drug eluting from the balloon consists on a fast inflation with a subsequent drug high dosage release. Then, inflation must last between 30s and a minute for drug release onto the artery walls; however, lesion dilatation may be optimized by means of larger inflation times according to the user's discretion.

Between the luer hub and the catheter shaft, there is the strain relief tube for avoiding catheter kinking at this point.

The connector hub is Y-shaped and presents two entry ports:

- The side port, which is directly connected to the smallest lumen of the catheter shaft for **luminor 35** and to outer lumen for **luminor 14m& 18**, allows contrast liquid passage inside it to dilate the balloon.
- The straight port is intended for guidewire insertion. The guidewire is inserted into the connector entry port and it comes out of the tip of the catheter to reach the lesion by crossing the biggest lumen (**luminor 35**) or the internal lumen (**luminor 14m & 18**).

To delimitate the catheter length, radiopaque markers are placed inside of the balloon that indicate the catheter position to the user, when it is inside the patient and it is viewed on a X-ray screen.

At the distal end is the tip, which is rounded and atraumatic in shape in order to avoid damaging the arteries while it is being advanced.

To inflate the balloon, the connector inflation port must be connected to an inflation device with a manometer. The balloon is inflated to predictable diameters with the pressure shown by the manometer. The balloon is intended to reach different diameters at different pressures (according to the compliance curve shown on the primary packaging).

The balloon catheter is supplied in different diameters and lengths of balloon, which are properly labelled in the primary and secondary packaging. The useful length of the catheter depends on the model: all references of **luminor 35** are available for 80 cm and 140 cm; all references of **luminor 14m** are available for 100 or 150cm and references of **luminor 18** are available for 100cm and 140/150cm depending on the item reference. All available product references are included in the tables below.

Variants

In order to distinguish the different available configurations of the product, they are referenced.

Reference code is divided in five parts.

- The first one informs about the kind of product; in this case **BP**: Peripheral Balloon.
- The second part: **DPCxx** indicates that it is a *Coaxial Drug-eluting balloon for (Pre)dilatation*. Where xx relates to the maximum guidewire size to be used in combination (14 indicates guidewire of 0.014 inches (**luminor 14m**) 18 indicates 0.018 inches (**luminor 18**) and 35 indicates 0.035 inches (**luminor 35**)).
- The rest describe the product sizes: useful catheter length (cm), balloon diameter (hundredths of mm) and stent length (in mm, 3 digits).

All available product reference items are included in the following tables.

1. References of **luminor 14M** with catheter length 100 cm

Ø mm	LENGTH (mm)					
	40	60	80	120	150	200
1.5	BP DPC14M 100 150 040	BP DPC14M 100 150 060	BP DPC14M 100 150 080	N.A	N.A	N.A
2	BP DPC14M 100 200 040	BP DPC14M 100 200 060	BP DPC14M 100 200 080	BP DPC14M 100 200 120	BP DPC14M 100 200 150	BP DPC14M 100 200 200
2.5	BP DPC14M 100 250 040	BP DPC14M 100 250 060	BP DPC14M 100 250 080	BP DPC14M 100 250 120	BP DPC14M 100 250 150	BP DPC14M 100 250 200
3	BP DPC14M 100 300 040	BP DPC14M 100 300 060	BP DPC14M 100 300 080	BP DPC14M 100 300 120	BP DPC14M 100 300 150	BP DPC14M 100 300 200
3.5	BP DPC14M 100 350 040	BP DPC14M 100 350 060	BP DPC14M 100 350 080	BP DPC14M 100 350 120	BP DPC14M 100 350 150	BP DPC14M 100 350 200
4	BP DPC14M 100 400 040	BP DPC14M 100 400 060	BP DPC14M 100 400 080	BP DPC14M 100 400 120	N.A	N.A

N.A.= Not Available

2. References of luminor 14M with catheter length 150cm

Ø mm	LENGTH (mm)					
	40	60	80	120	150	200
1.5	BP DPC14M 150 150 040	BP DPC14M 150 150 060	BP DPC14M 150 150 080	N.A	N.A	N.A
2	BP DPC14M 150 200 040	BP DPC14M 150 200 060	BP DPC14M 150 200 080	BP DPC14M 150 200 120	BP DPC14M 150 200 150	BP DPC14M 150 200 200
2.5	BP DPC14M 150 250 040	BP DPC14M 150 250 060	BP DPC14M 150 250 080	BP DPC14M 150 250 120	BP DPC14M 150 250 150	BP DPC14M 150 250 200
3	BP DPC14M 150 300 040	BP DPC14M 150 300 060	BP DPC14M 150 300 080	BP DPC14M 150 300 120	BP DPC14M 150 300 150	BP DPC14M 150 300 200
3.5	BP DPC14M 150 350 040	BP DPC14M 150 350 060	BP DPC14M 150 350 080	BP DPC14M 150 350 120	BP DPC14M 150 350 150	BP DPC14M 150 350 200
4	BP DPC14M 150 400 040	BP DPC14M 150 400 060	BP DPC14M 150 400 080	BP DPC14M 150 400 120	N.A	N.A

N.A = Not Available

3. References of luminor 18 with catheter length 100 cm

Ø mm	LENGTH (mm)						
	20	40	60	80	120	150	200
2	BP DPC18 100 200 020	BP DPC18 100 200 040	BP DPC18 100 200 060	BP DPC18 100 200 080	BP DPC18 100 200 120	BP DPC18 100 200 150	BP DPC18 100 200 200
2.5	BP DPC18 100 250 020	BP DPC18 100 250 040	BP DPC18 100 250 060	BP DPC18 100 250 080	BP DPC18 100 250 120	BP DPC18 100 250 150	BP DPC18 100 250 200
3	BP DPC18 100 300 020	BP DPC18 100 300 040	BP DPC18 100 300 060	BP DPC18 100 300 080	BP DPC18 100 300 120	BP DPC18 100 300 150	BP DPC18 100 300 200
3.5	BP DPC18 100 350 020	BP DPC18 100 350 040	BP DPC18 100 350 060	BP DPC18 100 350 080	BP DPC18 100 350 120	BP DPC18 100 350 150	BP DPC18 100 350 200
4	BP DPC18 100 400 020	BP DPC18 100 400 040	BP DPC18 100 400 060	BP DPC18 100 400 080	BP DPC18 100 400 120	BP DPC18 100 400 150	BP DPC18 100 400 200
5	BP DPC18 100 500 020	BP DPC18 100 500 040	BP DPC18 100 500 060	BP DPC18 100 500 080	BP DPC18 100 500 120	BP DPC18 100 500 150	BP DPC18 100 500 200
6	BP DPC18 100 600 020	BP DPC18 100 600 040	BP DPC18 100 600 060	BP DPC18 100 600 080	BP DPC18 100 600 120	BP DPC18 100 600 150	BP DPC18 100 600 200
7	BP DPC18 100 700 020	BP DPC18 100 700 040	BP DPC18 100 700 060	BP DPC18 100 700 080	BP DPC18 100 700 120	N.A	N.A
8	BP DPC18 100 800 020	BP DPC18 100 800 040	BP DPC18 100 800 060	BP DPC18 100 800 080	N.A	N.A	N.A

N.A = Not Available

4. References of luminor 18 with catheter length 140 cm or 150 cm

Ø mm	LENGTH (mm)						
	20	40	60	80	120	150	200
2	BP DPC18 150 200 020	BP DPC18 150 200 040	BP DPC18 150 200 060	BP DPC18 150 200 080	BP DPC18 150 200 120	BP DPC18 150 200 150	BP DPC18 150 200 200
2.5	BP DPC18 150 250 020	BP DPC18 150 250 040	BP DPC18 150 250 060	BP DPC18 150 250 080	BP DPC18 150 250 120	BP DPC18 150 250 150	BP DPC18 150 250 200
3	BP DPC18 150 300 020	BP DPC18 150 300 040	BP DPC18 150 300 060	BP DPC18 150 300 080	BP DPC18 150 300 120	BP DPC18 150 300 150	BP DPC18 150 300 200
3.5	BP DPC18 150 350 020	BP DPC18 150 350 040	BP DPC18 150 350 060	BP DPC18 150 350 080	BP DPC18 150 350 120	BP DPC18 150 350 150	BP DPC18 150 350 200
4	BP DPC18 150 400 020	BP DPC18 150 400 040	BP DPC18 150 400 060	BP DPC18 150 400 080	BP DPC18 150 400 120	BP DPC18 150 400 150	BP DPC18 150 400 200
5	BP DPC18 140 500 020	BP DPC18 140 500 040	BP DPC18 140 500 060	BP DPC18 140 500 080	BP DPC18 140 500 120	BP DPC18 140 500 150	BP DPC18 140 500 200
6	BP DPC18 140 600 020	BP DPC18 140 600 040	BP DPC18 140 600 060	BP DPC18 140 600 080	BP DPC18 140 600 120	BP DPC18 140 600 150	BP DPC18 140 600 200
7	BP DPC18 140 700 020	BP DPC18 140 700 040	BP DPC18 140 700 060	BP DPC18 140 700 080	BP DPC18 140 700 120	N.A.	N.A.
8	BP DPC18 140 800 020	BP DPC18 140 800 040	BP DPC18 140 800 060	BP DPC18 140 800 080	N.A.	N.A.	N.A.

N.A = Not Available

5. References of luminor 35 with catheter length 80 cm

Ø mm	LENGTH (mm)					
	20	40	60	80	120	150
5	BP DPC35 080 500 020	BP DPC35 080 500 040	BP DPC35 080 500 060	BP DPC35 080 500 080	BP DPC35 080 500 120	BP DPC35 080 500 150
6	BP DPC35 080 600 020	BP DPC35 080 600 040	BP DPC35 080 600 060	BP DPC35 080 600 080	BP DPC35 080 600 120	BP DPC35 080 600 150
7	BP DPC35 080 700 020	BP DPC35 080 700 040	BP DPC35 080 700 060	BP DPC35 080 700 080	BP DPC35 080 700 120	N.A

N.A = Not Available

6. *References of luminor 35 with catheter length 140cm*

Ø mm	LENGTH (mm)					
	20	40	60	80	120	150
5	BP DPC35 140 500 020	BP DPC35 140 500 040	BP DPC35 140 500 060	BP DPC35 140 500 080	BP DPC35 140 500 120	BP DPC35 140 500 150
6	BP DPC35 140 600 020	BP DPC35 140 600 040	BP DPC35 140 600 060	BP DPC35 140 600 080	BP DPC35 140 600 120	BP DPC35 140 600 150
7	BP DPC35 140 700 020	BP DPC35 140 700 040	BP DPC35 140 700 060	BP DPC35 140 700 080	BP DPC35 140 700 120	N.A

N.A.= Not Available

5. **INTENDED USE**

Paclitaxel eluting PTA balloon dilatation catheter iVascular Luminor is indicated for dilatation or opening of stenosis located in the peripheral arteries (the iliac, femoral, iliofemoral, popliteal, infrapopliteal and renal arteries), as well as for the treatment of obstructive lesions of arteriovenous fistulas, whether original or artificial, with a reference diameter balloon between 1.5 and 4.0 mm and lengths from 20 to 200 mm (*luminor 14m*), from 2 to 8 mm and lengths from 20 to 200 mm (*luminor 18*) and between 5 and 7 mm and lengths from 20 to 150 mm (*luminor 35*). It is also indicated for stent postdilation in the peripheral vascular system.


Figure 2.

Prevalent population: patients with peripheral arterial disease (PAD).

Peripheral arterial disease, also known as peripheral vascular disease (PVD), is induced by artery atherosclerosis (hardening). Atherosclerosis provokes peripheral artery narrowing (stenosis) and in Peripheral Artery Disease (PAD) located in lower limbs, the legs blood flow is constraint. The main interesting anatomical parts in lower extremity PAD include below the knee (BTK) arteries the superficial femoral artery (SFA) and the femoropopliteal segment.

The most advanced form of Peripheral Artery Disease (PAD) is Critical Limb Ischemia (CLI) which manifests pain while resting, gangrene and ulceration. A viability reduction in the limb can be experimented with Critical Limb Ischemia (CLI), with potential risk of progression to amputation, particularly in diabetic patients. By the way, Critical Limb Ischemia (CLI) treatment is required in order to prevent amputation. These treatments involve the closed artery dilatation

6. EXPECTED USES, USERS, PATIENTS, SORROUNDINGS

The product is designed for its use in Hemodynamic/Radiology Rooms by Radiologists or vascular surgeons. The procedure is carried out with local anaesthesia at the place of the insertion. The advancement of the catheter through the arteries is observed by Fluoroscopic X-ray on a screen. The user will see on the screen the catheter radiopaque markers and he will monitor its advance.


Figure 3. Hemodynamic labs


Figure 4. *X-ray image*

The technique for insertion consists of:


The access can be from the femoral (leg) or the radial (wrist) artery.


Figure 5. *Insertion technique (access)*


A puncture needle is inserted into the patient's groin (accessing to the superficial femoral artery) or wrist (accessing to the radial artery).

A guidewire is placed through this needle.
The needle is removed while the guidewire remains inside.

An introducer with dilator is introduced through the guidewire.

The dilator and the guidewire are removed and the introducer remains ready for the introduction of the wished number of catheters.

The guidewire is then advanced until crossing the lesion to dilate.

The balloon catheter, which is folded by the moment, is introduced on the guidewire by the lumen guide.

In case of complex or occlusive lesions, preliminary pre-dilatation with a conventional balloon without drug is recommended. Once the lesion is prepared (if it is required), drug eluting balloon catheter is placed at the injured area. The catheter preparation before its insertion involves attaching a three-way stopcock to the coaxial luer connector of the inflation port (Optionally with an extension between both ports, to improve the device flexibility). The three-way stopcock remains connected then: on one hand, to the side luer hub of the balloon coaxial luer hub and, on the other hand, to a luer hub connected to an inflation device, and the middle connection allows the purge of the balloon with a syringe.

The balloon catheter is advanced to the lesion. The pressure inside the balloon is raised with the inflation device until it gets inflated and it dilates the artery. Once the wished pressure is reached, it is maintained during 30-60s for drug transfer to the artery wall.

Later, the balloon is deflated and removed through the introducer.

The elements that contact the balloon catheter during the implantation are:

- The **arterial introducer**, which, as mentioned before, remains along the intervention allowing access to the artery of the different devices. The introducer has a backstop valve on his proximal end which allows the passage of the devices and avoids the loss of the patient's blood. Introducers are manufactured in different French gauges to allow the passage of different diameter devices. The French gauge is determined by the internal diameter of the introducer, as the different devices go inside, and the equivalence in mm is obtained by the formula:

$$\text{mm} = F/3$$

The used introducers for *luminor* are those of 4, 5 and 6 French:

French	4	5	6
Internal Ø of the introducer (mm)	1.33	1.67	2.00

Likewise, the introducer has a colour code in order to identify its French size:

Red: 4F, Grey: 5F, Green: 6F, Orange: 7F, Blue: 8F.

- **0.014, 0.018 or 0.035 inches guidewire**, in order to cross the balloon catheter's guide lumen.

A guidewire is a stainless steel wire thread with long length (at least 180 cm), whose distal part is very flexible to avoid damage to the arteries. The guidewire allows the advancement of

	INSTRUCTION FOR USE Vascular luminor	RUSSIA/ Edition 2
	PACLITAXEL ELUTING PTA BALLOON DILATATION CATHETERS	DATE: July 2019
		Page 12 of 32

treatment to catheters, balloon catheters, stents, etc. Its distal extreme flexibility is due to the following reasons:

- It must cross lesions that sometimes are very tortuous.
- It is the first device to go out by the introducer and to go into by itself inside the channel to treat.

The guidewire function is to cross lesions and remain in a place as far as possible to allow the canalization of the different devices used in the treatment.

The guidewire size is expressed as its outer diameter in inches, so, a 0.035" guidewire has an outer diameter of 0.89 mm.

- **three-way stopcocks / extensions** connecting different accessories, and
- **Inflation device** to dilate the balloon.
- Several standard **10-20cc** syringes with saline solution for washing the system.
- **Heparinized normal saline solution.**
- **Contrast medium** (use contrast media suitable for intravascular use).
- **Appropriate medication to the procedure**

7. CONTRAINDICATIONS

- Inability to cross the lesion with a guidewire.
- It must not be used in pregnant or lactating women or patients with diagnosed Paclitaxel hypersensitivity.

8. WARNINGS

- The device should only be used by experienced doctors familiar in the percutaneous transluminal angioplasty (PTA) technique.
- The percutaneous procedure must be carried out in hospitals equipped with emergency facilities to conduct open heart surgery or hospitals with easy access to other hospital where this type of operation can be performed.
- The product is for single use. It must not be re-sterilised or re-used once the procedure has been completed. Re-using the product in another patient may lead to cross contamination, infection or transmission of infectious diseases from one patient to another. Re-use of the product may cause alterations to it and affect its effectiveness.
- The product is supplied sterile. Check the expiry date and do not use products which are past this date.
- Inspect the packaging before opening it. If there are any defects or the packaging is damaged, do not use the product.
- Use aseptic techniques when the product is taken out of the primary packaging.

- Select the balloon diameter according to the diameter of the artery to be treated. Do not insert a balloon with a diameter larger than the artery to be treated.
- Do not touch the balloon directly or dry it with gauze, since doing so may cause the drug coating to come off.
- Do not expose the product to organic solvents such as alcohol, or to contrast media not indicated for intravascular use.
- Do not use air or gaseous media to inflate the balloon.
- Advance the product on the guide wire using fluoroscopy. Do not allow the product to advance without the guide wire inside it.
- Do not manipulate, advance or retract either the catheter or the wire guide while the balloon is inflated.
- In case a post-dilatation is necessary after use of product, it must be made with a conventional dilatation balloon, one without the drug, since the same vessel segment must not be dilated with more than one balloon with drug.
- Drug eluting stents must not be implanted in the same segment that has been previously treated with the balloon with drug, since it is possible for overdose or interaction among the active principles to occur.
- Published in vitro and in vivo studies prove the mutagenic effect of Paclitaxel. The carcinogenesis risk potential has not been studied. The published in vivo studies prove the teratogenic effect, therefore it is not recommended for pregnant women to use this product.

9. PRECAUTIONS

- Administer suitable medical therapy to the patient: anticoagulants, vasodilators, etc., in accordance with the procedure for inserting intravascular catheters.
- Use with standard 0.014, 0.018 or 0.035 inches guidewires and standard introducers. If the introducer length goes above 25cm or the braided type is used, it may be necessary to increase its size to avoid catheter friction.
- Before use, check the size and suitability of the catheter for its intended use.
- Use the product with precautions in complex lesions, since the abrasive nature of said lesions may alter the drug coating. In these cases it is recommended to predilate the lesion previously with a conventional, drugless balloon.
- Extreme care must be exercised when the device is taken off from the pack and passed through the hemostatic valve to make sure that the drug coating is not damaged.
- Proceed with extreme care so as not to damage the catheter while advancing it.
- The balloon must be inflated with a mixture of saline solution and contrast medium (preferably 50/50, though the proportion of contrast medium can be reduced in large balloons, like 4mm diameter or longer ones).
- If any resistance is noticed during advance, stop and determine its cause before

continuing.

- Do not exceed the rated burst pressure (RBP) as the balloon could burst.
- If any resistance to withdrawing the catheter is noted, it is recommended that the balloon catheter, the guide wire and introducer be removed as a single unit.
- Store in a cool, dry place away from direct sunlight.
- After use, this product may represent a biohazard. Handle and dispose it in accordance with accepted medical practices and pertinent local, state or federal laws and regulation

10. POSSIBLE ADVERSE EFFECTS / COMPLICATIONS

The possible adverse effects and/or complications which might arise before, during or after the procedure include the following:

- Death
- Ictus / embolism / thrombosis
- Haemodynamic deterioration
- Acute elastic recoil
- Arrhythmia
- Aneurysm or pseudoaneurysm
- Infections
- Total occlusion of the artery
- Reocclusion of the area treated: restenosis
- Spasm
- Perforation or dissection of the area treated
- Arteriovenous fistula
- Local bleeding with haematoma at the access site
- Hypo/hypertension
- Allergic or immune reaction to the drug.
- Histological **changes in the arterial wall.**

11. DESCRIPTION OF THE DRUG INCORPORATED ON THE DEVICE

11.1. Additional incorporated active substance

The action ancillary medicinal substance incorporated on the paclitaxel-eluting peripheral balloon dilatation catheters *iVascular Luminor* is Paclitaxel.


Paclitaxel is an organic compound and the primary active substance of Taxane substance category, which have shown antineoplastic properties. Paclitaxel as active substance was isolated from the bark of the "Pacific yew" tree, *Taxus brevifolia*, in 1967. WANI and compa

successfully identified its structural formula and demonstrated its cytotoxic effects in different murine tumour.

Paclitaxel is chemically known as 5 β , 20-Epoxy-1,7-dihydroxy-9-oxotax-11-ene-2 α ,4,10 β ,13 α -tetrayl 4,10-diacetate 2-benzoate 13-[(2R,3S)-3-(benzoylamino)-2-hydroxy-3-phenylpropanoate].

Its molecular formula is C₄₇H₅₁NO₁₄, whose molecular weight is 853.9 g/mol.

Its structure is:

Paclitaxel is a crystalline powder whose colour is between white and almost white, which is partially soluble in water, soluble in methanol and ethanol, and very soluble in methylene chloride.

11.2. Justification of the medicinal substance inclusion

The aim of the treatment of peripheral artery disease is to decrease cardiovascular death and risk of lower limb amputation in patients with critical limb ischemia (CLI), as well as it improves its patient's quality life¹.

From the origin, chosen method for treatment of injured femoral arteries by revascularization techniques has been balloon angioplasty². Afterwards, stents have been used as an alternative to decrease induced restenosis by balloon angioplasty. Theory on that pharmacological treatment of neointimal hyperplasia decreases restenosis has led to drug eluting stents (DES) development. However, DES may not be suitable for use in complex lesions (tortuous or small vessels, or large or diffuse calcified lesions)³. Therefore, their use in PAD, which is usually located in femoropopliteal and below the knee (BTK) arteries, may not be suitable⁴.

Obtained restenosis inhibition results with a Paclitaxel single-dosage^{4,5} have brought to the incorporation of drug eluting technologies alternatively to the stent. Within these technologies are situated drug eluting balloons (DEB), which may be used in femoropopliteal and below the knee (BTK) arteries, where stenting may be an unsuitable option⁴.

Drug eluting from the balloon consists on a fast inflation so that subsequent high dosage of drug is released within a very short period of time without needing a permanent implantation⁶ and thus, avoiding neointimal hyperplasia.

In some drug eluting stents, late thrombosis has been associated to the drug eluting polymer use, which results in a need of prolonged dual antiplatelet therapy for its avoidance^{7,8}. Drug eluting balloons do not require prolonged antiplatelet therapy, which is an additional advantage, due to it does not need any polymer for drug release.

Therefore, Paclitaxel coated balloons use in the treatment of peripheral artery disease may meet drug eluting stents' objectives, where stenting is an unsuitable treatment option.

	INSTRUCTION FOR USE Vascular luminor	RUSSIA/ Edition 2
	PACLITAXEL ELUTING PTA BALLOON DILATATION CATHETERS	DATE: July 2019
		Page 16 of 32

Bibliographic references

- ¹ Shammas NW. Epidemiology, classification, and modifiable risk factors of peripheral arterial disease. *Vasc Health Risk Manag.* 2007;3(2):229-34
- ² Gibbs JM, Pena CS, Benenati JF. Treating the diseased superficial femoral artery. *Tech Vasc Interv Radiol.* 2010 Mar;13(1):37-42.
- ³ Waksman R, Pakala R. Drug-eluting balloon: the comeback kid? *Circ Cardiovasc Interv.* 2009 Aug;2(4):352-8.
- ⁴ Tepe G, Schmitmeier S, Speck U, Schnorr B, Kelsch B, Scheller B. Advances on drug-coated balloons. *J Cardiovasc Surg (Torino).* 2010 Feb;51(1):125-43.
- ⁵ Tepe G, Zeller T, Albrecht T, Heller S, Schwarzwald U, Beregi JP, et al. Local delivery of paclitaxel to inhibit restenosis during angioplasty of the leg. *N Engl J Med.* 2008 Feb 14;358(7):689-99.
- ⁶ Banning AP, Lim CC. Drug-eluting balloons: what is their place on the interventionalist's shelf? *Heart* 2010. Aug;96(16):1257-8.
- ⁷ Bode C, Zehender M. The use of antiplatelet agents following percutaneous coronary intervention: focus on late stent thrombosis. *European Heart Journal Supplements* 2007; (Supplement D), D10-D19: D10-D9.
- ⁸ National Institute for Health and Clinical Excellence. TA152 - Drug-eluting stents for the treatment of coronary artery disease. 2008.

11.3. Medicinal substance quantity incorporated to the device

iVascular luminor peripheral balloon catheter involves a Paclitaxel dosage of 3 µg/mm² balloon surface.

11.4. Type of active substance. Quality, safety and utility

☒ Active substance known for the established indications

11.4.1. Data on quality of the used active substance

Its empirical formula is C₄₇H₅₁NO₁₄ and its molecular weight is 853,9 g/mol. It consists on a white or almost white crystalline powder, soluble in methanol, ethanol, methylene chloride, tetrahydrofuran and insoluble in water, that is obtained by a semi-synthetic process.

Paclitaxel control has been reported in EP and USP Pharmacopoeias. Paclitaxel incorporated in Luminor meets the Pharmacopoeia:

- The product is characterized by IR, meeting the standard material.
- The product is characterized by HPLC, retention time of the main peak agrees with the standard. Product purity, according to the test, is between 98,0 and 102,0% (recalculated on anhydrous substance).
- Most common impurities are identified: 0.2% is not exceeded in any case, and when summing them up 1.0% is neither exceeded.

11.4.2. Data on safety of the used active substance

Product safety has been demonstrated on animal preclinical trials, where this product has been compared with the standard balloon (without drug) and a marketed equivalent balloon with Paclitaxel.

On the other hand, safety of balloon with Paclitaxel, based on its drug content, will be equivalent to marketed balloons with Paclitaxel, due to its equivalence and thus, its equivalent drug content. Its safety is demonstrated in the Relevant Scientific Literature Report.

11.4.3. Data on utility of the used active substance

Paclitaxel is an anti-microtubule agent that induces microtubule formation from tubulin dimers and stabilizes them avoiding depolymerisation. This stability provokes the inhibition of normal dynamic reorganization of the microtubule network, which is essential for cell vital functions during interphase and mitosis. Therefore, it inhibits cell proliferation by a sustained mitosis block during metaphase-anaphase stage⁹.

It has antitumor activity and it is indicated for first and second line chemotherapy, as well as the treatment of some kinds of cancer (advanced ovary carcinoma, lung cancer, throat cancer) up to a 175 mg/m² body surface maximum dosage that are administrated by intravenous infusion during three hours every three weeks.

In vitro studies with vascular straight muscular cells have demonstrated inhibition of proliferation during at least 12 days as a consequence of 3 min incubation with Paclitaxel¹⁰. These features have been applied to Paclitaxel use in retarded release coating formulations on stents to reduce neointimal proliferation. Different animal studies have confirmed neointimal hyperplasia inhibition after over-dilatation and Paclitaxel eluting stenting in coronary arteries^{11,12,13}.

In turn, they have also been applied to drug eluting balloons. Studies where Paclitaxel was added to the contrast medium during angiographic control, or as a part of balloon coating, demonstrated its efficacy. In case of balloon, there was proven that with a very short contact time (60s) between drug eluting balloon and the artery walls, an equivalent or higher efficacy than with drug eluting stent use was achieved¹⁴.

Regarding to *iVascular luminor* Paclitaxel eluting balloon, it has been demonstrated its drug use and efficacy by means of preclinical trials in pigs and the Relevant Scientific Literature Report.

Bibliographic references

⁹ Gotlieb AI, May LM, Subrahmanyam L, Kalnins VI. Distribution of microtubule organizing centers in migrating sheets of endothelial cells. *J. Cell Biol.* 1981; 91:589-594

¹⁰ Axel DI, Kunert W, Goggelmann C, Oberhoff M, Herdeg C, Kuttner A et al. Paclitaxel inhibits arterial smooth muscle cell proliferation and migration in vitro and in vivo using local drug delivery. *Circulation* 1997; 96:636-645

¹¹ Heldman AW, Cheng L, Jenkins GM, Heller PF, Kim DW, Ware M, Jr. et al. Paclitaxel stent coating inhibits neointimal hyperplasia at 4 weeks in a porcine model of coronary restenosis. *Circulation* 2001; 103:2289-95

Drachman DE, Edelman ER, Seifert P, Groothuis AR, Bornstein DA, Kamath KR et al. Neointimal thickening after stent delivery of paclitaxel: change in composition and arrest of growth over six months. *J. Am. Coll. Cardiol.* 2000; 36:2325-32

	INSTRUCTION FOR USE Vascular luminor	RUSSIA/ Edition 2
	PACLITAXEL ELUTING PTA BALLOON DILATATION CATHETERS	DATE: July 2019
		Page 18 of 32

¹³ Farb A, Heller PF, Shroff S, Cheng L, Lolodgie FD, Carter AJ et al. Pathological analysis of local delivery of paclitaxel via a polymer-coated stent. *Circulation* 2001; 104:473-479

¹⁴ Scheller B, Speck U, Abramjuk C, Bernhardt U, Bohm M, Nickenig G. Paclitaxel balloon coating: a novel method for prevention and therapy of restenosis. *Circulation* 2004; 110:810-814

11.5. Monograph contained at the Pharmacopoeia of the active substance

Has the additional active substance got a monograph at the European Pharmacopoeia or, in its absence, at any National Pharmacopoeia of an EU member country or, in its absence, of national monograph of a third country?

☐ NO

☒ YES

✓ EUROPEAN PHARMACOPOEIA 6.3

- 01/2017:1794
- PACLITAXEL
- Paclitaxelum

11.6. Compatibility data with other medicinal substances

Possible Paclitaxel interactions with concurrent administered drugs have not been studied. Interactions with used level in antineoplastic treatment are detailed in Paclitaxel technical sheet; but, taking into account that balloon surface Paclitaxel quantity corresponds approximately to some hundredths of the usually quantity for antineoplastic treatment, interactions with other drugs are considered to be unlikely.

However, precaution must be taken when concurrent CYP3A4 or CYP2C8 substrate (terfenadine, ciclosporin, lovastatin, midazolam and ondansetron) or drugs with high PPI (Plasma protein binding) are administered.

Regarding to product compatibility with usually administered drugs in percutaneous transluminal angioplasty procedure, it has been evaluated and demonstrated in conducted animal (pig) preclinical trials, given that the followed protocol establishes habitual use medication in human practice (heparin, contrast liquid, antiplatelet therapy, etc).

12. BIOLOGICAL SAFETY

Not applicable. The paclitaxel-eluting peripheral balloon dilatation catheters *iVascular luminor* do not contain any kind of biological materials or related: no animal or human tissues, neither recombinant or fermentation products.

13. PACKAGING

The product is single supplied. Each primary packaging contains a Paclitaxel eluting balloon catheter for peripheral dilatation.

The catheter is presented with a balloon protective sheath which helps keeping it fold and a stylet protecting the guidewire lumen. The entire set is inserted into a dispenser to avoid damages on the catheter.

Figure 6. Final protection stylet (Stainless steel) and sheath (HDPE)

The balloon catheter is inserted into a dispenser in order to avoid damages on it while transport and storage


Figure 7. Balloon catheter inserted in a dispenser (HDPE: High density polyethylene)

The primary packaging comprises a bag containing the device to preserve product sterility and the compliance curve that shows the achieved balloon diameter in each increase of pressure.


Figure 8. Sterile bag (LLDPE: Linear low-density polyethylene) and compliance curve (glow cardboard)

The label for product identification is placed on the sterile bag. The bag contains the sterility indicators.


Figure 9. Label (semigloss couche paper) placed on the sterile bag (LLDPE)

This bag is inserted into a box with the instructions for use, to avoid damages in transport and storage. A label with all the data of the product and the manufacturer (the same label that is on the bag) is stacked on one side of the box. On the other one, a round, transparent, inviolable sticker is placed to indicate that the packaging has not been opened before its use.


Figure 10. Frontal and posterior view of the carton box with informative label provided by the manufacturer


Figure 11. Tamper-proof label

Finally, the individual packaging are placed in a shipping box manufactured from white, high resistance corrugated fireboard with capacity up to 25 boxes. This shipping box avoids damage of packaging and products during transport and storage.

Contents

- A paclitaxel-releasing dilatation balloon catheter OTW, covered by a protector over the balloon and a guide lumen protection stylet. All the set is introduced in a dispenser to avoid damaging the catheter, and it is packed into a sterile bag.
- One card with the compliance curve showing the nominal inflation pressure and the recommended maximum pressure.
- One leaflet with instructions for use.

14. HOW TO USE

14.1. *iVascular luminor 14m*

Equipment Required

- Heparinized normal saline solution
- Contrast medium (use contrast media suitable for intravascular use)
- 0.014" guide wire; do not use any other size of guide wire
- Introducer with haemostatic valve of the size specified on the label. Do not use smaller size as these might damage the catheter.
- Three-way stopcock
- Line extensions
- Inflation device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system

Preparing the Catheter

- In case of occlusion or very calcified lesions, it is recommended to predilate the lesion with a conventional balloon, without drug.
- Remove the catheter from the protective dispenser. Check that it is the right size. Gently remove the protective pod and stylet from the balloon.
- Connect a 10ml syringe containing sterile saline solution to the port for the guide wire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the system of air:
 - Attach a three-way stopcock to the contrast liquid port (the side part of the connector). An extension line could be optionally placed between the balloon connector and the three-way stopcock.
 - Close the air passage through the balloon.
 - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
 - Open the passage in the three-way stopcock between the syringe and the catheter.
 - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
 - When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

DANGER: If bubbles do not stop coming from the catheter into the syringe after 10 minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

Technique for Insertion / Treatment

- Position the 0.014" guide wire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.
- Insert the proximal end of the guide wire into the distal tip of the balloon catheter.
- Advance carefully the balloon catheter along the guide wire, through the guide catheter and with the help of fluoroscopy, until the area to be treated is reached. Extreme care must be exercised when the balloon catheter is passed through the hemostatic valve to make sure that the drug coating is not damaged. In case the balloon catheter chosen cannot be placed inside the lesion, choose a smaller size of a drugless conventional balloon in order to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Keep the pressure for 30-60 seconds for an optimal release of the drug. Most of the drug is released in the first 30 seconds of inflation but, in order to optimize the dilation of the lesion, it is possible to use longer inflation times greater than 1 minute at the operator's discretion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, carefully retract the catheter until it can be removed through the introducer.
- Inject contrast liquid through the introducer to check the treated lesion before removing the guide wire.
- Double antiplatelet therapy must be administered for a minimum of four weeks after the intervention, and the treatment may be prolonged at the operator's discretion.

14.2. iVascular luminor 18

Equipment Required

- Heparinized normal saline solution
- Contrast medium (use contrast media suitable for intravascular use)
- 0.018" guidewire (secondarily 0.014")
- Introducer with haemostatic valve of the size specified on the label. Do not use smaller sizes as these might damage the catheter.
- Three-way stopcock
- Line extensions
- Inflation device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system

Preparing the Catheter

	INSTRUCTION FOR USE Vascular luminor	RUSSIA/ Edition 2
	PACLITAXEL ELUTING PTA BALLOON DILATATION CATHETERS	DATE: July 2019
		Page 24 of 32

- In case of occlusion or very calcified lesions, it is recommended to predilate the lesion with a conventional balloon, without drug.
- Remove the catheter from the protective dispenser. Check that it is the right size. Gently remove the protective pod and stylet from the balloon.
- Connect a 10ml syringe containing sterile saline solution to the port for the guidewire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the system of air:
 - Attach a three-way stopcock to the contrast liquid port (the side part of the connector). An extension line could be optionally placed between the balloon connector and the three-way stopcock.
 - Close the air passage through the balloon.
 - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
 - Open the passage in the three-way stopcock between the syringe and the catheter.
 - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.
 - When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

DANGER: If bubbles do not stop coming from the catheter into the syringe after minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

Technique for Insertion / Treatment

- Position the 0.018" guidewire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times. If a guidewire measuring 0.014" was placed previously it is not necessary to remove it since the catheter is also compatible with a guidewire of this size.
- Insert the proximal end of the guidewire into the distal tip of the balloon catheter.
- Advance carefully the balloon catheter along the guidewire, through the guide catheter and with the help of fluoroscopy, until the area to be treated is reached. Extreme care must be exercised when the balloon catheter is passed through the hemostatic valve to make sure that the drug coating is not damaged. In case the balloon catheter chosen cannot be placed inside the lesion, choose a smaller size of a drugless conventional balloon in order to predilate the lesion.
- Connect the inflation device (1/2-1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.

- Keep the pressure for 30-60 seconds for an optimal release of the drug. Most of the drug is released in the first 30 seconds of inflation but, in order to optimize the dilation of the lesion, it is possible to use longer inflation times greater than 1 minute at the operator's discretion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guidewire in position, carefully retract the catheter until it can be removed through the introducer.
- Inject contrast liquid through the introducer to check the treated lesion before removing the guidewire.
- Double antiplatelet therapy must be administered for a minimum of four weeks after the intervention, and the treatment may be prolonged at the operator's discretion.

14.3. *iVascular luminor 35*

Equipment Required

- Heparinized normal saline solution
- Contrast medium (use contrast media suitable for intravascular use)
- 0.035" guide wire; do not use any other size of guide wire
- Introducer with a haemostatic valve of the size specified on the label. Do not use smaller sizes as these might damage the catheter. If an introducer more than 25cm long or of the braided type is used, it may be necessary to increase its size to avoid rubbing on the catheter.
- Three-way stopcock
- Inflation device (manual pump with built-in manometer)
- Several standard 10-20cc syringes with saline solution for washing the system

Preparing the Catheter

- In case of occlusion or very calcified lesions, it is recommended to predilate the lesion with a conventional balloon, without drug.
- Remove the catheter from the protective dispenser. Check that it is the right size. Gently remove the protective pod and stylet from the balloon.
- Connect a 10ml syringe containing sterile saline solution to the port for the guide wire (straight part of the connector) and irrigate the lumen until the liquid comes out of the tip.
- Purge the system of air:
 - Attach a three-way stopcock to the contrast liquid port (the side part of the connector).
 - Close the air passage through the balloon.
 - Attach a 10-20ml syringe with a third part of saline solution to the three-way stopcock.
 - Open the passage in the three-way stopcock between the syringe and the catheter.
 - With the syringe in a vertical position, withdraw the plunger upwards allowing the air bubbles to exit to the liquid.

	INSTRUCTION FOR USE Vascular luminor	RUSSIA/ Edition 2
	PACLITAXEL ELUTING PTA BALLOON DILATATION CATHETERS	DATE: July 2019
		Page 26 of 32

- When bubbles stop coming in, close the three-way stopcock at the catheter end and remove the syringe.

DANGER: If bubbles do not stop coming from the catheter into the syringe after minutes of negative pressure, this may be a clear indication that the balloon catheter has leaks, is broken or the connections between the syringe and the three-way stopcock are not properly closed. If bubbles continue to be seen after checking the connections, do not use the device. Return it to the manufacturer or distributor for inspection.

Technique for Insertion / Treatment

- Position the 0.035" guide wire across the lesion using PTA techniques, using fluoroscopy to determine its position at all times.
- Insert the proximal end of the guide wire into the distal tip of the balloon catheter.
- Advance carefully the balloon catheter along the guide wire, through the guide catheter and with the help of fluoroscopy, until the area to be treated is reached. Extreme care must be exercised when the balloon catheter is passed through the hemostatic valve to make sure that the drug coating is not damaged. In case the balloon catheter chosen cannot be placed inside the lesion, choose a smaller size of a drugless conventional balloon in order to predilate the lesion.
- Connect the inflation device (1/3 of its capacity of contrast liquid mixture: saline solution) to the three-way stopcock and open the passage between the balloon and the inflation device.
- Apply pressure to the inflation device so that the balloon inflates. Do not exceed the rated burst pressure (RBP) on the label and the compliance curve.
- Keep the pressure for 30-60 seconds for an optimal release of the drug. Most of the drug is released in the first 30 seconds of inflation but, in order to optimize the dilation of the lesion, it is possible to use longer inflation times greater than 1 minute at the operator's discretion.
- Draw the plunger on the inflation device back to deflate the balloon. Keep up the negative pressure for between 20 and 60 seconds depending on the size of the balloon. Make sure the balloon is fully deflated (by means of fluoroscopy) before moving the catheter.
- With negative pressure in the inflation device and with the guide wire in position, carefully retract the catheter until it can be removed through the introducer.
- Remove the guide wire.
- Double antiplatelet therapy must be administered for a minimum of four weeks after the intervention, and the treatment may be prolonged at the operator's discretion.

15. TECHNICAL SPECIFICATIONS

<u>General features</u>	<i>luminor 14M</i>	<i>luminor 18</i>	<i>luminor 35</i>
Balloon	Semi-compliant (10-15% compliance from nominal pressure to RBP)		
Balloon coating	Paclitaxel (3 µg/mm ²)		
Outer diameter (±0.05 mm)	3.9F/3.8F/2.9F (1.31 mm/1.28 mm/0.98 mm)	4.0F/3.9F/3.2F (1.33 mm/1.31 mm/1.07 mm)	5F (1.64 mm)
Nominal pressure	7 atm	7 atm	6-7 atm
Rated Burst Pressure (RBP)	16 atm (14 atm for diameter 8mm)		
Radiopaque markers	2 radiopaque markers of Pt / Ir, one at each end of the balloon, width 0.9±0.1 mm		
Introducer compatibility	4F (1.33 mm) - all diameters	4F (1.33 mm) - Ø2-5mm, 5F (1.64 mm) - Ø 6-8mm	5F (1.64 mm) - Ø 5 mm 6F (2 mm) - Ø 5-7 mm
Recommended guidewire	0.014" (0.36 mm)	0.018" (0.46 mm) 0.014" (0.36 mm) (optional)	0.035" (0.89 mm)
Deflation time	Maximum 10s	Maximum 8s	Maximum 8s
Available diameters (±0.5 mm)	1.5 to 4 mm	2 to 8 mm	5 to 7 mm
Available lengths (± 1.5 mm)	40 to 200 mm	20 to 200 mm	20 to 150 mm

16. STERILIZATION CONDITIONS

The devices are provided sterile by ethylene oxide (EO) terminal sterilization. The devices are non-pyrogenic. Note the product "Use By" date specified on the package. DO NOT use the device if the "Use by" date has passed.

The devices are intended for single-use only; do not reuse. Do not resterilize, as this can compromise device performance and increase the risk of cross contamination due to inappropriate reprocessing. Inspect all product prior to use. Do not use if the package is open or damaged, or if product is damaged.

	INSTRUCTION FOR USE Vascular luminor	RUSSIA/ Edition 2
	PACLITAXEL ELUTING PTA BALLOON DILATATION CATHETERS	DATE: July 2019
		Page 28 of 32

17. TRANSPORTATION, STORAGE

Handle with care. Store in original packaging. Keep dry. Keep away from sunlight.

Devices packaged by manufacturer may be transported by all modes of transport in covered vehicle in accordance with the rules of transportation of goods for each type of transport.

Products should be stored at +15°C to +30°C in an enclosed space in the original manufacturer's packaging for the labeled shelf life.

18. SHelf LIFE

The shelf life of Paclitaxel eluting PTA balloon dilatation catheter iVascular luminor is 36 months.

19. ENVIRONMENT PROTECTION REQUIREMENTS

The devices do not adversely impact the environment.

20. DISPOSAL

Products, which had contact with blood and/or other biological liquids of patients should be disposed and eliminated in accordance with the requirements for handling medical waste, established by SanPiN 2.1.7.2790-2010 for class B and other normative legal acts of the Russian Federation.

Unused products, including the ones which expired (which had no contact with the biological liquids of patients) should be disposed in accordance with the procedure of disposal of solid domestic waste established by the legislation of the Russian Federation for class A by SanPiN 2.1.7.2790.

21. DESIGN VERIFICATION TESTS

21.1. Methodology

Design verification protocol has been carried, providing specific reports for each standard. These tests, specified in the corresponding report, have been conducted according to the methodology established in the current standards. The obtained results are detailed, confirming the achievement of acceptance criteria, set for each test in accordance with the current standards.

21.2. Design verification tests

The following table includes the different tests performed in order to comply with the standard classified depending on the section of the standard that applies. Tested references and batch numbers, acceptance criteria, results and indications about whether they comply with the requirement or not are included as well.

Vascular luminor
Physical and mechanical bench testing according to current international standards

TEST	STANDARD	ACCEPTANCE CRITERIA	RESULTS	Ok/no OK
Biocompatibility	ISO 10555-1	The balloon catheter should be biocompatible	LUMINOR are biocompatible.	OK
Outer surface		Absence of strange elements under 2.5x Catheter surface free of drug	No observations of foreign material or drug moieties in the entire catheter surface.	OK
Resistance to corrosion		Not applicable. The device does not have metallic elements that may be susceptible of corrosion		
Tensile strength		All the components and joints must exceed values of tensile strength established by UNE EN ISO 10555-1 standards	All tubes and joints exceed the minimum established strength by UNE EN ISO 10555-1 standards	OK
Leaks absence		Liquid leak absence under pressure. Airways inwards the cone absence by aspiration	Samples have been kept more than 30s under 3000kPa and is noted: – Absence of leaks or drops – No bubbles nor bubble generation is observed	OK
Nominal size		Outer diameter of balloon in mm Effective length of catheter in cm Effective balloon length, in mm	The results are according to labelling of references tested	OK

<i>iVascular luminor</i> Physical and mechanical bench testing according to current international standards				
TEST	STANDARD	ACCEPTANCE CRITERIA	RESULTS	Ok/no OK
Radiodetectability	ISO 10555-4	The balloon is visible under fluoroscopy	Pt/ 10% Ir radiopaque markers, are fluoroscopically viewable	OK
Distal tip configuration		The tip must be soft, round and conic-shaped	The tip is made of a soft materials and is conic and atraumatic in shape end	OK
Leak/damage absence after fatigue		Nor leak nor damage after fatigue test	Once 10 cycles at 16atm are conducted, nor damages nor leaks are observed	OK
Taper 6%	20594-1	Both coaxial luers cones must be 6% tapered	Peripheral balloon coaxial luer cones meets 6% taper	OK
Control drug loading	Adaptation ISO 25539-2	Load within the expected range ($3,0 \pm 0,5 \mu\text{g}/\text{mm}^2$)	Loads are within the range established	OK
Drug stability		Percentage of impurities within the limits of acceptance of the pharmacopeia	The percentage of impurities is within the limits of acceptance of the pharmacopeia	OK
Mass loss	ISO 10993-13	Mass loss should be less than 2% at each established period	The mass loss is less than 2%	OK
Bioburden	ISO 11737-1	<50 cfu/g	<9 cfu/g	OK
Sterility	ISO 11737-2	The sample must be sterile	Sterile	OK
Ethylene chlorohydrin	ISO 10993-7	<9 mg/sample	<2 mg	OK
Ethylene oxide residuals	ISO 10993-7	<4 mg/sample	<0.1 mg	OK

Vascular luminor

Physical and mechanical bench testing according to current international standards

TEST	STANDARD	ACCEPTANCE CRITERIA	RESULTS	Ok/no OK
Endotoxins	LAL Chromogenical test	< 20 EU/sample	Endotoxins: < 0.125 EU/sample	OK
Shelf life (accelerated and real aging)	ISO 10555-1 ISO 10555-4 ISO 11737-2	Meet the acceptance criteria according to tests of 10551-1 and ISO 10555-4 at 3 years Maintenance of sterility at 3 years	Luminor meet the requirements at 3 years real aging.	OK

21.3. Conclusions

Paclitaxel-eluting peripheral balloon dilatation catheters, *iVascular luminor* meet the required tests by the applicable standards for its verification, and thus it is considered to be passed.

22. WARRANTY

The device and all its components have been designed, manufactured, tested and packed with a maximum level of precaution. LVD Biotech warrants the device until its expiration date, provided the packaging is not broken, tampered with or damaged.

STANDARDIZED SYMBOLS FOR USE IN LABELLING MEDICAL DEVICES

(EN ISO 15223-1: 2016 STANDARD)

 Catalogue reference	 Batch no.	 Expiry date	 Keep away from sunlight
 Manufacturer	 Sterilized with ethylene oxide	 Limitation in storage temperature	 Keep dry

 Consult instructions for use	 Do not re-use	 Do not re-sterilize	 Do not use if the pack is damaged
 Not pyrogenic			

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