www.vacutestkima.it

"APPROVED"
by the Director of
Vacutest Kima S.r.l.

VACUTEST KIMA S.r.l. _____ Monica Vianello

Via dell'Industria, 12
33020 ARZERGRANDE (UD)
August 2019
Partita IVA 03450130288

INSTRUCTION FOR USE FOR MEDICAL PRODUCT

Vacutest Urine Containers, sterile

(Revision 2)

2019

www.vacutestkima.it

MEDICAL PRODUCT NAME

Vacutest Urine Containers, sterile (hereafter referred to as the Container(s)).

DESIGN VARIANTS

1. Vacutest Urine Container, 120 ml, with screw cap and integrated collection device for urine collection with vacuum tube
2. Vacutest Urine Container, 120 ml, with screw cap and integrated collection device for urine collection with vacuum tube, single pack
3. Vacutest Urine Container, 60 ml, with screw cap and integrated collection device for urine collection with vacuum tube, single pack
4. Vacutest Urine Container, 60 ml, with screw cap and integrated collection device for urine collection with vacuum tube
5. Vacutest Urine Container, 3 l, with screw cap and integrated transfer straw

INTENDED USE

Containers are designed for urine collection, storage and/or transportation of the urine sample for further laboratory analysis as a supporting aid for in-vitro diagnostic.

FIELD OF APPLICATION

Containers are designed for using in medical treatment and preventive care institutions.

INDICATIONS FOR USE

To be used according to prescription for all groups of patients for urine collection, storage and/or transportation for further analysis.

LIST OF POSSIBLE ANALYSIS WHERE CONTAINERS COULD BE USED

The list (including, but not limited) is presented in the Table 1.

Table 1

Name	Ref	Type of analysis
Vacutest Urine Container, 120 ml, with screw cap and integrated collection device for urine collection with vacuum tube	21820	biochemical analysis, clinical analysis, Nechiporenko test
Vacutest Urine Container, 120 ml, with screw cap and integrated collection device for urine collection with vacuum tube, single pack	21800	
Vacutest Urine Container, 60 ml, with screw cap and integrated collection device for urine collection with vacuum tube	21890	
Vacutest Urine Container, 60 ml, with screw cap and integrated collection device for urine collection with vacuum tube,	21883	

Name	Ref	Type of analysis
single pack		
Vacutest Urine Container, 3 l, with screw cap and integrated transfer straw	90476S	biochemical analysis, Zemnitskiy test

CONTRAINDICATIONS

- Do not use when damaging to package integrity (for containers in single packs).
- Do not use when damaging to hermetic state of the container: if there are signs of the container's opening, of damage to the protective label placed on the cap.

POTENTIAL USERS

Medical personnel.

SAFETY MEASURES

- Do not use after expiry date.
- Do not use if package integrity is damaged (for containers in single packs).
- Do not use if hermetic state of the container is damaged.
- Do not reuse. Single use only!
- Do not use if there is a sharp edge, protrusion or roughness on the surface of the container that can cut, pierce or tear gloves, or the user's finger.

USE GUIDELINES

- Determine the integrity of consumer packaging by visual inspection. Verify the best before date of the container. Before sampling, open the consumer packaging (if any), open the cap, fill the container with biomaterial, tightly close the cap.
Do not use if package (if any) integrity is damaged.
- Do not use when damaging to hermetic state of the container: if there are signs of the container's opening, of damage to the protective label placed on the cap.
- Use containers only for their intended purpose.

CAUTION! Containers are single use products.

OPERATING CONDITIONS

The product should be used at the temperature from +10°C to +35°C and relative are humidity 75% at 15°C and 100% at 25°C.

COMPATIBILITY

Containers are compatible with vacuum tubes for urine collection used for safe transfer of the sample from container.

Compatible products are listed in Table 2.

Table 2

No.	Ref. No.	Description	Manufacturer	Registration data
Products for urine collection, transportation and analysis, single-use BD Vacutainer in kits and single packs			Becton Dickinson and Company, USA	Registration Certificate No. ΦC3 2007/00888 dd. 28.01.2016
1.	368500	Vacutainer Urine Tube, without reagent, 4 ml, 13x75 mm		
2.	368501	Vacutainer Urine Tube, without reagent, 6 ml, 13x100 mm		
3.	364992	Vacutainer Urine Tube, with preservative, 8 ml, 16x100 mm		
4.	365017	Vacutainer Urine Tube, with preservative, 8 ml, 16x100 mm		
5.	365000	Vacutainer Urine Tube, without reagent, 9,5 ml, 16x100 mm		
6.	364938	Vacutainer Urine Tube, without reagent, 10 ml, 16x100 mm		
7.	364915	Vacutainer Urine Tube, without reagent, 11 ml, 16x100 mm		
8.	364958	Vacutainer Urine Tube, with preservative, 4 ml, 13x75 mm		
9.	364955	Vacutainer Urine Tube, with preservative, 10 ml, 16x100 mm		
Vacuette system for vacuum urine collection, with accessories			Greiner Bio-One GmbH, Austria	Registration Certificate No. ΦC3 2011/09315 dd. 23.01.2015
10.	456007	Vacuette Urine Tube, without reagent, 6,5 ml, 13x100 mm		
11.	455028	Vacuette Urine Tube, without reagent, 9,5 ml, 16x100 mm		
12.	455007	Vacuette Urine Tube, without reagent, 10 ml, 16x100 mm		
13.	456009	Vacuette Urine Tube, with preservative, 6,5 ml, 13x100 mm		
14.	454239	Vacuette Urine Tube, with preservative, 4 ml, 13x75 mm		
15.	455003	Vacuette Urine Tube, with preservative, 10 ml, 16x100 mm		
16.	455029	Vacuette Urine Tube, with preservative, 9,5 ml, 16x100 mm		
17.	455049	Vacuette Urine Tube, with preservative, 9,5 ml, 16x100 mm		
18.	455048	Vacuette Urine Tube, with preservative, 10 ml, 16x100 mm		
Polymer products for laboratory analysis			Ningbo Greetmed Medical Instruments Co., LTD., China	Registration Certificate No. ΦC3 2012/11857 dd. 28.03.2012
19.	12007050	Urine tube without reagent 9 ml 16x100 mm		
20.	12007051	Urine Tube, with preservative 9 ml 16x100 mm		
21.	12007052	Urine tube without reagent 9 ml 16x100 mm		

A. Use of containers with volume of 60 ml and 120 ml


- Clean the hands and genitals thoroughly before using the container.
- After clean the genitals thoroughly collect in a container the first morning micturition (as more concentrated), paying attention to avoid the first spout.
 - Urine should be transferred to medical staff to the laboratory within 1.5-2 hours after collection. Recommended storage temperature during this period is from +2 to +8 °C.

1. Open the cap by unscrewing anti-clockwise.

2. Lay the cap upside down on a clean surface.

3. Do not touch internal surfaces of the container and cap.


4. Collect the urine sample. Fill the container up to $\frac{3}{4}$ of the capacity.

5. Turn the cap tightly in a clockwise direction to seal.

Urine transfer to the tube

6. Gently shake the sample before transferring it to the tube.

7. Partially raise the protective label (do not remove it completely).

8. Insert the tube. Gently apply pressure to pierce of the tube cap. Keep the tube immersed until it is full. (end of flow).

9. Remove the tube and fully restick the protective label.

B. Use of containers with volume of 3l

Urine collection

- Clean the hands and genitals thoroughly before using the container. The first morning spout should be avoided. But the time of micturition is fixed.
- Further all the urine urinated during 24 hours from the time of the first micturition to the same hour in a day;

1. Open a big cap by screwing anti-clockwise, lay the cap upside down on a clean surface. Do not touch internal surfaces of the cap
2. Micturition shall be made directly to the container, its wide throat. After each micturition the cap should be screw tightly. Urine should be kept in a closed container in the refrigerator at bottom shelf, avoiding freezing at temperature from +2 to +8°C.
3. When finishing the collection (the last micturition at time as the first micturition was, but in a day), the urine should be transferred to the medical personnel or to the laboratory during 1,5-2 hours after collection completed. Recommended storage temperature from +2 to +8°C.

Urine transfer to the tube

1. Before transfer of urine to the tube gently homogenize the urine sample (mix the urine in the closed container by circular motions).
2. Lift the label placed on a small cap with integrated transfer straw and avoid to detach it completely.
3. Insert the tube. Gently press the tube to allow the cap piercing. Wait for the complete filling of the tube (stop of the flow).
4. Pull out the tube and restick the protective label.

General view of the products

Containers are hollow cylindrical vessels with screw cap. Containers are provided with graduation.

Technical specifications and description

- Containers do not contain additives.
- Containers assembled shall be hermetically sealed.
- A single commercial pack of the containers (if any) shall be hermetically sealed.
- On the surface of the container there should be no tears, dents, barbs, traces from the press mold with a height of more than 0.5 mm and foreign impurities protruding above the surface. No more than three foreign impurities with an area of no more than 0.25 mm² each are allowed.
- Container's cannula should be tightly fixed in the cap.
- Metal parts should be manufactured of corrosion-resistant materials.
- A window for information and graduation should be applied on the side surface of the container by injection casting.
- The shelf life must be at least 3 years from the date of sterilization.
- The container must remain sterile throughout its shelf life.
- During operation, transportation and storage the products should be resistant to mechanical and climate effects.
- Only materials the safety of which is duly confirmed should be used for manufacture.
- Products should be resistant to the effects of biological fluids and secretions of body tissues they come into contact during operation.

Materials used during manufacture

Materials used during manufacture are presented in table 3.

Table 3

Purpose	Material
Container	Polypropylene 100% (CAS No. 9003-07-0)
Cap	Polyethylene (LPNP) CLEARFLEX - FLEXIRENE
Integrated collection device for vacuum tube, Cannula (20G)	Cannula (20 G): stainless steel AISI 304 Base: Polypropylene (CAS No. 9003-07-0) Cover glove: isoprene rubber -100% isoprene rubber Glue: epoxy compound LOCTITE ABLESTIK 927-10E(ECCOBOND 927-10 E) Lubricant: silicone fluid Dow Corning® 360
Note - dyes: 1. Amber - AMBER PE 9/F090-TP 2. Yellow -YELLOW PE 3/R067-AT	

Dimensions and weight of the products are presented in tables 4 and 5.

Table 4

Name	Reference	Container's weight, g	Weight of cap with urine collection device, g	Weight of a bigger cap, g	Total weight, g
Urine container with screw cap and integrated transfer straw, 3 l	90476S	170,5±0,5	6,5±0,5	10,0±0,5	187,0 ±0,5
Vacutest Urine Container with screw cap and integrated collection device for urine collection with vacuum tube, 120 ml	21820	11,5±0,5	9,0±0,5	-	20,5 ± 0,5
Vacutest Urine Container with screw cap and integrated collection device for urine collection with vacuum tube, 120 ml, single pack	21800				
Vacutest Urine Container with screw cap and collection device for urine collection with vacuum tube, 60 ml	21890	6,5±0,5	6,0±0,5	-	12,5 ± 0,5
Vacutest Urine Container with screw cap and collection device for urine collection with vacuum tube, 60 ml, single pack	21883				

Table 5

Parameter	Reference No.		
	21800/21820	21890/21883	90476S
Diameter of base, mm	49	39	132
Height, mm	82	65	276
Nominal volume, ml	120	60	3000
Graduation, ml	20-40-60-80-100-120	10-20-30-40-50-60	100-500-1000-1500-2000-2500-3000
Indication minimum level, ml	50	50	100
Indication maximum level, ml	100	60	3000
Diameter of cap, mm	66,5	54	Large cap – 83 Small cap – 35
Diameter of hole in the cap for vaccum tube, mm	19	19	19
Height of cap, mm	16,1	15,2	Large cap – 15 Small cap – 12
Diameter of cannula (±0,1), mm	0,9		
Lenth of cannula, mm	18		
Note – allowed deviation from dimensions ±1 mm			

COMPLETENESS OF SETS

1. Product in single pack or without pack
2. User's manual, 1 pc. for a shipping case

STERILIZATION

Containers are single use products and released sterile. Internal volume of the container shall be sterile.

Way of sterilization: irradiation. Radiation dose: 15 kGy.

SHELF LIFE

Shelf life: 3 years from date of sterilization.

CONDITIONS OF STORAGE AND TRANSPORTATION

- Transportation and storage can be carried out by any type of transport in closed transport facilities in accordance with the rules applicable to each type of transport and approved in the prescribed manner. Transportation conditions: from -50 to +50 °C.
- Containers should be stored in the manufacturer's packaging in a dry heated room at temperature from +5 to +40 °C and with relative humidity of not more than 80% at a distance of at least 1 m from heating appliances. Avoid exposure to direct sunlight and significant temperature fluctuations during storage.
- After transportation in conditions of below-freezing temperatures, the products should be kept in a shipping case under storage conditions for at least 12 hours prior to use.

DISPOSAL

- Used products shall be disposed in accordance with current legislation.
- Used products relate to the class of epidemiological waste and shall be disposed in accordance with current legislation.
- Products not used for their intended purpose due to expiration of storage or other reasons relate to the class of epidemiological safe waste, which is close in composition to solid household waste, and disposed of as household waste.

WARRANTIES OF THE MANUFACTURER

- The manufacturer warrants that container conforms to the requirements (sterility, integrity of individual pack) if the user observes conditions of transportation and storage, as well as operating rules.
- Warranty period – 3 years from date of sterilization.

TECHNICAL SERVICE AND REPAIR

Containers are single use products. They are not subject to technical service and repair.

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AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE RF

Limited Liability Company "BioSystems" ("BioSystems" LLC),

197101, Russia, Saint Petersburg, 3 Pinsky lane, let. A

Tel.: (812) 320-49-49

E-mail: main@bioline.ru

LIST OF STANDARDS APPLIED BY THE MANUFACTURER OF MEDICAL PRODUCTS

UNI EN ISO 9001:2015 certificate no. 4264/4 issued by ICIM spa

UNI CEI EN ISO 13485:2016 certificate no. 4265/4 issued by ICIM spa

CE: quality guarantee system through the issuing of CE declaration of conformity according to CE Directive 98/79/CE

UNI CEI EN ISO 15223-1 Symbols to be used on medical devices labels.

UNI EN 14971 risk management application on medical devices

UNI EN 556-1 requirements for medical devices indicated as sterile

UNI EN ISO 11137-1 sterilization of sanitary products – Radiation – Part 1

UNI EN ISO 11737-2 microbiological methods – sterility tests made during sterilization process validation.

EN 14254 In vitro diagnostic medical devices - Single-use receptacles for the collection of specimens, other than blood, from humans

EN ISO 11607-1 Packaging of medical products subject to final sterilization. Part 1. Requirements for materials, sterilization barrier systems and packaging systems.

PACKAGE

Containers 21800, 21883 shall be individually packed in a transparent bag.

Containers in single packs and without them shall be placed in a shipping case.

Data on packaging materials are presented below in table 6.

Table 6

Packaging / Ref. No.	21800	21820	21890	21883	90476S
Single pack	Polypropylene F203T (CAS 9003-07-0)	-	-	Polypropylene F203T (CAS 9003-07-0)	-
Shipping case	200 pcs. Material: cardboard	400 pcs. Material: cardboard	500 pcs. Material: cardboard	500 pcs. Material: cardboard	30 pcs. Material: cardboard

MARKING

Individual marking of containers 21820, 21800, 21883, 21890:

Individual marking of containers 90476S:

VACUTEST[®] kima

www.vacutestkima.it

There is an additional marking (on single pack) for marking of containers in single pack (21800, 21883) in the RF:

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
For marking of shipping case in the RF an additional package (on transport case) shall be applied:

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IMPORTANT MARKING SYMBOLS

Symbol	Meaning
	Manufacturer
	Best before...
	Lot number
	Reference number
	Irradiation sterilization
	In vitro diagnostic medical device
	Obligatory mark of conformity of products placed at the EEA market
	Single use
	Manufacture date

Symbol	Meaning
	Keep dry
	Avoid sunlight
	Top
	Bottom
	Side
	Take care! See user's manual
	Temperature range
	Fragile. Handle carefully
	Country of origin: Italy
	Keep away from children
	do not use if packaging is damaged

CAMERA DI COMMERCIO, INDUSTRIA, ARTIGIANATO E AGRICOLTURA DI PADOVA
PADOVA CHAMBER OF COMMERCE, INDUSTRY, HANDCRAFT AND AGRICULTURE

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In relation to the instance on record made by Mrs. MONICA VIANELLO, proprietor/legal representative/attorney of the Company VACUTEST KIMA SRL, we certify that the said applicant is a duly authorized signatory for any documents and acts to be used in trade with foreign countries.

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