


Dentsply Sirona Switzerland
Group of Companies
Chemin du Verger 3
1338 Ballaigues
Switzerland

info@dentsplysirona.com
www.dentsplysirona.com

DIRECTIONS FOR USE
PROTAPER GOLD ENDODONTIC STERILE INSTRUMENTS OF DIFFERENT SIZES
FOR DENTAL USE ONLY

Manufactured by
Maillefer Instruments Holding Sàrl
Chemin du Verger 3
Ballaigues, CH-1338
Switzerland

Luc Etienne	Antonio Santacroce
	
Vice President QA/RA Maillefer Instruments Holding Sàrl Ballaigues	Finance Director Maillefer Instruments Holding Sàrl Ballaigues
01.10.2019	01.10.2019
	

Content

1. DESCRIPTION	3
2. PHOTOS.....	
3. COMPOSITION	3
4. INDICATIONS FOR USE	3
5. CONTRAINDICATIONS	3
6. WARNINGS.....	3
7. PRECAUTIONS	3
8. ADVERSE REACTIONS.....	4
9. STEP BY STEP INSTRUCTIONS FOR PROTAPER GOLD FILES.....	4
10. DISINFECTION, CLEANING AND STERILIZATION.....	5
11. TERMS OF STORAGE, TRANSPORTATION AND USAGE	8
12. INFORMATION ABOUT DISPOSAL AND DESTRUCTION.....	8
13. EXPECTED LIFE OF THE PRODUCTS	9
14. SHELF-LIFE.....	9
Appendix to the direction for use of Protaper Gold endodontic sterile instruments of different sizes..	10

1. DESCRIPTION
Protaper Gold endo
Protaper Gold e
Forms of del
the direc

1. DESCRIPTION.

Protaper Gold endodontic sterile instruments Shaping Files, sizes SX, S1, S2; length 19, 21, 25, 31 mm.

Protaper Gold endodontic sterile instruments Finishing Files, sizes F1, F2, F3, F4, F5; length 21, 25, 31 mm.

Forms of delivery of Protaper Gold endodontic sterile instruments of different sizes see in the appendix to the direction for use.

2. PHOTOS

Protaper Gold instruments Shaping Files SX	
Protaper Gold instruments Shaping Files S1	
Protaper Gold instruments Shaping Files S2	
Protaper Gold instruments Finishing Files F1	
Protaper Gold instruments Finishing Files F2	
Protaper Gold instruments Finishing Files F3	
Protaper Gold instruments Finishing Files F4	
Protaper Gold instruments Finishing Files F5	

3. COMPOSITION

The cutting part of these instruments is made of a nickel-titanium alloy.

Instruments are sterile. Sterilization method: Gamma Radiation.

4. INDICATIONS FOR USE

These instruments are to be used only in a clinical or hospital environment, by qualified users.

Scope: Dentistry. Application field: for the removal of dentin and shaping of the root canal.

5. CONTRAINDICATIONS

As with all mechanically driven root canal instruments, Protaper Gold files should not be used in cases of severe and sudden apical curvatures due to heightened risk of separation.

6. WARNINGS

This product contains nickel and should not be used for individuals with known allergic sensitivity to this material.

7. PRECAUTIONS

- Straight-line access is a prerequisite for proper root canal treatment, Protaper Gold files are no exception.
- Multiple use disinfection and re-sterilization cycles may lead to increased risk of file separation.
- These instruments should not be immersed in a sodium hypochlorite solution.
- Instrument reprocessing: follow the reprocessing instructions on part 10.

- Clean the flutes frequently during instrumentation, inspecting for signs of distortion or wear, uneven flutes, dull spots.
- Frequently irrigate, recapitulate and irrigate the canal throughout the procedure, minimally after each file.
- Protaper Gold files should only be used in regions of the canal that have a confirmed and reproducible glide path. Establish a reproducible glide path using hand files, at least an ISO 015 size.
- Use the Shaping Files (S1, S2 and SX) with a brushing action on the withdrawal stroke in order to create straight line radicular access.
- Use the Finishing Files (F1, F2, F3, F4 and F5) with no brushing action.
- Use the appropriate finishing files to passively follow the canal to the working length, and then withdraw immediately.
- Protaper Gold files are manufactured with a process that results in a file that has a Gold appearance. Due to this proprietary processing, ProTaper Gold files may appear slightly curved. This is not a manufacturing defect. While the file can be easily straightened using only your fingers, it is not necessary to straighten the file prior to use. Once inside the canal, the Protaper Gold file will follow the anatomy.
- Always use minimal apical pressure. Never force the files down the canal.
- For optimal usage, torque control devices are recommended.
- The Protaper Gold rotary files can be used at motor speeds between 250 rpm and 350 rpm.

Recommended motor settings:

Protaper Gold		
File Size	Speed [rpm]	Torque [N•cm]
Protaper Gold S1 and SX	300	5.10
Protaper Gold S2 and F1	300	1.50
Protaper Gold F2, F3, F4, F5	300	3.10

The speed and torque settings indicated in the above table are for example only and may vary according to each user preferences and motor capabilities.

8. ADVERSE REACTIONS

As with all mechanically driven root canal instruments, Protaper Gold files should not be used in cases of severe and sudden apical curvatures due to heightened risk of separation.

9. STEP BY STEP INSTRUCTIONS FOR PROTAPER GOLD FILES

6.1 Radiographic Evaluation

Review different horizontally angulated radiographs to diagnostically determine the width, length, and curvature of any given root and canal.

6.2 Access Preparation

Create straight-line access to the canal orifice(s) with emphasis on flaring, flattening, and finishing the internal axial walls.

6.3 Protaper Gold SHAPING TECHNIQUE

The crown down technique is the technique of choice for rotary instruments.

- Create straight-line access to canal orifice.
- In the presence of a viscous chelator (such as Glyde File Prep root canal conditioner) passively scout the coronal 2/3 with 10 and 15 hand files. Gently work these instruments until a smooth, reproducible glide path is confirmed. Alternatively, mechanized glide path files (such as Proglider or Pathfiles) may be used after a 10 hand file.
- In the presence of NaOCl, "float" the S1 in the canal and passively "follow" the glide path. Before light resistance is encountered, laterally "brush" and cut dentin on the outstroke to improve straight-line access and apical progression. Always brush away from the furcation.

- Continue shaping with S1 as described until the depth of the 15 hand file is reached.
- Use the S2, exactly as described for the S1, until the depth of the 15 hand file is reached.
- In the presence of a viscous chelator or NaOCl, scout the apical 1/3 with 10 and 15 hand files and gently work them until they are loose at length.
- Establish working length, confirm patency and verify the presence of a smooth reproducible glide path in the apical 1/3.
- Use the S1, with a brushing action, until working length is reached.
- Use the S2, with a brushing action, until working length is reached.
- Reconfirm working length, irrigate, recapitulate and re-irrigate, especially in more curved canals.
- Use Finishing File F1, in a "non-brushing" action, with each insertion deeper than the previous insertion until working length is reached. Do not leave the file at working length for longer than one second.
- Gauge the foramen with a 20 hand file. If the instrument is snug at length, the canal is shaped and ready to be obturated.
- If the 20 hand file is loose at length, proceed to the F2 and, when necessary the F3, F4 and F5, with the same non-brushing motion to working length, gauging after each Finishing file with 25, 30, 40 or 50 hand files respectively.
- If necessary, use the SX with a brushing motion to move the coronal aspect of the canal away from furcal concavities and/or to create more coronal shape. SX can also be used to optimally shape canals in shorter roots.
- The Protaper Gold sequence is the same regardless of the length, diameter or curvature of the canal.

10. DISINFECTION, CLEANING AND STERILIZATION

Reprocessing procedure for dental instruments.

I - FOREWORD

Devices that are marked as "sterile" do not require any specific treatment before the first use.

For those devices that are not labelled "single use", re-processing of the devices should be carried out as per this DFU. For hygiene and sanitary safety purposes, these instruments must be cleaned and sterilized before each re-use to prevent any contamination.

II - GENERAL RECOMMENDATION

- 1) Use only a detergent solution, with disinfecting effect, which is approved for its efficacy (VAH/DGHH-listing, CE marking, FDA approval) and in accordance with the DFU of the detergent solution manufacturer. For all metal devices, it is recommended to use anticorrosion disinfecting and cleaning agents.
- 2) For your own safety, please wear personal protective equipment (gloves, glasses, mask).
- 3) The user is responsible for the cleaning and sterilization of the product for the first cycle and each further usage as well as for the usage of damaged or dirty devices where applicable after sterilization.
- 4) It is safest for the practitioner to use our devices only once. Should our devices be reused, we recommend that they should not be used more than 5 times. After each processing they should be carefully inspected before use: the appearance of defects such as deformations (bent, unwound), breakage, corrosion, loss of color coding or marking, indicate that the devices are not able to fulfil the intended use with the required safety level and must therefore be discarded.

For our root canal shaping instruments we recommend not to exceed the following maximum number of uses:

Type of canal	Stainless Steel instruments with a Diameter $\leq$ ISO 015	Stainless Steel instruments with a Diameter $>$ ISO 015	NiTi instruments
Extremely curved ($>30^\circ$) or S-shaped	1 canal max.	2 canals max.	2 canals max.

canals			
Moderately curved canals (10° to 30°)	1 canal max.	4 canals max.	4 canals max.
Slightly curved (<10°) or straight canals	1 canal max.	8 canals max.	8 canals max.

- 5) For the final rinsing step deionised water use is mandatory, whether using an automated washer-disinfector or a manual cleaning method. Tap water is permissible for the other rinsing steps.
- 6) NiTi instruments should not be used with Hydrogen Peroxide (H₂O₂) solution which is known to degrade them.
- 7) Only the active part of the NiTi instrument, which is in contact with the patient should be immersed in a NaOCl solution concentrate at NOT more than 5%.
- 8) Avoid device to dry out, prior to, or during pre-disinfection, or cleaning. Dried biological material can be difficult to remove.
- 9) Use only device appropriated support for reprocessing.
- 10) Do not use label systems or identification markers directly on the device.

III - STEP-BY-STEP INSTRUCTIONS

	Operation	Activities	Warning and remarks
1.	Disassembling	- Disassemble the device, if applicable.	- Remove and discard silicone stops.
2.	Pre-Disinfection	- Soak all devices immediately after use in a disinfection solution (We recommend the use of Prolystica 2X Concentrate Enzymatic Presoak and Cleaner at 0.4% for a minimum of 15 minutes). Use a tray made from high density polyethylene or stainless steel.	- Follow instructions and respect concentrations and immersion times given by the manufacturer (an excessive concentration may cause corrosion or others defects on devices). - The pre-disinfection solution should be a specific solution targeted by the supplier for pre-disinfection. It should be used at the dilution specified by the supplier. It should contain, or be combined with a proteolytic enzyme. - The pre-disinfection solution should be aldehyde free (to avoid blood impurities fixation) and without di- or triethanolamines as corrosion inhibitor. Change the pre-disinfection solution regularly i.e. When it becomes soiled, or when efficacy is diminished due to exposure to microbial loads. - Do not use pre-disinfecting solutions containing Phenol or any products, which are not compatible with the devices. - For visible impurities observed on instruments a pre-cleaning is recommended with a soft brush (made from either nylon, polypropylene, acrylic). Manually brush the device until visible impurities are removed.
3.	Rinsing	- Abundant rinsing (at least 1 min) under running water (ambient temperature).	- Use tap water for rinsing. - If a pre-disinfectant solution contains a corrosion inhibitor, it is recommended to do the rinsing step just before starting the cleaning step.
4a.	Automated Cleaning with washer- disinfector	- Place the devices in a kit, support, or container (made from stainless steel or titanium) to avoid any contact between devices or posts. - Place the devices in the washer-disinfector and execute the defined cycle (Ao value > 3000 or, at least 5 min at 90°C (194°F)). - Use a detergent solution with cleaning properties (we recommend Neodisher Mediclean Forte at 0.4%).	- Discard any devices with defects (broken, bent,...). - Avoid any contact between instruments or posts when placing in the washer-disinfector use kits, supports or containers. - Follow instructions and concentrations given by the manufacturer of the detergent solution. - Follow the instructions of the washer-disinfector and verify the success criteria after each cycle have been met as stated by the manufacturer. - The final rinse step should be with deionized water. For other steps follow the water quality defined by the manufacturer. - Use only approved washer-disinfector according to EN ISO 15883, maintained and validated

			regularly. - It is recommended to use an alkaline detergent with tensides, which has grease removal, disinfection (against bacteria/ fungi) and corrosion inhibition properties. The detergent should be approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and used in accordance with its DFU. The detergent should be aldehyde free and without di- or triethanolamines as corrosion inhibitor.
OR			
4b.i	Manual Cleaning assisted by an ultrasonic device	- Place the devices in a kit, support or container (made from stainless steel, polypropylene or titanium) to avoid any contact between devices. - Immerse in the detergent solution with cleaning properties (we recommend Neodisher Mediclean Forte at 2%), assisted by an ultrasonic device if suitable for at least 15 min.	- No visible impurities should be observed on the devices. - If visible impurities are observed on the devices, the device must be manually brushed with a soft brush (made from either nylon, polypropylene, acrylic) until visible impurities are removed. - Discard any devices with defects (broken, bent, and unwound). - Follow instructions, observe water quality, concentrations and cleaning time stated by the manufacturer of the cleaning solution. - It is recommended to use an alkaline detergent with tensides, which has grease removal, disinfection (against bacteria/ fungi) and corrosion inhibition properties. The detergent should be approved for its efficacy (VAH/DGHM-listing, CE marking, FDA approval) and used in accordance with the DFU of the detergent solution manufacturer). - The detergent should be aldehyde free and without di- or triethanolamines as corrosion inhibitor.
4b.ii	Rinsing	- Abundant rinsing (at least 1 min) under running water (ambient temperature).	- Use deionised water for rinsing. - If the previously used cleaning solution contains a corrosion inhibitor, it is recommended to do the rinsing step just before starting the autoclaving.
4b.iii	Drying	- Devices should be thoroughly dried before inspection and packaging.	- Dry with a single use non-woven cloth. - Devices should be dried until visual traces of moisture are eliminated. - Particular attention has to be paid to effectively dry joints or cavities within a device.
5.	Inspection	- If applicable assemble the devices (including the placement of new silicon stops). - Inspect the devices functionality. - Visually inspect devices with naked eye under appropriate lighting (min 500 lux) and sort out those with defects.	- Dirty devices must be cleaned again. - Do not re-use silicon stops. - Discard devices, which show any defect as described in the General Recommendation above (point 4).
6.	Packaging	- Place the devices in a kit, support or container to avoid any contact between instruments or posts and pack the devices in "Sterilization pouches".	- Device must be double-packaged using paper-plastic pouches for steam sterilization prior sterilization. Ensure that the pouches are suitable for steam sterilization and were validated and manufactured as per ISO 11607 and EN 868-5. - Use an appropriate packaging, moist-heat resistant (141°C, 286°F) and compliant with ISO 11607. - Avoid any contact between instruments or posts during sterilization. Use kits, supports or containers. - For sharp devices that are not contained within a box, silicon tubes should be placed around the devices to prevent packaging piercing. - Seal the pouches according to the

			recommendation of the pouch manufacturer. If a thermo-sealer is used, the process must be validated and the thermosealer must be calibrated and qualified. - Check the validity period of the pouch given by the pouch manufacturer to determine the shelf life.
7.	Sterilization	We recommend a steam sterilization at 134°C / 273.2°F during 18 minutes for the purpose of de-activating potential prions.	- The instruments and posts must be sterilized according to the packaging labelling. - When sterilizing multiple instruments in one autoclave cycle ensure that the sterilizer's maximum load is not exceeded. - Place the pouches in the steam sterilizer according to the recommendation given by the sterilizer manufacturer. - Use only Pre-Vacuum air Removal steam sterilizer that are matching the requirements of EN 13060 (class B, small sterilizer) and EN 285 (full size sterilizer), with saturated steam. - Use a validated sterilization procedure according to ISO 17665 with a minimum drying time of 20 min. - Respecting the maintenance procedure of the sterilizer is under the responsibility of the owner and should be performed following the requirements for medical devices sterilization (examples: planning of maintenance, qualification, acceptance criteria of condensate and water as per EN 285, annex 2). - Control the efficiency and acceptance criteria of the sterilization procedure (packaging integrity, no humidity, no colour change of packaging, positive physico-chemical indicators, conformity of actual cycle parameters, to reference cycle parameters). A special attention should be paid to the packaging integrity if the sterilization cycle 134°C (273.2°F), 18 minutes was used. - Store traceability records and define shelf-life according to packaging manufacturer guidelines. - Shorter sterilization cycles according to local regulations are possible but are not guaranteed to de-activate prions.
8.	Storage	- Keep devices in sterilization packaging in a clean environment, away from sources of moisture and direct sunlight. Store at temperature 0°C - +30°C; humidity 50 – 60%.	- After sterilization, the product should be manipulated with care in order to keep the integrity of the packaging (sterile barrier). - Sterility cannot be guaranteed if packaging is open, damaged or wet. - Check the packaging and the medical devices before using them (packaging integrity, no humidity and use by date). In case of damage, a complete rework should be performed.

11. TERMS OF STORAGE, TRANSPORTATION AND USAGE

Storage terms	Temperature 0°C - +30°C; humidity 50 – 60%, keep away from sunlight
Transportation terms	Temperature -20°C - +60°C, humidity 50 – 60 %, keep away from sunlight.
Usage terms	These instruments are to be used only in a clinical or hospital environment, by qualified users. Temperature 0°C - +30°C; Humidity 50 – 60%, keep away from sunlight. Instruments are properly operated inside the oral cavity, in the presence of biological fluids in it at nominal temperatures: +32 °C - +42 °C

12. INFORMATION ABOUT DISPOSAL AND DESTRUCTION

Hazard Class Waste - Class B.

13. EXPECTED
the expected life was
Extremely
Moderately
Slightly


13. EXPECTED LIFE OF THE PRODUCTS

The expected life was defined depending on the type of canals to be treated

- Extremely curved (>30°) or S-shaped canals: 2 canals maximum;
- Moderately curved canals (10° to 30°): 4 canals maximum;
- Slightly curved (<10°) or straight canals: 8 canals maximum.

14. SHELF-LIFE

Shelf-life of the product: 5 years.

Symbols	EN
	Manufacturer
	Reference number
	Lot number
	Expiry date
	Date of manufacture
	QR-Code
	CE mark
	Sterilized by irradiation
	Sterilizable in a steam sterilizer (autoclave) at the specified temperature
	Do not use if package is damaged
	Clockwise rotation
	Handle right angle RA
	Nickel Titanium
	Silicone
	Not returnable if seal is broken

Consult instructions for use

Manufacturer:

Maillefer Instruments Holding Sàrl
Chemin de Verger, 3,
CH-1338 Ballaigues
Switzerland
www.dentsplymaillefer.com

Authorized Representative in the Russian Federation:

«Dentsply Sirona» LLC, 115432, Moscow, Andropov prospect 18, building 6, offices 1 39 40.
Telephone: +7 495 725 10 87.

Appendix to the direction for use of Protaper Gold endodontic sterile instruments of different sizes

Forms of delivery of Protaper Gold endodontic sterile instruments of different size:

1. Protaper Gold endodontic sterile instruments Assortment in blisters 6 pcs. SX - F3 21 mm: Protaper Gold shaping instruments (Shaping Files), length 19 mm, size SX – 1 pc.; Protaper Gold shaping instruments (Shaping Files), length 21 mm, size S1 – 1 pc., S2 – 1 pc.; Protaper Gold finishing instruments (Finishing Files), length 21 mm, size F1 – 1 pc., F2 – 1 pc., F3 – 1 pc.
2. Protaper Gold endodontic sterile instruments Assortment in blisters 6 pcs. SX - F3 25 mm: Protaper Gold shaping instruments (Shaping Files), length 19 mm, size SX – 1 pc.; Protaper Gold shaping instruments (Shaping Files), length 25 mm, size S1 – 1 pc., S2 – 1 pc.; Protaper Gold finishing instruments (Finishing Files), length 25 mm, size F1 – 1 pc., F2 – 1 pc., F3 – 1 pc.
3. Protaper Gold endodontic sterile instruments Assortment in blisters 6 pcs. SX - F3 31 mm: Protaper Gold shaping instruments (Shaping Files), length 19 mm, size SX – 1 pc.; Protaper Gold shaping instruments (Shaping Files), length 31 mm, size S1 – 1 pc., S2 – 1 pc.; Protaper Gold finishing instruments (Finishing Files), length 31 mm, size F1 – 1 pc., F2 – 1 pc., F3 – 1 pc.
4. Protaper Gold endodontic sterile shaping instruments (Shaping Files) in blisters 6 pcs., length 19 mm, size SX.
5. Protaper Gold endodontic sterile shaping instruments (Shaping Files) in blisters 6 pcs., length 21 mm, size S1.
6. Protaper Gold endodontic sterile shaping instruments (Shaping Files) in blisters 6 pcs., length 21 mm, size S2.
7. Protaper Gold endodontic sterile shaping instruments (Shaping Files) in blisters 6 pcs., length 25 mm, size S1.
8. Protaper Gold endodontic sterile shaping instruments (Shaping Files) in blisters 6 pcs., length 25 mm, size S2.
9. Protaper Gold endodontic sterile shaping instruments (Shaping Files) in blisters 6 pcs., length 31 mm, size S1.
10. Protaper Gold endodontic sterile shaping instruments (Shaping Files) in blisters 6 pcs., length 31 mm, size S2.
11. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 21 mm, size F1.
12. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 21 mm, size F2.

13. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 21 mm, size F3.
14. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 21 mm, size F4.
15. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 21 mm, size F5.
16. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 25 mm, size F1.
17. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 25 mm, size F2.
18. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 25 mm, size F3.
19. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 25 mm, size F4.
20. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 25 mm, size F5.
21. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 31 mm, size F1.
22. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 31 mm, size F2.
23. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 31 mm, size F3.
24. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 31 mm, size F4.
25. Protaper Gold endodontic sterile finishing instruments (Finishing Files) in blisters 6 pcs., length 31 mm, size F5.

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