SECRETARY OF STATE
STATE OF INDIANA

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Connie Lawson
Secretary of State


APOSTILLE

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1. Country: United States Of America
2. This public document has been signed by *Sasha L. Brown*
3. Acting in the capacity of notary public in & for *Greene County*
4. And bears the seal/stamp of notary public in & for the State of Indiana

CERTIFIED

5. at Indianapolis, Indiana
6. This *Eighteenth* day of *June*, 2019
7. by the Secretary of State of Indiana
8. No. A2019-0618110221
9. Seal/Stamp:
10. Signature:

Connie Lawson

Connie Lawson
Indiana Secretary of State


*Effective May 1, 2003 all apostilles from the Indiana Secretary of State will have an electronically printed seal.
This document was revised March 17, 2012.*

COOK INCORPORATED
750 DANIELS WAY
BLOOMINGTON, IN 47404 U.S.A.
PHONE: 812.339.2235 TOLL FREE: 800.457.4500
WWW.COOKMEDICAL.COM

CERTIFICATION STATEMENT

As a responsible representative of Cook Incorporated, I confirm this is an accurate version of the Instruction for Use for "Stents to restore and maintain patency of the vessels" Submission.

Janice Mercado
Manager, Regulatory Affairs
Cook Incorporated


11 JUN 2019
Date

United States of America)
State of Indiana) ss:
County of Greene)

Subscribed and sworn to before me this 11 day of June 2019.

My commission expires on October 06, 2022.

Sasha L Brown
Notary Public, State of Indiana
Greene County
Commission# 658997

COOK INCORPORATED
750 DANIELS WAY
BLOOMINGTON, IN 47404 U.S.A.
PHONE 812.339.2235 TOLL FREE 800.457.4500
WWW.COOKMEDICAL.COM

INSTRUCTION FOR USE

Stents to restore and maintain patency of the vessels

Formula 418 Vascular Stent

DEVICE DESCRIPTION

The Formula 418 Vascular Stent is a balloon-expandable stent made of 316L stainless steel with a slotted tube configuration. It is pre-mounted on an over-the-wire balloon catheter, which serves as the delivery system (Formula 418 Vascular Balloon-Expandable Stent System). The stent is positioned between two radiopaque marker bands, which are located inside the balloon at the proximal and distal tapers of the balloon. The cannula design of the stent provides a low outside diameter profile, which permits use with a 5.0 French sheath or 6.0 French guiding catheter for 3-6 mm diameter sizes and a 6.0 French sheath and 7 French guiding catheter for 7-8 mm diameter sizes. The Formula 418 Vascular Balloon-Expandable Stent System is pre-mounted on 80 cm and 135 cm over-the-wire balloon catheter delivery systems. The stent is available in nominal expanded diameters of 3, 4, 5, 6, 7, and 8 mm with lengths of 12, 16, 20, 24 and 30 mm each. (Table 1)

Table 1: Device Specifications (diameters and lengths)

Order Number	Catheter Length (cm)	Stent Diameter (mm)	Stent Length (mm)	Minimum Catheter/Sheath* Diameter (inches)	Guiding Inside	Guiding Catheter/Sheath (French size)
FOV418-18-80-3-12	80	3	12	0.068		6.0/5.0
FOV418-18-80-3-16	80	3	16	0.068		6.0/5.0
FOV418-18-80-3-20	80	3	20	0.068		6.0/5.0
FOV418-18-80-3-24	80	3	24	0.068		6.0/5.0
FOV418-18-80-3-30	80	3	30	0.068		6.0/5.0
FOV418-18-80-4-12	80	4	12	0.068		6.0/5.0
FOV418-18-80-4-16	80	4	16	0.068		6.0/5.0
FOV418-18-80-4-20	80	4	20	0.068		6.0/5.0
FOV418-18-80-4-24	80	4	24	0.068		6.0/5.0
FOV418-18-80-4-30	80	4	30	0.068		6.0/5.0
FOV418-18-80-5-12	80	5	12	0.068		6.0/5.0
FOV418-18-80-5-16	80	5	16	0.068		6.0/5.0
FOV418-18-80-5-20	80	5	20	0.068		6.0/5.0
FOV418-18-80-5-24	80	5	24	0.068		6.0/5.0
FOV418-18-80-5-30	80	5	30	0.068		6.0/5.0
FOV418-18-80-6-12	80	6	12	0.068		6.0/5.0
FOV418-18-80-6-16	80	6	16	0.068		6.0/5.0
FOV418-18-80-6-20	80	6	20	0.068		6.0/5.0
FOV418-18-80-6-24	80	6	24	0.068		6.0/5.0
FOV418-18-80-6-30	80	6	30	0.068		6.0/5.0
FOV418-18-80-7-12	80	7	12	0.078		7.0/6.0
FOV418-18-80-7-16	80	7	16	0.078		7.0/6.0
FOV418-18-80-7-20	80	7	20	0.078		7.0/6.0
FOV418-18-80-7-24	80	7	24	0.078		7.0/6.0
FOV418-18-80-7-30	80	7	30	0.078		7.0/6.0
FOV418-18-80-8-12	80	8	12	0.078		7.0/6.0
FOV418-18-80-8-16	80	8	16	0.078		7.0/6.0
FOV418-18-80-8-20	80	8	20	0.078		7.0/6.0
FOV418-18-80-8-24	80	8	24	0.078		7.0/6.0
FOV418-18-80-8-30	80	8	30	0.078		7.0/6.0
FOV418-18-135-3-12	135	3	12	0.068		6.0/5.0
FOV418-18-135-3-16	135	3	16	0.068		6.0/5.0
FOV418-18-135-3-20	135	3	20	0.068		6.0/5.0
FOV418-18-135-3-24	135	3	24	0.068		6.0/5.0
FOV418-18-135-3-30	135	3	30	0.068		6.0/5.0
FOV418-18-135-4-12	135	4	12	0.068		6.0/5.0

FOV418-18-135-4-16	135	4	16	0.068	6.0/5.0
FOV418-18-135-4-20	135	4	20	0.068	6.0/5.0
FOV418-18-135-4-24	135	4	24	0.068	6.0/5.0
FOV418-18-135-4-30	135	4	30	0.068	6.0/5.0
FOV418-18-135-5-12	135	5	12	0.068	6.0/5.0
FOV418-18-135-5-16	135	5	16	0.068	6.0/5.0
FOV418-18-135-5-20	135	5	20	0.068	6.0/5.0
FOV418-18-135-5-24	135	5	24	0.068	6.0/5.0
FOV418-18-135-5-30	135	5	30	0.068	6.0/5.0
FOV418-18-135-6-12	135	6	12	0.068	6.0/5.0
FOV418-18-135-6-16	135	6	16	0.068	6.0/5.0
FOV418-18-135-6-20	135	6	20	0.068	6.0/5.0
FOV418-18-135-6-24	135	6	24	0.068	6.0/5.0
FOV418-18-135-6-30	135	6	30	0.068	6.0/5.0
FOV418-18-135-7-12	135	7	12	0.078	7.0/6.0
FOV418-18-135-7-16	135	7	16	0.078	7.0/6.0
FOV418-18-135-7-20	135	7	20	0.078	7.0/6.0
FOV418-18-135-7-24	135	7	24	0.078	7.0/6.0
FOV418-18-135-7-30	135	7	30	0.078	7.0/6.0
FOV418-18-135-8-12	135	8	12	0.078	7.0/6.0
FOV418-18-135-8-16	135	8	16	0.078	7.0/6.0
FOV418-18-135-8-20	135	8	20	0.078	7.0/6.0
FOV418-18-135-8-24	135	8	24	0.078	7.0/6.0
FOV418-18-135-8-30	135	8	30	0.078	7.0/6.0

*See manufacturer specifications for French (Fr) equivalent.

INTENDED USE

The Formula 418 Vascular Balloon-Expandable Stent System is intended for stenting of peripheral arteries.

CONTRAINDICATIONS

The Formula 418 Vascular Balloon-Expandable Stent System is contraindicated for use in:

- Patients for whom antiplatelet and/or anticoagulant therapy is contraindicated.
- Patients who have a lesion that cannot be crossed with a wire or a balloon angioplasty catheter.
- Patients who have stenoses that cannot be dilated to permit passage of the stent.
- Patients with bleeding disorders.
- Stenting of an arterial vessel where leakage from the artery could be exacerbated by placement of a stent.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.

WARNINGS

- The use of this device carries the associated risks of subacute thrombosis, vascular complications and/or bleeding events. Judicious selection of patients is necessary.
- Persons allergic to 316L stainless steel may suffer an allergic reaction to this implant.
- The stent should not be placed in a location with significant external compression which could crush the stent.
- The stent should not be deployed in areas subject to significant flexure or bending to reduce the risk of stent deformations or fracture.

PRECAUTIONS

General

- Implantation of the stent should be performed only by physicians who have received adequate training. Adequate training consists of specific training in vascular stenting techniques and procedures under the supervision of an interventional physician experienced in vascular stent implantation.

- Stent placement should only be performed at hospitals where emergency artery bypass graft surgery can be readily performed.
- Subsequent restenosis may require redilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of an endothelialized Formula 418 Vascular Stent is unknown at present.

Stent Placement

- Manipulation of the Formula 418 Vascular Balloon-Expandable Stent System requires fluoroscopic control.
- The stent is preloaded on its delivery balloon and should not be removed. Do not attempt to place stent on another balloon catheter for deployment.
- Deploy the stent over a wire guide not exceeding 0.018 inch diameter.
- DO NOT wipe or clean stent or catheter with organic solvent (i.e. isopropyl alcohol).
- DO NOT expose catheter to temperatures above 54°C (130°F).
- DO NOT use power injection systems with the delivery system.
- DO NOT rotate any part of the system during deployment.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is particularly important during removal of the catheter from packaging, placement over the wire guide and advancement through the large-bore Tuohy-Borst "Y" adapter and guiding catheter hub.
- Do not pre-inflate balloon prior to stent deployment (See Preparation of Balloon Catheter section in these Instructions for Use).
- Use only the appropriate balloon inflation medium. The standard inflation medium is a 1:1 mixture of contrast medium and normal saline. Do not use air or any gaseous substance as a balloon inflation medium.
- Expansion of the stent should not be undertaken if the stent is not properly positioned in the vessel (See Removal of Unexpanded Stent section in these Instructions for Use).
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the product label. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible intimal damage and dissection.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stent and may cause acute closure of the vessel, requiring additional intervention (further dilatation or placement of additional stents).
- Use prior to the expiration date on the package label.
- Intended for one-time use (Do not re-sterilize).
- The use of mechanical atherectomy devices within an implanted stent may cause damage to the stent.
- Repositioning of the device after deployment is not possible.

MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Formula 418 Vascular Stent is **MR Conditional** according to ASTM F2503, *Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment*. A patient with this stent can be scanned safely anytime after placement under the following conditions.

- Static magnetic field of 3-Tesla or 1.5-Tesla
- Maximum spatial magnetic gradient of 720 Gauss/cm or less
- Product of the spatial gradient and the static magnetic field of 21.6 T²/m or less
- MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg or less for 15 minutes of scanning

Non-clinical evaluation was conducted in an MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter

in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual). In non-clinical testing, a single and two overlapped Formula 7 mm x 20 mm stents produced a maximum temperature rise of 2.6 °C and 3.3 °C, respectively, during 15 minutes of MRI (i.e., for one scanning sequence) performed in a 3 Tesla MR system (3 Tesla/128 MHz, Excite, Software G3.0-052B, General Electric Healthcare) at an MR system reported whole-body-averaged SAR of 3.0 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.8 W/kg).

The effect of heating in the MRI environment for stents with fractured struts is unknown.

Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 15 mm of the Formula 7 mm x 20 mm stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

MedicAlert Foundation

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted at:

Mail: MedicAlert Foundation International

2323 Colorado Avenue

Turlock, CA 95382

Phone: +1-888-633-4298 (toll free)

+1-209-668-3333 (from outside the US)

Fax: +1-209-669-2450

Web: www.medicalert.org

POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of a vascular stent are presented in the following alphabetical list and include, but are not limited to:

- Abscess
- Allergic reaction to stainless steel or contrast agents
- Arrhythmias (ventricular fibrillation, ventricular tachycardia, other)
- Arteriovenous fistula
- Bowel infarct
- Death
- Dialysis
- Dissection
- Drug reaction to antiplatelet agents
- Drug reaction, allergic reaction to contrast media
- Emboli (air, tissue, or thrombotic emboli) resulting in tissue ischemia/infarction
- Emergency surgery to correct vascular complications
- Emergent renal artery bypass surgery
- Extremity ischemia/amputation
- Fever
- Gastrointestinal symptoms from anticoagulation/antiplatelet medication
- Hematoma at vascular access site
- Hemorrhage requiring transfusion
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection and pain at vascular access site

- Intimal tear
- Kidney infarct
- Myocardial infarction
- Myocardial ischemia
- Nephrectomy
- Peripheral neuropathy
- Pseudoaneurysm at vascular access site
- Pseudoaneurysm formation
- Renal artery thrombosis, aneurysm, rupture, perforation, occlusion, spasm, or restenosis
- Renal insufficiency or failure
- Stent migration or embolization
- Stent misplacement
- Stroke/cerebral vascular accident
- Tissue necrosis or ulceration

PRODUCT RECOMMENDATIONS

Wire Guide Use and Selection

The Formula 418 Vascular Balloon-Expandable Stent System is compatible with 0.018 inch wire guides.

Guiding Catheter/Introducer Selection

Correct guiding catheter/sheath selection and technique are necessary for use of the stent. Ensure that the inside lumen of the guiding catheter/sheath is of sufficient size to allow unobstructed passage of the Formula 418 Vascular Balloon-Expandable Stent System.

NOTE: Minimum inside diameter requirements for the guiding catheter/sheath are shown in Table 1.

Standard vascular catheter curves must be selected to provide adequate guiding catheter/sheath “back-up support” to achieve successful stent placement.

Stent Size Selection

The stent selected should have an expanded diameter approximately equal to, or slightly larger than, the estimated reference diameter of the artery to be stented. The length of the stent should be chosen to adequately cover the length of the lesion. A compliance card is provided with the product to facilitate the selection of an appropriate size of stent.

INSTRUCTIONS FOR USE

The Formula 418 Vascular Balloon-Expandable Stent System is used in conjunction with equipment required for a conventional PTA procedure including, but not limited to, a vascular access set, arterial sheath, guiding catheter, wire guide and inflation device.

Lesion/Vessel Preparation

Stent placement must take into account proximal atherosclerotic plaque, which may inhibit advancement of the stent, as well as atherosclerotic plaque beyond the lesion which may prevent advancement of the device across the primary lesion. In preparing the vessel, optimal vasodilation, anticoagulant and antiplatelet therapy are essential. This device is intended for treatment following suboptimal angioplasty.

Preparation of Balloon Catheter

1. Remove the Formula 418 Vascular Stent System from the package, remove protective sleeve from the distal tip of the catheter and inspect the stent to ensure it has not been damaged.
2. Prepare the balloon lumen with a standard 1:1 contrast-saline mixture as follows:

WARNING: Do not attempt pre-inflation technique to purge the balloon lumen. Do not use air or any gaseous medium to inflate the balloon.

- a. Using a 20 ml syringe containing 5 ml of contrast-saline mixture, apply negative pressure for 20-30 seconds.
- b. Release pressure, allowing negative pressure to draw mixture into balloon lumen.
- c. Detach syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
- d. Prepare inflation device in standard manner and purge to remove all air from syringe and tubing.
- e. Attach inflation device to balloon lumen directly, ensuring no air bubbles remain at connection.
- f. Pull negative pressure on inflation device.

CAUTION: Significant amounts of air in the balloon may cause uneven expansion of the stent and difficulty in deployment of the stent.

3. Moisten the stent with heparinized saline.

CAUTION: Do not use gauze sponges, as fibers may disrupt stent.

4. Flush the wire guide lumen of the balloon catheter in standard fashion to purge air.
5. Before inserting the Formula 418 Vascular Balloon-Expandable Stent System into the introducer or guiding catheter, be sure that a large-bore Tuohy-Borst "Y" adapter or introducer valve is on the guiding catheter/sheath.
6. Advance a 0.018 inch wire guide of appropriate length across the target lesion.
7. Insert the appropriate guiding catheter/sheath (see **Table 1**).
8. Advance the pre-mounted Formula 418 Vascular Balloon-Expandable Stent System over the wire into either the introducer valve or Tuohy-Borst "Y" adapter, making sure that the wire guide exits the wire guide port proximal to the balloon.
 - a. If using an introducer with a valve, make sure the flared end of the insertion tool (provided in the package) is loaded over the pre-mounted stent on the balloon catheter. Pass the insertion tool loaded with the premounted stent through the introducer valve. Push the Formula 418 Vascular Balloon-Expandable Stent System into the body of the introducer. Slide the insertion tool proximally up the catheter shaft away from the guiding catheter/sheath. A slight contact of the stent with the guiding catheter/sheath may be felt, but there must be no resistance.
 - b. If using a Tuohy-Borst "Y" adapter, advance the pre-mounted Formula 418 Vascular Balloon-Expandable Stent System over the wire and into the fully opened large-bore Tuohy-Borst "Y" adapter. Gently advance the Formula 418 Vascular Balloon-Expandable Stent System completely through the Tuohy-Borst "Y" adapter, taking care to maintain wire guide position. Then advance the Formula 418 Vascular Balloon-Expandable Stent System into the guiding catheter/sheath. A slight contact of the stent with the guiding catheter/sheath may be felt, but there must be no resistance.

WARNING: If resistance is encountered, do not force passage. Resistance may indicate damage to stent.

Positioning the Stent

1. Ensure guiding catheter/sheath stability before advancing the stent delivery balloon into the vascular artery.

WARNING: If initial guiding catheter/sheath position is lost, avoid pulling or pushing the guiding catheter/sheath over the stent. If this is done, the distal end of the guiding catheter/sheath may damage the stent.

2. Position the Formula Vascular Stent across the lesion, using both the distal and proximal balloon markers as reference points. Optimal placement requires the ends of the stent to extend beyond the margin of the lesion to be stented.

WARNING: If the Formula 418 Vascular Balloon-Expandable Stent System does not readily advance through the vessel, do not force. If the stent will not advance in spite of good guiding catheter/sheath support, consider dilating proximal obstructing plaque, or changing the wire guide or the guiding catheter/sheath. (Refer to instructions for Removal of Unexpanded Stent.)

NOTE: When stenting long lesions, be sure to cover the distal portion of the lesion first. It is very important not to leave distal lesions uncovered. Frequent injections of contrast around the uninflated Formula 418 Vascular Balloon-Expandable Stent System will allow visualization of the extent of the lesion and facilitate accurate placement of the stent.

Balloon Expansion/Stent Deployment

Prior to stent expansion, utilize high-resolution fluoroscopy to verify that the stent has not been dislodged during positioning.

1. To expand the stent, inflate the balloon with a 1:1 contrast-saline mixture to the recommended expansion pressure indicated on product label. Do not move stent system during deployment.
2. Complete expansion and apposition of the stent against the vessel wall is necessary for clinical success. Do not exceed rated burst pressure of the balloon as indicated on the product label.
3. Once the stent has been deployed, post-deployment inflation is at the discretion of the operator to achieve optimum angiographic appearance. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter.

Balloon Deflation and Removal

1. Completely deflate the balloon, by pulling negative pressure with the inflation device or a 20 ml syringe. This usually requires 10 seconds or less.

WARNING: Allow enough time for the balloon to fully deflate prior to removal.

2. Slowly withdraw the balloon catheter from the stent while maintaining negative pressure on the balloon. Maintain position of the guiding catheter/sheath to prevent it from being drawn into the vessel. Observe under fluoroscopy to ensure that the balloon disengages from the stent.

Removal of Unexpanded Stent

Do not attempt to pull an unexpanded stent back into the guiding catheter/sheath. The Formula 418 Vascular Balloon-Expandable Stent System should be withdrawn until the proximal end of the stent is aligned with the distal tip of the guiding catheter/sheath. Withdraw the guiding catheter/sheath and stent delivery system as a single unit, leaving the wire guide in place.

WARNING: If stent is removed, do not attempt to reuse the device. Damage to the stent may occur upon removal.

Stent-Assisted, High-Pressure Balloon Angioplasty

A deployed Formula 418 Vascular Balloon-Expandable Stent may be further expanded using a low-profile, high-pressure balloon catheter. The stent should have been well embedded into the vessel wall by the deployment balloon before the lesion is recrossed. If recrossing is necessary, the position of the stent should be monitored and care should be taken to avoid displacing the stent with subsequent manipulation of wire guides and balloons. An appropriate balloon catheter should be used with a very flexible wire guide.

If a Formula 418 Vascular Balloon-Expandable Stent is dilated with a high-pressure balloon, the balloon should be appropriately sized for the vessel diameter. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter. Use of a new balloon for post-deployment dilatation is suggested. Avoid balloon catheters with stiff tips that might displace the stent during advancement. Avoid balloons that "wing out" and/or re-wrap inconsistently, which might displace the stent during withdrawal of the balloon from the stent.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

FORMULA 535® VASCULAR BALLOON-EXPANDABLE STENT SYSTEM

DEVICE DESCRIPTION

The Formula 535 Vascular Stent is a balloon-expandable stent made of 316L stainless steel with a slotted tube configuration. It is pre-mounted on an over-the-wire balloon catheter, which serves as the delivery system (Formula 535 Vascular Balloon-Expandable Stent System). The stent is positioned between two radiopaque marker bands, which are located inside the balloon at the proximal and distal tapers of the balloon. The cannula design of the stent provides a low outside diameter profile. Refer to **Table 1** for guide catheter and sheath compatibility.

Table 1: Device Specifications (diameters and lengths)

Order Number	Catheter Length (cm)	Stent Diameter (mm)	Stent Length (mm)	Minimum Catheter/Sheath* Inside Diameter (inches)	Guiding Catheter/Sheath (French size)
FOV535-35-80-4-12	80	4.0	12	0.074	6.0/5.0
FOV535-35-80-4-16	80	4.0	16	0.074	7.0/5.0
FOV535-35-80-4-20	80	4.0	20	0.074	7.0/5.0
FOV535-35-80-4-24	80	4.0	24	0.074	7.0/5.0
FOV535-35-80-4-30	80	4.0	30	0.074	7.0/5.0
FOV535-35-80-4-40	80	4.0	40	0.074	7.0/5.0
FOV535-35-80-4-60	80	4.0	60	0.074	7.0/5.0
FOV535-35-80-5-12	80	5.0	12	0.087	7.0/6.0
FOV535-35-80-5-16	80	5.0	16	0.087	7.0/6.0
FOV535-35-80-5-20	80	5.0	20	0.087	8.0/6.0
FOV535-35-80-5-24	80	5.0	24	0.087	8.0/6.0
FOV535-35-80-5-30	80	5.0	30	0.087	8.0/6.0
FOV535-35-80-5-40	80	5.0	40	0.087	8.0/6.0
FOV535-35-80-5-60	80	5.0	60	0.087	8.0/6.0
FOV535-35-80-6-12	80	6.0	12	0.087	8.0/6.0
FOV535-35-80-6-16	80	6.0	16	0.087	8.0/6.0
FOV535-35-80-6-20	80	6.0	20	0.087	8.0/6.0
FOV535-35-80-6-24	80	6.0	24	0.087	8.0/6.0
FOV535-35-80-6-30	80	6.0	30	0.087	8.0/6.0
FOV535-35-80-6-40	80	6.0	40	0.087	8.0/6.0
FOV535-35-80-6-60	80	6.0	60	0.087	8.0/6.0
FOV535-35-80-7-12	80	7.0	12	0.087	8.0/6.0
FOV535-35-80-7-16	80	7.0	16	0.087	8.0/6.0
FOV535-35-80-7-20	80	7.0	20	0.087	8.0/6.0
FOV535-35-80-7-24	80	7.0	24	0.087	8.0/6.0
FOV535-35-80-7-30	80	7.0	30	0.087	8.0/6.0
FOV535-35-80-7-40	80	7.0	40	0.087	8.0/6.0
FOV535-35-80-7-60	80	7.0	60	0.087	8.0/6.0
FOV535-35-80-8-12	80	8.0	12	0.087	8.0/6.0
FOV535-35-80-8-16	80	8.0	16	0.087	8.0/6.0
FOV535-35-80-8-20	80	8.0	20	0.087	8.0/6.0
FOV535-35-80-8-24	80	8.0	24	0.087	8.0/6.0
FOV535-35-80-8-30	80	8.0	30	0.087	8.0/6.0
FOV535-35-80-8-40	80	8.0	40	0.087	8.0/6.0
FOV535-35-80-8-60	80	8.0	60	0.087	8.0/6.0
FOV535-35-80-9-20	80	9.0	20	0.100	9.0/7.0
FOV535-35-80-9-30	80	9.0	30	0.100	9.0/7.0
FOV535-35-80-9-40	80	9.0	40	0.100	9.0/7.0
FOV535-35-80-9-60	80	9.0	60	0.100	9.0/7.0
FOV535-35-80-10-20	80	10.0	20	0.100	9.0/7.0
FOV535-35-80-10-30	80	10.0	30	0.100	9.0/7.0
FOV535-35-80-10-40	80	10.0	40	0.100	9.0/7.0

FOV535-35-80-10-60	80	10.0	60	0.100	9.0/7.0
FOV535-35-135-4-12	135	4.0	12	0.074	6.0/5.0
FOV535-35-135-4-16	135	4.0	16	0.074	7.0/5.0
FOV535-35-135-4-20	135	4.0	20	0.074	7.0/5.0
FOV535-35-135-4-24	135	4.0	24	0.074	7.0/5.0
FOV535-35-135-4-30	135	4.0	30	0.074	7.0/5.0
FOV535-35-135-4-40	135	4.0	40	0.074	7.0/5.0
FOV535-35-135-4-60	135	4.0	60	0.074	7.0/5.0
FOV535-35-135-5-12	135	5.0	12	0.087	7.0/6.0
FOV535-35-135-5-16	135	5.0	16	0.087	7.0/6.0
FOV535-35-135-5-20	135	5.0	20	0.087	8.0/6.0
FOV535-35-135-5-24	135	5.0	24	0.087	8.0/6.0
FOV535-35-135-5-30	135	5.0	30	0.087	8.0/6.0
FOV535-35-135-5-40	135	5.0	40	0.087	8.0/6.0
FOV535-35-135-5-60	135	5.0	60	0.087	8.0/6.0
FOV535-35-135-6-12	135	6.0	12	0.087	8.0/6.0
FOV535-35-135-6-16	135	6.0	16	0.087	8.0/6.0
FOV535-35-135-6-20	135	6.0	20	0.087	8.0/6.0
FOV535-35-135-6-24	135	6.0	24	0.087	8.0/6.0
FOV535-35-135-6-30	135	6.0	30	0.087	8.0/6.0
FOV535-35-135-6-40	135	6.0	40	0.087	8.0/6.0
FOV535-35-135-6-60	135	6.0	60	0.087	8.0/6.0
FOV535-35-135-7-12	135	7.0	12	0.087	8.0/6.0
FOV535-35-135-7-16	135	7.0	16	0.087	8.0/6.0
FOV535-35-135-7-20	135	7.0	20	0.087	8.0/6.0
FOV535-35-135-7-24	135	7.0	24	0.087	8.0/6.0
FOV535-35-135-7-30	135	7.0	30	0.087	8.0/6.0
FOV535-35-135-7-40	135	7.0	40	0.087	8.0/6.0
FOV535-35-135-7-60	135	7.0	60	0.087	8.0/6.0
FOV535-35-135-8-12	135	8.0	12	0.087	8.0/6.0
FOV535-35-135-8-16	135	8.0	16	0.087	8.0/6.0
FOV535-35-135-8-20	135	8.0	20	0.087	8.0/6.0
FOV535-35-135-8-24	135	8.0	24	0.087	8.0/6.0
FOV535-35-135-8-30	135	8.0	30	0.087	8.0/6.0
FOV535-35-135-8-40	135	8.0	40	0.087	8.0/6.0
FOV535-35-135-8-60	135	8.0	60	0.087	8.0/6.0
FOV535-35-135-9-20	135	9.0	20	0.100	9.0/7.0
FOV535-35-135-9-30	135	9.0	30	0.100	9.0/7.0
FOV535-35-135-9-40	135	9.0	40	0.100	9.0/7.0
FOV535-35-135-9-60	135	9.0	60	0.100	9.0/7.0
FOV535-35-135-10-20	135	10.0	20	0.100	9.0/7.0
FOV535-35-135-10-30	135	10.0	30	0.100	9.0/7.0
FOV535-35-135-10-40	135	10.0	40	0.100	9.0/7.0
FOV535-35-135-10-60	135	10.0	60	0.100	9.0/7.0

*See manufacturer specifications for French (Fr) equivalent.

INTENDED USE

The Formula 535 Vascular Balloon-Expandable Stent System is intended for stenting of peripheral arteries.

CONTRAINDICATIONS

The Formula 535 Vascular Balloon-Expandable Stent System is contraindicated for use in:

- Patients for whom antiplatelet and/or anticoagulant therapy is contraindicated.
- Patients who have a lesion that cannot be crossed with a wire or a balloon angioplasty catheter.
- Patients who have stenoses that cannot be dilated to permit passage of the stent.
- Patients with bleeding disorders.
- Stenting of an arterial vessel where leakage from the artery could be exacerbated by placement of a stent.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.

WARNINGS

- The use of this device carries the associated risks of subacute thrombosis, vascular complications and/or bleeding events. Judicious selection of patients is necessary.

- Persons allergic to 316L stainless steel may suffer an allergic reaction to this implant.
- The stent should not be placed in a location with significant external compression which could crush the stent.
- The stent should not be deployed in areas subject to significant flexure or bending to reduce the risk of stent deformations or fracture.

PRECAUTIONS

General

- Implantation of the stent should be performed only by physicians who have received adequate training. Adequate training consists of specific training in vascular stenting techniques and procedures under the supervision of an interventional physician experienced in vascular stent implantation.
- Stent placement should only be performed at hospitals where emergency artery bypass graft surgery can be readily performed.
- Subsequent restenosis may require redilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of an endothelialized Formula 535 Vascular Stent is unknown at present.

Stent Placement

- Manipulation of the Formula 535 Vascular Balloon-Expandable Stent System requires fluoroscopic control.
- The stent is preloaded on its delivery balloon and should not be removed. Do not attempt to place stent on another balloon catheter for deployment.
- Deploy the stent over a wire guide not exceeding 0.035 inch diameter.
- DO NOT wipe or clean stent or catheter with organic solvent (i.e. isopropyl alcohol).
- DO NOT expose catheter to temperatures above 54°C (130°F).
- DO NOT use power injection systems with the delivery system.
- DO NOT rotate any part of the system during deployment.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is particularly important during removal of the catheter from packaging, placement over the wire guide and advancement through the large-bore Tuohy-Borst “Y” adapter and guiding catheter hub.
- Do not pre-inflate balloon prior to stent deployment (See Preparation of Balloon Catheter section in these Instructions for Use).
- Use only the appropriate balloon inflation medium. The standard inflation medium is a 1:1 mixture of contrast medium and normal saline. Do not use air or any gaseous substance as a balloon inflation medium.
- Expansion of the stent should not be undertaken if the stent is not properly positioned in the vessel (See Removal of Unexpanded Stent section in these Instructions for Use).
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the product label. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible intimal damage and dissection.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stent and may cause acute closure of the vessel, requiring additional intervention (further dilatation or placement of additional stents).
- Use prior to the expiration date on the package label.
- Intended for one-time use (Do not re-sterilize).
- The use of mechanical atherectomy devices within an implanted stent may cause damage to the stent.
- Repositioning of the device after deployment is not possible.

MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Formula 535 Vascular Stent is **MR Conditional** according to ASTM F2503, *Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment*. A patient with this stent can be scanned safely anytime after placement under the following conditions.

- Static magnetic field of 3-Tesla or 1.5-Tesla
- Maximum spatial magnetic gradient of 720 Gauss/cm or less
- Product of the spatial gradient and the static magnetic field of 21.6 T2/m or less
- MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg or less for 15 minutes of scanning

Non-clinical evaluation was conducted in an MR system (Excite, General Electric Healthcare) with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gaussmeter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual). In non-clinical testing, a single and two overlapped Formula 10 mm x 60 mm stents produced a maximum temperature rise of 2.5 °C and 3.7 °C, respectively, during 15 minutes of MRI (i.e., for one scanning sequence) performed in a 3 Tesla MR system (3 Tesla/128 MHz, Excite, Software G3.0-052B, General Electric Healthcare) at an MR system reported whole-body-averaged SAR of 3.0 W/kg (associated with a calorimetry measured whole-body-averaged value of 2.8 W/kg).

The effect of heating in the MRI environment for stents with fractured struts is unknown.

Image Artifact

MR image quality may be compromised if the area of interest is within the lumen or within approximately 15 mm of the Formula 10 mm x 60 mm stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

MedicAlert Foundation

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted at:

Mail: MedicAlert Foundation International

2323 Colorado Avenue

Turlock, CA 95382

Phone: 888-633-4298 (toll free)

209-668-3333 (from outside the US)

Fax: 209-669-2450

Web: www.medicalert.org

POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of a vascular stent are presented in the following alphabetical list and include, but are not limited to:

- Abscess
- Allergic reaction to stainless steel or contrast agents
- Arrhythmias (ventricular fibrillation, ventricular tachycardia, other)
- Arteriovenous fistula
- Bowel infarct
- Death
- Dialysis
- Dissection
- Drug reaction to antiplatelet agents
- Drug reaction, allergic reaction to contrast media

- Emboli (air, tissue, or thrombotic emboli) resulting in tissue ischemia/infarction
- Emergency surgery to correct vascular complications
- Emergent renal artery bypass surgery
- Extremity ischemia/amputation
- Fever
- Gastrointestinal symptoms from anticoagulation/antiplatelet medication
- Hematoma at vascular access site
- Hemorrhage requiring transfusion
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection and pain at vascular access site
- Intimal tear
- Kidney infarct
- Myocardial infarction
- Myocardial ischemia
- Nephrectomy
- Peripheral neuropathy
- Pseudoaneurysm at vascular access site
- Pseudoaneurysm formation
- Renal artery thrombosis, aneurysm, rupture, perforation, occlusion, spasm, or restenosis
- Renal insufficiency or failure
- Stent migration or embolization
- Stent misplacement
- Stroke/cerebral vascular accident
- Tissue necrosis or ulceration

PRODUCT RECOMMENDATIONS

Wire Guide Use and Selection

The Formula 535 Vascular Balloon-Expandable Stent System is compatible with 0.035 inch wire guides.

Guiding Catheter/Introducer Selection

Correct guiding catheter/sheath selection and technique are necessary for use of the stent. Ensure that the inside lumen of the guiding catheter/sheath is of sufficient size to allow unobstructed passage of the Formula 535 Vascular Balloon-Expandable Stent System.

NOTE: Minimum inside diameter requirements for the guiding catheter/sheath are shown in Table 1.

Standard vascular catheter curves must be selected to provide adequate guiding catheter/sheath “back-up support” to achieve successful stent placement.

Stent Size Selection

The stent selected should have an expanded diameter approximately equal to, or slightly larger than, the estimated reference diameter of the artery to be stented. The length of the stent should be chosen to adequately cover the length of the lesion. A compliance card is provided with the product to facilitate the selection of an appropriate size of stent.

INSTRUCTIONS FOR USE

The Formula 535 Vascular Balloon-Expandable Stent System is used in conjunction with equipment required for a conventional PTA procedure including, but not limited to, a vascular access set, arterial sheath, guiding catheter, wire guide and inflation device.

Lesion/Vessel Preparation

Stent placement must take into account proximal atherosclerotic plaque, which may inhibit advancement of the stent, as well as atherosclerotic plaque beyond the lesion which may prevent advancement of the device across the primary lesion. In preparing the vessel, optimal renal vasodilation, anticoagulant and antiplatelet therapy are essential. This device is intended for treatment following suboptimal angioplasty.

Preparation of Balloon Catheter

1. Remove the Formula 535 Vascular Stent System from the package, remove protective sleeve from the distal tip of the catheter and inspect the stent to ensure it has not been damaged.
2. Prepare the balloon lumen with a standard 1:1 contrast-saline mixture as follows:

WARNING: Do not attempt pre-inflation technique to purge the balloon lumen. Do not use air or any gaseous medium to inflate the balloon.

- a. Using a 20 ml syringe containing 5 ml of contrast-saline mixture, apply negative pressure for 20-30 seconds.
- b. Release pressure, allowing negative pressure to draw mixture into balloon lumen.
- c. Detach syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
- d. Prepare inflation device in standard manner and purge to remove all air from syringe and tubing.
- e. Attach inflation device to balloon lumen directly, ensuring no air bubbles remain at connection.
- f. Pull negative pressure on inflation device.

CAUTION: Significant amounts of air in the balloon may cause uneven expansion of the stent and difficulty in deployment of the stent.

3. Moisten the stent with heparinized saline.

CAUTION: Do not use gauze sponges, as fibers may disrupt stent.

4. Flush the wire guide lumen of the balloon catheter in standard fashion to purge air.
5. Before inserting the Formula 535 Vascular Balloon-Expandable Stent System into the introducer or guiding catheter, be sure that a large-bore Tuohy-Borst "Y" adapter or introducer valve is on the guiding catheter/sheath.
6. Advance a 0.035 inch wire guide of appropriate length across the target lesion.
7. Insert the appropriate guiding catheter/sheath (see **Table 1**).
8. Advance the pre-mounted Formula 535 Vascular Balloon-Expandable Stent System over the wire into either the introducer valve or Tuohy-Borst "Y" adapter, making sure that the wire guide exits the wire guide port proximal to the balloon.
 - a. If using an introducer with a valve, make sure the flared end of the insertion tool (provided in the package) is loaded over the pre-mounted stent on the balloon catheter. Pass the insertion tool loaded with the premounted stent through the introducer valve. Push the Formula 535 Vascular Balloon-Expandable Stent System into the body of the introducer. Slide the insertion tool proximally up the catheter shaft away from the guiding catheter/sheath. A slight contact of the stent with the guiding catheter/sheath may be felt, but there must be no resistance.
 - b. If using a Tuohy-Borst "Y" adapter, advance the pre-mounted Formula 535 Vascular Balloon-Expandable Stent System over the wire and into the fully opened large-bore Tuohy-Borst "Y" adapter. Gently advance the Formula 535 Vascular Balloon-Expandable Stent System completely through the Tuohy-Borst "Y" adapter, taking care to maintain wire guide position. Then advance the Formula 535 Vascular Balloon-Expandable Stent System into the guiding catheter/sheath. A slight contact of the stent with the guiding catheter/sheath may be felt, but there must be no resistance.

WARNING: If resistance is encountered, do not force passage. Resistance may indicate damage to stent.

Positioning the Stent

1. Ensure guiding catheter/sheath stability before advancing the stent delivery balloon into the vascular artery.

WARNING: If initial guiding catheter/sheath position is lost, avoid pulling or pushing the guiding catheter/sheath over the stent. If this is done, the distal end of the guiding catheter/sheath may damage the stent.

2. Position the Formula Vascular Stent across the lesion, using both the distal and proximal balloon markers as reference points. Optimal placement requires the ends of the stent to extend beyond the margin of the lesion to be stented.

WARNING: If the Formula 535 Vascular Balloon-Expandable Stent System does not readily advance through the vessel, do not force. If the stent will not advance in spite of good guiding catheter/sheath support, consider dilating proximal obstructing plaque, or changing the wire guide or the guiding catheter/sheath. (Refer to instructions for Removal of Unexpanded Stent.)

NOTE: When stenting long lesions, be sure to cover the distal portion of the lesion first. It is very important not to leave distal lesions uncovered. Frequent injections of contrast around the uninflated Formula 535 Vascular Balloon-Expandable Stent System will allow visualization of the extent of the lesion and facilitate accurate placement of the stent.

Balloon Expansion/Stent Deployment

Prior to stent expansion, utilize high-resolution fluoroscopy to verify that the stent has not been dislodged during positioning.

1. To expand the stent, inflate the balloon with a 1:1 contrast-saline mixture to the recommended expansion pressure indicated on product label. Do not move stent system during deployment.
2. Complete expansion and apposition of the stent against the vessel wall is necessary for clinical success. Do not exceed rated burst pressure of the balloon as indicated on the product label.
3. Once the stent has been deployed, post-deployment inflation is at the discretion of the operator to achieve optimum angiographic appearance. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter.

Balloon Deflation and Removal

1. Completely deflate the balloon, by pulling negative pressure with the inflation device or a 20 ml syringe.

WARNING: Allow enough time for the balloon to fully deflate prior to removal.

2. Slowly withdraw the balloon catheter from the stent while maintaining negative pressure on the balloon. Maintain position of the guiding catheter/sheath to prevent it from being drawn into the vessel. Observe under fluoroscopy to ensure that the balloon disengages from the stent.

Removal of Unexpanded Stent

Do not attempt to pull an unexpanded stent back into the guiding catheter/sheath. The Formula 535 Vascular Balloon-Expandable Stent System should be withdrawn until the proximal end of the stent is aligned with the distal tip of the guiding catheter/sheath. Withdraw the guiding catheter/sheath and stent delivery system as a single unit, leaving the wire guide in place.

WARNING: If stent is removed, do not attempt to reuse the device. Damage to the stent may occur upon removal.

Stent-Assisted, High-Pressure Balloon Angioplasty

A deployed Formula 535 Vascular Balloon-Expandable Stent may be further expanded using a low-profile, high-pressure balloon catheter. The stent should have been well embedded into the vessel wall by the deployment balloon before the lesion is recrossed. If recrossing is necessary, the position of the stent should be monitored and care should be taken to avoid displacing the

stent with subsequent manipulation of wire guides and balloons. An appropriate balloon catheter should be used with a very flexible wire guide.

If a Formula 535 Vascular Balloon-Expandable Stent is dilated with a high-pressure balloon, the balloon should be appropriately sized for the vessel diameter. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter. Use of a new balloon for post-deployment dilatation is suggested. Avoid balloon catheters with stiff tips that might displace the stent during advancement. Avoid balloons that “wing out” and/or re-wrap inconsistently, which might displace the stent during withdrawal of the balloon from the stent.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

FORMULA 414, 414 RX VASCULAR BALLOON-EXPANDABLE STENT SYSTEM

DEVICE DESCRIPTION

The Formula 414 Vascular Stent is a balloon-expandable stent made of 316L stainless steel with a slotted tube configuration. It is premounted on a rapid exchange balloon catheter, which serves as the delivery system (Formula 414 RX Vascular Balloon-Expandable Stent System). The stent is positioned between two radiopaque marker bands, which are located inside the balloon at the proximal and distal tapers of the balloon. The cannula design of the stent provides a low outside diameter profile, which permits use with a 5.0 French sheath or 6.0 French guiding catheter. Formula 414 RX Vascular Balloon-Expandable Stent System is pre-mounted on 80 cm and 135 cm rapid-exchange balloon catheter delivery systems. The stent is available in nominal expanded diameters of 4-7 mm with lengths of 12, 16, 20 and 24 mm each. (**Table 1**)

Table 1: Device Specifications (diameters and lengths)

Order Number	Catheter Length (cm)	Stent Diameter (mm)	Stent Length (mm)	Minimum Guiding Catheter/ Sheath* Inside Diameter (inches/mm)	Guiding Catheter/ Sheath (French Size)
FOVX414-14-80-4-12	80	4	12	0.068/1.73	6.0/5.0
FOVX414-14-80-4-16	80	4	16	0.068/1.73	6.0/5.0
FOVX414-14-80-4-20	80	4	20	0.068/1.73	6.0/5.0
FOVX414-14-80-4-24	80	4	24	0.068/1.73	6.0/5.0
FOVX414-14-80-5-12	80	5	12	0.068/1.73	6.0/5.0
FOVX414-14-80-5-16	80	5	16	0.068/1.73	6.0/5.0
FOVX414-14-80-5-20	80	5	20	0.068/1.73	6.0/5.0
FOVX414-14-80-5-24	80	5	24	0.068/1.73	6.0/5.0
FOVX414-14-80-6-12	80	6	12	0.068/1.73	6.0/5.0
FOVX414-14-80-6-16	80	6	16	0.068/1.73	6.0/5.0
FOVX414-14-80-6-20	80	6	20	0.068/1.73	6.0/5.0
FOVX414-14-80-6-24	80	6	24	0.068/1.73	6.0/5.0
FOVX414-14-80-7-12	80	7	12	0.068/1.73	6.0/5.0
FOVX414-14-80-7-16	80	7	16	0.068/1.73	6.0/5.0

FOVX414-14-80-7-20	80	7	20	0.068/1.73	6.0/5.0
FOVX414-14-80-7-24	80	7	24	0.068/1.73	6.0/5.0
FOVX414-14-135-4-12	135	4	12	0.068/1.73	6.0/5.0
FOVX414-14-135-4-16	135	4	16	0.068/1.73	6.0/5.0
FOVX414-14-135-4-20	135	4	20	0.068/1.73	6.0/5.0
FOVX414-14-135-4-24	135	4	24	0.068/1.73	6.0/5.0
FOVX414-14-135-5-12	135	5	12	0.068/1.73	6.0/5.0
FOVX414-14-135-5-16	135	5	16	0.068/1.73	6.0/5.0
FOVX414-14-135-5-20	135	5	20	0.068/1.73	6.0/5.0
FOVX414-14-135-5-24	135	5	24	0.068/1.73	6.0/5.0
FOVX414-14-135-6-12	135	6	12	0.068/1.73	6.0/5.0
FOVX414-14-135-6-16	135	6	16	0.068/1.73	6.0/5.0
FOVX414-14-135-6-20	135	6	20	0.068/1.73	6.0/5.0
FOVX414-14-135-6-24	135	6	24	0.068/1.73	6.0/5.0
FOVX414-14-135-7-12	135	7	12	0.068/1.73	6.0/5.0
FOVX414-14-135-7-16	135	7	16	0.068/1.73	6.0/5.0
FOVX414-14-135-7-20	135	7	20	0.068/1.73	6.0/5.0
FOVX414-14-135-7-24	135	7	24	0.068/1.73	6.0/5.0
FOVX414-14-80-4-20	80	4	20	0.068/1.73	6.0/5.0

*See manufacturer specifications for French (FR) equivalent

INTENDED USE

The Formula 414 RX Vascular Balloon-Expandable Stent System is intended for stenting of peripheral arteries.

CONTRAINDICATIONS

The Formula 414 RX Vascular Balloon-Expandable Stent System is contraindicated for use in:

- Patients for whom antiplatelet and/or anticoagulant therapy is contraindicated.
- Patients who have a lesion that cannot be crossed with a wire or a balloon angioplasty catheter.
- Patients who have stenoses that cannot be dilated to permit passage of the stent.
- Patients with bleeding disorders.
- Stenting of an arterial vessel where leakage from the artery could be exacerbated by placement of a stent.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.

WARNINGS

- The use of this device carries the associated risks of subacute thrombosis, vascular complications and/or bleeding events. Judicious selection of patients is necessary.
- Persons allergic to 316L stainless steel may suffer an allergic reaction to this implant.
- The stent should not be placed in a location with significant external compression which could crush the stent.
- The stent should not be deployed in areas subject to significant flexure or bending to reduce the risk of stent deformations or fracture.

PRECAUTIONS

General

- Implantation of the stent should be performed only by physicians who have received adequate training. Adequate training consists of specific training in vascular stenting techniques and procedures under the supervision of an interventional physician experienced in vascular stent implantation.
- Stent placement should only be performed at hospitals where emergency artery bypass graft surgery can be readily performed.
- Subsequent restenosis may require redilatation of the arterial segment containing the stent. The long-term outcome following repeat dilatation of an endothelialized Formula 414 Vascular Stent is unknown at present.

Stent Placement

- Manipulation of the Formula 414 RX Vascular Balloon-Expandable Stent System requires fluoroscopic control.

- The stent is preloaded on its delivery balloon and should not be removed. Do not attempt to place stent on another balloon catheter for deployment.
- Deploy the stent over a wire guide not exceeding 0.014 inch diameter.
- DO NOT wipe or clean stent or catheter with organic solvent (e.g. isopropyl alcohol).
- DO NOT expose catheter to temperatures above 54°C (130°F).
- DO NOT use power injection systems with the delivery system.
- DO NOT rotate any part of the system during deployment.
- Special care must be taken not to handle or in any way disrupt the stent on the balloon. This is particularly important during removal of the catheter from packaging, placement over the wire guide and advancement through the large-bore Tuohy-Borst "Y" adapter and guiding catheter hub.
- Do not pre-inflate balloon prior to stent deployment (See **Preparation of Balloon Catheter** section in these Instructions for Use).
- Use only the appropriate balloon inflation medium. The standard inflation medium is a 1:1 mixture of contrast medium and normal saline. Do not use air or any gaseous substance as a balloon inflation medium.
- Expansion of the stent should not be undertaken if the stent is not properly positioned in the vessel (See **Removal of Unexpanded Stent** section in these Instructions for Use).
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the product label. Use of pressures higher than specified on the product label may result in a ruptured balloon, with possible intimal damage and dissection.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stent and may cause acute closure of the vessel, requiring additional intervention (further dilatation or placement of additional stents).
- Use prior to the expiration date on the package label.
- Intended for one-time use (Do not re-sterilize).
- The use of mechanical atherectomy devices within an implanted stent may cause damage to the stent.
- Repositioning of the device after deployment is not possible.

MAGNETIC RESONANCE IMAGING

Non-clinical testing has demonstrated that the Formula 414 Vascular Stent is MR Conditional according to ASTM F2503, Standard Practice for Marking Medical Devices and Other items for Safety in the Magnetic Resonance Environment. A patient with this stent can be scanned safely anytime after placement under the following conditions:

- Static magnetic field of 3 Tesla or 1.5 Tesla.
- Maximum spatial gradient magnetic field of 720 Gauss/cm or less.
- Product of the spatial gradient and the static magnetic field of 21.6 T²/m or less.
- MR system reported whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg or less for 15 minutes of scanning.

POTENTIAL ADVERSE EVENTS

Adverse events that may be associated with the use of a vascular stent are presented in the following alphabetical list and include, but are not limited to:

- Allergic reaction to stainless steel or contrast agents
- Arrhythmias (VF, VT, other)
- Arteriovenous fistula
- Artery thrombosis, aneurysm, rupture, perforation, occlusion, spasm, or restenosis
- Bowel infarct
- Death
- Dialysis

- Dissection
- Emboli (air, tissue or thrombotic emboli), resulting in tissue ischemia/infarction
- Emergent artery bypass surgery to correct vascular complications
- Extremity ischemia/amputation
- Gastrointestinal syndroms from anticoagulation/antiplatelet medication
- Hematoma at vascular access site
- Hemorrhage requiring transfusion
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection and pain at vascular access site
- Intimal tear
- Myocardial infarction
- Myocardial ischemia
- Pseudoaneurysm at vascular access site
- Stent migration or embolization
- Stroke/CVA

PRODUCT RECOMMENDATIONS

Wire Guide Selection

The Formula 414 RX Vascular Balloon-Expandable Stent System is compatible with 0.014 inch wire guides.

Guiding Catheter/Sheath Selection

Correct guiding catheter/sheath selection and technique are necessary for use of the stent. Ensure that the inside lumen of the guiding catheter/sheath is of sufficient size to allow unobstructed passage of the Formula 414 RX Vascular Balloon-Expandable Stent System.

NOTE: Minimum inside diameter requirements for the guiding catheter/sheath are shown in Table 1.

Standard vascular catheter curves must be selected to provide adequate guiding catheter/sheath “back-up support” to achieve successful stent placement.

Stent Size Selection

The stent selected should have an expanded diameter approximately equal to, or slightly larger than, the estimated reference diameter of the artery to be stented. The length of the stent should be chosen to adequately cover the length of the lesion. A compliance card is provided with the product to facilitate the selection of an appropriate size of stent.

INSTRUCTIONS FOR USE

The Formula 414 RX Vascular Balloon-Expandable Stent System is used in conjunction with equipment required for a conventional PTA procedure including, but not limited to, a vascular access set, arterial sheath, guiding catheter, wire guide and inflation device.

Lesion/Vessel Preparation

Stent placement must take into account proximal atherosclerotic plaque, which may inhibit advancement of the stent, as well as atherosclerotic plaque beyond the lesion which may prevent advancement of the device across the primary lesion. In preparing the vessel, optimal vasodilation, anticoagulant and antiplatelet therapy are essential. This device is intended for treatment following suboptimal angioplasty.

Preparation of Balloon Catheter

1. Remove the Formula 414 RX Vascular Stent System from the package, remove protective sleeve from the distal tip of the catheter and inspect the stent to ensure it has not been damaged.
2. Prepare the balloon lumen with a standard 1:1 contrast-saline mixture as follows:

WARNING: Do not attempt pre-inflation technique to purge the balloon lumen. Do not use air or any gaseous medium to inflate the balloon.

- a. Using a 20 ml syringe containing 5 ml of contrast-saline mixture, apply negative pressure for 20-30 seconds.
- b. Release pressure, allowing negative pressure to draw mixture into balloon lumen.
- c. Detach syringe, leaving a meniscus of mixture on the hub of the balloon lumen.
- d. Prepare inflation device in standard manner and purge to remove all air from syringe and tubing.
- e. Attach inflation device to balloon lumen directly, ensuring no air bubbles remain at connection.
- f. Pull negative pressure on inflation device.

caution : Significant amounts of air in the balloon may cause uneven expansion of the stent and difficulty in deployment of the stent.

3. Moisten the stent with heparinized saline.

caution : Do not use gauze sponges, as fibers may disrupt stent.

4. Flush the wire guide lumen of the balloon catheter as follows:

- a. Remove the protective tube from the flushing sheath.
- b. Attach a syringe filled with heparinized normal saline to the flushing tool, insert the flushing tool into the distal end of the catheter and inject this saline into the lumen.
- c. Flush solution should be seen coming out of the wire guide port proximal to the balloon. Follow this procedure for any additional flushing.

5. Before inserting the Formula 414 RX Vascular Balloon-Expandable Stent System into the introducer or guiding catheter, be sure that

a large-bore Tuohy-Borst "Y" adapter or introducer valve is on the guiding catheter/sheath.

6. Advance a 0.014 inch wire guide of appropriate length across the target lesion.

7. Insert the appropriate guiding catheter/sheath (see **Table 1**).

8. Advance the premounted Formula 414 RX Vascular Balloon-Expandable Stent System over the wire into either the introducer valve or Tuohy-Borst "Y" adapter, making sure that the wire guide exits the wire guide port proximal to the balloon.

- a. If using an introducer with a valve, make sure the flared end of the insertion tool (provided in the package) is loaded over the premounted stent on the balloon catheter. Pass the insertion tool loaded with the premounted stent through the introducer valve. Push the Formula 414 RX Vascular Balloon-Expandable Stent System into the body of the introducer. Slide the insertion tool proximally up the catheter shaft away from the guiding catheter/sheath. A slight contact of the stent with the introducer may be felt, but there must be no resistance.

- b. If using a Tuohy-Borst "Y" adapter, advance the premounted Formula 414 RX Vascular Balloon-Expandable Stent System over the wire and into the fully opened large-bore Tuohy-Borst "Y" adapter. Gently advance the Formula 414 RX Vascular Balloon-Expandable Stent System completely through the Tuohy-Borst "Y" adapter, taking care to maintain wire guide position. Then advance the Formula 414 RX Vascular Balloon-Expandable Stent System into the guiding catheter/sheath. A slight contact of the stent with the introducer sheath/guiding catheter may be felt, but there must be no resistance.

WARNING: If resistance is encountered, do not force passage. Resistance may indicate damage to stent.

Positioning the Stent

1. Ensure guiding catheter/sheath stability before advancing the stent delivery balloon into the vascular artery.

WARNING: If initial guiding catheter/sheath position is lost, avoid pulling or pushing the guiding catheter/sheath over the stent. If this is done, the distal end of the guiding catheter/sheath may damage the stent.

2. Position the Formula Vascular Stent across the lesion, using both the distal and proximal balloon markers as reference points.

Optimal placement requires the ends of the stent to extend beyond the margin of the lesion to be stented.

WARNING: If the Formula 414 RX Vascular Balloon-Expandable Stent System does not readily advance through the vessel, do not force. If the stent will not advance in spite of good guiding catheter/sheath support, consider dilating proximal obstructing plaque, or changing the wire guide or the guiding catheter/sheath. (Refer to instructions for Removal of Unexpanded Stent.)

NOTE: When stenting long lesions, be sure to cover the distal portion of the dissection first. It is very important not to leave distal lesions uncovered. Frequent injections of contrast around the uninflated Formula 414 RX Vascular Balloon-Expandable Stent System will allow visualization of the extent of the dissection and facilitate accurate placement of the stent.

Balloon Expansion/Stent Deployment

Prior to stent expansion, utilize high-resolution fluoroscopy to verify that the stent has not been dislodged during positioning.

1. To expand the stent, inflate the balloon with a 1:1 contrast-saline mixture to the recommended expansion pressure indicated on product label. Do not move stent system during deployment.
2. Complete expansion and apposition of the stent against the vessel wall is necessary for clinical success. Do not exceed rated burst pressure of the balloon as indicated on the product label.
3. Once the stent has been deployed, post-deployment inflation is at the discretion of the operator to achieve optimum angiographic appearance. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter.

Balloon Deflation and Removal

1. Completely deflate the balloon, by pulling negative pressure with the inflation device or a 20 ml syringe.

WARNING: Allow enough time for the balloon to fully deflate prior to removal.

2. Slowly withdraw the balloon catheter from the stent while maintaining negative pressure on the balloon. Maintain position of the guiding catheter/sheath to prevent it from being drawn into the vessel. Observe under fluoroscopy to ensure that the balloon disengages from the stent.

Removal of Unexpanded Stent

Do not attempt to pull an unexpanded stent back into the introducer guiding catheter/sheath. The Formula 414 RX Vascular Balloon-Expandable Stent System should be withdrawn until the proximal end of the stent is aligned with the distal tip of the guiding catheter/sheath. Withdraw the guiding catheter/sheath and stent delivery system as a single unit, leaving the wire guide in place.

WARNING: If stent is removed, do not attempt to reuse the device. Damage to the stent may occur upon removal.

Stent-Assisted, High-Pressure Balloon Angioplasty

A deployed Formula 414 RX Vascular Balloon-Expandable Stent may be further expanded using a low-profile, high-pressure balloon catheter. The stent should have been well embedded into the vessel wall by the deployment balloon before the lesion is recrossed.

If recrossing is necessary, the position of the stent should be monitored and care should be taken to avoid displacing the stent with subsequent manipulation of wire guides and balloons. An appropriate balloon catheter should be used with a very flexible wire guide.

If a Formula 414 RX Vascular Balloon-Expandable Stent is dilated with a high-pressure balloon, the balloon should be appropriately sized for the vessel diameter. The stent may be post-dilated up to 1 mm beyond nominal labeled stent diameter. Use of a new balloon for post-deployment

dilatation is suggested. Avoid balloon catheters with stiff tips that might displace the stent during advancement.

Avoid balloons that “wing out” and/or re-wrap inconsistently, which might displace the stent during withdrawal of the balloon from the stent.

HOW SUPPLIED

Supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

ZILVER VENA VENOUS SELF-EXPANDING STENT

U.S CAUTION: INVESTIGATIONAL DEVICE Limited by federal (or United States) law to investigational use.

Do not re-sterilize.

DEVICE DESCRIPTION

The Zilver® Vena™ Venous Stent is a self-expanding, flexible, slotted-tube nitinol stent. Post-deployment, the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region.

Upon deployment, the stent provides support, while maintaining flexibility in the vessel.

The Zilver® Vena™ Venous Stent is available in the following sizes:

7 French (2.3 mm) Delivery System							
Stent Length (mm)		60		100		140	
Delivery System (cm)		80	120	80	120	80	120
Stent Inner Diameter (mm)	14	x	x	x	x	x	x
	16	x	x	x	x	x	x

The stent comes preloaded in a 7Fr (2.3 mm) delivery catheter. Handloading of the stent is not possible. Stent deployment is controlled by retraction of the handle while holding the metal cannula stationary.

INTENDED USE

The Zilver® Vena™ Venous Stent is intended for use in the iliofemoral veins for the treatment of symptomatic venous outflow obstruction.

CONTRAINDICATIONS

Patients with total venous occlusion that cannot be dilated to allow passage of the introducer system or wire guide.

WARNINGS

Nitinol (nickel titanium) may cause allergic reactions in some patients.

PRECAUTIONS

- This product should only be used by physicians trained and experienced in diagnostic and interventional vascular techniques. Standard techniques for interventional vascular procedures should be employed.
- Manipulation of the Zilver® Vena™ Venous Stent requires high-resolution fluoroscopic control.
- Do not use power injection systems with the delivery system.

- Prior to the procedure, the patient's underlying condition should be assessed for compatibility with anticipated procedural and postprocedural antiplatelet/anticoagulation therapy.
- Use in patients with a history of contrast sensitivity is not recommended unless the patient can be adequately premedicated.
- Arterial use is not recommended.

Stent Handling

- Do not attempt to remove the stent from the delivery system before use.
- Do not expose any part of the delivery system to organic solvents (e.g., alcohol).
- Carefully inspect the sterile package and stent system prior to use to verify that neither has been damaged during shipment.
- Use the stent system prior to the expiration date specified on the package.
- The device is intended for single use. Attempts to reprocess, re-sterilize, and/or reuse may lead to device failure and/or transmission of disease.

Stent Placement

- Ensure that the safety lock is not inadvertently removed prior to stent release.
- Do not rotate any part of the system during deployment.
- Repositioning of the device after deployment is not possible since the introducer catheter cannot be re-advanced over the stent once deployment begins.
- Once stent deployment has begun, the stent must be fully deployed.

Stent/System Removal

- Do not advance sheath after stent has been deployed. Delivery system can be removed without the need to recapture tip.

Post Implant

- Antiplatelet/anticoagulant therapy should be administered during and after procedure according to institutional standard of care.
- Use caution when re-crossing a stent to avoid stent damage or migration.

MRI INFORMATION

Non-clinical testing has demonstrated that the Zilver® Vena™ Venous stent is MR Conditional according to ASTM F2503. A patient with this stent may be scanned safely under the following conditions.

Static Magnetic Field

- Static magnetic field of 3.0 Tesla or less
- Spatial magnetic gradient field of 1580 Gauss/cm or less

The static magnetic field for comparison to the above limits is the static magnetic field that is pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual).

MRI-Related Heating

- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 3.0 W/kg for 15 minutes of scanning (i.e., per scanning sequence);
- 1.5 Tesla Temperature Rise

In non-clinical testing, the Zilver® Vena™ Venous stent produced maximum temperature rises of 2.6 and 3.8 °C (for a single stent and pair of overlapping stents) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 1.5 Tesla System (Siemens

Magnetom) at an MR system reported whole-body-averaged SAR of 3.8 W/kg (associated with a calorimetry measured whole body averaged value of 3.1 W/kg).

• **3.0 Tesla Temperature Rise**

In non-clinical testing, the Zilver® Vena™ Venous stent produced a maximum temperature rises of 2.2 and 4.6 °C (for a single stent and pair of overlapping stents) during 15 minutes of MR imaging (i.e., for one scanning sequence) performed in a MR 3 Tesla System (General Electric Excite) at an MR system reported whole-body-averaged SAR of 3.0 W/kg (associated with a calorimetry measured whole body averaged value of 2.8 W/kg).

Image Artifact

MR image quality may be compromised if the area of interest is within approximately 7 mm of the position of Zilver® Vena™ Venous stent as found during non-clinical testing using T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic stent.

For US Patients Only

Cook recommends that the patient register the MR conditions disclosed in this IFU with the MedicAlert Foundation. The MedicAlert Foundation can be contacted in the following manners:

Mail: MedicAlert Foundation International

2323 Colorado Avenue

Turlock, CA 95382

Phone: 888-633-4298 (toll free)

209-668-3333 from outside the US

Fax: 209-669-2450

Web: www.medicalert.org

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to, the following:

- Abrupt stent closure
- Allergic reaction to anticoagulant and/or antithrombotic therapy or contrast medium
- Allergic reaction to Nitinol
- Arrhythmia
- Arteriovenous fistula
- Death
- Embolism
- Fever
- Hematoma/hemorrhage at access site
- Hypersensitivity reactions
- Hypertension
- Hypotension, nausea or symptoms of a vasovagal response
- Infection/abscess formation at access site
- Intimal injury/dissection
- Myocardial infarction (MI)
- Pseudoaneurysm formation
- Pulmonary embolism
- Renal failure
- Restenosis or thrombosis of the stented vein
- Retinal, retroperitoneal, gastrointestinal, or genitourinary bleeding associated with anticoagulant therapy
- Septicemia/bacteremia
- Spasm

- Stent malapposition
- Stent migration
- Stent strut fracture
- Stroke
- Tissue necrosis
- Vessel perforation/rupture
- Worsened pain

PRODUCT RECOMMENDATIONS

Venous Access

For venous access, the use of an access set that accepts a 7.0 French (2.3 mm) introducer catheter is recommended.

Wire Guide Selection

The use of a .035-inch (0.89 mm) wire guide is recommended.

PTA Balloon Selection

For pre- and post-dilation, an appropriately sized balloon catheter is recommended.

Stent Selection

The chosen stent diameter should be at least 1 mm larger than the estimated diameter of the target vein.

HOW SUPPLIED

The Zilver® Vena™ Venous Stent is supplied sterilized by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from package, inspect the product to ensure no damage has occurred.

Zilver® Vena™ Venous Stent (Fig. 1)


Fig. 1

- a. Handle
- b. Hub
- c. Safety Lock
- d. Delivery System: Outer Sheath
- e. Tip of Delivery System Inner Catheter
- f. Side-arm Flushing Port

- g. Metal Cannula
- h. Radiopaque Marker on the Delivery System
- i. Gold Radiopaque Markers

The stent is supplied preloaded in a 7.0 French (2.3 mm) introducer catheter. A radiopaque marker (h) on the distal tip of the sheath is used to visualize deployment of the stent. Stent deployment is controlled by means of a hand-held device.

Multiple Stent Placement

If placement of more than one stent is required, the following recommendations should be considered:

- When more than one stent is required, resulting in stent-to-stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.
- In relation to the operator, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximally to the previously placed stent).
- Stents placed in tandem must overlap to allow for complete coverage of the lesion.

INSTRUCTIONS FOR USE

Stent Sizing

1. Determine the proper stent size after complete diagnostic evaluation. The stent diameter should be at least 1 mm larger than the estimated diameter of the target vein. Measure the total length of the obstructed vein to determine the length of the segment that requires the stent and the length of the stent needed. Stent from the most distal end of the lesion to the most proximal end of the lesion in a continuous fashion.

NOTE: If multiple stents are required to cover the length of the lesion, please refer to the Multiple Stent Placement section of these instructions for use for further recommendations.

Introduction of the Stent

1. Gain access at the appropriate site utilizing an introducer sheath or guiding catheter. Refer to product labeling for required introducer sheath or guiding catheter size.
2. Introduce an .035 inch (0.89 mm) wire guide through the introducer sheath or guiding catheter across the distal segment of the target lesion.
3. Predilatation before stent placement is at the discretion of the physician. Remove the balloon catheter, leaving the wire guide in place.


Fig. 2

4. Immediately before placing the delivery system into the body, use a 1 ml syringe to flush the delivery system with saline through the side-arm flushing port. Flush until a few drops of saline exit the distal tip, between the delivery system outer sheath (d) and inner catheter (e). (**Fig. 2**)
5. Use the 1 ml syringe to flush the wire guide lumen of the stent delivery system with saline through the hub. (**Fig. 3**)

Fig. 3

6. Insert the delivery system over the wire guide.
7. Under fluoroscopy, advance the delivery system past the most distal segment of the lesion.


Fig. 4

8. Pull back on the stent delivery system under fluoroscopy until the radiopaque markers on the stent (i) are at the desired position. The stent is now ready to be deployed. (**Fig. 4**)

Deployment of the Stent

1. Before deployment, it is important to straighten the proximal part of the delivery system as much as possible, remove any slack in the delivery system and to keep the handle in a stable position.
2. The stent expansion must be performed under fluoroscopic control.


Fig. 5

3. Hold the hub (b) on the metal cannula (g) steady. To deploy the stent, remove the red safety lock (c). (**Fig. 5**)


Fig. 6

4. Hold the hub end stationary. The stent will be deployed as you pull the handle (a) toward the hub (b). (Fig. 6) The radiopaque marker on the delivery system (h) indicates the progress of the deployment.

NOTE: Full deployment of the stent length will occur when the distal end of the sheath has been retracted past the proximal part of the stent.


Fig. 7

5. As deployment occurs, continue sliding the handle (a) toward the hub (b) in a slow, smooth and consistent fashion. (Fig. 7)

NOTE: Once stent deployment has begun, the stent must be fully deployed. Repositioning of the SilverR Vena Venous Stent is not possible since the delivery system's outer sheath cannot be re-advanced over the stent once deployment begins. Refer to the Multiple Stent Placement section of these instructions for use for information on missed lesions.


Fig. 8

6. The stent is fully deployed when the handle (a) reaches the hub (b). (Fig. 8)

Post Stent Deployment

1. Delivery System Removal - Do not advance sheath after the stent has been deployed. Delivery system can be removed without the need to recapture tip. Check for delivery system integrity after removal from the patient.

2. Perform a venogram to ensure full deployment of the device.

NOTE: If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation can be performed at the discretion of the physician. Select an appropriate size balloon catheter and dilate the lesion with conventional technique. The inflation diameter of the balloon used for post dilatation should approximate the diameter of the reference vessel. Remove the balloon from the patient.

3. Introducer sheath and wire guide may be removed at this point.

4. Close the entry wound as appropriate.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

ZILVER FLEX 35 VASCULAR STENT

CAUTION: U.S. federal law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Do not re-sterilize.

DEVICE DESCRIPTION

The ZILVER FLEX™ 35 Vascular Stent (ZFV) is a self-expandable stent made of nitinol. The special laser-cut pattern of the nitinol tube provides a construction with a strong radial force and high flexibility.


The stent has radiopaque markers at both ends (i).

The stent is supplied preloaded in a 6.0 French (2.0mm) introducer catheter. A radiopaque marker (h) on the distal tip of the sheath is used to visualize deployment of the stent. Stent deployment is controlled by means of a hand-held device. Refer to Stent Size Selection Table below.

Stent length (mm)	20		30		40		60		80	
Delivery system (cm)	80	125	80	125	80	125	80	125	80	125
Stent outer diameter (mm)	5	x	x	x	x	x	x	x	x	x
	6	x	x	x	x	x	x	x	x	x
	7	x	x	x	x	x	x	x	x	x
	8	x	x	x	x	x	x	x	x	x
	9	x	x	x	x	x	x	x	x	x
	10	x	x	x	x	x	x	x	x	x

Stent length (mm)	100		120		140		170		200	
Delivery system (cm)	80	125	80	125	80	125	80	125	80	125
Stent outer diameter (mm)	5	x	x	x	x	x	x	x	x	x
	6	x	x	x	x	x	x	x	x	x
	7	x	x	x	x	x	x	x	x	x

	8	x	x	x	x	x	x	x	x	x	x
	9	x	x	x	x	x	x	x			
	10	x	x	x	x	x	x	x			

- a. Handle
- b. Hub
- c. Safety Lock
- d. Delivery System: Outer Sheath
- e. Tip of Delivery System Inner Catheter
- f. Side-arm Flushing Port
- g. Metal Cannula
- h. Radiopaque Marker on the Delivery System
- i. Gold Radiopaque Markers

INTENDED USE

The product is intended for use in the iliac, superficial femoral artery (SFA) and above-the-knee popliteal artery for the following treatments:

- Arteriosclerotic stenosis
- Total occlusions that have been recanalized

The product provides mechanical support to maintain constant blood flow of the vessel.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

Placement of the stent requires advanced skills in interventional procedures. The following instructions will give technical guidance but do not obviate formal training in the use of the device.

Standard techniques for placement of an arterial access sheath and wire guide should be employed during use of the vascular stent.

CONTRAINDICATIONS

Patients with bleeding diathesis or other disorder, e.g. peptic ulceration or recent cerebrovascular accident, limiting the use of antiplatelet and anticoagulant therapy.

Patients with total arterial occlusion.

Coronary use is not recommended.

WARNINGS

- Possible allergic reactions should be avoided.

PRECAUTIONS

- Manipulation of the product requires high-resolution fluoroscopic control.
- The lesion and vessel can be predilated before placement of the stent.
- Do not try to remove the stent from the introducer system before use.
- Ensure that the safety lock is not inadvertently removed prior to stent release.

MRI INFORMATION

Non-clinical testing has demonstrated that the ZILVER FLEX™ 35 Vascular Stent is MR conditional. A patient with this stent can be scanned safely after placement under the following conditions.

Static Magnetic Field

- Static magnetic field of 3-Tesla or less
- Highest spatial magnetic gradient of 720-gauss/cm

Non-clinical evaluation was conducted in a 3-Tesla MR system with a maximum spatial magnetic gradient field of 720 Gauss/cm as measured with a gauss meter in the position of the static magnetic field pertinent to the patient (i.e., outside of scanner covering, accessible to a patient or individual) and resulted in a deflection angle of 3 degrees. An implant that displays a deflection angle of ≤ 45 degrees is considered to pose no harm to a patient or individual.

MRI-Related Heating

1.5 and 3.0 Tesla Systems

- Maximum MR system reported, whole-body-averaged specific absorption rate (SAR) of 2.9 W/kg for 15 minutes of scanning (i.e., per pulse sequence)

Non-clinical evaluation, including analysis of the effects of blood flow, was conducted for single and two overlapped versions of the ZILVER FLEX™ 35 Vascular Stents (lengths = 140 and 200 mm) producing a maximum temperature change of 3.1°C during 15 minutes of MRI performed at 1.5-Tesla (Magnetom, Siemens Medical Solutions, Malvern, PA) and 3-Tesla (Excite, General Electric Healthcare, Milwaukee, WI) under the following minimum conditions: 140 mm and 200 mm lengths:

	<u>1.5-Tesla</u>	<u>3-Tesla</u>
MR system reported, whole body averaged SAR	2.9-W/kg	3.0-W/kg
Calorimetry measured values, whole body averaged SAR	2.1-W/kg	2.8-W/kg

Image Artifacts

MR image quality may be compromised if the area of interest is in the same area or within approximately 5 mm of the position of the ZILVER FLEX™ 35 Vascular Stent, as found during non-clinical testing using the sequence: T1-weighted, spin echo and gradient echo pulse sequence in a 3.0 Tesla MR system (Excite, General Electric Healthcare, Milwaukee, WI). Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant. Heating in the MRI environment for stents with fractured struts is not known.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to:

- Abrupt stent closure
- Allergic reaction to nitinol
- Amputation
- Angina/coronary ischemia
- Arrhythmia
- Arterial aneurysm
- Arterial rupture
- Arteriovenous fistula
- Atheroembolization (Blue Toe Syndrome)
- Death
- Embolism
- Fever
- Hematoma/hemorrhage
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection/abscess formation at access site

- Intimal injury/dissection
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Myocardial infarction
- Pseudoaneurism formation
- Pulmonary embolism
- Renal failure
- Restenosis of the stented artery
- Septicemia/bacteraemia
- Stent malapposition
- Stent migration
- Stent strut fracture
- Stroke
- Spasm
- Tissue necrosis
- Worsened claudication/rest pain

PRODUCT RECOMMENDATIONS

Arterial access

For arterial access, we recommend the use of an access set that accepts a 6.0 French (2.0mm) introducer catheter.

Wire Guide Selection

The use of a .035-inch (0.89mm) wire guide is recommended.

PTA Balloon Selection

For pre- and post-dilation, an appropriately sized balloon catheter is recommended.

Stent Selection


The chosen stent diameter should be at least 1 mm larger than the diameter of the reference vessel.

PROCEDURE

Access

1. Gain access at the appropriate site utilizing an introducer sheath or guiding catheter. Refer to product labeling for required introducer sheath or guiding catheter size.
2. Introduce the extra or ultra stiff .035" (0.89mm) wire guide through the introducer sheath or guiding catheter across the distal segment of the target lesion.
3. Predilate if necessary. Remove the balloon catheter, leaving the wire guide in place.
4. Determine the proper stent size after complete diagnostic evaluation. The stent diameter should be at least 1 mm larger than the diameter of the reference vessel. Measure the length of the target lesion to determine the length of the stent required. Allow for the proximal and distal aspects of the stent to cover the entire target area. **NOTE:** If multiple stents are required to cover the length of the lesion, please refer to the Multiple Stent Placement section of these instructions for use for further recommendations.
5. Use the 1 ml syringe to flush the wire guide lumen with saline through the hub.
6. Immediately before placing the delivery system into the body, use a 1 ml syringe to flush the delivery system with saline through the side-arm flushing port. Flush until a few drops of saline exit the distal tip, between the delivery system outer sheath (d) and inner catheter (e).
7. Insert the delivery system over the wire guide.
8. Under fluoroscopy, advance the delivery system past the lesion.
9. Pull back on the stent delivery system under fluoroscopy until the radiopaque markers on the stent (i) are at the desired position. The stent is now ready to be deployed. (**Fig. 1**)


Fig. 1

Deployment of the Stent


1. Before deployment, it is important to straighten the proximal part of the delivery system as much as possible and to keep the handle in a stable position.
2. The stent expansion must be performed under fluoroscopic control.
3. Hold the hub (b) on the metal cannula (g) steady. To deploy the stent, remove the red safety lock (c). (**Fig. 2**)

Fig. 2


4. Hold the hub end stationary. The stent will be deployed as you pull the handle (a) toward the hub (b). (**Fig. 3**) The radiopaque marker on the delivery system (h) indicates the progress of the deployment. **NOTE:** Full deployment of the stent length will occur when the distal end of the sheath has been retracted past the proximal part of the stent.

Fig. 3


5. As deployment occurs, continue sliding the handle (a) toward the hub (b) in a slow, smooth and consistent fashion. (Fig. 4) **NOTE: Once stent deployment has begun, the stent must be fully deployed. Repositioning of the ZILVER FLEX™ 35 Vascular Stent is not possible since the delivery system's outer sheath cannot be readvanced over the stent once deployment begins. Refer to the Multiple Stent Placement section of these instructions for use for information on missed lesions.**

Fig. 4

6. The stent is fully deployed when the handle (a) reaches the hub (b). (Fig. 5)

Fig. 5

7. Perform an arterial angiogram to ensure proper placement of the device. **NOTE:** If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation (standard PTA) can be performed at the discretion of the physician. Select an appropriate size PTA balloon catheter and dilate the lesion with conventional technique. The inflation diameter of the PTA balloon used for post dilatation should approximate the diameter of the reference vessel. Remove the PTA balloon from the patient.

8. Delivery System Removal - Do not advance sheath after the stent has been deployed. **Delivery System can be removed without recapturing tip.** Check for delivery system integrity post removal from the patient.

9. Introducer sheath and wire guide may be removed at this point.

Multiple Stent Placement

If placement of more than one stent is required, the following recommendations should be considered:

- When more than one stent is required, resulting in stent-to-stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.
- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations (i.e., a second stent should be placed proximal to the previously placed stent).
- Stents placed in tandem should slightly overlap.

HOW SUPPLIED

Supplied sterile by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

ZILVER 518, ZILVER 635 STENT VASCULAR DESIGN

CAUTION: U.S. Federal Law restricts this device to sale by or on the order of a physician (or properly licensed practitioner).

Do not re-sterilize.

DEVICE DESCRIPTION

The Zilver® Stent Vascular Design (ZIV) is a self-expandable stent made of nitinol. It is a flexible, slotted tube that is designed to provide support while maintaining flexibility in the vessel upon deployment. Post-deployment, the stent is designed to impart an outward radial force upon the inner lumen of the vessel, establishing patency in the stented region.

The stent has radiopaque markers at both ends (h).


The stent is supplied preloaded in a 5.0 French (1.67 mm) or 6.0 French (2.0 mm) introducer catheter. A radiopaque marker (i) on the distal tip of the sheath is applied to visualize deployment of the stent. Stent deployment is controlled by means of a hand held device.

Refer to Stent Size Selection Tables below.

5.0 French (1.67 mm)												
Stent length (mm)	20		30		40		60		70		80	
Delivery system (mm)	80	125	80	125	80	125	80	125	80	125	80	125

Stent outer diameter (mm)	4	X	X	X	X	X	X	X	X			X	X
	5	X	X	X	X	X	X	X	X			X	X
	6	X	X	X	X	X	X	X	X			X	X
	7	X	X	X	X	X	X	X	X			X	X
	8	X	X	X	X	X	X	X	X			X	X
	9	X	X	X	X	X	X	X	X			X	X
	10	X	X	X	X	X	X	X	X		X	X	X

6.0 French (2.0 mm)													
Stent length (mm)		20		30		40		60		70		80	
Delivery system (mm)		80	125	80	125	80	125	80	125	80	125	80	125
Stent outer diameter (mm)	4												
	5	X	X	X	X	X	X	X	X	X		X	X
	6	X	X	X	X	X	X	X	X	X		X	X
	7	X	X	X	X	X	X	X	X	X		X	X
	8	X	X	X	X	X	X	X	X	X		X	X
	9	X	X	X	X	X	X	X	X	X		X	X
	10	X	X	X	X	X	X	X	X	X		X	X
	12			X	X	X	X	X	X	X		X	X
	14			X	X	X	X	X	X	X		X	X

- a. Handle
- b. Hub
- c. Safety Lock
- d. Introducer catheter
- e. Introducer tip
- f. Flushing Port
- g. Metal Cannula
- h. Radiopaque Markers
- i. Radiopaque Marker

INTENDED USE

The product is intended for use in the iliac artery for the following treatments:

- Arteriosclerotic stenosis
- Total occlusions that have been recanalized

The product provides mechanical support to maintain constant blood flow of the vessel.

The product is intended for use by physicians trained and experienced in diagnostic and interventional techniques.

Placement of the stent requires advanced skills in interventional procedures. The following instructions will give technical guidance but do not obviate formal training in the use of the device.

Standard techniques for placement of an arterial access sheath and wire guide should be employed during use of the vascular stent.

CONTRAINDICATIONS

Patients with bleeding diathesis or other disorder, e.g. peptic ulceration or recent cerebrovascular accident, limiting the use of antiplatelet and anticoagulant therapy.

Patients with total arterial occlusion.

Coronary and below the knee use is not recommended.

WARNINGS

Persons allergic to nitinol(nickel titanium) may suffer an allergic reaction to this implant.

PRECAUTIONS

- Manipulation of the product requires high resolution fluoroscopic control.
- The lesion and vessel can be predilated before placement of the stent.
- Do not try to remove the stent from the introducer system before use.
- Ensure that the safety lock is not inadvertently removed prior to stent release.

MR CONDITIONS

Non-clinical testing has demonstrated that the Zilver Vascular Stent is MR Conditional. It can be scanned safely under the following conditions:

- Static magnetic field of 3 Tesla or less
- Spatial gradient field of 720 Gauss/cm or less
- Whole body averaged specific absorption rate (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes of scanning (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla), respectively.

In non-clinical testing, the Zilver Vascular Stent produced maximum temperature rises of 0.1, 3.8, 0.8, and 0.1 degrees C (for a single stent at 1.5 Tesla, a pair of overlapping stents at 1.5 Tesla, a single stent at 3 Tesla, and a pair of overlapping stents at 3 Tesla, respectively) at whole body averaged specific absorption rates (SAR) of 1.5 W/kg (for a single stent at 1.5 Tesla) and 3 W/kg (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) for 20 minutes (for a single stent at 1.5 Tesla) and 15 minutes (for a single stent at 3 Tesla and a pair of overlapping stents at 1.5 and 3 Tesla) of MR scanning in a 1.5 Tesla/64 MHz General Electric MR scanner, a 1.5 Tesla Magnetom Siemens Medical Solutions MR Scanner (to evaluate a pair of overlapping stents), and a 3 Tesla Excite General Electric MR scanner. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the Zilver Vascular Stent. Therefore, it may be necessary to optimize MR imaging parameters for the presence of this metallic implant. Heating in the MRI environment for stents with fractured struts is not known.

POTENTIAL ADVERSE EVENTS

Potential adverse events that may occur include, but are not limited to:

- Abrupt stent closure
- Allergic reaction to nitinol
- Amputation
- Angina/coronary ischemia
- Arrhythmia
- Arterial aneurysm
- Arterial rupture
- Arteriovenous fistula
- Atheroembolization (Blue Toe Syndrome)
- Death
- Embolism

- Fever
- Hematoma/hemorrhage
- Hypersensitivity reactions
- Hypotension/hypertension
- Infection/abscess formation at access site
- Intimal injury/dissection
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Myocardial infarction
- Occlusion
- Pseudoaneurism formation
- Renal failure
- Restenosis of the stented artery
- Septicemia/bacteraemia
- Stent malapposition
- Stent migration
- Stent strut fracture
- Stroke
- Spasm
- Trombosis
- Tissue necrosis
- Vessel trauma/vessel rupture
- Worsened claudication/rest pain

CHEMICAL TOLERANCES

- The stent must not come in contact with metal other than nitinol. Please pay attention to this if the patient has an arterial stent already.

PRODUCT RECOMMENDATIONS

Arterial access

For arterial access, we recommend the use of an access set that accepts the introducer catheter (5.0 French [1.67 mm] or 6.0 French [2.0 mm]).

Wire Guide Selection

For vascular introduction a .018-inch (0.46 mm) wire guide is recommended for the 5.0 French (1.67 mm), and a .035-inch (0.89 mm) wire guide for the 6.0 French (2.0 mm) introducer catheter.

PTA Balloon Selection

For pre- and postdilation, an appropriately sized balloon catheter is recommended.

Stent Selection

The chosen stent diameter should be at least 1 mm larger than the diameter of the reference vessel. Determine the proper stent size after complete diagnostic evaluation. The stent deployment must be performed under fluoroscopic control. Measure the length of the target lesion to determine the length of the stent required. Allow for the proximal and distal aspects of the stent to cover the entire are.

PROCEDURE

Access

1. Upon removal from the package, inspect the product to ensure no damage has occurred.
2. Introduce the recommended access catheter using standard percutaneous technique.
3. Introduce the wire guide through the access catheter across the distal segment of the target lesion.
4. Remove the access catheter, leaving the wire guide in place.

5. Predilate if necessary. Remove the balloon catheter leaving the wire guide in place.
 6. Determine proper stent size.
 7. Prior to insertion, flush the delivery system with saline through the flushing port.
 8. Insert the delivery system over the wire guide.
 9. Place the delivery system under fluoroscopy to determine that the radiopaque markers on the delivery system (h) are at the desired position. The stent is now ready to be deployed (Fig. 2).
- NOTE: If multiple stents are required to cover the length of the lesion, please refer to the Multiple Stent Placement section of these instructions for use for further recommendations.


Fig. 2

Deployment of the stent

1. Before deployment, it is important to straighten the proximal part of the delivery system as much as possible and to keep the handle in a stable position.
2. The stent expansion must be performed under fluoroscopic control.

Fig. 3a



Fig. 3b

3. Remove the safety lock (Fig. 3a). Gently slide the handle towards the hub (Fig. 3b) and the stent will deploy. The radiopaque markers on the delivery system indicate the progress of the deployment.
4. Relocation of the stent during deployment is only possible to a more proximal position; carefully pull the system towards you (Fig. 4). When proper location has been achieved, continue deployment.

Fig. 4



Fig. 5

When the handle reaches the hub, the stent has been fully deployed (Fig. 5). Post stenting dilation can now be performed if required. Remove the wire guide.

5. Perform completion angiogram to ensure proper placement of the device, remove the introducer catheter.

6. Delivery System Removal - Do not advance sheath after the stent has been deployed. Delivery system can be removed without the need to recapture tip. Check for delivery system integrity post removal from the patient.

Multiple Stent Placement

If placement of more than one stent is required, the following recommendations should be considered:

- The stent must not be in contact with metals other than nitinol. Please take this into consideration if the patient already has a metal stent.
- In relation to the lesion site, the distal area of narrowing should be stented first, followed by the proximal locations, i.e. a second stent should be placed proximal to the previously placed stent.
- Stents placed in tandem should slightly overlap.

HOW SUPPLIED

Supplied sterile by ethylene oxide gas in peel-open packages. Intended for one-time use. Sterile if package is unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. Store in a dark, dry, cool place. Avoid extended exposure to light. Upon removal from the package, inspect the product to ensure no damage has occurred.

REFERENCES

These instructions for use are based on experience from physicians and (or) their published literature. Refer to your local Cook sales representative for information on available literature.

COMPLIANCE WITH RELEVANT REGULATIONS AND TECHNICAL STANDARDS

Stents to restore and maintain patency of the vessels are in compliance with the Council Directive (MDD) 93/42/EEC.

Cook Incorporated, Cook Ireland Ltd and William Cook Europe manufactures all medical products in accordance with Medical Devices Directive 93/42/EEC.

Cook Incorporated manufacture all medical products in accordance with the United States Food and Drug Administration (FDA) Quality System Regulations, 21 CFR Part 820, Quality System Regulations.

All of Cook facilities are also certified to ISO 13485 (2003) Quality Management System.

WARRANTY

An authorized representative of the manufacturer in the territory of the Russian Federation is the ZAO "SHAG".

Russia, 119002, Moscow, Karmanitskiy per., 9.office 501a.

phone: (495) 956-13-09;

Fax: (495) 956-1310,

e-mail: info@schag.ru

Manufacturer's warranty binding on quality products.

SHELF LIFE

The "use by" date as denoted by symbols on the labelling. The shelf life is shown by indication the time of manufacture (factory symbol) to the stated expiration date (hourglass symbol). The difference between both dates is the shelf life.

The Zilver[®] Vena[™] Venous Self-Expanding Stents (ZVT7-), Self-Expanding Vascular Stents (ZIV5-, ZIV6-, ZFV6-), Formula 418 Balloon-Expandable Stent (FOV418), Renal (FORX4-) and Vascular Stent (FOVX414-), and the Formula 535 Balloon-Expandable Stent (FOV 535) are labeled with a three year shelf-life.

STERILIZATION DETAILS

These devices are supplied sterile in peel open packages utilizing an EtO (Ethylene Oxide) validated method to a sterility assurance level (SAL) of 10⁻⁶. These devices are intended for one-time use, sterile if unopened or undamaged. Do not use the product if there is doubt as to whether the product is sterile. These devices should be stored in a dark, dry, cool place; avoid extended exposure to light.

STORAGE RULES

These devices should be stored in a dry location, away from temperature extremes.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation at the temperature of from 0 till + 40°C and humidity of from 60 till 95%.

REPAIR/MAINTENANCE


These devices are single use only, do not attempt to repair, carry out maintenance or reuse.

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with Council Directive 93/42/EEC

For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms N2.1.7.2790-10 of December 12, 2010 "Sanitary-epidemiological requirements for the treatment of medical waste". Risk class of medical waste: class B.

Regulatory Affairs Manager

These devices should be stored in a dry location, away from temperature extremes.

TRANSPORTATION RULES

Transportation shall account for the type of transportation employed, the destination, the vendor used for transportation and an estimation of the frequency and harshness of any handling required during transportation at the temperature of from 0 till + 40°C and humidity of from 60 till 95%.

REPAIR/MAINTENANCE


These devices are single use only, do not attempt to repair, carry out maintenance or reuse.

DISPOSAL RULES

Upon completion of procedure, the devices are to be disposed of as per institutional guidelines for biohazardous medical waste.

For EU countries -in accordance with Council Directive 93/42/EEC

For Russian Federation – in accordance with national requirements, regulated by sanitary rules and norms N2.1.7.2790-10 of December 12, 2010 "Sanitary-epidemiological requirements for the treatment of medical waste". Risk class of medical waste: class B.

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